

"Build today, then strong and sure, With a firm and ample base; And ascending and secure. Shall tomorrow find its place." Longfellow, Henry W. (1807-1882)








Where there is a successful business, it is because someone has taken a courageous decision.

Peter Drucker

Located at Km. 585.1 of the A-2 highway, at the exit Abrera (Barcelona), covering an area of 21000 m², and factory build 14,800m². Situated 28 Km from the city of Barcelona, with excellent links by motorway AP-7, exit 25 (Martorell) or by the A-2 highway, Km 585.1 Barcelona-Madrid, industrial estate "Can Amat" next to the production centre of SEAT.







A-2 highway, Km. 585,1 - 08630 ABRERA (Barcelona-SPAIN) - P. +34 93 770 08 77* (6 lines) - F. +34 93 770 23 62 - www.jpselecta.es - selecta@jpselecta.es



J.P. SELECTA S.A.

Scientific Laboratory Equipment Manufacturer



COMECTA S.A.

Exclusive Representation



AQUISEL S.L.

manufacturer of Disposable Clinical Test Kits





Own Brands:



ISO 9002:2000

Analytical Division

International Patents

Tubes for Blood Analysis

Pipettes for Determination of Erythrocyte Sedimentation Rate E.S.R.

Sealed Containers for Urine Samples

Fixation and concentration of fecal parasites

Greeting



1940-2015

The publication of this new General Catalogue meets 66 years since I undertook my activity as a producer in a simple workshop. After this time the Selecta Group has now reached a significant volume of transactions, filling a market segment so broad that we have been projected to the level of the leading European manufacturers.

Our impact in the international markets, consolidated in the last years, and the domain in the Spanish area, is the response to the determination of implementing standardization, design, research and innovation to our production processes.

But I firmly believe that we must add to this success, a tool that we make available to our customers and that has been confirmed as essential to our project: the General Catalogue. I must confess that I have dedicated a very personal care, investing much of my time writing and editing more than our 30 catalogues and thousand individual brochures we have published during this time.

I have always concerned about the accuracy of the data that the catalogue provides, the perfection of images and the rationalization of texts. This new catalogue 2015-2017, with nearly 400 pages and 25 different equipment specialties, and bringing together more than 1000 equipments and 3000 different codes, must be a good example of this concern.

I trust that you can take advantage of the usefulness of this GENERAL CATALOGUE, and I remind you that, in a society as computerized as the present, I could not miss a last generation web page with complete information which can be printed and presented in several different languages: www.jpselecta.es

Jupice but

Sincerely

Joaquín Puig Alemany President of Selecta Group



When we launch a new catalogue, we try to evaluate the effectiveness of the former, in this case the one published in 2011. To sum up, during the period 2011-2013, despite the stagnation of the world's economy, the growth of Selecta Group shows that our model of business communication, by using the General Catalogue, transmits the product quality and benefits that has made us highly competitive both nationally and internationally.

That's the reason why we have respected the design that has given this character to our publications. Our new catalogue responds to this philosophy.

In this Catalogue 2015-2017 you will find the widest range of laboratory products a manufacturer may exhibit, and you will check the different engineering innovations achieved in our company.

The restructuring process of our $\mathbb{R} + \mathbb{D} + I$ department to better accommodate to technological trends has been one of the proposed goals and thus we can offer important new features to our users:

- Implementation of a TFT (Thin Film Transistor) touchscreen which regulates all parameters and establishes programs and graphs on a large number of our devices.
- Expansion of our products offer with several latest equipments, identified with the "NEW" icon.

We want these changes and improvements to be interpreted as a sign of vitality. We trust we could expand our offer in the coming years, and confirm our competitiveness and take advantage of the global economic balance which is being announced.

We are available to complete this General Catalogue 2015-2017 interpretation since this is our commitment to strategic communication plan.

Your faithfully,

Jordi Puig Herrera Managing Director of Selecta Group



MANAGEMENT TEAM

From left to right: DAVID PECANINS, Quality control manager - PEDRO CARAFI, Purchasing manager - JOAQUÍN ARQUE, Finance LUIS TREMOLEDA, Advisor - TONY EJARQUE, Commercial director - JOAQUÍN PUIG, President - JORDI PUIG, Managing director ANTONIO BERGÉ, Production manager - CARLOS BARBERO, Logistic manager - FRANCISCO SANTIAGO, Quality control RICARD CARDÚS, Techincal director - ANTONIO ALARCÓN, Comecta, s.a. director - MANEL SOTELO, Aquisel, s.l director

J.P. Selecta World wide



INFORMATION

- The apparatus presented in this catalogue is designed specifically for the application intended.
 Contact us for confirmation of suitable use in any other distinct application.
- 2. In accordance with existing standards. All tests on our equipment are completed at an ambient temperature of 20 °C.
- 3. For correct, installation, use and maintenance of our equipment read all the instructions before use.
- 4. J. P. Selecta SA reserves the right to modify some of the parameters or characteristics of the equipments described in this catalog, in order to improve their benefits or adapt them to changing current regulations at all times.

VOLTAGES AND FREQUENCY

- 5. Unless stated otherwise in the catalogue, our equipment is manufactured to be connected to the mains as:
 - Single phase 230-V ±10% 50 Hz.
 - Three phase 230/400V \pm 10% 50 Hz.
 - Single phase 110-V \pm 10% 60 Hz. (on demand).
- 6. If required other voltages and frequencies can be supplied. In some application the power can be modified.
- 7. Heating equipment without motors can operate on both 50 and 60 Hz.
- 8. In all cases where different voltages, frequencies and plugs type are requested. Written confirmation will be required before order acceptance, otherwise standard power and Schuko plug will be supplied.

CALIBRATION

If required certificates of calibration can be supplied that are traceable to known certified standards.

ACCESSORIES AND SPARES

To maintain the equipment to the original specification and design, it is important that genuine J.P. Selecta spares be used. In the case that original spares are not used, J.P. Selecta are not responsible for any damage or anomaly caused.

J.P. SELECTA HAVE OPENED A WINDOW TO THE WORLD

Discover all our new developments by connecting to our Web site:

http://www.jpselecta.es

You can see and match to your requirements all of our equipment and products.



OUR GUARANTEE

We guarantee the supply of spare parts for the life time of the production of the model manufactured.



Image showing our technical services department still repairing equipment manufactured in the 1970's



RS-232 communication software

SUITABLE FOR RECORDING DATA, REMOTE CONTROL AND ATOMISATION.

Suitable for all equipment with



J.P. SELECTA, s.a.



Comunication software for laboratory equipment





J.P. SELECTA,s.a.

Autovía A-2. Km 585,1 08630 Abrera (Bacelona) Spain www.jpselecta.es jpselecta@jpselecta.es

Many of the **J.P.SELECTA**, **s.a.** manufactured equipment comes with this software as standard. Installed on a PC with a direct communication link to the equipment for data transfer or control. The software function depends on the equipment interface but in general (ovens, baths, thermostats etc.) can be connected for:


- Logging of data: The PC can constantly receive temperature data for ovens, baths, etc.

- Capture of data: Accepts data for printing or incorporation into a report.

- Remote control: Communication direct to the microprocessor control panels of, ovens, baths, etc. with automated control direct from

the PC to the instrument..

J.P.SELECTA, s.a. can manufacture software to order for automated laboratory processes, e.g. temperature cycling and time etc.







Quality Control Systems

J.P. Selecta complies with the obligatory controls and regulations related to the safety of electrical equipment, the measurement and control within the laboratory, and the completion of the standards and CE directives that affect the manufacturer of our equipment.

ELECTRICAL SAFETY 73-23 CE. ELECTROMAGNETIC COMPATIBILITY 89-336 CE. PRESSURE VESSEL 97-23 CE

Details of Standards and directives:

UNE-EN 61010-1 Part 1: General requirements.

UNE-EN 61010-2-010 Part 2-010: Required for the particular use of heated laboratory equipment

and materials.

UNE-EN 61010-2-020 Part 2-020: Particular requirement for laboratory centrifuges.

UNE-EN 61010-2-041 Part 2-041: Requirement for Autoclaves that use steam vapour for the treatment

of materials used within medicine and laboratory processes.

UNE-EN 61010-2-051 Part 2-051: Particular requirement for shaking and mixing equipment.


UNE-EN 61326 Requires electromagnetic compatibility.






Certificate of registration AENOR for the design, production, commercialisation and post sales service of :ovens, centrifuges, stirrers, thermostatic immersion baths, ultrasonic baths, hotplates and blocks, digesters, Kjeldahl steam distillations, fat, fibre and alcohol extractors, refrigerators, freezers and environmental cabinets.











THE INTERNATIONAL CERTIFICATION NETWORK

CERTIFICATE

IQNet and AENOR hereby certify that the organization

J.P. SELECTA, S.A.

CR N-II, KM 585,1 08630 - ABRERA(BARCELONA) España

for the following field of activities

The design, production, commercialization and aftersales of:
oven, centrifuges, autoclaves, stirrers, inmersion thermostats, thermostatic baths, distillator, ultrasonic cleaners, hotplates, heating mantles,
digestion blocs, kjeldahl distilling unit, fat extractor, fiber extractor, alcohol distilling unit, blood bank refrigerator cabinets, refrigerated
cabinets, chest freezer.

The commercialization of in Vitro diagnostic medical devices

has implemented and maintains a

Quality Management System

which fulfills the requirements of the following standard

ISO 9001:2008

First issued on: 1999-08-01

Last issued: 2014-02-22

Validity date: 2017-02-22

spañola de AENOR

Registration Number: ES-0109/1996

- I Net

Michael Drechsel President of IONet Yormalización y Certificación

Aveling BRITO
Chief Executive Officer

IQNet Partners*:

AENOR Spain AFNOR Certification France AIB-Vincotte International Belgium ANCE Mexico APCER Portugal CCC Cyprus CISQ Italy CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil FONDONORMA Venezuela ICONTEC Colombia IMNC Mexico Inspecta Certification Finland IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding CmbH and NSAI Inc.

* The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under www.iqnet-certification.com



On request, we can supply calibration certificates of the working conditions of the instruments before they leave the factory. These certificates can be used for OQ, IQ validations and for equipment that needs to comply to a ISO9000 quality system.

These certificates can be produced for Ovens, Incubators, Refrigerated products, Autoclaves, Baths, Inmersion thermostats, furnaces, etc. instrument that require a temperature measurement by realising several temperature sensors and readings to obtain the probable temporal and spatial uncertainty.

Note: Consult for other types of equipment.

EXAMPLE OF CALIBRATION CERTIFICATE OF BACTERIOLOGICAL AND CULTURES OVEN:


Part No. 1001738 Calibration certificate, 2 points.

The calibration certificate has been completed based on two measuring points. This means to the user that the requested certificate should indicate two temperatures at which the effective calibration should be completed. For more than two measuring point please consult J.P. Selecta.

When requesting the certificate full details of the company and address can be included in the certificate information.

Note: When asking for a certificate, please indicate data (company name and address) for its validity.


















Selecta Group









Pressurized and acclimatized room with sophisticated serial production automatisms.





OVENS WITH TOUCH SCREEN Pages. 139 to 156



MUFFLE FURNACES WITH TOUCH SCREEN Pages. 162 to 163



CENTRIFUGES WITH TOUCH SCREEN Pages. 115 to 133



STIRRERS WITH TOUCH SCREEN Page. 31



OVERHEAD STIRRER Page. 41



AUTOCLAVE ELEVATOR Page. 70



IMMERSION THERMOSTAT WITH TOUCH SCREEN

Pages. 101 to 105



DISTILLERS Pages. 263 to 265



ELECTRONIC REGULATOR WITH TOUCH SCREEN Page. 293



REACTOR Page. 298 and 299



SHAKERS FOR PLATELET BAGS Page. 79



PLASMA BAGS DEFROST Page. 83



PORTABLE FOCULATOR Page. 262

AAAA.



ORBITAL SHAKER INCUBATORS Pages. 42 and 43



INCUBATION CHAMBER Page. 79



LABORATORY FUME CUPBOARDS Page. 284



INCUBATOR FOR PETRI CAPSULES Page. 154



FIXATION AND CONCENTRATION OF FAECAL PARASITES
Page. 309




GENERAL INDEX BY INSTRUMENT TYPE

COLOUR IDENTIFICATION













Pages 25 to 45





















Pages 55 to 70











Pages 71 to 77

Refrigerated Cabinets, Freezers,...











Pages 78 to 83













Pages 84 to 97

Water, Oil, Reciprocal Shaking, Dry-Blocks, Ultrasonic Cleaners,













Pages 98 to 110

Recirculation, Refrigerated, Cooling Coils,...













Pages 111 to 133

Digital Control, Brushless, Microprocessor Control, Refrigerated, Blood Bags,...













Pages 134 to 164

Drying, Bacteriological, Vacuum, CO2, Low Temperature,...











Pages 165 to 211

Clamps, Stands, Spatulas, Racks, Burners, Stirrer Bars, Porcelain,...











Pages 212 to 224













Pages 225 to 236













237 to 248














Pages 249 to 265











266 to 282

Microscopes, Refractometers, Melting Point, Polarimeter,...











283 to 301


Safety Cabinets, Viscometers, Monitor and Control, Reactos, Filters,...













302 to 313

Medical devices, Petri plates

ALPHABETIC INDEX

Α				G	
Alcohol burner	195	Colorimeters	243 to 244	Glass drying oven	138
Alcohol distillation	259 to 260	Conductivity meter	248	Glass slides	215
ANALYTICAL CONSUMABLES	302 to 311	Container cryo cooler	190	Glass thermometers	296-297
ANALYTICAL TECHNIQUE	237 to 248	Containers sterilization	194	Gloves	200
AUTOCLAVES	55 to 70	Control thermometers	296 to 297	_	
Automatic Kjeldahl unit	253 to 256	Controller electronic sensoterm	293	Н	
Automatic time switch	295	Cooling bath	218	HAEMATOLOGY	79 to 83
Automatic sealer	83	Cooling cabinet plate	217	Hand operated centrifuge	114
		Cooling immersion unit	105 to 106	Hand microtome	216
В		Crane for autoclaves	70	Heating blocks	108 to 110
BACTERIOLOGICAL INCUBATORS	151 to 157	Crucible lids	182	HEATING MANTLES	225 to 231
Bag sealer	68	Crucibles stainless steel, Ni & Zirc	o182	Heating tapes	231
BALANCES	267 to 269	Cryothermostats	104	Hidraulic press	191
Baskets stainless steel	181 to 184	Cuvette quartz and glass	241	High temperature oven	147
Bath paraffin wax	217	Chest freezer	77	Histology Paraffin	92
Bath reciprocating shaker	93	Chromatography tank	219	Homogeneizers	44 to 45
Bath viscometer	94	Chronometer digital		HOTPLATES	232-236
BATHS		_		Hotplates, ceramic	232
Baths for water and oil	86 to 92	D		Humidity meter	294
Baths ultrasonic cleaning	95 to 97	Density measurement	268	_	
Beakers stainless steel		Desiccator hidroscopic		1	
Biological oxygen demand COD	261	Desiccator vacuum		Illuminator cold light	278
Block digester Kjeldahl		Detergent biosel	63	Immersion cooler	
Blood banks instrumentation		Diamond point pencil		Immersion thermostats	100 to 109
Blood banks cabinets	81 to 82	Digestion & distilling tubes		Incubation chamber	157
Bossheads	166	Digestion system		Incubator COD	261
Box cuvette storage	189	Digital thermometer		Incubator CO ₂	
Box for cryovials 15 / 2 ml		Dispenser wax		Incubator orbital shaker	
Box slide		Dispensing pipettes		Incubators	
Brushes		Dissecting kit		Incubator for Petri capsules	154
Bucket with handle	181 to 184	Distillation still		Infra-red incinerator	
Burners	195 to 198	Dry block heater	108 to 110	INSTRUMENTAL	165 to 211
		Drying rack		INSTR. BIOLOGY AND HISTOLOGY	212 to 224
C		Drying and sterilization ovens			
CO ₂ Incubator	158			K	
Condenser		E		Kidney dish	183
Callipers	201	Electrical contact thermometers	293	KJELDAHL DISTILLING UNITS	251 to 256
Caps, stoppers serology		Extractor fat and oils		Kolle handle loop	214
Cell counter		Extractor fiber unit	258	_	
Centrifugal pumps	53	_		L	
CENTRIFUGES		F		LABORATORY ACCESSORIES	199 to 211
Centrifuge blood bags		Faecal parasites fixation/concentra	ation309	LABORATORY-HOSPITAL	
Centrifuge, refrigerated12		Flake ice maker		Laboratory jack	
Ceramic hotplates		Floculators		Laboratory lamp	
Circulation pump		FILTRATION		Laboratory fume cupboards	
Clamps and stands		Forceps dissection		Laminar flow cabinets	
Cold plate histology		FURNACE MUFFLE		Liquid transfer	
Colony counter digital		Fusion spoon flat bottom		LOW TEMPERATURE71 to	
,g minimum					

M					
Magnifiers	201	Pumps circulation	50 to 53	Thermometers	296 to 297
Mantles bank for Kjeldahl	230	Pumps vacuum	50 to 53	Thermoreactor COD	261
Melting point	282			THERMOSTATIC BATHS	98 to 110
Membrane vacuum	50 to 52	R		Thermostat dry blocks	108 to 110
Messuring range	201	"		Thermostat viscometer	290
Metal blocks	108 to 110	Racks	185 to 189	Tilting mixer shaker	35-109 to 110
Microcentrifugues	114	Refractometers	279 to 281	Timer alarm	295
Micrometer palmer	201	Reactors		Tong for flasks	170
Microplate reader	245	Refrigerated cabinets	159 to 161	Tripod base stand	174
Microplate washer	245	Refrigerators	73 - 77 - 80 - 83	Tubes with screw cap	310
Microtiter	244	Rocking mixer	36		
MICROSCOPES	270 to 279	Roller & tilt mixer	35	U	
MONITOR AND CONTROL	293 to 295	Rotator	35	Ultrasonic cleaning baths	95 to 97
MUFFLE FURNACES	162 to 164	Rotatory evaporator	47 - 48	Ultrasonic homogeneizers	44
Multi-parameter water quality	y meter244	S		Ultraviolet sterilizer	224
N		Sand bath	87	V	
Nutritional analysers and water	er 249 to 265	Scalpel blades		Vacuum dry oven	148 to 149
matritional analysons and water	01	Scissors universal		Vacuum pump	
0		Scoop with handle		Vacuum sealer	
Oenological distiller unit	250 to 260	Sealers		Video microcamera	
Orbital shakers		Self-sealing parafilm		Viscometer bath	
OVENS		Shaker, orbital and horizontal		Viscometer flow cup	
Oven,vacuum		Slide boxes		Viscometer rotational	
Ovens, refrigerated		Slides		Viscometers glass	
Ovens, sterilization		Spatulas		viocomotoro giuco illinilinilini	
Ovens, universal		Spectacles		W	
Overhead stirrer		Spectrophotomers		Water circulating unit	107
0.00.000		Sponge scrubber		Water distillers	
P		Stabilising ring		Wax dispenser	
Palmer micrometer	201	Staining tank		Weights for balances	
Paraffin wax cabinet		Stands		Westergreen	
Paraffin wax dispenser		Steam distillation Kjeldahl		Wire baskets	
Parasites faecal fixation/conc		Stereomicroscopes			
Peristaltic pumps		Sterile stoppers			
Petri dishes		Stirring shafts			
Petri capsules incubator	154	Stirrer magnetic			
Ph electrodes		Storage cabinets for anatomy			
pH meters	246 to 247	Strip reader			
Photometers		Suctioner for pipettes			
Pipettes	220 to 222	Support arm for electrodes			
Pipettes rack		Support westergreen			
Polarimeter		_			
Polycarbonate tubes		T			
Polypropylene tubes		Table drain rack	193		
Poupinel sterilizer		Tanks,immersion thermostats			
Porcelain		Test tube racks			
Pumps centrifugal		Thermocycler			
- -		-			

Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page
022113224	I 1000442233	1000725234	1000985287	1001129199	1001265297	1001467171	1001597194
0106006265	1000443233	100078063-194	1000986287	1001130199	1001266297	1001468171-213	1001598194
0203681162	1000444233	1000781231	1000988287	1001131202-222	1001267297	1001469171-213	1001599181
0203692162	1000481193	1000782231	1000990287	1001132202-222	1001268297	1001470171-213	1001600181
068620352	100048733-204	1000783231	1000995286	1001133202-222	1001296105-106	1001471171-213	1001601181
068621052	1000488214	1000784231	1000996288	1001138218	1001298297	1001472171-213	1001602181
068621652	100049563-194	100078670	1000998286	1001139218	1001299296	1001478206	100161034
0703052265 0928292131	100049663-194 1000500295	100078833-204 100078933-204	1000999286 1001000286	1001140218 1001141218	100130674 100130774	1001479206 1001480206	1001611287 1001612287
100000133	1000504175	100079933-204	1001000286	1001142218	100130774	1001481206	1001613287
1000007193-230	100053633-204	100079133-204	100100336-94	1001143202	1001314192	1001482206	1001616286
100000933-204	100053733-204	1000801171-124	100100732	1001144202	1001330170	1001483206	1001617286
1000011172	100053833-204	1000802214	1001009234-293	1001145202	1001331170	1001484206	1001619159
100001233-204	100053933-204	1000803214	1001013289	1001146202	1001332199	1001485206	1001620159
1000013214	1000544100-106	1000805180	1001015289	1001147202	1001333202	1001486206	1001623287
1000014214	1000545100	1000806171-213	1001016289	1001148202	1001334202	1001487206	1001624287
1000015193-259	1000546100	1000807171-213	1001019169	1001149202	1001335202	100150182	1001625288
1000016199 100001933-204	1000547100 1000548100	1000809213 1000818175	1001030286 1001031286	1001150202 1001151202	1001336202 1001337171	100150282 100150382	1001626288 1001627288
100001933-204	1000548182	1000825175	1001031286	1001152213	1001337171	100150382	1001628288
100002563-96	1000564182	1000826175	1001033286	1001161177	1001340188	1001518178	1001629259
100002663-96	1000565182	1000827175	1001034286	1001162219	1001341188	1001523188	1001630259
100003233-204	1000566182	1000829194	1001035286	1001163219	1001342192	1001524189	1001631259
1000033214	1000570182	1000830187	1001036287	1001164219	1001346175	1001525189	1001632259
1000037167	1000571182	1000831186	1001037287	1001168199	1001347175	1001526189	1001633259
100003933-204	1000572182	1000832187	1001038287	1001169177	1001352218	1001527189	1001634259
100004533-204	1000573182	100085263-96	1001039287	1001170177	1001354180	1001528188	1001635259
1000049172	1000574182 1000575182	100085363-96	100104449	1001171177	1001359213	1001529189	1001636259 1001640192
1000050172 1000063171	1000576182	1000855178 1000856179	100106133-204 100106233-204	1001189171-124 1001190171-124	1001360213 1001361213	1001530180 100153374	1001640192 1001660311
1000064171	1000577182	1000857179	100106333-204	1001191171-124	1001362213	100153674	1001661311
1000067171	1000580171-124	1000858179	100106433-204	1001192213	1001363213	1001540107	1001662311
1000069166	1000595199	1000859180	1001070287	100119991	1001364213	1001541196	1001663311
1000093224	1000596194	1000860180	1001071287	100120091	1001365213	1001542189	1001664311
1000096224	1000603186	1000861180	1001072287	100120191	1001366215-278	1001543189	1001665311
1000100296	1000604186	1000862180	1001073287	100120291-105	1001367215	1001544189	1001666311
1000101172	1000605186	1000863193	1001074287	100120391-105	1001368216	1001545189	1001675158
1000104230	1000606186	1000864171-124	1001075287 1001090177	100120491-105	1001369216	1001552174	1001699192 1001702193-222
1000105171-215 1000106230	1000607186 1000608186	1000865171-124 1000866214	1001090177	100120591 1001212185	1001370216 1001371216	1001556114 1001557117-122	1001702193-222
1000100172	1000609186	1000867214	100109417	1001215193-220	1001374106	1001558117-122	100171234
100012086-296	1000610186	1000869216	1001095169	100121763-194	1001394131	1001559133	1001713189
1000123289	1000611186	100087291-296	1001096169	100121863-194	1001403150	1001562133	1001714189
1000206176	1000612186	1000873296	1001097169	1001221169	1001422259	1001563133-185	1001715186-310
1000208176	1000620182	1000875196	1001098169	100122263-194	1001423169	1001565133	1001716186-310
1000216180	1000621182	1000876196	1001099169	1001223193	1001434195	1001566133	1001717216
1000217180	1000622182	1000877197	1001100169	1001224287	1001435202	1001567133	100171834
1000236176 1000238176	1000623182 1000624182	1000878197 1000893106-293	1001101169 1001102172	1001225287 100123034	1001436203 1001437203	1001568133 1001570133	100171934 1001720169
1000251176	1000625182	1000895293	1001104177	100123134	1001437203	1001574133	1001733199
1000259194	1000635258	100089733-204	1001105177	100123394	1001439203	1001581178	1001734177
1000275213	1000641261	100089833-204	1001114197-215	100123494	1001440203	1001582178	1001735177
1000283170	1000642261	100090033-204	1001115188	100123594	100145394-290	1001583178	1001736177
1000284170	1000646259	100090133-204	1001116188	1001251187	100145494-290	1001584179	1001737177
1000285169	1000661171	100090233-204	1001117188	1001252187	100145594-290	1001585179	100173811
1000294177	1000662171	1000903214	1001118188	1001254192	100145694-290	1001586179	1001744299
1000347289 1000394106	1000684176-197	1000906181 1000907181	1001119188	1001255171-200 1001256195	100145794-290	1001587179	1001745299
1000395106	1000685176-197 1000686176-197	1000907181	1001120188 1001121188	1001256195	100145894-290 100145994-290	1001588179 1001589179	1001746298 1001747298
1000396106	1000687170-137	1000909181	1001121188	1001258297	1001459166	1001569163-170	1001748298
1000397100	1000688176	1000912224	1001123188	1001259297	1001461166	1001591163-170	1001749298
1000412180	1000689176	1000973157	1001124188	1001260297	1001462166	1001592180	1001750298
1000413180	1000705289	100097836-94	1001125188	1001261297	1001463168	1001593180	1001751298
1000431213	1000707289	100097936-94	1001126194	1001262297	1001464169	1001594181	1001752299
1000433187	1000716169	100098036-94	1001127194	1001263297	1001465169	1001595181	1001753299
1000434187	1000718233	100098136-94	1001128194	1001264297	1001466192	1001596181	1001754298

Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page
1001755298	1009860182	1026112185	1125115194	1581704305	2000239154	2101504160	300086997
1001757298	1010280179	1026124185	1125121194	1581705305	2000240154	2101505160	3000886294
1001758298	1010313178	1026148185	1125124194	1581814305	2000368163	2101506160	3000897234-293
1001765298	1010318178	1026412185	1191100181	1581942305	2000369163	2101507161	300091097
1001773299	1010321178	1026424185	1196120178	1581944305	2000376162	2101508160-161	3000914262
1001780201	1010420170	1026448185	1196150178	1581945305	2000377162	210150980	3000915262
1001781201	1010430170	1027680178	1198302182	1581952305	2000381138	210151080	300091686
1001801160-161 1001804160	1010440170 1010450170	1028130170 1028140170	1198402182 1198502182	1581954305 1581962305	2000787137 2000911137	2101618159 210162174	300091734-86 3000919232
1001805160	1012126181	1028230170	1198602182	1581964305	2000911157	210162280-81	3000920232
1001806160-161	1012128181	1028240170	1276120178	1581965305	2000964156	2101623140-156	3000921232
100181077	1012131207	1028301207	1276150178	1581974305	2000965156	2101624140-156	300095740-157
100181177	1012132207	1028302207	1276210178	1581975305	2000994151	2101625140-156	300095879
100181277	1012133207	1028303207	1277100178	1581984305	2001251142	2101626140-156	300097439
100181373	1012134207	1028304207	1277150178	1600230308	2001252142	2200850164	300097539
100181473-74	1012135207	1028305207	1277210178	1630233308	2001253142	2200851164	300097639
100181573	1012201207	1040913169	1278140178	1640234308	2001254142	2200852164	300097739
1001926169 1002010176	1012202207 1012203207	1040919169 1040929169	1394855172 1501103305	1650000305 1660223308	2001255142 2001256152	2200853164 2200854164	300098079 300098179
1002012176	1015161207	1040925169	1501106305	1681154305	2001257152	2200855164	300100152
1002014176	1015162207	1042044230	1501113305	1740100307	2001258152	3000156234	300100252
1002126181	1015163207	1042066230	1501116305	1750224307	2001259152	300035635	300100494
1002128181	1015164207	1051837218	1501123305	1760200307	2001260152	3000389100	300101039
1002280234	1015165207	1052010219	1501126305	1770001307	2001406147	300040086	300101140
1002290234	1017181207	1052012187	1501133305	1770002307	2001477163	300040186	3001172157
1002310176	1017182207	1052014219	1501136305	1770003307	2001615151	300040287	300120896
1002312176	1017183207	1052016187	1501206305	1770004307	2002371154	3000416224	3001214107
1002314176	1017184207	1052022187	1501223305	1780225307	2002372154	300043538	300122136
100241863-194 100241963-194	1017185207 1019160181	1052026187 1052050216	1501226305-307 1501306305	1800400309 1800401309	2002961146 2002962146	3000438216 300044535	300125083 300125183
100242063-194	1019170181	1052100216	1501326305	1800402309	2002971146	3000447229	3001295105
100242670	1019180181	1087102176-197	1501406305	1800404309	2002972146	3000448229	300130049
1002429169	1019201207	1087122176-197	1501416305	2000002138-140	2003721146	3000449229	3001373105
1002521184	1019202207	1087142176-197	1501606305	2000003138-140	2003741146	3000450229	300142394-290
1002531100-186	1019203207	1087162176-197	1501706305	2000007146	2005163144	300045992-217	3001452290
1002532100-186	1019204207	1087182176-197	1501806305	2000009138-140	2005164144	300051296	300148895
1002801170	1019205207	1087202176-197	1501816305	2000011141-151	2005165144	300051396	300172335
1002802170	1019260181	1090100296	1501946305	2000012141-151	2005166144	300051496	3001752298
1002874175 1002875175	1019270181 1019830182	1090150296 1090250296	1501956305 1501963305	2000013141-151 2000015141-151	2005167144 2005168144	300051596 3000543102-288	300182035 300221535
1002876175	1019840182	1090360296	1501966305	2000015141-131	2005169144	3000544102-288	3002301234
1002932169	1019850182	1091100181	1501976305	2000021141-151	2005170144	3000545103-298	3002302234
1004012186	1019860182	1092050296	1501986305	2000022141-151	210104077	3000546103-288	300251235
1004024186	1021210181	1092100181	1561401305	2000023141-151	210104177	3000547103-288	300301135
1004540169	1021212181	1093100296	1561521305	2000024141-151	210118177	3000548104	3003141227
1004713199-209	1021215181	1095060296	1561811305	2000025141-151	210118277	3000549104	3003142227
1004715 199-209	1024132193	1095100296	1561941305	2000030148	210127073	3000550105	3003145227
100477163-194 100477263	1024148193 1024188193	1095150296 1095210179	1561951305 1561961305	2000063154 2000064154	210127173 210127273	3000551105 300058787	3003178231 3004192229
1005207219	1024261207	1095250296	1561971305	2000071147	210127373	300058887	3004195229
1005208219	1024262207	1095360296	1561981305	2000072144	210127474	300061796	300525035
1005209219	1024263207	1096050296	1581104305	2000073144	210127574	300068396	300610035
1008110175	1024264207	1096100296	1581105305	2000074144	210127875	3000718234	3009100106-296
1008112175	1024265207	1096120178	1581114305	2000075144	210127975	3000778107	3009200296
1008114175	1025115194	1096150178	1581124305	2000081138	210128075	3000788198-214	3031410227
1008116175	1025121194	1096190179	1581134305	2000091138	210128276	3000833262	3031420227
1009100296 1009160181	1025124194 1025412187	1097050296 1098302182	1581204305 1581205305	2000095144 2000200141	210128377 210128477	3000834262 3000835262	3031430227 3031450227
1009170181	1025424187	1098402182	1581209305-307	2000200141	210128576	3000836262	3031470228
1009180181	1025500185	1102100179	1581304305	2000205151	210149781	300083796	3031471228
1009260181	1025512185	1102120179	1581305305	2000206151	210149881	300083997	3031472228
1009270181	1025524185-310	1102140179	1581404305	2000207151	210149981	300085035	3031473228
1009280181	1025548185	1102160179	1581405305	2000208141	210150081	300086597	3031474228
1009830182	1025812185-291	1102200179	1581524305	2000209141	210150174	300086697	3031710231
1009840182	1025824185	1102240179	1581525305	2000210141	2101502160	300086797	3031715230
1009850182	1025848185	1102280179	1581604305	2000238154	2101503160	300086897	3041910229

Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page
3041920229	400175659	4120127247	5000051285	505092141-205	5150503300	5303125300	5310506201
3049110229	400175759	4120129248	5000052285	505093241-205	5151001300	5303150300	531201037
3049120229	400175860	412013066	5000750198	505093341-205	5151002300	5303180300	531201137
3049150229	400213660	412013160	5001613284	505093441-205	519541048	5304070300	531201237
3314100227	400241664	4120160248	5001614284	505093541-205	519542048	5304090300	531201337
3314200227	400242064	4120161248	5001615284	505093641-205	519542148	5304110300	531201437
3352104297	400242170	4120162248	5001616284	505093741-205	519542248	5304125300	531201537
3352204297	400242270	4120169247	5001617284	505093841-205	525203754	5304150300	531201637
3352604297	400242770	4120174247	5001618284	505093941-205	525203854	5304180300	531201737
3491100229	400242870	4120177247	5001619284	505095041-205	525203954	5305070300	531209442
3491200229	4002430255	4120183244	5001620284	505095141-205	525204054	5305090300	531209542
3521015297	400251662	4120184244	5001621284	505095241-205	525330052	5305110300	531209643
3522015297	400251762	4120185244	5001622284		525340052	5305125300	531209743
				505095341-205			
3523015297	400251862	4120186244	5001623284	505095441-205	525350052	5305150300	531210541-42
4000200231	4002627253	4120187244	5001624284	5057022166-174	525360052	5305180300	531210641-42
4000474150	4002628158	4120189244	5001625284	505960044	525390052	5306070300	531210741-42
4000475150	4002629154	4120190244	5001626284	505960144	5291047301	5306090300	531210841-42
4000476150	4002841257	4120191244	5001627284	505960244	5291070301	5306110300	531210941-42
4000490217	4002842257	4120192244	5001628284	505960344	5291090301	5306125300	531211042
4000501217	4002851254	4120193244	5004000310	5100014241	5291110301	5306150300	531211142
4000507252-261	4005071261	4120194244	5004001310	5100015241	5291125301	5306180300	531211242
4000508252-261	4005081261	4120220248	5004002310	5100016241	5291150301	5310001209	531211341-42
4000509252-261	4005091261	4120222248	5004003310	5100017241	5293047301	5310002209	531212042
4000600258	4042300252-253	4120223248	5004004310	5100018241	5293070301	5310003209	5313000278
4000601258	410004168-211	4120300246	5004005310	5100019241	5293090301	5310004209	5313003272
4000624258	410004468-211	4120400246	5004006310	5100020241	5293110301	5310005209	5313004272-273
4000629251	410004568	4120401246	5004007310	5100021241-260	5293125301	5310006209	5313007272-273
4000630251	410004668	4120410248	5004008310	5100022241	5293150301	5310007209	5313009272
4000631251	4100050211	4120411248	5004009310	5100023241-260	5297001301	5310008209	5313013270
4000632158	4100051211	4120500246	5004010310	5100024241	5297002301	5310009209	5313015270-271
4000633199-252	4100057211	4120510248	5004011310	5100025241	5297003301	5310003209	5313016271-278
4000634233-258	410007568-211	4120600247	5004012310	5109000238	5297004301	5310011209	5313017273
4000635158	410010068-211	4120700244	5004013310	5109001238	5297005301	5310012209	5313018273
4000638261	410020068-211	4120800299	5004014310	5109002238	5297006301	5310013209	5313019273
4000639261	4100519218	4120801299	5004015310	5109003238	5297007301	5310014209	5313021275
4000640261	410171268-211	4120802299	5004016310	5109004238	5297008301	5310015209	5313022270-278
4000643261	4120009243	4120803299	5004017310	5109100110	5297009301	5310016209	5313026274
4000644261	4120010243	412905049-299	5004018310	5109101110	5297010301	5310017209	5313027272
4000645261	4120013244	4312004239	5004019310	5109102110	5297011301	5310018209	5313028272
4000648259	4120018241	4312006240	5004020310	5109103110	5297012301	5310019209	5313029276
4001045252-253	4120019260	4312007239	5004021310	5109104110	5297013301	5310020209	5313034272-273
4001047251	4120020240	4312015243	5004022310	5109105110	5297014301	5310021209	5313035272-273
4001048251	4120021240	4313030244	5004023310	5109106110	5297015301	5310022209	5313044271
4001049251	4120025239	4313040260	5004024310	5109107110	5297016301	5310023209	5313045271
4001050252	4120026239	4512009243-256	5004025310	5109108110	5297017301	5310024209	5313046273
4001051252	4120029256	4903001263	5004026310	5109200109	5297018301	5310025209	5313055275
4001052252	4120030241	4903002264	5004027310	5109997244	5297019301	5310026209	5313060275
400121963	4120031241	4903004263	5004028310	5109999244	5297020301	5310027209	5313061275
4001425258	4120032241	4903005264	501000082-267	5110000245	5297021301	5310028209	5313062275
4001473259	4120032260	4903006263	501000182-267	5110020242	5297022301	5310029209	5313063275
4001489149							
	4120034260	4903007264	501000282-267	5110021242	5297023301	5310030209	5313064276-277
4001490148	4120101247	4903008264	501000382-267	5110022242	5301070300	5310031209	5313065276
4001538261	4120102247	4903010265	501000482-267	5110023242	5301090300	5310032209	5313066276
4001610230	4120103247	4903011265	5011148198	5110024242	5301110300	5310033209	5313067276-277
400161153-230	4120104247	4903012265	5014415198	5110025242	5301125300	5310034209	5313068276-277
400161253-230	4120107247	4903013265	5014449198	5110026242	5301150300	5310035209	5313069277
4001613259	4120108247	4903014265	502010073-191	5110027242	5301180300	5310036209	5313070275
4001626259	4120109247	4905000224	502010173-191	5110028242	5302070300	5310100252-253	5313071275
4001676158	412011364-295	4905001224	502010373-191	5110029242	5302090300	5310101252-253	5313072275
4001724252-253	4120114247	5000011197	5028192200	5110033242	5302110300	5310102202	5313073275
400172766	4120117241	5000041200	5028193200	5120210281	5302125300	5310500201	5313074275
400172963-265	4120120247	5000042163-200	5028194200	5120230288	5302150300	5310501201	5313075272
400173063-265	4120121247	5000044200	5028195200	5130501300	5302180300	5310502201	5313076272
400173467	412012366	5000045200	5028196200	5130502300	5303070300	5310503201	5313077272
400173567	4120125247	5000046200	5050111298	5150501300	5303090300	5310504201	5313078271
400174558	4120126247	5000050200	5050903298	5150502300	5303110300	5310505201	5313081271

Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page
5313082271	Part No Page 5600022292	Part No Page 5609511285	Part No Page 581000345	Part No Page 5838137269	Part No Page 5900499133-295	Part No Page 591010050	600061940
5313083271	5600023292	5609512285	581000445	5838138269	5900500295	591010150	600062040
5313084271	5600024292	5609513285	581000545	5838139269	5900501295	591010651	600070887
5313089271	5600025292	5715202125-132	581000645	5838140269	5900600166-174	591010751	600070987
5313090270	5600026292	5715203132	581000745	5838141269	5900601201	591010851	600071091
5313091270	5600027292	5726000208-268	581010045	5838142269	5900602201	591011051	600071191
5313092270	5600028292	5726001208-268	581010145	5838143269	5900603201	591012051	600071291
5313093270	5600029292	5726004208-268	581010245	5838144269	590061053	591012151	600071394
5313094270	5600030292	5726005208-268	581011045	5838145269	590061153	591012251	600071491
5313095270	5600031292	5726006208-268	5810200191	5838146269	590061253	591012351	6000719234
5313096271	5600032292	5726007208-268	5810201191	5838147269	590061353	591012451	6000720234
5313097271	5600033292	5726008208-268	5810202191	5838148269	590062050	592904154	6000778174
5313098271	5600034292	5726009208-268	5810203191	5838150269	590062150-150	599120352	6000779174
5313099271 5313111272	5600035292	5726010208-268	5810204191	5838151269	590062450	599120452	6000798174
5313113270	5600036292 5600037292	5726011208-268 5726012208-268	5830024198 5830025267	5838152269 5838153269	590062550 590065096	599120552 599120652	6000799174 6000879176
5313115273	5600038292	5726013208	5830026267	5838154269	5900700222	599120752	6000880176
5313119272	5600039292	5726016208	5830027267	5838155269	5900701222	599120852	600089691-167
5313222275-277	5600040292	5726019208	5830031267	5838156269	5900800222	6000008230	6000901173
5313308271	5600041292	5726025208	5830032267	5838157269	5900801222	6000009230	6000902173
5313309271	5600042292	5726113208	5830033267	5838158269	5901000220	600001192	6000903173
5313311271	5600043292	5726116208	5830039267	5838159269	5901001280	6000016172	6000904215
5435030271-278	5600044292	5726119208	5831302269	5838160269	5901002280	6000026172	600099886
5435032278	5600045292	5726125208	5831303269	5838161269	5901003280	600003386	6001092106
5435034272	5600046292	5726213208	5831304269	5838162269	5901006280	6000036167	6001093100-106
5435035278	5600047292	5726216208	5831822269	5838163269	5901007279	6000046173	6001100172
5435036278	5600050291	5726219208	5832302269	5838164269	5901008279	600005587	6001103173
5435038278	5600051291	5726225208	5832303269	5838165269	5901012281	6000070199	6001150172
5490034220	5600052291	5726238208	5832304269	5838166269	5901013281	6000090173	6001160173
5490035220	5600053291	5726413208	5832822269	5838167269	5901016278 5901017278	6000093180 6000094180	600117393 600118689
5490036220 5490037220	5600054291 5600055291	5726416208 5726419208	5838101268 5838102269	5838168269 5838169269	5901020280	6000117167	600118789
5490039220	5600056291	5726425208	5838103269	5838170269	5901021280	6000120172	600118889
5490040220	5600057291	5726513208	5838104269	5838171269	5901022280	600012391	600119486
5490041220	5600058291	5726516208	5838105269	5838172269	5901023280	6000128171-124	600119589
5490042220	5600059291	5726519208	5838106269	5838173269	5901024280	6000135172	600119689
5490044220	5600060291	5726525208	5838107269	5838174269	5901025280	600013788	600119789
5490047220	5600061291	5752395132	5838108269	5838175269	5901981274	600013888	600119889
5490050221	5600065291	5752396132	5838109269	5838176269	5901982274	600014088	6001200172
5490053221	5600066291	5752397132	5838110269	5838177269	5901983274	600014188	600120996
5490054221	5600067291	5752398132	5838111269	5838178269	5901984275	600014288	600123689
5490055221	5600068291	5752399132	5838112269	5838179269	5901985275	6000170172	600123789
5490056221	5600069291	5752401132	5838113269	5838180269	5901986275	6000190173	600123889
5490057221	5600070291	5752421132	5838114269	5838181269	5901988276	6000220172	600123989
5490058221 5490059221	5600071291 5600072291	5752561132 5752571132	5838115269 5838116269	5838182269 5838183269	5901990277 5902008275	6000235172 6000270172-227	6001278180 6001300172
5490060221	5600073291	5752951132	5838117269	5838184269	5902014277	6000278172-227	600140091-106
551503050	5600074291	5752955132	5838118269	5838185269	5903023294	6000291182	6001421107
551503150	5600075291	5753201132	5838119269	5838186269	5903025294	600038588	6001422107-299
5600000268	5600076291	5753205132	5838120269	5838187269	5903029295	600038688	6001423107
5600001291	5600396215-278	5753241132	5838121269	5838188269	5903030294-299	600038788	600146292
5600002291	5600397215-278	5753243132	5838122269	5838189269	5903049265	600038888	600147292
5600003291	5600398215-278	5753401132	5838123269	5838190269	5904009220	600038988	600148292
5600004291	5600399215-278	5753411132	5838124269	5862301269	5904010220	6000390100-106	6001491107-299
5600005291	5600400215-278	5753501132	5838125269	5862302269	5904011220	6000391100-288	600160391
5600006291	5601050190	5753512132	5838126269	5862303269	5904012220	6000392100-106	6001800106
5600007291	5601051190	5800038210	5838127269	5862350269	5904013133-220	6000393100-106	6001901173
5600008291	5601052190	5800042210	5838128269	5862401269	5904014133-220	600041041-174	6001902173
5600009291	5601053190	5800044210	5838129269	5862402269	5904015220	6000411174	6001903173
5600010291 5600011291	5601054190	5800047210 5800049210	5838130269 5838131269	5862403269 5862450269	5904017133-220	6000482173 6000534173	6002100172 6002150172
5600012291	5601055190 5609401282	580004983	5838132269	5862450269 5900010220	5904018220 5905006220	6000535173	6002200172
5600013291	5609402282	580006183	5838133269	5900019216	5905027221	6000567182	6002268174
5600014291	5609404282	5810000145	5838134269	5900020220	5905028221	6000568182	6002278174
5600020292	5609405282	581000145	5838135269	5900100220	5905029221	6000569182	6002300172
5600021292	5609510285	581000245	5838136269	5900200220	591000951	6000570182	6002314236

Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page	Part No Page
6002316236	603201193	7000532170	700083731	7001393131	7001704116	7002576119	8029050184
6002868175	603201290	7000536176	700083831	7001395128	7001705120	7003208108	8080207183-215
600289141-205	603201390	7000538176	7000846199	7001396128-131	7001706121-122	7003611197	8080208183
600289241-205	6051622218	700054086	7000848199	7001397128-131	7001707120-122	7003612196-197	8080801183
600289341-205	6051635218	700054186-192	7000856199	7001398128-131	7001708120-122	7003613197	8080802183
600289441-205	6051645218	7000550118	7000858199	7001399128-131	7001709120-122	7003614197	8080803183
600289541-205	6051690218	7000552128	7000910132	7001400128-131	7001710120-122	700400186	8080804183
600289641-205	6077812174	7000560126-127	700091335	7001401128-131	7001711120	700400286-192	8080904184
600289741-205	6099120179	7000561127	7000929116	7001403132	700172335	7004031169	8081810183
600300047	6099150179	700057286	7000930120-122	7001407122	700172535	7005000132	8081812183
600301047	6099180179	700058987	7000931120-122	7001408121	700173127	7005010176	8081814183
600323191	6099210179	7000600196	7000934125	7001443125	7001768121	7005012176	8081815183
600323291	6100120179	7000609119-123	7000935125	7001444124	7001769121	7005014176	8081817183
600323491	6100150179	7000657196	7000937195	7001445124-125	7001770125	7005211167	8085301183
600323691-96	6100180179	7000658196	7000938127	7001446124-125	7001771125	7005310176	8085302183
6003294230	6100210179	7000661196	7000939127	7001447124-125	7001772125	7005312176	8085303183
6003296230	6121312178	7000663196	700095132	7001448124-125	7001773125	7005314176	8085504183
6003901106	6121315178	7000664196	700095232	7001449124-125	7001774125	7008010175-197	8085505183
6003921106	6155100235	7000666196	700095332	7001450124	7001775125	7008012175-197	8085514183
6003931106	6156100235	7000700117-122	700095432	7001451128	7001776125	7008133166-172	8085515183
600459192-217	6157100235	7000706289	7001000132	7001474108	7001777125	7008253166-172	8085530184
600512196	6200255173	7000714108	700100527	7001475108	7001778125	7008412199	8085531184
600512296	7000004199	7000715108	700100627	7001476124	7001779120-122	7008512199	8085532184
600512396	7000004233	7000716108	7001017289	7001477124	7001780120-122	7009079123-125	8092902184
600512496-97	7000014114	7000717108	7001018289	7001485126	7001781120-122	7009133125	8092905184
600513196	7000015114	7000731123-125	7001019289	7001486126	7001782120-122	7009223196	8129010184
600513296	7000015114	7000731123-123	7001013289	7001490130	7001783125-128	7012019195	8129020184
600513396	7000017114	7000733125	7001020289	7001496289	7001784125-128	7022810175	8129030184
600513496-97	7000018114	7000734119-123	7001021289	7001500114-132	7001785125-128	7022812175	8129040184
600514196	7000025167	7000735119-123	7001065118	700151130	7001786125-128	7022814175	8129050184
600514296	70000244195	7000736119-123	7001066118	7001539196	7001787125-128	7031220108	8929010184
600514396	7000130168	7000737119-123	7001067118	7001546196	7001788125-128	7031612108	8929020184
600514496-97	7000131168	7000738119-123	7001068118	7001547196	7001789125-128	7041921229	8929030184
600515196	7000134195	7000739119-123	7001070118	7001551195	7001790125-128	7041922229	8929040184
600515396	7000140168	7000740119-123	7001072118	7001561109-243	7001791125-128	7041951229	8929050184
600683196	7000141168	7000741124-125	7001073120-122	7001568260	7001792125-128	7041952229	899120052
600683296	7000160168	7000742124-125	7001074120-122	7001569109	7001793125-128	7051720218	899120152
600683396	7000161168	7000743124	7001075128	7001570109	7001794125-128	7051725218	899120252
6007712174	7000181167	7000745127	7001088132	7001573118	7001795125-128	7051735218	
6007912174	7000208108	7000746127	7001106195	700160930	7001796125-128	7051740218	
600837196	7000225196	7000747127	7001107195	700161031	7001797125-128	7094519195	
600837396	700024228	7000748127	7001108195	7001618108	7001798125-128	7121431223	
6009730182	700024329	7000749127	7001109195	700163827	7001799125-128	7204055166	
6009740182	7000251167	7000750127	7001110195	7001645125	7001800125-128	7419101229	
6009750182	7000256108	7000751127	7001111195	7001646125	7001801122	7419102229	
6009760182	7000292199	7000752127	7001112195	7001647119-123	7001803128	7419201229	
6010140184	7000293199-258	7000753127	7001113195	7001648119-123	7001804128	7419202229	
6010252184	7000295177	7000754127	7001224108	7001650126-127	7001805128	7462200108	
6014204230	7000306108	7000755127	7001239289	7001651126-127	7001806128	7471200108	
6014206230	7000346108	7000756127	7001256108	7001652126-127	7001807128	8000140184	
6019730182	700038436	7000763128	7001315120-122	7001653126-127	7001808128	8000250184	
6019740182	700038540	7000764128	7001316120-122	7001654127	7001809128	8000252184	
6019750182	7000403168	7000767132	7001317120-122	7001655128	7001810128	8000500184	
6019760182	7000404169	7000768132	7001318120-122	7001656128	7001811128	8001000184	
6021312178	7000405169	7000769132	7001319120-122	7001657127	7001812128	8002000184	
6021315178	7000406169	7000770132	7001320120-122	7001659127	7001813118	8002902184	
6022610174	7000408166-174	7000771132	7001321120-122	7001660127	7001814118	8002905184	
6022612174	700044629-293	7000772132	7001322120-122	7001661127	7001815118	8003000184	
6022613174	7000483216	7000773132	7001323121-122	7001662127	7001816118	8004000184	
6022710174	7000489196	7000793120-122	7001324121-122	7001670121-122	7002240114	8005000184	
6022712174	7000497215-235	7000794120-122	7001325121-122	7001671121	7002356117	8012902184	
6023141236	7000498215-235	7000795120-122	7001326120	7001672121	7002357116	8012905184	
6023161236	7000506176	7000796120-122	7001387127	7001673121	700243129	8029010184	
6028610175	7000508176	7000797120-122	7001388114	7001700120-122	7002510234	8029020184	
6028612175	7000521167	7000798120-122	7001391131	7001702121	7002556118	8029030184	
6028614175	7000531166	7000808175-197	7001392131	7001703116	7002560117-133	8029040184	



Magnetic stirrers 27 to 34 pages **Rotator, tilting and vortex mixer shakers 35** page **Rocking Mixer shakers** 36 page **Orbital and linear motion shakers** 37 to 40 pages **Overhead stirrer** 41 page **Orbital shaker incubators** 42 to 43 pages **Ultrasonic homogenizer** 44 page **Grinding and mixer homogeneizer** 45 page



STIRRERS Stirrer and shaker summary table

QUICK OVERVIEW

								φυισι σ	VERVIEV	
				M	AGNETICS					
Single Direction	Part No.	Controller	Volume max.	Heater	Controller	Range	Range	Rapid	Plate Ø	
Stirring			litres		Туре	temperature °C	r.p.m.	Overview	cm	
AGIMIN	7001731	Electronic	1	No	-	-	40 / 2000	ABS	10	
AGIMICRO	7001638	Analogue	2	No	-	-	40 / 1700	ABS	12	
AGIMATIC-S	7000242	Analogue	10	No	-	-	60 / 1600	INOX	14.5	
AGIMATIC-N	7000243	Analogue	10	Yes	Electronic/Energy	350	60 / 1600	INOX	14.5	
AGIMATIC-E	7002431	Analogue	10	Yes	Electronic/Temperature	350	60 / 1600	INOX	14.5	
AGIMAN	7000446	Analogue	2	Yes	Electronic/Energy	350	60 / 1600	For Heatin	g Mantle	
AGIMATIC-ED	7001511	Digital	10	Yes	Electronic/Temperature	350	60 / 1600	INOX	14.5	
AGIMATIC-ED-C	7001609	Digital	10	Yes	Electronic/Temperature	350	60 / 1600	CERAMICS	14,5	
Large volume stirring										
AGIMATIC-HS	7001005	Digital	25	No	-	-	50 / 980	INOX	28 x 34	
AGIMATIC-HL	7001006	Digital	60	No	-	-	50 / 980	INOX	44 x 38	
Reversible stirring action										
AGIMATIC REV-TF1	Г 7000837	TFT touch scre	en 10	Yes	Electronic/Temperature	350	50 / 1500	INOX	14.5	
AGIMATIC REV-TFT-	C 7000838	TFT touch scre	en 10	Yes	Electronic/Temperature	350	50 / 1500	CERAMICS	14,5	
Multi- positions										
MULTIMATIC 5-S	7000951	Analogue	5 x 0.6	No	-	-	1600	INOX	28 x 33	
MULTIMATIC 9-S	7000952	Analogue	9 x 0.3	No	-	-	1600	INOX	28 x 33	
MULTIMATIC 5-N	7000953	Analogue	5 x 0.6	Yes	Electronic/Energy	300	1600	INOX	25 x 25	
MULTIMATIC 9-N	7000954	Analogue	9 x 0.3	Yes	Electronic/Energy	300	1600	INOX	25 x 25	
				ROTARY 8	TILTING MIXERS					
Tube Mixers	Pa	rt No.	Vol. m	ax. litres	Heaters	Ra	Range r.p.m.		Capacity	
ORBIT	300	00445	5 x	0.25	No		12		Selection	
MOVIL ROD	700)1723		-	No	10	1 / 80	35 mm (
MOVIL TUB	700	00913		-	No		15	12 Tul	jes	
VORTEX-VIB	700	1725		-	No	2	800	See Tube S	election	
				VIBRA	ATING SHAKER					
Oscillating	Part No	o. Ma	x. Load	Heater		Capacity		08	scil. min.	
VIBROMATIC	70003	84	1 Kg	No	8 flasks of 50 mm Ø ma	x. (optional +	Erlenmeyer 25	0 ml) 10	0 / 950	
					O - FRO SHAKERS					
To & Fro	Part N	0.	Heater			nge	Range	PI	atforms	
& orbital movement						ature °C	r.p.m.		cm	
ROTATERM	30004		Yes	Electr	onic/Digital	.250	20/230		11x31	
ROTABIT	30009	974	No		-	-	20/230	According to a	ccessorie	





Shaker for Petri plates "AG-200", shaker for bottles "AG-200-A", and shaker for Erlenmeyer flasks "AG-200-B".

Thermo shaker "TRM-4" for Microtiter.

Incubated / refrigerated high volume orbital shakers "D-110", "D-210", "D-1102" and "D-2102".

Ultrasonic homogeneizers "CY-500".



A new generation of stirrer for the most exacting work



Magnetic Stirrer "Agimin" without heating

MICROPROCESSOR CONTROLLED, ELECTRONIC CONTROL, PORTABLE.

FEATURES

Maximum stir volume 1 litre. Speed controllable from 40 to 2000 r.p.m. with storage of the ultimate speed setting.

Double power source: 4 x 1.5V AA batteries or with a mains adapter.

Moulded ABS plastic case.

Excellent resistance to chemical agents.

CONTROL PANEL

Luminous On/Off push button. Luminous r.p.m. push button. Increase/ decrease speed push button. Low battery indicator lamp.

MODEL

Part No.	Height	Ø	Depth	Power	Weight
	cm	cm	cm	W	Kg
7001731	5	10	15	0.9	0.35





Comes complete with support bar and clamp for 3 electrodes (1 x 13 mm \emptyset and 2 x 16 mm \emptyset), 4 batteries 1.5V AA, mains adapters and 2 PTFE coated stir bars.



Magnetic stirrer "Agimicro" without heating

FEATURES

Maximum stir volume: 2 Litres.
Adjustable stir speed from 40 to 1400 r.p.m.
Durable moulded ABS plastic body.
Excellent chemical resistance.

MODEL

Part No.	Height	Ø Plate	Power	Weight
	cm	cm	W	Kg
7001638	5	12	2	0.5

Comes complete with a 8 \emptyset x 32 mm P.T.F.E. coated stir bar.

CONTROL PANEL

Mains switch with luminous "on" indicator. Analogue control of speed in

Analogue control r.p.m.

ACCESSORIES

Support bar and electrode clamp (1 of 13 and 2 of 16 mm Ø). Used with pH and conductivity meters.
Part No. 7001639







Powerful magnetic stirrer "Agimatic-H"without heating

FOR LOW VISCOSITY VOLUMES FROM 1 TO 60 LITRES.
USING REDUCED VOLUMES THE UNIT CAN STIR MEDIUM TO HIGH VISCOSITY LIQUIDS.

Maximum stirring power

FEATURES

Electronic speed control from 50 to 980 r.p.m, will maintain a constant stir speed irrespective of the liquid viscosity.

Digital speed display.

AISI 304 stainless steel top plate with an epoxy coated metal case.

CONTROL PANEL

Main switch with a luminous "on" indicator. Electronic speed control. Digital speed display in r.p.m.

MODEL

	Part No.	Speed control r.p.m.	Maximum stir volume in litres	Height/Width/Depth (exterior) cm	Power W	Weight Kg	
AGIMATIC-HS	7001005	50 to 980	up to 25	12 28 34	35	8	
AGIMATIC-HL	7001006	50 to 980	up to 60	14 44 38	40	12	

Model HS comes complete with a 16 Ø x 127 mm stir bar and model HL with a 19 Ø x 150 mm stir bar.






Magnetic Stirrers "Agimatic" and "Agiman"



A new generation design of stirrers at the forefront of technology

COMMON FEATURES

Unique functional design, with epoxy covered alloy case. Choice of stainless steel or ceramic top plate (model dependent). An internal case acts as a spill guard and insulator, made from AlSl304 stainless steel. Tilted control panel for easy viewing with splash protection. Location hole at the back for a retort support rod, supplied as an accessory. Electronic speed control from 60 to 1600 r.p.m. and from 50 to 1500 r.p.m. for models "Agimatic Rev" with constant speed irrespective of liquid viscosity. Easy to service and to replace modular electronics. All heated models have a connection at the back to connect an electronic contact thermometer "Sensoterm II", electronic controller "Electemp-TFT" or Pt100 temperature probe, model dependent. Agimatic and Agiman exterior measures: Height 14.5, Width 19 and Depth 26 cm. Weight: 3,2 Kg. All hotplates have a safety over temperature cut out, set to the maximum working temperature.





Magnetic Stirrers "Agimatic" and "Agiman" analogue control



FOR SPEEDS FROM 60 TO 1600 rpm.. Maximum Stir Volume: 10 Litres.

PLATE MADE FROM AISI 304 STAINLESS STEEL WITH REMOVABLE SAFETY RING. PLATE Ø 14.5 cm.



"Agimatic-S" without heating

CONTROL PANEL

- 1. Mains switch illuminated.
- 2. Electronic stir speed control.



MODEL

	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
AGIMATIC-S	7000242	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.





CONTROL PANEL

- 1. Switch on button with indicator lamp for heating and for stirring speed.
- 2. Electronic stir speed control.
- 3. Electronic heater energy control from 0-100%, synchronised with an indicator lamp that functions when the plate is hot.
- 4. Illuminated indicator when the heater is "on".



MODEL

	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
AGIMATIC-N	7000243	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.





"Agimatic-E" with heating
ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER. FROM 50 °C TO 350 °C PRECISION ±2 °C.

CONTROL PANEL

- 1. Switch on button with indicator lamp for heating and for stirring speed.
- 2. Electronic stir speed control.
- 3. Electronic temperature control of the heater plate from 50 °C to 350 °C synchronised with an indicator lamp that functions when the plate is hot.
- 4. Illuminated heater indicator when the heater is "on". Precision ±2 °C.





MODEL	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
AGIMATIC-E	7002431	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.





Magnetic stirrer "Agiman"

ELECTRONIC CONTROL OF SPEED AND POWER TO THE HEATER. WITH A MECHANISM FOR ATTACHING HEATING MANTLES.

CONTROL PANEL

Switch on button with indicator lamp for heating and for stirring speed.

Electronic stir speed control.

Electronic heater energy control from 0-100%, synchronised with an indicator lamp that goes on when the mantle is hot.

Illuminated heater indicator when the mantle is "on".

ACCESSORIES

Attachable heating mantle for the "Agiman" magnetic stirrer.

Part No. Capacity flask ml 3000447 250 3000448 500 3000449 1000 3000450 2000

Complements for Agimatic-S-N-E and Agiman. See accessories on page 34 and 35.

MNNFI

MODEL	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Weight Kg	
AGIMAN	7000446	un to 2	14.5 19 26	3.5	

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.





Magnetic Stirrers "Agimatic-ED" and "Agimatic-ED-C" digital control



FOR SPEEDS FROM 60 TO 1600 RPM.

FOR ADJUSTABLE TEMPERATURES FROM 50 °C TO 350 °C. PRECISION ± 2 °C.

MAXIMUM STIR VOLUME: 10 LITRES.

PLATE Ø IS 14,5 cm.

WITH RENOVABLE SAFETY RING AGAINST SPILLAGE.



"Agimatic-ED" with heating


ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.

SILUMIN UPPER PLATE COVERED BY AISI 304 STAINLESS STEEL.

CONTROL PANEL

- 1. Digital display indicating stir speed or temperature.
- 2. Push button, parameter selector, Stir speed or
- 3. Push button, increase parameter value.
- 4. Push button, decrease parameter value.
- 5. Indicator lamp, safety thermostat "on"
- 6. Push button Stop-Start

Precision ±2 °C.





MODEL

Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
AGIMATIC-ED 7002431	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.





"Agimatic-ED-C" with heating

ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.

SILUMIN UPPER PLATE WITH CERAMIC SURFACE.

CONTROL PANEL

- 1. Digital display indicating stir speed or temperature.
- 2. Push button, parameter selector, Stir speed or Temperature.
- 3. Push button, increase parameter value.
- 4. Push button, decrease parameter value.
- 5. Indicator lamp, safety thermostat "on"
- 6. Push button Stop-Start

Precision ±2 °C.



MODEL

	Part No.	Maximum stir volume litres			Weight Kg
AGIMATIC-ED-C	7002431	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

Ceramic surface:

Excelent conducitity and heat distribution. Unalterable surface.

Higher acceleration to reach temperature.

Extreme resistance against chemical agents and scratches.





Complements for Agimatic-ED and Agimatic ED-C. See accessories on pages 34 and 35.



Reversible magnetic Stirrers "Agimatic-Rev-TFT" and "Agimatic-Rev-TFT-C"



Forefront technology with TFT touch screen

FOR SPEEDS FROM 50 TO 1500 rpm.

FOR ADJUSTABLE TEMPERATURES FROM 50 °C TO 350 °C. PRECISION ± 2 °C.

MAXIMUM STIR VOLUME: 10 LITRES. PLATE Ø IS 14,5 cm.

MAGNETIC REVERSIBLE ACTION STIRRERS, WITH ADJUSTABLE TIME STIRRING PERIODS AND PRECISE SPEED CONTROL, PROVIDES THE FOLLOWING ADVANTAGES FOR REPEAT PROCESSES:

A.- FASTER MIXING, REDUCED STIRRING TIME. B.- CONTROLLABLE VORTEX EFFECT. C.- QUICKLY ELIMINATES AIR FROM WITHIN THE SAMPLE MIXTURE, DEGASSING. D.- MINIMUM SPLASHING, NO LOSS OF SAMPLE.

SELECT MOVEMENT FOR ONE DIRECTION OR REVERSE.

ELECTRONIC CONTROL OF THE FREQUENCY INVERTING PERIOD, FROM 10-30 SECONDS.

MEMORY UP TO 10 WORKING PROGRAMS. AUTOMATIC CYCLES UP TO 6 SEGMENSTS WITH TEMPERATURE PROGRAMMING, RPM AND TIME.

REAL TIME GRAPHICAL REPRESENTATION OF AUTOMATIC CYCLES.



"Agimatic-Rev-TFT" with heating

ELECTRONIC DIGITAL CONTROL OF STIR SPEED AND HEATER TEMPERATURE. WITH REMOVABLE SAFETY RING AGAINST SPILLAGE.

SILUMIN UPPER PLATE COVERED BY AISI 304 STAINLESS STEEL.

CONTROL PANEL

- 1. Probe position indicator.
- 2. rpm, temperature and time reading.
- 3. Start up.
- **4.** Stop.
- 5. Settings.
- 6. Backward.
- 7. Increase/decrease parameters value.
- 8. Cursor.

MODEL





	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
Agimatic-Rev-TFT	7000837	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.



"Agimatic-Rev-TFT-C" with heating

ELECTRONIC DIGITAL CONTROL OF STIR SPEED AND HEATER TEMPERATURE. WITH REMOVABLE SAFETY RING AGAINST SPILLAGE.

SILUMIN UPPER PLATE WITH CERAMIC SURFACE.

CONTROL PANEL

- 1. Probe position indicator.
- 2. rpm, temperature and time reading.
- 3. Start up.
- 4. Stop.
- 5. Settings.
- 6. Backward.
- 7. Increase/decrease parameters value.
- 8. Cursor



	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
Agimatic-Rev-TFT-C	7000838	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

Ceramic surface:


Excelent conducitity and heat distribution.

Unalterable surface.

Higher acceleration to reach temperature.

Extreme resistance against chemical agents and scratches.









Complements for Agimatic-Rev-TFT and Rev-ED-C see accesories on pages 34 and 35.



Multiple position magnetic stirrers "Multimatic"



COMMON FEATURES

External plate made from AISI 304 stainless steel with engraved stirrer positions.

Independent stirring speed control for each position.

Front panel protected against splashes.

Models 5-S and 5-N: Suitable for 5 flasks, 600 ml max. (max. Ø 90 mm).

9-S and 9-N: Suitable for 9 flasks, 300 ml max. (max. Ø 75 mm).

Epoxy coated external case.

APPLICATIONS

Specially recommended for bio-pharmaceutical diagnostic assays.



Stirrers "Multimatic" 5-S and 9-S without heating

FOR 5 AND 9 POSITIONS.



CONTROL PANEL

Illuminated mains switch.
Individual position stir speed control.

MODELS

	Part No.	Height / Width / D (exterior) cm		Speed control r.p.m.	Maximum stir volume ml	Power W	Weight Kg
MULTIMATIC 5-S	7000951	11 28 3	5 28 33	Up to 1600	5 x 600	16	7
MULTIMATIC 9-S	7000952	11 28 3	5 28 33	Up to 1600	9 x 300	27	9

Comes complete with a 5 or 9 mm stir bar, model dependent covered in PTFE of 8 Ø x 32 mm long.

ACCESSORIES

Transparent methacrylate bath.

Fits onto the non heating stir plate of the "Multimatic" 5-S y 9-S

The bath can withstand temperatures up to 60 °C. See section Thermostats and immersions (page 98).

Part No.	Capacity litres	Height / Width / Depth (exterior) cm	Weight Kg	
1001007	14	16 29 35	4	





Hotplate stirrers "Multimatic" 5-N y 9-N with heating

FOR 5 AND 9 POSITIONS, OVER TEMPERATURE SAFETY THERMOSTAT FITTED.




CONTROL PANEL

Mains illuminated power switch.
Individual speed control for each position.
Electronic heater energy control from 0-100% synchronised with an indicator lamp when the heater is on.

	Part No.	Height / Width/ Depth (exterior) cm	Width / Depth (plate) cm	Maximum temperature °C	Speed control r.p.m.	Maximum stir volume ml	Heater Power W	Total Power W	Weight Kg
MULTIMATIC 5-N	7000953	16 28 35	25 25	300	Up to 1600	5 x 600	680	700	12
MULTIMATIC 9-N	7000954	16 28 35	25 25	300	Up to 1600	9 x 300	680	700	14

Comes complete with a 5 or 9 mm stir bar, model dependent covered in PTFE of 8 \emptyset x 32 mm long.

ACCESSORIES FOR MAGNETIC STIRRERS



Stir-bars followers for magnetic stirrers.

Encased in PTFE with pivot ring, temperature resistant to +275 °C.

Part No.

1000009	6 Ø x 10 mm long.	without pivot ring.
1000012	4 Ø x 12 mm long.	without pivot ring.
1000019	7 Ø x 20 mm long.	with pivot ring.
1000032	8 Ø x 32 mm long.	with pivot ring.
1000039	8 Ø x 40 mm long.	with pivot ring.
1000045	10 Ø x 45 mm long.	with pivot ring.
1001061*	12 Ø x 75 mm long.	with pivot ring.
1001062*	16 Ø x 102 mm long.	with pivot ring.
1001063*	16 Ø x 127 mm long.	with pivot ring.
1001064*	19 Ø x 150 mm long.	with pivot ring.
* With pivo	t ring for large volumes.	



Concentric stir-bar with beaker cover "CENTRIMAN".

Encased in P.T.F.E. with an adjustable axis bar made of AISI 304 stainless steel. Suitable for applications with volumes of 100 to 2000 ml.

Soft start stirring and adjustable height prevents the stir-bar jumping inside the vessel.


Conical fitting lid of polypropylene with excellent chemical and thermal properties that resists temperatures up to 160 °C. The lid has an auto-lubricating central axis hole with two lateral air displacement holes.

Part No.

1000900 Concentric stir-bar with conical fitting lid for beakers volumes of 100 ml. to 500 ml. Length 32 mm. 1000901 Concentric stir-bar with conical fitting lid for beakers volumes of 500 ml. to 1000 ml. Length 50 mm.

1000902

Concentric stir-bar with conical fitting lid for beakers volumes of 1000 ml. to 2000 ml. Length 55 mm.



Triangular stir-bar "TRIMAG".

Encased in heat resistant PTFE of temperatures up to 275 °C. useful for "dissolving solids" and stirring silts. Liquids flow upwards as well as sideways. Quiet with excellent stability. Suitable for flat bottomed vessels.

Part No.

1000536	9 side x 35 mm long
1000537	14 side x 40 mm long
1000538	12 side x 50 mm long
1000539	14 side x 55 mm long





Starhead stir bar.

Circular, encased in heat resistant PTFE for temperatures up to 275 °C. Useful for tubes, flasks and round bottom vessels.

Part No.

25 Ø x 15 mm height. 1000790 1000791 40 Ø x 17 mm height.



Box of 17 assorted stir bars.

Encased in temperature resistant PTFE for temperatures up to 275 °C.

Contents:

1 stir-bar with pivot ring 6x10 mm.

1 stir-bar with pivot ring 4 x12 mm.

1 stir-bar with pivot ring 7x20 mm.

2 stir-bars with pivot ring 8x32 mm.

2 stir-bars with pivot ring 8x40 mm.

2 stir-bars with pivot ring 10x45 mm.

1 triangular stir-bar of 9x35 mm.

1 triangular stir-bar of 14x40 mm.

1 triangular stir-bar of 12x50 mm.

1 triangular stir-bar of 14x55 mm.

1 stir-bar Rotor-Pat of 11x35 mm.

1 oval stir-bar of 16x30 mm.

1 oval stir-har of 20x40 mm.

1 star head stir-bar of 20x15 mm.

1 star head stir-bar of 40x17 mm.

Part No. 1000487



Stir Bar "ROTOR-PAT". Encased in PTFE useable up to temperatures of 275 °C.

Excellent stability with strong vortexing turbulence, useful for round bottom or convex bottom flasks. Part No.

1000001 10 Ø max. x 32 mm long.



Oval stir bar, coated with PTFE, resistant to high temperatures of up to +275 °C. Can be used for all types, including round bottom vessels.

1000788 16 Ø x 30 mm long. 20 Ø x 40 mm long.



Floating magnetic stir-bar "IMANFLOTER".

Encased completely in PTFE Dismountable structure, autoclavable. Designed specifically for low constant speed stirring. Recommended for tissue culture applications. The floating conception of the Imamfloter minimises the grinding effect of cell damage as there is no friction on the bottom of the flask.

Medium speed levels prevent the IMANFLOTER from jumping within the vessel.

Part No. Total length Stir bar length 39 mm 1000897 50 mm 60 mm 1000898 45 mm



Stir-bar retriever. Flexible, encased in PTFE. Part No.

1000020 Length: 280 mm.

ACCESSORIES FOR MAGNETIC STIRRERS AND HOTPLATE STIRRERS



Support rod with thread of 12 mm \emptyset x 500 mm high, stirrer dependent, attaches to the main body of the stirrer.

Part No. 1000725



Clamp with bosshead , "AUTOMATIC" Chrome coated alloy with screw closure.

For Ø from 0 to 25mm.

Part No. 7000251



Electronic contact thermometer with digital readout "Sensoterm II".

Complete with AISI 304 stainless steel sheathed Pt1000 temperature sensor.

FEATURES:

Conceived and designed as external sensor probe, where the controller is away from the measuring temperature zone of heat and vapour. The advantages over a contact thermometer means additional safety, no broken glassware or mercury spills.

SPECIFICATION

Temperature range from: -50 °C to +300 °C. Temperature sensor : 3.5 Ø x 250 mm immersion.

Stability: 1 °C ±1 digit Resolution: 1 °C.

Temperature display: LCD 8 mm.

Working environment: 0 to 60 $^{\circ}\mathrm{C}$ and relative humidity

up to 80 %.

Dimensions: 130 high x 71 wide x 24 mm deep (wi-

thout probe)



Identification ICON

Part No. 1001009



Adapters for round flasks of 250, 500 and 1000 ml. Made of Aluminium alloy. Suitable for the stirrers with a Ø greater than 14.5 cm.

Maximum adapter temperature 250 °C.

Part No.

1001712 250 ml flasks. 1001718 500 ml flasks. 1001719 1000 ml flasks.



Temperature controller Electemp-TFT.

FEATURES:

P.I.D. digital temperature controller. User configurable microprocessor of function parameters.

Temperature range: -150 °C to 1300 °C

Sensor options: Pt100 or Thermo resistor K type probe. See page 293.

SPECIFICATION:

Stability: ±0.1 °C.

Resolution: ±0.1 °C Pt 100 /±1 °C K Type probe.

Operation: ON-OFF/P.I.D.

AUTOTUNING. Adapts automatically to the working me-

dium (Solid, liquid or gas.)

Dimensions: 10.5 high x 13.5 wide x 21cm deep.

Maximum power (at 230 V): 3100 W.

Weight: 1.250 Kg



Identification ICON

Part No. 3000897

ACCESSORIES:

Probe Pt100 Part No. **1000893** Probe K Type. Part No. **1000895**



Aluminium PTFE coated bath.

Suitable for stirrers and hotplate stirrers. "Agimatic N, E. ND. ED and Rev-E".

Suitable as water or oil baths for temperatures up to 200 °C with magnetic stirring from 150 to 1500 r.p.m. For rotary evaporators etc.

Part No.

1001230 capacity 2 litres. 11 cm high x 18 cm \emptyset . 1001231 capacity 4 litres.13 cm high x 24 cm \emptyset .

ACCESSORIES:

Suitable only for tank Part No. 1001230. Set of 4 reduction ring lid, stainless steel AISI304, with location hole for a thermometer.

Part No. 3000917



Temperature sensor Pt100. Controls the sample temperature in the digital hotplate stirrers:

Part Nos: 7001511, 7001609, 7000837 and 7000838.

Immersion Probe: 4mm Ø x 130 mm long. Cable connector 150 cm long.



Identification ICON

Part No. 1001610



Rotator "Orbit"

FEATURES

Rotation Speed: 12 r.p.m.

Adjustable rotation angle axis from vertical to horizontal. An elastic shaft drive and clip supports permit the flasks and tubes to be loaded and unloaded while the motor is running. Epoxy coated outer casing.

Recommended for extractions of cell cultures, blood tubes and other haematological applications.

MODEL

Part No.	Speed r.p.m.		Width terior)	/ Depth cm	Power W	Weight Kg
3000445	Fixed 12	19	23	20	15	6



ACCESSORIES

Erlenmeyer flask adapters, clips made from stainless steel.

Part No. Capacity

3000850 8 x 50 ml Erlenmeyer flasks. 3006100 6 x 100ml Erlenmeyer flasks.

3005250 5 x 250ml Erlenmeyer flasks.

Tube adapters, made from Aluminium, with nickel/chrome clips for different tube \emptyset .

Part No. Capacity

3000356 35 tubes from 6 to 9 mm Ø. **3003011** 30 tubes from 10 to 13 mm Ø.

3002512 25 tubes from 14 to 16 mm Ø.

Part No. Capacity

3002215 22 tubes from 17 to 20 mm Ø. **3001820** 18 tubes from 21 to 24 mm Ø. **3001723** 17 tubes from 23 to 28 mm Ø.



Tilting Mixer Shaker "Movil-Tub"

FEATURES

Fixed speed tilting shaker mixer with removable silicon, anti slide tray. Recommended especially for haematology, chemistry and microbiology applications.

MODEL

Part No.	Tube capacity		Tube ength mm	Speed r.p.m.		/ Widt xterior	h / Depth) cm	Tilt angle	Power W	Weight Kg
7000913	12	20	160	15	12	34	20	25°	6	2,5



Roller and tilt mixer "Movil-Rod"

WITH VARIABLE SPEED CONTROL FROM 10 TO 80 r.p.m.

FEATURES

Roll and tilt mixer with 7 rollers, double action, gently rolls the tubes while rocking, produces a sinusoidal movement, providing excellent sample homogenisation. Recommended for suspensions of biological and haematological fluids.

MODEL

Part No.	For tube	Speed	Length	Height / Width / Depth	Power	Weight
	Ø up to mm	r.p.m.	(rollers) cm	(exterior) cm	W	Kg
7001723	35	10 - 80	24	9.8 34.9 26.5	12	3.6





Mixer shaker for tubes "Vortex-Vib"

SHAKE BY MIXER HEAD.

Strong and economic.

FEATURES

With stirring selector by pressure, or fixed in continuous. Starting switch, 4 interchangeable platforms for multiple tubes and 50 ml flasks. Base with 4 suction cups for full adherence.

MODEL

Part No.	Frequency	Height /Width /Depth	Power	Weight
	r.p.m.	(exterior) cm	W	Kg
7001725	2800	120 170 170	50	4





Optimum shaking frequency: 950 oscillations per minute

APPLICATIONS

Dissolution, extractions and homogenisation.

FEATURES

Electronic control from 100 to 950 oscillations per minute. (1900 vibrations) and timer from 1 to 999 minutes or continuous operation.

Oscillation amplitude: 8 mm.

Easily detachable arms equipped with 8 clamps that allow all types of tubes and flask with diameters of between 5 to 50 mm can be held in position.

Maximum load: 4 kg.

CONTROL PANEL

- 1. Digital display indicator of osc./min. and time.
- 2. Parameter display selector.
- 3. Parameter selector, speed or time.
- 4. Push button increase value.
- 5. Push button decrease value.
- 6. Push button start-stop.



MODEL

Part No.	Oscillations p.m.	Maximum capacity Kg	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7000384	100 - 950	4	23 77 27	34	11



ACCESSORIES



Platform support for Erlenmeyer flasks, made from AISI 304 stainless steel with clip adapter position holes for 50 ml to 1000 ml flasks. Platform dimensions: 38 cm long x 19 cm wide.

Weight: 2 Kg.

Platform capacity:

- 18 Erlenmeyer flasks of 50 ml.
- 12 Erlenmeyer flasks of 100 ml.
- 8 Erlenmeyer flasks of 250 ml.
- 6 Erlenmeyer flasks of 500 ml.
- 3 Erlenmeyer flasks of 1000 ml.

Part No. 3001221



Adapters for Erlenmeyer flasks.

Made from tempered stainless steel.

i ai t ivo.		
1001003	Adapter for	50 ml.
1000978	Adapter for	100 ml.
1000979	Adapter for	250 ml.
1000980	Adapter for	500 ml.
1000981	Adapter for	1000 ml.

2 in 1. Platform support for Erlenmeyer flasks.



Model Vibromatic with adapter platform.



Orbital shaker "AG-200", "AG-200-A", "AG-200-B"





COMMON FEATURES

External ABS case with adapter platform for the different models. Innovative single button to select start, shaking speed, time and stop, Digital display of time and shaking speed.

Maintenance free induction drive motor which produces a uniform movement according to the rotation speed program.

Acoustic notification at the end of the cycle.



Main switch.

Single button selector for shaking speed and time. Display indicating the shaking speed. Display indicating the programmed time.

Measures: Alto 10 x Ancho 28 x fondo 26

TECHNICAL INFORMATION

Rotation frequency: 50-250 r.p.m. Timer from 1 min to 99h 59 min Rotation amplitude: 20mm

Power: 50 W.

Power supply: 220V/50Hz Load capacity: till 2,5Kg

Measures cm: Height 10 x Width 28 x Depth 26



Shaker for Petri plates "AG-200"

With non-slip rubber platform of 24cm width and 20cm depth. Capacity for 4 Petri plates.

Useful for slides and little flasks shaking.

Weight (Kg): 7. Part no: 5312010



Shaker for bottles "AG-200-A"

With non-slip rubber platform of 26cm width and 22cm depth, with 3 adjustable rubber retaining bars suitable for any kind of glasses: Erlenmeyer, flasks, bottles, etc.

Weight (Kg): 8 Part no: 5312011



Shaker for Erlenmeyer flasks "AG-200-B"

Tempered aluminium base of 27cm width and 22cm depth, with multi-adapter location holes to support different capacities flasks from 50ml to 2000ml. Weight (Kg): 8.

Part no: 5312012

Erlenmeyer flasks adapters

Made of stainless steel with upper spring. Part no. **5312013** Adapter for 50 ml. Platform capacity: 12 units

Part no. 5312014 Adapter for 100 ml.

Platform capacity: 8 units Part no. 5312015 Adapter for 250 ml.

Platform capacity: 6 units Part no. 5312016 Adapter for 500 ml.

Platform capacity: 4 units

Part no. 5312017 Adapter for 1000 ml.

Platform capacity: 2 units







Orbital and linear motion shaker "Rotaterm"





WITH HEATED PLATFORM TOP PLATE. ELECTRONIC DIGITAL SPEED AND TEMPERATURE DISPLAY. INCLUDES A SAFETY OVER TEMPERATURE DEVICE. TEMPERATURE STABILITY ±1 °C.

2 in 1 System suitable for orbital or linear motion

FEATURES

Adjustable temperature range from ambient +5 °C. to 250 °C.

Electronic 3 digit temperature display, monitored by a Pt 100 temperature probe.

Resolution: 1 digit.

Electronic digital speed control from 20 to 230 oscillations per minute and timer from 1 to 999 minutes or continuous operation.

Adjustable oscillation amplitude: 15 and 20 mm.

Choice of: orbital or linear motion.

Epoxy covered external casing.

Aluminium top plate platform with evenly distributed heater elements.

The upper platform is made of AISI 304 stainless steel with adjustable tensile spring clips that allow all kinds of glassware, flasks, Erlenmeyer flasks, racks, etc. to be held in position.

Digital calibration temperature circuit.



Part No 3000435 with platform and tensile spring clips.


CONTROL PANEL

Temperature control.

- 1. Digital temperature display.
- 2. Push button increase temperature.
- 3. Push button reduce temperature.
- 4. Push button start-stop heater.
- 5. Over temperature alarm.

Temperature and speed control.

- 6. Digital display of oscillation per minute and time (in minutes).
- 7. Displayed parameter.
- 8. Push button select displayed parameter, speed or time function.
- 9. Push button increase parameter value.
- 10. Push button decrease parameter value.
- 11. Push button start-stop shaker.



11 7 10

MODEL

Part No.	Platform	Width / Depth (platform) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3000435	Spring clip	41 31	30 46 38	2300	32

Diagram of Temperature Vs Time.

For 500ml Erlenmeyer flasks.

A. Plate to 250 °C: Silicon 150 °C. **B.** Plate to 200 °C: $\rm H_2O$ 97 °C.

C. Plate to 150 °C: H₂²0 82 °C.

D. Plate to 100 °C: H₂O 61 °C.

E. Plate to 50 °C: H₂O 37 °C.



Orbital and horizontal shaker "Rotabit"

ELECTRONIC DIGITAL CONTROL OF SPEED AND RUN TIME.





2 in 1 system, versatile orbital and horizontal shaker

FEATURES

Speed control from 20 to 230 oscillations per minute with timer from 1 to 999 minutes or continuous operation. Shaking amplitude of: 15 and 20 mm.

Choice of orbital movement or a horizontal "To and Fro" motion.

Load capacity: till 30 Kg

Easy to interchange platforms, see accessories for platform selection.

Epoxy covered outer casing.

CONTROL PANEL

- 1. Digital display showing oscillations and time.
- 2. Displayed parameter.
- 3. Push button, select parameter function, speed or time.
- 4. Push button increase parameter value.
- **5.** Push button reduce parameter value.
- **6.** Push button Start- Stop.



53:00

MODEL

Part No.	Control oscillations p.m.	•	/ Width / Depth kterior) cm	Power W	Load capacity Kg	Weight Kg
3000974	20 - 230	20	42 43	80	30	13.5

This shaker can be coupled to the "Boxcult" incubation chamber with platform Part No. 3000975 (see page 41).

ACCESSORIES

Universal platforms with non-slip surface with four adjustable silicon covered retaining bars suitable for all types of glassware, Erlenmeyer flasks, etc. to be held in position.

Part No.	Height / Width/ Depth (exterior) cm			Weight Kg
3000975	9	42	35	5
3000976	9	70	35	7

NOTE: This shaker has multiple options and accessories, a flow chart for easy selection of the accessories can be found at the end of this section.



Platform Part No. 3000975 Platform Part No. 3000976

Base for Erlenmeyer flask clips.

Made from AISI 304 stainless steel with multi-adapter location holes to support Erlenmeyer flask clips. Capacities 50, 100, 250, 500 and 1000 ml.

Part No.	For platform number Part No.	Width / Depth (platform) cm	50 ml	Capa 100 ml	city, No. of p 250 ml	oositions 500 ml	1000 ml	Weight Kg
3000977	3000975	40 34	36	25	16	9	5	1.5
3001010	3000976	69 34	66	45	21	18	10	2,2



Erlenmeyer flasks spring clips.

Made with tempered stainless steel.

Part No.

100 1003 Adapter for 50 ml. **1000978** Adapter for 100 ml.

 1000979
 Adapter for 250 ml.

 1000980
 Adapter for 500 ml.

 1000981
 Adapter for 1000 ml.



ACCESSORIES



Orbital To & Fro shaker "Rotabit" complete with accessory. "Boxcult" incubation chamber.

Incubation chamber "Boxcult"

The "Boxcult" can be fitted onto the "Rotabit", Part No. 3000974, with the universal platform, (SMALL) Part No. 3000975.

Forced air fan circulation.

Adjustable temperature, ambient +5 °C to 57 °C.

Homogeneity ±2.5%

Stability ±0.5%

Built of transparent Perspex, access is throught a wide hinged front door which provides easy access for loading and unloading. Samples can be viewed without opening the incubator maintaining a constant homogeneous temperature.

Part No.	•	/ Width nterior)	/ Depth cm	•	/ Width xterior)	n / Depth cm	Power W	Weight Kg
3000957	50	47	47	61	51	51	220	11

See incubator section for more details (page 134).

Set of 5 reduction rings adapters for separating flasks of 50 to 100 ml. Part No. 6000620

Separating funnel device.

Two positions, for mixing Horizontal to Vertical for easy decanting. Capacity: 5 separating funnels of 250 to 500 ml or 2 funnels of 1000 ml. Part No. 6000619

This accessory can only be fitted to the platform Part No. 3000976. The "Rotabit" Part No. 3000974 and platform need to be ordered separately.



Double platform part no 3001011, with shaker Part No. 3000974

Double platform.

Two levels.

Fits onto the "Rotabit" Part No. 3000974.

Complete with anti slide mats and four support bars that support the upper level. Each platform has four silicon clamping bars.


Part No.	Width / Depth (platform) cm	Total combined height cm	
3001011	55 35	40	



Separating funnel decanting accessory Part No. 6000619 Fits onto platform Part No. 3000976.



Flow diagram of the selection of accessories for the "Rotabit" To & Fro Orbital Shaker



SUITABLE FOR ERLENMEYER							
Capacity	Part No.	A	В				
50 ml	1001003	36	48				
100 ml	1000978	25	30				
250 ml	1000979	16	20				
500 ml	1000980	9	12				
1000 ml	1000981	5	8				



Overhead stirrer "SE-100"



BRUSHLESS MOTORS FREE OF MAINTENANCE AND SPARKS.

PARTIALLY SEALED ENCAPSULATED CASES AGAINST SUBSTANCES AND AGRESSIVE VAPOURS. THERMAL PROTECTION FOR OVERHEATING WITH AUTO POWER OFF AND LED INDICATOR.

HEIGHT ADJUSTABLE STIRRING ROD.



Excellent blending results in an accurate and repeatable control.

Stirring by electronic control of speed, from 60 to 1500 rpm.

Specially designed for those samples in which its viscosity changes during the process.

Max stirring volume: **Up to 40 L** (in H20)
Max stirring viscosity: **Up to 100.000 mPas.**Speed digital reading and pre-selection (rpm).
Adaptable spindle to stirring rods up to 9.5 mmØ.

MODEL

Part No.	Height/Width/Depth	Power	Weight
	(external) cm	W	Kg
7000385	17.5 7.5 13.5	100	2.8







Shaft with AISI 304 stainless steel rotor. Shaft: 7mm Ø.

Part No.	Rotor type	Ø shaft mm	Total length mm	Features
6002891	1	90	400	Double fly vane
6002892	2	60	400	Fly vane
6002893	3	50	400	Fixed blade
6002894	4	50	400	Double cross vane
6002895	5	58	400	Helix
6002896	6	70	450	Six hole paddle
6002897	7	70	450	Three hole paddle
5050921	8	46	400	Radial turbine
5050932	9	45	400	Turbo helix 45 mm Ø
5050933	10	66	400	Turbo helix 66 mm Ø



Bosshead "MAXIM"

For large diameters. Chrome plated with one opening to hold rods up to 30 mm \emptyset . Second opening up to 20 mm \emptyset . Part No. 7000408





SHAFTS WITH PTFE PADDLES

Shaft of 8 mm Ø.

Onait of O III	III D.					
Part No.	Rotor type	Ø shaft mm	Total Length mm	Features	Speed range r.p.m.	
5050934	11	50	500	Helix PTFE	500	
5050935	12	80	500	Anchor PTFE	500	



ADJUSTABLE HEIGHT SHAFT STIRRERS

Shafts made of AISI 304 stainless steel coated in PTFE and PTFE stirrer paddles. Shaft: 8 mm Ø. Adjustable height paddles using a locking thread. Can be used with differing paddles on the same shaft to create a

can be used with differing paddles on the same shart to cre multiple mixing action.

Part No.	Rotor type	Ø shaft mm	Total Length mm	Features			
5050936	13	70	400	Shaft PTFE			
5050937	14	40	400	Helix PTFE			
5050938	5050938 Shaft PTFE Spare or additional paddle for multiple stir effect						
5050939	Helix PT	FF Snare or a	additional naddle	for multiple stir effect			



Rod of AISI 304 stainless steel of 20 Ø x 800 mm high. Designed to support weights up to 20 Kg at 300 mm height from the base. Useful for stirrers.

Weight: 7 Kg. Part No. **6000410** NEW



High volume orbital shaker incubator "1102" and 2102"



MODEL 1102 ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C TO 60 °C.
MODEL 2102 REFRIGERATED FOR ADJUSTABLE TEMPERATURE FROM +5 °C TO 60 °C.
FAN AIR CIRCULATED, WITH OR WITHOUT REFRIGERATION.

DIGITAL ELECTRONIC CONTROL OF SPEED, TEMPERATURE AND TIME. FOR LONG OPERATING PERIODS.

SAFETY:

TRIPLE PROTECTION MOTOR DRIVE:

OVER TEMPERATURE, IRREGULAR MOVEMENTS, AUTOMATIC STOP WHEN THE DOOR IS OPENED, TEMPERATURE ALARM, THERMAL CUTOUT.

FFATURES

Metallic external case epoxy-coated.

Made with two doors framework type and double glass windows which allow visibility of the material exposed in the inside

Two removable upper trays for high volume Erlenmeyer flasks accommodation in the lower tray.

Adjustable outlet hole for the air in the inside of the chamber.

With wheels for movements and adjustable support for a stabled fixing.

Current and set digital reading of temperature, time and speed parameters.

Acoustic alarm and illuminated indicator at the end of the cycle and also to indicate temperature and / or time programs errors.

Automatic switch off mechanism when the shaker door is opened.

Motor over load protection, with automatic switch off activation.

Induction drive motor, where the drive mechanism has an anti vibration system. Refrigerated models have a hermetically sealed compressor unit with a ventilated condenser.

Internal case made of stainless steel AISI 316 with polished finish and glazed door for a handy visibility of the samples process.

Internal platform with multi-adapter positions for differing accessories, such as: conical flask adapters, microtitre plate holders and universal tray with adjustable tension springs for other non standard sized vessels.

MODEL 1102 2102 Part No. 5312094 5312095 Rotation amplitude 26 mm 26 mm Controllable speed range from 30 to 300 r.p.m. from 30 to 300 r.p.m. Controllable speed in steps of 1 r.p.m. 1 r.p.m. Refrigeration No Yes Heating Yes Yes Controllable temperature range from ambient +5 °C from 5 °C to 60 °C to 60 °C to 60 °C Controllable temperature in steps of 0.1 °C 0.1 °C	
Rotation amplitude 26 mm Controllable speed range from 30 to 300 r.p.m. from 30 to 300 r.p.m. Controllable speed in steps of 1 r.p.m. 1 r.p.m. Refrigeration No Yes Heating Yes From ambient +5 °C to 60 °C Controllable temperature range Controllable temperature in steps of 0.1 °C Refrigeration 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 26 mm 6 rom 30 to 300 r.p.n. 7 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. 2 res	
Controllable speed range from 30 to 300 r.p.m. from 30 to 300 r.p.m. 1 r.p.m. 1 r.p.m. Yes Heating Yes from ambient +5 °C to 60 °C Controllable temperature range Controllable temperature in steps of 1 r.p.m. 1 r.p.m. 1 r.p.m. Yes Yes 1 r.p.m. Yes Yes 1 r.p.m. Yes Yes 1 r.p.m. 1 r.p.m. 1 r.p.m. 1 r.p.m. Yes Yes Yes 1 r.p.m. Yes Yes 1 r.p.m. Yes Yes 1 r.p.m. Yes	
Controllable speed in steps of 1 r.p.m. 1 r.p.m. Refrigeration No Yes Heating Yes Yes Controllable temperature range from ambient +5 °C to 60 °C Controllable temperature in steps of 0.1 °C 0.1 °C	
Refrigeration No Yes Heating Yes Yes Controllable temperature range from ambient +5 °C to 60 °C to 60 °C Controllable temperature in steps of 0.1 °C Ves from 5 °C to 60 °C 0.1 °C	n.
Heating Yes Yes Controllable temperature range from ambient +5 °C from 5 °C to 60 °C to 60 °C Controllable temperature in steps of 0.1 °C 0.1 °C	
Controllable temperature range from ambient +5 °C to 60 °C Controllable temperature in steps of from ambient +5 °C to 60 °C to 60 °C O.1 °C	
Controllable temperature range to 60 °C to 60 °C Controllable temperature in steps of 0.1 °C 0.1 °C	
to 60 °C to 60 °C Controllable temperature in steps of 0.1 °C 0.1 °C	
on the state of th	
Chamber temperature uniformity $\pm 1-2$ °C $\pm 1-2$ °C $\pm 1-2$ °C	
Timer from 0 to 500 hours from 0 to 500 hours	S
Usable platform dimensions 734 x 458 mm 734 x 458 mm	
Platform number 2 2	
Maximum Erlenmeyers capacity *4 x 5000 ml, or *8 x 3000 ml, or *8 x 2000n	nl
(see accessories) or 24 x 1000 ml, or 44 x 500 ml, or 56 x 250	
or 104 x 100 ml or 104 x 50 ml	,
Dimensions Height x Width x Depth 1440 x 950 x 700 mm 1440 x 950 x 700 mm	m
Weight 265 Kg 272 Kg	
Power 930 W 1200 W	


*Note: An upper tray divided into two removable pieces which allow high volume Erlenmeyer flasks accommodation, from 2000 to 5000 ml. in the lower tray or the Universal platform which is supplied as an accessory.

230v/50-60Hz





CONTROL PANEL



- O. LCD Display.
- 1. Heater "ON" indicator.
- 2. Modify and confirm push button.
- 3. Time function push button.
- 4. Temperature push button.
- 5. Increase value push button.
- **6.** Reduce value push button.
- 7. Rotation push button.
- 8. Speed push button.
- 9. Start / Stop Push button.
- 10. Alarm indicator.

Refrigerated model also carry:

Refrigeration compressor "ON" indicator.
Adjust temperature push button.

ACCESSORIES Erlenmeyer and flask adapters.

Made from hardened sprung stainless steel.



Part No. 5312105 Adapter for 50 ml
Part No. 5312106 Adapter for 100 ml
Part No. 5312107 Adapter for 250 ml
Part No. 5312108 Adapter for 500 ml
Part No. 5312109 Adapter for 1000 ml
Part No. 5312110 Adapter for 2000 ml
Part No. 5312111 Adapter for 3000 ml
Part No. 5312112 Adapter for 5000 ml



Adapter for microtiter holder From 85 x 130 mm., allows up to 6 plates. Part No. 5312113



Universal platform tray with elastic tension clips that hold in position any type of vessel, flasks, beakers, racks,

Dimensions 734 x 458 x 75 mm Part No. **5312120**

Power supply



Orbital shaker incubators "D-1102" y "D-2102"



MODEL D-1102 FOR ADJUSTABLE TEMPERATURES FROM AMBIENT +4°C UP TO 60°C. MODEL D-2102 REFRIGERATED FOR ADJUSTABLE TEMPERATURES FROM 4°C UP TO 60°C. FAN AIR CIRCULATED, WITH OR WITHOUT REFRIGERATION. DIGITAL ELECTRONIC CONTROL OF SPEED, TEMPERATURE AND TIME. FOR LONG OPERATING PERIODS.

APPLICATIONS

For use in cell and bacterial cultures, fermentation technology, etc...

It is widely used in research institutions, universities, epidemiological prevention centres, doping tests, agriculture and molecular biology.

FEATURES

Metallic external case epoxy-coated.

Current and set digital reading of temperature, time and speed parameters. Acoustic alarm and illuminated indicator of end of cycle and also of malfunction in programmed temperatures &/or time. Automatic switch off protection mechanism in case the door is opened.


Motor overload protection with automatic switch off and activation.

Compact unit with brushless motor assembly built-in over an anti-vibration mount and hermetic compressor group with aired condenser in models with refrigeration.

Case with methacrylate lid for an easy viewing of sample processes. Internal platform with multiadapter positions to different accessories, such as: adapters for Erlenmeyer, microtiters plates or universal tray with adjustable tension springs for vessels of different sizes.

MODEL	D-1102	D-2102			
Part No.	5312096	5312097			
Rotation amplitude	25 mm	25 mm			
Controllable speed range	from 25 to 300 r.p.m.	from 25 to 300 r.p.m.			
Controllable speed in steps of	1 r.p.m.	1 r.p.m.			
Refrigeration	No	Yes			
Heating	Yes	Yes			
Cantuallable temperature renge	From ambient +4 °C	From ambient 4°C			
Controllable temperature range	up to 60 °C	up to 60 °C			
Controllable temperature in steps of	0.1 °C	0.1 °C			
Chamber temperature uniformity	±1 °C	±1 °C			
Timer	From 0 to 999 hours	From 0 to 999 hours			
Usable platforms dimensions	460 x 460 mm	460 x 460 mm			
Platforms number	1	1			
Erlenmeyers maximum capacity	9 x 1000 ml or 9 x 50	00 ml or 16 x 250 ml,			
(see accessories)		l or 25 x 50 ml			
Height x Width x Depth	620 x 630 x 840 mm	620 x 630 x 840 mm			
Weight	76 Kg	98 Kg			
Power	800 W	800 W			
Power supply	220v/50-60Hz				





- 1. Left push button.
- 2 Right push button.
- 3. Menu selection push button.
- 4. Enter push button.
- 5. Alphanumeric screen.
- 6. Increase value push button.
- 7. Reduce value push button.
- 8. Start/Stop push button.

ACCESSORIES Erlenmeyer flasks adapters.

Built-in stainless steel with upper spring.



Part No. 5312105 Adapter for 50 ml Part No. **5312106** Adapter for 100 ml Part No. **5312107** Adapter for 250 ml Part No. **5312108** Adapter for 500 ml Part No. **5312109** Adapter for 1000 ml



Adapter for microtiter holder From 85 x 130 mm, allows up to 6 plates. Part No. 5312113



COMECTA Ultrasonic homogenizers "CY-500"

INTENSE CAVITY SHAKING. FOR 10 TO 500ML VOLUMES.



OPERATION

The ultrasonic generator transforms the electric current from 220V 50Hz into ultrasonic energy of 20Hz, to feed the ultrasonic transducers. The transducer produces an elastic distortion by following the alternate voltage, which traduces in a longitudinal mechanical vibration that produces a cavitation effect in the titanium alloy probe submerged in the solution. This generates a multitude of micro bubbles which releases a considerable energy in the probe end that allow a sample intense shaking and vibration.

APPLICATIONS

Biology, chemistry, environment, general analysis laboratory.

It allows disintegration of cells and biological tissue, DNA protein extraction, ether hydrolysis RNA, reactions output acceleration and increase, earth and sediments treatment according to EPA SW methods which simplifies the Soxhet method in time, quality control, R&D, sample gas remove, dissolution and homogenization, emulsion, dispersion, liposomes formation, proteins microencapsulation, etc.



FEATURES

The equipment consists of a control and power unit and an ultrasonic transducer (homogenizer).

Control unit: Epoxy coated external metallic case and ABS front panel. It consists of a LCD display, cycle time programming push buttons, impulse intervals, maximum power and temperature with 10 programs configurable by the user. Probe and connector for the sample temperature control.

Ultrasonic transducer consists of a piezoelectric converter with a ¼ inches titanium alloy probe (5,6mmØ and 60mm height).

It is hold by means of a retort stand and an adjustable clamp, max. 44 cm height.

ACCESSORY Anti-noise cabin

Made in double isolated layer that protects from high noises emited by the ultrasound effects. It comes with a transparent door.

Part No. 5059601



CONTROL PANEL

On/off switch.

Parameters indicator LCD display.

Ultrasonic power adjustable control (1-99%).

Push button from 0 to 9 for parameters selection: Time (1" - 2hours).

Maximum impulses and temperature.

SET push button: select parameters.

CLEAR push button: cancel parameters.

ENTER push button: validate parameters.

START/STOP push button: cycle beginning and end.

PAUSE push button: pause function.

RECALL-SAVE push button: memory to save 10 programs.

MODEL

Part					Power	Weight
No.	Frequency	(e:	(exterior) cm		W	Kg
5059600	20	22	19	33	500	6,7

It is supplied with a support-rod adjustable in height.



V solid mixer homogeneizer "VS-5", "VS-8" and "VS-50"

CAPACITY: FROM 5 TO 50 KG. (INDUSTRIAL DESIGN)



APPLICATIONS

For the use in pharmaceutical, chemistry and alimentary industry.

For powder, granulated and small concentrated batches that are sent to the production plant after being investigated in laboratory.

FEATURES

Made of stainless steel AISI 304.

It have two loading entrances in V and an unloading hole which is protected with isolated covers with silicon gaskets.

Engine group with shaking device at a constant speed of 18 rpm.

CONTROL PANEL

Start button. Stop button. Container positioner button. Timer from 1 to 99,9h.

MODELS

Model	Part no.	Volume Liters	Load capacity (50%) Liters	RPM	ŭ,	/ Width . xternal c		Power W	Weight Kg
VS-5	5810100	5	2,5	18	58	72	37	550	47
VS-8	5810101	8	4	18	64	81	37	550	52
VS-50	5810102	50	20	15	130	150	50	750	250



SPARE PART

Silicone gasket for VS-5, VS-8 and VS-50 Part No. 5810110 (Unit).



Multifunctional grinding mill "TR-20"

FOR CRASHING AND SHEARING (INDUSTRIAL DESIGN)

APLICATIONS

For a wide range of products grinding, solid, fragile, soft and semi-hard ones. For the use in pharmaceutical, chemistry and alimentary industry in general.

FEATURES

Made of stainless steel AISI 304.

Device for different kind of easily changeable sieves. From 0,12 mm Ø to 0,80 mm Ø. It is supplied with a sieve of $0,18 \text{ mm } \emptyset$.

Only suitable for grinding samples with relative humidity below 20%.

Model	Part No.	Capacity Kg / h	Db	RPM	Measures (cm) Height./ Width. / Depth.		Power Kw	Weight Kg	
TR-20	5810000	10 a 30	85	4200	70	40	80	2,2	85



SPARE PARTS Shear Part no. 5810006





Device for an easy sieve insertion.

VS-50 Model







EVAPORATORS



PERISTALTIC DOSING PUMP



MEMBRANE VACUUM AND PRESSURE PUMP "VACUM-PRES"



MEMBRANE VACUUM PUMPS



ROTARY VEIN VACUUM PUMPS



WATER RECIRCULATION VACUUM PUMP



LIQUID TRANSFER PUMP, ELECTRIC



LIQUID TRANSFER PUMP, MANUAL



CENTRIFUGAL PUMPS

Evaporators
Peristaltic dosing pumps
Rotary vein pump
Membrane vacuum pumps
Vaccum pumps
Water recirculation pumps
Centrifugal pumps
Liquid transfer pumps

page	46 and 47
page	49
page	50
pages	50 to 52
page	52 and 53
page	53
page	53
page	54



Rotatory evaporator "RS 3000-V" FOR 1000, 2000 AND 3000 ML. FLASKS.

Quality and safety. Rotary speed and bath temperature regulation

FEATURES

Compact rotatory evaporator composed of stainless steel AISI 304 case. Manual elevator system and bath made of dur-al alloy chemically treated. Rotatory speed: from 20 to 120 rpm.

Bath temperature: up to 200 °C.

External probe temperature: up to 100 °C.

Vacuum breaker with air inlet filter.

CONTROL PANEL



- 1. Starting main switch ON/OFF.
- 2. Digital reading indicator of external probe temperature.
- 3. Bath temperature digital regulator.
- 4. Rotator pump switch.
- 6. Water recirculation pump switch.
- 7. Regulator knob of rotating speed. Vacuometer built-in the case.

5. Heater bath switch. RS 3000-V

It is completely supplied with a tap for solution, a PT100 probe, a coolant, one 3000 ml. flask, one 2000 ml flask, one 1000 ml. flask and its accessories.


6003000 0 to 120

Regulation t.p.m.

Part no



GLASS AND ITS ACCESSORIES DIAGRAM SPARE PARTS



ACCESSORY

Height / Width / Depth

(exterior) cm

46

26

Weight

Kg

28

W

1010

Vacuum regulator for ambient air entrainment. Part No. 6003010





High capacity rotatory macro evaporator "COM-1010" and "COM-1020"

CAPACIDAD COM-1010: 10 LITERS CAPACIDAD COM-1020: 20 LITERS

Supplemented with a wide range of optional accessories.

APPLICATIONS

The rotatory macro evaporator is mainly used in the level pilot and in biological, pharmaceutical, chemical, food and other product productions.

FEATURES

With a large diameter evaporator flask that allows the diffusion of the solution and a quick evaporation.

The macro evaporator allows the addition of multiple accessories to form a complete distillation equipment.





MODELS	COM-1010	COM-1020				
Part No.	5195410	5195420				
Evaporator flask	10 liters	20 liters				
Receptor flask	5 liters 10 liters					
Speed regulation	Continuous speed regulation					
Rotary speed	10-130 r/min					
Heating bath	Stainless steel					
Temperature control	Digital dis	splay, 0-99° C .				
Achievable vacuum	399	9.9 P a				
Evaporator capacity	Water >3.2L/h Alcohol>6.5L/h	Water >3.2L/h Alcohol>6.5L/h				
Lift	Powe	red lift				
Power supply	220 V/50 Hz	380 V/50 Hz III				
Power	4.8 Kw	6.3 Kw				
Dimensions (cm)	170x55x92 (HeightxWidthxDepth)	210x60x120 (HeightxWidthxDepth)				

ACCESSORIES 1. Membrane pump COM-401. Part No. 5195421 2. Coolant circulation COM-30L. Part No. 5195422 Vacuum pump. (See page. 50) Part No. 5900621 Vacuum trap. (See page. 50) Part No. 4001611



Peristaltic pump "Percom N-M"

ELECTRONIC SPEED CONTROL WITH DIGITAL DISPLAY,
WITH MULTI HEAD ADAPTER FOR TUBES FROM 4 TO 12 mm Ø (ext.). MAXIMUM WALL THICKNESS 1,5 mm.
ADJUSTABLE FLOW RATES FROM 1080 TO 73800 ml/hr.

FEATURES

Small footprint bench top model Quiet with minimal vibrations. Speed control 20 to 200 r.p.m.

CONTROL PANEL

Push button Stop/Start/switch.
Digital speed display.
Push button increase speed control.
Push button decrease speed control.


BACK

External control connection using a foot switch for Start / Stop, supplied as an accessory.

MODEL

Part No.	Height / Width / Depth	Pump rate of H₂O	Power	Weight
	(exterior) cm	ml/hr	W	Kg
3001300	16 12 18	1080 to 73800	12	2.6

Graph of pump rate/r.p.m. dependent on the internal tube Ø.



Graph showing with liquids: Distilled water. These characteristics vary according to the density of the liquid.





Peristaltic dosing pump "Peristaltic PR-2003"

ELECTRONIC SPEED CONTROL AND DIGITAL DISPLAY.
MULTI TUBE ADAPTER FROM 4 TO 12 mm Ø (ext.). MAXIMUM WALL THICKNESS 1,6mm.
ADJUSTABLE DOSING RATE PROGRAMMABLE FROM 10 ML TO 10 LITRES.
MAXIMUM PUMP RATE 1500 ml / minute.

Main On/Off switch.

FEATURES

Small foot print, bench top. Silent and free from vibrations.

Pump tube held by a simple clamp

Glass dosing tip "PECIPUMP" provides precise dispen-

Dispense rate: The pump can dispense set volumes over a set time period.

Flow rate: The pump dispenses at a constant rate over a set time or continuous operation.

MODEL


BACK


Part No. Height/Width/Depth (exterior) cm		Pump rate of H ₂ O	Power	Weight			
		ml/hour	W	Kg			
412905	0	25	19	23	90,000	65	4

Foot switch connection Start - Stop, for operation from

a distance, supplied as an accessory.

Graph of pump rate/r.p.m. dependent on the internal tube \emptyset .





Graph showing with liquids: Distilled water. These characteristics vary according to the density of the liquid.

ACCESSORIES FOR PERISTALTIC PUMPS "PERCOM N-M" AND "PR2003"

Standard silicon tubing. (see page: 52)

START STOP foot pedal controller. Part No. 1001044

CONTROL PANEL

Numeric keypad to enter values. LCD display of 16 x 2 characters.

Calibration memory.

Change pump direction, control from the panel.

Rotation direction indicator

Foot switch connection for Start/Stop.

TECHNICAL SPECIFICATIONS

TECHNICAL OF ECHTICATIONS
Speed (r.p.m.):
Resolution (ml): 1 or 0.1
Dispensing range (ml): 10 - 9999
Pump rate See table
Nº of programs:
Nº of repeat doses:9999
Repetition interval:1 - 9
Maximum tube wall (mm):
Maximum external tube diameter (mm):12
Note: The precision depends on the tube \emptyset and the speed.



COMECTIVE Membrane vacuum pumps "MV 50 and MV 100"

FOR AIR, GASES AND STEAM.



Free of lubrication • Free of maintenance • Quiet

APPLICATIONS

Vacuum filtration.
Vacuum distillation.
Rotating evaporators.
Gas extraction and transport.
Vacuum drying.

COMUN FEATURES

Constructed in aluminium.

Membrane of polifluoropreno, resistant to the corrosion and long-life. Refrigeration system that allows a continuous operation of 24 hours.



MODEL "MV 100"

CONTROL PANEL

Vacuum and pressure gauges. Adjustable pressure and vacuum. Tensile handle for easy transport. Button On-Off.

CONTROL PANEL

Vacuum indicator. Vacuum regulator. Button On-Off.

MODELS	Part No.	Vaccum Rate m³/h	Vaccum Limit mbar	Pressure bar	Heigh / Width / Depth (exterior) cm	Power W	Power supply	Weight Kg
MV 50	5910100	1,8	200	2	23,5 16 21	160	220V/50Hz	8
MV 100	5910101	3,6	200	2,75	23,5 16 30	160	220V/50Hz	9

They are supplied with a 1 meter of silicone tube \emptyset 8 x 12mm.



COMECTA Vacuum pump

ROTARY VEIN PUMP WITH ANTI RETURN VALVE PREVENTS OIL FLOW BACK, SUITABLE FOR GENERAL LABORATORY APPLICATIONS. OVER TEMPERATURE MOTOR PROTECTION CUT-OUT AND MAIN ON/OFF SWITCH.

FEATURES

Asbestos veins and internal joints.

Aspiration connection 16 mm Ø in model 5900620 and 16 mm Ø in model 5900621.

Analogue vacuum gauge.

High oil volume and forced lubrication.

Oil reservoir charge capacity:

25 cc in model 5900620 and 55 cc in model 5900621.

Exhaust filter and ballast.

Shock absorber mounted.

Low noise level (62db)

Maximum working temperature 60 °C.

Portable with lifting handle.





MODELS

MIONET9						
Part No.	Pump rate m³/hr	Vacuum limit mbar	Height / Width / Depth (exterior) cm	r.p.m.	Power W	Weight Kg
5900620	3,6	0,1	23 27 12	1400	120	7,6
5900621	3.6	0.05	27 35 14	1400	180	11

ACCESSORY

High pressure **rubber tubes** of 18 x 26 mm Ø. Part No 5900624 High pressure **rubber tubes** of 8 x 15 mm Ø. For Vacuum oven connection "VACIOTEM". Part No 5900625

Pump oil for vacuum pumps, 1 litre bottle. Part No. 5515030

SPARE Vacuum gauge (Fits both models). Part No. 5515031





Membrane vacuum pump

SUITABLE FOR AIR, GAS AND VAPOURS.

FEATURES

Transfer, evacuation and compression, free from contamination. The models N-86 and N-022 of 2.5 and 4 bar of relative pressure, can be used also as a compressor.

Can be operated dry.

No lubrication required.

Maintenance free.

Silent.

The double head models are made with linear two phase membranes.

All Models (except model N-022), are made with a patented system of diaphragms manufactured with optimum tension, that results in a smaller more durable pump.

Approximate loss 6 x 10-3 mbar x L/s

Complete with thermal breaker and fuse.

Medium to high resistance to chemicals, (model dependant). We advise that you check the compatibility of the pump performance with the materials that will be in contact with the pump parts (see table).









MODELS / SPECIFICATIONS	N-86	N-022	N-816	N-820	N-840
Part No.	5910106	5910107	5910108	5910009	5910110
Vacuum pump rate (to 1 Atm, of pressure)	5.5 l/min	13 l/min	16 I/min	20 l/min	34 I/min
vacuum pump rate (to 1 Aun, or pressure)	0.33 m³/h	0.78 m³/h	0.96 m³/h	1.2 m³/h	2.04 m³/h
Vacuum limit (mbar) absolute	160 mbar	100 mbar	20 mbar	8 mbar	8 mbar
Pressure (bar) relative	2.5 bar	4 bar	0.5 bar	1 bar	1 bar
Tube connection, Ø interior	4 mm	6 mm	6 mm	10 mm	10 mm
Number of heads and phases	1	1	2	2	2
Head material	PPS	Aluminium	PPS	PTFE	PTFE
Membrane material (diaphragms)	PTFE	PTFE	PTFE	PTFE	PTFE
Valve material	FFPM	Stainless steel	FFPM	FFPM	FFPM
Motor protection level	IP-20	IP-20	IP-20	IP-44	IP-44
Gas operating temperature and permissible ambient			5 to 40 °C		
Current (A) / Power (W)	0.63 / 65	0.7 / 100	0.6 / 100	0.7 / 120	1.5 / 245
Dimensions height x width x depth (cm)	14 x 9 x 17	19 x 15 x 20	14 x 9 x 36	21 x 15 x 31	23 x 17 x 34
Weight (Kg)	1.9	4	3.95	9.3	12.6

APPLICATIONS AND CHEMICAL RESISTANCE	N-86	N-022	N-816	N-820	N-840
Chemical resistance (depends on medium being pumped.)	medium-high	medium-high	medium-high	medium-high	medium-high
For gasses and vapours very aggressive or corrosive	-	-	-	✓	✓
For gasses and vapours slightly aggressive or corrosive	√	✓	1	1	✓
For liquid transfer and evacuation	✓	1	✓	1	
For liquid transfer, evacuation and as a compressor	✓	1	-	-	-
Ideal for vacuum filtration	✓	✓	✓	1	✓
Ideal for vacuum filtration, evaporation, drying and gel drying	-	-	-	✓	✓

SPARE PARTS	N-86	N-022	N-816	N-820	N-840
Spares Kit: 1 membrane +2 valves + 2 gaskets	5910120	-	-	-	-
Spares Kit: 1 membrane +2 valves + 1 gasket	-	5910121	-	-	-
Spares Kit: 2 membrane +4 valves + 4 gaskets	-	-	5910122	5910123	5910124

COMECTA

Silicon tube for peristaltic and membrane

pumps

MODELS

MODET9	Part No.	interior Ø mm	exterior Ø mm
Standard tube*	8991202*	4	7
	8991201*	6	9
	8991200*	8	11



	Part No.	interior Ø mm	exterior Ø mm	
	5253300*	0.8	4	
	5253400*	1.7	4.9	
	0686203	3	7	
Calibrated	5253600*	3.1	6.3	
tubing*	5253500	4.8	9.8	
	0686210	5	9	
	5253900	6.3	11.3	
	0686216	8	12	

- (*) Tubing supplied in 1 meter lengths.
- (*) Part Numbers with asterisk* are adaptables to peristaltic pumps "Percom N-M" and "PR-2003". See page 49



Sink weights and glass connections

FOR CENTRIFUGAL AND PERISTALTIC PUMPS.

FEATURES

Made of glass.

Adapters for centrifugal and peristaltic pumps. For better aspiration and dispensing.



A. Sink weight "ASPIPUMP" Surface olive that fixes to the tube. The sink weight maintains the tube at the bottom of the vessel being used for the aspirant liquid.

5991204 4 mm Ø interior. 5991206 6 mm Ø interior. 5991208 8 mm Ø interior.

Part No. Tube size




B. Spout "PECIPUMP"

Surface olive is pushed onto the inside of the dispensing end of the tube, liquids can then be dispensed easily and accurately.

Part No. Tube size

5991203 4 mm Ø interior. 5991205 6 mm Ø interior. 5991207 8 mm Ø interior.



A. Sink weight "ASPIPUMP" an efficient way to aspirate to the last drop.

B. Spout "PECIPUMP" an exact way to dispense.



Membrane vacuum pump "Vacum-Sel"

FEATURES

Small footprint, bench top. Quiet and vibration free.

Maximum working temperature: 40 °C.

CONTROL PANEL

Mains power switch.
Vacuum gauge, scale 0 to -1 bar.
Vacuum in.
Air out.



Part No.	Pump Rate litres/min.	Vacuum limit	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3001001	1.08	-0.6 bar ≈ 440 mm Hg	19 14 19	35	3





Membrane vacuum and pressure pump "Vacum-Pres"

FEATURES

Small footprint, bench top. Quiet and vibration free.

Maximum working temperature: 40 °C.

CONTROL PANEL

Mains power switch.
Vacuum gauge, scale 0 to -1 bar.
Pressure gauge, scale 0 to 2.5 bar.
Control for pressure/vacuum.
Vacuum in.

Pressure out.

MODEL

Part No.		ACUUM	Pres				/ Depth	Power	Weight
	Pump Rate /mir	n. Vacuum limit	Pressure Rate /min.	Max. Pressure	(ex	terior)	cm	W	Kg
3001002	1.08	-0.6 bar=440 mm Hg	1.3	1.2 bar	19	14	19	35	3.2





Water recirculation vacuum pump

DESCRIPTION

Made from a centrifugal pump that pulses water at a high velocity through suction aspirators, creating a pump rate of 10 litres/min from each pump connection. The vacuum achieved, 5 to 40mbar is dependent on the water temperature in the reservoir. Aspirated gasses from both inputs mix in the 10 litre reservoir. Although the interior of the pump is made from stainless steel and polypropylene it is a requirement to change the water in the pump to limit any attack of any chemicals or vapours absorbed in the water. For better protection of equipment it is recommended to connect a "Scrubber" to wash and neutralise any corrosive gasses. The pump includes an exit valve that allows the water to be circulated externally. This can be used to empty and change the water in the reservoir to form an external re-circulation of water.

FEATURES

Reduced water consumption.

2 independent connections, aspiration pump vacuum rate 10 litres/min.

Vacuum limit: 5 and 40 mbar (dependent on water temperature).

Capacity: 10 litres.

Low noise level (<65dBA).

Valve connections for entrance and exit of vacuum and external water recirculation.

CONTROL PANEL

2 vacuum connections.

2 analogue vacuum gauges.

Illuminated mains switch.

ACCESSORY

Extraction system and neutralisation of vapours.

Made up of a "Scrubber" unit that blocks the passage of condensed acid vapours and neutralises the gas vapour in NaOH, and a recirculating water pump that provides an adequate pump rate to aspirate vapours.

It is supplied with 3 Kg. acid vapours neutralizer solution.

"Scrubber" unit

Part No.	. 3	Height / Width / Depth (exterior) cm		
4001611	32 3	31 16	2	

Spare part:

3 Kg. acid vapours neutralizer solution. Part No. 4001610



MODEL

Part No.	Height / Width / Depth		th Vacuum limit	Pump rate	Power	Weight	
	(ex	terior)	cm	mbar	litre/minute	W	Kg
4001612	44	39	28	betwen 5 and 40	10	280	10



Complete unit image for a Kjeldahl distillation block.

COMECTA Centrifugal pumps "MP"



FEATURES

Centrifugal pump for re-circulating and transferring of liquids.

A selection of models with pump rates of 8 to 60 litres/min.

The motor power to the turbine is transmitted via a magnetic coupling.

All parts in contact with liquids are chemically resistant materials.


Manual priming.

SERIES MP

Dimensions per model in mm.

MODEL	Depth	Height	Width	Ø
MP-15R	129	82	74	60
MP-20R	210	108	116	92
MP-40RX	277	130	122	102
MP-100R	322	175	150	110

SEISUN



TECHNICAL DATA	MP-15R	MP-20R	MP-40Rx	MP-100R		
Part No.	5900610	5900611	5900612	5900613		
Pump rate (max) P (pressure) =0 (litres/min)	8	17	47	60		
Precision (min) Q (pump rate)=0 (bar)	0.21	0.37	0.40	0.50		
Tube connection (aspiration)	14 mm Ø	18 mm Ø	26 mm Ø	26 mm Ø		
Tube connection (output)	14 mm Ø	18 mm Ø	26 mm Ø	26 mm Ø		
Normal working temperature (°C)	0 - 80	0 - 80	0 - 80	0 - 80		
Function cycle	continuous	continuous	continuous	continuous		
Power consumption (W)	10	20	65	250		
Pump material: Body	Po	lypropylene with	a mixture of fibre	glass		
Turbine material:	Po	lypropylene with	a mixture of fibre	glass		
Pump spindle/axle material:	ceramic aluminium oxide					
Seals and junctions material:	EPDM	EPDM	EPDM	EPDM		
Voltage:	220-240V 50-60Hz					



COMECTA Liquid transfer pump, manual "RP-90P"

MADE FROM CHEMICALLY AND CORROSION RESISTANT POLYPROPYLENE. ROTARY MOVEMENT.

FEATURES

Manual rotary suction movement that extracts liquid or sucks up liquid from a drum, composes of a pump handle, dip tube 1.2 M long by 1 inch \emptyset , with adjustable clamp that fastens to the barrel and an outlet tube of 3/4". Easy to use, light weight and very robust.

MODEL

Part No.	Pump rate litres/min.	Weight Kg	
5252040	20	1.6	





COMECTA Liquid transfer pump, electric "FA-100" SUITABLE FOR NON CORROSIVE AND NON FLAMMABLE LIQUIDS.

FEATURES

All metal with electric air displacement motor drive, with steel suction tube 1 m long and 38 mm Ø. Maximum height 13 m.

MODEL

Part No.	Pump rate litres/min.	R.p.m.	Power W	Weight Kg
5929041	86	10000	120	6.6





SPARES

1. Drum faucet: 2" in diameter.


Made of brass. Part No. 5252038

2. Drum Plug key for opening drum, standard model.

Made of brass. Part No. 5252039

3. Air bleed drum cap 2" in diameter.

Made of brass. Part No. 5252037





Vertical desktop autoclaves	Pages	58 to 59
Pressure regulator		
Vertical stand autoclaves.	Pages	60 to 63
Microprocessor control and atmospheric purge		
Water distiller for autoclaves	Pages	63
Horizontal stand autoclave	Page	64
Microprocessor control		
Horizontal stand autoclaves	Pages	65 to 67
Microprocessor control and fractional vacuum purge		
Polypropylene and polyester bag sealers	Pages	68
Autoclaves for preserved food	Pages	69 to 70



ELECTRICAL AUTOCLAVES FOR STERILIZATION

Wide range of electrically heated bench top and top loading autoclaves. Pressure range from 0.62 to 2 bar, Capacity: 1.5, 8, 12, 18, 20, 23, 25, 50, 80 and 150 litres.

Complies with the pressure and autoclave directives



PRINCIPAL OF STEAM STERILIZATION:

STERILIZATION IS THE PROCESS OF THE DESTRUCTION AND ELIMINATION OF ALL FORMS OF MICROBIAL LIFE, INCLUDING BACTERIA AND ITS SUB-FORMS.

THE LABORATORY AUTOCLAVE CONSISTS OF A VESSEL IN WHICH STEAM FLOWS AT A HIGH PRESSURE WITH A STERILIZATION PURPOSE. SATURATED STEAM WILL IMMEDIATELY CONDENSE AND THUS RELEASING ITS HIGH EVAPORATION HEAT, WHICH IS IN FACT THE DESTRUCTIVE AGENT. IN ORDER THAT STERILIZATION WOULD BE EFFECTIVE, STEAM MUST BE SATURATED AND, CONSEQUENTLY, ALL THE AIR INSIDE THE VESSEL MUST BE PURGED BEFORE INITIALIZING THE PROCESS.

IN THE FOLLOWING PAGES J.P. SELECTA OFFERS AN EXTENSIVE RANGE OF AUTOCLAVES WITH A CHOICE OF ATMOSPHERIC PURGING OR VACUUM PURGING SYSTEMS. A CHOICE OF MICROPROCESSOR CONTROLLED SYSTEMS ARE OFFERED INCLUDING DRYING CYCLES AT THE END OF THE STERILIZATION CYCLE.



Electrical Autoclaves for sterilization









MICROPROCESSOR CONTROL OF TEMPERATURE AND STERILIZATION CYCLE,

ATMOSPHERIC AIR PURGE: AUTESTER P 1,5 (mod. vertical)

PRESOCLAVE III 50 (mod. vertical)

PRESOCLAVE III 80 (mod. vertical)



MICROPROCESSOR CONTROL OF TEMPERATURE AND STERILIZATION CYCLE, **PULSED VACUUM AIR PURGE:**

CONTROL BY PRESSURE WITH ATMOSPHERIC PURGE: MICRO 8 MED 12 MED 20

AUTESTER ST DRY PV-B 12 (mod. horizontal) AUTESTER ST DRY PV-B 18 (mod. horizontal) AUTESTER ST DRY PV-B 23 (mod. horizontal)











RS-232 and USB interface output of parameters to a computer, printer or USB adapter.

The printer can be built-in the autoclaves together with the control panel:

AUTESTER ST DRY PV-III 50 (mod. vertical) **AUTESTER ST DRY PV-III 80 (mod. vertical)** AUTESTER ST DRY PV-III 150 (mod. vertical) PRESOCLAVE III 50 (mod. vertical) PRESOCLAVE III 80 (mod. vertical)



SUMMARY TABLE OF MODEL RANGE

OUICK OVERVIEW

VOID OF HOUSE IN MICE IN MICE IN MICE									SOLOW OFFICE
MODELS	Part No.	System DRY	Capacity litres	Temperature °C	Pressure bar	Power W	Control	Purge control	Туре
MICRO 8	4001745	NO	8	121 and 134	1 and 2	1000	Pressure thermostat	Manual	Vertical bench top
MED 12	4001756	NO	12	121 and 134	1 and 2	1800	Pressure thermostat	Atmospheric	Vertical bench top
MED 20	4001757	NO	20	121 and 134	1 to 2	2500	Pressure thermostat	Atmospheric	Vertical
PRESOCLAVE III 50	4001758	NO	50	115 to 134	0.62 to 2	3600	Electronic controller	Atmospheric	Vertical
PRESOCLAVE III 80	4002136	NO	80	115 to 134	0.62 to 2	3600	Electronic controller	Atmospheric	Vertical
AUTESTER ST DRY PV-B 12	4001727	YES	12	121 and 134	1 and 2	1800	Microprocessor	Purged vacuum	Horizontal bench top
AUTESTER ST DRY PV-B 18	4001734	YES	18	121 and 134	1 and 2	1500	Microprocessor	Purged vacuum	Horizontal bench top
AUTESTER ST DRY PV-B 23	4001735	YES	23	121 y 134	1 and 2	1700	Microprocessor	Purged vacuum	Horizontal bench top
AUTESTER ST DRY PV-III 25	4002420	YES	25	105 to 134	0.21 to 2	2510	Microprocessor	Purged vacuum	Horizontal bench top
AUTESTER ST DRY PV-III 50	4002516	YES	50	105 to 134	0.21 to 2	3600	Microprocessor	Purged vacuum	Vertical
AUTESTER ST DRY PV-III 80	4002517	YES	80	105 to 134	0.21 to 2	3600	Microprocessor	Purged vacuum	Vertical
AUTESTER ST DRY PV-III 15	04002518	YES	150	105 to 134	0.21 to 2	7500	Microprocessor	Purged vacuum	Vertical
STERIL-FOOD COM 20	4002421	NO	20	70 a 121	1	2500	Microprocesador	Automatic	Vertical
STERIL-FOOD MAN 80	4002427	NO	80	70 to 121	1	3600	Microprocessor	Automatic	Vertical
STERIL-FOOD COM 80	4002422	NO	80	70 to 121	1	3600	Microprocessor	Automatic	Vertical
STERIL-FOOD MAN 150	4002428	NO	150	70 to 121	1	7500	Microprocessor	Automatic	Vertical
STERIL-FOOD COM 150	4002423	NO	150	70 to 121	1	7500	Microprocessor	Automatic	Vertical



Autoclaves for sterilization "MICRO 8", "MED 12" and "MED 20"

PRESSURE REGULATOR THERMOSTAT OF 1 AND 2 bar. TEMPERATURES 121 °C AND 134 °C. VERTICAL DESKTOP MODELS. "MICRO 8" CAPACITY: 8 LITRES. "MED 12" CAPACITY: 12 LITRES. VERTICAL STAND MODEL "MED 20" CAPACITY: 20 LITRES.

An economic alternative. Semiautomatics.

SAFETY:

SAFETY VALVE: PREVENTS THE PRESSURE WITHIN THE AUTOCLAVE OVERPASSING THE MAXIMUM SET LIMIT. SAFETY PRESSURE: DISCONNECTS POWER TO THE HEATER IN CASE OF OVERPASSING THE PRESSURE SET LIMIT. SAFETY THERMOSTAT: DISCONNECTS THE HEATING IN CASE OF OVERTEMPERATURE OR LACK OF WATER IN THE TANK.

THERMALLY INSULATED LID MADE OF MOLDED THERMOPLASTIC.

COMMON FEATURES

AISI 304 stainless steel external case and top.

Stainless steel chamber and lid, with easy opening.

INCOLOY special stainless steel heater element, which provides long life, low corrosion and high temperatures resistance.

Manual air purge when starting ("MICRO 8" model).

Automatic atmospheric purge ("MED 12" and "MED 20" models).

Safety device to prevent lid opening.

Starting main switch.




REAR PART Mains source.

Two fuse carriers. Safety thermostat.



SIDE PART Drain outlet.

Steam and drainage selective valve.



- 1. Pressure gauge.
- 2. Over temperature alarm light indicator.
- 3. Door open light indicator.
- 4. Sterilization phase light indicator.
- 5. Sterilization end of cycle light indicator.
- 6. Pressure selector 121°C or 134°C.
- 7. Light indicator at 121°C.
- 8. Light indicator at 134°C.
- 9. Timing from 0 to 60'.



Autoclave for sterilization "MICRO 8"

DESKTOP MODEL



MODEL

Part	Capacity	Ø Height	Height/Width/Depth	Power	Weight
No.	litres	(usable) cm	(exterior) cm	W	Kg
4001745	8	24 15	40 38 38	1000	13

It comes complete with a stainless steel wire basket of Ø 22x14cm high. Part no. 1004771.



Autoclave for sterilization "MED 12"

DESKTOP MODEL

Features: see page 58







MODEL

Part	Capacity	Ø Height	Height/Width/Depth	Power	Weight
No.	litres	(usable) cm	(exterior) cm	W	Kg
4001756	12	24 26	46 38 38	1800	18

It comes complete with a stainless steel wire basket of Ø 22x24cm high. Part no. 1004772.



Autoclave for sterilization "MED 20"

VERTICAL STAND MODEL.
CONFORMS TO APPARATUS PRESSURE REGULATION DIRECTIVES.







Part	Capacity	Ø Height	Height/Width/Depth	Power	Weight
No.	litres	(usable) cm	(exterior) cm	W	Kg
4001757	20	24 41	70 38 38	2500	41

It comes complete with two stainless steel wire baskets of Ø 22x24cm and Ø 22x14cm high. Part no. 1004772 and 1004771.





Autoclaves for sterilization "Presoclave III" 50 and 80



MICROPROCESSOR TEMPERATURE AND CYCLE CONTROL AND ATMOSPHERIC PURGE.
ELECTRONIC CONTROL OF TIME AND TEMPERATURE RANGE FROM 115 °C TO 134 °C. (0.62 TO 2 bar).
CONFORMS TO APPARATUS PRESSURE REGULATION DIRECTIVES.
VERTICAL MODEL: 50 AND 80 LITRES.

SAFETY

SAFETY OVER PRESSURE VALVE: PREVENTS OVER PRESSURE ABOVE THE MAXIMUM SET LIMIT.

SAFETY PRESSURE SWITCH: DISCONNECTS POWER TO THE HEATER IF THE UNIT EXCEEDS THE MAXIMUM PRESSURE LIMIT.

SAFETY THERMOSTAT AND WATER LEVEL: DISCONNECT POWER TO THE HEATER IF INSUFFICIENT WATER OR OVER TEMPERATURE.

SAFETY LID: PREVENTS THE LID BEING OPENED IF THE CHAMBER IS UNDER PRESSURE.

SAFETY INSULATION: THERMALLY INSULATED LID.

An economic alternative, semi-automatic

FEATURES

External case made from stainless steel.

INCOLOY stainless steel heater elements, hard waring, resistant to corrosion and high temperatures.

Sterilization chamber and door made from stainless steel.

Automatic purge to atmosphere

Lid.6 point closure.

CONTROL PANEL



- 1. Pressure gauge.
- 2. Illuminated alarm indicator.
- 3. Illuminated indicator end of cycle.
- 4. Illuminated indicator door open
- 5. Temperature display.
- 6. Time display
- 7. Temperature set pushbutton (range from 115 °C to 134 °C).
- 8. Pushbutton increase value.
- 9. Push button START
- 10. Push Button STOP
- 11. Push button Reduce value.
- 12. Push button select time, (range from 3' to 99').
- 13. Push button open lid

MODELS	Part No.	Capacity litres	Ø Height (usables cm)	Height / Width / Depth (exterior cm)	Power W	Weight Kg
PRESOCLAVE III 50	4001758	50	30 70	117 48 62	3000	75
PRESOCLAVE III 80	4002136	80	40 60	107 58 72	3600	100

Does not include baskets and drums, order separately.

ACCESSORIES: AISI 304 Stainless steel wire baskets and drums.

For Presoclave III 50. Part. no: 4001725.

Capacity, 3 baskets Ø 25 and 20 cm. height. Part No. $1000495\,$

For Presoclave III 80, Part. no: 4002136.

Capacity, 2 baskets Ø 36 and 28 cm. height. Part No. 1000496

Sterilization drum made from AISI 304 stainless steel.

For Presoclave III 50, Part. no: 4001725.

Capacity, 3 drum Ø 25 and 18 cm. height. Part No. 1002418

For Presoclave III 80, Part. no: 4002136.

Capacity, 2 drum Ø 35 and 27 cm. height. Part No. 1002419

Deteraent Bio-Sel.

Water distillation unit especially for autoclaves.

USB adapter model. Pen-Drive included (Memory board) for data storage.

Part No: **4120131**







BACK VIEW:

- 1. Mains power.
- 2. Safety valve.
- 3. Steam empty valve.
- 4. Safety thermostat
- 5. Drain selector valve.

COMPLEMENTS That needs to be factory installed:

Pt100 second sensor only temperature reader. Part No. 4001219 (see page 63). Printer of continuous paper roll with print out of temperature, pressure, time and mode. Part No. 4002421 (see page 62).

Loading and unloading crane. Part No. 1002426 (see page 70).



Autoclaves for sterilization "Autester ST DRY PV III"



MICROPROCESSOR PROCESS CONTROL.

AUTOMATIC VACUUM DRY AND PURGE "PREVACUUM"

PROGRAMMABLE TEMPERATURE FROM 105 °C TO 134 °C (0.21 TO 2 bar). CONFORMS TO APPARATUS PRESSURE REGULATION DIRECTIVES.

VERTICAL FLOOR STANDING MODEL: 50, 80 AND 150 LITRES.



CONFORMITY

SAFETY FOR EQUIPMENT FOR LABORATORY USE EN 61010-1

EN 61010-2-040 SAFETY FOR STERILIZER EQUIPMENT FOR LABORATORY USE

EN 61326 ELECTROMAGNETIC COMPATIBILITY FOR EQUIPMENT FOR LABORATORY USE.

FEATURES

Stainless steel exterior and top.

Stainless steel lid.

Lid, six closure point system.

Pulsed vacuum system and atmosphere purge.

Incorporated steam generator, connected to the main chamber.

Water reservoir for clean water for the steam generator boiler.

External water connection if required.

6 function modes.

RS 232 for PC parameters printing (back side).

Printer output.

USB output for PEN-DRIVE data logging and Computer connection.

Connection for a second temperature sensor.

Memory for 9 customer set programs.


SAFETY

Safety: The lid cannot be opened if the chamber is under pressure.

Safety valve: Prevents over pressure over the maximum limit.

Safety pressure switch: switches the heater power off to if over pressure.

Safety Lid: insulated lid protection.



CONTROL PANEL

TFT touch screen with virtual pushbuttons.



- 1. Graphical representation of real time process
- 2. Start/Stop virtual pushbutton.
- 3. Running program.
- 4. Press indicator.
- Temperature indicator.
- 6. Running program parameters.



BACK VIEW

- 1. Mains input.
- 2. External water supply.
- 3. Safety valve.
- 4. Vacuum vapour valve.
- 5. Air in filter valve.
- 6. Safety thermostat.
- 7. Safety thermostat. 8. Empty tank valve.
- 9. Drain tank valve.
- 10. RS232 connector.
- 11. 2nd Temperature sensor con nection.



MICROPROCESSOR FEATURES

The microprocessor controls the performing of the air extraction process, sterilizing and drying, executing the selected procedures and commands by the user, adjusting the following parameters:

- Sterilization temperature from 105°C up to 134°C.
- Sterilization time from 3' to 59'.
- Drying time from 20' to 60'.
- Up to 10 user-configurable programs.
- Six operating modes:
- Solid 1 Solid 2 fractions. Solid 3 fractions -Drying - Liquids - Liquids 2 probes
- Verification test
- Vacuum Bowie Dick.
- Settings:
- Units languages (Spanish, English, French).
- USB or printer module. Clock: date/time.
- External water: on/off.
- Password

- Calibration
- Temperature probe 1 Temperature probe 2
- Water dosing.
- Audible and visual alarm for overtemperature and insufficient temperature.
- Network failures detection and memorization
- Safety therrmostat by software
- Sterilization cycles data logging with capacity for over 500 cycles.
- Internal cleaning program
- Software for PC
- Temperature correcction
- Summarized user manual on screen
- Graphical representation of real time process

Error messages of functions and program execution are displayed and an acoustic alarm sounds. The vapour generator automatically cuts out.

MESSAGES:

1. E 2PROM
2. RTC FAULT
3. SENSOR PRECISION
4. HEATER PROBE
5. CHAMBER PROBE
6. OVER PRESSURE
7. VACUUM
10. TEMPERATURE LOW FAULT
11. INPUT ERROR
12. STERILIZATION
13. LEVEL 1
14. LEVEL 2
15. DOOR
7. VACUUM
16. PRINTER NO PAPER

7. VACUUM 16. PRINTER NO PAPER 8. PRESSURE 17 PAPER JAM 9. BOILER TEMPERATURE 18. ELECTRICAL FAULT

99	6		
			3
-		HIDHARINI (BIRININI)	

MODELS	Part No.	Dry	Capacity litres	Ø / Height (usables) cm	Height / Width / Depth (exterior) cm	Temp. ° C	Water Reservoir H ₂ O capacity.	Basket size	Power W	Weight Kg
AUTESTER ST DRY PV III 50	4002516	Yes	50	30 63	117 48 62	105 to 134	10	3	3000	108
AUTESTER ST DRY PV III 80	4002517	Yes	80	40 60	107 58 72	105 to 134	10	2	4850	136
AUTESTER ST DRY PV III 150*	4002518	Yes	150	50 70	118 80 95	105 to 134	20	3	7500*	250

Comes without baskets and drums.

COMPLEMENTS: Factory fitted options.



Printer of continuous paper roll with print out of temperature, pressure, time and mode. Part No. 4002417

Pt100 Sensor (See page 63). Loading and unloading crane. Part No. 1002426 (see page 70).

ACCESSORIES: AISI 304 Stainless steel wire baskets and drums.

For Autester ST DRY PV-II 50, Part No: 4002515.

Capacity, 3 baskets Ø 25,5 and 20 cm. height. Part No: 1000495

For Autester ST DRY PV-II 80, Part No: 4002425.

Capacity, 2 baskets Ø 36 and 28 cm. height. Part No: 1000496

For Autester ST DRY PV-II 150, Part No: 4002426.

Capacity, 3 baskets Ø 44,5 and 22 cm. height. Part No: 1000780

Sterilization drum made from AISI 304 stainless steel.

For Autester ST DRY PV-II 50, Part No: 4002515.

Capacity, 3 drum Ø 25 and 18 cm. height. Part No. 1002418

For Autester ST DRY PV-II 80, Part No: 4002425.

Capacity, 2 drum Ø 35 and 27 cm. height. Part No. 1002419

Detergent Bio-Sel. (see page 63)

Water distillation unit especially for autoclaves. (see page 63).

PARAMETERS AND PROGRAMS PRE INSTALLED IN THE AUTOCLAVE MICROPROCESSOR (Note: All these parameters are adjustable by the user)

Nº Programme	Description	Temperature Sterilization	Time Sterilization	Time Dry	Purge
0	Free	-	-	-	-
1	Unwrapped instruments	134 °C	4 min	15 min	2 (1)
2	Wrapped instruments	134 °C	12 min	30 min	3 (2)
3	CJD	134 °C	20 min	20 min	3 (2)
4	Delicate unwrapped	121 °C	15 min	15 min	2 (1)
5	Delicate wrapped	121 °C	30 min	30 min	3 (2)
6	Disinfecting	105 °C	25 min	15 min	1 (3)
7	Liquids	121 °C	30 min	-	1 (3)
8	Solids	121 °C	15 min	30 min	1 (3)
9	Solids	134 °C	4 min	30 min	1 (3)

- (1) Double Vacuum purge.
- (2) Triple Vacuum purge.
- (3) Atmospheric purge (No vacuum).

^{*} The part number 4002426 due to its power rating is only manufactured for three phase 230V or 400V. Specify preference when ordering.

OPTIONS: for vertical autoclaves, need to be factory fitted.



Pt 100 Sensor. 2nd Probe for reading sample temperature.


To guarantee in all circumstances that the system has completed the cycle to the set time and temperature parameters, the Autester ST DRY PV II autoclave can be installed with a temperature probe that can be placed within the culture flask, residual packets etc.

When the probe is connected to the controller of the autoclave, it is automatically recognised by the microprocessor, whereby the real time temperature values detected by the probe, are displayed and printed. For autoclaves PRESOCLAVE III this Pt100 sensor only reads the temperature with an external thermometer. Part No. 4001219

Loading and unloading crane. Part No. 1002426 (see page 70).


ACCESSORIES FOR VERTICAL AUTOCLAVES





Part No.	Ø ci	Height n	Suitable for autoclaves Part No.	Number of baskets
1004771	22	14	4001745	1
1004772	22	24	4001756	1
1004771	22	14	4001757*	1
1004772	22	24	4001757 *	i
1000495	27	20	4001758 - 4002516	3
1000496	36	28	4002136 - 4002517	2
1000780	44	22	4002518	3

Note:* For autoclave Part No. 4001757 supplied with 1 basket Part No. 1004771 + 1 basket Part No. 1004772



AISI 316 stainless steel wire dividers for autoclave baskets. Part No.

1001217 4 divider compartments for basket part no. 1000495.

1001218 6 divider compartments for basket part no. 1000496.

1001222 7 divider compartments for basket part no. 1000780.



Sterilization Drums.

Part No.	Ø Height cm	Suitable for autoclaves Part No.	Capacity Drums
1002418	25 18	4001758 - 4002516	3
1002419	35 27	4002136 - 4002517	2
1002420	28 18	Not suitable for our a	autoclaves.



Water distiller, specially for autoclaves "DEST-4".

Distillation capacity: 1.5 litres/hour. Internal capacity 4 litres.

Case made from thermo resistant plastic.

Internal stainless steel.

water required Cooling by forced air though a condenser. Without glass heater elements. Easy to use. Safety cut out system if over temperature due to insufficient water. Comes complete with a 4 litre three way dispenser plastic bottle with filter. Note: for daily use it is recommended to change the filter every 3 months.

Part No.	Capacity	Conductivity	Ø Height	Power	Weight
	Distillation litre/Hour	µs/cm	cm	W	Kg
4001729	1.5	5	29 39	750	3.5

Spare: Dispenser spout with filter. Part No. 4001730



Detergent "BIO-SEL"

We recommend the use of BIO-Sel detergent for cleaning, it is also biodegradable and free from Phosphates.

Part No.

1000025 Capacity 1 litre. 1000026 Capacity 4 litres.

1000852 Box of 12 x 1 litre flasks.

1000853 Box of 4 x 4 litre flasks.

additional



Autoclave for liquids and solids sterilization "Autester ST DRY PV III"



CAPACITY: 25 LITRES. SPACE FOR 12 BOTTLES OF 250 ML, 9 OF 500 ML, 3 OF 1000 ML AND 2 OF 2000 ML. CONFORMS TO THE REGULATION - 97/23/CEE EQUIPMENT UNDER PRESSURE.

Complete selection of programs. Displays in English - Spanish - French.

FEATURES

Sterilization chamber and lid made from AISI 316 stainless steel.

External case epoxy coated.

Digital control of all autoclave operations.

Large back lit LCD 2 x 24 character display.

Air filter, purifies air before entering in the sterilization chamber.

Separate clean water reservoir and waste water.

RS-232 Interface output of parameters to a computer or for printer.

USB output for PEN-DRIVE data logging and Computer connection.

SAFETY

Safety valve.

Prevents the door being opened if under pressure or high temperature is produced. Temperature safety thermostat.

Door closed correctly detector

Thermally protected lid.

MICROPROCESSOR FUNCTIONS

Controlled parameters are:

- Sterilization temperature from 105°C up to 134°C.
- Sterilization time from 3' to 59'.
- Drying time from 20' to 60'.
- Up to 10 user-configurable programs.
- · Six operating modes:
- Solid 1 Solid 2 fractions. Solid 3 fractions Drying Liquids Liquids 2 probes
- Verification test
- Vacuum Bowie Dick.
- · Settings:
- Units languages (Spanish, English, French). USB or printer module. Clock: date/time. - External water: on/off.
- Password
- Calibration
- Temperature probe 1 Temperature probe 2 Water dosing.
- Audible and visual alarm for overtemperature and insufficient temperature.
- Network failures detection and memorization
- Safety therrmostat by software
- Sterilization cycles data logging with capacity for over 500 cycles.
- Internal cleaning program
- · Software for PC
- Temperature correcction
- Summarized user manual on screen
- Graphical representation of real time process

CONTROL PANEL

TFT touch screen with virtual pushbuttons.



- 1. Graphical representation of real time process
- 2. Start/Stop virtual pushbutton.
- 3. Running program.
- 4. Press indicator.
- 5. Temperature indicator.
- 6. Running program parameters.









MODEL	Part No.	Capacity litres	Ø/De (usable		Height/ (ext			Power W	Weight Kg
AUTESTER ST DRY PV III 25	400242	0 25	30	35	58	46	59	2510	52



Comes complete with three trays, extractor forceps, and additional support for instruments and sterile wrapping

PARAMETERS AND PROGRAMS PRE INSTALLED IN THE AUTOCLAVE MICROPROCESSOR (Note: All these parameters are adjustable by the user) Description Temperature Time Time Purge

Programme		Sterilization	Sterilization	Dry	. 3
0	Free	-	-	-	-
1	Unwrapped instruments	134 °C	4 min	15 min	2 (1)
2	Wrapped instruments	134 °C	12 min	30 min	3 (2)
3	CJD	134 °C	20 min	20 min	3 (2)
4	Delicate unwrapped	121 °C	15 min	15 min	2 (1)
5	Delicate wrapped	121 °C	30 min	30 min	3 (2)
6	Disinfecting	105 °C	25 min	15 min	1 (3)
7	Liquids	121 °C	30 min	-	1 (3)
8	Solids	121 °C	15 min	30 min	1 (3)
9	Solids	134 °C	4 min	30 min	1 (3)

- (1) Double Vacuum purge.
- (2) Triple Vacuum purge.
- (3) Atmospheric purge (No vacuum).

MESSAGES:

WILOGAULG.	
1. E 2PROM	10. TEMPERATURE LOW FAULT
2. RTC FAULT	11. INPUT ERROR
3. SENSOR PRECISION	12. STERILIZATION
4. HEATER PROBE	13. LEVEL 1
5. CHAMBER PROBE	14. LEVEL 2
6. OVER PRESSURE	15. DOOR
7. VACUUM	16. PRINTER NO PAPER
8. PRESSURE	17 Paper Jam
9. BOILER TEMPERATURE	18. ELECTRICAL FAULT

ACCESSORIES

Ink Printer (no thermal).

Paper size 2 1/4" (56 mm) wide.

Includes communications and power cable. Model

Height/Width/Depth Part No. Power Weight Kg 4120113 4 16 10 DC 5V 3A 1



Water distiller especially for Autoclaves. (See page 63). Bag sealers (See page 68).

Special tray for bottles of 0.25, 0.5, 1 or 2 litres. Part No. 4002416



Autoclave for sterilization "Autester ST DRY PV" Class B

CO REST

AUTOMATIC VACUUM DRY AND PURGE.
PROGRAMS TO 121 °C AND 134 °C (1 AND 2 bar).
MODEL HORIZONTAL. CAPACITY: 12, 18, AND 23 LITRES.

APPLICATIONS

"Class B" autoclaves are used within the field of general medicine, dentistry, veterinary, tattoo parlours and sterilization applications in general, to sterilize any material that could have been in contact with body fluids including blood.

The general capabilities of the instrument are regulated by EN 13060 that guarantees the differing sterilization applications.

The normal processes include: Solids with and without wrapping, porous material with and without wrapping and utensils that have cavities and holes (e.g. dentistry materials.)

To guarantee sterilization of these materials the "AUTESTER ST DRY PV" Class B, autoclaves, have an efficient staged vacuum, that performs a complete extraction of air that allows penetration of the steam vapour.

The same vacuum system is used to obtain the dry phase, which is required for perfect storage of sterilized materials.

The "AUTOCLAVES ST DRY PV" Class B includes verification and validation programs.

CONFORMITY:

EN 13060 SMALL GENERAL PURPOSE STEAM STERILIZERS. EN 61010-1 EQUIPMENT SAFE FOR LABORATORY USE.

EN 61010-2-040 SAFETY FOR STERILIZER EQUIPMENT FOR LABORATORY USE

EN 61326 ELECTROMAGNETIC COMPATIBILITY FOR EQUIPMENT USED IN THE LABORATORY.

SAFETY:

SAFETY VALVE: PREVENTS THE PRESSURE EXCEEDING THE MAXIMUM SET LIMIT.


SAFETY THERMOSTAT: DISCONNECTS THE HEATER AGAINST OVER TEMPERATURE OR LOW WATER IN THE RESERVOIR.

THERMALLY INSULATED LID.









18 and 23 litres

Staged fractional vacuum: Extracts air from the load that permits absolute steam vapour penetration.

Sterilization stage: Sterilization within the time and temperature selected.

Dry cycle: Extracts the humidity from the load leaving it suitable for storage.



Autoclave for sterilization Autester "ST DRY PV" 12 Class B

CAPACITY: "ST DRY PV": 12 LITRES.





According to directive: 97/23/CEE equipment under presure, 3.3.

FEATURES

Stainless steel internal chamber.

Staged vacuum purge.

Dry by vacuum

Incorporated steam generator.

Water reservoir to supply the chamber

High vapour penetration of porous material

Used in digestive studies, ophthalmology, operating theatres etc .

Printer interface.

USB adapter for parameters memorization (optional).


Safety valve. Prevents over pressure over the maximum working limit.

Safety thermostat. Disconnects power to the heater if over temperature or low water in the reservoir.

Thermally protected lid.



CONTROL PANEL



TEMPERATURE CONTROL AND AUTOMATIC MICROPROCESSOR CYCLE.

- 1. Textile materials
- **2.** 121 °C.
- 3. 134 °C.
- 4. Wrapped solids.
- **5.** Solids (Not wrapped).
- $\textbf{6.} \ Low \ level \ water \ indicator.$
- 7. Autoclave ready for operation.8. Sterilization cycle in process
- **9.** Program selection.

- 10. Select wrapped or non wrapped.
- **11.** Test.
- 12. Fast sterilization.
- 13. "maintain hot" the chamber.
- 14. Preheat.
- 15. Pressure indicator (bar).
- **16.** Temperature indicator (°C).
- 17. Time Indicator (min.).

CYCLES	134 °C Wrapped	134 °C Not Wrapped	121 °C Wrapped	121 °C Not Wrapped	Textiles	Test Bowie Dick	Test Vacuum
Sterilization Temperature	134 °C	134 °C	121 °C	121 °C	134 °C	134 °C	-
Sterilization pressure	2.1 bar	2.1 bar	1.1 bar	1.1 bar	2.1 bar	-	2.1 bar
Vacuum purge time	8 min.	4 min.	8 min.	4 min.	8 min.	-	-
Sterilization time	12 min.	4 min.	30 min.	20 min.	12 min.	3.5 min.	10 min.
Dry time	10 min.	10 min.	10 min.	10 min.	15 min.	1 min.	-

MODEL	Part No.	Capacity litres	Ø / Depth cm	External dimensions cm Height Width Depth	Trays dimensions cm Height Width Depth	Power W	Weight kg
AUTESTER ST DRY PV B 12	4001727	12	20 35	39,2 48,7 60	1,8 14 29,5	1800	42,2



Comes complete with:

spare joint, door adjustment key, support with three trays, extractor forceps, empty water filter and waste water out.

ACCESSORIES

USB adapter model.

Pen-Drive included (Memory board) for data storage.

Part No: 4120130
Thermal Printer.

Paper size 2 1/4" (56 mm) wide.

Includes communications and power cable.

Power: DC 5V 3A Part No: 4120123

Water distiller especially for Autoclaves. (See page 63).

Bag sealers (See page 68).





Autoclave for sterilization "Autester ST DRY PV" 18 Class B "ST DRY PV" 23 Class B





CAPACITY: "ST DRY PV" 18: 18 LITRES. "ST DRY PV" 23: 23 LITRES..

According to directive: 97/23/CEE equipment under presure, 3.3.

FEATURES

Digital electronic control. Stainless steel internal chamber.

Staged vacuum purge.

Dry by vacuum

Incorporated steam generator.

Water reservoir to supply the chamber

High vapour penetration of porous material

Used in digestive studies, ophthalmology, operating thea-

Automatic or manual filling of the water tank.

Printer interface.

USB adapter for parameters memorization (optional).

SAFETY

Mechanism door lock, prevents the door being opened accidentally.

Safety valve. Prevents over pressure over the maximum working limit.

Safety thermostat. Disconnect power to the heaters if over temperature or low water in the reservoir are produced.


Thermally insulated lid.

Dirty water tank full alarm.

CONTROL PANEL

- 1. LCD screen.
- 2. Push button increase parameters.
- 3. Push button decrease parameters.
- 4. Push button menu.
- 5. Push button Start-Stop.





CYCLES	134 °C Emergency	134 °C Solid	134 °C Hollow	134 °C Sanguinary	121 °C Solid	121 °C Hollow	134 °C Glass	121 °C Tampon	134 °C Bowie Dick	134 °C Helix	134 °C User-defined	Test Vacuum
Starilization Tomporature	134 °C	134 °C	134 °C	134 °C	121 °C	121 °C	134 °C	121 °C	134 °C	134 °C	134 °C	vacuum
Sterilization Temperature											134 G	-
Sterilization pressure	2,1 bar	2,1 bar	2,1 bar	2,1 bar	1,1 bar	1,1 bar	1,1 bar	1,1 bar	2,1 bar	2,1 bar	-	-
Vacuum purge number	1	1	4	3	1	3	1	3	3	3	3	-
Sterilization time	4 min.	4 min.	4 min.	18 min.	20 min.	20 min.	4 min.	30 min.	3,5 min.	3,5 min.	4 min.	-
Dry time	8 min.	4 min.	10 min.	10 min.	8 min.	10 min.	4 min.	18 min.	3min.	3 min.	3 min.	-

MODEL	Part No.	Capacity litres	Ø / depth cm	External dimensions cm Height Width Depth	Trays dimensions cm Height Width Depth	Power W	Weight kg
AUTESTER ST DRY PV 18	4001734	18	24,5 36	40 47 83	2 20,5 28	1500	46,6
AUTESTER ST DRY PV 23	4001735	23	24,5 47	40 47 94	2 20,5 38	1700	48



Comes complete with:

spare joint, door adjustment key, support with three trays, extractor forceps, empty water filter and waste water out.

ACCESSORIES

USB adapter model.

Pen-Drive included (Memory board) for data storage.

Part No: 4120130

Ink Printer (no thermal).

Paper size 2 1/4" (56 mm) wide.

Includes communications and power cable.

Power: DC 5V 3A Part No: **4120123**

Water distiller especially for Autoclaves. (See page 63).

Bag sealers (See page 68).





Bag sealer "S-400"

SUITABLE FOR POLYPROPYLENE-POLYESTER STERILISING BAGS WITH A PAPER BASE.

FEATURES

Metal epoxy covered case.

Pressure sealer. Manual cut.

Stainless steel AISI 304 reel support at the back.

CONTROL PANEL

Illuminated mains switch.

Lamp on when the thermostat has reached the pre-set working temperature.

MODEL

Part No.	Seal Width	Maximum seal	Height / Width / Depth	Power	Weight
	mm	length mm	(exterior) cm	W	Kg
4100041	9	270	12 35 21	200	5.3





Bag sealer "Seal-1000"

FOR POLYPROPYLENE AND POLYESTER BAGS WITH A STERI-LE PAPER BASE.

FEATURES

Metal construction with stainless steel tray.

Pressure sealer unit with lateral leaver with two bilateral cutting edges.

Constant sealing temperature (from 180 $^{\circ}\text{C}$ and 200 $^{\circ}\text{C}).$

External support for rolls of up to 240mm wide, with adjustable separator sides to accommodate different roll sizes.

CONTROL PANEL

Illuminated start switch.

Indicator, optimum seal temperature.

Part No.	Seal width	Maximum seal	Height / Width / Depth	Power	Weight
	mm	length mm	(exterior) cm	W	Kg
4101712	10	260	30 36 43	400	6.5





Automatic bag sealer "Seal-Com 2"

IN ACCORDANCE WITH DIN 58953 SECTION 7 STANDARD. FOR SEALING STERILE PAPER, POLYPROPYLENE-POLYESTER BAGS. CONTINUOUS SEALING, AUTOMATIC ADVANCE.

Continuous Thermosealer

FEATURES

Exterior made from stainless steel. Heating time: 4 minutes.

Seal width: 10 mm.

Sealing speed: 6 m/min.

Digital control and display of tempera-

ture and time.

Acustic sign of cycle starting.

CONTROL PANEL

Luminous "ON" switch. Luminous "Manual or Automatic"

switch.



MODEL

Part No.	Height / Width / Depth (sealer) cm	Power W	Weight Kg	
4100044	16 48 26	400	13	

SEALING MACHINE ACCESSORIES



Rolls of paper for sterilization "Polisteril"

Paper base polypropylene and polyester bags showing differing ink prints of sterilized items.

Part No.	Bag width cm	Bag length m
4100075	7,5	100
4100100	10	100
4100200	20	100



Sterile paper "Policort" trimmer

Completely made of metal with back and reel support, for a capacity of 2 x 200mm or 4 x 100 mm wide reels. Circular cutting knife edge.

Measurements: 13 high x 51 wide x 32 cm deep.

Weight: 2.8 Kg. Part No. **4100045**



Date stamping machine for sterile bags.

Part No. **4100046**

STERILIZED 30 OCTOBER 2015 Valid for 3 Months.

OPTIONS: Ultrasonic cleaning baths for dentistry, "Ultrasons-Digit" (see page 96). Dry sterilizer "DRYTIME" poupinel oven (see page 137).



Autoclaves for preserved food pasteurization/sterilization "STERIL-FOOD" 20, 80 and 150

SELECTABLE TEMPERATURE BETWEEN 70 °C & 121°C.

CAPACITY: 20, 80 & 150 L.

AUTOCLAVES "MAN" MODEL NEED A COMPRESSED AIR EXTERNAL INSTALATION.

AUTOCLAYES "COM" MODEL INCLUDE COMPRESOR AND DO NOT NEED COMPRESSED AIR INSTALLATION.

WATER SUPPLY NEEDED WITH LOW-CAL CONTENT AND DRAINAGE. (WATER PRESURE: 3 BAR MINIMUM).



NEW

CONFORMITY EN 61010-1 SAFETY FOR EQUIPMENTS FOR LABORATORY USE.

EN 61010-2-040 SAFETY FOR STERILIZER EQUIPMENTS. EN 61326 ELECTROMAGNETIC COMPATIBILITY. EN/97/23/EC SAFETY FOR PRESSURE EQUIPMENTS.

APPLICATIONS

Preserved food sterilization.

Very suitable for small productions and wide variety.

FEATURES

External case, top, tank and lid made of stainless steel.

Tank fill by external water inlet.

Temperature control probe inside the sample.

Counter pressure system to prevent the cans breakage.

Temperature, time and pressure registered in USB.

SAFETY

Device that avoids opening the lid when the chamber is under pressure.

Safety valve: Prevents overpressure over the maximum limit.

Safety pressure switch: Disconnects heater when there's overpressure.

Thermally insulated lid.

Door properly closed detector.

Double filter retention

CONTROL PANEL

- 1. Parameters LCD display.
- 2. Microprocessor single button.
- 3. Operation Led.
- 4. Start/stop push button.
- 5. USB adapter.
- 6. Open lid knob.
- 7. USB function light indicator.



MICROPROCESSOR FEATURES

Microprocessor controls sterilization and cooling process:

- Pasteurization temperature from 70°C to 100°C.
- Sterilization temperature from 100°c to 121°C.
- Sterilization time from 1' to 99':
- Selectable cooling temperature.
- Programs from 0 to 9 user configurable.
- Configuration: Units. Languages (Spanish, English, French).
 - Timer: date/hour. External water: on/off.
- Password
- Calibration Control probe temperature.
 - Sample probe temperature

Any error message of functions or execution is displayed in the panel and an acoustic alarm sounds, automatically cutting the steam generator out.

BENEFITS

Water quick cooling.

High efficiency in heat treatment

Respect the organoleptic properties of food.

Easy and quick changing of heat treatment parameters.

10 user defined programs.

Models with built-in compressor.





- 1. Mains input.
- 2. External water supply.
- 3. Safety valve.
- 4. Vacuum steam valve.
- 5. Safety thermostat.
- 6. RS232 connector.
- 7. Second probe connector.
- 8. Compressed air inlet.

MODELS	Part No.	Drying	Capacity litres	Ø / Height (usable) cm	Height/Width/ Depth (external) cm	Temp. ° C	Capacity baskets	Power W	Weight Kg
STERIL-FOOD COM 20	4002421	NO	20	24 41	70 38 38	70 to 121	2	2500	55
STERIL-FOOD MAN 80	4002427	NO	80	40 60	107 58 72	70 to 121	2	3600	136
STERIL-FOOD COM 80	4002422	NO	80	40 60	107 58 72	70 to 121	2	3600	136
STERIL-FOOD MAN 150*	4002428	NO	150	50 70	118 80 95	70 to 121	3	7500*	250
STERIL-FOOD COM 150*	4002423	NO	150	50 70	118 80 95	70 to 121	3	7500*	250

It is supplied without baskets or drums. *Part numbers 4002428 and 4002423, due to their power consumption are manufactured for three phase 230V or 400V. Specify preference when ordering.

COMPLEMENTS FOR "PRESOCLAVE III 50 - 80", "AUTESTER ST DRY PV III 50 - 80 - 150" AND "STERIL-FOOD 80 - 150"

Loading and unloading crane

Baskets elevator for loading and unloading of autoclaves of 50, 80 and 150L With adapter for handling baskets without deforming.

Axial rotation for placing baskets to a trolley. Total lift height: 2.4 m (autoclave assembled)

Outside completely stainless steel.

Maximum load: 50Kg. Power: 230V 300W Must be factory installed.

Control panel: ON / OFF button + UP / DOWN pushbutton integrated into a sin-

gle command. Part No: 1002426



ACCESORIES FOR STERIL-FOOD

Decalcifier "C-3"

APPLICATIONS

Pre-treatment of water with high content of lime, for general laboratory use.

Recommended for distillation units use, from 25° French water hardness.

TECHNICAL DATA

AISI 304 stainless steel metallic tank with manual device of salt regeneration by two-way wrenches.

Resin capacity: 12 litres. Regeneration by salt: 2 Kg. Cyclical output by regeneration: 1200 litres at 60° French degrees / 4800 litres at 35° French degrees.

Cyclical output by regeneration for distillation: from 300 to 800 litres.

Maximum mains inlet pressure: 4,5 Kg/cm2.

Maximum input permissible hardness: 60° French degrees. Output hardness: 1° French degree.

External measurements: 62 height x 19 cm Ø. Weight: 20 Kg. Part Number 0703052



It comes complete with 34 in and out hoses.

AISI 304 stainless steel wire baskets. Maximum load capacity: 8 kg.

Part No.	For autoclave	Ø / Height	Capacity	Bottle capacity	
	litres	(Usable) cm	baskets	Ø72x85mm	Ø82x115mm
1000496	80	36 x 28	2	55	55
1000780	150	45 x 22	3	162	150
1000786	150	47 x 15	5	150	95

Note: To optimize the load in the autoclave, we manufacture baskets on demand suitable to the size and weight of your containers.







Flake ice maker **73** page Refrigerated cabinets "Medilow" S, M, L & LG **73** page Refrigerated medical cabinets, S, L & M 74 page **Upright freezers "Templow" S, M, L, P and Z** 75 and 76 page Storage refrigerators "Stocklow" LS & GS **77** page **Chest freezers "Combatemp" 77** page



LOW TEMPERATURE

Model selection table

QUICK OVERVIEW

Equipment type	Model	Forced Air	Regulation		CAPACITY		
Refrigerated Cabinets				180 litres	236 litres	544 litres	544 litres Glass
+1 +40 °C	Medilow	Yes	Digital	2101270	2101271	2101272	2101273
			J				
Storage refrigerators				236 litres	352 litres	319 litres	
+0 +15 °C	Pharmalow	Yes	Digital	2101274	2101275	-	
+0 +15 °C	Pharmalow	Yes	Digital	-	-	2101501	
Storage freezers				513 litres			
-1425 °C	Templow-M	No	Digital	2101279			
				304 litres			
-1525 °C	Templow-S	No	Digital	2101278			
				600 litres			
-1035 °C	Templow-L	Yes	Digital	2101280			
				345 litres (Two compartments)			
+330 °C	Templow-P	Yes	Digital	2101285			
				252 litres (Two compartments)			
+240 °C	Templow-Z	Yes	Digital	2101282			
Storage refrigerators				600 litres	600 litres Glass		
+1 +15 °C	Stocklow LS	Yes	Digital	2101283	_		
+1 +15 °C	Stocklow GS	Yes	Digital	-	2101284		
	ottomon do	100	2.9				
Chest Freezers				206 litres	215 litres	343 litres	352 litres
-1428 °C	Combatemp	No	Digital	2101040	-	2101041	-
-1040 °C	Combatemp	No	Digital	-	2101181	-	2101182
Flake ice maker				40 litres	85 litres	150 litres	
. m. o. oo maroi	"IMS 45"	Compressor	-	5020100	-	-	
	"IMS 85"	Compressor	-	-	5020101	_	
	"IMS 150B"	Compressor	-	_	-	5020103	
	11113 1000	001114110001				3020100	

OTHER LOW TEMPERATURE DEVICES

Thermocycler for thermal cycles PCR-K96



(See page 238)

Shaking Bath Unitronic



(See page 93)

Cryo circulators Frigiterm-TFT



(See page 103)

Bath immersion coolers Digit-cool



(See page 105)

Refrigerated cabinet Refricub



(See page 218)

Cold Plate Plac-center



(See page 217)

Cooling Coils for baths Frigedor



(See page 107)

Microtube, table top cooler boxes for microtubes of 0.5 to 2 ml / 15 ml.



(See page 190)

Precision Peltier cooled incubator / cabinets **Prebatem**



(See page 155 and 156)



cabinets Hotcold



Incubator cabinets **Blood bank**



Flake ice maker "IMS" 40, "IMS" 85 and "IMS 150" "IMS" 40 MODEL: PRODUCTION 40 Kg-24 HOURS.

"IMS" 85 MODEL: PRODUCTION 85 Kg-24 HOURS. "IMS" 150B MODEL: PRODUCTION 150 Kg-24 HOURS.

MODELS	Part No.	24h.(kg) Production	Tank Capacity (Kg)	Height /	Width / Depth cm	Power W	Weight Kg
IMS-40	5020100	40	8	72	38 54	280	40
IMS-85	5020101	85	13	88	55 61	420	57
IMS-150B	5020103	150	35	98	71 69	695	70

Features: See page. 191.





Refrigerated cabinets "Medilow" S, M, L & LG with forced air fan circulation.

WITH DIGITAL ELECTRONIC TEMPERATURE CONTROL.

TEMPERATURE RANGE CONTROLLABLE FROM +2 °C TO +40 °C. WITH A STABILITY OF ±1.5 °C. RESOLUTION 1 °C. **CAPACITY: 180, 236 AND 544 LITRES.**

SAFETY: CONFORMS TO EN61010.1. OVER TEMPERATURE SAFETY THERMOSTAT INSTALLED.



Models Medilow M and L. Medilow M. Capacity: DBO 3 blocks



Model Medilow LG, double glass door.



Model Medilow S. Capacity: DBO 2 blocks

APPLICATIONS

Determination of BOD at 20°C. Measurement of enzyme activity at 25°C. Storage and Preservation of samples, Incubation, Cosmetics, Botany, Food Products, etc.

TEMPERATURE CONTROL

Electronic temperature control with digital display. Pt 100 temperature probe. Probe failure alarm.

FEATURES

Epoxy covered external case. Easy to clean, strong, molded ABS Interior. Vibration and CFC free hermetically sealed compressor.

Turbine fan air circulation, enclosed within the unit. External hole for connecting inside probes, etc. Internal power connection to a maximum of 500 W. Adjustable level front supports.

Door can be hung to open left or right. Door lock.

CONTROL PANEL

Temperature control. Live-time temperature display. Pre-selected temperature display. Increase value, push button. Decrease value, push button. Program parameter selection, push button. Configure push button.

MODELS	Part No.	Capacity litres	Door glass	Height / (int	Width / terior) c		Heig	ht / W (exter		/ Depht cm	Number of shelves	Internal power points	Power W	Weight Kg
MEDILOW-S	2101270	180	-	70	51	44	8	5	60	60	3	2	250	36
MEDILOW-M	2101271	236	-	106	47	44	12	5	60	61	3	2	330	47
MEDILOW-L	2101272	544	-	145	60	56	16	4	75	73	4	6	330	73
MEDILOW-LG*	2101273*	544	YES	145	60	56	16	4	75	73	4	6	330	82

Comes complete with steel shelves, PVC laminated.

104	, ,		-	U	000	0_
SPARE PARTS	S She	lves				
For Medilow S/M/	/L/LG	Part No.	2101270	2101271	2101272/21	01273
Shelves			1001813	1001814	10018	15

^{*}Anti mist double glazed insulated door, with fluorescent interior light.



Refrigerated medical cabinet "Pharmalow" S, L and M

WITH FAN ASSISTED CIRCULATION. TEMPERATURES RANGE FROM +0 °C UP TO +15 °C.

APPLICATIONS

Products that require a constant below atmosphere temperature, preservation of pharmaceutical products, medications, micro-tubes, etc.

COMMON FEATURES

External case and metallic door.

White coated and interior made of anti-toxic white

Double glazed door and interior light for model M. Hermetically sealed compressor.

Homogeneous, turbine fan assisted air circulation.

Automatic de-frost with water evaporation system. Adjustable level, front base supports. Door can be mounted to open left or right. Door lock.

CONTROL PANEL

Temperature thermostat regulator. Digital thermometer display for models S and L Electronic digital temperature control model M.

Refrigerator "Pharmalow" S and L

MODELS	Part No.	Capacity litres	Height / ! (into	Width erior) (Height / Width / Depht (exterior)cm		Number of shelves	Power W	Weight Kg	
Pharmalow S	2101274	236	106	47	44	125	60	61	4	170	47
Pharmalow L	2101275	352	145	47	44	164	60	61	4	180	56

removable methacrylate partitions.

Frame support, partitioned metallic drawers.

partitions in total, (maximum capacity: 6 drawers).

Comes complete with steel shelves, PVC laminated.

ACCESSORIES Need to be factory fitted.

Model Pharmalow S and L

with frame support and drawers.

Frame support for Pharmalow L (Part No. 2101275), with 5 drawers and 60 partitions in total, (maximum capacity: 10 drawers).

Modular system with perforated base, sliding drawers, each drawer is divided into three fixed compartments, each compartment with four

Frame support for Pharmalow S (Part No. 2101274), with 4 drawers and 48

Part No. 1001307

Part No. 1001306

OPTIONAL

Drawer with 3 compartments and 12 partitions. Part No.1001308



Details of frame support, drawers, dividers and partitions.

Storage refrigerator with double glazed anti mist door "Pharmalow" M WITH DIGITAL ELECTRONIC TEMPERATURE CONTROL

MODEL

Part No.	Capacity litres		/ Width nterior)	n / Depht cm		Ancho A	/ Fondo cm	Power W	Power HP	Number of shelves	Weight Kg
2101501	319	139	48	45	198	60	64	360	1/4	6	98

Comes complete with steel shelves, PVC laminated.

ACCESSORY "Pharmalow" S - L and M. Must be factory fitted.

Circular graphic register independently powered.

Temperature range: -10 °C up to 40 °C. Precision ±2% at 12 °C. Chart speed selection of 24 hours or 7 days. Disc Ø 125mm Graphic register dimensions:

18 cm high x 22 cm wide x 35 cm deep.

Powered by an LR-6 1.5V battery. Battery life 1 year.

Part No. 2101621





Model Pharmalow S with 4 shelves.

SPARE PARTS Shelves	
For Pharmalow S/L. Part No.	2101274/2101275
Shelves	1001814



SPARE PARTS

Part No

Brackets (4 units).



Upright freezers "Templow" S, M,L,P and Z.



APPLICATIONS

Preservation and storage of products that require to be preserved at low temperatures.

Upright freezer "Templow" S

DIGITAL TEMPERATURE CONTROL, TEMPERATURE RANGE FROM -15 °C TO -25 °C. 8 DRAWERS WITH INDIVIDUAL FREEZER COMPARTMENTS.

FEATURES

Epoxy covered metallic external case. Strong, molded, chemically resistant ABS interior. Manual defrost with water evaporation system. Door can be hung to open left or right.

CONTROL PANEL

Mains power ON/ OFF switch.
Temperature thermostat regulator.
Temperature indicator lamp.
Audible and visual alarms.
Push button alarm reset/stop.
Push button with quick freeze indicator.



Upright freezer "Templow" M

TEMPERATURE RANGE FROM -14 °C TO -25 °C.
7 INDIVIDUAL FREEZER COMPARTMENTS WITH 14 BASKETS.

FEATURES

Epoxy covered external case. Strong, chemically resistant, ABS molded case. Specially coated "Rilsan", wire mesh baskets. Manual defrost, with water evaporation system. The door can be hung to open left or right. Lockable door.

CONTROL PANEL

Temperature thermostat regulator. Digital temperature display.



Upright freezer "Templow" L

WITH FAN ASSISTED AIR CIRCULATION. TEMPERATURE RANGE FROM -10 °C UP TO -35 °C. COMES WITH 4 SHELVES.

FEATURES

External and internal case made completely of AISI 304 stainless steel, with water sump and drain. Chemically resistant.

Acclimatized hermetically sealed CFC free compressor. Automatic defrost with electrical water evaporation system.

Regulable shelves.

Adjustable front level supports.
The door can be hung to open, left or right.
Door lock.

CONTROL PANEL

Mains power ON/ OFF switch. Electronic temperature regulator. Digital temperature display.



Comes with steel PVC laminated shelves and baskets.



Refrigerated cabinets with double compressor "Templow" P

FOR TEMPERATURES FROM 3°C UP TO 8°C IN UPPER CASE. FOR TEMPERATURES FROM -9°C UP TO -30°C IN LOWER CASE. WITH 4 TRAYS IN UPPER CASE AND 3 DRAWERS IN LOWER CASE.

FEATURES

External case and metallic door.

White coated and interior made of anti-toxic white sheet.

Double compressor which offers maximum temperature stability.

Gradual adjustment of the trays for a versatile use of the interior.

Door with lock and key.

CONTROL PANEL

Digital electronic control.

Alarms of low and high temperature, door open and power failure.

Parameters configuration.

Upper case temperature indicator display.

Lower case temperature indicator display. Parameters

value increase/decrease.

Last 3 warning signs pre-setting. 10 mm diameter hole for PT100 probe.



Refrigerated cabinet "Templow" Z

FOR TEMPERATURES FROM 2 °C UP TO 9°C IN UPPER CASE. FOR TEMPERATURES FROM -10 °C UP TO -40°C IN LOWER CASE. WITH 3 TRAYS AND 1 DRAWER IN UPPER CASE. WITH 4 DRAWERS IN LOWER CASE.

FEATURES

External case and interior completely made of ABS resistant to chemical agents.

Hermetically sealed tropicalized compressor.

Automatic defrosting with water evaporation by temporized electrical defrost.

Supports for height adjustable shelves.

Door lock.

CONTROL PANEL

High precision electronic control with digital indicator. Password lock.

Temperature printout (last 7days).

Parameters configuration.

Upper case temperature indicator display.

Lowe case temperature indicator display.

Parameters value increase/decrease.

All adjusted parameters activation.



MODELS	Code	Capacity litres	Range °C	Height/Width/Depth (interior) cm	Height/Width/Depth (exterior) cm	Power W	Weight Kg
TEMPLOW P	2101285	Upper 240 Lower 105	Upper 3 to 8 Lower -9 to -30	Upper 110 44 43 Lower 59 43 43	200 60 61	345	85
TEMPLOW Z	2101282	Upper 153 Lower 99	Upper 2 to 9 Lower -10 to -40	Upper 80 47 46 Lower 41 39 36	177 62 60	400	78

Templow P: Comes complete with 4 PVC steel coated trays and 3 plastic drawers.

Templow Z: Comes complete with 3 glass trays and 5 plastic drawers.



Storage refrigerators "Stocklow" LS and GS with fan Assisted Circulation.



MODEL L, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C. MODEL G. WITH DOUBLE GLAZED ANTI MIST DOOR.TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C.



APPLICATIONS

Preservation and storage of products that require below ambient temperatures.

FEATURES

Chemically resistant.

Tropical hermetic compressor.

Automatic defrost with electrical water evaporation system.

Support rack and shelves.

Adjustable front level supports.

The door can be hung to open left or right.

Door lock.

There are two versions:

Stocklow LS which is completely made of metal, both interior and exterior, AISI 304 stainless steel door.

Stocklow GS inside made of AISI 304 stainless steel and epoxy white external case with glazed door that enables internal view.

CONTROL PANEL

Mains power, ON/OFF switch. Electronic temperature regulator. Digital temperature display. Interior light only (Double glazed door model).







Model Stocklow GS, with metallic door.

MODELS	Part No.	Capacity litres	Door	Height / Width / Depht (interior) cm	Height / Width / Depht (exterior) cm	Number of shelves	Power W	Weight Kg
STOCKLOW LS	2101283	600	Metallic	138 62 68	210 73 84	3	250	120
STOCKLOW GS	2101284	600	Glass	138 62 68	210 73 84	3	250	125

Comes complete with steel shelves, PVC laminated.



Chest Freezers "Combatemp"

CONTROLLABLE TEMPERATURE RANGE FROM -14 °C TO -28 °C. CONTROLLABLE TEMPERATURE RANGE FROM -10 °C TO -45 °C.

APPLICATIONS

Suitable for preservation and freezing of sample.

FEATURES

Epoxy covered external case. Interior of aluminium. with quide supports for baskets or travs. The interior of the door is made of easy to clean. strong, chemically resistant, moulded ABS. Drop down lockable door with compensating support hinges. Hermetically sealed compressor, silent and vibration free

Automatic interior light when the door is opened. (all models except Combatemp

Digital electronic temperature control. N.T.C. Sensor.

CONTROL PANEL

Temperature control. Digital live and preset temperature display.

Push-hutton temperature selector.



ione and vibration moo.			r don button temperature ecretic.				-					
MODELS	Part No.	Capacity litres	Freezing capacity in 24hrs Kg	Temperature up to °C		nt / Width (interior)		•	/ Width xterior)	/ Depht cm	Power W	Weight Kg
COMBATEMP S-30	2101040	206	28	-28	74	63	53	92	75	71	160	38
COMBATEMP L-30	2101041	343	38	-28	74	100	53	92	113	71	200	56
COMBATEMP S-40	2101181	215	30	-40	63	89	41	92	113	76	230	57
COMBATEMP L-40	2101182	352	40	-40	65	117	50	92	137	81	400	72

Note: Part No. 2101040 is supplied with 1 basket. Part No. 2101041 is supplied with 2 baskets.

SPARE PARTS Baskets			
For Combatemp. Part No.	2101040 / 2101041	2101181	2101182
Baskets	1001810	1001811	1001812



Shakers for platelet bags 79 page **Incubation chamber 79** page **Platelet bags incubator cabinets** 80 page **Blood bank refrigerated cabinets** 81 to 82 pages **Shaking monitor 82** page Plasma bags defrost baths page 83 **Automatic sealer and plasma extractor** 83 page



Horizontal shakers for platelet bags "APL-54" and "APL-108"



MODEL "APL-54": CAPACITY FOR 54 PLATELET BAGS.
MODEL "APL-108": CAPACITY FOR 108 PLATELET BAGS.



NEW

Platelet bags constant shaking for long periods of time.





FEATURES

Epoxy covered case with treated aluminium platform with AISI 304 stainless steel racks for bags.

Stirring automatic restart in 30 seconds and pause device for bags handling.

Acoustic alarm in case of power failure.

LCD display informing of the events log: it counts the number of power failures and their duration.

When power is restored, the device automatically restarts the stirring.

CONTROL PANEL

Single button to adjust:

- -Date and time programmation.
- -Shaking stop button for 30 seconds.
- -Alarms reset.

- -Moving to review the number of incidences.
- -Shaker stop time programmation.

MODELS	Part No.	№ of Shelves	Capacity bags	Height / Width / Depth Cm.	Power W	Voltage V	Weight Kg
APL-54	3000980	9	54	40 44 39	50	230	22
APL-108	3000981	18	108	40 84 38	50	230	35

It can be supplied for 110 V- 60 Hz on demand.



Incubation chamber "PLX-IN"

FOR PLATELET BAGS "APL-54" SHAKING, CAPACITY 54 PLATELET BAGS. FOR REGULABLE TEMPERATURES FROM 18 °C UP TO 60 °C. RESOLUTION 0,1 °C. STABILITY \pm 0,1 °C. HOMOGENEITY \pm 0,5 °C. SET ERROR \pm 0,3 °C.

Peltier effect. No compressor group. Quiet, vibrations free, high accuracy and low consumption.

FEATURES

Made of transparent methacrylate with swinging door that allows an easy access to the working area and an inside chamber view during bags shaking.

AISI 304 stainless steel lower base.

Auto shut-off protective device in case of door opening and automatic restart in 30 seconds.

CONTROL PANEL

Start-up main switch with digital temperature regulation.

MODEL

Part No.	Height / (ex	Width terior)		Height / Width (interior)			Power w	Weight Kg
3000958	44	63	65	33	60	50	360	38

It can be supplied for 110 V- 60 Hz on demand.

Note: For a correct bag preservation, it is recommended to keep a temperature between 20 and 24 °C.



Incubation chamber "PLX-IN" Part No. 3000958 with "APL-54" shaker Part No. 3000980



Platelet bags incubated cabinets "Blood Bank" E and F

EXCLUSIVE DESIGN ADAPTED TO ACCOMODATE PLATELET BAGS STIRRERS "APL". FORCED AIR CIRCULATION.







SAFETY:

LOW TEMPERATURE THERMOSTAT LIMITER.

VISUAL AND ACOUSTIC INDICATORS: EXCESSIVE LOW TEMPERATURE (+18 °C), EXCESSIVE HIGH TEMPERATURE (+24 °C), OPEN DOOR, MAINS FAILURE AND BATTERY FAILURE.

TEMPERATURE GRAPHIC REGISTER. 7 DAYS RECORDER, AUTONOMY 6 MONTHS. ACOUSTIC ALARM ON ALL FAIL CONDITIONS. ACOUSTIC INDICATOR.









APPLICATIONS

Storage cabinet specially designed for platelet bags preservation.

COMMON FEATURES

AISI 304 stainless steel external cabinet and internal lining.

Double glazed door, reversible lockable door, easy to change, with automatic feedback mechanism and manual locking.

Hermetically sealed compressor mounted on antivibration settlement with forced type evaporator and ventilated type condenser.

Shelves and supports in steel PVC laminated.

2 power outlet base for model "Blood-Bank E" max. power 500 W 6 power outlet base for model "Blood-Bank F" max. power 500 W

Inner air circulation for a more temperature homogeneity.

Refrigerant R134a.

Temperature sensors located in similar conditions to platelet bags. Circular graphic register built-in.

Microprocessor temperature control that allows configuration and visualization of the following parameters by means of the display:

- Entry code.
- External alarm signal.
- SET adjustment.
- Alarm adjustment.
- Calibration.

CONTROL PANEL

- 1. Digital display
- 2. Battery life indicator.
- 3. Open door indicator.
- 4. Mains failure indicator.
- 5. Increase value button.
- 6. Decrease value button. 7. Validate / RESET alarm.
- 8. Seven days graphic register.



OPTIONAL COMPLEMENT

Alarm messages transmision

Automatic transmission system of alarm or incident SMS messages by using a GSM module SIM card activated. Ability to connect up to 2 phones.

GSM Module. Part No: 2101622



MODELS	Part No.	Capacity litres	Capacity stirrers	Height / Width / Depth (interior) cm		Height / Width / Depth (exterior) cm		Power W	Power HP	Number of shelves	Weight Kg		
Blood Bank E	2101509	600	2 x APL-54	136	58	69	207	70	82	350	3/8	2	115
Blood Bank F	2101510	1300	4xAPL-54 o 2xAPL-108	136	128	69	208	140	82	1100	3/8	2	190

They are supplied with trays.



Blood bank refrigerated cabinets "Blood Bank" A, B, C and D

FORCED AIR CIRCULATION. FIXED TEMPERATURE +4 °C. STABILITY ±1.5 °C. RESOLUTION 1 °C. SEVEN DAYS GRAPHIC TEMPERATURE REGISTER. CAPACITIES FROM 90 TO 480 BLOOD BAGS.



SAFETY:

ALARMS SYSTEMS VISUAL AND ACOUSTIC: LOW TEMPERATURE $(+2^{\circ}C)$, OVER TEMPERATURE $(+6^{\circ}C)$, DOOR OPEN, MAINS FAILURE AND BATTERY FAILURE. INDEPENDENT BATTERY BACKUP GRAPHIC TEMPERATURE REGISTER. 7 DAY RECORDER, AUTONOMY 6 MONTHS. ACOUSTIC ALARM ON ALL FAIL CONDITIONS. REMOTE ACOUSTIC ALARM.







Models A and B with optional drawer rack.

Model C

Model D

APPLICATIONS

AISI 304 storage cabinets specifically designed and manufactured for blood bag storage.

COMMON FEATURES

Stainless steel. AISI 304. external cabinet and interior lining. Double glazed antimist, reversible lockable door, easy to change, self close door if left ajar. Hermetically sealed compressor mounted on shock absorbers with evaporator and forced air condenser.

Fan air circulation maintaining homogeneous temperature distribution. Refrigerant R134a.

Temperature sensors are located and calibrated to read simulated blood bag conditions.

Independently powered circular graphic register.

Microprocessor controlled temperature control via a display that can configure the following parameters:

- Open code.
- Activate external alarm.
- Adiust SET.
- Adjust alarm.
- Calibration.

CONTROL PANEL

- 1. Digital display
- 2. Battery life indicator.
- 3. Open door indicator.
- 4. Mains failure indicator.
- 5. Increase value button.
- 6. Decrease value button. 7. Validate / RESET alarm.
- 8. Seven days graphic register.



OPTIONAL COMPLEMENT

Alarm messages transmision

Automatic transmission system of alarm or incident SMS messages by using a GSM module SIM card activated. Ability to connect up to 2 phones.

GSM Module. Part No: 2101622



MODELS	Part No.	Capacity litres	Capacity bags	Hei	jht / Widt (interior)	h / Depht cm	Не		/ Width xterior)	/ Depht cm	Power W	Power HP	Number of shelves	Weight Kg
Blood Bank A	2101497	319	110	139	48	45	19	98	60	64	360	1/4	5	98
Blood Bank B	2101498	442	150	126	58	55	19	38	70	71	360	1/4	5	111
Blood Bank C	2101499	600	200	136	58	69	2)7	70	82	350	3/8	5	115
Blood Bank D	2101500	1300	480	136	128	69	21	18	140	82	1100	3/8	10	190

Comes complete with steel shelves, PVC laminated.



OPTIONAL ACCESSORIES,

Metallic drawers, to be ordered separately.

Modular rack system of drawers, divided into compartments and subdivided by plastic dividers for different bag classification, bags are placed upright longitudinal.

System of 5 drawers of 18 bags per drawer. Capacity 90 bags. Part No 1001501

Blood Bank B,

System of 5 drawers of 30 bags per drawer. Capacity 150 bags. Part No. 1001502

Blood Bank C,

System of 5 drawers for 40 bags per drawer. Capacity 200 bags. Part No. 1001503

Blood Bank D,

System of 12 drawers (6 in each compartment) for 40 bags per drawer. Capacity 480 bags. Part No. 1001504



Detailed image of the drawers with dividers.

SPARE PARTS

Brackets and shelves

į.	2.40.1010 4114 01101100			
ĺ	For Bloodbank code	2101497	2101498	2101499 / 2101500
l	Brackets - 4 units	1001801	1001801	1001801
î.	Shelves	1001804	1001805	1001806
t	Fach shalf requires A hrankets			

Shaking monitor for blood collection "HEDO2"

FOR BAGS BLOOD COLLECTION PROCESSES, SWINGING PLATFORM,

FEATURES

HEDO2 is a digital blood collection monitor, created to control and monitoring blood collection processes. It has been designed to perform simple, precise and accurate operations.

It is composed of:

- -Easily removable tray.
- -Wired or wireless interface networking.
- -Automatic clamp closing at the end of the donation.
- -Automatic tare calculation.
- -Real time monitoring of the donation showed on the digital display.
- -Stand alone mode (operating with rechargeable battery). 8 hours life.
- -Tube sensor presence.

Power supply: 110-220V / 50-60 Hz.

MODEL

Part No.	Max weight	Precision	Resolution	Height / Width / Depth	Power	Weight
	Capacity g	g	g	(external) cm	W	Kg
5010000	999	± 3	1	19 32 32	60	5

It is supplied with a transport suitcase.



ACCESSORIES

Wireless data transfer Hedo-2. Part No. 5010001 Wireless connection module Hedo-2, Part No. 5010002 Barcode reader Hedo-2, Part No. 5010003 Alarm lamp pool Hedo-2. Part No. 5010004



Plasma bags defrost baths and temperature controllers for red cells packed "Unitronic-OR"

NEW

"UNITRONIC-OR-4" CAPACITY FOR 4 BAGS.

"UNITRONIC-OR-12" CAPACITY FOR 12 BAGS.

SHAKING SPEED, TIME AND TEMPERATURE DIGITAL ELECTRONIC REGULATION.

FOR REGULABLE TEMPERATURES FROM AMBIENT +5 UP TO 36 °C.

STABILITY \pm 0,05 °C. Homogeneity \pm 0,1 °C. Set temperature error \pm 1 °C. Resolution 0,1 °C.



SAFETY:

NORMA DIN 12879.0. SAFETY THERMOSTAT ADJUSTABLE WITH MANUAL RESET FOR SWITCHING THE BATH OFF IN CASE OF OVERTEMPERATURE.

Regulable reciprocal shaking from 10 to 70 oscillations per minute

Controllable time from 1 minute up to 99 hours 59 minutes. End of cycle acoustic alarm



MODELS	Part No.	Bags capacity		Height / Width / Depth (exterior) cm		Power W	Weight Kg
OR-4	3001250	4	40	54	44	1500	25
OR-12	3001251	12	40	74	44	2000	34

FEATURES

Linear motion shaking tray from 10 to 100 oscillations per minute, with a variable stroke length of 32 or 46 mm.

Digital selection and reading of temperature, speed and time.

Temperature sensor: Pt 100 thermo-resistant.

RS-232 outlet for parameter's reading and printing in a computer.


External case of AISI 304 stainless steel.

Internal tank with an AISI 304 stainless steel rack for the correct accommodation of plasma bags.

Drain tap built-in.

CONTROL PANEL

- 1. Main switch selector.
- 2. LCD screen that shows:
- Temperature, oscillations and time.
- Alarms: probe and temperature failure.
- Configuration: °C or °F, calibration and working limits.
- 3. Start / Stop button.
- 4. Parameter value decrease button.
- 5. Parameter value increase button.
- Temperature, oscillations and time selector button.
- Configuration button.



MODEL

Blood bags tube automatic sealer "Sealcom-B"

FOR TUBES FROM Ø 3 TO 6MM.



Part No.	Sealing	Height / Width / Depth	Power	Weight
	width mm Ø	(exterior) cm	sealing W	Kg
5800060	3 to 6	15 16 34	400	7

FEATURES

Metal epoxy covered external case.

AISI 304 stainless steel sealing head.

A microprocessor provides an automatic sealing by dielectric heating when inserting part of the tube in the piston.

Sealing time from 1 to 2 seconds.

CONTROL PANEL

Main switch, with light indicator.

Light indicating operation, sealing and end of sealing.



COMECTA Manual plasma extractor "PR-B"

EASY TO USE

MANUAL SYSTEM THAT ACCEPTS ALL KINDS OF PLASMA BAGS. TRANSPARENT PLATFORM FOR VISUAL CONTROL OF THE PLASMA. STAINLESS STEEL AISI 304 CABINET.

MODEL

Part No.	Height / Width / Depth (exterior) cm	Weight Kg		
5800061	25 17 30	2,6		







pages page	86 to 88 87
page	87
pages	88 to 91
page	92
page	92
page	93 and 94
page	94
page	95 to 97
	page page page page page page page



BATHS; WATER AND OIL, SAND, AND ULTRASONIC Summary table of the different models QUICK OVERVIEW



BATHS WATER / OIL	l					CA	.PACITY / PART I	ln			
lange	Model		Control	Safety thermost	at 1.5 litres	U.A.	IFAUITT / FARTT	vu.		Stability	RS-232
+30 110 °C	Aquater	m Th	nermostat	No	3000916					±2°C	No
					1 litre						
+30 110 °C	Baher	Th	nermostat	No	7000540					±3°C	No
					4 litres						
+40 80 °C	Termofi	n Th	nermostat	No	3000459	(Paraffin)				±1 °C	No
					2.5 litres	· ·					
+5 110 °C	Univeb	a Th	nermostat	Yes	3000400					±1 °C	No
	Omvob	u ''	ioiiiiootat	100						±1 0	140
+5 200 °C	Univeb	. TI		Voc	2.5 litres					±2 ° C	No
+3 200 °G	Unived	a II	nermostat	Yes	3000401					±2 ℃	No
					11 litres						
+5 110 °C	Termob	at Th	nermostat	Yes	3000587					±2 °C	No
					11 litres						
+5 200 °C	Termob	at Th	nermostat	Yes	3000588					±3 °C	No
					2 litres	5 litres	12 litres	20 litres	45 litres		
+5 110 °C	Preciste	rm Th	nermostat	Yes	6000137	6000138	6000140	6000141	6000142	±1 °C	No
+5 200 °C	Preciste	rm Th	nermostat	Yes	6000385	6000386	6000387	6000388	6000389	±2 °C	No
+5 99.9 ° C	Precisd	ig I	Digital <mark>µ</mark>	Yes	-	6001195	6001196	6001197	6001198	±0.1 °C	No
+5 200 °C	Precisd	•	Digital µ	Yes	-	6001236	6001237	6001238	6001239	±0.5 °C	No
+5 80 °C	Agibat-2	20	Digital	Yes	-	-	-	6032012	-	±0.5 °C	No
					2x2 litres	4x2 litres	6x2 litres				
+5 200 ° C	Precisb	at Th	nermostat	No	6001462	6001472	6001482	-	-	±2°C	No
Baths with stirring o	or shaking				27 litres						
+5 99.9 °C	Unitronic Recip	orocating	Digital µ	Yes	6032011	(Shaker bath)	-	-	-	±0.05 °C	Yes
0 99.9 °C	Unitronic Recipi	rocating C	Digital <mark>µ</mark>	Yes	6001173	(Shaker bath)	-	-	-	±0.05 °C	Yes
aths for viscometer	rs										
+ 5 100 °C : with micropro	VB-142 cessor.		Digital Ites ambier	Yes It temperature	-+5 °c.	-	-	3001423	-	±0.05 °C	No
and baths		0	11 . 1 . //45	II. / B I. II I.	. / ME II . / D . II	Ultrasonic cleanin	•	•			
	Model	Control	Height / Wi (usabl	dht / Depth Height e) cm	t / widnt / Depth (usable) cm		Part No.	Capacity litres		Widht / Depth ables) cm	
+400 ° C	Combi-Sand	Electroni	c 5 18	38 5	30 30		3000512	1	6	14 14	
Part No.			60007		3000709		3000683	2,6	9	23 13	
Ultrasonic cleaning Part No.	baths with heating ar Capacity	nd with digital e Height / Wic	•		er heater		3001208	4	10	30 14	
	litres	(usable	e) cm	. 5	W		3000513	6		30 14	
3000865	5	15 2			75		3000837	6,5		24 Ø	
3000866	6	15 30			75		3000617	9		50 14	
3000867	10	15 30			225		3000514	9		30 23	
3000868	22	15 50			375		3000515	20	15	50 30	
3000869	40	20 50			000		-	1			1
BATH TEMP	FRATURE		Silio	on oil Ethylene	glycol 50%						
LIQUIDS GU				-tinyitene	Water						
						Silicon oil	50 cps				
	°C	-40	-20	0 20	0 40	60 80	100	120	140 16	50 180	200



Water bath "Aquaterm"

ADJUSTABLE TEMPERATURE FROM 30 °C TO 110 °C. STABILITY ±2 °C.

Economic option

FEATURES

AISI 304 stainless steel tank with epoxy coated base.

Heater element embedded in high temperature, corrosive

resistant INCOLOY.

CONTROL PANEL

Temperature regulator thermostat switch, synchronised with an indicator lamp. MODEL

MODEL						
Part No.	Capacity litres	Ø (usable)	Height cm	Ø Height (exterior) cm	Power W	Weight Kg
3000916	1.5	20	11.5	22 20.5	640	1.3



ACCESSORIES IN AISI 304. STAINLESS STEEL

Perforated shelf. Part No. 6000998

Flat Lid with 5 concentric ring set.

Part No. 3000917



Boiling bath "Baher"

lacksquare adjustable temperature from 30 °C to 110 °C. Stability ± 3 °C.

FEATURES

Rapid boiling (6-7 minutes).

AISI 304 stainless steel tank and epoxy coated base.

Heater element embedded in INCOLOY stainless steel.

CONTROL PANEL

Temperature regulator thermostat switch, synchronised with an indicator lamp.

MUDEL	
Part No.	Capacity

Part No.	Capacity	Ø	Height	Ø Height	Power	Weight
	litres	(usable)	cm	(exterior) cm	W	Kg
7000540	1	14	11	15 19	640	1



ACCESSORIES MADE OF STAINLESS STEEL

Lid with bakelite handle. Rack 24 haemolysis tubes up to 13 mm diameter. Part No. 7000572 Part No. 7004002

Rack for 24 tubes up to 18 mm diameter. Part No. 7000541

Perforated shelf. Part No. 6001194



Water and oil baths "Univeba"

WATER BATH: ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C TO 110 °C STABILITY ±2 °C. OIL BATH: ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C TO 200 °C STABILITY ±3 °C.

SAFETY:

CONFORMS TO THE DIN STANDARD 12877 BY INCORPORATING A MANUALLY RESET OVER TEMPERATURE SAFETY CUT-OUT.

APPLICATIONS

To be used with a rotary evaporator or as a stand alone Bains Marie bath.

FEATURES

Tank interior and exterior case made of AISI 304 stainless steel. The heater element is embedded in a protective INCOLOY stainless steel shield.

CONTROL PANEL

Illuminated main power switch indicator lamp. Hydraulic temperature thermostat regulator synchronised with an indicator lamp. Safety thermostat cut out indicator lamp.



MODELS

Part No.	Maximum Temperature °C	Capacity litres	Ø (usable)	Height cm	Height / Widht / Depth (exterior) cm	Power W	Weight Kg
3000400	110	2.5	18	13	19 22 27	390	4
3000401	200	2.5	18	13	19 22 27	640	4

ACCESSORIES MADE OF STAINLESS STEEL

- 1- Heating element cover shelf. Part No. 6000033
- 2- Reduction ring cover lid of 110, 80, 60 and 37 mm Ø. With a thermometer location hole in the lid. Part No. 7004001
- **3- Tube rack** for 24 haemolysis tubes up to 13 mm \emptyset . Part No. 7004002
- 4- Tube support rack for 24 tubes up to 18 mm Ø. Part No. 7000541
- 5- Thermometer. Part No.

1000120 up to 120 °C. 1000873 up to 220 °C.



Water and oil baths "Termobat"

ADJUSTABLE TEMPERATURE FROM AMBIENT+5 °C TO 110 °C. STABILITY ± 2 °C TO 200 °C. STABILITY ± 3 °C.

SAFETY:

CONFORMS TO THE DIN STANDARD 12877 BY INCORPORATING A MANUALLY RESET OVER TEMPERATURE SAFETY CUT-OUT.

Large capacity

APPLICATIONS

To be used with a rotary evaporator for flasks up to 10 litres or as a stand alone bains marie.

FEATURES

Tank interior made of AISI 304 stainless steel, with an epoxy coated base.

The heater element is embedded in a protective INCO-LOY stainless steel shield.

CONTROL PANEL

Main power switch, with luminous indicator lamp.

A hydraulic temperature thermostat regulator is synchronised with an indicator lamp.

MODELO

MODELO						
Part No.	Maximum Temperature °C	Capacity litres	Ø Height (usable) cm	Ø Height (exterior) cm	Power W	Weight Kg
3000587	110	11	29 17	29 25	2100	5.5
3000588	200	11	29 17	29 25	2100	5.5

ACCESSORIES MADE OF AISI 304 STAINLESS STEEL

Reduction ring cover lid of 210, 180, 145, 115 and 80 mm Ø. With a location hole for a thermometer.





Heater element cover shelf.

Part No. 6000055



Thermometer.

Part No.

1000120 up to 120 °C. 1000873 up to 220 °C.



Sand baths "Combiplac-Sand"

ADJUSTABLE PLATE TEMPERATURE UP TO 400 °C.



FEATURES

Electronic energy heater controller. The steel hotplate surface has an evenly distributed heater element providing homogeneous heat distribution.

The non warping AISI 304 stainless steel sand bath can be easily removed for cleaning. The Hotplate can be used independently without the sand bath. The hotplate casing is made of AISI 304 stainless steel.



Mains power switch with luminous "On" indicator lamp.

Electronic heater power control between 0 to 100 %.

Heater operation indicator lamp.

The back of the unit is fitted for an electronic thermometer "Sensotherm" or digital electronic controller "Electemp-TFT" described under regulators and controllers (see page 293).



MODELS

Part No.	Height / Widht/ Depth (usable) cm	Height / Widht / Depth (exterior) cm	Power W	Weight Kg
6000708	5 18 38	18 22 41	2300	13
6000709	5 30 30	21 32 35	2800	15



Circulation pump "Ciragit"

FEATURES

Metal construction, with all components in contact with liquid made of AISI 304 stainless steel.

The unit can be adjusted to fit most tanks or baths, can be submerged and held in position by an adjustable height support clamp.

Pump circulation with adjustable flow rates allowing external circulation of liquids.

ON/ OFF illuminated power switch with light when on. Interchangeable flow output connector:

- position A: Internal tank circulation.
- position B: External circulation.

SPECIFICATIONS

OUTPUT water:

20 °C , 12 I/ min.- 150 mbar **OUTPUT** silicon oil:

50 °C 4 I/min.- 130 mbar

Dimensions: 9.5 x 12.5 x 23 cm

Power: 30 W

Weight: 2.5 Kg.

Part No. 3000402





Water and oil baths "Precisterm"

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT +5 °C UP TO 110 °C, STABILITY ± 1 °C. Adjustable temperature range from ambient +5 °C up to 200 °C, stability ± 2 °C.

Extensive range of sizes and accessories

Stainless steel interior and external case



FEATURES

Capacity; from 2 litres up to 45 litres.

Double skinned tank of AISI 304 stainless steel with an exterior case of AISI 304 stainless steel. Heating elements embedded in INCOLOY heat and corrosion resistant stainless steel.

CONTROL PANEL

Illuminated mains power switch.

Hydraulic thermostat temperature control, synchronised with the "ON" indicator lamp. Safety thermostat cut out indicator lamp.

MODELS

Part No.	Maximum Temperature °C	Capacity litres	•	/ Widht usable) (Height / (ex	Widht (terior) (Power W	Weight Kg
6000137	110	2	15	15	13	20	20	27	390	3.5
6000385	200	2	15	15	13	20	20	27	640	3.5
6000138	110	5	15	30	14	20	37	27	600	6
6000386	200	5	15	30	14	20	37	27	1110	6
6000140	110	12	15	30	30	20	40	42	1100	8
6000387	200	12	15	30	30	20	40	42	1400	8
6000141	110	20	15	48	30	20	58	42	1500	10
6000388	200	20	15	48	30	20	58	42	2100	10
6000142*	110	45	15	62	50	22	69	62	2400	18
6000389*	200	45	15	62	50	22	69	62	3800	18

SAFETY:

DIN 12877 STANDARD.
THE UNIT WILL CUT OUT IF THERE IS
INSUFFICIENT LIQUID OR IF THE BATH
TEMPERATURE EXCEEDS THE SET
TEMPERATURE.
A MANUAL RESET ADJUSTABLE

A MANUAL RESET ADJUSTABLE
TEMPERATURE SAFETY THERMOSTAT IS
PLACED AT THE BACK OF THE UNIT.



ACCESSORIES. See page 91.

^{*} With built in drain tap.



Water baths "Precisdig" with or without anti-evaporation folding lid

DIGITAL CONTROL AND DISPLAY OF TEMPERATURE AND TIME.

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT+5 °C TO 99.9 °C, STABILITY ±0.1 °C, HOMOGENEITY ±1 °C. ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT+5 °C TO 200 °C, STABILITY±0.5 °C, HOMOGENEITY ±2.5 °C.

SAFETY: DIN 12877 STANDARD. THE UNIT WILL CUT OUT IF THERE IS INSUFFICIENT LIQUID OR IF THE BATH TEMPERATURE EXCEEDS THE SET TEMPERATURE. A MANUAL RESET ADJUSTABLE TEMPERATURE SAFETY THERMOSTAT IS PLACED AT THE BACK OF THE UNIT.

Programmable temperature and time. Interior and exterior made from Stainless steel





MODELS								- 8		
Part No.	Temperature maximum °C	Capacity litres	•	/ Widh usable)	nt / Depth) cm		/ Widh xterior	it / Depth) cm	Power W	Weight Kg
6001195	99,9	5	15	30	14	20	37	27	600	6
6001236	200	5	15	30	14	20	37	27	1100	6
6001196	99,9	12	15	30	30	20	40	42	1100	8
6001237	200	12	15	30	30	20	40	42	1400	8
6001197	99,9	20	15	48	30	20	58	42	1500	10
6001238	200	20	15	48	30	20	58	42	2100	10
6001198	99,9	45	15	62	50	20	69	62	2400	18
6001239	200	45	15	62	50	20	69	62	3800	18

All baths with heating shelf and drainage. Faucet included.

MODELS WITH MAKROLON® FOLDING LID

Part No.	Temperature maximum °C	Capacity litres	Height /			•		t / Depth id) cm	Power W	Weight Kg
6001186	99,9	12	15	30	30	40	40	42	1100	8
6001187	99,9	20	15	48	30	40	58	42	1500	10
6001188	99,9	45	15	62	50	40	69	62	2400	18


All baths with heating shelf and drainage. Faucet included.

FEATURES

Double skinned tank of AISI 304 stainless steel with an external case of AISI 304 stainless steel.

Heating elements embedded in INCOLOY heat and corrosion resistant stainless steel.

Comes complete with drain tap.



CONTROL PANEL

- 1. Illuminated mains switch.
- 2. Temperature mode indicator.
- 3. Time mode indicator.
- 4. Temperature and time indicator.
- **5.** Start time indicator.
- 6. Waiting / delay time indicator.
- **7.** Temperature set button.
- 8. Time set button.
- 9. Increase, push button.
- 10. Decrease, push button.
- 11. Stop and start time.
- 12. Read Temperature.
- 13. Set start up time: from 1 min. to 9 hours 59 min, or up to 99.9 hours once the set temperature has been reached.
- 14. Set waiting / delay time: from 1 to 24 hours.
- 15. Safety thermostat indicator lamp.

MAKROLON® FOLDING LID PROPERTIES

Due to its transparency, you can visualize the bath's inside. It improves temperature stability, and restricts evaporation and pollution

ACCESSORIES. See page 91.



Water bath with inner stirring "Agibat-20" with digital control and display of temperature and time

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT +5°C UP TO 80°C. STABILITY ±0.1°C. SET TEMPERATURE ACCURACY ±1 °C, WORKING TEMPERATURE PRECISION 0.1 °C.

P.Selecta original model. More space. More accuracy. Economical.







FEATURES

Double skinned tank, inner case made of AISI 304 stainless steel. External case of AISI 304 steel. Submerged heating elements in INCOLOY special alloy. It comes complete with drain tap.

STIRRING SYSTEM

These baths have an aspiration and outflow internal circulation pump which allow a constant flow stirring.

Flow: 8 litres per minute.

CONTROL PANEL

See page 8

Three positions switch: only heater, heater + stirring and start/stop.

MODEL WITH MAKROLON® FOLDING LID

Due to its transparency, you can visualize the bath's inside. It improves temperature stability, and restricts evaporation and pollution.

MODELS

Part No.	Lid Makrolon®	Maximum Temp °C	Litres Capacity	Height/Width/Depth (usable) cm	Height/Width/Depth (exterior) cm	Power W	Weight Kg
6032012	Yes	80	20	15 48 30	40 58 42	2010	11
6032013	No	80	20	15 48 30	40 58 42	2010	12

All baths with heating shelf and drainage. Faucet included.

ACCESSORIES. See page 91.

SAFETY:

DIN 12877 STANDARD. THE UNIT WILL CUT OUT IF THERE IS INSUFFICIENT LIQUID OR IF THE BATH TEMPERATURE EXCEEDS THE SET TEMPERATURE. A MANUAL RESET ADJUSTABLE TEMPERATURE SAFETY THERMOSTAT IS PLACED AT THE BACK OF THE UNIT.





ACCESSORIES FOR"PRECISTERM", "PRECIDSDIG" AND "AGIBAT-20"

Suitable for all "PRECISTERM", "PECISDIG" and "AGIBAT-20" models.



Removable constant liquid level.
Part No. 6001400
Note: Not compatible with gale and Makrolon® lid.



Fixed, constant liquid level.
Part No. 6000123

Únicamente para baños modelo "PRECISTERM"



Thermometer support with locking screw (up to 40 mm).
Part No. 6000896



- 1. Gable lid with thermometer location position hole.
- 2. Flat lid with concentric reduction rings of 105, 80, 60 and 37 mm Ø. With thermometer location position hole.

Note: The perforated lid avoids the lock of the Makrolon $\ensuremath{\mathfrak{B}}$ folding lid.

- 3. Support frame for tube racks.
- 4. Tube Racks various tube diameters available.




Capacity litres	Model code AGIBAT-20	Model code PRECISTERM	Model code PRECISDIG	1 Gable lid	2 Lid with concentric reduction rings	3 Support frame
2		6000137 - 600038	ō	-	6003231 1 place	-
5		6000138 - 600038	6 - 6001195*- 6001236*	6000710	6003232 2 places	1001199
12		6000140 - 600038	7 - 6001196*- 6001237*	6000711	6003234 4 places	1001200
20	6032012*6032013*	6000141 - 600038	8 - 6001197*- 6001238*	6000712*	6003236 6 places	1001201
45		6000142*- 600038	9*- 6001198*- 6001239*	6000714	6001603 12 places	1001205

^{*} Comes complete with drain tap. * For all baths except those carrying Makrolon® lid.

Capacity litres	for PRECISDIG baths with Makrolon® lid	2 Lid with concentric reduction rings	3 Support frame
12	6001186	6003234 4 places	1001200
20	6001187	6003236 6 places	1001201
45	6001188	6001603 12 places	1001205

	4 Tube racks (requires support frame)											
Part No.	Tube Ø mm	Number of tubes per rack	Rack Dimensions Height / Widht / Depth (usable) cm	Number of racks by bath capacity litres 5 12 20 45								
1001202	13	36	8 8.2 23.6	1	3	5	10					
1001203	16	24	8 8.2 23.6	1	3	5	10					
1001204	20	24	8 8.2 23.6	1	3	5	10					





Stabilising rings "circum" prevents flasks from floating when placed in a water bath. Suitable for all types of flasks,erlenmeyers, etc. See page. 177.



Histology paraffin section mounting bath "Termofin"

ADJUSTABLE TEMPERATURE FROM 40 °C TO 80 °C. STABILITY ±1 °C.

APPLICATIONS

Suitable for the fusion and preservation of paraffin wax in a liquid state that enables the user to see the floated sections.

FEATURES

Aluminium inner tank and lid with a P.T.F.E. matt black coating. Thermometer location hole in the lid to read temperature.

CONTROL PANEL

Mains power ON/OFF switch. Hydraulic temperature thermostat regulator synchronised with an indicator lamp.

MODEL

Part No.	Capacity	Ø	Height	Ø Height	Power	Weight
	litres	(usable)	cm	(exterior) cm	W	Kg
3000459	4	25	10	34 14	430	2.3

Comes complete with lid and thermometer.



ACCESSORY

Thermometer support (can only be used without the lid).

Part No. 6004591

NOTE: Wax dispenser (see page 217).



Bank of oil and water baths "Precisbat"

ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 200 °C, STABILITY \pm 2 °C. INDEPENDENT TEMPERATURE CONTROL FOR EACH POSITION.

FEATURES

Inner vessel made of AISI 304 stainless steel with the external case made of AISI 304 stainless steel.

The heater element is embedded in stainless steel INCO-LOY alloy, that is corrosion and heat resistant.

Each position has a lid with concentric reduction rings of; 105, 80, 60 and 40 mm in Ø.

An adjustable support frame is mounted at the back of the unit to facilitate the mounting and support of flasks, condensers and other apparatus, etc.

Capacity: 2 litres per position.

CONTROL PANEL

Independent control for each position.

Mains ON/OFF with luminous "ON" indicator lamp for each position.

Hydraulic thermostat temperature control synchronised with the "ON" indicator lamp.



MODELS

Part No.	Temperature maximum °C	Number of places	Height / Widht / Depth (usable each position) cm	Height / Widht / Depth (exterior) cm	Power W	Weight Kg
6001462	200	2	12 14 12	19 35 29	1300	5.5
6001472	200	4	12 14 12	19 67 29	2560	10
6001482	200	6	12 14 12	19 100 29	3800	15

ACCESSORIES



Perforated shelf (One per tank) made of AISI 304 stainless steel. Part No. 6000011



Clamp for supporting flasks and condensers. Part No. 7000403



Thermometer up to 220 °C. Part No.1000873



Bosshead for supporting clamps. Part No. 1000069



Precise reciprocal shaking baths controled by microprocessor "Unitronic Reciprocating Shaking Bath" and with refrigeration "Unitronic Reciprocating Shaking Bath C"



UNITRONIC RECIPROCATING SHAKING BATH LINEAR: TEMPERATURE RANGE FROM AMBIENT +5 °C to 99.9 °C. Unitronic reciprocating shaking bath c linear: Temperature range from 0 °C to 99.9 °C. Stability \pm 0.05 °C homogeneity \pm 0.1 °C. Set temp. Error \pm 1 °C. Resolution 0.1 °C.

SAFETY:

DIN 12877 STANDARD. THE UNIT WILL CUT OUT IF THERE IS INSUFFICIENT LIQUID OR IF THE BATH TEMPERATURE EXCEEDS THE SET TEMPERATURE. A MANUAL RESET ADJUSTABLE TEMPERATURE SAFETY THERMOSTAT IS PLACED AT THE BACK OF THE UNIT.

Reciprocal linear agitation, controllable from 10 to 100 oscillations per minute

Controllable time from 1 minute to 99 hours 59 minutes

MICROPROCESSOR CONTROL WITH DIGITAL DISPLAY AND PUSH BUTTON TEMPERATURE, SPEED AND TIME SELECTION. RS-232 INTERFACE FOR PRINTER OR COMPUTER OUTPUT.



APPLICATIONS

Cell culture and biological samples, yeasts and general applications that require stable and reproducible temperatures.

COMMON FEATURES

Linear motion sliding tray from 10 to 100 oscillations per minute, with a variable stroke length of 32 or 46 mm.

Digital display and selection of temperature, oscillations and time.

Temperature sensor: thermo-resistant Pt100. Rack frame support made of AISI 304 stainless steel incorporated in the tank that can accommodate a selection of flasks and racks.

Internal tank made of AISI 304 stainless steel with an external case of AISI 304 stainless steel.

Built in drain tap.

Ventilated hermetically sealed compressor for refrigeration in the "Unitronic - C".

CONTROL PANEL

- 1. Mains power switch.
- 2. LCD screen that shows:
 - Temperature, oscillations per minute and time.
 - Alarms: set temperature and probe failure.
 - Configuration: °C or °F, temperature calibration and set working temperature limits.
- 3. Push button Start/ Stop.
- 4. Push button to decrease value.
- 5. Push button to increase value.
- **6.** Push button for parameter selection: temperature, strokes per min. and time.
- 7. Push button for configuration.



"Unitronic Reciprocating Shaking Bath"



"Unitronic Reciprocating Shaking Bath C" Refrigerated Bath



MODELS Part No.	Capacity litres	Height / Width / Depth (bath) cm	Height / Width / Depth (support carriage) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
UNITRONIC RECIPROCATING 6032011	27	18 48 29	17 35 24	30 74 40.5	1500	25
UNITRONIC RECIPROCATING-C 6001173	27	18 48 29	17 35 24	30 74 63.5	1750	43

ACCESSORIES FOR UNITRONIC RECIPROCATING SERIES



Support Platform for Erlenmeyer flasks.

Made of AISI 304 stainless steel with location holes for various configurations of Erlenmeyer flask clip supports. Sizes: 50, 100, 250, 500 and 1000 ml.

Part No.	Width / (platfor	- 1	50 ml	No. of flash 100 ml	ks and capa 250 ml		atform 1000 ml	Weight Kg
3001004	23.5	34.5	20	12	6	6	3	1



Racks: made of AISI 304 stainless steel.

Capacity: 2 racks.

Part No.	Tube Ø mm	No tubes per rack	Height / \ (ra	Width ack) c	
1001233	13	70	8.5	11	34
1001234	16	48	8.5	11	34
1001235	20	48	8.5	11	34



Adapter clips for Erlenmeyer flasks.

Made of tempered stainless steel.

Part No.

1001003	Adapter for	50 ml.
1000978	Adapter for	100 ml
1000979	Adapter for	250 ml.
1000980	Adapter for	500 ml.
1000981	Adapter for	1000 ml.



Gable Lid: made of AISI 304 stainless steel.

Part No. 6000713



Precision Viscometer Bath "VB-1423"

DIGITAL TEMPERATURE CONTROL FROM AMBIENT +5 °C TO 100 °C. STABILITY ± 0.05 °C. Homogeneity ± 0.05 °C. Set error value ± 0.09 °C. Resolution 0.1 °C.

SAFETY:

SAFETY THERMOSTAT CONFORMS TO THE DIN 12876. MANUAL RESET.

FEATURES

Temperature sensor; Pt100 thermo-resistor, stainless steel AISI 304 lid with three viscometer locations ports, three independent lids and an additional location port for the control thermometer. The main body of the bath is made of a 20 litre borosilicate glass tank. A white plate is located at the back to help optimize and read the viscometers.

CONTROL PANEL

- 1. Main power switch with luminous ON.
- 2. Temperature regulator:
 - 3. Real time temperature display.
 - **4**. Push button increase value.
 - 5. Push button decrease value.
 - **6.** Push button configure operation.
- 7. Safety thermostat safety lamp.

Precise constant temperature Kinematic Viscometer Bath

Manufactured for the calibration of viscometers, conforms to UNE 400313, ISO3105, ASTM 445 & 2515 specifications





MODEL

MODEL						
Part No.	Temperature control range °C	Capacity litres	Height / Ø (tank) cm	Height / Ø (total) cm	Power W	Weight Kg
3001423	amb.+5 up to 100	20	32 30	47 30	1000	8

ACCESSORIES

Universal viscometer support. made from PTFE with stainless steel AISI 304 support. Suitable for the following viscometers :



- Cannon-Fenske for transparent liquids.
- Cannon-Fenske for opaque liquids.
- Ubbelohde.
- Ostwald.
- BS U Tube.
- -Cannon-Manning semi-micro.
- Ubbelohde type BS/IP/SL, BS/IP/SL(S) & type BS/IP/MSL.
- DIN Ubbelohde.

Part No. 1001453

Calibration Chronometers (see page 295).

Thermometers for viscometer baths.

Part No.

1001454	Thermometer ASTM	120C at	38.6 to	41.4 °C divisions of	0.05 °C.
1001455	Thermometer ASTM	121C at	98.6 to	101.4 °C divisions of	0.05 °C.
1001456	Thermometer ASTM	91C at	20.0 to	50.0 °C divisions of	0.1 °C.
1001457	Thermometer ASTM	92C at	40.0 to	70.0 °C divisions of	0.1 °C.
1001458	Thermometer ASTM	93C at	60.0 to	90.0 °C divisions of	0.1 °C.
1001459	Thermometer ASTM	94C at	80.0 to	110.0 °C divisions of	0.1 °C.



Ultrasonic Cleaning Baths "Ultrasons",

"Ultrasons UB-1488", "Ultrasons-H", "Ultrasons-HD" and "Ultrasons-P"

FUNDAMENTAL THEORY

The principal of ultrasonic cleaning consists of the use of high frequency sound waves (40 kHz), produced by a generator through a transducer, which propagates them mechanically inside the tank, this produces a cavitation effect which leads to the formation of millions of low pressure microscopical bubbles which carry out molecular cleaning, eliminating impurities, polluting agents and dirt from the parts or material which must be cleaned.

Complex parts can be cleaned without disassembly since the cavitation penetrates wherever the cleaning solution is in contact with the surface.

The generator of these baths is completely transistorised and incorporated, working through the "Sweeping frequency" system, which ensures uniform cavitation at all points of the resonance tank.

COMMON FEATURES

Double bodied tank, made entirely of stainless steel.

Emptying drain, anti-parasite filter, adjustable timer and heating, model dependent.

APPLICATIONS

Laboratories in general: acceleration of chemical reactions, degassing of liquids, cell disruption, cleaning of sieves, pipettes, micro-pipettes, cuvettes, trays, viscometers, decomposition of radioactive substances, etc.

Dentistry: Cleaning of prosthesis, instruments, etc.

Optics: Cleaning of frames, lenses, and contact lenses, etc.

SUMMARY TABLE

Medicine: Cleaning of instruments in general, forceps, probes, scalpels etc. **Jewellery and Watch making:** Cleaning of all kind of watches and jewels. **Industry:** Cleaning of all electronic components, printed circuit boards, semiconductors, petrol injector sieves and filters, etc.





Ultrasonic bath, Automatic Lift and Dry "Ultrasons-UB-1488"

Selecta original model. Quick dry, without having to control the instrument



MODEL

Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Generator power	Power	Weight
	litres	(usable) cm	(exterior) cm	W	total W	Kg
3001488	6	15 30 14	30 40 34	150	2150	12

Comes complete with basket.

FEATURES

Bath with automatic basket elevator that operates at the end of the time cleaning cycle.

Hot air drying by two-turbine air fans that switch on at the start of the elevation and finalise the complete operation. Waste drain at the front of the instrument. Removable support basket.



CONTROL PANEL

- 1. Cleaning time control: 0-15 minutes.
- 2. Push button drying time: 0-15 minutes.
- 3. Push button increase value.
- 4. Push button reduce value.
- **5**. Push button stop / start.
- **6**. Digital timer display.



Ultrasonic cleaning baths "Ultrasons" without heating





CONTROL PANEL

Timer switch for 0 -15 minutes or continuous operation. Indicator lamp when the unit is in operation.

MODELS

Part No.		Capacity litres		/ Widtl sable)	n / Depth cm	Ŭ	nt / Widtl (exterior)	n / Depth cm	Generator power W	Timer 0-15'	Waste drain	Weight Kg
3000512	1	1	6	14	14	20	20	18	50	Yes	No	2,5
3000683	2	2,6	9	23	13	22	29	19	100	Yes	No	4
3001208	8	4	10	30	14	24	33	18	100	Yes	Yes	7
3000513	3	6	15	30	14	30	36	20	150	Yes	Yes	9,5
3000837	5	6,5	14	24	Ø	30	27,5	Ø	150	Yes	No	8
3000617	4	9	15	50	14	30	56	20	200	Yes	Yes	10,7
3000514	6	9	15	30	23	30	36	30	200	Yes	Yes	13
3000515	7	20	15	50	30	30	56	36	400	Yes	Yes	22

COMPLEMENT



Detergent "BIO-SEL"

Bio-degradable, phosphate free. Excellent cleaning power just by mixing only 20ml per litre of water. BIO-SEL, in a 2% solution has a 9.7 pH level. Low toxicity. This detergent disperses, emulsifies, segregates, dissolves, suspends and breaks down dirtiness.

1000025 Capacity 1 litre. 1000026 Capacity 4 litres.

1000852 Box of 12 x 1 litre bottles.

1000853 Box of 4 x 4 litres bottles.





Lids are not compatible with baskets.



Lid with concentric reduction ring



ACCESSORIES

Baskets and Lids.

Made of AISI 304 stainless steel.

Cleaner model Part No.	BASKET Part No.	Heigh	BASKE t / Widt	:1S LI th / Depth (cm)	D WITH CONCEN Part No.	IRIC REDUCTION RING Number of positions	Part No.
3000512	6005123	5	11	11	6005122	1	6005121
3000683	6006833	7	20	12	6006832	1	6006831
3001208	6001209	7	26	13	-	-	6005131
3000513	6005133	12	26	13	6005132	2	6005131
3000837	6008373	13	23	Ø	-	-	6008371
3000617	6005144	12	46	12	6005134	3	6005124
3000514	6005143	13	29	23	6005142	2	6005141
3000515	6005153	13	49	29	6003236	6	6005151



ltrasonic Bath "Ultrasons-Digit"

5 CYCLES ELECTRONIC DIGITAL TIME CONTROL **CAPACITY 1.4 LITRES**

FEATURES

Internal tank made from stainless steel, with minimum and maximum level indicator.

External case, basket and lid made from PVC. Maximum cleaning power, frequency 42 KHz

CONTROL PANEL

Pushbutton start/ stop

Timer 5 programs selection in seconds: 90" - 180" - 280" - 380" and 480".

Wide application range from jewellery, watches, optical components. Dentistry and general laboratory use.

MODEL Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Power	Weight
	litres	(usable) cm	(exterior) cm	W	Kg
5900650	1.4	6 13 16	15 17 22	70	1.75

Comes complete with basket and see through transparent PVC lid.





Heating cleaning bath "Ultrasons-HD"

DIGITAL ELECTRONIC REGULATION OF TEMPERATURE AND TIME FOR TEMPERATURES FROM AMBIENT +5°C TILL 80°C

Optimum cleaning efficiency



It comes complete with a universal basket, lid with handle and drain device.

Bath Part. No. 3000865 (5 Litres capacity) is also supplied with basket for small items.

FEATURES

Double body container, completely made of stainless steel and with waste drain device.

CONTROL PANEL

Illuminated main switch.

Start-stop push button.

Increase/decrease time push button in minutes (from 0' to 60').

Increase/decrease temperature push button in °C.

MODEL

Part No.	Capacity litres	•	'Width able)			t/Wid terior)	th/Dept cm	h Power Generator W	Power Heater W	Weight Kg
300086	5 5	15	24	14	27	31	18	120	75	3,8
300086	6 6	15	30	15	27	35	18	180	75	5,8
300086	7 10	15	30	24	35	38	28	240	225	8,2
300086	8 22	15	50	30	27	58	35	600	375	15
300086	9 40	20	50	40	40	65	50	900	1000	25

It comes complete with a universal basket, small basket for small items, lid with handle and drain device.



Heated Ultrasonic baths "Ultrasons-H"

AJUSTABLE TEMPERATURE FROM AMBIENT +5 °C HASTA 75 °C.

APPLICATIONS

For long pieces cleaning (especially for pipettes).



Note: For cleaning long objects whose length may exceed that of the tank, such as tubes, pipettes, etc, it has been proved that they are submitted to the same cavitation effect by introducing them in a test tube filled with liquid and placed vertically in the cleaner.

FFATURES

Evenly distributed heating element, adhered to the tank for maximum heat transfer. CONTROL PANEL

Adjustable thermostat heater.

Heater "on" indicator lamp.

Timer switch from 0-15 minutes or continuous operation.

In operation on indicator lamp.

MODEL

Part No.	Capacity litres		Width able)			t/Wid terior)		oth Power Generator W	Power Heater W	Weight Kg
3000839	9	15	50	14	30	56	20	200	800	11

ACCESSORIES

Basket, Height 12/Width 46/Depth12 cm. Part No 6005144

3 Lid with concentric reduction ring. Part No ${\bf 6005134}$

Lid with handle. Part No 6005124

See page: 91





Ultrasonic pipette cleaner "Ultrasons-P"

FEATURES

Designed for automatically cleaning and rinsing pipettes, mounted inside a basket, pipettes are ultrasonically cleaned and rinsed.

Minimum rinsing time by using a water siphon action through the inside of the pipettes.

Metallic areas in contact with liquid are made of AISI 316 stainless steel.

The outer casing is made of AISI 304 stainless steel.

Cylindrical vessel in transparent acrylic, with lid to avoid splashing.

Connector(1") water inlet and siphon-drain tube.

Basket with handle in polypropylene for the easy handling and carrying of the pipettes.

CONTROL PANEL

Timer from 0-15 minutes or continuous operation.

Indicator lamp showing when the bath is in operation.

MODE

Part No.	Capacity	Ø / Height	Height / Widht/ Depth	Power	Weight
	litres	(usable) cm	(exterior) cm	W	Kg
3000910	10	15 64	84 34 40	250	15



Immersion thermostats. Analogue control	page	100
Immersion thermostats. Digital & Touch screen control	pages	101 to 106
Circulating ultrathermostat	page	103
Circulating criothermostat	page	103
Bath circulators with thermostatic equipment	pages	104 to 105
Digital & Touch screen control		
Immersion thermostat for baths	pages	105 to 106
Refrigerated units for baths	page	107
Water recirculators	page	107
Thermostat dry-blocks	page	108 to 109
Thermo shakers	page	109 to 110



THERMOSTAT BATHS AND IMMERSION THERMOSTATS

Summary table of the different models EQUIPMENT FOR STABLE AND PRECISE TEMPERATURE CONTROL.


QUICK OVERVIEW

- /		
- ~		-
	RS232	7

EQUIFINIENI FU	II O IADLL A	AND I HEOLO	E I E IVIII E I II	I OIL OUI	TITIOL.					
Range		Model	Pari	t No.	Control		Volume máx. litres	Stability °C	Set error °C	RS-232
+5 100 °C	Ter	motronic II	3000	389	Analogue		20	±0,2	-	No
+5 150 °C	Digi	term-S-150	3000)543	Digital µ		20	to 150 °C ±0,05	±1	Yes
+5 200 °C	Digito	erm-TFT-200	3000)544	TFT Touch scree	en µ	20	to 200 °C ±0,1	±2	Yes
								to 100 °C ±0,05		
0 100 °C	D	igit-Cool	300	1373	Digital µ		20	to 100 °C ±0,05	±1	Yes
. F 4F0 9 0	Dinita	C 450 0	0 000	DE 40	Distrikat u		20	to 0 °C ±0,05		Yes
+5 150 °C	•	rm-S-150-2)548	Digital µ			±0,05	-	
+5 150 °C	•	erm-S-150-2)549	Digital µ		27	±0,05	-	Yes
+5 200 °C	•	m-TFT-200 2			TFT Touch scree	•	20	±0,1	±2	Yes
+5 200 °C	Digiter	rm-TFT-200	27 3000)551	TFT Touch scree	en µ	27	±0,1	±2	Yes
+5 200 °C	III	term-TFT-20	0 200	DE 4 E	TET Touch core	en	0	to 100 °C ±0,05	±1	Yes
+3 200 6	Ultrat	leriii-1 - 1-20	U 3001)545	TFT Touch scre	enµ	8	to 200 °C ±0,1	±2	Yes
								to 100 °C ±0,05		
-10 100 °C	Frigi	term-TFT-10	3000)546	TFT Touch scree	enµ	8	to 0 °C ±0,05	±1	Yes
								to -10 °C ±0,1		
								to 100 °C ±0,05		
-30 100 °C	Frigit	term-TFT-30	3000)547	TFT Touch scree	en µ	8	to 0 °C ±0,05	±1	Yes
								to -10 °C ±0,2		
+5 up to 90	°C In	terheater	600	1491	Digital		20 l/min.	-	±1	No
Refrigerated coils										
-20 20 °C		Frigedor	300	0778	No controlle	r	20	-	-	No
-20 20 °C	Fri	gedor-Reg	300	1214	Digital		20	±1,5	±1,5	No
THERMOSTATS DRY- Range	BLOCKS Model	Control	Blocks-1	Blocks-3	Blocks-4	Stability				
30 200 °C	Tembloc	Digital µ	7462200	-	-	±0,75	°C	With multiple adapters.		
30 200 °C N	/lultiplaces	Digital µ	-	7471200	-	±0,75	°C	With multiple adapters.		
FIXED 37 °C	Clinic-bloc	-	7001569	-	-	±0,5 °	C	Capacity 20 tubes of 75	x 13mm. VAC.	
FIXED 37 °C	Bio-bath	-	7001561	-	-	±0,5 °	C	Capacity 15 cuvettes of	10 x 10mm.	
+5 60 °C	TRM-4	Digital µ	-	-	5109200	±0,5 °		For 4 microtiter plates.		
+5 100 °C	TR100-6	Digital µ	5109100	_	-	±0,5 °		WITH SHAKING. With m	ultiple adapters.	
: with micronroc									-p	

 $\ensuremath{\mu}\xspace$ with microprocessor.

^{+5:} ambient +5 °C.





SAFETY:

STANDARD DIN 12876. Safety over temperature thermostat with signal lamp alarm

FEATURES

Electronic temperature control.

Pt100 temperature sensor.

Maximum volume to maintain maximum temperate: 20 litres. Heater elements made of INCOLoY stainless ste-

Fitted with a stirring circulating pump.

All parts in contact with liquids made of stainless steel AISI 304

Fixing nut on the side.

Minimum required tank depth: 14 cm

CONTROL PANEL

Illuminated mains switch. Temperature control selector. Heater "on" indicator lamp.



MODEL

Part No.	Temp. Range	Stability	Height / Width / Depth	Pump	Power	Weight
	°C	°C	(exterior) cm	Pressure / Flowrate H ₂ O	W	Kg
3000389	Amb.+5 to 100	±0.2	32* 11 16	150 mbar 5 l/min.	1000	2.5

^{*} Height measured with the thermometer included.

It is supplied with a thermometer.

ACCESSORIES

Methacrylate tank "Clinic-Term"

Part No.	Capacity	Temp. max.	Height / Width / Depth	Capacity	Weight
	litres	°C	(usable) cm	Rack	Kg
1000544	8	60	14 16 39	4	2

Racks for the "Clinic-Term" made of AISI 304 stainless steel.

Racks for 50 micro-tubes of 11mm Ø. Part No. 1000546

Racks for 50 haemolysis tubes and coagulation tests with a diameter of 13 mm. Part No. 1000545

Racks for 14 tubes with a diameter of 18mm. Part No. 1000547

Racks for 8 cuvettes for spectrophotometers, 10 mm square, standard for monotest flask of different sizes. Part No. 1000548

Note: The methacrylate tank should not be used over 60 °C.

Methacrylate transparent tank.

Designed to be used with the immersion thermostat "Termotronic II".

Part No.	Capacity	Temp. max.	Height / Width / Depth	Capacity	Weight
	litres	°C	(usable) cm	Rack	Kg
1000397	7	60	10 20 38	4	2

Rack for 24 haemolysis tubes up to 13 mm Ø, made from AISI304 S.S. Part No. 1002532 Rack for 14 test tubes up to 16 mm Ø, made from AISI 304 S.S. Part No. 1002531

Metal baths. Manufactured with a double external AISI 304 stainless steel skin and an AISI 304 stainless steel pressed bath interior.

p								
Part No.	Capacity litres	Height / W (usab	/idth / ble) cı		Height / (ex	Width terior)	•	Weight Kg
6000390	9	15 2	29	22	20	34	28	3.5
6000391	12	15 3	30	31	20	37	35	4.3
6000392*	20	15 4	48	30	20	55	35	6.6
6000393*	27	20 4	48	30	25	55	35	7.6
6001093*	45	15 6	62	50	23	70	56	11







* With drain tap and lifting handles.



IMMERSION THERMOSTATS WITH DIGITAL AND TOUCH SCREEN CONTROL

The most precise solution for maintaining constant liquid temperatures

SAFETY:

STANDARD DIN 12876.

ADJUSTABLE OVER TEMPERATURE SAFETY THERMOSTAT WITH MANUAL RESET AND LOW LIQUID LEVEL FLOAT SENSOR.

COMMON FEATURES

Temperature probe Pt 100.

Heating element made of stainless steel INCOLOY resistant to high temperatures and corrosive environments.

Circulation Pump for internal or external temperature control for:

(Polarimeters, refractometers and viscometers, etc).

Cooling coil to maintain ambient temperatures using tap water.

External temperature control connection using a Pt 100 probe (see page 106).

External metal case coated in a corrosive resistant epoxy.

All parts in contact with liquid are made of AISI 304 stainless steel.

PARTS THAT ARE IN CONTACT WITH LIQUID

- -Safety over temperature thermostat.
- -Liquid level float.
- -Cooling coil.
- -Heating element made of stainless steel INCOLOY that can withstand high temperatures and is corrosion resistant.
- -Refrigeration coil connections.
- -Pt 100 temperature probe.
- -Circulation pump.
- -Pump outlet for external circulation.
- -Circulating pump outlet.





Outlet nut of the fluid flow fixed in A or B interchangeable:

- A. Internal stirring within the tank.
- B. External circulation.

H₂O FLOW OUTPUT

Flow I/m

Cuvette inner circulation level Up to 7,05

External circulation level raised to a meter. Up to 5,16

Tests made with silicone tubes of 8 mm internal Ø.

TWO DIFERENT CONTROL SYSTEMS:



Digital display and



TFT touch screen.

DIGITERM-S-150

Specific functions of models with digital screen:

- Real time temperature reading.
- Display resolution 0.1 °.
- Connection device for external temperature Pt100 probe.
- Permanent view of set temperature.

- Low liquid level indication.
- Heating indicator.
- Monitoring temperature: high limit / low limit.
- Automatic calibration of a point.



DIGITERM TFT-200, ULTRATERM-TFT-200, FRIGITERM-TFT-10, FRIGITERM-TFT-30

Specific functions of models with touch screen:

- Acoustic and visual alarm.
- Clock and calendar.
- Single or cyclic On / Off programming.
- Up to 10 work programs.
- Up to 6 segments per program.
- Stability time in each segment (from 1 min to 99h).
- Power control program (between 10% and 100%).
- Alarms and events storage.
- Probe error detection.
- Self Diagnostics.
- Network failure detection and saving.
- Over temperature and low temperature alarms and memorization (date, start time, end time and temperature).

- Safety thermostat (TS) by software.
- Mechanic safety thermostat (TS).
- USB and RS -232 output.
- PC software.
- User manual on screen .
- Configurable parameters: Date / time, temperature correction , data collection interval, language (English, Spanish and French) , °C / °F selection , over temperature and low temperature limit.
- Temperature control auto-tuning.
- Circulation pump speed control.
- 10 points automatic calibration.







Immersion thermostat with digital control "Digiterm S-150"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5°C TO 100°C. With external refrigeration from -20 to 100.

High precision temperature control

SAFETY:

STANDARD DIN 12876. 1. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET 2. OVER TEMPERATURE ALARM - 3. LOW LIQUID LEVEL FLOAT

FEATURES

Maximum working volume to maintain maximum temperature: 20 litres.

A clamp and finger screw attach the unit to the tank, alternatively the extendible support bridge accessory can be used.

Minimum fixing depth to the tank: 14 cm.

RS-232 interface for computer control or printing data.

CONTROL PANEL


- 1. Bath temperature display indicator.
- 2. Set temperature display indicator.
- 3. Heater function light indicator.
- 4. Set value decrease button.
- 5. Set value increase button.











Accessories (See page. 106).



Immersion thermostat with TFT touch screen "Digiterm TFT-200"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5°C TO 200°C.

With external refrigeration from -20 to 200.

High precision temperature control

SAFETY:

STANDARD DIN 12876. 1. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET 2. OVER TEMPERATURE ALARM - 3. LOW LIQUID LEVEL FLOAT

FEATURES

Maximum working volume to maintain maximum temperature: 20 litres.

A clamp and finger screw attach the unit to the tank, alternatively the extendible


support bridge accessory can be used.

Minimum fixing depth to the tank: 14 cm.

RS-232, USB and Ethernet output.

CONTROL PANEL

- 1. Set/ current temperature indicator.
- 2. Pump flow indicator (%).
- 3. Stop button.
- 4. Start button.
- 5. Push up parameter.
- 6. Push down parameter.
- 7. Ramps / cycles button.
- 8. Set / current time indicator in hours / minuts.
- 9. Return button.









MODELS	Part No.	Control range °C	Stability °C	Height / Width / Depth (exterior) cm	Pump Pressure / Flowrate H ₂ O	Power W	Weight Kg
DIGITERM TFT-200	3000544	amb.+5 to 200	at 100 °C ±0,05 at 200 °C ±0,1	28 18 19	150 mbar-5,6 l/min.	2060	4,5

Accessories (See page. 106).



Recirculation thermostat "Ultraterm-TFT-200"

ADJUSTABLE TEMPERATURES FROM AMBIENT+5 °C TO 200 °C.





SAFETY:

STANDARD DIN 12876. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET LOW LIQUID LEVEL SENSOR PROTECTION.

FEATURES

Maximum working volume to maintain maximum temperature: 8 litres.
AISI 304 stainless steel top and lid with an inner tank of AISI 304 stainless steel.
Complete with a pump for internal and external liquid circulation with inlet and outlet connectors. A cooling coil to maintain ambient temperatures using tap water is also provided.

A RS-232 interface for print out or computer control comes as standard. A drain tap is also fitted.



MODEL

Part No.	Capacity litres	Control range °C	Stability °C	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Pump Pressure / Flowrate H ₂ O	Power W	Weight Kg
3000545	8	From amb. +5 to 200 °C With external refrigeration	to 100 °C ±0.05 to 200 °C ±0.1	15 20 14	36 28 36	150 mbar-12 l/min	2060	9.5
		from -20 to 200						



Refrigerated recirculation baths "Frigiterm-TFT-10" and "Frigiterm-TFT-30" ADJUSTABLE TEMPERATURE FROM -10 °C TO 100 °C. STABILITY: 0 °C TO 100 °C, ± 0.05 °C / -10 °C, ± 0.1 °C.

ADJUSTABLE TEMPERATURE FROM -10 °C TO 100 °C. STABILITY: 0 °C TO 100 °C, ± 0.05 °C / -10 °C, ± 0.1 °C. ADJUSTABLE TEMPERATURE FROM -30 °C TO 100 °C. STABILITY: 0 °C TO 100 °C, ± 0.05 °C / -30 °C, ± 0.2 °C. SET ERROR : ± 1 °C. RESOLUTION: 0.1 °C.

SAFETY:

STANDARD DIN 12876. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET. LOW LIQUID LEVEL SENSOR PROTECTION.

FEATURES

Maximum working volume to maintain maximum temperature: 8 litres.
AISI 304 stainless steel top and lid with an inner tank of AISI 304 stainless steel.
Complete with a pump for internal and external liquid circulation with inlet and outlet connectors. A drain tap is placed at the side of the unit.
A hermetically sealed compressor is mounted on anti-vibration mounts.

A RS-232 interface for connection to a printer or computer.



MODELS

	Part No.	Capacity litres	Control range		Width / Depth able) cm	Height / Width (exterior)		Pump Pressure / Flowrate H ₂ O	Power W	Compresser power H.P.	Weight Kg
FRIGITERM-10	3000546	8	-10 + 100	12	20 14	36 66	44	150 mbar-12 l/min	1150	1/8	28
FRIGITERM-30	3000547	8	-30 + 100	12	20 14	36 66	44	150 mbar-12 l/min	1460	3/8	30



PRECISE CIRCULATING BATHS Choice of baths and immersion circulators "Digiterm-S-150" and "Digiterm-TFT-200"




CONTROLLABLE TEMPERATURES FROM AMBIENT +5 °C TO 150 OR 200 °C. CAPACITY FROM 20 TO 27 LITRES.

APPLICATIONS

Incubation processes, Enzymatic reactions, Fermentation, laboratory cultures in general.

COMMON FEATURES

Double skinned baths with an AISI 304 stainless steel internal tank with an AISI 304 stainless steel external, for convenience, handles at the side and drain tap are provided as standard. A detachable immersion thermostat for controlling stirring and temperature is fitted. A cooling coil with the use of tap water or cooled re-circulator can be used to maintain ambient temperatures. Internal or external circulation pump (Polarimeters, refractometers and viscometers, etc)



Graph of temperature and time for the Tectron Bio 100 and Digiterm-S-100 and Digiterm-TFT-200, units with capacity of 20 litres.

Tests completed with water up to 100°C:

- **A.** Digiterm-TFT-200 with lid at 100 °C.
- B. Digiterm-TFT-200 no lid at 90 °C.
- C. Digiterm-S-100 with lid at 100 °C.
- D. Digiterm-S-100 with no lid at 90 °C.

Graph of temperature and time for the Digiterm 100 and 200, units with capacity of 20 litres.

Tests completed with Silicon oil up to 200°C:

- A. Digiterm-TFT-200 with lid at 200 °C.
- B. Digiterm-TFT-200 with lid at 100 °C.

27 litres capacity baths take 20% more time to reach the same temperature.



Bath circulators "Digiterm-S-150"

CUVETTES SET OF 20 AND 27 L. CAPACITY, WITH THERMOSTATIC EQUIPMENT INCLUDED. FORT TEMPERATURES FROM AMB+5 TO 150 °C.

STABILITY: ±0.05 °C. RESOLUTION: 0.1 °C.



SAFETY STANDARD DIN 12876. SAFETY ADJUSTABLE OVER TEMPERATURE CUT OFF WITH MANUAL RESET. Low Liquid Level Cut off.




MODELS	Part No.	Temperature Range °C	Cuvettes Capacity litres		: / Widt nterior)	h / Depth cm	Не		/ Width exterior)	n / Depth) cm	Pressure mbar	ump Pump rate I/min	Power W	Weight Kg
DIGITERM-S-150-20	3000548	150	20	15	30	29	3	36	55	35	150	12	1060	9.2
DIGITERM-S-150-27	3000549	150	27	20	30	29	4	11	55	35	150	12	1060	10.8

ACCESSORIES see page 106.



Bath circulators "Digiterm-TFT-200"









SAFETY: CONFORMS TO THE DIN 12876. STANDARD. ADJUSTABLE OVER TEMPERATURE SAFETY CUT OFF WITH ADJUSTABLE MANUAL RESET. LOW LIQUID LEVEL CUT OFF. ELECTRICAL FAULT ALARM.

FEATURES

Temperature sensor probe Pt 100.

Heating element made of stainless steel INCOLOY, resistant to high temperatures and corrosive environments.

Circulation pump with flow rate, suitable for external temperature control for: (Polarimeters, refractometers, viscometers, etc).

Cooling coil to maintain ambient temperatures, using tap water.

External temperature control connection using a Pt 100 probe (see page106).

External metal case with a corrosion resistant epoxy coating.

All parts in contact with liquid are made of AISI 304 stainless steel.

Handles at the side and drain tap.



MODELS	Part No.	Temperatura Range °C	Cuvettes Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Pump Pressure Pump rate mbar I/min	Power W	Weight Kg
TFT-200 20	3000550	200	20	15 30 29	36 55 35	150 12	2060	10
TFT-200 27	3000551	200	27	20 30 29	41 55 35	150 12	2060	12

ACCESSORIES see page 106.

ACCESSORIES for the DIGITERM-S-150 and DIGITERM-TFT-200. Made of AISI 304 stainless steel.



- 1. Gable lid. Part No. 3001295
- 2. Lifting rack support. Capacity 3 tube racks. Part No. 1001296
- 3. Tube racks for lifting rack support.

Part No.	For Ø tubes mm	Capacity Rack tubes		t/Widtl xterior)	n/Depth cm
1001202	13	36	8	8.2	23.6
1001203	16	24	8	8.2	23.6
1001204	20	24	8	8.2	23.6



Immersion thermostat for baths "Digit-Cool"

DIGITAL SELECTOR AND TEMPERATURE CONTROL ADJUSTABLE FROM 0 °C TO 100 °C. STABILITY: 100 °C \pm 0.1 °C, and at 0 °C \pm 0.05 °C. Set error: \pm 0.1 °C. resolution: 0.1 °C.



SAFETY: STANDARD DIN 12876. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET. LOW LIQUID LEVEL SENSOR PROTECTION.

FEATURES

Maximum volume to achieve the working temperature: 20 litres.

Minimum depth of the elements with contact with the liquid: 15 cm.

Hermetically sealed compressor.

RS-232 out put to a computer control or printer.

Portable control and stirring unit with lifting handle.

CONTROL PANEL

Main power switch.

LCD Display showing working parameters.

Push button configuration selector.

Push button parameter selection.

Push button increase/decrease set value.

Start / Stop, push button.

MODEL

MODEL						
Part No.	Control Range	Height / Width / Depth	Pump	Power	Power	Weight
	O°	(exterior) cm	Pressure / Flow H ₂ O	W	H.P.	Kg
3001373	0 +100	45 21 51	150 mbar / 12 l/min	1460	1/5	22

ACCESSORIES See page 106.

	Autonomous unit with cooling system	
--	-------------------------------------	--

Compressor unit built-in





Digitcool thermostat with bath and 4 place ring set lid.

ACCESSORIES for thermostatic unit digit-cool.

Stainless steel baths.

Double skinned with an AISI 304 stainless steel exterior with a bath of AISI 304 stainless steel interior. Supplied with a drain tap and lifting handles.





1. Lifting rack support. . Capacity 3 tube racks. Part No. 1001296

2.	Tube	racks	for	lifting	rack	SI	ıpp	ort.	
_		-	α.		_		_		٠.

Part No.	For Ø tubes mm	Capacity Rack tubes	•	t / Widt exterior	h / Depth) cm
1001202	13	36	8	8.2	23.6
1001203	16	24	8	8.2	23.6
1001204	20	24	8	8.2	23.6



Part No.	Capacity litres	•	ht / Width / Depth Height / (usable) cm (ex		Width terior)		Weight Kg	
6000392	20	15	48	30	20	55	35	6.6
6000393	27	20	48	30	25	55	35	7.6



Ring set lid 105, 80, 60 and 37 mm Ø reduction rings and a check thermometer location hole in the lid. Capacity 4 places.

Part No. 1001374

Digitcool thermostat with bath and 4 place ring set lid.

ACCESSORIES FOR IMMERSION THERMOSTATS





Pt100 temperature sensor with handle

4 mm \emptyset x 135 mm long submersion. Supplied with a 150cm cable and connector. Part No. 1000893



Made of AISI 304 stainless steel, adjustable. to fit the immersion thermostat models "Tectron and "Digiterm".

Support for apertures from 22 to 44 cm across. Part No. 6001092



Example showing the support bridge. Part No. 6001092.

Stainless steel baths.

Double skinned with an AISI 304 stainless steel exterior and internal bath of AISI 304 stainless steel.

MODELS

Part No.	Capacity litres	•	Width sable)	/ Depth cm	Height / (ex	Width terior)		Weight Kg	
6000390	9	15	29	22	20	34	28	3.5	
6000391	12	15	30	31	20	37	35	4.3	
6000392	* 20	15	48	30	20	55	35	6.6	
6000393	* 27	20	48	30	25	55	35	7.6	
6001093	* 45	15	62	50	23	70	56	11	

^{*} With drain tap and lifting handle.

Stainless steel baths, thermally insulated. Ideal for low temperatures. Suitable for J.P. Selecta refrigerated units. Similar to our other bath but with additional internal expanded foam insulation that prevents external condensation or heat loss, working temperature range from -40°C to +90°C.

MODELS

MODELO									
Part No.	Capacity litres	•	Height / Width / Depth (usable) cm		Height / Width / Depth (exterior) cm		Weight Kg		
6003901	9	15	29	22	24	40	34	6	
6003921*	20	15	48	30	24	61	41	9	
6003931*	27	20	48	30	29	61	41	11	

^{*} With drain tap and lifting handles.



Removable constant liquid level. Suitable for any type of bath. Adjustable height to obtain the required liquid level.

Part No. 6001400 Bath depth up to 15 cm. 6001800 Bath depth up to 20 cm.



Thermometer support. Stainless steel, suitable for thermometer. Fixation aperture up to 40 mm. Part No. 6000896



Methacrylate baths. Maximum temperature 60 °C.

MODELS

MIODELO					
Part No.	Capacity litres	Height / (u:	Thickness mm		
1000394	8	15	18	31	6
1000544	9,5	14,5	16,5	39	6
1000395	13	15	23	40	8
1000396	30	20	30	50	8



Thermometer.



Suitable for "Tectron-Bio" immersion thermostat. Scale -20 to 100 °C. Part No. 3009100



Refrigerated cooling coils for baths "Frigedor" and "Frigedor-Reg"

TEMPERATURES FROM -20 °C TO +20 °C.

APPLICATIONS

Designed for bath and tank applications that require below ambient temperatures.


COMMON FEATURES

The unit is bench mountable and contains within the epoxy coated case a CFC free hermetically sealed compressor with condenser and evaporator, the cooling coil is made of AISI 304 stainless steel.

COIL DIMENSIONS

Refrigerated immersion length: 900 mm.

Cooling coil Ø: 45 mm. Coil length: 150 mm.



Graph showing the cooling performance of the 8 litre H_2O "Frigedor" With insulated stirring tank.

MODEL FRIGEDOR

No temperature controller incorporated. Continually operates the compressor.

CONTROL PANEL

Mains switch with "ON" indicator lamp.

MODEL FRIGEDOR REG

Equipped with a temperature controller with digital control and display. Includes a Pt 100 temperature probe. Resolution: 1 digit.

CONTROL PANEL

Mains switch with indicator "ON" lamp. Digital temperature controller with push buttons, connector for the Pt 100 temperature probe. (See page 106).





MODELS	Part No.	Temperature range °C	Stability °C	Height / Width / Depth (exterior) cm	Cooling potential	Power W	Weight Kg
FRIGEDOR	3000778	-20 to +20	-	41 21 34	to -20 °C = 50 W	285	14
FRIGEDOR-REG	3001214	-20 to +20	±1.5	41 21 34	to -20 °C = 50 W	285	14

We recommend our thermally insulated baths (see page 106).



Water recirculators "Interheater"

FROM AMBIENT +5°C TO 90°C.

SAFETY: STANDARD DIN 12876. OVER TEMPERATURE SAFETY ADJUSTABLE THERMOSTAT WITH MANUAL RESET.

DESCRIPTION

Water circulation unit, designed to feed closed water loop circuits for condensers, distillation columns, reactor jackets, viscometers, electrophoresis baths etc. designed specifically for:

- Constant temperature.
- Constant pump rate of refrigerated fluid.
- Closed loop circuit, avoids the build up of scale in cooling coils,and cooling circuits.

FEATURES

Digital electronic temperature regulation.

Circulation flow: 20 L/min.

Inox covered steel case.

RS-232, USB and Ethernet output for recording tem-

peratures.

Hose connectors \emptyset 12mm internal.

230V single phase.

INTERHERATER Heating power 2 or 4 KW control panel selectable.

CONTROL PANEL

Mains switch.

Digital temperature control.

Pressure gauge.

Water level in the reservoir.

"Autotunning" function for optimizing the temperature control of each facility.

Supports 9 programmable temperature profiles including temperature ramps and steps.

MODEL

MODEL							
Part No.	Part No. Height / Width / Depth (exterior) cm			Regulation range °C	Power W	Weight Ka	
	(0/1	101101	OIII	•	••	ng -	
6001491	70	43	56	5 to 90	4100	22	

ACCESSORY

Roll of 25 m of reinforced silicone hose.

Part number 1001540



WITH DIGITAL ELECTRONIC CONTROL OF TEMPERATURE AND TIME

SAFETY: CONFORMS TO EN 61010 DIRECTIVE, OVER TEMPERATURE CUT OUT SYSTEM.

High precision. Free from contamination.

APPLICATIONS

Clinical, biochemical and chemical laboratories. Incubation of DNA, enzyme assays, residual testing in dairy, incubation and fusion of AGAR, cell culture. Thermo control to dry, boiling, evaporation, concentration, hydrolysis, digestion etc.

COMMON FEATURES

Dry block systems have several advantages over wet contact methods such as: no risk of evaporation of heating liquid, higher working temperatures, better long term temperature stability and no sample cross contamination or germination of bacteria in the bath.

Digital temperature calibration. Timer from 1 to 999 mi

nutes or continuous operation. Over temperature alarm. Heating elements are distributed evenly across the contact surface maintaining an even distribution of heat.

Easy to use digital temperature and time control with a digital display (3 Digits).

Built in temperature sensor Pt 100.

Interchangeable blocks made of anodised DURAL alloy to accommodate a variety of tube diameters.

(Customer specific sizes also available.)

External case made of epoxy covered steel with an AISI 304 stainless steel top.



CONTROL PANEL

- 1. Digital display showing temperature and time in minutes.
- **2.** Indicator show displayed parameter (Time or temperature.).
- 3. Push button selector to show time or temperature.
- 4. Push button increase displayed value.
- 5. Push button decrease displayed value.
- **6.** Push button Start / Stop.
- 7. Over temperature alarm indicator.



Dry block heater for tubes "Tembloc"



ACCESSORIES

Metal blocks, size: 105 Ø x 55 mm high.

Part No.

7000346 for 34 x 6 mm Ø tubes

7000714 for 18 x 1.5 ml. Eppendorf tubes

7001224 for 24 x 12 mm Ø tubes.

7001618 for 18 x 16 mm Ø tubes.

7000208 for 8 x 20 mm Ø tubes.

7000256 for 6 x 25 mm Ø tubes.

7000715 Blank block with no holes for customized demand.

MODEL


Part No.	Temperature	Stability	Homogeneity	Set error	Resolution	Height / Width / Depth	Power	Weight
	°C	°C	%	%	°C	(exterior) cm	W	Kg
7462200	30 - 200	±0.75	±1.5	±2	1	11 19 29	350	3.7

Supplied complete with thermometer 0-200 °C and extraction tool for blocks.



Dry Block heater for tubes "Multiplaces"

CAPACITY: THREE BLOCKS.



ACCESSORIES

Metal Blocks: 55 high x 95 wide x 75 mm deep.

Part No.

7000306 For 30 x 6 mm \emptyset tubes (total 3 blocks = 90 tubes).

7000716 For 20 x 1.5ml, Eppendorf tubes. (total 3 blocks = 60 tubes.)

7031220 For 20 x 12 mm \emptyset tubes (total 3 blocks = 60 tubes).

7031612 For 12 x 16 mm \emptyset tubes (total 3 blocks = 36 tubes).

7003208 For 8 x 20 mm \emptyset tubes (total 3 blocks = 24 tubes).

7001256 For 6 x 25 mm Ø tubes (total 3 blocks = 18 tubes).

7000717 Blank block without holes, can be customised to customer requirements.

7001474 Block for 15 cuvettes 10 x 10. Capacity: 3 blocks.

7001475 Block for 96 well microplates of 1.2 ml. Capacity: 1 block.

MODEL

Part No.	Temperature °C	Stability °C	Homogeneity	Set error	Resolution °C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7471200	30 - 200	±0.75	±1.5	±2	1	11 31 31	700	7.4

Supplied complete with thermometer 0-200 °C and extraction tool for blocks.



Metallic thermostat dry block "Clinic-Bloc" and "Bio-Bath"

FIXED TEMPERATURE OF 37 °C WITH CHECK THERMOMETER.

SAFETY: NORM EN 61010. OVER TEMPERTURE CUT-OUT FITTED.

"CLINIC-BLOC"

Capacity 20 tubes (5x4) of 75x13 mm. VAC. Part No.7001569

"CLINIC-BLOC" 1537

Capacity 15 cuvettes (5x3) of 10x10 mm. Part No. 7001570

MODEL

Part No.	•	/ Widt terior)	h / Depth cm	Stability °C	Power W	Weight Kg
7001569	11	18	28	±0.5	10	2.4
7001570	11	18	28	±0.5	10	2.4

"BIO-BATH"

Capacity 8 tubes (5x4) of 75x13 mm VAC.

MODEL

Part No.	Height / Width / Depth	Stability	Power	Weight
	(exterior) cm	°C	W	Kg
7001561	8.8 9.6 7	±0.5	8	0.5









Thermo shaker "TRM-4" for Microtiter plates

WITH ELECTRONIC CONTROL AND DIGITAL DISPLAY OF TEMPERATURE, SHAKING SPEED AND TIME. FOR TEMPERATURES FROM AMBIENT +5 °C. TO 60 °C. STABILITY \pm 0,5 °C.



APPLICATIONS

For clinical analysis, immunology and nutritious quality control use.

FEATURES

External ABS case with a platform for heating and shaking 2-4 96-well Microtiter plates

Upper lid with aluminium plate coated inside, and simultaneous heating system to the main platform.

Digital display of current and set point temperature, time and shaking speed parameters.

Acoustic and luminous indications at the end of the cycle and due to malfunction in temperature or time programmed.

A built-in motor generates a uniform movement depending on the speed programmed, protected to over temperature with auto switch off and switch on activation. Temperature calibration to meet the user's needs.

Protection device for over temperature that provides safety and reliability.

CONTROL PANEL

LCD display indicating temperature, speed in r.p.m. and time.

Push button SET to program temperature, speed in r.p.m. and time.

Push button ▲ increase value.

Push button ▼ decrease value.

Push button START / STOP for starting and stopping the shaking system.

Push button OFF for any function.

SPECIFICATIONS

Temperature: Ambient +5°C to 60°C. Heating time: ≤23 min (to 60°C). Shaking range: 100-1200rpm. Orbital rotation: 2mm.

Timer: From 1 min. to 99h. 59 min.





MODEL

Part No.	Temperature	Stability	Block	Height / Width	Height / Width / Depth	Power	Weight5
	°C	°C	homogeneity °C	(platform) cm	(exterior) cm	W	Kg
5109200	ambient +5 to 60	±0.5	±0.5	21.5 30	18 40 39	132	9.5



Metallic dry block thermo shaker "TR100-G"



WITH ELECTRONIC CONTROL AND DIGITAL DISPLAY OF TEMPERATURE, SHAKING SPEED AND TIME. FOR ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C. UP TO 100 °C.

APPLICATIONS

Wide applications of TR100-G for DNA analysis, lipids and other cellular components extraction, DNA libraries creation, DNA amplification, electrophoresis pre-denaturalization, serum solidification, etc.

FEATURES

External ABS case with multiadapter device for different blocks depending on the analysis to be made.

Digital display of current and set point temperature, time and shaking speed parameters. Acoustic and luminous notification at the end of the cycle and due to malfunction in temperature or time programmed.

A built-in motor generates a uniform movement depending on the speed programmed, protected to over temperature with auto switch off and switch on activation. Temperature calibration to meet the user's needs.

The thermostat includes a heater platform to adapt multiple blocks, depending on the analysis to be made.

Protection device for over temperature that provides safety and reliability.

CONTROL PANEL

LCD display indicating temperature, speed in r.p.m. and time.

Push button SET to program temperature, speed in r.p.m. and time.

Push button A increase value.

Push button ▼ decrease value.

Push button START / STOP for starting and stopping the shaking system.

Push button OFF for any function.



SPECIFICATIONS

Temperature: Ambient +5°C to 100°C. Heating time: ≤23 min (to 100°C). Shaking range: 200-1500rpm.

Orbital rotation: 2mm.

Timer: From 1 min. to 99h. 59 min.

MODEL

Part No.	Temperature	Stability	Block	Height / Width / Depth	Power	Weight
	°C	°C	homogeneity °C	(exterior) cm	W	Kg
5109100	ambient +5 to 100	±0,5	±0,5	18 21 30	88	9.5

ACCESSORIES

ABS coated metallic blocks adapted to TR100-G thermostat by means of fixing screws to the block. Easy to clean and autoclavable.



For 96 microtubes of 0,2 ml. with polypropylene lid. Part No. 5109104

For 54 microtubes of 0,2 ml. with polypropylene lid. Part No. 5109105

For 15 microtubes of 0,5 ml. and 20 microtubes of 1,5 ml. with polypropylene lid. Part No. 5109106 For 35 microtubes of 1,5 ml. with polypropylene lid. Part No. 5109107

For 35 microtubes of 2 ml. with polypropylene lid. Part No. 5109108

Thermocycler for thermal cycles "B960" See page 238.



Hand operated centrifuges page 114
Microtubes and vials centrifuge page 114
Centrifuge angle rotor for 15 ml. page 114
Digital electronic control centrifuges pages 115 to 117
High speed centrifuges and refrigerated pages 118 to 131
TFT Touch screen control

Centrifuges for flasks and blood bags page 123 to 131
Tubes of glass and plastic and buckets page 132
Kit for blood and plasma centrifuges page 133



Centrifuge comparison by model and tube capacity

Description

Capacity

Dimensions

Description

Capacity

Reaction

vial

Ø 6,3/7

47/49 mm

0.25/0.5 ml 1.5/2.2 ml

Micro-

haematocrit

Ø 1.4

75 mm

Reaction

vial

Ø 10.7/13

42 mm

Tube

RIA

5 ml

Ø 12

Hemolisis

5 ml

Ø 12.5

75 mm

Tapval or

VAC tube

5 ml

Ø 13 X

Glass

7 ml

Ø 12

100 mm



	To the second				Difficusions			55/75 mm	75/80 mm		
SUMMARY TABLE OF CENTRIFUGES BY MODEL IN RELATION TO VOLUME, CAPACITY SPEED AND MAXIMUM RCF (G). QUICK OVERVIEW											
Model	Part No.	Max. Volume ml	Max. Volume of tube ml	R.p.m. max.	R.c.f. max. (xg)						
	7000014	30	15	1500	360						
Manual	7000015	60	15	1500	360						
	7000016	2 x butyrometer	-	1500	370						
Cencom II	7002240	90	15	4000	1760						
Centro 4-BL	7002560	20	5	5000	2907				4		
Centro 8-BL	7002356	120	15	4500	2515				8	8	
Centrolit II-BL	7002357	40,5	2,2	12000	14103	24	18	18			
Centromix-BLT	7002556	120	25	12000	14200				8	4	
Mixtasel-BLT	7002576	400	100	4200	2901				24 16	24	
Centronic-BLT	7001768	400	100	13500	19048	24	40	30	24 16	24	
Centrofriger-BLT	7001769	400	100	15300	23550	24	40	30	24 16	24	
Meditronic-BLT	7001476	800	200	11000	12619			30 48	48 28	48	
Medifriger-BLT	7001477	000	200	15000	23470			00 40	40 20	40	
Macrotronic-BLT	7001485	2000	500	10000	12296			80	80 48	96	
Macrofriger-BLT	7001486	2000	300	12000	17706			00	00 40	30	
Macrofriger-BLT-Blo	od 700 1490	9600	1000	4070	5370					144	

See page 132 for tubes and accessories. Tubes supplied as an accessory unless stated.

Plastic with cap	Glass	Glass	Plastic with cap	Glass	Falcon	Glass	Bottle	Glass	Bottle	Bottle	Butyrometer	ASTM
10 ml	15 ml	25 ml	30 ml	50 ml	50 ml	100 ml	200 ml	250 ml	500 ml	1000 ml		
Ø 16.1 x 81.1 mm	Ø 16 x 100 mm	Ø 24 x 100 mm	Ø 25.3 x 92 mm	Ø 34 x 100 mm	Ø 29.5/35 x 117 mm	Ø 44/45 x 100 mm	Ø 56.5 x 104 mm	Ø 56 x 147 mm	Ø 69 x 159 mm	Ø 98/100 x 160-175 mr	n	Ø 58 x 158 mm
											THE	M M
Tapval or VAC tube	Falcon				Plastic with cap				Blood bag	Microtiter		ASTM
10 ml	15 ml				80/85 ml				500/600 ml	6x96		100 ml
Ø 16 X 100 mm	Ø 16/23 x 120 mm				Ø 38 x 112 mm					Ø 85 x 130 mm		Ø 37 X 200 mm
	2 4											
	6 6										2	
	8											
4	8	4										
16	18 16	8	8	8	4	4						
10 16	18 16	8	8	8	4	4						
10 16	18 16	8	8	8	4	4						
20 28	28 12	8	6	8	4 6	6	4			6x96		
20 68	68 16	36	20	16	8 6	4	4	4	4	12x96		4 4
	102	54		24		12	6		6 12	6		



Hand operated centrifuges

FEATURES

For 15 ml and butyrometer tubes.

Comes with a table holding clamp.

Adaptable heads for 2 or 4 x 15 ml tubes or for 2 butyrometer tubes.

MODELS

Part No.	Comes with head for	Max. speed r.p.m.	Max. R.c.f.	Weight Kg
7000014	2 tubes	1500	360	1.1
7000015	4 tubes	1500	360	1.2
7000016	2 butyrometers	1500	370	1.2

Comes with head and metal buckets.



ADDITIONAL ACCESSORIES

Head for two tubes. Part No. 7000017 Head for four tubes. Part No. 7001388

Bucket for 15 ml tubes.
Part No. 7001500
Head for two butyrometers.
Part No. 7000018





COMECTA Microcentrifuge "Cencom I"

FOR MICROTUBES OF 2,0 - 1,5 - 0,5 AND 0,2 ML.

DIGITAL ELECTRONIC CONTROL OF SPEED AND TIME. ADJUSTABLE SPEED UP TO 14000 r.p.m.

FEATURES

ABS external case.

Maintenance free induction drive motor.

Rubber suspended motor which reduces the vibrations. High speed microcentrifuge that minimizes the samples over temperature.

Extremely quiet operation.

For 12 microtubes with identification code.

Polycarbonate upper lid.

Peephole on the lid to visualize the centrifuge inside.

CONTROL PANEL

Start switch.

Push button Start/Stop.

LCD display indicating time and speed.

Time selector push button.

Speed selector push button.

Push button ▲ increase different parameters.

Push button ▼ decrease different parameters.

Push button motor drive rotation for quick accele-

ration Short/Spin.

Push button open lid when cycle has ended.



MODEL

Part No.	Tubes capacity	Height/ Width/ Depth (exterior) cm	Max. speed r.p.m.	R.c.f.(xg)	Power W	Time selected	Weight Kg
5022000	12	15 21 25	14500	14000	105	till 99' 99"	4

Comes complete with 12 adapters of 1.5-0.5 and 0.2ml.



Centrifuge Angle rotor "Cencom II"

WITH ANGLE ROTOR FOR 6 X 15 ml tubes and 6 additional adapters for 15 ml falcon tubes. Digital electronic control of speed and time. Adjustable speed up to 4000 r.p.m. (1790 r.c.f.)

SAFETY:

SAFETY ACCORDING TO EN. 61010.2.020 SAFETY MICRO-SWITCH THAT DISCONNECTS THE POWER TO THE MOTOR WHEN THE LID IS OPENED.

PROTECTIVE FUSE LOCATED AT THE BACK OF THE INSTRUMENT.

APPLICATIONS

General laboratory use and teaching.

FEATURES

External case made from coated ABS plastic.

Rubber suspended motor, reduces any motor vibrations.

Fixed angle 6 place rotor with identification code.

Transparent polycarbonate lid that allows the user to see inside the centrifuge. Speed selection from 1000 to 4000 r.p.m. in 500 r.p.m. steps.

Time selection from 1' to 60' in steps of 1'.

CONTROL PANEL

Digital display of speed and time.

Push button Start-Stop.

Push button selector for speed and time.

Push button increase value. Push button decrease value. Cross section view of the angle rotor showing bucket and bucket support.





MODEL

Part No.	Tube capacity	Height/Ø	Max. speed	Max.	Power	Weight
	ml	(exterior) cm	r.p.m.	R.c.f.	W	Kg
7002240	6 x 15	28 26	4000	1760	80	4.5

Supplied complete with 6 x 15 ml tubes with adapters and 6 adapters for 15 ml Falcon tubes.

ACCESSORY

Adapter. Suitable for 13 \emptyset x 75 mm VAC tubes. Rotor capacity, 6 adapters.

Part No. 1001556

Digital microprocessor controlled Centrifuges,

Series BL and BLT



INDUCTION DRIVE, MAINTENANCE FREE (NO CARBON BRUSHES).

- -SPEED CONTROL BY VARIABLE FREQUENCY SYSTEM.
- -RELIABLE FUNCTIONS AND REPEATABILITY.
- -LOGICAL EASY TO USE CONTROL PANEL.
- -SUSPENDED, DAMPENED VIBRATION FREE MOTOR DRIVE, GUARANTEES SMOOTH OPERATION.
- -QUIET OPERATION, INTERNAL INSULATION WHICH ABSORBS NOISE AND VIBRATIONS.
- -SUCTION SUPPORT FEET THAT PREVENTS BENCH CREEP.

SAFETY:

AUTOMATICALLY STOPS IF OUT OF BALANCE.
LID CANNOT BE OPENED IF THE MOTOR IS RUNNING.
SYSTEM CANNOT BE STARTED IF THE LID IS OPEN.
SPEED LIMITED BY EACH ROTOR.

TWO DIFERENT CONTROL SYSTEMS:



Digital display and



TFT touch screen

CENTROLIT II-BL - CENTRO 4-BL - CENTRO 8-BL



Specific functions of models with digital screen:

- -Electronic control of speed.
- -Electronic control of time.

CENTROMIX BLT - MIXTASEL BLT - CENTRONIC BLT - CENTROFRIGER BLT - MEDITRONIC BLT - MEDIFRIGER BLT - MACROTRONIC BLT - Macrofriger blt - Macrofriger blt-blood



Specific functions of models with touch screen:

- Acoustic and visual alarm
- 10 working programs with 6 parameters for each program: acceleration, braking, rpm, time, type of rotor.
- Network failure detection and storage.
- Safety thermostat (refrigerated models only).
- USB output.

- Software for PC
- Temperature correction (refrigerated models only).
- Languages: Spanish, English and French.
- Overtemperature and Low temperature alarm.
- Manual user on screen.
- Rotors list available on screen.
- Door open alarm.
- Centrifugation cycles with up to 4 time-speed steps.
- Rotor and motor running hours log.



CENTRIFUGES ELECTRONIC DIGITAL CONTROL "Centrolit II-BL", "Centro-4-BL", "Centro-8-BL", With induction drive maintenance free motor.



UNE-EN 61010.2.020. ELECTROMAGNETIC LID LOCK. PREVENTS THE OPENING OF THE CENTRIFUGE IF THE ROTOR IS IN MOTION. PREVENTS THE STARTING OF THE MOTOR IF THE LID IS OPEN. ROTOR AND BUCKETS ARE IDENTIFIED BY LASER SIGNAL. AUTOCLAVABLE.

COMMON FEATURES

Electronic control with digital display of speed and time, with acceleration and brake ramp.

Timer controllable from 1' to 60' or continuous operation.

Rubber suspended motor that guarantees a smooth operation.

Digital speed display.

Pulse button, that drives the motor while pressed.


Automatic stop if out of balance, with indicator light.

Suction base supports that adhere to the work bench surface.

Quiet operation, insulated interior that absorbs noise and vibrations. In the case of mains failure, manual open override aperture.

Noise level 50-60dBA.

Rotors and buckets with laser identification. Autoclavable.



CONTROL PANEL

- 1. Speed display.
- 2. Push button increase speed.
- 3. Push button reduce speed.
- 4. Run Time display.
- 5. Push button increase run time.
- 6. Push Button decrease run time.
- 7. Push button pulse start (RUN/STOP).
- 8. Push button start.
- 9. Push button stop.
- 10. Push button open lid.



Centrifuge "Centrolit II-BL"SUITABLE FOR HAEMATOCRIT AND SMALL VIALS OF 0.25/0.5 ml and tubes of 1.5/2.2 ml. FREE MAINTENANCE INDUCTION DRIVE MOTOR.





ACCESSORIES

Interchangeable rotors made of hardened anodised aluminium. Autoclavable.

Part No. head	Rotor type	Volume rotor	Tube capacity maximum	Max. R.p.m.	Max. R.c.f. (xg)
7001704	Micro-haematocrit	-	24 x 1.4 mm Ø	12000	13143
7001703	Angle+lid	40.5 ml	18 x 0.25/0.5 ml + 18 x 1.5/2.2 ml	12000	14103
7001702	Angle PCR	6,4 ml	32 x 0.2 mm Ø	12000	12500

Angle rotors supplied without tubes.

Micro-haematocrit rotor supplied with one box of 100 tubes.



Micro-haematocrit reader. Suitable for micro-haematocrit rotor 7001704. Part No. 7000929



Centrifuge "Centro-4-BL"WITH SWINGOUT HEAD FOR 4 TUBES OF 75 X 13 mm VAC. FREE MAINTENANCE INDUCTION DRIVE MOTOR.



Control panel features (see page 116).

Centrifuge kit for separating plasma rich platelets from whole blood.





MODEL

Part No.	Height / Width / Depht	Power	Max.	Max.	Weight
	(exterior) cm	W	R.p.m.	R.c.f. (xg)	Kg
7002560	20 28 34	170	5000	2907	6

Rotor comes without tubes.

ABS plastic moulded external case.



Centrifuge "Centro-8-BL"

WITH ANGLE ROTOR FOR 7-15 ml Tubes and 75 X 13 mm vac (with adapters). FREE MAINTENANCE INDUCTION DRIVE MOTOR.



Control panel features (see page 116).



MODEL

Part No.	Max. Volume ml.	Tube capacity		Width / Depth terior) cm	Power W	Max. R.p.m.	Max. R.c.f. (xg)	Weight Kg
7002356	120	8 x 15 ml	20	28 34	170	4500	2515	6.5

Supplied complete with 8 x 15 ml glass tubes.



Bucket adapter for 7 ml tubes. Rotor capacity: 8 adapters. Part No. **7000700**

Bucket adapter for 75 x 13 mm VAC. tubes.

Rotor capacity: 8 adapters. Part No. 1001557

SPARES. Tubes (see page 132).




MICROPROCESSOR CONTROLLED CENTRIFUGES "Centromix BLT", "Mixtasel BLT", "Centronic BLT" and "Centrofriger BLT"

COMMON FEATURES

Robust all metal construction:

- -External case made from DUR-Al alloy.
- -Internal bowl and top plate all stainless steel
- -Internal steel safety chamber within the case.



Rotor identification head (only models "Centronic BLT" and "Centrofriger BLT") when the centrifuge lid is closed the rotor communicates all the data of the functions of it. This prevents, before use, any risk of failure of using the wrong rotor in a program or exceeding the rotors running parameters. The advantage of this system is that we can add additional rotors to the range without changing the centrifuge model.

Electronic circuits, microprocessor digital control of all parameters and functions: speed, acceleration, brake, RCF, temperature, timer, 10 program storage, automatic rotor recognition and alarm system.

Induction drive maintenance free motor.

Low noise level 50-60 dBa.

Rotors, buckets and adapters are identified with an inalterable engraved laser code. Air circulation system, to limit any excessive rise in temperature within the centrifuge, all air circulated within the chamber is channelled out of the back and is not re-circulated, this means that the centrifuges internal temperature will remain low and constant.

Hermetically sealed compressor, mounted on anti-vibration mounts with re-circulating evaporator around the centrifuge chamber. (only applies to refrigerated centrifuges).

CONTROL PANEL



- 1. Enter parameters key.
- 2. Speed analog indicator.
- 3. RPM reading kev.
- 4. RCF reading kev.
- 5. Acceleration ramp programmed.
- **6.** Temperature programmed. (refrigerated models only).
- 7. Reading the real time temperature. (refrigerated models only).
- 8. Rotor model in use.
- 9. Door opening key.
- 10. RPM or RCF programmed.
- 11. RPM or RCF real-time reading
- 12. Start process key.
- 13. Stop process key.
- 14. Programmed braking ramp.
- **15.** Time set-up in process.
- 16 Elapsed time reading.



Centrifuge "Centromix BLT"

WITH A CHOICE OF ROTORS AND INTERCHANGEABLE ADAPTERS. FREE MAINTENANCE INDUCTION DRIVE MOTOR.



External DUR-AL case, with AISI 304 stainless steel bowl.



HEAD AND BUCKET REDUCERS

Reducers for bucket. Part No. 7001066: 7001813 For 1 tube 15 ml (16,5 Ø x 100 mm). 7001814 For 1 tube 10 ml (16 Ø x 80 mm). 7001815 For 1 tube 7 ml (12 Ø x 100 mm). 7001816 For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for bucket. Part No. 7001067: 7000700 For 1 tube 7 ml (12 Ø x 100 mm). 1001557 For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for bead. Part No. 7000555: 7000700 For 1 tube 7 ml (12 Ø x 100 mm). 1001557 For 1 tube 7 ml (12 Ø x 100 mm).

Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depht (exterior) cm	Power W	Weight Kg
7002556	120	4 x 25 ml	23 30 35	180	12

CCESSORIES Rotors and interchangeable buckets and adapters.										
Part No. head	Rotor type	Part No. bucket	Totel max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg			
7000555	Angle	-	150 ml	10 x 15 ml	16.5 x 100	4400	2063			
7001065	Swing out	-	-	-	-	-	-			
	Swing out	7001066	100 ml	1 x 25 ml	24 x 100	4400	2601			
	Swing out	7001067	60 ml	1 x 15 ml	16.5 x 100	4400	2601			
7001065	Swing out	7001068	28 ml	1 x 7 ml	12 x 100	4400	2601			
	Swing out	7001573	40 ml	2 x 5 ml VAC	13 x 75	4400	2163			
	Swing out	7001070	60 ml	VAC	16 x 100	4400	2601			
	Swing out	7001072	-	undrilled	-	4400	2601			
7001702	Angle PCR	-	6,4 ml	32 x 0.2 ml	-	12000	12500			
7001704N	/licrohaemato)	-	24 x 1.4 mm Ø	-	12000	13143			
7001703	Angle+lid	-	40.5 ml	18 x 0.25/0.5 ml	-	12000	14103			



Polyvalent with rotors and interchangeable adapters

Control panel features (see page 118).





External case of DUR-AL alloy with ABS plastic coated front AISI 304 stainless steel bowl.

7002576	400	4 x 100 ml	34 41 54	210	22
Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depht (exterior) cm	Power W	Weight Kg
MUDEL					

ACCESSORIES

Swing out rotor with interchangeable sealed buckets and adapters, each with identification code. Autoclavable.





















7000731

7000734

7000736

7001647

7000737

7000738 7001648

Adapters

7000609

Part No. head	Rotor type	Bucket open	Sealed bucket	Part No. adapter	Totel max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	Swing out	7000731	7000732	7000734	200 ml	1 x 50 ml	34 x 100	5000	4137
	Swing out	7000731	7000732	7000735	200 ml	2 x 25 ml	24 x 100	5000	4137
	Swing out	7000731	7000732	7000736	420 ml	7 x 15 ml and VAC*	16.5 x 100	5000	4137
7009079	Swing out	7000731	7000732	7001647	240 ml	6 x 10 ml	17 x 110	5000	4137
	Swing out	7000731	7000732	7000737	336 ml	12 x 7 ml	12 x 100	5000	4137
	Swing out	7000731	7000732	7000738	240 ml	12 x 5 ml	12.5 x 75	5000	4137
	Swing out	7000731	7000732	7001648	140 ml	7 x 5 ml and VAC*	13 x 75	5000	4137
	Swing out	7000731	7000732	7000739	144 ml	12 x 3 ml	11 x 100	5000	4137
	Swing out	7000731	7000732	7000740	-	undrilled	-	5000	4137
7000609	Swing out	-	-	-	-	6 x 96 microtiter	-	4200	2702

^{*}VAC: Blood tube with stopper.

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.

^{**} Bucket part nos. 7000731 and 7000732 accommodate 1 x 200 ml polypropylene bottle with cap without the need of an adapter.

ACCESSORIES

Interchangeable angle rotors with lids, made of hardened anodised aluminium. Rotors come with identification system. Autoclavable rotors and buckets.

Part No. head	Rotor type	Max. volume	Tube capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
7001705	Angle+lid	270 ml	18 x 15 ml	16.5 x 100	4200	2669
7001711	Angle+lid	400 ml	8 x 50 ml	34 x 100	4200	1921





Swing out rotor, interchangeable buckets and adapters each with identification code. Autoclavable.

BUCKETS



7001326























7001700

Sealed bucket with hermetically

sealed lid.

7000795

ADAPTERS













7001315

7001316

7001317 7001318

7001320 7001319

7001322 7001321

Part No. Rotor	Rotor type	Part No. bucket	Sealed bucket	Part No. adapter	Totel max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	Swing out	7001700	-	-	400 ml	1 x 100 ml	44 x 104	4200	2901
	Swing out	7000793	-	-	200 ml	1 x 50 ml	34 x 100	4200	2901
	Swing out	7000794	-	-	100 ml	1 x 25 ml	24 x 100	4200	2901
	Swing out	7000795	-	-	240 ml	4 x 15 ml	16.5 x 100	4200	2901
	Swing out	7000796	-	-	140 ml	5 x 7 ml	12 x 100	4200	2901
	Swing out	7000797	-	-	100 ml	5 x 5 ml	12.5 x 75	4200	2441
	Swing out	7000930	-	-	240 ml	4 x 15 ml Falcon	16.5 x 120	4200	3030
	Swing out	7000931	-	-	200 ml	1 x 50 ml Falcon	30 x 117	4200	3030
	Swing out	7001073	-	-	16 VAC	4 VAC*	16 x 100	4200	2901
7001326	Swing out	7001074	-	-	16 VAC	4 VAC*	13 x 75	4200	2441
	Swing out	7000798	-	-	undrilled	-	-	4200	2901
	Swing out	-	7001315	100 ml glass tube**	400 ml	1 x 100 ml	44 x 104	4200	2901
	Swing out	-	7001315	7001316	200 ml	1 x 50 ml	34 x 100	4200	2901
	Swing out	-	7001315	7001317	100 ml	1 x 25 ml	24 x 100	4200	2901
	Swing out	-	7001315	7001318	240 ml	4 x 15 ml	16.5 x 100	4200	2901
	Swing out	-	7001315	7001319	140 ml	5 x 7 ml	12 x 100	4200	2901
	Swing out	-	7001315	7001320	100 ml	5 x 5 ml	12.5 x 75	4200	2901
	Swing out	-	7001315	7001321	64 ml	4 VAC*	13 x 75	4200	2901
	Swing out	-	7001315	7001322	undrilled	-	-	4200	2901

^{*}VAC(r): Blood tube with cap (bung).

HEAD AND BUCKET REDUCERS

Reducers for head. Part No. 7001711: **7001707** For 1 tube 50 ml (29 Ø x 103 mm). **7001708** For 1 tube 25 ml (24 Ø x 100 mm).

7001709 For 1 tube 30 ml (25 Ø x 92 mm). **7001710** For 1 tube 15 ml (16,5 Ø x 100 mm).

7001780 For 1 tube 10 ml (16 Ø x 80 mm). 7001781 For 4 tubes 7 ml (12 Ø x 100 mm). **7001782** For 4 tubes 5 ml (12,5 Ø x 75 mm). Reducers for bucket. Part No. 7001700: $\textcolor{red}{\textbf{7001779}}$ for 1 tube 85 ml (38 Ø x 112 mm).

^{**}Sealed bucket part no. 7001315 accommodates a 100 ml glass tube without the need of an adapter. Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.



High speed centrifuges with microprocessor control "Centronic BLT" and "Centrofriger BLT"



Versatile. Large range of rotors and adapters. Maintenance free induction motor drive.



Model "Centronic BLT"

Control panel features (see page 115).

External case of DUR-AL alloy with ABS plastic coated front AISI 304 stainless



Model "Centrofriger BLT" refrigerated, temperature range -12 °C to 40 °C.

MODEL

Part No.	Max. Volume ml	Max. number of tubes	Height / Width / Depht (exterior) cm	Power W	Weight Kg
7001768	400	4 x 100 ml	34 41 54	440	22

MODEL

Part No.	Max. Volume ml	Max. number of tubes	Height / Width / Depht (exterior) cm	Power W	Weight Kg
7001769	400	4 x 100 ml	34 41 71	960	39

Note: The rotor being used limits the actual temperature in the chamber due to the speed and ambient temperature. Internal temperature indicated on the LCD display.

ACCESSORIES


Angle rotors with lid, made of anodised hardened aluminium with identification system. All rotors are autoclavable.



1001020		1001323	7001071	7001072	700107	J	7001400 """	outoutator rouder
Part No. head	Rotor type	Max. volume	Tubes capacity	Ø x Length mm	CENTRO Max. R.p.m	INIC BL-II Max. R.c.f. (xg)	CENTROFI Max. R.p.m	RIGER BL-II Max. R.c.f. (xg)
7001670	Angle+lid	270 ml	18 x 15 ml	16.5 x 100	4200	2669	4200	2669
7001706	Angle+lid	400 ml	8 x 50 ml	34 x 100	6000	3920	6000	3920
7001325	Angle+lid	500 ml	6 x 50 ml	29 x 103	6000	18844	6000	24204
7001324	Angle+lid	180 ml	6 x 30 ml	25 x 92	13500	15890	15300	20410
7001323	Angle+lid	100 ml	10 x 10 ml	16 x 80	13500	18335	15300	23550
7001671	Angle+lid	40.5 ml	18x0.25/0.5 ml + 18x1.5/2.2 ml	-	13500	15992	15300	20541
7001672	Angle+lid	60 ml	40 x 0.25/0.5 ml	-	13500	13298	15000	20120
7001673	Angle+lid	45 ml	30 x 1.5/2.2 ml	-	13500	19048	15000	23716
7001408	Angle+lid	Microhaematocrit	24	1.4 x 75	12000	15935	12000	15935
7001702	Angle PCR	6,4 ml	32 x 0.2 ml	-	12000	12500	12000	12500
Note: All buck	ets and adapters o	come without tubes un	lless stated, see page 132 for tul	bes and accessories.				

ACCESSORIES

Swing out rotor with interchangeable sealed buckets and adapters, each with identification code. Autoclavable.























7001407



Sealed bucket with hermetically sealed lid.

















7001315	7001316	7001317	7001318	7001319	7001320
Rotor	Part No.	Sealed		Part No.	To

Part No. head	Rotor Type	Part No. bucket	Sealed bucket	Part No. adapter	Totel max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	Swing out	7001700	-	-	400 ml	1 x 100 ml	44 x 104	5000	4137
	Swing out	7000793	-	-	200 ml	1 x 50 ml	34 x 100	5000	4137
	Swing out	7000794	-	-	100 ml	1 x 25 ml	24 x 100	5000	4137
	Swing out	7000795	-	-	240 ml	4 x 15 ml	16.5 x 100	5000	4137
	Swing out	7000796	-	-	140 ml	5 x 7 ml	12 x 100	5000	4137
	Swing out	7000797	-	-	100 ml	5 x 5 ml	12.5 x 75	5000	3318
	Swing out	7000930	-	-	240 ml	4 x 15 ml Falcon	16.5 x 120	5000	4109
	Swing out	7000931	-	-	200 ml	1 x 50 ml Falcon	30 x 117	5000	4109
	Swing out	7001073	-	-	16 VAC	4 VAC*	16 x 100	5000	4137
7001407	Swing out	7001074	-	-	16 VAC	4 VAC*	13 x 75	5000	3318
	Swing out	7000798	-	-	undrilled	-	-	5000	4137
	Swing out	-	7001315	100 ml glass tube**	400 ml	1 x 100 ml	44 x 104	5000	4137
	Swing out	-	7001315	7001316	200 ml	1 x 50 ml	34 x 100	5000	4137
	Swing out	-	7001315	7001317	100 ml	1 x 25 ml	24 x 100	5000	4137
	Swing out	-	7001315	7001318	240 ml	4 x 15 ml	16.5 x 100	5000	4137
	Swing out	-	7001315	7001319	140 ml	5 x 7 ml	12 x 100	5000	4137
	Swing out	-	7001315	7001320	100 ml	5 x 5 ml	12.5 x 75	5000	4137
	Swing out	-	7001315	7001321	64 ml	4 VAC*	13 x 75	5000	4137
	Swing out	-	7001315	7001322	undrilled	-	-	5000	4137

^{*}VAC(r): Blood tube with cap.

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.

HEAD AND BUCKET REDUCERS


Reducers for head. Part No. 7001706: 7001707 For 1 tube 50 ml (29 Ø x 103 mm). **7001708** For 1 tube 25 ml (24 Ø x 100 mm). **7001709** For 1 tube 30 ml (25 Ø x 92 mm). **7001710** For 1 tube 15 ml (16,5 Ø x 100 mm). **7001780** For 1 tube 10 ml (16 Ø x 80 mm). 7001781 For 4 tube 7 ml (12 Ø x 100 mm). **7001782** For 4 tube 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001325:

 For 1 tube 25 ml (24 Ø x 100 mm). For 1 tube 30 ml (25 Ø x 92 mm). For 1 tube 15 ml (16,5 Ø x 100 mm). For 1 tube 10 ml (16 Ø x 80 mm). For 2 tubes 7 ml (12 Ø x 100 mm). For 2 tubes 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001324: For 1 tube 15 ml (16,5 Ø x 100 mm). For 1 tube 10 ml (16 Ø x 80 mm).

 For 1 tube 7 ml (12 Ø x 100 mm). For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001670: For 1 tube 7 ml (12 Ø x 100 mm). For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001323: For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for bucket. Part No. 7001700: For 1 tube 85 ml (38 Ø x 112 mm).

^{**}Sealed bucket part no. 7001315 accommodates a 100 ml glass tube without the need of an adapter. Micro-haematocrit rotor comes with a box of 100 tubes and calculus graph.

ACCESSORIES



^{*}VAC: Blood tube with stopper.

Swing out

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.



7000609

Microprocessor controlled centrifuges "Meditronic BLT", "Medifriger BLT", "Macrotronic BLT" and "Macrofriger BLT"



2702

4200

COMMON FEATURES

Induction drive maintenance free motor.

Low noise level 50-60dBa.

Robust all metal construction:

- -External case made from DUR-Al alloy.
- -Internal bowl and top plate all stainless steel
- -Internal steel safety chamber within the case.

Electronic circuits, microprocessor digital control of all parameters and functions: speed, acceleration, brake, RCF, temperature, timer, 10 program storage, automatic rotor recognition and alarm system.



Rotor identification head, when the centrifuge lid is closed the rotor communicates all the data of the functions of it. This prevents, before use, any risk of failure of using the wrong rotor in a program or exceeding the rotors running parameters. The advantage of this system is that we can add additional rotors to the range without changing the centrifuge model.

All Rotors, buckets and adapters are identified with an inalterable engraved laser code.

Air circulation system, for centrifuges without refrigeration to limit any excessive rise in temperature within the centrifuge, all air circulated within the chamber is channelled out of the back and is not re-circulated, this means that the centrifuges internal temperature will remain low and constant.

Hermetically sealed compressor, mounted on anti-vibration mounts with re-circulating evaporator around the centrifuge chamber. (only applies to refrigerated centrifuges).

CONTROL PANEL

Touch screen



6 x 96 microtiter

- 1. Enter parameters key.
- 2. Speed analog indicator .
- 3. RPM reading key.
- 4. RCF reading key.
- 5. Acceleration ramp programmed.
- **6.** Temperature programmed. (refrigerated models only).
- 7. Reading the real time temperatu-
- re. (refrigerated models only).

- 8. Rotor model in use.
- 9. Door opening key.
- 10. RPM or RCF programmed.
- 11. RPM or RCF real-time reading
- 12. Start process key.
- Stop process key.
- 4. Programmed braking ramp.
- Time set-up in process.
- 6 Elapsed time reading.

^{**} Bucket part nos. 7000731 and 7000732 accommodate 1 x 200 ml polypropylene bottle with cap without the need of an adapter.



High speed centrifuges with microprocessor control "Meditronic-BLT" and "Medifriger-BLT"



Capacity up to 800 ml. Maintenance free induction motor drive.



Model "Meditronic-BLT"

Model "Medifriger-BLT" refrigerated, temperature range -12 °C to 40 °C.

MODEL

Part No.	Max. Volume	Max. Number	Height / Width / Depth	Power	Weight
	ml	of tubes	(exterior) cm	W	Kg
7001476	800	4 x 200 ml	38 46 52	520	44

MODEL

Part No.	Max. Volume	Max. Number	Height / Width / Depth	Power	Weight
	ml	of tubes	(exterior) cm	W	Kg
7001477	800	4 x 200 ml	38 68 52	1230	70

Note: The rotor being used limits the actual temperature in the chamber due to the speed and ambient temperature. Internal temperature indicated on the LCD display.

ACCESSORIES

Angle rotors with lid, made of anodised hardened aluminium with identification system. All rotors are autoclavable.



70014	149	7001448		7001447	7001446		7001445	70	001444
Part No. head	Rotor type	Max. volume	C	Tube apacity	Ø x Length mm	MEDI Max. R.p.m	TRONIC BL-S Max. R.c.f. (xg)	MEDIFR Max. R.p.m	IGER BL-S Max. R.c.f. (xg)
7000741	Angle+lid	600 ml	6	x 100 ml	44 x 104	4500	3170	4500	3170
7000742	Angle+lid	400 ml	8	x 50 ml	34 x 100	4600	3075	4600	3075
7000743	Angle+lid	360 ml	24	x 15 ml	16.5 x 100	4800	3761	4800	3761
7001449	Angle+lid	510 ml	6	x 85 ml	38 x 105	8000	6868	12000	15455
7001448	Angle+lid	400 ml	8	x 50 ml	29 x 103	8000	7802	12000	15938
7001447	Angle+lid	180 ml	6	x 30 ml	25 x 92	9000	7424	14000	17965
7001446	Angle+lid	200 ml	20	x 10 ml	16 x 80	9000	8964	11000	13390
7001445	Angle+lid	100 ml	10	x 10 ml	16 x 80	9000	7293	15000	20124
7001444	Angle+lid	45 ml	30	x 1.5/2.2 ml	-	11000	12619	15000	23470
Note: All bucket	Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.								

ACCESSORIES

Swing out rotor with interchangeable sealed buckets and adapters, each with identification code. Autoclavable.

Swing out rotor Sealed bucket with hermetically sealed lid. Open bucket Adapters

7000732 7000731 5715202 7000733 7000734 7000935 7000

7009079 1000732 1000731 1000734 1000734 1000935 1000935 100093

7000738

7001648

7000739

Part No. head	Rotor type	Bucket open	Sealed bucket	Part No. adapter	Totel max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	Swing out	7000731	7000732	Bottle 200 ml 5715202**	800 ml	1 x 200 ml	56.6 x 104	5000	4137
	Swing out	7000731	7000732	7000733	400 ml	1 x 100 ml	44 x 104	5000	4137
	Swing out	7000731	7000732	7000734	200 ml	1 x 50 ml	34 x 100	5000	4137
	Swing out	7000731	7000732	7000935	200 ml	1 x 50 ml Falcon	30 x 117	5000	4137
	Swing out	7000731	7000732	7001645	240 ml	2 x 30 ml	25 x 92	5000	4137
	Swing out	7000731	7000732	7000735	200 ml	2 x 25 ml	24 x 100	5000	4137
	Swing out	7000731	7000732	7000736	420 ml	7 x 15 ml and VAC*	16.5 x 100	5000	4137
7009079	Swing out	7000731	7000732	7000934	180 ml	3 x 15 ml Falcon	16.5 x 120	5000	4137
	Swing out	7000731	7000732	7001646	-	4 x 17 Ø x 110 mm	-	5000	4137
	Swing out	7000731	7000732	7001647	240 ml	6 x 10 ml	17 x 110	5000	4137
	Swing out	7000731	7000732	7000737	336 ml	12 x 7 ml	12 x 100	5000	4137
	Swing out	7000731	7000732	7000738	240 ml	12 x 5 ml	12.5 x 75	5000	4137
	Swing out	7000731	7000732	7001648	140 ml	7 x 5 ml and VAC*	13 x 75	5000	4137
	Swing out	7000731	7000732	7000739	144 ml	12 x 3 ml	11 x 100	5000	4137
	Swing out	7000731	7000732	7000740	-	undrilled	-	5000	4137
7000609	Swing out	-	-	-	-	6 x 96 microtiter	-	4200	2702

^{*}VAC: Blood tube with stopper.

7000736

7000934

7001646

7001647

7000737

HEAD REDUCERS Reducers for head Part No. 7000741:

HEUUCEIS IOI HEUU. I AIT NO. 1000141.
7001770 For 1 tube 85 ml (38 Ø x 105 mm).
7001771 For 1 tube 50 ml (34 Ø x 100 mm).
7001772 For 1 tube 50 ml (29 Ø x 103 mm).
7001773 For 1 tube 25 ml (24 Ø x 100 mm).
7001774 For 1 tube 30 ml (25 Ø x 100 mm).
7001775 For 3 tubes 15 ml (16,5 Ø x 100 mm).
7001776 For 3 tubes 10 ml (16 Ø x 80 mm).
7001777 For 7 tubes 7 ml (12 Ø x 100 mm).
7001778 For 7 tubes 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7000742:
7001707 For 1 tube 50 ml (29 Ø x 103 mm).
7001708 For 1 tube 25 ml (24 Ø x 100 mm).
7001709 For 1 tube 30 ml (25 Ø x 92 mm).
7001710 For 1 tube 15 ml (16,5 Ø x 100 mm).

7001780 For 1 tube 10 ml (16 Ø x 80 mm).
7001781 For 4 tubes 7 ml (12 Ø x 100 mm).
7001782 For 4 tubes 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7001449:
7001783 For 1 tube 50 ml (34 Ø x 100 mm).
7001784 For 1 tube 50 ml (29 Ø x 103 mm).
7001785 For 1 tube 25 ml (24 Ø x 100 mm).
7001786 For 1 tube 30 ml (25 Ø x 92 mm).
7001787 For 2 tubes 15 ml (16,5 Ø x 100 mm).
7001788 For 2 tubes 10 ml (16 Ø x 100 mm).
7001789 For 5 tubes 7 ml (12 Ø x 100 mm).
7001790 For 5 tubes 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7001448:
7001791 For 1 tube 25 ml (24 Ø x 100 mm).
7001792 For 1 tube 30 ml (25 Ø x 100 mm).

7001793 For 1 tube 15 ml (16,5 Ø x 100 mm).
7001794 For 1 tube 10 ml (16 Ø x 80 mm).
7001795 For 2 tubes 7 ml (12 Ø x 100 mm).
7001796 For 2 tubes 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7001447:
7001797 For 1 tube 15 ml (16,5 Ø x 100 mm).
7001798 For 1 tube 10 ml (16 Ø x 80 mm).
7001799 For 1 tube 7 ml (12 Ø x 100 mm).
7001800 For 1 tube 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7001445/7001446:
7001801 For 1 tube 5 ml (12,5 Ø x 75 mm).
Reducers for head. Part No. 7000743:
7000700 For 1 tube 7 ml (12 Ø x 100 mm).
1001557 For 1 tube 5 ml (12,5 Ø x 75 mm).

7000609

7000740

^{**} Bucket part nos. 7000731 and 7000732 accommodate 1 x 200 ml polypropylene bottle with cap without the need of an adapter. Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.



Microprocessor controlled Centrifuge "Macrotronic-BLT" and "Macrofriger-BLT"



Multi-capability. Capacity from 1.5/2.2 ml vials up to tubes of 500 ml and 500ml blood bags.

Maintenance free induction motor drive.

Control panel features (see page 123). External case made of DUR-AL alloy with an AISI 304 stainless steel bowl.





Model "Macrofriger-BLT" refrigerated, temperature range -12 °C to 40 °C.

Model "Macrotronic-BLT"

MUDEL					
Part No.	Max. volume ml	Max. No. of bags	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7001/185	2000	/I v 500 ml	45 BO BB	720	76

MODEL

Part No.	Max. Max. No.		Height / Width / Depth	Power	Weight
	volume ml of bags		(exterior) cm	W	Kg
7001486	2000	4 x 500 ml	97 60 66	1230	156

Note: The rotor being used limits the actual temperature in the chamber due to the speed and ambient temperature. Internal temperature Indicated on the LCD display.



Part No. head	Rotor type	Max. volume	Tube capacity	Ø x Length mm	MACRO Max. R.p.m	TRONIC BL Max. R.c.f. (xg)	MACROI Max. R.p.m	FRIGER BL Max. R.c.f. (xg)
7000561	Anale+lid	600 ml	12 x 50 ml	34 x 100	4530	3785	4530	3785
7000560	Angle+lid	840 ml	56 x 15 ml	16.5 x 100	4030	3268	4030	3268
7001650	Angle+lid	1500 ml	6 x 250 ml	62 x 122	8000	10373	9000	13129
7001651	Angle+lid	510 ml	6 x 85 ml	38 x 105	10000	10731	12000	15452
7001652	Angle+lid	400 ml	8 x 50 ml	29 x 103	10000	11066	12000	15935
7001653	Angle+lid	360 ml	12 x 30 ml	25 x 92	10000	12296	12000	17706
7001654	Angle+lid	200 ml	20 x 10 ml	16 x 80	10000	11066	12000	15935

Swing out rotor with interchangeable sealed buckets and adapters, each with identification code. Autoclavable. Open bucket Swing out rotor Sealed bucket with hermetically sealed lid. 7000745 7000746 7000756 **Adapters** 7001387 7001657 7000747 7001658 7000748 7000939 7001659 7001660 7000749 7000750 7000938 7001661 7001662 7000751 7000752 7000753 7000754 7000755 Part No. Rotor Bucket Sealed Part No. Totel max. Volume Bucket/Adapter Ø x Length Max. Max. bucket R.c.f. (xg) head adapter (4 adapters) capacity mm R.p.m. 4615 Swing out 7000746 7000756 7001387 800 ml 1 x 200 ml 56.6 x 104 5000 Swing out 500 ml 1 x 125 ml 5000 4615 7000746 7000756 7001657 51 x 115 Swing out 7000746 7000756 7000747 400 ml 1 x 100 ml 44 x 104 5000 4615 x 105 Swing out 7000746 7000756 7001658 320 ml 1 x 80 / 85 ml 38 5000 4615 Swing out 7000746 7000756 7000748 400 ml 2 x 50 ml 34 x 100 5000 4615 Swing out 7000746 7000756 7000939 400 ml 2 x 50 ml Falcon 30 x 117 5000 4615 Swing out 7000746 7000756 7001659 400 ml 2 x 50 ml 29 x 103 5000 4615 **Swing out** 4615 7000746 7000756 7001660 600 ml 5 x 30 ml 25 x 92 5000 Swing out 7000746 7000756 7000749 500 ml 5 x 25 ml 24 x 100 5000 4615 7000745 Swing out 7000746 7000756 7000750 720 ml 12 x 15 ml 16.5 x 100 5000 4615 12 VAC* Swing out 7000746 7000756 7000750 720 ml 16.5 x 100 5000 4615 Swing out 7000746 7000756 7000938 240 ml 4 x 15 ml Falcon 16.5 x 120 5000 4615 Swing out 7000746 7000756 7001661 12 VAC* 13 x 75 5000 4615 Swing out 7000746 7000756 7001662 480 ml 12 x 10 ml 16.8 x 85 5000 4615 Swing out 7000746 7000756 7000751 560 ml 20 x 7 ml 12 x 100 5000 4615 **Swing out** 7000746 7000756 7000752 400 ml 20 x 5 ml 12.5 x 75 5000 4615 Swing out 7000746 240 ml 20 x 3 ml 11 x 100 5000 4615 7000756 7000753 Swing out 7000746 120 ml 20 de 1.5 / 2.2 ml 4615 7000756 7000754 5000 7000746 5000 Swing out 7000756 7000755 undrilled 4615

*VAC: Blood tube with stopper.

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.

ACCESSORIES

Swing out rotor





Open bucket





Adapters















Bucket for 3 x 96 well plates

Part No. head	Rotor type	Bucket open	Part No. adapter	Total max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	swing out	7001451	7001395	2000 ml	1 x 500 ml	70 x 165	3200	2358
	swing out	7001451	7001655	1000 ml	1 x 250 ml	62 x 122	3200	2358
	swing out	7001451	7001396	800 ml	1 x 200 ml	56.6 x 104	3200	2358
	swing out	7001451	7001397	800 ml	2 x 100 ml	44 x 104	3200	2358
7001450	swing out	7001451	7001398	800 ml	4 x 50 ml	34 x 100	3200	2358
7001450	swing out	7001451	7001399	900 ml	9 x 25 ml	24 x 100	3200	2358
	swing out	7001451	7001400	1020 ml	17 x 15 ml	16.5 x 100	3200	2358
	swing out	7001451	7001401	672 ml	24 x 7 ml	12 x 100	3200	2358
	swing out	7001656	_	_	3 plates x 96 microtiter	_	3200	2117

Swing out rotor










158 Height x 58 mm Ø max. for buckets Part No. 7000764





Part No. head	Rotor type	Open Bucket	Total max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	swing out	7000763	1000 ml	1 x 250 ml	56 x 147	3520	3567
7000552	swing out	7001075	4 x 1 ASTM	1 x ASTM cylindrical	36 x 200	3720	3567
	swing out	7000764	4 x 1 ASTM	1 x ASTM pear shaped	58 x 158	3720	3567

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.

HEAD REDUCERS

Reducers for head. Part No. 7001650: For 1 tube 100 ml (44 Ø x 104 mm). 7001804 For 1 tube 85 ml (38 Ø x 105 mm). For 1 tube 50 ml (34 Ø x 100 mm). For 2 tubes 50 ml (29 Ø x 103 mm). For 1 tube 50 ml (30 Ø x 117 mm). 7001808 For 3 tubes 30 ml (25 Ø x 92 mm). For 7 tubes 15 ml (16,5 Ø x 100 mm). For 7 tubes 10 ml (16 Ø x 100 mm). For 12 tubes 7 ml (12 Ø x 100 mm). For 12 tubes 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7001651:

 For 1 tube 50 ml (34 Ø x 100 mm). For 1 tube 50 ml (29 Ø x 103 mm). For 1 tube 25 ml (24 Ø x 100 mm). For 1 tube 30 ml (25 Ø x 92 mm). For 2 tubes 15 ml (16,5 Ø x 100 mm). For 2 tubes 10 ml (16 Ø x 100 mm). For 5 tubes 7 ml (12 Ø x 100 mm). 7001790 For 5 tubes 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001652: For 1 tube 25 ml (24 Ø x 100 mm). For 1 tube 30 ml (25 Ø x 100 mm).

 For 1 tube 15 ml (16.5 Ø x 100 mm). For 1 tube 10 ml (16 Ø x 80 mm). For 2 tubes 7 ml (12 Ø x 100 mm). For 2 tubes 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7001653: For 1 tube 15 ml (16,5 Ø x 100 mm). For 1 tube 10 ml (16 Ø x 80 mm). For 1 tube 7 ml (12 Ø x 100 mm). For 1 tube 5 ml (12,5 Ø x 75 mm). Reducers for head. Part No. 7000560: For 1 tube 7 ml (12 Ø x 100 mm). For 1 tube 5 ml (12,5 Ø x 75 mm).



Refrigerated high capacity microprocessor controlled centrifuge suitable for tubes, flasks and blood bags "Macrofriger-BLT-Blood"




CONTROLLABLE TEMPERATURES FROM 0 TO 35 °C.

CAPACITY FOR FLASKS UP TO 1000 ml and especially for triple and quad blood bags from 500-800 ml. Maximum capacity 12 bags.




Original J.P. Selecta design. With maintenance free induction drive motor



APPLICATIONS


Specially designed for blood bank applications, where repeatability, precision and reliability of: speed, time, temperature, acceleration and deceleration are significant factors in obtaining perfect separations of the distinct blood fractions.



Serology investigatios.

Molecular biology.

Serum separation, pharmaceutical applications.



Cross-section of the motor, rotor drive mechanism, showing the windbreak protector with lid that reduces the air resistance and maintains a regulated temperature.

FEATURES

Case: internal structure of tubular steel with lateral Dur-al supports.

Internal Stainless steel bowl .

Aluminium windbreak cover with detector and alarm sensor that prevents the centrifuge operation without it.

Safety shield: 8 mm thick steel.

Swing-out rotor, buckets and adapters made of high-density material, engraved identification unalterable.

Suspended motor with silent-block and shock absorbers that guarantee smooth, quiet operation without vibrations.

Drive force: equipped with an induction drive maintenance free and highly efficient pollution free motor.

Semi-hermetically sealed, compressor pollution free with no hydrocarbons or fluorocarbons.

Mounted on an anti vibration adjustable rubber feet.

Microprocessor digital control of all functions: speed, acceleration, brake, R.C.F., temperature, time, 20 program memory, rotor identification, operation fault alarm. Digital indicator of the rotor and blood bag buckets operating hours.

SAFETY:

SAFETY, CONFORMS TO UNE-EN 61010-2-020.
PROTECTIVE INTERNAL CYLINDER BETWEEN THE CASE AND BOWL.
THE LID CANNOT BE OPENED WHILE THE ROTOR IS SPINNING.
THE MOTOR CANNOT SPIN WHILE THE LID IS OPEN.

IMBALANCE DETECTOR, AUTOMATICALLY STOPS IF THERE IS ANY IMBALANCE. DOUBLE ELECTROMAGNETIC LID LOCK.

IF THERE IS A POWER FAILURE, A FAIL SAFE SYSTEM CAN BE USED TO OPEN THE LID.

EMERGENCY STOP BUTTON.

CONTROL PANEL

Touch screen



- 1. Enter parameters key.
- 2. Speed analog indicator .
- 3. RPM reading key.
- 4. RCF reading key.
- **5.** Acceleration ramp programmed.
- **6.** Temperature programmed. (refrigerated models only).
- **7.** Reading the real time temperature. (refrigerated models only).
- 8. Rotor model in use.
- 9. Door opening key.
- 10. RPM or RCF programmed.
- 11. RPM or RCF real-time reading
- 12. Start process key.
- 13. Stop process key.
- 14. Programmed braking ramp.
- 15. Time set-up in process.
- 16 Elapsed time reading.





Image showing the rotor and double buckets for quad and triple blood bags of 500 and 800 ml.

TECHNICAL SPECIFICATION


Maximum volume:9600 ml	
Maximum bucket volume:	
Speed:)70 r.p.m.
Accuracy:± 5 r.p.m.	
Resolution r.p.m.:	
R.C.F.:	′O
Imbalance:±5 g	
Accuracy:4 g	
Temperature range:	35 °C
Temperature stability:	
Resolution:1°C	
Temperature range:	<u>.</u>
Centrifuge radius:	
Acceleration profile:	
Brake profile:	lratic .
Voltage:	
Power:	1

MODEL

Part No.	Max. Volume ml	Max. Capacity		Height / Width / Depth (exterior) cm		Power W	Weight Kg
7001490	9600	6 x 1600 ml	105	96	96	5500	420

ACCESSORIES

Swing out rotor 6 place, treated



Part No. 7001391

Vasos v reductores intercambiables con señal de identificación. Autoclavables.



Buckets with two compartments for blood bags of 500 to 800 ml. For triple and quad bags Part No.**7001392**



Single bucket. Capacity: 1 blood bag of 500 to 800 ml or 1 flask of 1 litre or with tube adapters. Part No. 7001393

Bags holder set for double cup 7001392 Part No. 0928292

Adapters for single bucket Part No. 7001393.

Bottle

















7001398

7001399

7001400

Double Bucket	Single Bucket	Part No. adapter	Total max. Volume (6 adapters)	Capacity: bags, tubes by adapter	Ø x Length mm	R.p.m. Min./ Max.	R.c.f. (xg) Min./Max.
7001392	-	-	9600 ml	2 bags 500 or 800 ml	-	300/4070	29/5370
-		Bottle 1000 ml 1001394	6000 ml	1 bottle 1000 ml or 1 bag	98 x 175	300/4070	29/5370
-	7001393	7001395	3000 ml	1 bottle 500 ml	70 x 165	300/4070	29/5370
-	Capacity	7001396	1200 ml	1 bottle 200 ml	56.6 x 104	300/4070	29/5370
-	without adapter:	7001397	1200 ml	2 tubes 100 ml	44 x 104	300/4070	29/5370
-	1 bag of 500 to 800 ml	7001398	1200 ml	4 tubes 50 ml	34 x 100	300/4070	29/5370
-	or 1x bottle	7001399	1350 ml	9 tubes 25 ml	24 x 100	300/4070	29/5370
-	1000 ml	7001400	1530 ml	17 tubes 15 ml	16.5 x 100	300/4070	29/5370
-	1000 1111	7001401	1008 ml	24 tubes 7 ml	12 x 100	300/4070	29/5370

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories. Other adapters can be manufactured to the customers specification.

COMPLEMENTS



Blood bank refrigerated cabinets "Blood Bank" A, B, C and D

FORCED AIR CIRCULATION.

FIXED TEMPERATURE +4 °C. STABILITY ±1.5 °C. RESOLUTION 1 °C. SEVEN DAYS GRAPHIC TEMPERATURE REGISTER. CAPACITIES FROM 90 TO 480 BLOOD BAGS.

FEATURES

See page. 81





MODELS	Part No.	Capacity litres	Capacity bags		t / Width interior)	r / Depht cm		/ Width xterior)	/ Depht cm	Power W	Power HP	Number of shelves	Weight Kg
Blood Bank A	2101497	319	110	139	48	45	198	60	64	360	1/4	5	98
Blood Bank B	2101498	442	150	126	58	55	198	70	71	360	1/4	5	111
Blood Bank C	2101499	600	200	136	58	69	207	70	82	350	3/8	5	115
Blood Bank D	2101500	1300	480	136	128	69	208	140	82	1100	3/8	10	190
Comes complete w	ith steel shelves	, PVC laminated											

COMMON ACCESSORIES. Tubes of glass and plastic.

Plastic tubes and flasks with screw cap lids.

L							
	Volume ml	Dimensions mm	Туре	Part No.	Material	Max. (xg) range	Temperature °C
			Comecta	5752399	Polypropylene	50000	4 121
	10	16 Ø x 80	Herolab	5752421	Polypropylene	57000	-180 14 5
			Herolab	5752401	Polycarbonate	57000	-150 121
			Comecta	5752398	Polypropylene	50000	4 121
	30	25 Ø x 92	Herolab	5752955	Polypropylene	64000	-180 14 5
			Herolab	5752951	Polycarbonate	64000	-150 121
			Comecta	5752397	Polypropylene	50000	4 121
	50	29 Ø x 103	Herolab	5753205	Polypropylene	69000	-180 14 5
			Herolab	5753201	Polycarbonate	69000	-150 121
			Comecta	5752396	Polypropylene	40000	4 121
	80	38 Ø x 105	Herolab	5753243	Polypropylene	50000	-180 14 5
			Herolab	5753241	Polycarbonate	50000	-150 121
	200	56.6 Ø x 104	Herolab	5715202	Polypropylene	10000	-180 14 5
	200	30.0 Ø X 104	Herolab	5715203	Polycarbonate	10000	-150 121
			Comecta	5752395	Polypropylene	13200	4 121
	250	62 Ø x 122	Herolab	5753411	Polypropylene	33000	-180 14 5
			Herolab	5753401	Polycarbonate	33000	-150 121
	500	70 Ø x 165	Herolab	5753512	Polypropylene	14000	-180 14 5
	300	7U W X 100	Herolab	5753501	Polycarbonate	14000	-150 121
	1000	98 Ø x 175	Herolab	5752571	Polypropylene	4936	-180 14 5
	1000	30 N X 1/3	Herolab	5752561	Polycarbonate	4936	-150 121
100							



Type HEROLAB
polycarbonate tubes
transparent Polypropylene
tubes translucent.
Autoclavable.

Type COMECTA polypropylene translucent. Autoclavable.





Glass Tubes.		
Volume	Dimensions mm	Part No.
Micro-haematocrit*	1.4 Ø x 75	7000910
7 ml	12 Ø x 100	7000773
15 ml	16.5 Ø x 100	7000772
25 ml	24 Ø x 100	7000771
50 ml	34 Ø x 100	7000770
100 ml	44 Ø x 104	7000769
250 ml	56 Ø x 147	7000768
ASTM	Conical	7001088
ASTM	Pear shape	7000767
*Box of 100 tubes		

SPARE									
Buckets for angle roto	Buckets for angle rotors.								
Suitable for model	Part No. 7001403 15 ml plastic	Part No. 7001500 15 ml metal	Part No. 7005000 50 ml metal	Part No. 700 1000 100 ml metal					
Cencom II Part No. 7002240	Yes								
Centro-8 Part No. 7002356		Yes							
Centromix II-BL		Rotor							
Part No. 7002555		7000550							
Meditronic BL-S		Rotor	Rotor	Rotor					
Part No. 7001376		7000743	7000742	7000741					
Medifriger BL-S		Rotor	Rotor	Rotor					
Part No. 7001377		7000743	7000742	7000741					
Macrotronic BL		Rotor	Rotor						
Part No. 7001385		7000560	7000561						
Macrofriger BL		Rotor	Rotor						
Part No. 7001386		7000560	7000561						

ALL THE ROTORS, BUCKETS AND ADAPTERS COME WITH LASER ETCHED IDENTIFICATION.





Centrifuge kit for separating plasma rich platelets from whole blood

MADE UP OF:





Centrifuge "CENTRO-4-BL"

Swing-out rotor (not an angle rotor), that performs the correct separation of the distinct blood fractions.

Capacity: 4 tubes, 75 x 13 mm VAC tubes.

Electronic control with digital display. Maximum speed of 5000 r.p.m. Maximum time 15 minutes or continuous operation. Fixed acceleration and deceleration rate.

Strong moulded external case of ABS plastic, extremely light weight, (6 Kg). Part No. 7002560





Metallic thermostat dry block "BIO-BATH"

Small size with a fixed temperature of 37 °C with check thermometer, 8 tube positions for 75x13 VAC tubes. Note: Reduces the waiting time for the patient.

MODEL

Part No.	Height / Width / Depht	Stability	Power	Weight
	(exterior) cm	°C	W	Kg
7001561	8.8 9.6 7	0.5	8	0.5





Part No. 1001562 Support rack made from stainless steel for 3 pipettes

Part No. 1001563 Tube rack made from stainless steel for 24 tubes (6x4) of 75x13 mm VAC tubes.
Part No. 1001574 Flask containing 10% calcium chloride, quantity 30 units. 10 ml. per flask.

Part No. 1001568 Strap bracelets of pure silicon "Esmarch". 3 units.





1 Part No. 1001558 Pipette tips, sterile with filters packed as individual tips, size 5-200 µl, and quantity 250 units.*

1 Part No. 1001559 Pipette tips, sterile with filters packed as individual tips, size 100-1000 ul, and quantity 250 units.*

2 Part No. 1001565 Tubes VE Sodium citrate 75x13mm VAC tubes, blue 4 ml sterile tubes. Box of 100 units.

2 Part No. 1001566 Tubes 75x13 mm sterile,100 units in bags of 4 units.

3 Part No. 1001567 Micro extractor for the extraction of blood by vacuum. Box of 50 units. Complete with two adapters.

Note: * All above components are individually sterilized.

If aseptic conditions are required a laminar flow cabinet should be used for sample manipulation.



Part No. **5900499** Timer alarm bench top "RELAVIS" 0-60 minutes.



J.P. SELECTA has designed an effective small centrifuge and thermostat dry block that is suitable for easy transportation for consultants, emergency services, reception labs and clinical trials etc. To maintain the complete kit, a robust **metallic case** with moulded interior is provided. The complete kit comprises of:

1 centrifuge "Centro-4", 1 metallic thermostat dry block "Bio-Bath", mains cables, 3 fixed volume pipettes "Dispen-Fix", 1 rack 24 tubes, 1 timer "Relavis", and 2 boxes for consumables, tubes, tips, micropercusors, strips and flasks. In the lid are clip positions for easy access for the pipette rack and 3 pipettes. Total weight with all items: 11.5Kg.

Part No. 1001570 Metallic case with moulded compartments Armaflex, two perspex boxes for tubes and consumable Clips for pipette support.





Poupinel sterilizers
Drying and sterilization ovens
Universal precision ovens
High temperature ovens
Vacuum ovens
Desiccators
Bacteriological incubators
Cooled low temperature incubator
Incubator chamber
CO2 incubator
Cooled incubators
Electric muffle furnaces

137 page 138 to 150 pages 143 to 144 pages 147 page 148 to 149 pages 150 page 151 to 154 pages pages 155 to 156 157 page 158 page 159 to 161 pages 162 to 164 pages



LABORATORY OVENS AND INCUBATORS

Leading edge technology. Maximum Precision.



Drying and sterilization ovens:

Drying: For all drying processes of diverse laboratory material or glass material in general, printed circuits, granule and powder, etc. **Sterilization:** They guarantee microorganisms destruction, either pathogen ones or not, which will be over or inside the material. They allow sterilization of powder and non-volatile viscose substances. For a good sterilization, a temperature between 160 °C and 180 °C and 2 hours exposition is usually required.

Vacuum drying ovens:

They are developed for applications of thermal and drying treatments of heat-sensitive products.

Bacteriological culture ovens:

For microorganisms or culture incubation in clinical diagnosis, in sanitary or nutritious industry. The samples are preserved at a determined temperature and period of time.

Cooled low temperature ovens:

For microorganisms or culture incubation, in clinical diagnosis, in sanitary or nutritious industry. The samples are preserved at a determined temperature and period of time.

Anaerobic cell and tissue cultures for CO2 ovens:

Essential element in laboratories for research, cell biology, molecular biology, different cancer sorts and general pharmaceutical laboratories.

Conforms to the international directives for safety and precision. MODEL RANGE:

- Drying and sterilization.
- Universal: Programmed for cultures and sterilization.
 - Vacuum drying.
 - Bacteriological cultures.
 - Low temperature- High Precision Peltier systems.
 - CO₂ Incubators.
 - Precise refrigerators and cooled incubators.
- More than 70 models with capacities from 19 to 720 litres.
- Controllable temperatures from -10 to 250 °C and 400 °C.
- Analogue or digital control through a microprocessor for temperature and time.
 - Wide range of accessories for varying applications.



OVENS, INCUBATORS AND FURNACES Summary table of the different models



							ŲUIGK (JVERVIE	/V	
MODEL RANGE Sterilizers poupinel	MODELS	CONTROL	2.5 litres	19 litres	CAPACITY			Safety	RS-232	US
100 250 °C	DRYTIME	ANALOGUE	2000911	-				YES	-	-
60 250 °C	DRYTERM	ANALOGUE	-	2000787				YES	-	
GLASS DRYING			126 litres							
40 170 °C	DRYGLASS	ANALOGUE	2000381					YES	-	
CONVECTION NATURAL			19 litres	36 litres	52 litres	80 litres	150 litres			
40 250 °C	CONTERM	ANALOGUE	2000208	2000209	2000200	2000210	2000201	YES		
Ambient+5 250 °C	DIGITHEAT-TFT	TFT Touch screen µ	2001251	2001252	2001253	2001254	2001255	YES	YES	γ
FORCED AIR, FAN CONVECTION	BENCH TOP	·		33 litres	47 litres	76 litres	145 litres			
Ambient+5 250 °C	DIGITRONIC-TFT	TFT Touch screen µ		2005163	2005165	2005167	2005169	YES	YES	γ
	DIGITRONIC-TFT glass door	•		2005164	2005166	2005168	2005170	YES	YES	γ
FORCED AIR, FAN CONVECTION I		<u> </u>	216 litres	288 litres	400 litres	720 litres				
Ambient+5 250 °C	DRYBIG 230/400V III PHAS	SES DIGITAL µ	2002961	2002971	2003721	2003741		YES	YES	γ
Ambient+5 250 °C	DRYBIG 230V I PHASE	•	2002962	2002972	_	-		YES	YES	γ
HIGH TEMPERATURE			80 litres							
60 400 °C	HIGHTEMP 230/400V III PHA	SES DIGITAL u	2001406					YES	-	
VACUUM OVEN			3 litres	47 litres						
35 200 °C	VACIOTEM T	DIGITAL µ	-	4001489				YES	YES	γ
35 200 °C	VACIOTEM TV	DIGITAL µ	-	4001490				YES	YES	γ
Ambient+5 170 °C	VACUO-TEMP	DIGITAL	4000474					YES	-	
DESICCATOR			55 litres							
			1001403					YES	-	
INCUBATION CHAMBER			110 litres							
Ambient+5 57 °C	BOXCULT	DIGITAL	3000957					YES	-	
INCUBATORS BENCH TOP			19 litres	36 litres	52 litres	80 litres	150 litres			
Ambient+5 80 °C	INCUBAT	ANALOGUE	2000205	2000206	2001615	2000207	2000994	YES	-	
Ambient+5 80 °C	INCUDIGIT-TFT	TFT Touch screen µ	2001256	2001257	2001258	2001259	2001260	YES	YES	γ
INCUBATORS LARGE AND FLOOI	R STANDING	<u> </u>	288 litres	400 litres	720 litres					
Ambient+5 80 °C	INCUBIG-TFT	TFT Touch screen µ	2000238	2000239	2000240			YES	YES	γ
LOW TEMPERATURE CABINETS		·	36 litres	80 litres	150 litres					
5 60 °C	PREBATEM-TFT	TFT Touch screen µ	2000963	2000964	2000965			YES	YES	γ
CO2 INCUBATOR		·			150 litres					
Ambient+5 50 °C	INCUBATOR CO2	DIGITAL µ			4002628			YES	YES	γ
WITH REFRIGERATION			160 litres	319 litres	442 litres	600 litres				
+5 65 ℃	HOTCOLD-S	DIGITAL µ	2101618	-	-	-		YES	-	
0 50 °C	HOTCOLD A-B-C	DIGITAL µ	-	2101502	2101503	2101504		YES	YES	γ
-10 50 °C	HOTCOLD UB - UC	DIGITAL µ		-	2101505	2101506		YES	YES	γ
5 50 °C	HOTCOLD GL	DIGITAL µ	-	-	-	2101507		YES	YES	γ
MUFFLE FURNACE			3 litres	3.6 litres	8 litres	9 litres				
Up to 1150 °C	SELECT-HORN-TFT	TFT Touch screen	-	2000366	_	2000367		YES	YES	γ
Up to 1100 °C	R	TFT Touch screen	2000368	-	2000369	-		YES	-	-
with micronrocessor										



Poupinel dry heat sterilizer "Drytime" ADJUSTABLE TEMPERATURES FROM 100 °C UP TO 250 °C.

ADJUSTABLE TEMPERATURES FROM 100 °C UP TO 250 °C STABILITY: ± 6 °C.

APPLICATIONS

For quick surgical sterilization of diverse instruments surgical odontological, etc.

FEATURES

Heating by shielded elements in the base which provide a rapid temperature rise.

Flap door.

Inner chamber in AISI 304 stainless steel. Removable tank with extraction clamps.

Epoxy-coated outer casing.

SAFETY

Over temperature cut out incorporated. EN.61010 Standard.

CONTROL PANEL

Mains switch.

Mains indicator lamp.

Hydraulic thermostat for temperature control.

Timer 0 to 120 min. with automatic off.

Analogue temperature reading thermometer.





Part No. 00009096

MODEL

Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Power	Weight
	litres	(interior) cm	(exterior) cm	W	Kg
2000911	2.5	5 30 16	17 40 32	430	8



Poupinel dry heat sterilizer "Dryterm"

ADJUSTABLE TEMPERATURES FROM 60 °C UP TO 250 °C. STABILITY: ±10 °C.

APPLICATIONS

For surgical sterilization of diverse instruments surgical odontological,etc.

FEATURES

Heating by shielded elements in the base that provides a rapid rise in temperature.

Flap door.

Inner chamber made of AISI 304 stainless steel, complete with a heater cover, three shelf runners and two perforated shelves 10 mm high.

Epoxy-coated outer casing.

SAFETY

Over temperature cut out incorporated. EN.61010 Standard.

CONTROL PANEL

Hydraulic thermostat temperature control. Locking device for thermostat knob.

Timer 0 to 120 min. with automatic switch off.


Heater "ON" indicator.

Analogue temperature reading thermometer.

MODEL

Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Power	Weight
	litres	(interior) cm	(exterior) cm	W	Kg
2000787	19	25 32 23	37 54 34	770	19







SAFETY:

EN.61012 STANDARD OVER TEMPERATURE SAFETY CUT OUT FITTED.
DIN 12880. STANDARD (CLASS 2 AND 3.1) ADJUSTABLE SAFETY THERMOSTAT FITTED.

FEATURES

Hydraulic thermostat for temperature control.

Air circulation by turbo fan.

Inner chamber made of AISI 304 stainless steel with shelf runners.

Removable tempered glass sliding doors.

Ventilation port for steam.

Epoxy coated external case.

STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

CONTROL PANEL

Dual heating power selector switch.

Mains indicator lamp.

Hydraulic thermostat for temperature control.

Locking system of thermostat knop.

Heater "ON" operation indicator lamp.

Analogue thermometer.

Adjustable over temperature safety thermostat, that cuts off the power if the control thermostat fails, manual reset with "on" indicator lamp.



MODEL

Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Shelf	Power	Weight
	litres	(interior) cm	(exterior) cm	Positions	W	Kg
2000381	126	45 70 40	66 94 54	8	3000	65

ACCESSORIES

Accessories must be factory installed.



Part No.

2000002 Timer switch 0-120 minutes. **2000003 Timer switch** 0-12 hours.

2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes.

SPARES

Part No.

2000081 Shelf guides x 2.

2000091 Shelf.

Each shelf requires 2 guides.



OVENS AND INCUBATORS PREMIER RANGE

MODELS:

- NATURAL AIR CONVECTION. DRYING AND STERILIZATION.
- FAN ASSISTED CIRCULATION. UNIVERSAL APPLICATIONS.
- NATURAL AIR CONVECTION. BACTERIOLOGY AND INCUBATION.

CONTROL: ANALOGUE OR DIGITAL MICROPROCESSOR CONTROL OF TEMPERATURE AND TIME, MODEL DEPENDENT. COMPLIES WITH THE STANDARDS: DIN 50011 - DIN 58945. REQUIRED FOR HEATING, STABILITY AND HOMOGENEITY.



SAFETY:

STANDARD EN.61010. INCORPORATED FIXED OVER TEMPERATURE DEVICE .
STANDARD DIN 12880. (CLASS 2 AND 3.1)SAFETY THERMOSTAT CONTROLLER FITTED.

Leading edge technology





Detailed longitudinal cross section.

COMMON FEATURES

Construction.

- 1. External case treated with a corrosive resistant epoxy coating.
- 2. Internal part: Easy to clean AISI 304 stainless steel double chamber, self adjusting door seal and adjustable shelves and guides.
- **3.** Control panel: independent insulated control panel to facilitate all types of instruments, controls and regulators.
- 4. Adjustable air inlet.
- **5.** Flexible floating door seal, self adjusting that maintains the best possible seal.

Technical Properties.

- **6.** Excellent thermal qualities of the insulation has the optimum performance according to heater capacity and power consumption, with minimal external temperature loss.
- 7. Independent heating chamber for the heating elements to obtain an even heat distribution and rapid temperature equilibrium and stabilization.

Fan assisted convection models have a turbo fan. All incubators for bacteriology and cell culture have a second inner door of tempered glass.

Technology from J. P. Selecta:

- 8. Adjustable guide and shelf positions.
- Double seal around the chamber to provide a gentle but effective seal.
- 10. Floating spring door that adjusts the pressure and absorbs the thermal expansion.
- 11. Adjustable door pressure system closure. Internal tempered glass door.

NOTE:

For all models, the values for stability and homogeneity shown are based on temperature conditions with the ventilation closed.

The optimum homogenization of temperature within the chamber is based on a reasonable load that does not surpass more than 70% of the volume of the chamber. The graphic results shown for temperature for each model are based on the above criteria.

CONTROL PANELS

Models with Analogue control.

- 1. Main switch.
- 2. "On" indicator lamp.
- 3. Temperature control thermostat.
- 4. Heating "ON" indicator lamp.
- 5. Analogue thermometer temperature indicator.
- 6. Vacant positions for additional accessories.
- 7. Controllable safety thermostat that disconnects power to the heater in case of a fault in the main thermostat, manual reset (Directive DIN12880.2 class 2 and 3.1) and function signal lamp.



Models with 4.3 inches TFT touch screen. 1. Main switch.

2. TFT touch screen:

Visual audible alarm.

Clock calendar.

Single or cyclic On / Off programming.

Up to 10 work programs.

Up to 6 segments per program.

Stability time in each segment (from 1 min to 99h).

Alarms and events storage.

Probe error detection.

Self Diagnostics.

Ramps between segments.

Door open alarm.

Network failure detection and saving.

Over temperature and low temperature alarms and memorization (date, start time, end time and temperature).

Safety thermostat (TS) by software.

Mechanic safety thermostat (TS).

PC software.

User manual on screen.

Temperature control auto-tuning.

Configurable parameters: Date / time, temperature correction, data collection interval, language (English, Spanish and French), °C / °F selection, over temperature and low temperature limit.

- 3. RS-232 output.
- 4. USB output.
- 5. Security thermostat.
- 6. Ethernet output para for LAN connection.



MODEL SUMMARY TABLE

MODEL COMMENT	IDEL				
Models	CONTERM	DIGITHEAT	DIGITRONIC	INCUBAT	INCUDIGIT
TYPE	Drying Oven	Drying Oven	Universal	Bacteriological Incubator	Bacteriological Incubator
CONTROL	Temperature	Temperature + time	Temperature + time	Temperature	Temperature + time
DISPLAY	Analogue	Digital	Digital	Analogue	Digital
AIR	Convection	Convection	Fan assisted	Convection	Convection
CIRCULATION	natural	natural		natural	natural
CAPACITY LITRES	19 - 36 - 52 - 80 - 150	19 - 36 - 52 - 80 - 150	33 - 47 - 76 - 145	19 - 36 - 52 - 80 - 150	19 - 36 - 52 - 80 - 150

ACCESSORIES



2000002 Timer switch 0-120 minutes. Suitable for CONTERM.

2000003 Timer switch 0-12 hours. Suitable for CONTERM and INCUBAT.

2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes. Suitable for **CONTERM** and **INCUBAT**.



Part No.

2000016 Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours. Suitable for DIGITHEAT, DIGITRONIC and INCUDIGIT.



Optional communication modules

Part No. 2101623 Module for Wifi network.

Part No. 2101624 Module for Bluetooth.

Part No. 2101625 Module RF.

Part No. 2101626 RS-232 to RS-485 converter.

Suitable for DIGITHEAT, DIGITRONIC and INCUDIGIT.



DRYING AND STERILIZATION OVENS







NATURAL CONVECTION.

TEMPERATURE THERMOSTAT CONTROL WITH ANALOGUE THERMOMETER. FOR ADJUSTABLE TEMPERATURES FROM 40 °C UP TO 250 °C. STABILITY: ± 0.3 °C UP TO 150 °C. Homogeneity: $\pm 3,25$ °C UP to 150 °C.

SAFETY

STANDARD EN.61010. INCORPORATED FIXED OVER TEMPERATURE DEVICE .
STANDARD DIN 12880. (CLASS 2 AND 3.1) SAFETY THERMOSTAT CONTROLLER FITTED.

FEATURES, CONTROL PANEL, SAFETY, STANDARD AND ACCESSORIES (see pages 139 and 140).





Models Conterm, Part No. 2000208, 2000209 and 2000210.


Model Conterm type Poupinel, Part No. 2000200 and 2000201.

STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODELS

Part No.	Capacity litres	•	Width terior)	/ Depth cm	•	Width terior)	/ Depth cm	Shelf Positions	Power W	Weight Kg
2000208	19	30	25	25	50	60	44	5	640	27
2000209	36	40	30	30	60	65	49	7	950	35
2000200	52	33	47	33	53	82	52	5	1075	44
2000210	80	50	40	40	70	74	59	8	1230	54
2000201	150	50	60	50	70	95	68	8	2150	76



Performance graph of temperature and time.

A. Set at 250 °C: 1 h 30'. B. Set at 180 °C: 1 h 12'. C. Set at 100 °C: 1 h.

ACCESSORIES Accessories must

Accessories must be installed in the factory.



Part No.
2000002 Timer switch 0-120 minutes.
2000003 Timer switch 0-12 hours.



2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes.

SPARES

Shelves and guides.

Oven Part No.	2000208	2000209	2000200	2000210	2000201			
Guides set (2 units)	2000011	2000012	2000012	2000013	2000015			
Shelves	2000021	2000022	2000024	2000023	2000025			
Each shelf requires two quides (one set)								



Drying and sterilization ovens "Digitheat-TFT"

NATURAL CONVECTION.

MICROPROCESSOR CONTROL WITH TFT TOUCH SCREEN.

ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C UP TO 250 °C. STABILITY: ± 0.3 °C, UP TO 150 °C. Homogeneity: $\pm 3,25$ °C, UP TO 150 °C. SET ERROR: ± 2 % of the working temperature. Resolution: 1 °C.





SAFETY:

STANDARD EN.61010. INCORPORATED FIXED OVER TEMPERATURE DEVICE .
STANDARD DIN 12880. (CLASE 2 AND 3.1)CONTROLLABLE SAFETY THERMOSTAT FITTED.



Reaches working temperature with minimum delay

FEATURES, CONTROL PANEL, SAFETY, STANDARD AND ACCESSORIES (see pages 139 and 140).






Model Digitheat, Part No. 2001251, 2001252 and 2001254.

STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODELS

Part No.	Capacity litres	•	(Width Iterior)	om / Depth	•	Width terior)	n / Depth cm	Shelf Positions	Power W	Weight Kg
2001251	19	30	25	25	50	60	44	5	600	24
2001252	36	40	30	30	60	65	49	7	900	35
2001253	52	33	47	33	53	82	52	5	1000	44
2001254	80	50	40	40	70	74	59	8	1200	59
2001255	150	50	60	50	70	95	68	8	2100	73



Performance graph of temperature and time.

- A. Set at 250 °C: 60'.
- B. Set at 180 °C: 54'.
- C. Set at 100 °C: 48'.

ACCESSORIES

2000016 Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.



SPARES

Shelves and guides.

Oven Part No.	2001251	2001252	2001253	2001254	2001255			
Guides set (2 units)	2000011	2000012	2000012	2000013	2000015			
Shelves	2000021	2000022	2000024	2000023	2000025			
Each shelf requires two guides (one set).								



Universal precision ovens "Digitronic-TFT"

FAN ASSISTED CIRCULATION.

BACTERIOLOGICAL ASSAYS, DRYING PROCESSES AND STERILIZATION.
MICROPROCESSOR CONTROL WITH TFT TOUCH SCREEN.
ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 250 °C.

STABILITY: ± 0.3 °C, up to 100 °C. Homogeneity: ± 1 °C, up to 100 °C. Set error: $\pm 2\%$ of the working temperature. Resolution: 1 °C.







SAFETY:

EN.61012 STANDARD OVER TEMPERATURE SAFETY CUT OUT FITTED.
ADJUSTABLE OVER TEMPERATURE SAFETY THERMOSTAT DIN 12880. (CLASS 2 AND 3.1) FITTED.

Multipurpose oven. Fast response and recuperation of temperature.


FEATURES

- 1. TFT touch screen.
- 2. Inner chamber made of AISI 304 stainless steel.
- **3.** Pre-mixing chamber made of AISI 304 stainless steel.
- **4.** Homogeneously distributed shielded heating elements with complete air circulation throughout.
- **5.** Low external temperature due to excellent thermal insulation.
- **6.** Flexible silicon door gasket around the entrance of the chamber.
- Excellent door seal due to the floating inner door that adjusts and absorbs the thermal expansion.
- **8.** Turbo fan made of AISI 304 stainless steel that makes to circulate the air at the working temperature.
- **9.** Diagram showing the air flow from the pre-mixing chamber around the heating elements prior to entry to the oven's chamber.
- 10. Independent insulated control box.
- 11. Epoxy coated outer case.
- Ventilator with adjustable outlet (access at the back of the unit).
- 13. Adjustable height positions for guides and shelves.
- 14. Shelves made of AISI 304 stainless steel.


27,3 °C

15. Toughened double safety glass door for viewing the contents of the oven without having to open the door. (Model dependent).

CONTROL PANEL, SAFETY, STANDARD AND ACCESSORIES (see pages 139 and 140).







Model Digitronic with solid metal door. Part No. 2005163 and 2005167. (With toughened glass window door. Part No. 2005164 and 2005168).









Model Digitronic type Poupinel, door with toughened double glass window Part No. 2005166 and 2005170. Model Digitronic type Poupinel, Part No. 2005165 and 2005169.


STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODELS

INIODELO												
Part No.	Capacity litres	Door Type	Heating rate to 100 °C minutes	Recovery time* minutes	Complete air exchange per hour		t / Widt (interior	th / Depth) cm	Height / Width / Depth (exterior) cm	Shelf Positions	Power W	Weight Kg
2005163	33	metal	15	7	16	40	28	30	60 65 55	7	1200	38
2005164	33	glass	15	7	16	40	28	30	60 65 55	7	1200	40
2005165	47	metal	16	7	16	33	45	32	53 81 58	5	1200	46
2005166	47	glass	16	7	16	33	45	32	53 81 58	5	1200	50
2005167	76	metal	17	9	14	50	38	40	70 75 65	8	1600	58
2005168	76	glass	17	9	14	50	38	40	70 75 65	8	1600	64
2005169	145	metal	17	10	12	50	58	50	70 95 72	8	2000	74
2005170	145	glass	17	10	12	50	58	50	70 95 72	8	2000	79

^{*} Recovery time: the door was opened for 1 minute. After that, this is the time to recover the set temperature to 100 °C.



SPARES

Shelves and guides.

Oven Part No.	2005163 2005164	2005165 2005166	2005167 2005168	2005169 2005170					
Guides (2) (Set)	2000012	2000033	2000013	2000015					
Shelves	2000072	2000073	2000074	2000075					
Each shelf requires two guides i.e. one set.									

ACCESSORIES



Preparation of furnaces for drying moisture saturated samples. (Arids, muds, sands ...)
When adding turbine, the number of renewals of the air inside the furnace per hour multiplies by 10.

Must be factory installed Part No. 2000095



Drying and sterilization ovens "Dry-Big"

FAN ASSISTED CIRCULATION.
DIGITAL CONTROL AND DISPLAY OF TE

DIGITAL CONTROL AND DISPLAY OF TEMPERATURE AND TIME ADJUSTABLE TEMPERATURES FROM 40 °C UP TO 250 °C

STABILITY: ± 0.4 °C, up to 100 °C. Homogeneity: ± 2 °C, up to 100 °C set error : $\pm 2.5\%$ of the working temperature. Resolution: 1 °C



SAFETY:

STANDARD EN.61010. FIXED OVER TEMPERATURE DEVICE FITTED.
STANDARD DIN 12880. (CLASS 2 AND 3.1)ADJUSTABLE SAFETY THERMOSTAT FITTED.


Fast working and recovery temperature

FEATURES

- 1. Microprocessor controlled with digital display of temperature and time, pre-programmable time start and run time once the set temperature has been achieved through the Pt100 temperature sensor.
- 2. Inner chamber made of AISI 304 stainless steel.
- 3. Pre mixing chamber made of AISI 304 stainless steel.
- 4. Shielded heating elements with complete air circulation, homogeneously distributed throughout.
- **5.** Low external temperature due to excellent thermal insulation.
- **6.** Ventilation fan to force the air to circulate in the oven.
- 7. Diagram showing the air flow from the premixing chamber round the heating elements to the oven chamber.
- 8. Independent insulated control box.
- 9. Air inlet.
- 10. Ventilator with adjustable outlet of 120 Ø mm.
- 11. Shelf guides.
- 12. Shelves made of AISI304 stainless steel.
- **13.** Flexible silicon door gasket around the entrance of the chamber.
- 14. Epoxy coated outer case.

CONTROL PANEL

- 1. Illuminated mains switch.
- 2. Temperature mode indicator.
- 3. Time mode indicator.
- 4. Display for temperature and time.
- 5. Operating, Status mode.
- 6. Delay time state indicator.
- 7. Push button temperature selector.
- 8. Push button time selector.
- 9. Push button "increase" value or parameter.
- 10. Push button "decrease" value or parameter.
- 11. Push button Stop/Start.
- 12. Set temperature.
- **13.** Set run time: time period from 1 minute to 9 hours 59 minutes, or up to 99.9 hours, once the set temperature value has been reached.
- **14.** Set wait time before starting the run, time period from: 1 to 24 hours.
- **15.** RS-232 Interface output to a computer, for printer or USB adapter.
- Adjustable safety thermostat that overrides the microprocessor in case of failure, with manual reset and indicator lamp.









STANDARD EQUIPMENT


2 Shelves.

MODELS

MODELS										
Part No.	Voltage	Capacity litres	Heating rate to reach 100 °C, minutes	Recovery time* minutes	Air exchanges per hour	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Number of shelf positions	Power W	Weight Kg
2002961	230 / 400 three phase 230 single phase	216	16	10	12	60 60 60	87 112 84	6	4000	150
2002971	230 / 400 three phase 230 single phase	288	18	10	11	80 60 60	107 112 84	8	5000	161
DOUBLE [OOR CABINET									
2003721	230 / 400 three phase	400	18	13	6	100 80 50	128 132 74	10	5250	200
2003741	230 / 400 three phase	720	19	13	6	120 100 60	150 152 80	12	6000	264

Energy saving, three phase units are recommended.

^{*}Recovery time, the door was opened for 60 seconds, time taken to recover to the set temperature of 100 °C.



Performance graph of temperature and time.

A. Set at 250 °C: 1 h 6'.

B. Set at 180 °C: 42'. C. Set at 100 °C: 24'.

ACCESSORIES

4120131 USB adapter model. Pen-Drive included (Memory board) for

Accessories that must be installed in the factory:





with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.

2000007 Digital programmable microprocessor. Capacity: 10 programs of 100 segments. Programmable timer: up to 99 hours 59' 59". Program repetition: up to 99 times. Programs can also be linked for up to 4 stages.

RS-232 interface for data download to a printer or computer.

SPARES

Shelves.				
Oven Part No.	2002961/62	2002971/72	2003721	2003741
Sholvos	วกกกกรว	วกกกคว	วบบบบยว	2000064



High temperature oven "Hightemp"

FAN CONVECTION.

DIGITAL CONTROL AND DISPLAY OF TEMPERATURE AND TIME.

ADJUSTABLE TEMPERATURES FROM 60 °C UP TO 400 °C.

STABILITY: ±1 °C. UP TO 300 °C. HOMOGENEITY: ±3 °C. UP TO 300 °C.

SET ERROR: ±2 % OF THE WORKING TEMPERATURE.

SAFETY:

STANDARD DIN 12880. ADJUSTABLE OVER TEMPERATURE THERMOSTAT FITTED.

FEATURE

Digital electronic temperature control. Independent control box chamber thermaly insulated.

Shielded heating elements.

Fan circulation motor with thermal cut out, motor operates independently from the heating elements, the motor can be activated during the cooling cycle.

Inner chamber in AISI 304 heat resistant stainless steel with a high tolerance against corrosion and high temperatures.

Fixed position shelf guides.

Ventilation device with adjustable outlet.

Epoxy-coated outer casing.

STANDARD EQUIPMENT

2 shelves made of AISI 304 stainless stee

CONTROL PANEL

Main switch.

Mains indicator lamp.


Heater switch.

Heater operation indicator lamp.

Digital electronic temperature control.



Electronic safety thermostat with a K type probe that cuts off power to the heating elements in case of a controller fault. (standard to DIN 12.880 class 2).



Performance graph of temperature and time.

A. Set at 400 °c: 1h 50'.

B. Set at 200 °c: 1h.

ACCESSORIES. Shelves made of AISI 304 stainless steel. Part No. 2000071



Part No.	Voltage	Capacity litres	Height / (in	Width terior) of shelf		•	t / Width exterior)	/ Depth cm	Nº of shelf positions	Power W	Weight Kg
2001406	230 / 400 three phase	80	50	40	40	80	120	61	4	4000	158

ACCESSORIES

Accessories that must be installed in factory.

Part No.



2000002 Timer switch 0-120 minutes.

2000003 Timer switch 0-12 hours.



2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes.



Vacuum drying oven "Vaciotem-TV"

DIGITAL TEMPERATURE CONTROL, ELECTRONIC VACUUM PRESSURE DISPLAY AND TIMER. CONTROLLABLE TEMPERATURE FROM 35 °C TO 200 °C STABILITY ± 1 °C, UP TO 100 °C. HOMOGENEITY ± 3 °C, UP TO 100 °C. SET ERROR ± 2 °C. RESOLUTION 1 °C.





SAFETY:

OVER TEMPERATURE CUT OUT FITTED IN ACORDANCE WITH THE EN.61010 STANDARD. DIN 12880, STANDARD ADJUSTABLE SAFETY THERMOSTAT FITTED.

FEATURE

Digital electronic control of: temperature, vacuum pressure and pre-selected programmable timer.

Temperature sensor Pt100

Automatic air inlet at the end of the operation cycle.

Heating element placed evenly around the chamber.

Chamber made of AISI 304 stainless steel.

Trays made of anodised aluminium.

Door with hardened glass window, which sits on to a silicon gasket that absorbs any contractions and expansions that may occur.

Vacuum port with bleed valve.

Air valve at the front.

Vacuum pump connection at the back.

Epoxy covered outer case.

RS-232 Interface output for parameters to a computer, printer or USB adapter.



CONTROL PANEL

- 1. RS232 interface.
- 2. Air inlet.
- 3. Air inlet valve.
- 5. Vacuum pressure indicator lamp.
- 6. Air inlet valve indicator lamp, end of cycle.
- 7. Running indicator lamp.
- 8. Under vacuum indicator lamp.
- 9. Digital vacuum display in mbar.
- 10. Push button to select vacuum.
- 11. Push button to select electronic valve at the end of the cycle.
- 12. Push button to increase value.
- 13. Push button to decrease value.
- 14. Push button to STOP/START.
- 15. Indicator of mode temperature.
- 16. Indicator of mode time.
- 17. Indicator of operating.
- 18. Indicator of mode waiting time.
- 19. Digital display of temperature or time.

- 20. Push button to select temperature.
- 21. Push button to select time.
- 22. Push button to increase value.
- 23. Push button to decrease value.
- 24. Push button to STOP/START.
- 25. Mains switch.
- **26.** Safety thermostat in operation.

28 14 29 10 13 11 12 18 19 15 24 16 23 20 21 22 26 25

BACK

- 27. Air inlet.
- **28.** Pump power connection.
- 29. Vacuum connection.
- **30.** Adjustable safety thermostat

MODEL

Part No.	Vacuum Max.	Capacity litres	Ø / Depth (interior) cm	Height / Width / Depth (exterior) cm	Shelves	Power W	Weight Kg
4001490	10 ⁻² mm Hg	47	34 52	54 76 70	2	2000	73

Note: To obtain the optimum homogeneity at the set temperature, the load should not surpass more than 70 % of the volume of the chamber.

ACCESSORIES



4120131 USB adapter model. Pen-Drive included (Memory board) for data storage.

Accessories see page 150


SPARE PARTS

Shelves. (2) Part No. 2000030



Vacuum oven "Vaciotem-T"

DIGITAL TEMPERATURE AND TIMER CONTROL. Controllable temperature from 35 °C to 200 °C.





SAFETY:

OVER TEMPERATURE CUT OUT FITTED IN ACORDANCE WITH THE EN.61010 STANDARD.

DIN 12880. STANDARD ADJUSTABLE SAFETY THERMOSTAT FITTED.

FEATURE

Digital electronic control of temperature and pre-selected programmable timer.

Running time range: from 1 minute to 9hrs 59 min. or 99.9 hrs.

Pre-program start time, (wait time range): 1 hr to 24 hrs

Temperature sensor Pt100

Heating element placed evenly around the chamber.

Chamber made from AISI 304 stainless steel.

Trays made from anodised aluminium.

Door with hardened glass window, which sits on to a silicon gasket that absorbs any contractions and expansions that may occur.

Vacuum port with bleed valve

Air valve at the front

Vacuum pump connection at the back.

Epoxy covered outer case.

RS-232 Interface output of parameters for a computer, printer or USB adapter.




CONTROL PANEL

- 1. RS 232 connector.
- 2. Air inlet.
- 3. Air inlet valve.
- 4. Vacuum gauge.
- 15. Temperature mode indicator.
- 16. Time mode indicator.
- 17. Operation indicator.
- 18. Waiting time indicator.
- 19. Time and temperature digital display.
- **20.** Push button to select temperature.
- 21. Push button to select time.
- 22. Push button to increase value.
- 23. Push button to reduce value.
- 24. Push button to STOP/START.
- 25. Mains switch.
- **26.** Safety thermostat indicator lamp.
- 27. Vacuum pump control switch.

BACK

- 28. Vacuum pump power connection.
- 29. Vacuum connection.
- 30. Air inlet.
- **31.** Adjustable safety thermostat.



MODEL

Part No.	Vacuum Max.	Capacity litres	Ø / Depth (interior) cm	Height / Width / Depth (exterior) cm	Shelves	Power W	Weight Kg
4001489	10 ⁻² mm Hg	47	34 52	54 76 70	2	2000	73

Note: To obtain the optimum homogeneity at the set temperature, the load should not surpass more than 70 % of the volume of the chamber.

ACCESSORIES



4120131 USB adapter model. Pen-Drive included (Memory board) for data storage.

Accessories see page 150

SPARE PARTS

Shelves. (2)

Part No. **2000030**

VACUUM EQUIPMENT ACCESSORIES FOR VACIOTEM-T AND VACIOTEM-TV



COMECTA Vacuum pump "VACUM-10 Pa"

ROTARY VEIN PUMP WITH ANTI RETURN VALVE PREVENTS OIL FLOW BACK, SUITABLE FOR GENERAL LABORATORY APPLICATIONS. OVER TEMPERATURE MOTOR PROTECTION CUT-OUT AND MAIN ON/OFF SWITCH. RECOMMENDED FOR THE "VACIOTEM T AND TV" AND THE DESICCATOR "VACUO-TEMP".

FEATURE

Heat resistant veins and internal joints Aspiration inlet flange: 16 mm Ø. High oil volume and forced lubrication. Exhaust filter and ballast. Shock absorber mounted.
Free from vibrations
Low noise level (62db).
Maximum working temperature 60 °C.
Portable, with lifting handle included.

MODEL

Part No.	Pump rate m³/h	Vacuum limit mbar	Height / Width / Depth (exterior) cm	r.p.m.	Power W	Weight Kg
5900621	3,6	0.06	27 35 14	1400	180	11





Heated vacuum desiccator "Vacuo-Temp"

WITH TEMPERATURE THERMIC LIMITER.
TIME AND TEMPERATURE DIGITAL ELECTRONIC CONTROL.
ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C TO 170 °C.
STABILITY: ±1 °C. RESOLUTION: 1 °C.
TIME FROM 1' TO 999', OR CONTINUOUS.

FEATURES

AISI 304 stainless steel outer casing.

Polished aluminium alloy flat surface plate with an effective vacuum seal.

Tempered glass bell jar with silicon gasket seal. Shielded heating element.

Pt 100 temperature probe.

Vacuum pump connection at the back of the unit.
Vacuum bleed valve.

CONTROL PANEL

Main switch.

Analogue vacuum gauge.

Digital time & temperature display.

Overheating alarm.

Visualized parameter indicator.

Push button for the visualized parameter.

Push button to increase the parameter.

Push button to decrease the parameter.

Button On-Off.



MODEL

Part No.	Vacuum Max	Usable volume litres	Ø heating plate cm	Heigth / Width / (exterior) c			Power W	Weight Kg
4000474	10 ⁻² mm Hg	3	23.5	17	28	34	540	9

Supplied complete with bell jar and silicon seal.

SPARES

Tempered glass bell 15 cm high and 23 cm \emptyset . Part No. 4000475

Silicon seal. Part No. 4000476



Desiccator for materials

WITH HYDROMETER CONTROL.

APPLICATIONS

Cabinet with protection against humidity and dust for anhydrous, biological and chemical preservation of samples.

FFATIIRE

Made of robust transparent 12mm thick methacrylate. The door has a silicon seal and magnetic catch.

Volume: 55 Litres.

Dimensions 50 cm high x 38 cm wide x 29 cm deep. Supplied complete with three perforated shelves and a stainless steel AISI 304 tray to hold desiccating material.

Part No. 1001403





BACTERIOLOGICAL INCUBATORS





Bacteriological incubators "Incubat"

NATURAL CONVECTION.

TEMPERATURE THERMOSTAT CONTROL WITH ANALOGUE THERMOMETER. ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 80 °C. STABILITY: ± 0.1 °C, UP TO 37 °C. HOMOGENEITY: ± 0.5 °C, UP TO 37 °C INTERNAL GLASS DOOR.

FEATURES, CONTROL PANEL, STANDARD AND ACCESSORIES (see pages 139 and 140).

SAFETY:

OVER TEMPERATURE CUT OUT INCORPORATED ACCORDING TO THE EN.61010 STANDARD. ADJUSTABLE SAFETY THERMOSTAT DIN 12880. FITTED.





°C 70 60 50 40 30 20 10 20 40 60 80 120 Minutes.

STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODELS

Part No.	Capacity litres	•	/ Width terior) o	r / Depth cm	•	/ Width terior)	/ Depth cm	Shelves positions	Power W	Weight Kg
2000205	19	30	25	25	51	57	49	5	165	26
2000206	36	40	30	30	60	62	54	7	245	36
2001615	52	33	47	33	53	79	57	5	275	46
2000207	80	50	40	40	70	72	64	8	315	54
2000994	150	50	60	50	70	92	74	8	535	78

SPARES

Shelves and guides.

Oven Part No.	2000205	2000206	2001615	2000207	2000994					
Set guides (2 units)	2000011	2000012	2000012	2000013	2000015					
Shelves	2000021	2000022	2000024	2000023	2000025					
Each shelve requires two quides (one set).										

Performance graph of temperature and time.

A. Set at 80 °C: 1 h 54'.

B. Set at 56 °C: 1 h 46'.

C. Set at 37 °C: 1 h 18'.

ACCESSORIES

Accessories must be factory installed.



Part No. 2000003 Timer switch 0-12 hours.



2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes.



Digital bacteriological incubators "Incudigit-TFT"

NATURAL CONVECTION.

DIGITAL CONTROL AND DISPLAY OF TEMPERATURE AND TIME.

ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C UP TO 80 °C.

STABILITY: ±0.1 °C, UP TO 37 °C. HOMOGENEITY: ±0.5 °C, UP TO 37 °C.

SET ERROR: ±2% OF THE WORKING TEMPERATURE, RESOLUTION 0.1 °C

INTERNAL TEMPERED GLASS DOOR.

DOUBLE CHAMBER, MINIMUM RISK OF SAMPLE CONTAMINATION.
INSIDE WITHOUT OPENINGS AND WITH ROUNDED CORNERS. EASY TO CLEAN.

FEATURES, CONTROL PANEL, STANDARD AND ACCESSORIES (see pages 139 and 140).



OVER TEMPERATURE CUT OUT INCORPORATED ACCORDING TO THE EN.61010 STANDARD.

ADJUSTABLE SAFETY THERMOSTAT DIN 12880. FITTED.


















Performance graph of temperature and time.

- A. Set at 80 °C: 1 h 12'.
- B. Set at 56 °C: 54'.
- C. Set at 37 °C: 48'.

STANDARD EQUIPMENT
2 shelves and 4 shelf guides.

MODELS

MODELO										
Part No.	Capacity litres	•	/ Width terior) o	n / Depth cm	•	/ Width terior)	n / Depth cm	Shelves positions	Power W	Weight Kg
2001256	19	30	25	25	51	57	49	5	150	26
2001257	36	40	30	30	60	62	54	7	225	36
2001258	52	33	47	33	53	79	57	5	250	46
2001259	80	50	40	40	70	72	64	8	300	54
2001260	150	50	60	50	70	92	74	8	525	75

SPARES

Shelves and guides.

Oven Part No.	2001256	2001257	2001258	2001259	2001260
Guides (2) (Set)	2000011	2000012	2000012	2000013	2000015
Shelves	2000021	2000022	2000024	2000023	2000025
Each self requires two guides	i.e. one set.				

ACCESSORIES

must be factory installed.



2000016 Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.



Incubators for bacteriology and cell culture "Incubig-TFT"

SAFETY: STANDARD EN.61010 OVER TEMPERATURE CUT OUT FITTED. STANDARD DIN 12880. ADJUSTABLE SAFETY THERMOSTAT FITTED.

NATURAL CONVECTION.



STABILITY: ± 0.1 °C, UP TO 37 °C. Homogeneity: ± 0.5 °C, UP to 37 °C. Set error: $\pm 2\%$ of the working temperature, resolution 0.1 °C

INTERNAL TEMPERED GLASS DOOR.







Capacities up to 720 litres

FEATURE

Microprocessor control and 4.3 inches TFT touch screen display.

Large surface area heating elements.

Inner chamber made of AISI 304 stainless steel.

Double door, interior door of tempered glass that allows the user to see the contents of the chamber without opening the door.

Adjustable air vent.


Epoxy covered external case.

STANDARD EQUIPMENT

For Part No. 2000238, 2 shelves and 4 shelf guides. For Part No. 2000239 and 2000240, 2 shelves.

Model Part No. 2000238.





Performance graph of temperature and time.

A. Set at 80 °C: 1 h 45'.

B. Set at 56 °C: 1 h 10'.

C. Set at 37 °C: 54'.

Note: To obtain the optimum homogeneity at the set temperature, the load should not surpass more than $70\,\%$ of the volume of the chamber.







CONTROL PANEL

4.3 inches TFT touch screen models:

1. Main switch.

2. TFT touch screen:

Visual audible alarm.

Clock calendar.

Single or cyclic On / Off programming.

Up to 10 work programs.

Up to 6 segments per program.

Stability time in each segment (from 1 min to 99h).

Alarms and events storage.

Probe error detection.

Self Diagnostics.

Ramps between segments.

Door open alarm.

Network failure detection and saving.

Over temperature and low temperature alarms and

memorization (date, start time, end time and temperature).

Safety thermostat (TS) by software.

Mechanic safety thermostat (TS).

USB and RS -232 output.

PC software.

User manual on screen.

Temperature control auto-tuning.

Configurable parameters: Date / time, temperature correction, data collection interval, language (English, Spanish and French), °C / °F selection,

over temperature and low temperature limit.

- 3. RS-232 output.
- 4. USB output.
- 5. Security thermostat.
- 6. Ethernet output para for LAN connection.



MODELS

Part No.	Туре	Capacity litres	•	/ Widtl nterior)		•	Width terior)		№ of shelf guides	Power W	Weight Kg
2000238	1 door	288	80	60	60	97	91	76	8	570	87
2000239	2 door	400	100	80	50	130	114	75	10	1100	160
2000240	2 door	720	120	100	60	152	134	85	12	1600	225



4120131 USB adapter model.

Pen-Drive included (Memory board) for data storage.

Accessories must be factory installed.



2000016 Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.

SPARES Shelves and guides.

0240
0064
•

Each self requires two guides i.e. one set.



Incubator for Petri capsules

NATURAL CONVECTION.

MICROPROCESSOR REGULATION AND TEMPERATURE DIGITAL CONTROL.

FOR ADJUSTABLE TEMPERATURES FROM AMBIENT +5°C TO 60°C.

STABILITY: ±0,1°C TO 37°C. HOMOGENEITY: ±0,1°C TO 37 °C. SETPOINT ERROR: ±0,1°C. RESOLUTION: 0,1°C.

Small size. Culture visual control. Transportable.

APPLICATIONS

Specially designed for bacteria and fungi cultures in Petri capsules at the same temperature of human body.

Culture surface 320 x 220 mm (Inner height: 20mm)

Culture visual monitoring.

Transparent cover.

Easy access to samples.

Approximate capacity: (single level) (mm)

15 Petri capsules of Ø55.

10 Petri capsules of Ø80.

7 Petri capsules of Ø90.

6 Petri capsules of Ø100.

3 Petri capsules of 120x120.

2 Petri capsules of Ø140.

MODEL

Part No.	Margin T ^a °C		/ Width terior) (/ Depth cm		/ Width terior)	n / Depth cm	Power W.	Weight Kg
4002629	Amb +5 to 60	7	45	25	2	32	22	300	3.5





Cooled low temperature incubator "Prebatem-TFT"

FORCED AIR FAN CIRCULATION.

MICROPROCESSOR CONTROLLED WITH DIGITAL DISPLAY

ADJUSTABLE TEMPERATURES FROM 5 °C UP TO 60 °C. RESOLUTION 0.1 °C SEMICONDUCTOR HEATING AND COOLING SYSTEM.

QUIET-STABLE - FREE FROM VIBRATIONS - VERY ACCURATE - LOW POWER CONSUMPTION.

INNER TEMPERED GLASS DOOR.











Leading edge technology, Peltier effect. No compressor.

SAFETY: CONFORMS TO THE DIN 50011 STANDARD FOR TEMPERATURE STABILITY AND HOMOGENEITY.

CONFORMS TO THE DIN 12880.STANDARD ADJUSTABLE SAFETY THERMOSTAT FITTED.


APPLICATIONS

Biotechnology, Bacteriology, Plasma fractionation, Biology, Enzymatic test, Research, Serum studies, metrology, Botany, Phytopharmacy, Cosmetics, Water analysis and Agricultural research, feeding, new techniques for protein crystallization.

FEATURE

- 1.4.3 inches TFT touch screen.
- 2. Inner chamber and elements made of AISI 304 stainless steel.
- 3. Premixing temperature chamber.
- Semiconductor- static radiator for heating and cooling.
- **5.** Excellent thermal insulation within the chamber.
- 6. Turbo fan to make the air circulate.
- **7.** Diagram showing the homogeneous air flow from the premixing chamber of the semiconductor cooling / heating system.
- 8. Independent insulated control box .
- 9. Support rack for trays.
- 10. Shelves of AISI 304 stainless steel.
- **11.** Epoxy coated outer case.

PERFORMANCE	Specific	ation
	at 10 °C	at 37 °C
Stability	±0.5 °C	±0.1 °C
Homogeneity	±0.1 °C	±0.3 °C
Set error	±0.25 °C	±0.20 °C



CONTROL PANEL

- 1. Main switch.
- 2. TFT touch screen:

Visual audible alarm.

Clock calendar.

Single or cyclic On / Off programming.

Up to 10 work programs.

Up to 6 segments per program.

Stability time in each segment (from 1 min to 99h).

Alarms and events storage.

Probe error detection.

Self Diagnostics.

Ramps between segments.

Door open alarm.

Network failure detection and saving.

Over temperature and low temperature alarms and memorization (date, start time, end time and temperature).

Safety thermostat (TS) by software. Mechanic safety thermostat (TS).

USB and RS -232 output.

PC software.

User manual on screen.

Configurable parameters: Date / time, temperature correction, data collection interval, language (English, Spanish and French), °C / °F selection, over temperature and low temperature limit.

- 3. RS-232 output.
- 4. USB output.
- **5.** Security thermostat.
- 6. Ethernet output para for LAN connection.



CONTROL PANEL

Main switch.

Mains indicator lamp.

Microprocessor control and digital temperature display.

Adjustable safety thermostat.




STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODELS

Part No.	Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Shelf guides	Power consumption W/hr. at 5 °C at 40 °C	Power W	Weight Kg
2000963	36	40 30 30	60 65 49	7	70 50	310	54
2000964	80	50 40 40	70 75 59	8	75 55	310	73
2000965	150	50 60 50	70 95 68	8	90 60	310	94



Performance graph of temperature and time.

A. Set at 50 °C: 40'.

B. Set at 0°C: 48'.

Note: To obtain the optimum homogeneity at the set temperature, the load should not surpass more than 70 % of the volume of the chamber.



ACCESSORIES

Accessories must be factory installed.



Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours. Part No. **2000016**



Optional communication modules

Part No. 2101623 Module for Wifi network. Part No. 2101624 Module for Bluetooth.

Part No. 2101625 Module RF.

Part No. 2101626 RS-232 to RS-485 converter.



Shelves and guides.

Oven Part No.	2000963	2000964	2000965
Guides (2) (Set)	2000012	2000013	2000015
Shelves	2000022	2000023	2000025

Each self requires two guides i.e. one set.



Incubation chamber "Boxcult"

FAN CONVECTION.

ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 57 °C. STABILITY: ± 0.25 °C, up to 37 °C homogeneity: ± 1 °C, up to 37 °C. Set error: ± 2 % of the working temperature. Resolution 0.1 °C.

SAFETY:

SAFETY STANDARD EN 61010. OVER TEMPERATURE SAFETY THERMOSTAT FITTED.

FEATURE

Made of transparent methacrylate that allows the user to see inside the incubator during operation. To facilitate the access to the working area the unit has a wide front door, and a removable base made of AISI 304 stainless steel. The fan convection circulation system ensures an even and rapid recovery of temperature.

A 30 mm \emptyset port at the rear can be used to connect power to apparatus inside the chamber.

Supplied as accessories, the removable base allows the Boxcult to be mounted on the "Rotabit" reciprocal / orbital shaker. (described in the stirrer section.)

The metallic top of the chamber includes the heating elements, air circulation fan and temperature control.

CONTROL PANEL

Main switch.

Digital electronic temperature control.



Incubation chamber "Boxcult" Part No. 3000957 with base Part No.3001172 and support rack with two shelves Part No. 1000973. Supplied as accessories.

MODEL

Part No.	Capacity	Height / Width / Depth	Height / Width / Depth	Power	Weight
	litres	(interior) cm	(exterior) cm	W	Kg
3000957	110	50 47 47	61 51 51	220	18

Supplied without bottom base, or stainless steel rack and shelves.



Orbital & shaker stirrer "Rotabit" part number 3000974 with incubation chamber "Boxcult"

For an easy handling, all control devices are out-

side the chamber enclosure

ACCESSORIES

Removable bottom base made of AISI 304 stainless steel. Part No. 3001172

Stainless steel rack with 4 shelves positions, each one separated by 9 cm. Comes complete with 2 removable shelves. Useful dim. 43 cm long and 41 cm wide. Part No. 1000973



CO2 Incubators for anaerobic cell and tissue cultures "Incubator CO2"

MICROPROCESSOR CONTROL WITH DIGITAL DISPLAY OF TEMPERATURE AND CO2.

ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C TO 50 °C

STABILITY: ±0.2 °C, UP TO 37 °C. HOMOGENEITY: ±0.5 °C, UP TO 37 °C. RESOLUTION: 0.1 °C.

ALARM RANGE: FROM AMBIENT+5 °C TO 50 °C. RESOLUTION: 0.1 °C. CO₂ RANGE: FROM 0 TO 20%. STABILITY: ±0.3%. RESOLUTION: 0.1%





SAFETY:

STANDARD DIN 12880. DOUBLE INDEPENDENT OVER TEMPERATURE SAFETY THERMOSTAT.

CO2 DEVIATION FROM SET VALUE. OPEN DOOR INDICATOR, ELECTRICAL FAULT INDICATOR. LOW CO2 PRESSURE.

Infrared CO2 sensor

Chamber sterilization function

FEATURE

External case of steel coated with epoxy with insulated chamber.

The chamber is made of stainless steel with removable shelf supports and easy clean system.

Two doors; one interior of tempered glass with silicon gasket and a heated external steel door with magnetic seal to prevent condensation on the glass door.

Smooth door action, to prevent jolts or vibrations disturbing the contents of the incubator.

The CO_2 input is by a metal tube of 6 mm Ø x 4 mm at the back of the unit.

RS-232 Interface output for a computer, printer or USB adapter.

CONTROL SYSTEM

Digital electronic control of temperature and CO_2 , by a single multilevel control button and LCD screen, that controls all functions within the chamber.

HUMIDITY CONTROL

The humidity level within the chamber is at a constant 98% RH level, that is produced directly by water evaporation previously introduced at the bottom of the chamber.





CONTROL PANEL

- 1. Visual alarm indicator.
- 2. LCD display of all parameters.
- 3. Button single control of multilevel functions.
- 4. Printer (Optional)
- 5. Main On switch.

MODEL

Part No.	Capacity litres	Height/Width/Depth (interior) cm		Height/Width/Depth (exterior) cm			Shelf guide positions	Power W	Weight Kg	
4002628	150	65	50	46	95	65	73	9	800	110

Comes with two shelves.

ADDITIONAL Shelves stainless steel. Part No. 1001675

ACCESSORIES



USB adapter model.

Pen-Drive included (Memory board) for data storage.

Part No. 4120131

Printer: temperature, CO2, time and status.

(Needs to be factory fitted.)
Part No. **4001676**



Fyrite CO₂ analyser.

Monitor for checking the CO_2 % concentration. The unit has a graduated scale of 0 to 20 %. Reagent valid for 300 analysis. Should not be used with explosive gasses.

Part No. 4000632

Reagent flask 64 ml. Part No. 4000635



PRECISE COOLED INCUBATORS HOTCOLD



HOTCOLD S HOTCOLD A-B-C HOTCOLD UB-UC

HOTCOLD GL

CONTROLLABLE TEMPERATURES FROM +5 °C TO 65 °C CONTROLLABLE TEMPERATURES FROM 0 °C TO 50 °C CONTROLLABLE TEMPERATURES FROM -10 °C TO 50 °C

CONTROLLABLE TEMPERATURES FROM 0 °C TO 50 °C (DEPENDING ON WORKING MODE)

SAFETY:

DIN STANDARD 12880.2

SAFETY THERMOSTAT FITTED THAT DISCONNECTS POWER TO THE HEATER IF THE CONTROLLER FAILS. MANUAL RESET.

APPLICATIONS

Enzymatic tests, serum and plasma fractions BOD tests, cosmetics, botany, pharmacy, industry, agriculture, bacteriology, biotechnology and research.



Refrigerated cabinet "Hotcold S"

FORCED AIR CIRCULATION. DIGITAL ELECTRONIC CONTROL OF TEMPERATURE AND TIME, ADJUSTABLE FROM +5 °C TO 65 °C. STABILITY ± 0.1 °C, up to 20 °C. Homogeneity ± 0.5 °C, up to 20 °C. Set error ± 2 °C. Resolution 0.1 °C.



Epoxy coated external case. Interior AISI304 stainless steel. Door with double glazed glass to maintain internal temperature. Illumination switch with internal fluorescent light. Side port for the introduction of external cables probes and tubes etc.

Cooling gas R134a.

4 wheels with brake.

Two safety power sockets.

CONTROL SYSTEM

Electronic digital controller for temperature and time. Timer and off programmable from 1' to 99 hrs 59'.

Programmable defrost.

High and low temperature alarm.

Temperature calibration.

CONTROL PANEL

- 1. Display for temperature / time.
- 2. Temperature indicator.
- 3. Time indicator.
- 4. Alarm indicator.
- 5. Heater functioning indicator.
- 6. Push button for set temperature.
- 7. Push button for set time.
- 8. Mains switch.
- 9. Push button to increase value.
- 10. Push button to decrease value.
- 11. Push button to confirm value.

STANDARD EQUIPMENT

2 shelves and 4 shelf guides.

MODEL

MUDEL				İ		even it	ıau uistiivutivii vi up	lu 10% uiiil voiuiile	
HOTCOLD	Part No.	Range °C	Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Number of shelves	Motor HP	Power W	Weight Kg
S	2101618	+5 +65	160	50 65 50	163 128 63	10	3/8	400	70
U	2101010	TU TUU	100	00 00 00	100 120 00	10	3/0	400	10


Guides (2) (Set).

Shelves.

1001619

1001620







Optimum temperature homogenization can be achieved with an even load distribution of up to 70% unit volume.



Precise refrigerated cabinets "Hotcold A-B-C-UB-UC"







DIGITAL ELECTRONIC CONTROL OF TEMPERATURE AND TIME.
HOTCOLD A-B-C ADJUSTABLE TEMPERATURE FROM 0 °C TO 50 °C.
HOTCOLD UB-UC ADJUSTABLE TEMPERATURE FROM -10 °C TO 50 °C.






SAFETY:

DIN STANDARD 12880.2 SAFETY THERMOSTAT FITTED THAT DISCONNECTS POWER TO THE HEATER IF THE CONTROLLER FAILS. MANUAL RESET.





PRECISION TABLE



CONTROL PANEL

- 1. Graphic display.
- 2. Push button increase value.
- 3. Push button decrease value.
- 4-5. Push button move cursor.
- 6. Push button, validate set value.
- 7. Push button set program.
- 8. Push button start.
- 9. Push button stop.
- **10.** Printer for time and temperature (Optional). Part No. 2101508, (needs to be factory fitted.)

Models A, B and UB.

FEATURES

Exterior case, door and interior made from AISI 304 stainless steel. steel shelves. PVC laminated.

Reversible door which can be opened by either side, with easy to change the lock and the joint, manual lockout.

Hermetically sealed compressor with anti vibration mounts with fan forced evaporation unit with ventilated condenser

Homogeneous internal temperature by forced circulating air.

Refrigerant R134 for models B and C.

Refrigerant R404 for model UB and UC.

Two safety power sockets.

Models C and UC.

Two external ports for external connections.

Temperature calibration.

Adjustable timer.

RS-232 Interface output for a computer, printer or USB adapter.

CONTROL SYSTEM

Programmable graphic display of temperature and time. All operation modes are programmable in up to 5 cycles and each cycle can be programmed at 1 hour intervals. These cycles can be repeated indefinitely or can be manually finished. Up to 5 different temperatures can be programmed into each cycle. All parameters can be stored. A temperature safety thermostat comes as standard.

MODELS

HOTCOLD	Part No.	Range °C	Capacity litres	•	: / Widtl iterior)	h / Depth cm	•	/ Width xterior)	/ Depth cm	Number of shelves	Included shelves	Power W	Power HP	Weight Kg
A	2101502	0 +50	319	139	48	45	198	60	64	14	3	236	1/5	78
В	2101503	0 +50	442	126	58	55	198	70	71	14	3	236	1/5	89
C	2101504	0 +50	600	136	58	69	207	70	82	14	4	351	3/8	100
UB	2101505	-10 +50	442	126	58	55	198	70	71	14	3	597	1/2	94
UC	2101506	-10 +50	600	138	58	69	207	70	82	14	4	593	1/2	110

NOTE: the HOTCOLD has internal power sockets that allows the use of a non-heating mixer shaker or stirrer or equipment for BOD assays to be powered internally. Alternatively power cables can be fed through external ports at each side of the unit. See chapter Mixers stirrers and shakers.

ACCESSORIES



USB adapter model.

Pen-Drive included (Memory board) for data storage.

Part No. 4120131

Printer shows temperature and time. Needs to be factory fitted.

Part No. 2101508

SPARES

Brackets and guides.

Cabinet Part No.	2101502	2101503 / 2101505	2101504 / 2101506						
Brackets set (4 units).	1001801	1001801	1001801						
Shelves	1001804	1001805	1001806						
Each self requires 4 brackets i.e. one set.									



Precision refrigerated cabinets "Hotcold GL"

FORCED AIR CIRCULATION.

DIGITAL ELECTRONIC CONTROL OF TEMPERATURE.TIME AND HUMIDITY. SUITABLE FOR TEMPERATURES FROM 5 °C TO 50 °C.

SAFETY:

SAFETY STANDARD: CONFORMS TO THE DIN 12880.2. ADJUSTABLE OVER TEMPERATURE CUT OUT FITTED THAT CUTS OFF HEATING IF OVER TEMPERATURE FAILS, MANUAL RESET.

APPLICATIONS

Refrigerated climate cabinet for botany testing of plants flowers, seed germination, photosynthesis agriculture etc. that require control of temperature, humidity and

Specifically designed with four function modes:

Mode A: Refrigerated incubator from 5 °C to 50 °C.

Mode B: Refrigerated incubator with illumination from 10 °C to 50°C

Mode C: Refrigerated incubator from 18 °C to 40 °C with an adjustable humidity

range from 50 to 98%.

Mode D: Refrigerated incubator with illumination from 18 °C to 40 °C with an

adjustable humidity range from 50 to 98%.

	Mode A	Mode B	Mode C	Mode D
Humidity	NO	NO	YES	YES
Illumination	NO	YES	NO	YES
Temperature range	5 / 50 °C	10 / 50 °C	18 / 40 °C	18 / 40°C
Stability	±0.5 °C	±0.5 °C	±0.5 °C	±0.5 °C
Homogeneity	±1.0 °C	±1.0 °C	±1.0 °C	±1.0 °C
Resolution	0.1 °C	0.1 °C	0.1 °C	0.1 °C
Set Error	±1.0 °C	±1.0 °C	±1.0 °C	±1.0 °C
Humidity range	-	-	50 / 98% Hr	50 / 80% Hr
Illumination range (Choice of 3)	-	0 / 4 K / 12 K lux	-	0 / 4 K / 12 K lux
Humidity resolution	-	-	1%	1%
Humidity precision 50 - 75%	-	-	±3%	±4%
75 - 80%	-	-	± 4 %	±5%

FEATURES

Exterior case, door and interior made from AISI 304 stainless steel. Reversible door can be fitted to open from either side, with automatic closing if left open.

The door interior supplies fluorescent illumination to the chamber, the power of which can be selected as 0 / 4K or 12 K Lux.

Hermetically sealed compressor with anti vibration mounts with fan forced evaporation unit with ventilated condenser.

Fan circulated homogeneous temperature.

All operation modes are programmable in up to 5 cycles of which each cycle can be programmed at 1 hour intervals. These cycles can be repeated indefinitely or can be manually terminated.

The humidity is constant during the program.

There are 10 Program storage memories of all parameters. The fan, temperature radiator and two thermal safety internal electrical sockets are located in the upper chamber.

Two external ports are located on each side for the introduction of tubes and cables for other diverse applications.

A humidity tray is located at the back of the unit for controlling humidity and is generated through eva-

RS-232 Interface output for a computer, printer or USB adapter.

CONTROL PANEL

- 1. Graphic display.
- 2. Push button increase value.
- 3. Push button decrease value.
- 4-5. Push button move cursor.
- 6. Push button, validate set value.
- 7. Push button set program.
- 8. Push button start.
- 9. Push button stop.
- **10.** Printer for time and temperature (Optional).

Part No. 2101508, (needs to be factory fitted).

MODEL

HOTCOL	.D Part No.	Range °C	Capacity litres		/ Widt terior)	- 1		Width terior)		Number of shelves	Motor HP	Power W	Weight Kg
GL	2101507	5 +50	600	138	58	69.5	208	75	115	14	3/8	950	265

NOTE: The HOTCOLD has internal power sockets that allow the use of a non-heating mixer shaker or stirrer or equipment for BOD assays to be powered internally. Alternatively power cables can be fed through external ports at each side of the unit. See chapter Mixers stirrers and shakers.



PRECISION TABLE

STABILITY	±0.5	°C
HOMOGENEITY	±1	°C
SET ERROR	±1	°C
RESOLUTION	0.1	°C





STANDARD EQUIPMENT

2 shelves and 8 brackets.



ACCESSORY



USB adapter model.

Pen-Drive included (Memory board) for data storage. Part No. 4120131

Printer shows temperature and time. Needs to be factory fitted.

Part No. 2101508

SPARES

Part No.

1001801 Brackets (4) (Set) **Shelves** 1001806

Each self requires 4 brackets i.e. one set.

Muffle Furnaces



Electric Muffle Furnaces "Select-Horn-TFT"

TEMPERATURE CONTROLLABLE UP TO 1150 °C.

SET ACCURACY: ±1 °C OF THE SET VALUE. RESOLUTION: 1 DIGIT.

DIGITAL ELECTRONIC CONTROLLER FOR TEMPERATURE AND TIME WITH TFT COLOUR TOUCH SCREEN.



PROBE BREAK DISCONNECTS THE POWER TO THE FURNACE AUTOMATICALLY. MICROSWITCH THAT DISCONNECTS THE POWER OF THE HEATER ELEMENTS WHEN THE DOOR IS OPEN. FLIP DOOR THAT CAN ALSO BE USED AS A SUPPORT TRAY AND USER PROTECTED FROM THE HOT INTERNAL SURFACE.







APPLICATIONS

Incineration processes, drying, degradation, re-heating, thermal treatments etc.

FEATURES

Interior chamber constructed from high quality lightweight refractory bricks, with a high alumina content with no asbestos or iron oxide.

Evenly distributed exceptional long life heating elements, annealed frequently at a high fusion point. Excellent thermal insulation made from Ceramic fibre of low density and thermal conductivity.

Low consumption with maximum performance.

Rapid temperature recovery after the door has been

Flap door with easy to change components.

Support tray made from special steel used as a base to support assay material.

USB and RS-232 output.



CONTROL PANEL

General

Main switch.

TFT touch screen 4.3".

Clock calendar.

Two working modes, normal or programming,

SPA - FRE - ENG menu.

Self-test on starting.

Temperature control auto-tuning.

°C/°F selection.

Type K probe.

Normal mode

Set point temperature selection

Up ramp or no ramp.

Stability time from 1 min to 99h or continuous.

Programming mode

10 programs capacity.

6 seaments per program.

Stability time in each segment from 1 min to 99h (or continuous in the last segment)

Up ramps between segments or no ramps

Daily - weekly On / Off programming.





Temperature ramps graph

Alarms

Network failure detection alarm.

Probe error detection alarm.

Over temperature and low temperature alarms.

Visual audible warning alarms.

Up to 100 alarms storage (date, start time, end time and alarm type).

Datalogging

Datalogging memory up to 15000 data. Logging interval from 5 seconds to 30 min. Data download via RS -232 or USB.

PC software for on-line registration (via RS-232).



Image of the flap door system in operation.

MODELS

MODELO					
Part No.	Capacity litres	Height / Width / Width (interior) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
2000376	3.6	10 15 24	52 54 56	2500	54
2000377	9	15 20 30	58 59 65	3000	70

Supplied complete with support tray, made from annealed steel.

SPARES

Support tray made from special steel used as a base to support assay material.

Code **0203681** for furnace Part No. 2000366 Code 0203692 for furnace Part No. 2000367





Electric muffle furnaces "R-3 L" and "R-8 L" 1100 °C



FOR TEMPERATURES ADJUSTABLE UP TO 1100 °C. MICROPROCESSOR CONTROL WITH TFT TOUCH SCREEN. PRECISION ±2 °C OF THE SET VALUE.

RESOLUTION: 1 DIGIT.



Metal external case with vent at the back of the unit. Interior and door made of ceramic fibre, resistant and durable (No asbestos). Heater situated at the side and bottom of the chamber.

USB and RS-232 output.



General

Main switch.

TFT touch screen 4.3".

Clock calendar.

Two working modes, normal or programming.

SPA - FRE - ENG menu.

Self-test on starting.

Temperature control auto-tuning.

°C/°F selection.

Type K probe.

Normal mode

Set point temperature selection

Up ramp or no ramp.

Stability time from 1 min to 99h or continuous.

Programming mode

10 programs capacity.

6 segments per program.

Stability time in each segment from 1 min to 99h (or continuous in the last segment)

Up ramps between segments or no ramps Daily - weekly On / Off programming.



Alarms

Network failure detection alarm.

Probe error detection alarm.

Over temperature and low temperature alarms. Visual audible warning alarms.

Up to 100 alarms storage (date, start time, end time and alarm type).





Temperature ramps graph

Datalogging

Datalogging memory up to 15000 data. Logging interval from 5 seconds to 30 min. Data download via RS -232 or USB. PC software for on-line registration (via RS-232).

MODEL	Part No.	Capacity litres		Width / terior) cı		Height / (ex	Width terior)		Power W	Weight Kg
N-3	2000368	3	11.5	12.5	20	43	34	47	1700	18
N-8	2000369	8	13	20	30	51	44	56	1800	28
Supplied complete with a refractory ceramic tray as a base and support for material to be assayed.										

ACCESSORIES FOR MUFFLE FURNACES

Adaptable only for "Select-Horn-TFT" furnaces Part No. 2000376 and 2000377 All accessories need to be fitted in the factory prior to delivery.





Exterior exhaust tube.

Located at the furnace back with a ventilator motor to extract gases and vapours. With an 80 mm Ø hat adapter.

Gases and Vapours can be extracted outside through the connecting tube.

Power consumption: 30 W.

Part No. 2001477

COMPLEMENTS



Gloves Thermal "Kevlar 800"

Conforms to EN 388, EN407 and EN420 standards. For use with temperatures up to 800 °C, Made from seamless terry knit, with double face fibres, high level of protection against heat and flame.

Length 36 cm, universal fit. Part No. 5000042



Crucible tongs.

With thermally protected plastic coated handles. With bow, curved tips. Part No. 1001590 Total length 220 mm. Part No. 1001591 Total length 330 mm.



Crucibles made of zirconium Zr. Crucibles made of pure nickel Ni. Crucibles made of glazed porcelain. Crucibles made of stainless steel. Crucibles made of quartz. (See page 182).



COMECTA: Electric muffle furnace "N-30 L" 1300 °C



FOR TEMPERATURES ADJUSTABLE UP TO 1300 °C. ELECTRONIC DIGITAL TEMPERATURE CONTROL. PRECISION ±2 °C OF THE SET VALUE. **RESOLUTION: 1 DIGIT.**


FEATURES

Metal external case with vent at the back of the unit. Interior and door made from ceramic fibre, resistant and durable (No asbestos). Heater situated at the side and bottom of the chamber.

CONTROL PANEL

Illuminated mains On/Off switch. Temperature control with digital display of both the set and actual temperature. Programmable in steps of 1 °C. Fitted with a type K probe.







MODEL	Part No.	Capacity litres	Height / Width / Depht (interior) cm	Height / Width / Depht (exterior) cm	Power W	Voltage V	Weight Kg
N-30 L	2200853	30	27.5 24 43	63 87 84	4600	230	120

Supplied complete with a refractory ceramic tray as a base and support for material to be assayed.



COMECTA Electric muffle furnaces "N-3 L", "N-8 L", "N-13 L", "N-22 L"and "N-80 L" 1100 °C

FOR TEMPERATURES ADJUSTABLE UP TO 1100 °C. ELECTRONIC DIGITAL TEMPERATURE CONTROL. PRECISION ±2 °C OF THE SET VALUE. **RESOLUTION: 1 DIGIT.**


FEATURES

Metal external case with vent at the back of the unit. Interior and door made from ceramic fibre, resistant and durable (No asbestos). Heater situated at the side and bottom of the chamber.

CONTROL PANEL

Illuminated mains On/Off switch. Temperature control with digital display of both the set and actual temperature. Programmable in steps of 1 °C. Fitted with a type K probe.






Temperature



MODELO	Part No.	Capacity litres	Height / Width / Depht (interior) cm	Height / Width / Depht (exterior) cm	Power W	Voltage V	Weight Kg
N-3 L	2200850	3	11,5 12,5 20°	43 34 47	1700	230	18
N-8 L	2200851	8,2	14 20 30	50 44 53	1800	220	33
N-13 L	2200852	13	18 22,5 36	55 50 70	1800	230	38
N-22 L	2200854	22	15,5 27,5 50	61 60 89	3000	230	58
N-80 L	2200855	80	48 40 40	157 94 98	7500	400 / 3 N	170

Supplied complete with a refractory ceramic tray as a base and support for material to be assayed.



Clamps and Stands 166 to 177 pages **Spatulas, scoops and vessels** 178 to 182 pages Stainless steel. Series Labor-Hospital + pages 183 to 184 Racks, baskets and storage for tubes, pipettes and flasks 185 to 194 pages **Flakes ice maker** 190 page **Burners and hot air gun** 195 to 198 pages **Diverse laboratory instruments** 199 to 211 pages



Clamps and Stands

CLAMPS MADE OF LIGHT WEIGHT "DURAL CROM" ALLOY EXTRAORDINARILY STRONG.
TONGS MADE OF AISI 304 STAINLESS STEEL.
RETORT STANDS - LABORATORY JACKS AND SCAFFOLDING.







BOSSHEAD

Chrome plated. For Ø from 5 to 17 mm.

Part No. 1000069



BOSSHEAD "MINI"

Chrome plated. For \emptyset from 5 to 12 mm. Part No. 1001460




BOSSHEAD "SUNTECH"

For retort bars of 5 to 22 mm \emptyset .

With screw clamp joints. Ideal for supporting overhead stirrers.

Part No. 5900600



BOSSHEAD "PASS"

Original Phywe.

For joining and clamping multipurpose, round or rectangular rods up to 12 mm across.

For cross T and parallel arrangements.

Part No. 7204055



BOSSHEAD XL

Maximum opening. For rods 13 to 32 mm \emptyset . Part No. **5057022**



BOSSHEAD "MAXIM"

For large diameters.

Chrome plated with one opening to hold rods up to 30 mm \emptyset . Second opening up to 20 mm \emptyset .

Part No. **7000408**



BOSSHEAD "GIROFIX"

Conforms to DIN standard 12895-E16.

Chrome plated with swivel toothed clamp for holding rods at any angle. For rod \emptyset from 5 to 15mm.

Part No. 7000531



BOSSHEAD "MEDI"

For Ø from 5 to 20 mm.

Part No. 1001461 Chrome plated. Part No. 1001462 Anodised black.



OIRIGINAL BOSSHEAD "TREVINGE"

Made of chrome plated "Dural" alloy.

Due to its strength, it is suitable for multilevel and three dimension constructions.

Part No. **7008133** For Ø of 9 to 13 mm. Part No. **7008253** For Ø of 19 to 25 mm.



BOSSHEAD FOR THERMOMETERS

Thermometer support for water and oil baths with locking screw (up to 40 mm). Thermometer rod diameters up to 7mm Ø.

Part No. **6000896**



THERMOMETER CLAMP

Swivel bosshead and press lock, made of AISI 304 stainless steel. For Ø from 6-15 mm. Part No. 6000036



THERMOMETER CLAMP


Hook at one end for a thermometer. Chrome plated.

Part No. 1000037





Model burette stand Part no. 6000117 and double burette part no. 7000181



BURETTE STAND

With AISI 304 stainless steel rod of 10 mm Ø x 600 mm high. Weight 2 Kg. Made specifically for the "DUPLEX" burette clamp with two plates at the base, black and white. Part No. 6000117

DOUBLE BURETTE CLAMP "DUPLEX"

Self adjusting clamp for 2 burettes, with central bosshead screw. Made of Chrome plated metal with PVC coated jaws. Enables complete visibility of the burette scale while being held. Part No. 7000181



BURETTE CLAMP "TRIPS"

Made of chrome plated metal for \emptyset from 4 to 20 mm. Part No. **7000521**



BOSSHEAD BURETTE AND TUBE CLAMP "TRIPS"

Made of chrome plated metal for \emptyset from 4 to 20 mm.

Part No. 7005211



SPRING CLAMP FOR BURETTES AND THERMOMETERS "AUTOMATIC"

Made of chrome plated metal with spring clamp jaw. For \emptyset from 0 to 25 mm. Designed for test tubes.

Part No. **7000025**



BOSSHEAD WITH "AUTOMATIC" CLAMP FOR BURETTES AND THERMOMETERS

Made of chrome plated metal with spring clamp jaw. For Ø from 0 to 25 mm.

Part No. 7000251



CLAMP FOR BURETTES AND TUBES

Made of chrome plated metal with two prong jaws coated with PVC. For Ø from 7 to 25 mm.

Part No. 7000130



CLAMP FOR BURETTES AND CONDENSERS

Made of chrome plated metal with PVC coated jaws. For \emptyset from 12 to 45 mm.


Part No. 7000131



CLAMP "TRIPLEX"

Three pronged PVC coated clamp made of chrome plated metal. For \emptyset from 0 to 35 mm.

Part No. 7000160



UNIVERSAL CLAMP "MINIMAX CROM"

Chrome plated with cork protected jaws. For \emptyset from 0 to 80 mm.

Part No. 1001463



BOSSHEAD WITH BURETTE AND TUBE CLAMP

Made of chrome plated metal with PVC coated double prong jaws. For \emptyset from 7 to 25 mm.

Part No. 7000140



SWIVEL CLAMP BOSSHEAD FOR CONDENSERS

Swivel clamp with double PVC lined jaws For \emptyset from 12 to 45 mm. Part No. 7000141



CLAMP "TRIPLEX"

Three pronged PVC coated clamp made of chrome plated metal. For \emptyset from 0 to 60 mm.

Part No. **7000161**



UNIVERSAL CLAMP "MINIMAX"

Made of Dur-aluminium with cork lined jaws.

For \emptyset from 0 to 80 mm.

Part No. 7000403



SWIVEL UNIVERSAL CLAMP "MINIMAX"


For Ø from 0 to 80 mm. Part No. **7000404**



UNIVERSAL CLAMPS"GIROMAX" WITH METAL ROTATION

With cork lined jaws. With a hinge joint of 180° rotation angle. For Ø from 0 to 80 mm.

Part No. 7004031



BEAKER CLAMP "MAXIM"

Made of chrome plated metal with clamp end covered with stabilized glass fibre. For \emptyset from 80 to 180 mm.

Part No. 7000406



CLAMP "MEDIMAX"

With chain and spring clamp made of AISI 304 stainless steel. For \emptyset from 90 to 145 mm.

Part No. 7000405



CONICAL JOINT CLIPS

Made from AISI 304 stainless steel. Part No.

1000716 joint 7.5/16 1001019 10/19 1001221 12/21 1001423 14/23 1001926 19/26 1002429 24/29 29/32 1002932 1004540 45/40

GROUND ROUND JOINT CLIPS

Made of chromed steel.

Part No.

 1040913
 joint
 round
 13 mm.

 1040919
 " 19 mm.

 1040929
 " 29 mm.

 1040935
 " 35 mm.



CONICAL JOINT CLIPS

Made of PMP (polymethyl pentane).

Part No. joint colour 1001095 10/19 black 1001096 14/23 yellow 1001097 19/26 blue 1001098 24/29 green 1001099 29/32 red 1001100 34/35 orange 45/40 1001101 brown



GROUND ROUND JOINT CLIPS

Made of AISI 304 stainless steel with ratchet fixed pin.

Part No. 1001720 round joint 19 mm. Part No. 1001464 round joint 29 mm. Part No. 1001465 round joint 35 mm.



TEST TUBE TONGS

Made of nickel plated steel. For tubes up to 30 mm Ø. Total length 145 mm. Part No. 1000285





CRUCIBLE TONGS

Original "J.P.Selecta" design with opposite equilibrium.

Made of Electro-bright finished AISI 304 stainless steel.

Part No. 1001330 Total length 400 mm.



CRUCIBLE TONGS

With bow, curved tips and flat hinge.

Electro-bright AISI 304 stainless steel finish.
Part No. 1010420 Total length 200 mm.
Part No. 1010430 Total length 300 mm.

Part No. 1010440 Total length 400 mm.
Part No. 1010450 Total length 500 mm.



CRUCIBLE TONGS

With bow, curved tips.

With thermally protected plastic coated handles.

Part No. 1001590 Total length 220 mm. Part No. 1001591 Total length 330 mm.




STERILIZATION TONGS

Made of AISI 304 electro-bright finished stainless steel. Total length 200 mm.

Part No. 1002801 Straight end. Part No. 1002802 Curved end.

Used for bandages, syringes, instruments, test tubes etc.



TONGS FOR FLASKS AND TUBES

Electro-bright finished AISI 304 stainless steel. For diameters from 15 to 60 mm.

Tong tips covered with stabilized glass fibre.

Part No. 1028130 Total length 300 mm.
Part No. 1028140 Total length 400 mm.



TONGS FOR FLASKS AND TEST TUBES


Sterilization Tongs. Made of AISI 304 electro-bright finished stainless steel. For diameters from 25 to 60 mm. Tong ends covered with stabilized glass fibre.

Part No. 1028230 Total length 300 mm.
Part No. 1028240 Total length 400 mm.



DISH TONGS

Made of eletro-bright finished AISI 304 stainless steel.
Part No. 1000283 Total length 250 mm.



BEAKER TONGS

Made of AISI 304 electro-bright finished stainless steel. Made for beakers from 150 to 1000 ml. For diameters from 60 to 100 mm.

Tong tips covered with stabilised glass fibre. Part No. 1000284 Total length 450 mm.



BEAKER TONGS "PROTECVAS"

Clamps beakers up to 2000 ml. With a span from 50 to 120 mm Ø. Jaws protected with rubber pads. Total length 300 mm. Part No. 7000532



BEAKER TONGS

Original "J.P.Selecta" design with opposite equilibrium. Made of Electrobright finished AISI 304 stainless steel. For beakers from 150 to 1000 ml.

For diameters from 60 to 100 mm.

Tong ends covered with a glass fibre sleeve.

Part No. 1001331 Total length 450 mm.



LABORATORY SCISSORS

AISI 410 stainless steel.

1 Scissors universal "Multi purpose" with a serrated/ dented blade. Suitable for cutting hard plastic, cardboard and soft wire etc.

Part No. 1001468 Cut length 50 mm. Total length 180 mm.

2 Long laboratory scissors.

Part No. 1001469 Cut length 100 mm. Total length 250 mm.

3 Straight multi purpose scissors. General use and bandage dressing. Part No. 1001470 Cut length 50 mm. Total length 140 mm.

4 Curved multi purpose scissors. General use and bandage dressing. Part No. 1000806 Cut length 40 mm. Total length 145 mm.


5 Dissection scissors

Part No. 1000807 Cut length 35 mm. Total length 125 mm.

6 Small dissection scissors

Part No. 1001471 Cut length 25 mm. Total length 110 mm.

7 Flat serrated scissors ideal for laminates, negatives, chromatography, etc. Part No. 1001472 Flat Length 40 mm. Total length 160 mm.



HOFFMAN CLAMPS / MOHR CLAMPS

1 Hoffman clip dual, open-close.

Part No. 1000063 Maximum opening: 15 mm.

2 Hoffman clip closed

Part No. 1000662 Maximum opening: 15 mm.

3 Hoffman clip open

Part No. 1000661 Maximum opening: 15 mm.

4 Hoffman clip dual medium, open-closed.
Part No. 1001467 Maximum opening: 25 mm.

5 Hoffman clip large, open and close.

Part No. 1000064 Maximum opening: 30 mm.

6 Mohr clamp, stainless steel AISI 304 electro bright finish.

Part No. 1000067 Total length: 70 mm.

7 Mohr clamp, nickel wire.

Part No. 1001337 Total length: 70 mm.



FORCEPS FOR LIFTING SLIDES

Made of AISI 304 stainless steel, with flat jaws. Length: 105 mm. Suitable for lifting, slides, watch glass or small flat object.

Part No. 1000105



DISSECTING FORCEPS

1 Dissecting Straight forceps round tips, made of AlSl304 stainless steel.
Part No. 1000580 Total length 140 mm.

2 Dissecting Straight forceps fine tips, made of AlSl304 stainless steel. Part No. 1000801 Total length 140 mm.

3 Dissecting curved forceps round tips, made of AISI 304 stainless steel.
Part No. 1000865 Total length 140 mm.

4 Dissecting curved forceps fine tips, made of AISI304 stainless steel. Part No. 1000864 Total length 140 mm.

MAXIMUM OPENING

5 Sprung straight forceps with curved fine tips,made of AISI 304 stainless steel. Part No. 1001190 Total length 160 mm.

6 Sprung straight forceps with fine tips, made of AISI 304 stainless steel. Part No. 1001189 Total length 160 mm.

7 Curved point, AISI 304 stainless steel forceps. Part No. 6000128 Total length 160 mm.

8 Fine point, AISI 304 stainless steel forceps. Part No. 1001191 Total length 120 mm.



HEAT AND SPILLS PROTECTION PAD GLOVE

Made of flexible silicon that fits all hand sizes.

Part No. 1001255





BOSSHEAD "P" TYPE

Nickel alloy bosshead. For fixing perpendicular structures. For 12 mm Ø bars. Part No. 1000011

ACCESSORY

Allen hex key. Suitable for "P" type bosshead. Part No. 1000111



ORIGINAL BOSSHEAD "TREVINGE"

Made of chrome plated "Dural" alloy.

Due to its strength, suitable for multilevel and three dimension constructions.

Part No. **7008133** For Ø of 9 to 13 mm. Part No. **7008253** For Ø of 19 to 25 mm.



RODS AND TUBES

12 mm Ø bars and rods. To be used with the "TREVIN" and "P" type bossheads.

			,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,		
Dur-aluminiu	ım rods.	AISI 304 stainless steel tube.			
Part No.	Length	Part No.	Length		
6000026	60 mm.	6000016	60 mm.		
6000220	200 mm.	6000120	200 mm.		
6000235	350 mm.	6000135	350 mm.		
6000270	700 mm.	6000170	700 mm.		
6002100	1000 mm.	6001100	1000 mm.		
6002150	1500 mm.	6001150	1500 mm.		
6002200	2000 mm.	6001200	2000 mm.		
6002300	3000 mm.	6001300 (rod)	3000 mm.		



UNIVERSAL MOUNTING BASE

Bench top multiple rod support.

Suitable for mounting and supporting vertical and horizontal bars.

Part No. 1000050



MOUNTING FOOT SUPPORT FOR RODS

For rod Ø of 12 mm.

Part No. 1000101 Rod position horizontal.
Part No. 1000049 Rod position vertical.



SIMULTANEOUS FOOT SUPPORT

For rod Ø of 12 mm.

Simultaneous positions, horizontal and vertical. Part No. 1001102



BENCH CLAMP "PASS"

Original Phywe. Clamps to the bench while holding from 4 to 12 mm diameter rod. Bench top width up to 40mm. Part No. 1394855

SUPPORT STAND KITS



700 high x 700 mm long.

Made up of:

8 rods 700 mm long.

2 rods 350 mm long.

16 "P" type bossheads.

4 foot supports.

1 Allen key.

Part No. 6000090 Dur-aluminium

Part No. 6000190 AISI 304 stainless steel



1000 high x 1000 mm long.

Made up of:

10 rods of 1000mm long.

2 rods 700 mm long.

3 rods 350 mm long.


34 "P" type bossheads.

6 foot supports.

1 Allen key.

Part No. 6000902 Dur-aluminium

Part No. 6001902 AISI 304 stainless steel



700 high x 1000 mm long.

Made up of:

4 rods 1000 mm long.


5 rods 750 mm long.

2 rods 350 mm long. 20 "P" type bossheads.

4 foot supports.

1 Allen key.

Part No. 6000901 Dur-aluminium Part No. 6001901 AISI 304 stainless steel



1000 high x 2000 mm long.

Made up of:

5 rods 2000 mm long.

7 rods 1000 mm long.

8 rods 200 mm long.

2 rods 350 mm long.

45 "P" type bossheads.

8 foot supports.

1 Allen key.

Part No. 6000903 Dur-aluminium

Part No. 6001903 AISI 304 stainless steel



SUPPORT "STAND"

Solid base, three vertical rods of 800 mm in height, two horizontal rods of 500 mm in length, six "P" type bossheads. Made of Dur-aluminium.

Weight: 6.5 Kg. Part No. 6000482



CENTRAL ROD RETORT STAND

Vertical rod of Dur-aluminium of 10 Ø x 600mm high. Weight: 3 Kg. Base dimensions 315 x 200 mm. Part No. 6000046



TRIPOD SUPPORT "PASS"

Original Phywe.

To support square and round rods from 6 to 14 mm dia. Weight: 2 Kg.

Base with levelling screws.

Part No. 6200255



DOUBLE BASE SUPPORT T "MINI"

Small foot support with screw clamp for removable rods from 8 to 12 mm \emptyset .

Weight: 1.2 Kg. Part No. 6001103



STANDARD RODS

Made from AISI 340 stanless steel. Suitable for double "T MINI" and Tripod "PASS"

Part No. 6001160 8 Ø x 500 mm length Part No. 6000534 10 Ø x 600 mm length Part No. 6000535 12 Ø x 700 mm length



U SHAPED STAND

AISI 304 stainless steel of 20 Ø x 800 mm high. Very stable. Weight: 7 Kg. Part No. 6000411



DOUBLE BASE SUPPORT T "MAXI"

Rod of AISI 304 stainless steel of 20 Ø x 800 mm high. Designed to support weights up to 20 Kg at 300 mm height from the base. Useful for stirrers. Weight: 7 Kg.

Part No. **6000410**



DOUBLE BASE SUPPORT T TELESCOPE

Stainless steel base AISI 304 of 32 mm Ø.
Telescopic elevation from 725 to 1025 mm.
Maximum support weight 20 Kg. Actual weight: 7.7 Kg
Compatible with Part No. 5057022 (See page 166).
Part No. 5057011



TRIPOD SUPPORT STAND

AISI 304 stainless steel rod.

Part No.	Ø x Height	Weight Kg
6002278	8 x 500	1.4
6022710	10 x 600	1.8
6022712	12 x 700	2.5



Rods made of Dur-aluminium

Part No.	Ø x Height	Weight Kg
6000798	8 x 500	1.2
6000799	10 x 600	1.5
6007912	12 x 700	2

UNIVERSAL ARM SUPPORT

With holder for up to 4 sensors or electrodes. Very stable heavy base.

Ideal for working at varying heights. Part No. 1001552



RETORT STAND

Conforms to the DIN 12892. Standard.

Rods made of AISI 304 stainless steel

	Rod	Plate size	
Part No.	Ø x height	length x width	Weight Kg
6002268	8 x 500 mm.	200 x 125 mm.	2
6022610	10 x 600 mm.	250 x 160 mm.	3
6022612	12 x 700 mm.	315 x 200 mm.	3.6
6022613	12 x 800 mm.	315 x 200 mm.	3.8

Rods of Dur-aluminium

	Rod	Plate size	
Part No.	Ø x height	length x width	Weight Kg
6000778	8 x 500 mm.	200 x 125 mm.	1.4
6000779	10 x 600 mm.	250 x 160 mm.	2.4
6007712	12 x 700 mm.	315 x 200 mm.	3
6077812	12 x 800 mm.	315 x 200 mm.	3.4



TRIANGLE SUPPORTS


Made of 4 mm Ø $\,$ AISI 304 stainless steel with electrobright finish.

Part No.	Side	Part No.	Side
1002874	40 mm.	1002877	70 mm.
1002875	50 mm.	1002878	80 mm.
1002876	60 mm.	1028710	100 mm.



REFRACTORY TUBE AND WIRE TRIANGLE SUPPORTS

Part No.	Side
1000825	50 mm.
1000826	60 mm.
1000827	75 mm.



TRIPOD STAND "TRISOP"

Made of galvanised steel.
Part No. Side Total height
1000524 130 mm. 200 mm.



TRIPOD SUPPORTS WITH ROUND TOP

	stainless	

Part No.	Ø ext. ring	lotal height
6002868	80 mm.	210 mm.
6028610	100 mm.	210 mm.
6028612	120 mm.	210 mm.
6028614	140 mm.	210 mm.

Nickel plated steel

illokol piutou otool			
Part No.	Ø ext. ring	Total height	
1000818	80 mm.	210 mm.	
1008110	100 mm.	210 mm.	
1008112	120 mm.	210 mm.	
1008114	140 mm.	210 mm.	
1008116	160 mm.	210 mm.	



TRIPOD STAND WITH RING TOP.

Galvanised steel.

Part No.	Ø ext. ring	Total height
7000808	80 mm.	210 mm.
7008010	100 mm.	210 mm.
7008012	120 mm.	210 mm.



TRIPOD STAND WITH RING TOP, ADJUSTABLE HEIGHT

Galvanised steel.

Part No. Ø ext. ring Adjustable height 7022810 100 mm. 170 to 250 mm. 7022812 140 mm. 170 to 250 mm. 7022814 180 mm. 170 to 250 mm.



1. DECANTING FLASK SUPPORT

Made of polypropylene. With support flasks up to 12 cm \emptyset . Part No. 1001346

2. DECANTING FLASK SUPPORT WITH BOSSHEAD

Made of polypropylene. With support flasks from 30 to 150 mm \emptyset . Complete with reduction ring. Part No. 1001347



CIRCULAR SUPPORT PLATFORMS WITH BOSSHEAD

To support flat bottom vessels and instruments. Chrome plated bosshead.

Part No. Plate Ø 6000879 80 mm. 150 mm.



OPEN RING BOSSHEAD SUPPORTS

Made of galvanised steel.


Part No.	Ring exterior Ø	Usable length A
7000506	60 mm.	100 mm.
7000508	80 mm.	130 mm.
7005010	100 mm.	160 mm.
7005012	120 mm.	170 mm.
7005014	140 mm.	180 mm.



CLOSED RING BOSSHEAD SUPPORTS

Made of galvanised steel.


Part No.	Ring exterior Ø	Usable length A
7000536	60 mm.	100 mm.
7000538	80 mm.	130 mm.
7005310	100 mm	160 mm.
7005312	120 mm.	170 mm.
7005314	140 mm	180 mm



OPEN RING RETORT SUPPORTS

Made of galvanised steel.

Part No.	Ring exterior Ø	Usable length A
1000206	60 mm.	180 mm.
1000208	80 mm.	195 mm.
1002010	100 mm.	200 mm.
1002012	120 mm.	215 mm.
1002014	140 mm.	225 mm.



CLOSED RING RETORT SUPPORTS

Made from galvanised steel.

Part No.	Ring exterior Ø	Usable length A
1000236	60 mm.	180 mm.
1000238	80 mm.	195 mm.
1002310	100 mm.	200 mm.
1002312	120 mm.	215 mm.
1002314	140 mm	225 mm



GAUZE IRON WIRE WITH CERAMIC FIBRE

Part No.

1087102	100 mm Side
1087122	120 mm Side
1087142	140 mm Side
1087162	160 mm Side
1087182	180 mm Side
1087202	200 mm Side



Heat conductivity Temperature °C 19 18 17 16 15 NEOCERAM 0 13 12 11 10 100 200 300 400 500 600



CERAMIC PLATES "NEOCERAM"

Resistant to temperatures from -200 °C to 700 °C.

Part No.

 1000684
 120 x 120 x 5 mm.

 1000685
 140 x 140 x 5 mm.

 1000686
 160 x 160 x 5 mm.

Resistant to temperatures from -200 °C to 700 °C. Excellent heat conductor.

20% better heat transfer compared to a conventional gauze ceramic.

High chemical resistance to aggressive compounds. Easy to clean and long lasting.

TRIPOD STAND FOR CERAMIC PLATES

With adjustable sides to prevent the plate from moving or tilting.

Total height 210 mm.

Part No. For plates of 1000687 120 x 120 mm. 1000688 140 x 140 mm. 1000689 160 x 160 mm.



LABORATORY JACK "LAB-ELE"

With anodised aluminium base.

Part No. 1000294 Usable platform: 115 x 140 mm.

Elevated height: 260 mm. Maximum support weight: 1 kg.

Top plate: 200 x 200 mm. Part No. 7000295

Maximum elevated height: 260 mm. Max. support weight: 1.7 kg. With hole and clamp screw for a support rod of 10 Ø x 600 mm high.



LABORATORY JACKS "LAB-STEEL"

Made completely of AISI 304 stainless steel.

Robust construction and stable body for supporting heavy weights. Models:

"LAB-STEEL" MINI

Platform Dims.: 145 x 145 mm. Elevation height Part No. 1001104

55 to 155 mm. Maximum support weight: 1.3 Kg.

"LAB-STEEL" MEDI

Part No. 1001161 Platform Dims.: 200 x 200 mm. Elevation height

65 to 270 mm. Maximum support weight: 2 Kg.

"LAB-STEEL" MAXI

Part No. 1001105

Platform Dims.: 250 x 250 mm. Elevation height 65 to 270 mm. Maximum support weight: 3 Kg.



TABLE TOP "ALMANDINE"

Made of polyester mineral resin.

Features:

Excellent chemical resistance, with resistance to class 1 and 2 chemicals. Resists abrasion, boiling water, dry heat, fire and scratches.

Stains easily removed with a detergent or abrasive cleaner.

Applications:

Protects bench top surfaces against aggressive chemicals, burns and scratches. Suitable for working with: Acids, alkali bases, solvents or dyes.

Preparation of horizontal agar gels.

Preparation of petri dishes, etc.

		,	
Part No.	Width	Thickness	Depth
	4.0		4.5

1001090 48 cm 0.6 cm 45 cm Without level With adjustable level supports 1001091 48 cm 0.6 cm 45 cm



P.V.C coated steel.

Applications:

The weight of the ring stabilises and prevents flasks from floating when placed in a water bath. Suitable for Conical, volumetric and all types of flasks.

Part No.	Colour	Mode	Ø interior mm.	Weight g.
1001169	red	open	50	207
1001170	red	open	70	340
1001171	red	open	100	470
1001734	blue	closed	48	270
1001735	blue	closed	70	360
1001736	blue	closed	90	440

Ring support holder, made of stainless steel.

Part No. 1001737



Spatulas, scoops and vessels.

STAINLESS STEEL AND NICKEL EVAPORATING BASINS.
STAINLESS STEEL, NICKEL AND ZIRCONIUM CRUCIBLES.





MICRO SPATULA WITH FLAT ENDS

Made od AISI 304 stainless steel.
Part No. Length Width
1027680 80 mm. 4 mm.
1276120 120 mm. 4 mm.
1276150 150 mm. 4 mm.
1276210 210 mm. 4 mm.




MICRO-SPATULA CURVED ENDS

Made of stainless steel AISI 304.
Part No. 1001581 Length 150 mm.



MICRO SPATULA WITH FLAT AND SPOON ENDS

Made of AISI 304 stainless steel.
Part No. Length Width
1277100 100 mm. 4 mm.
1277150 150 mm. 4 mm.
1277210 210 mm. 4 mm.





SPATULA FLAT AND GROOVED ENDS

Made of Stainless steel AISI 304.

Part No. 1000855 Total length 180 mm. Flat width 9 mm. Grooved width 8 mm.



SPATULA WITH CURVED ENDS

Made of stainless steel AISI 304. Part No. Length Width 1010313 130 mm. 11 mm. 1001518 150 mm. 11 mm. 1010318 180 mm. 11 mm. 1010321 210 mm. 11 mm.



SPATULA WITH SPOON AND FLAT ENDS

Part No. Length Width
1096120 120 mm. 18 mm. stainless steel AISI 304.
1096150 150 mm. 20 mm. stainless steel AISI 304.
1196120 120 mm. 18 mm. Pure nickel.
1196150 150 mm. 20 mm. Pure nickel.



CHISEL AND GROOVED SPOON END SPATULA

Part No. Length Width
6021312 120 mm. 18 mm. stainless steel AISI 304.
6021315 150 mm. 20 mm. stainless steel AISI 304.
6121312 120 mm. 18 mm. Pure nickel.
6121315 150 mm. 20 mm. Pure nickel.



FLAT SPATULA WITH SPOON

Made of Stainless steel AISI 304. Part No. 1001583 Length 180 mm.



SPATULA DOUBLE SPOON ENDS

Made of Stainless steel AISI 304.
Part No. 1001582 Length 150 mm.



DOUBLE SPATULA SEMI-MICRO

Made of Stainless steel AISI 304. Total length 140 mm Rectangular end: 35 long x 7 mm wide. Oval end: 35 long x 7mm wide. Part No. 1278140

4 mm 2 55 mm 2 7 mm 2 7 mm 3 20 mm 6 mm 40 mm 40 mm 10 mm 40 mm 21 mm 35 mm 38 mm

SPATULA WITH POCKET CLIP

Part No.

1. 1001585 Micro-spatula flat spoon end. Length 150 mm.

2. 1001586 Micro-spatula curved end. Length 150 mm.


3. 1001587 Spatula flat with grooved end. Length 150 mm.

4. 1001588 Spatula double end flat. Length 150 mm.

5. 1001589 Spatula flat and spoon ends. Length 150 mm.



Specially designed for lab coats with upper pocket.



FLAT AND CURVED END SPATULA

Made of AISI 304. stainless steel.

Part No.	Length	Width	
6099120	120 mm.	20 mm.	
6099150	150 mm.	24 mm.	
6099180	180 mm.	24 mm.	
6099210	210 mm.	24 mm.	



SPATULA DOUBLE CURVED ENDS

Made of AISI 304. stainless steel.

Part No.	Length	Width
6100120	120 mm.	15 mm.
6100150	150 mm.	25 mm.
6100180	180 mm.	25 mm.
6100210	210 mm.	25 mm.



GROOVED SPATULA

Grooved spoon, made of AISI 304 stainless steel.

Total length: 190 mm.

Designed to remove residues from flasks and bottles due to the flat sharp end and curved spoon.

Part No. 1096190



SPATULA SPOON FLAT-CONCAVE

Made Stainless steel AISI 304.

Part No. 1001584 Length 250 mm.



PASTE KNIFE WITH WOODEN HANDLE

Blade made from tempered Stainless steel.

Part No.	Blade length	Blade width	Part No.	Blade length	Blade width	
1010280	80 mm.	14 mm.	1102160	160 mm.	24 mm.	
1102100	100 mm.	16 mm.	1102200	200 mm.	29 mm.	
1102120	120 mm.	18 mm.	1102240	240 mm.	34 mm.	
1102140	140 mm.	22 mm.	1102280	280 mm.	34 mm.	




FLAT SPATULA WITH SPOON

Analytical type, made of AISI 304 Stainless steel. Total length: 210 mm. Rectangular end: 40 long x 9 mm wide.

Spoon end: 30 long x 18 mm wide.

Part No. 1095210



FLEXIBLE SPATULA

With stainless steel AISI 304 blade and handle.

Part No.	Total length	Blade length	Blade width
1000856	160 mm.	52 mm.	11 mm.
1000857	160 mm.	52 mm.	15 mm.
1000858	170 mm	60 mm	20 mm



MULTI PURPOSE KNIFE

With a black bakelite® handle. Designed for comfort and safety Made of tempered stainless steel.

Total length: 230 mm. Blade length: 95 mm. Part No. 1000862



SPATULA POCKET POUCH

Part No. **1001530**

Contents:

- 1 Flexible spatula. Part No. 1000856 . 1 Semi micro spatula. Part No. 1278140
- 1 Micro spatula flat spoon. Part No. 1277150. 1 Dissection scissors. Part No. 1000807
- 1 Forceps blunt end. Part No. 1000580



SODIUM SPOON

With brass basket. Ø 25 mm. Total length: 190 mm.

Part No. 1000805




POWDER SCOOPS

Length: 100 mm. Width: 26 mm.

Part No. 6000278 Stainless steel AISI 304.

Part No. 6001278 Pure nickel.



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Part No. 1001592 Medium Total length 320 mm.
Part No. 1001593 Big Total length 360 mm.





VIBRATING SPATULA AND BALANCE WEIGHT FORCEPS

Vibrating spatula.

With manual vibrating device. Made of AISI 304 stainless steel.

Part No. 1000413 Curve end. Total length 176 mm. Width 12 mm.
Part No. 1000412 Grooved end. Total length 176 mm. Width 12 mm.

Balance weight forceps.

Made of AISI 304 stainless, electrobright finish. Total length 93 mm.

Part No. 1000217

Forceps for balance weights.

Made of AISI 304 stainless, electrobright finish. Total length 100 mm.

Part No. 1000216



FUSION SPOONS

Made of AISI 304 stainless, wooden handle. Total length: 350 mm.

Flat base. Volume 90 ml. Part No. **6000093**

Round base. Volume 75 ml.

Part No. 6000094



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Made of etailion eteel Alei ee i.					
Part No.	Total length Usable length		Width	Capacity	
1000859	180 mm.	100 mm.	50 mm.	120 ml.	
1000860	260 mm.	150 mm.	75 mm.	450 ml.	
1000861	310 mm.	200 mm.	100 mm.	1500 ml.	
1001354	330 mm.	210 mm.	150 mm.	2000 ml.	



FUNNEL WITH HANDLE

Made of AISI304 stainless steel.

Part No. 1000906 100 mm Ø x 100 mm high x 10 mm exit.

Part No. 1000907 120 mm Ø x 112 mm high x 10 mm exit.

Part No. 1000908 150 mm Ø x 159 mm high x 20 mm exit.

Part No. 1000909 200 mm Ø x 183 mm high x 20 mm exit.



FLAT BOTTOM BASINS

Stainless steel AISI 304	Nickel Ni			
Part No.	Part No.	Ø	Height	Volume
1009160	1019160	60 mm.	25 mm.	40 ml.
1009170	1019170	70 mm.	25 mm.	60 ml.
1009180	1019180	80 mm.	30 mm.	90 ml.
1091100	1191100	100 mm.	35 mm.	180 ml.



PETRI DISH WITH LID

Stainless steel AISI 304	Nickel Ni			
Part No.	Part No.	Ø	Height	Volume
1002126	1012126	60 mm.	20 mm.	50 ml.
1002128	1012128	80 mm.	20 mm.	90 ml.
1021210	-	100 mm.	25 mm.	150 ml.
1021212	-	120 mm.	25 mm.	260 ml.
1021215	_	150 mm.	30 mm.	550 ml.



SPOON WITH HANDLE

Made of AISI304 stainless steel.

Part No. 1001594 Capacity 100 ml. Total length 240 mm.
Part No. 1001595 Capacity 200 ml. Total length 280 mm.
Part No. 1001596 Capacity 300 ml. Total length 360 mm.



EVAPORATION BASINS, ROUND BOTTOM

Stainless steel AISI 304	Nickel Ni			
Part No.	Part No.	Ø	Height	Volume
1009260	1019260	60 mm.	25 mm.	35 ml.
1009270	1019270	70 mm.	25 mm.	50 ml.
1009280	_	80 mm.	30 mm.	75 ml.
1092100	_	100 mm.	35 mm.	125 ml.



BUCKETS WITH HANDLES

Without graduation.


Part No. Height Ø Capacity 1001599 270 mm 270 mm 10 litres 1001600 300 mm 300 mm 15 litres

BUCKET LIDS

Part No. 1001601 Suitable for buckets 10 litres Part No. 1001599
Part No. 1001602 Suitable for buckets 15 Litres Part No. 1001600



DYE BEAKER
Stainless steel AISI 304.
Volume: 500 ml.
Part No. 6000291





CRUCIBLES

CRUCIBLES					
AISI 304 STA	AINLESS S	TEEL			Lids
Part No.	Ø	Height	Thickness	Volume	Part No.
1098302	30 mm.	32 mm.	2 mm.	20 ml.	6009730
1098402	40 mm.	42 mm.	2 mm.	45 ml.	6009740
1009830	30 mm	32 mm.	0.6 mm.	20 ml.	6009730
1009840	40 mm.	42 mm.	0.6 mm.	45 ml.	6009740
1009850	50 mm.	55 mm.	0.8 mm.	85 ml.	6009750
1009860	60 mm.	65 mm.	0.8 mm.	140 ml.	6009760
BURE NIOW	-1 11.				
PURE NICKE		Hataka	Third	W-1	Lids
Part No.	Ø	Height	Thickness	Volume	Part No.
1198302	30 mm.	32 mm.	2 mm.	20 ml.	6019730
1198402	40 mm.	42 mm.	2 mm.	45 ml.	6019740
1198502	50 mm.	55 mm.	2 mm.	85 ml.	6019750
1198602	60 mm.	65 mm.	2 mm.	140 ml.	6019760
1019830	30 mm.	32 mm.	0.6 mm.	20 ml.	6019730
1019840	40 mm.	42 mm.	0.6 mm.	45 ml.	6019740
1019850	50 mm.	55 mm.	0.6 mm.	85 ml.	6019750
1019860	60 mm.	65 mm.	0.6 mm.	140 ml.	6019760
ZIRCONIUM	l Zr				Lids
ZIRCONIUM Part No.	I Zr Ø	Height	Thickness	Volume	Lids Part No.
		Height 30 mm.	Thickness 1 mm.	Volume 20 ml.	
Part No.	Ø	•			Part No.
Part No. 1000563	Ø 33 mm.	30 mm.	1 mm.	20 ml.	Part No. 6000567
Part No. 1000563 1000564	Ø 33 mm. 47 mm.	30 mm. 43 mm.	1 mm. 1 mm.	20 ml. 50 ml.	Part No. 6000567 6000568
Part No. 1000563 1000564 1000565 1000566	Ø 33 mm. 47 mm. 52 mm. 59 mm.	30 mm. 43 mm. 48 mm.	1 mm. 1 mm. 1 mm.	20 ml. 50 ml. 75 ml.	Part No. 6000567 6000568 6000569 6000570
Part No. 1000563 1000564 1000565 1000566	Ø 33 mm. 47 mm. 52 mm. 59 mm.	30 mm. 43 mm. 48 mm. 51 mm.	1 mm. 1 mm. 1 mm. 1 mm.	20 ml. 50 ml. 75 ml. 100 ml.	Part No. 6000567 6000568 6000569 6000570 Lids
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No.	Ø 33 mm. 47 mm. 52 mm. 59 mm.	30 mm. 43 mm. 48 mm. 51 mm.	1 mm. 1 mm. 1 mm. 1 mm.	20 ml. 50 ml. 75 ml. 100 ml.	Part No. 6000567 6000568 6000569 6000570 Lids Part No.
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml.	Part No. 6000567 6000568 6000569 6000570 Lids Part No. 1000620
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571	Ø 33 mm. 47 mm. 52 mm. 59 mm. 48S Ø 45 mm. 50 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml.	Part No. 6000567 6000568 6000569 6000570 Lids Part No. 1000620 1000621
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572	Ø 33 mm. 47 mm. 52 mm. 59 mm. 4SS Ø 45 mm. 50 mm. 60 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml.	Part No. 6000567 6000568 6000570 Lids Part No. 1000620 1000621 1000622
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml.	Part No. 6000567 6000569 6000570 Lids Part No. 1000620 1000622 1000623
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573 1000574	Ø 33 mm. 47 mm. 52 mm. 59 mm. 48SS Ø 45 mm. 50 mm. 60 mm. 35 mm. 40 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm. 50 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml. 40 ml.	Part No. 6000567 6000569 6000570 Lids Part No. 1000620 1000621 1000623 1000624
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml.	Part No. 6000567 6000569 6000570 Lids Part No. 1000620 1000622 1000623
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573 1000574	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm. 40 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm. 50 mm. 62 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml. 40 ml.	Part No. 6000567 6000569 6000570 Lids Part No. 1000620 1000621 1000623 1000624
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573 1000574	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm. 40 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm. 50 mm. 62 mm.	1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml. 40 ml.	Part No. 6000568 6000569 6000570 Lids Part No. 1000620 1000621 1000622 1000623 1000624 1000625
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573 1000574 1000575 QUARTZ GLA	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm. 40 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm. 50 mm. 62 mm.	1 mm. 1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm.	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml. 40 ml. 75 ml.	Part No. 6000568 6000569 6000570 Lids Part No. 1000620 1000621 1000622 1000623 1000624 1000625 Lids
Part No. 1000563 1000564 1000565 1000566 QUARTZ GLA Part No. 1000570 1000571 1000572 1000573 1000574 1000575 QUARTZ GLA Part No. 1000575	Ø 33 mm. 47 mm. 52 mm. 59 mm. ASS Ø 45 mm. 50 mm. 60 mm. 35 mm. 40 mm. 55 mm.	30 mm. 43 mm. 48 mm. 51 mm. Height 36 mm. 40 mm. 48 mm. 44 mm. 50 mm. 62 mm.	1 mm. 1 mm. 1 mm. 1 mm. 1 mm. Thickness 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. 1,5 mm. Thickness	20 ml. 50 ml. 75 ml. 100 ml. Volume 38 ml. 49 ml. 90 ml. 24 ml. 40 ml. 75 ml.	Part No. 6000568 6000569 6000570 Lids Part No. 1000620 1000621 1000622 1000623 1000624 1000625 Lids Part No.

STAINLESS STEEL AISI 304- PURE NICKEL - ZIRCONIUM - QUARTZ

STAINLESS STEEL is usable at temperatures from -180 °C to +500 °C. Resistant to nitric acid, nitrate and potassium permanganate. Good resistance to concentrated alkaline solutions. Average resistance to acetic acid, hydrogen sulphide diluted, phosphoric low concentration, and solutions of sulphites, boric acid, organic acids and bromides.

NICKEL unalterable in the air and stainless up to temperatures of 500 °C, good resistance to molten alkali and seawater.

Excellent performance in the presence of wet chlorine or hydrochloric acid gas up to 500 °C. Slowly attacked by organic acids, hydrochloric and sulfuric acids. Low resistance to diluted nitric acid.

THE ZIRCONIUM unalterable in the air and clean, resistant up to 550 °C in an oxidizing atmosphere at 1500 °C in inert gas (argon or nitrogen), in many cases to replace platinum crucibles in all melting and in particular the sodium peroxide and carbonates. About 100 melting per crucible, making them very economical with reference to the others, either porcelain crucibles, steel or nickel. Resistant to most alkali solvents (Na, K, Li), carbonates, hydroxides, peroxides, borates, nitrates, chlorides, some fluorides and sulfuric acid up to 75% concentrate.

Hydrofluoric acid attacks it.

QUARTZ GLASS has a high temperature resistant (liquid glass transition 1130°C) as well as a low thermal expansion. It is considered a unique material for several scientific applications due to its purity (can only be attacked by hot phosphoric acid and hydroflouric acid) and thermal characteristics. Furthermore, quartz glass has extraordinary optical properties through its spectral transparency and high irradiation resistance.

Stainless steel receptacles. QUALITY AISI 304.

LABORATORY-HOSPITAL + SERIES





ROUND BASINS

Part No.	Ø exterior	Ø interior	Height	Volume
8081810	133 mm.	118 mm.	45 mm.	400 ml
8081812	175 mm.	157 mm.	58 mm.	850 ml
8081814	217 mm.	199 mm.	66 mm.	1650 ml
8081815	255 mm.	245 mm.	110 mm.	3400 ml
8081817	318 mm.	297 mm.	117 mm.	7000 ml



RECTANGULAR TRAYS

Conical shape

	Exterior	Exterior	Exterior	Interior	Interior
Part No.	Height	Width	Depth	Width	Depth
8080801	40 mm.	190 mm.	125 mm.	158 mm.	112 mm.
8080802	45 mm.	244 mm.	161 mm.	211 mm.	147 mm.
8080803	50 mm.	305 mm.	195 mm.	270 mm.	177 mm.
8080804	55 mm.	338 mm.	261 mm.	301 mm.	241 mm.

Long shape

	Exterior	Exterior	Exterior	Interior	Interior
Part No.	Height	Width	Depth	Width	Depth
8080208	18 mm.	175 mm.	62 mm.	145 mm.	49 mm
8080207	28 mm.	380 mm.	160 mm.	351 mm.	145 mm



STERILIZATION BOXES WITH LID

Part No. Height Width Depth 8085504 30 mm. 170 mm. 70 mm. 8085505 50 mm. 200 mm. 100 mm.



BOXES WITH ROUND LIDS

Part No.

8085514 140 Ø x 65 mm high. Volume 1000 ml. 8085515 180 Ø x 75 mm high. Volume 1500 ml.



KIDNEY DISHES

Part No.

8085301 32 high x 200 mm long.
8085302 37 high x 250 mm long.
8085303 42 high x 300 mm long.





TALL BEAKERS AND JUGS

Mirror finish stamped beakers and jugs with a curved lip.

Easily cleaned and resistant to large numbers of chemical agents.

1 Standard 2 With Spout 3 With handle and spout

Standard	2 with Spout	3 with nandle and	spout	
Part No.	Part No.	Part No.	Ø x Height	Volume
8002902	8012902	8092902	72 x 84 mm.	250 ml.
8002905	8012905	8092905	86 x 129 mm.	500 ml.
8029010	8129010	8929010	102 x 150 mm.	1000 ml.
8029020	8129020	8929020	132 x 184 mm.	2000 ml.
8029030	8129030	8929030	160 x 180 mm.	3000 ml.
8029040	8129040	8929040	180 x 200 mm.	4000 ml.
8029050	8129050	8929050	200 x 220 mm.	5000 ml.



LIDS FOR REAKERS AND JUGS

JEAKENO AND JUUO		
Beaker size	Part No.	Beaker size
250 ml.	8003000	3000 ml.
500 ml.	8004000	4000 ml.
1000 ml.	8005000	5000 ml.
2000 ml.		
	Beaker size 250 ml. 500 ml. 1000 ml.	Beaker size Part No. 250 ml. 8003000 500 ml. 8004000 1000 ml. 8005000



BUCKET WITH HANDLE

Buckets. Graduated from 1 to 12 litres. 320 \emptyset x 300 mm high. Volume 12 litres.

Part No. 8000252 Bucket lids Part No. 1002521



MEDICINE CUPS

Part No.

8085530 Graduated. Capacity 30 ml. **8085531** Graduated. Capacity 50 ml.

8085532 With handle, no graduations. Capacity 300 ml.



SERILISATION JAR 100 Ø x 155 mm high. Volume 1000 ml. Part No. **8080904**



ROUND BASIN370 Ø x 145 mm high. Volume 9000 ml.
Part No. 8000140



SUPPORT FOR ROUND BASIN
Portable with rubber wheels and brake.
Basin Part No. 8000140.
Part No. 6010140




SUPPORT FOR BUCKETSPortable with rubber wheels and brake.

For bucket Part No. 8000252. Part No. **6010252**



Racks, baskets and storage for tubes pipettes and flasks

TUBE RACKS








Part No.	Made of	Number	For tube	Dimensions ext.		r tube Dimensio	t. mm
		of tubes	Ø up to	Α	В	C	D
1001563	Stainless steel	24 (6 x 4)	13 mm VAC	95	70	70	60



Part No.	Made of	Number	For tube	Dime	nsio	ns e	kt. m	m
		of tubes	Ø up to	Α	В	C	D	
1001212	PVC coated wire	12 (4 x 3)	13 mm	65	50	65	50	
1025500	Stainless steel wire AISI 304	12 (4 x 3)	13 mm	65	50	65	50	



Part No.	Made of	Number	For tube	Dime	nsior	ıs ext	t. mm
		of tubes	Ø up to	Α	В	C	D
1026112	PVC coated wire	12 (6 x 2)	16 mm	120	45	80	65
1026412	PVC coated wire	12 (6 x 2)	25 mm	180	65	110	95
1025512	Stainless steel wire AISI 304	12 (6 x 2)	16 mm	120	45	80	65
1025812	Stainless steel wire AISI 304	12 (6 x 2)	25 mm	180	65	110	95



Made of	Number	For tube	Dime	ensio	ns ex	t. mm
	of tubes	Ø up to	Α	В	C	D
PVC coated wire	24 (8 x 3)	16 mm	160	65	80	65
PVC coated wire	24 (8 x 3)	25 mm	235	95	110	95
Stainless steel wire AISI 304	24 (8 x 3)	16 mm	160	65	80	65
Stainless steel wire AISI 304	24 (8 x 3)	25 mm	235	95	110	95
	PVC coated wire PVC coated wire Stainless steel wire AISI 304	PVC coated wire 24 (8 x 3) PVC coated wire 24 (8 x 3) Stainless steel wire AISI 304 24 (8 x 3)	PVC coated wire 24 (8 x 3) 16 mm PVC coated wire 24 (8 x 3) 25 mm Stainless steel wire AISI 304 24 (8 x 3) 16 mm	PVC coated wire 24 (8 x 3) 16 mm 160 PVC coated wire 24 (8 x 3) 25 mm 235 Stainless steel wire AISI 304 24 (8 x 3) 16 mm 160	Of tubes Ø up to A B PVC coated wire 24 (8 x 3) 16 mm 160 65 PVC coated wire 24 (8 x 3) 25 mm 235 95 Stainless steel wire AISI 304 24 (8 x 3) 16 mm 160 65	PVC coated wire 24 (8 x 3) 16 mm 160 65 80 PVC coated wire 24 (8 x 3) 25 mm 235 95 110 Stainless steel wire AISI 304 24 (8 x 3) 16 mm 160 65 80



Part No.	Made of	Number of tubes	For tube Ø up to		ensio B		t. mm D
		48 (8 x 6) 48 (8 x 6)	16 mm 25 mm 16 mm 25 mm	235 160	125 180 125 180	110 80	95 65




Part No.	Made of	Number	For tube	Dimensions ext.		t. mm	
		of tubes	Ø up to	Α	В	C	D
1000603	Stainless steel wire AISI 304	50 (10 x 5)	13 mm	162	84	35	35
1000605	Stainless steel wire AISI 304	50 (10 x 5)	13 mm	162	84	60	60
1000607	Stainless steel wire AISI 304	50 (10 x 5)	18 mm	210	110	40	40
1000609	Stainless steel wire AISI 304	50 (10 x 5)	18 mm	210	110	60	60
1000611	Stainless steel wire AISI 304	50 (10 x 5)	18 mm	210	110	85	85



Part No.	Made of	Number	For tube	Dimen	sions	s ex	t. mm
		of tubes	Ø up to	Α	В	C	D
1000604	Stainless steel wire AISI 304	100 (10x10)	13 mm	162 1	62	35	35
1000606	Stainless steel wire AISI 304	100 (10x10)	13 mm	162 1	62	60	60
1000608	Stainless steel wire AISI 304	100 (10x10)	18 mm	210 2	210	40	40
1000610	Stainless steel wire AISI 304	100 (10x10)	18 mm	210 2	210	60	60
1000612	Stainless steel wire AISI 304	100 (10x10)	18 mm	210 2	210	85	85



Suitable for Nessler tubes





Suitable for Falcon tubes with lids

Part No.	Made of	Number	For tube	Dime	nsior	ıs ex	t. mm
		of tubes	\emptyset up to	Α	В	C	D
1001715	Stainless steel wire AISI 304	15 (5 x 3)	16 mm	120	80	90	75
			Falcon 15 ml				
1001716	Stainless steel wire AISI 304	8 (4 x 2)	30 mm	145	80	90	75
			Falcon 50 mm				



With rack linking guides

Part No.	Made of	Number	For tube	Dimensions ext.		t. mm	
		of tubes	Ø up to	Α	В	C	D
1002532	Stainless steel AISI 304	24 (8 x 3)	12 mm	160	62	65	-
1002531	Stainless steel AISI 304	14 (7 x 2)	16 mm	160	62	65	-



Suitable for butyrometer

Part No.	Made of	Number of tubes	For tube	Dimensions ext. mm A B C D
	Stainless steel AISI 304 Stainless steel	12 (6 x 2)	Butyrometer	200 70 165 -
1023424	AISI 304	24 (8 x 3)	Butyrometer	285 95 165 -



With decanting aid

Part No.	Made of	Number	For tube	Dimensions ext. mm
		of tubes	Ø up to	A B C D
1000433	Anodised aluminium colour silver	120 (15 x 8)	11 mm	300 170 50 -
1000434	Anodised aluminium colour silver	120 (15 x 8)	14 mm	300 170 50 -



Stackable and numeric

Stackable racks

Made of anodised aluminium with alphanumeric from 1 to 10 and A to E.

Part No.	External dimensions mm.		Tube Capacity	For tubes of Ø up to	
	Height	Width	Depth		
1000830	60	165	95	50	11 mm
1000831	60	195	100	50	14 mm
1000832	60	220	110	50	16 mm




Transparent methacrylate racks

Part No.	Made of	Number	For tube	Dime	ension	s ext.	mm
		of tubes	Ø up to	Α	В	C	IJ
1052012	Transparent methacrylate	24 (8 x 3)	12 mm	170	90	50	50
1052016	Transparent methacrylate	24 (8 x 3)	16 mm	190	100	60	60
1052022	Transparent methacrylate	24 (8 x 3)	22 mm	255	150	80	80
1052026	Transparent methacrylate	24 (8 x 3)	26 mm	300	170	110	110



Part No.	Made of	Number	For tube	Dim	Dimensions ex		t. mm
		of tubes	Ø up to	Α	В	C	D
1001251	Methacrylate	25 (5 x 5)	11 mm	85	80	35	30
1001252	Methacrylate	25 (5 x 5)	12 mm	85	80	55	50



PROPISEL RACKS

Made of high density white polypropylene.

Easy to assemble. Stackable. Autoclavable. These racks do not float and are suitable for water baths for temperatures up to 140 °C.

Each rack has a cross reference numbering systems for easy tube identification. Easily cleaned with a disinfectant. Colour white.

Accessory rack identification cards are available in: yellow, blue and red.

Models					
Part No.	Number	For tubes	External dims. cr		
	of tubes	up toØ mm	Height	Width	Depth
1001115	90	13	7	10.5	25
1001116	60	16	7	10.5	25
1001117	40	20	7	10.5	25
1001118	24	25	7	10.5	25
1001119	21	30	7	10.5	25

ACCESSORY

Indentification cards, autoclavable.

Part No.

1001120 Yellow 1001121 Blue 1001122 Red

9999999999	
	PROPISEL-FIX RACKS
5555555555	Made of high density polypropylene. Comes with decanting aid that holds the tubes in place in the

rack while the complete rack is emptied in one ope-

ration. Features: Same feature as the "Propisel".

lmage	showing	the	decanting	aid	in
use an	d tubes si	till in	olace.		

Models					
Part No.	Number	For tubes	Exter	nal dim	s. cm
	tof tubes	up toØ mm	Height	Width	Depth
1001123	90	11 to 13	7	10.5	25
1001124	60	14 to 16	7	10.5	25
1001125	40	19 to 20	7	10.5	25

ACCESSORY

Indentification cards, autoclavable.

Part No.

1001120 Yellow 1001121 Blue 1001122 Red



RACK FOR 20 10 mm ✓ SQUARE CUVETTES

Made of polypropylene. Positioning cross-reference. Height 35 mm. Width 70 mm

Depth 210mm Colour: white

Part No. 1001528



CIRCULAR RACK

Made in polyethylene for 4 tubes of 25 mm Ø, 8 tubes of up to 19 mm Ø and 16 up side down peg tube supports. Dimensions: 95 high x 180 mm Ø.

Part No. 1001339



RACK FOR 1.5 AND 2 ml MICROTUBES

Made of Polycarbonate, stackable, numbered positions and autoclavable. Number of tubes: 24 microtubes. Diameter: up to 11 mm \emptyset .

Part No. 1001340



RACK FOR 1.5 AND 2 ml MICROTUBES

Made of solid polypropylene, ideal for refrigerator and freezer storage. Numbered positions and autoclavable. Capacity: up to 80 micro-tubes.

Part No. 1001341



RACK FOR 20 MICROTUBES OF 1.5 TO 2 ml

Made of polypropylene. Racks can be linked and fitted together. Numbered locations. Total length 21 cm. Part No. 1001523



THREE LEVEL MICROTUBE RACK FOR 0.2 - 0.5 AND 1.5 ml

Made from polypropylene and can be linked at either end to form longer racks. Numbered positions. Lower level 24 micro tubes of 0.2 ml. Middle level 14 micro tubes of 0.5 ml. Upper level 12 micro tubes of 1.5 ml.

Part No. 1001542

Standard colour yellow.



MULTI SUPPORT TUBE RACK TUBE SIZE Ø 10/13 AND 14/17 mm

Suitable for drying and storing of tubes. Autoclavable. Suitable for use in a water bath.

Standard colour blue.

Part No. 1001543 Group of 80 spikes for tubes

of 10-13 mm Ø.

Part No. 1001544 Group of 50 spikes for tubes

of 14-17 mm Ø.



REVERSIBLE TUBE RACK FOR MICROTUBES 0.2/0.5 AND 1.5/2 ml

Made from polypropylene with removable lid that fits on any side.

Alphabetic identification positions.

Face A for 168 micro tubes of 0.2 ml and 12 x 1.5 ml. Face B for 110 micro tubes of 0.5 ml and 12 x 1.5 ml. Standard colour yellow.

Part No. 1001545




REVERSIBLE TUBE RACKS FOR 0.5 - 1.5 AND 2 ml tubes

Face A for 96 micro tubes of 0.5 ml.
Face B for 96 micro tubes of 1.5 to 2 ml.
Made from polypropylene with removable lid.
Alphabetic tube position identification.

Height 26 mm, with lid 50 mm. Width 115 mm. Length 250 mm.

Standard colour yellow.

Part No. 1001524



BOX TO STORE 100 CRYO VIALS AND 1.5 TO 2 ml MICROTUBES

Made from polyethylene, with hinged, lid and pressure latch close.

Alphabetical grid tube location within the lid. Height with lid 55 mm. Width 140 mm. Length 140 mm. Standard neutral colour.

Part No. 1001525



BOX FOR 50 CRYO VIALS AND 1.5 TO 2 ml Microtubes

Made from polyethylene, with hinged, lid and pressure latch close.

Alphabetical grid tube location within the lid. Height with lid 55 mm. Width 80 mm. Length 140 mm. Standard neutral colour.

Part No. 1001526



STORAGE BOX FOR 96 TUBES OF 0.2 ml

Made from polypropylene, with hinged lid and pressure latch close.

Alphanumeric identification positions.

Height with lid 30 mm. Width 90 mm. Length 130 mm. Standard colour yellow.

Part No. 1001527



BOX FOR STORING 12 X 10 mm CUVETTES

Made from Polypropylene, with hinged lid and pressure latch close.

Removable base.

Cuvette heights up to 60 mm high.

Height with lid 60 mm. Width 115 mm. Length 115 mm. Standard colour neutral.

Part No. 1001529



FREEZER BOXES FOR VIAL STORAGE

Made from polycarbonate.

Suitable for freezer or liquid nitrogen storage for temperatures between -196 °C to +121 °C.

Alphanumeric identification, with ventilation holes and drain. Pressure lid.

Part No. 1001713 Box of 25 vials.

Height with lid 55mm. Width 75 mm. Length 75 mm. Part No. 1001714 Box of 100 vials, Lid with hinge. Height with lid 55mm, Width 130mm. Length 130 mm.

COMECTA

Chiller and cooler storage boxes




COMMON FEATURES

Specially designed for protecting and maintaining PCR* reagents, enzymes, bacteria, viruses and other biological samples cool or frozen during use. Also provides protection from temperature fluctuation, defrost cycles and freezer power failure.

Made from polycarbonate with a pre-filled chamber of non-toxic gel, complete with alphabetical identification grid in the lid.

Freezer temperature resistant down to -135 °C.

Tubes are isolated and have no contact with the cooling solution.





"QUICK FREEZE" CONTAINER

Quick freeze container with the capacity for 12 tubes 0.5 ml or 2 ml tubes.

Unique reproducible freezer system for biological samples using ethanol samples of DNA and RNA.

It will maintain the interior temperature for up to 45 minutes after removal from a freezer.

Freezing speed:

2 ml solution of alcohol at -20 °C, 3 minutes.

2 ml solution of water 0 °C, 4 minutes.

Size (with lid): 95 high x 100 depth x 140mm width.

Part No. 5601050



MINI FREEZER CONTAINERS "COOLER 1051" AND "COOLER 1052"

Cooler 1051

Freezer container at 0 °C with the capacity for 12 tubes of 0.5 to 2 ml.

Alphanumeric matrix on the lid.

It will maintain temperature for up to 2 to 5 hours

Size (with lid): 95 high x 100 depth x 140mm width.

Part No. 5601051

Cooler 1052

Freezer container at 0 °C to -20 °C with a capacity for 32 tubes of 0.5 to 2 ml.

Alphanumeric matrix on the lid.

Double insulation body and lid that can maintain the internal temperature for a period of up to 2 hours.

Size (with lid): 105 high x 130 depth x 270 mm width.

Part No. **5601052**



MINI FREEZER CONTAINERS "COOLER 1053"

Freezer container at 0 $^{\circ}$ C with the capacity for 12 tubes of 16/17 mm Ø and 15 ml centrifuge tubes. Alphanumeric matrix on the lid.

It will maintain temperature for up to 2 to 5 hours.

Size (with lid): 150 high x 140 depth x 190 mm width.

Part No. **5601053**



CONTAINER "CRYO 1 °C COOLER"

Freezer container with constant reproducible progressive cooling rate of 1 °C/min with a capacity for 18 cryo tubes of 1.2 to 2 ml

Operation: fill the container with alcohol isopropanol (C₃H₈O) and place tube rack inside replace the lid and put into a freezer.

The container can be placed in a water bath for thawing.

Capacity of C₃H₈O = 250ml

Size (with lid): 86 high x 117 mm Ø.

Part No. 5601055



"ICE BUCKET"

Container with lid, made of double-bodied polystyrene of 2 mm thickness with thermal insulation that is chemically resistant, made from expanded polystyrene.

Ideal for maintaining samples submerged in ice for short time periods. Capacity 3-5 litres.

External dimensions: 19 height x 27 width x 27 cm depth.

Internal dimensions: 17 height x 25 width x 25 cm depth.

Part No. 5601054



Flake ice maker "IMS-40", "IMS-85" and "IMS-150B" "IMS" 40 MODEL: PRODUCTION 40 kg-24 HOURS.



"IMS" 85 MODEL: PRODUCTION 85 kg-24 HOURS. "IMS" 150 MODEL: PRODUCTION 150 kg-24 HOURS.



For preservation and storage of analytic samples, biological tissues, cellular cultures and other applications that require maintenance of the temperature below ambient.

FEATURES

External case made of AISI 304 stainless steel.

Thermal isolated inside.

High capacity for flake ice production from running or purified water.

Water pressure between 1,5 and 3 bars.

Automatic switch off when there's water shortage and when the ice tank capacity is full.

CONTROL PANEL

Main power switch.

Luminous indicator for: Flake ice production.

Water shortage. Ice tank full.



■ © 013
100 miles 100 mi

MODELS	Part no.	24 h. (Kg) Production	Tank Capacity(Kg)	Height / Width / Depth cm	Power W	Power supply	Weight Kg
IMS-40	5020100	40	8	72 38 54	280	220V/50Hz	40
IMS-85	5020101	85	13	88 55 61	420	220V/50Hz	57
IMS-150B	5020103	150	35	94 63 61	685	220V/50Hz	76

COMPLEMENT "ICE BUCKET"

See page 190



Manual hidraulic press "PM-15"

MAXIMUM PRESSURE LOAD: 15 TONS

FEATURES

Compact and easy to use, designed for a maximum pressure load of 15 tons, especially for pressing in chemical research laboratories and materials for preparing ceramic composite. This machine has an integrated auto return hydraulic pump.

Tablets from 6.3mm to 19mm can be compacted with the MT matrix accessory, for KBR compression technique and fluorescence analysis works. Cylinder elevation of 0-25mm.

MODEL

Part No.	Max. Pressure tons	Height A	/ Widtl mm	n / Depth	Weight Kg	Max. Working area mm Height / Width / Depth
5810200	15	55	13	21	30	150 x 160 x 160

ACCESORIES

Part No. 5810201 MT Matrix 12.5mm Part No. 5810202 MT Matrix 6.3mm Part No. 5810203 MT Matrix 19mm Part No. 5810204 MT Matrix 3.2mm









FREEZING CONTAINER

For sample freezing with liquefy gas (liquid nitrogen), liquids and solutions.

Spherical inside. Exterior made of AISI 304 stainless steel. Double body with thermal isolated PRTFE joint. Hit-resistant.

External measures: 14 high x 24cm wide. Capacity: 1.5 L. Weight Kg: 2,9.

Part No. 1001640



FLOATING TUBE RACK 1.5 AND 2 ml MICRO-TUBES

Made of polypropylene, with a support to hold microtubes. Ideal for constant temperature application.

Capacity: 16 microtubes.

Dimensions: 4 high x 10 wide x 10 cm deep.


Part No. 1001342



THERMOMETER AND PIPETTE RACK

Made of polypropylene. No. of positions: 12 positions.

Part No. 1001314



PIPETTE RACK "CIRPIP"

Portable upright rack for pipettes and thermometers, made of PVC. No. of positions: 93 of various sizes. Stable heavy base with drain


Weight 0.4 Kg. Part No. 1001254



PIPETTE RACK "STABIL"

Pipette rack made from polyethylene for 28 pipettes or thermometers up to 14 mm \emptyset . Stable heavy base with drain and two adjustable height supports. Weight 1.5 Kg.

Part No. 1001466



SAMPLE CONTAINERS "CONTESEL"

Applications: designed for dispensing, storage or preservation of organic and inorganic samples.

Polypropylene container and quick fit screw cap lid.
Lid system with flap AQUISEL (international patent.)
That allows the manipulation of the sample without having to remove the lid and maintains a seal.
Dosing by means of an automatic or standard pipette.
Vessel graduated to: 20, 40, 60, 80 and 100 ml, frosted screen for writing.

Includes gullet for pouring.

Comes complete with violet coloured lid valve.

Part No.	colour	capacity	total	quantity
		ml	height cm	
1001699	violet	100	7	300

CIRCULAR RACKS

Made of AISI 304 stainless steel.

138 mm \emptyset x 120 mm high. Suitable for the "Baher" and "Univeba" baths (see BATHS for more details).

Part No. 7000541 for 24 tubes up to 18 mm in \emptyset . Part No. 7004002 for 24 tubes up to 13mm in \emptyset .



SUPPORT AND DRAIN RACK FOR PIPETES "SOP-PIP"

Made AISI 304 stainless steel. Holds up to 112 pipettes of different sizes. Part No. 1001223



PIPETTE RACK

Made of AISI 304 stainless steel. No. of places: 6 dosing pipettes of various sizes up to 28 cm in length. (See pages 215 - 216). Part No. 1001215




TABLE TOP OR WALL MOUNTABLE DRAIN RACK

Made of PCV coated wire. 640 high x 500 mm wide. With 79 positions, complete with drip tray. Part No. 1000481



TILTED SUPPORT RACK "SOP-INC"

Suitable for pipettes, hydrometers and thermometers rack, made of methacrylate. 6 positions. Part No. 1000015



SUPPORT RACK FOR PIPETTE AID "COM-PIP"

Made of AISI 304 stainless steel. No. of places: 6 dosing pipettes of various sizes. (See page 215-216). Part No. 1001702



KJELDAHL FLASK SUPPORT

Made of AISI 304 stainless steel with electrobright finish. Suitable for flasks up to 500 ml. Part No. 1000002



ROUND BOTTOM FLASK SUPPORT

Made of moulded rubber (E.P.D.M.). Temperature resistance up to 100 °C. Specially designed with stepped surface to support different flask diameters. From 100 ml up to 10 litre capacity.

Dimensions: 50 mm high x 160 mm Ø.

Part No. 1000863



DRAIN RACK

Made of PCV coated wire. 32 positions. 300 high x 400 mm wide. Part No. 1024132







88 positions. 620 high x 500 mm wide. Part No. 1024188



WIRE BASKETS

Made of AISI 304 stainless steel sheet with electro-bright finish. Designed to fit the Selecta autoclaves. (See page 55 to 70).

Part No.

1004771 220 Ø x 140 mm high. 1000596 220 Ø x 240 mm high.


1000495 265 Ø x 185 mm high.

Part No.

1000496 365 Ø x 265 mm high.

1000780 435 Ø x 205 mm high.

With cross dividers.



ACCESSORY

Basket dividers, stainless steel AISI304, electrobright finish.

Part No. For baskets, Part No.

1000495 1001217 4 compartment divider 1001218 6 compartment divider 1000496 1001222 7 compartment divider 1000780



WIRE BASKETS

Electropolished stainless steel AISI 304.

Part No. **Dimensions** 210 Ø x 180 mm high 1025115 150 Ø x 120 mm high 1025121 1025124

Part No. 1125115 1125121 240 Ø x 180 mm high 1125124



STERILIZATION DRUMS

Made from AISI 304 stainless steel. Designed to fit the Selecta autoclaves. (See page 55 to 70.)

Part No. Ø cm Height cm 1002418 28 18 1002419 35 27

1002420 28 18 Doesn't match Selecta autoclaves.



STERLIZER CAN FOR **PIPETTES**

Made from Stainless steel AISI 304.

Part No. height Ø mm mm 1000829 400 72 1001126 400 130

STERILIZATION CONTAINER **FOR PETRI DISHES**

Made of AISI 304 stainless steel. Complete with Petri dish rack for up to 120 mm Ø.

	~.	
Models	height / Ø	No. of
Part No.	mm	plates
1001127	260 / 130	10
1001128	400 / 130	18



PETRI DISH RACK

Made of AISI 304 stainless steel, with electrobright finish. For 10 dishes of 80 to

100 mm Ø. Part No. 1000259



PIPETTE BASKET

Made of AISI 304 stainless steel with electrobright finish. 90 mm Ø x 250 mm high.

Part No. 1000251



STERILIZATION CONTAINER **FOR PIPETTES**

Made from Stainless steel AISI 304. Lid interior and base lined with silicon to protect pipette end and tip.

Part No.	height mm	⊿ mm
1001597	290	65
1001598	340	65



Burners, alcohol, butane, propane and natural gas



ALCOHOL BURNER Screwcap jar with wick. Capacity 120 ml. Part No. 1001256




ALCOHOL BURNER Chrome plated top with wick height control. With lid and chain. Capacity 85 ml. Part No. 1001434



BUNSEN BURNER

With both gas and air regulator. Tube Ø exterior: 11 mm. Total height: 140 mm. Part No.

7000134 For butane/propane gas. For natural gas.



SIMPLE BUNSEN BURNER Air regulator only, nickel plated. External tube \emptyset : 11 mm. Total height: 130 mm. Part No.

7000044 For butane/propane gas. 7000937 For natural gas.



MEKER BURNER

With permanent flame. Gas and air control. External tube Ø at its widest point: 22 mm. Total height: 175 mm.

Part No.

7001106 For butane/propane gas. 7001107 For natural gas.



MEKER BURNERS

With gas and air regulator.

Model Small. External tube \emptyset at its maximum point: 22 mm. Total height: 175 mm.

Part No.

7001108 For butane/propane gas.

7001109 For natural gas.

Model Medium. External tube \emptyset at its maximum point: 27 mm. Total height: 185 mm.

Part No.

7001110 For butane/propane gas.

7001111 For natural gas.

Model Large. External tube \emptyset at its widest part: 32 mm. Total height: 190 mm.

Part No.

7001112 For butane/propane gas.

7001113 For natural gas.



TECLU BURNER

Part No.

With air and gas regulator. External tube Ø: 13 mm. Total height: 180 mm.

7012019 For butane/propane gas.

7094519 For natural gas.



NATURAL GAS BURNERS

With air and gas regulator.

With approved gas pressure connection 1/2" nut and olive, supplied as an accessory.

External tube Ø: 11 mm. Total height: 165 mm.

Part No. For gas Working Weight Normal pressure power Kg **7001551** natural 18 mbar 0.69 kw 1.6

ACCESSORY

Gas tube with connections (see page 196).



BUNSEN BURNER WITH SAFETY GAS CUT OUT

With gas and air regulator

Designed for maximum laboratory safety.

Equipped with a flame detector that will shut off the gas supply within 20 seconds if the flame is extinguished. With tube connection.

Tube Ø ext.: 11 mm. Total height: 140 mm.

Tube b ext	ii iiiiii. iotai	noigni. 140 mm	
Part No.	Gas type	Pressure	Nominal
			power
7001546	butane	28/37 mbar	0.69 kw
70015/17	natural	18 mhar	N 69 kw



BUNSEN BURNER WITH SAFETY GAS CUT OUT, NATURAL GAS

With gas and air regulator.

Designed for maximum laboratory safety.

Equipped with a flame detector that will shut off the gas supply within 20 seconds if the flame is extinguished. Tube connection of 1/2 inch. Supplied as an accessory. Tube Ø ext.: 11 mm. Total height: 165 mm.

Part No.	Gas type	Pressure	Nominal	Weigh
			power	Kg
7001539	natural	18 mbar	0.69 kw	1.6



BUNSEN BURNER

7000661

With air and gas regulators.

Ø ext. tube: 11 mm. Total height: 140 mm.						
Part No.	Gas type	Pressure	Nominal			
			power			
7000657	butane	28/37 mbar	0.69 kw			
7000658	natural	18 mbar	0.69 kw			
Tube Ø ext.	Tube Ø ext.: 13 mm. Total Height: 170 mm.					
Part No.	Gas type	Pressure	Nominal			
			power			
7000660	butane	28/37 mbar	1.14 kw			

18 mbar

1.14 kw

natural



BURNER FOR BUTANE GAS

Suitable for use with bottled gas "International Camping Gas". Approved following ACTA of the GENERAL LABORATORY ANALYSIS AND SERVICES FOR INDUSTRY. Manufactured in accordance with the "regulation for apparatus used with combustible gasses" decree 1651/1974. External tube Ø: 14mm. Direct regulated pressure. Nominal power 0.72 Kw at 1 bar.

Part No. **7000666**

ACCESSORY



Gas cable with connections, flexible coated stainless steel 100 cm long with 1/2 inch connection with quick fit and safety valve. Conforms to the UNE 60715 & BS669. Safety regulations.

Part No. 1001541



BUNSEN WITH PERMANENT PILOT FLAME

Permanent flame. Total height: 150 mm. Part No.

1000875 For butane/propane gas (45/57 mbar). **1000876** For natura gas (18/25 mbar).



BUNSEN BURNER WITH CONICAL TUBE TECLU

With air and gas regulators.

Total height: 150 mm.

Part No.	Gas type	Pressure	Nominal
7000663 7000664	butane natural	28/37 mbar 18 mbar	power 0.69 kw 0.69 kw

ACCESSORIES



Wide flame burner head made of tin. Suitable for the "Camping Gas" burner.

Part No. 7003612



Tripod support with flat ring, fitted with a burner, suitable for "Camping Gas International" bottle.

Exterior Ring \emptyset : 140 mm. Interior Ring \emptyset : 110 mm. Part No. 7000489

diameter 21 mm. Total height 170 mm. Part No. 7000225 Suitable for butane/propane. Part No. 7009223 Suitable for Natural gas.

With air and gas regulators. Tube \emptyset ext. at its maximum

HIGH POWER BURNER



BUNSEN J2000 BLUE UNIVERSAL

For Butane gas. Instant safety flame push button ignition incorporating a LOCK/ON mechanism that prevents accidental fires, Unit switches off when the users hand is removed. Gas fill by a simple lighter gas fuel as used common gas lighters. Long life flint. Power by a 1.5 V battery.



APPLICATIONS

Microbiology, bacteriology, dentistry and for general use where a fine hot flame is required for fast heat and fusion.

Part No. 1001114



BURNER "MICROTORCH"

APPLICATIONS:

Used for soldering, dentistry, general laboratory use. Suitable for plastic moulds, glass and metal tubes. etc. Butane gas.

Portable with a removable bench top base. Microflame instant piezo electric start. Flame temperature controllable to 1300 °C. Refill using common gas lighter cartridge. Gas capacity: 26 g. Total time before refill: 2 hours. Part No. 5000011



BLOW PIPES

Part No.

1000877 Butane gas. **1000878** Natural Gas.

ACCESSORIES

FLAME SPREADER BURNER HEAD



Straight form. Stainless steel AISI 304.

Part No. **7003611** For tube of Ø 11 mm. Part No. **7003613** For tube of Ø 13 mm.



GAUZE IRON WIRE WITH CERAMIC FIBRE

Part No.

 1087102
 100 mm side.

 1087122
 120 mm side.

 1087142
 140 mm side.

 1087162
 160 mm side.

 1087182
 180 mm side.

 1087202
 200 mm side.

FLAME SPREADER BURNER HEAD



Straight form. Made of tin.

Part No. **7003612**Part No. **7003614**

For tube of Ø 11 mm. For tube of Ø 13 mm.



TRIPOD STAND WITH RING TOP

Galvanised steel.

 Part No.
 Ø exterior ring
 Total height

 7000808
 80 mm.
 210 mm.

 7008010
 100 mm.
 210 mm.

 7008012
 120 mm.
 210 mm.


 (For variable height tripods see page 175).



CERAMIC PLATES "NEOCERAM"

Resistant to temperatures from -200 °C to 700 °C. Part No.

1000684 120 x 120 x 3 mm. 1000685 140 x 140 x 3 mm. 1000686 160 x 160 x 3 mm.



Excellent heat conductor.

20% better heat transfer compared to a conventional gauze ceramic.

High chemical resistance to aggressive compounds. Easy to clean and long lasting.



Infrared tip sterilizer "Sterilbio"

APPLICATIONS

Sterilization of culture tips, needles and small instruments. Incinerates organic material at 900 °C for 5/8 seconds.

FEATURES

Made from ceramic elements with a quality quartz tube interior.

Does not need purged oxygen and can be used in anaerobic environments.

Eliminates problems of aerosols and bacterial contamination.

Support at the side for 4 Kolle handles, needles, lancet or dissecting forceps etc.

MODEL

Part No.			/Depth	Ø/Le	J .	Power	Weight
	(e:	xterior)	cm	(tube inte	erior) cm	W	Kg
3000788	20	13	18	1.4	15	110	1.5





COMECTA Hot Air gun "Salki 1500"

WORKING TEMPERATURE 400 °C AND 590 °C.

APPLICATIONS

Suitable for heating liquids, solids and gasses:

- Dry
- Soft solder
- Fusions
- Evaporation
- Curing adhesives
- Plastic moulds
- Softening surfaces
- Fast defrost
- Film drying R.X.
- Chromatography applications drying paper and electrophoresis.
- Heating moulds and drying dental moulds etc.

FEATURES

Industrial model with vertical angle support.

Brushless motor.

Three positions:

- Position 0 stop.
- Position I temperature 400 °C.
- Position II temperature 590 °C.

MODEL

Part No.	Temperature °C	Power W	Weight Kg
5000750	until 590	1500	0.8



ACCESSORIES



Burner head 75 mm. Part No. **5011148**





Tube deflector. Surface 75 mm. Part No. **5014449**

Diverse laboratory instruments

COMECTA Pocket Balance "Pocket Scale"

FEATURES

Capacity up to 200 g, measuring in g, oz and qty. Digital LCD display.

Stabilization time ±3 seconds.

Tare by subtraction to maximum weight.

Piece counting function and automatic reset to zero.

Automatic calibration, external,

Batteries 3 x 1.5 V.

Lid to cover the weighing pan and operation menu.

MODEL

Part No.	Max. weight.	Precision	Height / Width / Depth	Width / Depth	Weight
	capacity g	±g	(exterior) mm	(plate) mm	Kg
5830024	200	0.1	24 74 120	74 71	0.13

Comes with protective cover.



Several Instruments





DIAMOND POINT PENCIL For glass engraving. Part No. 1001130



TUNGSTEN SCRIBE

For glass engraving, slides, ceramics, plastic etc. Retractable tip and pocket clip.

Part No. 1001674



GLASS TUBE CUTTER

For tubes with diameters up to 30 mm.

Part No. 1000016



CORK PRESS Model Crocodile shape.

Part No. 6000070





LABORATORY VIBRO-ENGRAVER SIGNOGRAPH-COM 25 DIAMOND

For engraving all types of glass, metal, plastic material, ceramic etc. Mains voltage 220-230 V. 50Hz. Power consumption: 25 W. Comes complete with a set of 5 engraving tips: Diamond, carbide, tempered steel and two for decorative work.

Part No. 1001129



WATER JET VACUUM PUMPS

Part No. Material Pressure Pump rate Aspiration 8 I/min 2.6 I/min 4000633 PVC 2 Kg 2 Kg Metal 5 I/min 1 I/min 7000293 7000292 0.5 I/min Metal 2 Kg 5 I/min



BLOW PIPE

Student model. Total length: 250 mm.

Part No. 1000595



Made of grey steel. Only suitable for dry products.

Ø Usable Height Part No. 130 mm 65 mm 1004713 1004715 150 mm 75 mm



CORK BORER SHARPENER

With a wooden handle, blade of tempered steel. Part No. 7000004

CORK BORER

Electrobrite finish. Electrobrite finish. Universal handle. Individual handle. Part No. Part No. 7000846 6 bores. 7000856

6 bores. 7000848 8 bores. 7000858 8 bores. 7008412 12 bores. 7008512 12 bores.



MINI CORK BORER

1001332 6 bores: 4, 4.9, 5.5, 6.4, 8 and 9.5 mm Ø.



PIPETTE FILLER, RUBBER BULB.

With three glass ball valves for filling, fine adjustment and empty. Accepts pipettes from 6 to 10 mm Ø.

Part No. 1001733



RICHARDSON RUBBER BULB

With net covered decompression bulb. For blowing out bench top automatic burettes.

Part No. 1001168





SPECTACLES "ULTRASPEC 2000"

Ideal for laboratory use. UV resistant.

Conforms to: DIN 4646-47 and DIN 58211; Spain MT 16/MT 17 regulations. Safety glasses made from one single lens, extremely light weight, anti-scratch coated. Panoramic vision due to a 180° single lens. One size fits all. Ventilation anti mist vents at the side of the legs.

Part No. 5028192



SPECTACLES "ULTRASPEC 1000"

Ideal for working in the laboratory, UV resistant.

Conforms to: DIN 4646-47 and DIN 58211; Spain MT 16/MT 17. Safety glasses made from one single lens, extremely lightweight, comfortable with a 180° panorama. Adjustable legs: length and angle. Detachable frame for easy cleaning.

Style: Sporty.
Part No. 5028194



2. GOGGLES FLEXIBLE "IVYMEN 3000"

Ideal for general laboratory use.

Conforms to: EN-166.

Panoramic view. Prescription spectacles can be used underneath the

Made from transparent flexible plastic that moulds to the face. Held on by an stretch band. Indirect vents at the side prevent misting.

Part No. 5028196

1. SPECTACLES "IVYMEN 1000"

Ideal for general laboratory use.

Conforms to: DIN 4646-47 and DIN 58211.

Safety Spectacles of one single lightweight lens, comfortable protection. With a wide panorama of 180°. Adjustable legs: length and angle. Removable frame for easy cleaning. Style: Sporty.

Part No. 5028195

3. SPECTACLES "IVYMEN 2000"

Ideal for general laboratory use.

Conforms to: DIN 4646-47 and DIN 58211

Safety spectacles of one single lightweight anti scratch lens. Perceptual panoramic vision of 180°. Can be use over normal prescription spectacles. Indirect vents at the side of the legs prevent misting.

Part No. 5028193

GLOVES ANTI-ACID "SOLVEX-PLUS"

Conforms to CE standards.

Effective protection against chemical products and solvents, manufactured with a thick coating of strong durable Nitrile. Flexible with anti bacterial and fungal protection for hygiene and against skin irritations.

Length 32 cm. Universal fit.

Part No. **5000044**

GLOVE HAND PROTECTOR AGAINST HEAT AND SPILLS

Adaptable silicon glove. Part No. 1001255



GLOVES THERMAL "KEVLAR 400"

Conforms to: EN-388 and EN407

These gloves can withstand temperatures up to 400 °C. Manufactured from Kevlar terry knit, seamless, with double face fibres, ambidextrous fit, high level of protection against heat and flame. Length 25 cm, universal fit.

Part No. 5000041



APRON ANTI-ACID "DELNEO 125"

Conforms to: EN 342 and EN 511 standards. For effective protection against chemicals and solvents, made from high resistant Neoprene.

Part No. **5000045**



GLOVES THERMAL "KEVLAR 800"

Conforms to EN 388 and EN407 standards. For use with temperatures up to 800 °C, Made from terry knit, seamless, with double face fibres, high level of protection against heat and flame. Length 36 cm, universal fit.

Part No. 5000042



PROTECTIVE MASK "P3-SL"

Conforms to DIN-EN 149:2001 standards. High range auto-filtering mask, for particles, dust, mists and smoke use. Made of cellulose and with a filtration valve.

Part. No 5000046





TABLE TOP MAGNIFIERS WITH ILLUMINATION

Magnifiers with high quality lenses and double lateral illumination, plastic body.



made from epoxy coated metal.

Lens size cm Magnification Power Weight Kg 15 height x 19 width 3 x 2 x 9 W 3.5 Part No. 5310500

Extendable and height adjustable angle poise, 89 cm With support base, large base with adjustable height and orientation.

> Lens size cm Magnification Power Weight Kg 16 height x 19 width 2 x 7 W 3 x Part No. 5310506



VERNIER CALLIPERS WITH DIGITAL DISPLAY DIN 862

5 digits display.

AISI 304 Stainless steel with double measuring scale. Reproducibility: 0.01 mm or 0.0005 inches

1.5 V battery power.

Battery life: 1 year of continuous use.

Comes complete with wooden case.

Part No.

5900601 Measuring range: 0 - 150 mm or 0 - 6". 5900603 Measuring range: 0 - 300 mm or 0 - 12".



BIFOCAL AND MONOFOCAL "MAGNIFIER"

Magnifying glass bifocal "Magnifier"

Completely made of high impact plastic.

Part No. Ø Lens Magnification **Bifocal** 5310501 65 mm 3 x 5 x 5310502 90 mm 2 x 5 x

Magnifying glass monofocal "Magnifier"

Lens made of high quality ground crystal.

Rim made from chromed steel with plastic handle. Part No. Ø Lens Magnification

5310503 50 mm 2 x 5310504 75 mm 2 x 2 x 5310505 100 mm



PALMER MICROMETER WITH DIGITAL DISPLAY DIN 862

5 diaits.

Measuring range: 0-25 mm or 0-1 inch. Resolution: 0.001 mm or 0.0001 inches. Precision: 0.002 mm. 1.5 V battery power.

Battery life: 1 year of continuous use. Comes complete with wooden case.

Part No. 5900602



AUXILIARY TABLES HEIGHT ADJUSTABLE FOR LABORATORY

Cabinet in a corrosive resistant epoxy coated.

Supplied with 4 lockable wheels for an easy transportation and 4 suction cups for a static clamping.

Specially designed to support different laboratory instruments (Ovens, furnaces, centrifuges, etc...).

5 adjustable heights: 72, 75, 78, 81 or 84 cm.

Measurements Part No. Load Weight Capacity Kg cm Kg 75 width x 60 depth 25 1001781 120 100 width x 70 depth 29 1001780 150

NEW



MEASURING RANGE


Made of polypropylene with silicon joints, for volumes up to 25 ml. Precise and safe, ergonomic design that allows filling and emptying without having to change the setting. Flexible thermoplastic seal against any type of burette. A button on the side empties the burette automatically. Easily dismantled and cleaned.

Part No.	Capacity	Identification colou
1001131	0 to 2 ml	Blue
1001132	0 to 10 ml	Green
1001133	0 to 25 ml	Red



PIPETTE RACK "COM-PIP"

Made from AISI 304 stainless steel.
Capacity 6 pipettes of differing sizes.
Part No. 1001702



REVERSIBLE POLYPROPYLENE UNIVERSAL TRAY

Suitable for 8 and 12 channel micro-pipettes, 5 ml capacity per channel or 1 single capacity of 50 ml, graduated at 10-25-40 and 50 ml.





BRUSHES "COMCEP"

Made with natural bristles. Applications:

- 1 Pipettes and capillaries of 3 to 5 mm Ø. (Brush)
- 2 Graduated pipettes and tubes from 5 to 10 mm Ø. (Brush)
- 3 Haemolysis tubes 8 to 10 mm Ø. (Brush)
- 4 Test tubes and graduated flasks of 10 to 16 mm Ø. (Brush)
- 5 Burettes of 25 to 50 ml and tubes of 15 to 20 mm Ø. (Brush).
- 6 Cylinders, flasks and test tubes etc. from 20 to 32 mm Ø. (Fan Brush).
- 7 Erlenmeyer flasks, test tubes, flasks etc. from 25 to 45 mm Ø. (Fan Brush).
 8 Flasks, Erlenmeyer flasks and test tubes etc. 35 to 80 mm Ø. (Fan Brush).
- 9 Beakers cylinders and Erlenmeyers etc. from 50 to 100 mm Ø. (Brush).
- 10 Test tubes of 16 to 25 mm Ø. (Brush).
- 11 Test tubes of 16 to 25 mm Ø. (Fan Brush).
- 12 Flasks, round bottom flasks of 0 to 100 mm Ø.
- 13 Beakers of 50 to 250 mm Ø. (brush with handle).

Part No.	Style	Brist	le mm	Length total	Bristle
		Ø	lenght	mm	type
1001143	1	4	90	400	Goat hair
1001144	2	10	100	700	Goat hair
1001145	3	10	90	240	Pig bristle
1001146	4	15	100	320	Pig bristle
1001147	5	25	100	700	Pig bristle
1001148	6	30	100	300	Pig bristle
1001149	7	35	100	440	Pig bristle
1001150	8	55	100	440	Horse hair
1001151	9	80	130	700	Horse hair
1001333	10	25	130	350	Pig bristle
1001334	11	25	100	320	Pig bristle
1001335	12	110	450	450	Nylon
1001336	13	40	250	400	Nylon



POLYCARBONATE CONICAL FLASKS

DURABLE, CAPACITY 250 ml

Rigid, transparent, autoclavable, unbreakable, reusable.

Density: 1.20 g/ml.

Temperature resistance from -135 $^{\circ}$ C to +135 $^{\circ}$ C.

Excellent chemical resistance including acids. Exceptions chemical that attack

polycarbonate. Two models with or without screw-top, PP caps included.


Part No.	Capacity	PP cap	Specific
	ml		weight
5310100	250	Not	1.2
5310101	250	Yes	1.2


Laboratory sponges and fiber scrubbers

Features

Designed specially for cleaning plastic and glass items.

Made from polyethylene foam that is durable and long lasting (20 times more than a sponge or brush.) Do not absorb water, easily cleaned with hot or cold water.







Models

 Used for all types of tubes and cylinders such as: test tubes distillation tubes, thermal flasks, cryogenic vessels, measuring cylinders etc.

N^o	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
1	1001435	70	100	480

2. Used for cleaning short vessels such as: jars, beakers, etc.

N^o	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
2	1001436	70	100	250

3. Scrubber suitable for all types of assay tube up to 25 mm Ø. Adjustable width.

Nº	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
3	1001437	22	170	220

4. Suitable for all types of tubes and vessels including Conical flasks, etc. Very flexible, can be used in vessels with small necks, such as separating funnels, by simply pressing down the scrubber that is sufficiently flexible to take the internal shape and also be easily extracted.

Nº	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
4	1001438	20	380	440

5. Fiber scrubber with brush polyethylene at the end, suitable for cleaning flasks, culture jars and vessels in general.

Nº	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
5	1001439	50	140	300


Fiber scrubber with brush polyethylene at the end, suitable for any kind of vessels cleaning, test-tubes, culture jars and flasks in general.

Nº	Part No.	Sponge width	Sponge length	Total length
		mm	mm	mm
6	1001440	50	180	500



Stirrer bars and propellers

TO CHOOSE THE CORRECT STIRRER BAR: THE LENGTH SHOULD BE APPROX. 2/3 OF THE DIAMETER OF THE VESSEL BEING USED.



STIR-BARS FOLLOWERS FOR MAGNETIC STIRRERS.

Encased in PTFE with pivot ring, temperature resistant to +275 °C.

Part No.		
1000009	6 Ø x 10 mm long	without pivot ring.
1000012	4 Ø x 12 mm long	without pivot ring.
1000019	7 Ø x 20 mm long	with pivot ring.
1000032	8 Ø x 32 mm long	with pivot ring.
1000039	8 Ø x 40 mm long	with pivot ring.
1000045	10 Ø x 45 mm long	with pivot ring.
1001061*	12 Ø x 75 mm long	with pivot ring.
1001062*	16 Ø x102 mm long	with pivot ring.
1001063*	16 Ø x127 mm long	with pivot ring.
1001064*	19 Ø x150 mm long	with pivot ring.

^{*} With pivot ring for large volumes.



CONCENTRIC STIR-BAR WITH BEAKER COVER "CENTRIMAN"

Encased in PTFE with an adjustable axis bar made of AISI 304 stainless steel. Suitable for applications with volumes of 100 to 2000 ml.

Soft start stirring and adjustable height prevents the stirbar jumping inside the vessel.

Conical fitting lid of polypropylene with excellent chemical and thermal properties that resists temperatures up to 160 °C. The lid has an auto-lubricating central axis hole with two lateral air displacement holes.

Part No.

1000900 Concentric stir-bar with conical fitting lid for beakers volumes of 100 ml to 500 ml. Length


32 mm.

1000901 Beaker capacities of 500 ml up to 1000 ml.

Stir-bar length 50 mm.

Beaker volumes of 1000 ml up to 2000 ml. 1000902

Stir-bar length 55 mm.



TRIANGULAR STIR-BAR "TRIMAG"

Encased in heat resistant PTFE of temperatures up to 275 °C. useful for "dissolving solids" and stirring silts. Liquids flow upwards as well as sideways. Quiet with excellent stability. Suitable for flat bottomed vessels.

Part No.

1000536 9 side x 35 mm long. 1000538 12 side x 50 mm long. 1000537 14 side x 40 mm long. 14 side x 55 mm long. 1000539



STARHEAD STIR BAR

Circular, encased in heat resistant PTFE for temperatures up to 275 °C. Useful for tubes.ma-

trix flasks and round bottom

Part No. 25 Ø x 15 mm high. 1000790 1000791 40 Ø x 17 mm high.



BOX OF 17 ASSORTED STIR BARS

Box with 17 assorted stir-bars encased in temperature resistant PTFE for temperatures up to 275 °C. Contents:

1 stir-bar with pivot ring 6x10 mm.

1 stir-bar with pivot ring 4 x12 mm.

1 stir-bar with pivot ring 7x20 mm.

2 stir-bars with pivot ring 8x32 mm.

2 stir-bars with pivot ring 8x40 mm.

2 stir-bars with pivot ring 10x45 mm.

1 triangular stir-bar of 9x35 mm.

1 triangular stir-bar of 14x40 mm.

1 triangular stir-bar of 12x50 mm.

1 triangular stir-bar of 14x55 mm.

1 stir-bar Rotor-Pat of 11x35 mm.

1 oval stir-bar of 16x30 mm.

1 oval stir-bar of 20x40 mm.

1 star head stir-bar of 20x15 mm. 1 star head stir-bar of 40x17 mm.

Part No. 1000487



STIR BAR "ROTOR-PAT"

Encased in PTFE useable up to temperatures of

Excellent stability with strong vortexing turbulence, useful for round bottom or convex bottom flasks. Part No.

1000001 10 Ø max. x 32 mm long.



OVAL STIR-BARS

Encased in PTFE for temperatures up to 275 °C. Can be used in all types of vessel and recommen-

ded for round bottom flasks. Part No.

1000788 16 Ø x 30 mm long. 1000789 20 Ø x 40 mm long.



FLOATING MAGNETIC STIR-BAR "IMANFLOTER"

Encased completely in PTFE. Dismountable structure, autoclavable. Designed specifically for low constant speed stirring. Recommended for tissue culture applications.

The floating conception of the Imamfloter minimises the grinding effect of cell damage as there is no friction on the bottom of the flask.

Medium speed levels prevent the IMANFLOTER from jumping within the vessel.

Stir bar length Part No. Total length 1000897 50 mm. 39 mm. 1000898 60 mm. 45 mm.



STIR-BAR RETRIEVER Flexible, encased in PTFE.

Part No.

1000020 Length: 280 mm.

OVERHEAD STIRRER ACCESSORIES "RZR"



STIRRING SHAFTS WITH ROTORS

Shaft with AISI 304 stainless steel rotor. Shaft: 7mm Ø.


Part No.	Rotor type	Ø shaft mm	Total length mm	Features
6002891	1	90	400	Double fly vane
6002892	2	60	400	Fly vane
6002893	3	50	400	Fixed blade
6002894	4	50	400	Double cross vane
6002895	5	58	400	Helix
6002896	6	70	450	Six hole paddle
6002897	7	70	450	Three hole paddle
5050921	8	46	400	Radial turbine
5050932	9	45	400	Turbo helix 45 mm Ø
5050933	10	66	400	Turbo helix 66 mm Ø



SHAFTS WITH PTFE PADDLES

Shaft of 8 mm Ø.

Part No.	Rotor type	Ø shaft mm	Total Length mm	Features	Speed range r.p.m.
5050934	11	50	500	Helix PTFE	500
5050935	12	80	500	Anchor PTFE	500



ADJUSTABLE HEIGHT SHAFT STIRRERS

Shafts made of AISI 304 stainless steel coated in PTFE and PTFE stirrer paddles. Shaft: 8 mm Ø. Adjustable height paddles using a locking thread.

Can be used with differing paddles on the same shaft to create a multiple mixing action.

Part No.	Rotor type	Ø shaft mm	Total Length mm	Features
5050936	13	70	400	Shaft PTFE
5050937	14	40	400	Helix PTFE
5050938	Shaft P	「FE Spare or	additional paddl	e for multiple stir effect
5050939	Helix PT	FE Spare or a	ndditional paddle	for multiple stir effect









- Shafts with rotors made from AISI 316Cb stainless steel. Shaft: 10 mm Ø.

0 1 0 T. I

- Shaft made from AISI 316Cb stainless steel with POM plastic rotor. Shaft: 10 mm Ø.

Part No.	Rotor type	Ø shaft mm	lotal Length mm	Features	Speed range r.p.m.	
505098	50 15-S	60	500	Doble cone-S Stainless	200 - 800	
505098	51 16-M	80	500	Triple cone-M Stainless	200 - 700	
505098	52 17-M	80	500	Triple cone-M Plastic	200 - 700	
505098	5 3 18-L	120	500	Triple cone-L Stainless	120 - 500	
505098	54 19-L	120	500	Triple cone-L Plastic	120 - 500	

VISCO-JET®

A suitable system for all types of stirring.


Dynamic stirring with turbulent flow on vessel walls and bottom.

Homogenised liquids and fluids of low viscosity.


The fluid is aspirated by the widest part of the cone and is expulsed by acceleration from the smaller \emptyset , creating in the middle a turbulent effect; prevents the formation of foam, clots, aeration and heating up of the sample, degasses the medium being stirred.

Reduces the vortex formed, for short stir periods.

Graphic effect of the incident fluid flow using differing paddle forms: paddle, helix, turbine and conical Visco-Jet®



Axial flow, diagram showing the agitation flow. Representative of rotors with blades angles of 90° in relation to the stirring axis. (Rotor type; Helix).



Radial flow, diagram showing the agitation within the liquid medium. Representative of rotors such as paddles and anchors that run parallel to the stirring axis



Radial flow, diagram showing the agitation within the liquid. Representative of radial turbine rotors that run parallel to the stirring axis.



Dynamic radial flow, diagram showing the principal of agitation through a cone within the liquid. Representative of the Visco-Jet*. Rotors. Excellent results for minimum power consumption.

Sterile stoppers



Features

Designed for microbiology and cell culture applications.

Made from very flexible cellulose, tapered to fit most sizes of: Erlenmeyer, matrix flasks and tubes.

Sterile to 200 °C.

Permeable to air.

Transpirable.

Hermetically adjustable by pressing on the top of the glass vessel.

RANGE OF SIZES 1:1



For Ø of 9 to 11.5 mm Part No. 1001478




For Ø of 12.5 to 14.5 mm Part No. 1001479



For Ø of 13.5 to 15.5 mm Part No. 1001480



For Ø of 17.5 to 18.5 mm Part No. 1001481



For Ø of 18 to 22 mm Part No. **1001482**



For Ø of 25 to 27 mm Part No. **1001483**




For Ø of 28 to 31 mm Part No. 1001484



For Ø of 31 to 35 mm Part No. **1001485**



For Ø of 36 to 42 mm Part No. **1001486**



For Ø of 44 to 48 mm Part No. 1001487



Caps for biological culture tubes "Sero-Tap"

Made of anodised aluminium with internal AISI 304 stainless steel spring inside the cap. Not hermetically sealed. Autoclavable at 1 bar.



Made of anodised aluminium with internal AISI 304 stainless steel spring inside the cap. Not hermetically sealed. Autoclayable at 1 bar.

motionity ocalea. Autobiavat	no at i bai.	
For tubes	Part No.	For tubes
Ø 12 - 13 mm		Ø 19 - 20 mm
Colour Blue	1019201	Colour Blue
Colour Silver	1019202	Colour Silver
Colour Gold	1019203	Colour Gold
Colour Garnet	1019204	Colour Garnet
Colour Black	1019205	Colour Black
For tubes	Part No.	For tubes
Ø 15 - 16 mm		Ø 24 - 26 mm
Colour Blue	1024261	Colour Blue
Colour Silver	1024262	Colour Silver
Colour Gold	1024263	Colour Gold
Colour Garnet	1024264	Colour Garnet
Colour Black	1024265	Colour Black
For tubes	Part No.	For tubes
Ø 17 - 18 mm		Ø 28 - 30 mm
Colour Blue	1028301	Colour Blue
Colour Silver	1028302	Colour Silver
Colour Gold	1028303	Colour Gold
Colour Garnet	1028304	Colour Garnet
Colour Black	1028305	Colour Black
	For tubes Ø 12 - 13 mm Colour Blue Colour Silver Colour Gold Colour Garnet Colour Black For tubes Ø 15 - 16 mm Colour Blue Colour Silver Colour Gold Colour Garnet Colour Black For tubes Ø 17 - 18 mm Colour Blue Colour Silver Colour Gold Colour Gold Colour Gold Colour Blue Colour Silver Colour Gold Colour Gold Colour Gold Colour Gold	Ø 12 - 13 mm Colour Blue 1019201 Colour Silver 1019202 Colour Gold 1019203 Colour Garnet 1019204 Colour Black 1019205 For tubes Part No. Ø 15 - 16 mm 1024261 Colour Blue 1024262 Colour Gold 1024263 Colour Gold 1024264 Colour Black 1024265 For tubes Part No. Ø 17 - 18 mm Colour Blue 1028301 Colour Gold 1028302 Colour Gold 1028303 Colour Garnet 1028304

Self-sealing parafilm ® M









APPLICATIONS

Parafilm [®] M, is ideal for laboratories because it is a thermoplastic, self-sealing film that enables to cling even to irregular surfaces, culture tubes, flasks, Erlenmeyers, bottles, Petri capsules and any kind of recipients with an adequate sealing, protected against moisture loss, spillage, pollution or evaporation.

FEATURES

Parafilm $^{\circledR}$ M, is a translucent waxy polyethylene film of 127 μ m width which slowly stretches to about three or four times its original length.

PROPERTIES

It is resistant to acids, gases and alkaline solutions as well as methanol and many common reagents. It resists operating temperatures from -45 $^{\circ}$ C to 50 $^{\circ}$ C. It can be sterilized by using radiation techniques as hydrogen peroxide.

MODELS

Part No.	Roll width mm.	Roll length m.
1012201	50	75
1012202	100	38
1012203	100	75

ACCESSORIES

Paraffin dispenser $^{\circledR}$ made of blue ABS plastic for 50 and 100 mm. width Parafilm $^{\circledR}$ rolls, with built-in cutter

Capacity: One or two 50 mm. rolls or one 100 mm. roll.

Part No: 1012200





POLYPROPYLENE TOPS FOR CULTURE TUBES "BACTI-CAPS"

General features:

Adaptable to the tubes by three flexible internal flaps that hold onto the edge of the tube, not a hermetic seal, for minimal sample evaporation.

Highly resistant to culture media, acids, alkalis, alcohol and ether.

Easy to clean, can be sterilized at 121 °C (1 atm, vapour pressure).

Not suitable for dry sterilization.

Come in 5 colours:

YELLOW, BLUE, WHITE, RED and GREEN.

Part No.	For tube Ø 12 - 13 mm	Part No.	For tube Ø 15 - 16 mm
5726013	Yellow Colour	5726016	Yellow Colour
5726113	Blue Colour	5726116	Blue Colour
5726213	White Colour	5726216	White Colour
5726413	Red Colour	5726416	Red Colour
5726513	Green Colour	5726516	Green Colour
Packet of 1	00 units.	Packet of 10	00 units.

Part No.	For tube Ø 18 - 19 mm
5726019	Yellow Colour
5726119	Blue Colour
5726219	White Colour
5726419	Red Colour
5726519	Green Colour
Packet of 10	00 units.

,	
Part No.	For tube
	Ø 24 - 25 mm
E70000E	Vallau Calaur
5726025	Yellow Colour
5726125	Blue Colour
FZOCOOF	White Octors
5726225	White Colour
5726425	Red Colour
EZOCEOE	Green Colour
5726525	areen Golour
Packet of 10	00 units.

Part No.	For tube Ø 37 - 38 mm	
5726238	White Colour	
Boxes of 25 units		



WEIGHING BOATS "ACCU-WEIGH"

Anti-static, Suitable for weighing with reproducible tare weight. Free from contaminants, Wide range of sizes.



General features:

Made from high impact polystyrene with a flat base, round edge, smooth surface, contamination free and hydrophilic.

Flexible and easy to manipulate, suitable for liquids and solids.

Supplied in two forms:

- Diamond, black or white colours.
- Rectangular, white colour.

Part No.	Shape	Capacity ml	Externa Height	al dimensi Width	ons mm Depth	Weight g	Colour	Quantity
5726004	Diamond	5	3.8	31	54	0.25	White	100 units
5726005	Diamond	5	3.8	31	54	0.25	Black	100 units
5726006	Diamond	30	13.5	55	85	1.00	White	100 units
5726007	Diamond	30	13.5	55	85	1.00	Black	100 units
5726008	Diamond	100	20	92	132	2.6	White	100 units
5726009	Diamond	100	20	92	132	2.6	Black	100 units
5726010	Rectangular	7	8.5	44	44	0.65	White	100 units
5726011	Rectangular	100	24	80	80	2.2	White	100 units
5726012	Rectangular	250	24	134	134	6.1	White	100 units





ALUMINIUM DISHES

Pressed with crenulated edge.

No static.

Withstands temperatures up to 450 °C.

Can be used as a lid by mounting one on top of another.

Ideal for sample testing of solids, humidity, mixtures of emulsions and moisture determinations by weight loss.

Suitable use as a weighing boat, cooling or evaporation of samples. Single use and disposable.

Part No.	Thickness mm	Height. / Ø (usable) mm	Quantity
5726000	0,2	5 42	packet of 112 units
5726001	0,2	10 57	packet of 110 units



GLAZED PORCELAIN, RESISTANT TO THERMAL SHOCK UP TO 1050 °C. BETTER CHEMICAL AND ACID RESISTANCE COMPARED TO GLASS. INCOMPATIBLE WITH HYDROFLUORIC ACID AND ALKALI SALTS.



GLAZED PORCELAIN TRAYS

Part No.	Height	Length	Base	Width	Base
	mm	mm	length mm	mm	width mm
5310001	15	50	42	28	17
5310002	15	60	52	30	20
5310003	18	87	80	60	50
5310004	18	120	110	60	45



GLAZED PORCELAIN EVAPORATION BASINS, ROUND BOTTOM WITH SPOUT

Part No.	Capacity	Length	Ø exterior
	ml	mm	mm
5310005	35	25	60
5310006	60	30	74
5310007	100	35	90
5310008	300	52	126
5310009	500	58	168



FLAT BOTTOM CRUCIBLE WITH HANDLE AND **SPOUT**

Part No.	Capacity	Height	Ø Top	Ø Bottom
	ml	mm	mm	mm
5310010	50	34	60	36
5310011	100	45	80	52
5310012	200	52	96	56
5310013	300	62	110	66
5310014	500	70	132	80



STAINING BEAKER INTERIOR AND BOTTOM **GLAZED**

Part No.	Capacity	Total	Ø	Support ring
	ml	Height mm	mm	Ø mm
5310028	500	145	90	100



CRUCIBLES MEDIUM HEIGHT WITH GLAZED LIDS CRUCIBLES, LOW FORM WITH LIDS,

CHUCIDL	ועשואו 6ם.	JIVI NEIG	ni wiin GL	AZEV LI
Part No.	Capacity	Height	Ø Bottom	Ø Top
	ml	mm	mm	mm
5310015	15	30	21	35
5310016	25	36	23	40
5310017	50	46	30	53
5310018	100	57	34	63



GLAZED

Part No.	Capacity	Height	Ø Bottom	Ø Top
	ml	mm	mm	mm
5310019	15	23	19	43
5310020	25	27	24	47
5310021	50	38	28	59



MORTAR AND PESTLE, WITH SPOUT, UNGLAZED INTERIOR AND PESTLE

Part No.	Capacity	Length	Ø Bottom	\emptyset superior	Handle	Ø Pestle
	ml	mm	mm	mm	length mm	mm
5310023	100	45	60	100	110	30
5310024	300	62	80	130	140	37
5310025	500	74	90	160	170	42
5310026	1000	90	114	216	200	48



SPATULA SPOON GLAZED Part No. Total Length mm 5310022 140

COMBUSTION BOATS, UNGLAZED

Part No.	Height	Width	Depth	Exterior
	mm	mm	Width mm	length ml
5310027	10	15	9	95

AGATE MORTAR AND PESTLE

Suitable for grinding and mixing, medium and hard substances.

Very resistant to abrasions and contamination.

Hardness: 7 Mohs scale.

Compression resistance: 110 Kg/mm² Breaking resistance: 210 Kg/mm²

Mortar and pestle polished externally and internally.

Models	Capacity	Ø/B	ase	Dimensions
Part No.	ml	(interio	r) mm	(exterior) mm
5310029	15	45	11	60
5310030	30	65	16	80
5310031	70	80	24	100
5310032	140	105	32	120



SPARES Pestle		
Part No.	For Mortar Part No.	Handle length mm
5310033	5310029	77
5310034	5310030	95
5310035	5310031	100
5310036	5310032	130

GREY IRON MORTAR WITH PESTLE

In grey iron.

Part No.	Ø usable	Height
1004713	130 mm	65 mm
1004715	150 mm	75 mm



COMECTA Sealer "305-C"

FEATURES

Sealer for any type of thermo-plastic material up to 0.4 mm thickness.

Easy to use: Pressure and release.

Manual cutting by slip cutting edge.

Energy saving system by a micro-switch, activated only when the sealing arm is in place.

Heat impulse for only 0.2 to 1.3 seconds.

Time controll according to the thickness of plastics to be seal.

A LED and audible beep are activated when the seal is complete.



MODEL

Part No.	Seal Width mm	Maximum seal length mm	Maximum plastic thickness mm	Seal time seconds	Height / Width / Depth (exterior) cm	Power W	Weight Kg
5800042	5	300	0.4	of 0.2 to 1.3	19 48 8.2	1100	5.2

Supplied complete with 2 heaters, 2 PTFE heater covers and cutting blade.

ADDITIONAL SPARES: Maintenance kit 305-C: 2 heaters, 2 PTFE heater covers and cutting blade Part No. Part No. 5800044



COMECTA Sealer "455"

HIGH THROUGHPUT SEALER.

FEATURES

Ideal for large size samples.

Sealer for any type of thermo-plastic material up to 0.8 mm thickness.

Easy to use by pressing a foot pedal which enables to have the hands free in order to manipulate the bag. Includes a platform of 45 cm wide x 12 cm depth to support the sample. This platform is adjustable in height up to 19.5 mm to adapt different volumes.

Energy saving system by a micro-switch, activated only when the sealing arm is in place.

Heat impulse for only 0.2 to 2.2 seconds.

Time controll according to the thickness of plastics to be seal.

A LED and audible beep are activated when the seal is complete.

Replaceable PTFE anti stick coated heaters.

The pressure of the pedal can be adjusted.

Includes a recoil mains cable at the back of the unit.

All metal robust construction.

Adjustable base level.





ΝИ		n	ы	
	ш	ш	ы	

Part No.	Seal Width mm	Maximum seal length mm	Maximum plastic thickness mm	Seal time, seconds	Height / Width / Depth (exterior) cm	Power W	Weight Kg
5800047	5	450	0.8	of 0.2 to 2.2	90 54 38	1000	29

Supplied complete with 2 heaters and a fuse.

ADDITIONAL SPARES: Maintenance kit 455: 2 heaters and fuse. Part No. 5800049



COMECTA Portable sealer "150-C"

FEATURES

Plastic sealer, PE, PP and thermoplastic material such as cellophane, Aluminium covered paper and craft paper up to 0.6 mm thickness.

Easy to use pressure clamp.

Parallel seal for both sides of the bag.

Superior heaters covered in PTFE that produces a parallel hermetic seal.

Two temperature selection switch 120 °C and 180 °C





MODEL

Part No.	Seal Width mm	Maximum seal length mm	Maximum plastic thickness mm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
5800038	14	150	0.6	9 17 18	50	0.6



Vacuum sealer "Sealcom-V"

OPTIMUM QUALITY SEAL AND PRESERVATION. SUITABLE FOR SEALING SOLIDS AND LIQUIDS.

APPLICATIONS

Vacuum seal for any type of thermoplastic material: polyethylene, polypropylene, polyester, aluminium etc. Ideal for sealing solids, liquids, powders, granules, seeds, fruit, food, medical apparatus, forensics, chemical products, sludge samples, mud, cores, etc. Prevents the oxidation, mould and decomposition. Prolongs the quality, durability and freshness of items to be stored.

FEATURES

Bench top. Stainless steel AISI 304 robust construction with transparent lid. Easy to use and program. After closing the door the complete process is automatic. Operates using a oil vacuum pump (type AV-30 or N32).

CONTROL PANEL

Vacuum gauge with scale 0 to -0.1 Mpa.

Sealer temperature adjustment, high, medium or low.

Push button emergency stop.

LED indicator of operation.

Digital selection and display of the following parameters:

- Vacuum time 0-99 seconds.
- Seal time 0 to 9.9 seconds
- Discharge time 0 to 9.9 seconds



MODEL	Part No.	Seal Width mm	Max. seal length mm	Pump rate m³/h	Vacuum limit Mpa	Height / Width / Depth (usable) cm	Height / Width / Depth o (exterior) cm	Power W	Weight Kg
SEALCOM-V	4100050	10	260	10	-0.1	5 28 32	28 32 52	370	35

Supplied with a maintenance kit that includes 2 heaters, 3 PTFE heater covers. 1 lid seal and a 250 ml bottle of oil.

ADDITIONAL SPARES

Maintenance kit for Sealcom-V: 2 heaters, 3 PTFE heater covers, 1 lid seal, Part No. 4100051 Vacuum pump oil: 1 litre bottle (Type AV-30). Part No. 4100052



Bag sealer "S-400"

SUITABLE FOR POLYPROPYLENE-POLYESTER STERILIZING BAGS WITH A PAPER BASE.

FEATURES

CONTROL PANELS

Metal epoxy covered case.

Illuminated mains switch.

Pressure sealer. Manual cut. Stainless steel AISI 304 reel support at Lamp "ON" when the thermostat has rea-

the back.

ched the pre-set working temperature.

MODEL

Part No.	Seal Width	Maximum seal	Height / Width / Depth	Power	Weight
	mm	length mm	(exterior) cm	W	Kg
4100041	9	270	12 35 21	200	5.3



SEALING MACHINE ACCESSORIES

Rolls of paper for sterilisation "Polisteril"

Paper base polypropylene and polyester bags showing differing ink prints of sterilized items.

Part No.	Bag width cm	Bag length
4100075	7,5	100
4100100	10	100
4100200	20	100



SEE A FULL RANGE OF BAG SEALERS FOR STERILIZATION AND ACCESSORIES ON PAGE 68

Automatic bag sealer "Seal-Com 2". Part No: 4100044



Bag sealer

"Seal-1000". Part No: 4101712





RECTANGULAR HOTPLATE "PLACFIN"



KOLLE MANDLES LOOPS Infrared sterilizers.



LABORATORY LIGHT



PIPETTE DISPENSERS AND AUTOMATIC DISPENSERS



PARAFFIN WAX FLOATING BATHS



BLOOD AND CELL COUNTERS



COLONY COUNTERS

Instruments for dissection and slide boxes
Magnifier screen
Hand microtome
Wax dispenser and paraffin wax floating
Cooling plates and baths
Stackable anatomy - Histology
Dispensing pipettes
Digital counter for blood and colony counter
Ultraviolet sterilizer and westergren supports

pages 213 to 216
page 216
page 217
page 217
pages 217 to 219
page 219
pages 220 to 222
pages 223 to 224
page 224



Made of AISI 410 stainless steel.

1 Scissors universal "Multi purpose" with a serrated/notched blade. Suitable for cutting hard plastic, cardboard and soft wire etc.

Cut length 50 mm. Part No. 1001468 Total length 180 mm.

2 Large laboratory scissors.

Part No. 1001469 Cut length 100 mm. Total length 250 mm.

3 Multipurpose straight scissors, general use and bandages.

Cut length 50 mm. Part No. 1001470 Total length 140 mm.

4 Multipurpose scissors, general use and bandages.

Cut length 40 mm. Part No. 1000806 Total length 145 mm. 5 Dissection scissors.

Part No. 1000807 Cut length 35 mm. Total length 125 mm.

6 Mini dissection scissors.

Part No. 1001471 Cut length 25 mm. Total length 110 mm.

7 Flat scissors with serrated edge, ideal for laminates,

negatives and chromatography.

Part No. 1001472 Flat length 40 mm. Total length 160 mm.



DISSECTION NEEDLES

Plastic handled, AISI 304 stainless steel needle.

Point Total length Part No. 140 mm 1000431 Lancet 1000275 140 mm Fine



DISSECTION NEEDLES

Hexagonal metallic handle with a AISI 304 stainless steel needle.

Part No. **Point Total length** 1001359 40 mm Lancet 1001360 Fine 140 mm



SCALPEL HANDLE No. 4

Made of AISI 304 stainless steel, comes with three sterile blades No. 20, 22 and 23,

Part No. 1001361

Spares:

Part No.

1001362 10 sterile blades nº 20 1001363 10 sterile blades nº 22

1001364 10 sterile blades nº 23



ANATOMY SCALPEL

Made of AISI 304 stainless steel. Blade length: 35 mm.

Part No. 1000809



AISI 304 stainless steel forceps for picking up swabs and textiles etc. Part No. 1001192





DISSECTION KIT

Model Student.

Contains: 1 fine point probe, 1 lancet probe, 1 blunt end forceps, 1 fine tip forceps, 1 scalpel handle no 4, 3 sterile blades for scalpel handle no 4, 1 scissors, 1 magnifying glass, 1 style probe, 1 scalpel.

Part No. 1001365



DISSECTION KIT

Wooden box.

Contents: 1 fine needle, 1 lancet needle, 1 dissection forceps blunt end, 1 forceps fine point, 1 scalpel handle no 4 with 3 sterile blades Nos. 20, 22 & 23, 1 scalpel and 1 scissors.

Part No. 1001152



DISSECTING FORCEPS AND GENERAL LABORATORY USE

Made of AISI 304 stainless steel with electro-bright finish.

Туре	Description	Total length	Part No.
1	Dissecting Straight forceps round tips	140 mm	1000580
2	Dissecting Straight forceps fine tips	140 mm	1000801
3	Dissecting curved forceps round tips	140 mm	1000865
4	Dissecting curved forceps fine tips	140 mm	1000864

5	Sprung straight forceps with curved fine tips with wood laminate sides	160 mm	1001190
6	Forceps straight serrated fine point, with wood laminate sides.	160 mm	1001189
7 8	Laboratory forceps curved point Laboratory forceps fine point	160 mm 120 mm	6000128 1001191




KOLLE HANDLES

Nickel metal handle with loop end of Ni-Cr. Handle length: 43 mm. Internal loop \emptyset : 3 mm. Part No.

1000033 Bakelite handle. Handle length: 200 mm.

1000488 "Pasteur", with metal handle. Handle length: 290 mm.



1. KOLLE HANDLE

Loop end.

Handle length: 43 mm. Loop Ø: 3 mm.

Part No.

1000013 Made of Ni-Cr.

1000014 Made of Platinum iridium (10% iridium).

2. NEEDLE FOR KOLLE HANDLE

Made of Ni-Cr.

Part No. 1000802

3. LANCET TIP FOR KOLLE HANDLE


Made of Ni-Cr. Total long: 50 mm.

Part No. 1000803



CTANI

Support rack for Kolle handles dissection needles made of hardened methacrylate. Capacity: 6 places.
Part No. 1000903



STERILE CULTURE LOOPS

Disposable flexible polystyrene culture handle, one loop and one pointed end.

Part No. Capacity Colour internal loop Ø In packets of 1000866 1 μ I green 1 mm 20 units 1000867 10 μ I blue 4 mm 20 units



INFRARED TIP STERILIZER "STERILBIO"

Applications

Sterilization of culture tips, needles and small instruments. Incinerates organic material at 900 $^{\circ}\text{C}$ for 5/8 seconds.

Features

Made from ceramic elements with a quality quartz tube interior.

Does not need purged oxygen and can be used in anaerobic environments.

Eliminates problems of aerosols and bacterial contamination.

Support at the side for 4 Kolle handles, needles, lancet or dissecting forceps etc.

Model

Part No.	Heigl	nt/Width	/ Depth	Ø/Le	ength	Power	Weight
	(exterior) cm		(tube int	erior) cm	W	Kg	
3000788	20	13	18	1.4	15	110	1.5




BURNER J-2000 UNIVERSAL

For Butane gas. Instant safety flame, push button ignition, incorporating a LOCK/ON mechanism that prevents accidental fires, the unit switches off when the users hand is removed. Gas filled, using simple lighter gas fuel, as used in common gas lighters. Long life flint. Power by a 1.5 V battery.

Applications: Microbiology, bacteriology, dentistry and for general use where a fine hot flame is required for fast heat and fusion.

Part No. 1001114



GLASS SLIDES

Slides, unground of 26 x 76 x 1.5 mm thickness.

Box of 72 units. Part No. 5600400


Slides of 26 x 76 x 1.5 mm with ground end of 26 x 20 mm for writing.

Box of 72 units. Part No. 5600399 Cover slide of 18 x 18 mm. Box of 100 units. Part No. 5600396

Cover glass of 22 x 22 mm. Box of 100 units.

Part No. 5600397 Cover glass of 24 x 24 mm.

Box of 100 units. Part No. 5600398



SLIDE DISPENSER

Made of ABS plastic with a transparent polycarbonate lid. With single slide manual extraction. Capacity 50 slides of 76 x 26 mm.

Part No. 1001366



RECTANGULAR HOTPLATE "PLACFIN"

Adjustable temperature from ambient+5 °C to 130 °C.

Stability ±3 °C. Not suitable for boiling.

Coated with chemically treated aluminium, enamelled outer case.

Part No. Width/Depth Height/Width/Depth Power Weight (plate) cm (exterior) cm 7000497 400 23 22 8 22 30

ACCESSORY: Slide support rack for Placfin


"Placfin" slide support adapter. Made of AISI 304 stainless steel. Capacity: up to 18 slide of 76 x 25 mm (in a horizontal or inclined position).

Part No. 7000498



Kg

3.7



SLIDE STAINING STAND


Extendible and made of AISI 304 stainless steel. Fits tanks and trays up to 50 cm in length. Capacity 12 to 14 slides inclined at an angle.

Part No. 6000904

ACCESSORY: Long rectangular instrument tray.

Made from AISI 304 stainless steel, 350mm long x 150mm wide.

Part No. 8080207



SLIDE SUPPORT RACK, FOR 26 x 76 mm SLIDES

Made from ABS plastic with transparent polycarbonate lid. Holds 100 slides in a vertical position, with numbered guides and index table.

Part No. 1001367



FORCEPS FOR LIFTING SLIDES

Made of AISI 304 stainless steel, with flat jaws. Length: 105 mm. Suitable for lifting, slides, watch glass or small flat object.

Part No. 1000105



LABORATORY LIGHT "LABOR-LAM"

Used for the classification and manipulation of blood images of; anti serum transparencies, albumin control, RH tests, colony counting, micro examination, colloidal rings, observation of precipitates, acid and oil filters etc. Swivel accessory for bench or wall mount bracket. Light: white light with diffusing methacrylate screen. Temperature 30 °C to 35 °C on the screen surface. Screen light: 300 mm long x 85 mm wide. External dimensions: 320 mm long x 130 mm wide x 55 mm high. Power: 8 W.

Part No. 3000438

ACCESSORY: Table top support. Swivel type. Part No. 7000483



STORAGE BOX FOR 26 x 76 mm SLIDES

Made of white ABS with a hinge lid and catch, numbered slots for slides and a cork base. An index for slide identification is placed inside the lid. Part No. 1001371 for 100 slides



SLIDE BOXES

Holds 5 slides, made of polypropylene. Ideal for personal use. Dimensions:82 mm long x 30 mm wide.

Part No. 1000869



Hand microtome

Cylindrical hand microtome, that can achieve a minimum cut of 10 $\mu m,$ using a manual advanced mechanism.

External milled knob with a micrometer scale of 50 divisions of 10 μm .

Maximum sample section of 20 mm.

Precise internal sample holder with a maximum aperture of 8 mm.

Part No. 5900019

Supplied complete with a carry case, dissection forceps, blade and a 5 mm thick glass base.




SLIDE BOX FOR STORING 26 x 76 mm SLIDES

PVC base, colour white, with numbered dividers.

Index in the lid interior.

Part No. 1052050 for 50 slides Part No. 1052100 for 100 slides



STORAGE BOX FOR SLIDES OF 26 x 76 mm

Made of polystyrene with numbered slots and index in the lid.

Part No. 1001368 for 25 slides
Part No. 1001369 for 50 slides
Part No. 1001370 for 100 slides



EPS SLIDE HOLDER FOR 5 SLIDES

Polystyrene slide holder designed for transporting and preservation of samples and the mailing of slides with safety, such as blood samples, serum and section samples. System of 5 divided sections to prevent rubbing against each other in transit. Up to 10 slides (75 x 24mm). A fitted lid holds the slides in place.

Dimensions: 53 mm height x 100 mm long x 63 mm deep.

Part No. 1001717 Packs of 25 units







Cooling plate "Plac-Center"

DIGITAL DISPLAY WITH ANALOGUE TEMPERATURE CONTROLLER.
ADJUSTABLE TEMPERATURES FROM -10 TO 5 °C.

APPLICATIONS

Designed to keep material in a solid state.

FEATURES

External case made of AISI 304 stainless steel. Cooling plate made of black anodised aluminium. Anti vibration, hermetically sealed CFC free compressor.

CONTROL PANEL

Illuminated Mains ON/OFF switch.

Analogue electronic temperature control with digital display (2 digits).

Resolution: ±1 digit. Temperature probe: PTC.

MODEL

Part No.	Width/Depth	Height/Width/Depth	Power	Power	Weight
	(plate) cm	(exterior) cm	HP	W	Kg
4000501	30 15	27 33 50	1/3	160	15





Wax dispenser "Dispenser"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5 TO 100 °C. STABILITY ±1.5 °C.

FEATURES

External metal case coated with epoxy. Internal tank made of aluminium coated with black indelible PTFE. Independently heated tap.

CONTROL PANEL

Illuminated Mains switch.

Indicator lamp shows when the heater is on.

Adjustable temperature thermostat controller.

Adjustable tap temperature controller.

MODEL

Part No.	Capacity	Ø / Height	Ø / Height	Power	Weight
	litres	(usable) cm	(exterior) cm	W	Kg
4000490	4	20 15	28 30	400	5





Paraffin wax floating out bath for histology samples "Termofin"

ADJUSTABLE TEMPERATURE FROM 40 °C TO 80 °C. STABILITY ±1 °C.

APPLICATIONS

Suitable for the fusion and preservation of paraffin wax in a liquid state that enables the users to see the floated sections.

FEATURES

Bath and lid made of black PTFE. coated aluminium with a location hole in the lid for a thermometer.

CONTROL PANEL

Mains ON/OFF switch. Adjustable temperature thermostat synchronised with an indicator lamp while heating.

MODEL

Part No.	Capacity	Ø / Height	Ø / Height	Power	Weight
	litres	(usable) cm	(exterior) cm	W	Kg
3000459	4	25 10	34 14	430	2.3


Supplied complete with lid and thermometer.

ACCESSORY

Thermometer support (can only be used without the lid).

Part No. 6004591

PTFE coated lid







STAINING TANK

O I A I I I I I I I I I I I I I I I I I	AIII						oon bormetically eagled lid
Part No.	Туре	ŭ	/ Width / cterior) m	•	Slide quantity	Applications	cap hermetically sealed lid.
1001138	1	85	105	70	10 of 76 x 26 mm	Tank with lid, ı	requires rack No. 2 that is placed inside.
1001139	2	30	90	48	10 of 76 x 26 mm	Tank rack, bot	tomless, submerge in No.1 tank.
1001140	3	105	60	55	16 of 76 x 26 mm	Hellendahl ver	tical tank with lid.
1001141	4	50	85	70	10 of 76 x 26 mm	Horizontal stai	ining tank with lid.
1001142	5	105	45	45	10 of 76 x 26 mm	Coplin, vertica	l staining tank with lid.
1001352	6	110	56	56	10 of 76 x 26 mm	Coplin, stainin	g tank made of polyethylene with lid.



1. ANGLE FORM MOULDS

For paraffin wax covers.

Supplied as two separate angular plates made of AISI 304 stainless steel.

Part No.	Height / Width / Depth	Part No.	Height / Width / Depth
	(usable mould) mm		(usable mould) mm
6051622	14 22 44	6051645	20 45 85
6051635	20 35 65	6051690	25 90 150

2. COMPACT MOULDS WITH LID

For tissue covered moulds.

Part No.	Height / Width / Depth	Part No.	Height / Width / Depth
	(usable mould) mm		(usable mould) mm
7051720	9 20 20	7051735	9 26 35
7051725	9 25 25	7051740	9 26 40

3. ROUND INSERT MOULDS

Made of perforated AISI 304 stainless steel with lid.

Part No.	Ø mm	Height mm	
1051837	37	8	



Cooling bath "Refricub"

ELECTRONIC CONTROL WITH DIGITAL DISPLAY ADJUSTABLE TEMPERATURE DOWN TO -25 °C.

APPLICATIONS

Suitable for maintaining paraffin wax and waxes in a solid form for long periods of time, etc.

FEATURES

Temperature homogeneity: ±2 °C.

Precision of temperature controller: 1 °C.

External case made of epoxy covered metal with AISI 304 stainless steel tank and top.

An anti vibration hermetically sealed CFC free compressor is mounted in the base of the unit, supplied with ventilated evaporator, regulating valve and condenser.

CONTROL PANEL

Illuminated mains indicator switch.

Electronic controller for temperature with two digit digital display with a resolution of: ± 1 digit. PTC probe.

Part No.	Capacity	Height/Width/Depth	Height/Width/Depth	Power	Power	Weight
	litres	(tank) cm	(exterior) cm	HP	W	Kg
4100519	18	20 30 31	80 46 44	3/8	350	22





Chromatography tanks

APPLICATIONS

For quantitative analysis. Useful for storing and protecting the plates by means of the hermetic upper lid.

FEATURES

Made of heavy thick glass, with a flat bottom, ground rims with rounded edges. Solid glass rectangular lid, smooth base with ground edges to make an air tight fit tank. Rectangular glass lid with handle, smooth base and ground surface for hermetically sealed closure with the tank.

MODELS	Part No.	Heigh/Width/Depth mm	N⁰ of plates	With Lid
MINICROM	1001162	139 135 80	1	Yes
CROMACROM	1001163	139 135 80	1	Yes
MULTICROM	1001164	240 236 118	5	Yes

The model MULTICROM has 5 plate slots to maintain the plate upright.





Stackable Anatomy - Histology sample storage chest

FEATURES

Manufactured in steel epoxy coated. Resistant to corrosive atmospheres.

STACKABLE MODULAR SYSTEM

Cabinet for storing slides.

Complete with 14 drawers for slides of 26 x 76 mm.

Drawer capacity: 350 slides.

Part No. 1052014

Cabinet for storing transparencies.

Complete with 7 drawers for transparencies of 50 x 60 mm.

Drawer capacity: 230 transparencies.

Part No. 1005207

Cabinet for storing paraffin wax mounted samples.

Complete with 10 drawers.

Part No. **1052010**


MODELS

Part No.	•	/ Width usable) o		Ŭ	/ Width xterior)	n / Depth cm	Weight Kg
1052014	9	2.6	45	14	48	48	9.5
1005207	7.2	5.8	47	12	48	48	6.5
1052010	2.6	43.5	47	34.5	48	48	18.5

Base and lid units

Mounted on the top or as a support base. Fits the individual modules to form distinct user definable cabinets.

Part No.		Ŭ	/ Width xternal)	/ Depth cm	Weight Kg
1005208	Base	5,5	48	48	2,5
1005209	Lid	5,5	48	48	2,5





Fixed volume dispensing pipettes "Dispen-Fix" with validation test according to serial number.

FEATURES

Fixed volume pipettes from 1 to 1000 µl, model dependent.

Simple and easy to use.

Discard tip ejection for all models.

Comes complete with calibration tool.

MODELS	Part No.	Capacity µI	Accuracy %	Precision %
Dispen-Fix	5904009	5	±2.0	< 0.8
Dispen-Fix	5904010	10	±1.6	< 0.5
Dispen-Fix	5904011	20	±1.0	< 0.3
Dispen-Fix	5904012	25	±1.0	< 0.3
Dispen-Fix	5904013	50	±0.8	< 0.3
Dispen-Fix	5904014	100	±0.8	< 0.25
Dispen-Fix	5904015	200	±0.8	< 0.15
Dispen-Fix	5904017	500	±0.8	< 0.2
Dispen-Fix	5904018	1000	±0.8	< 0.13



WITH VALIDATION TEST ACCORDING TO SERIAL NUMBER. **AUTOCLAVABLE UP TO 121 °C.**

FEATURES

Variable volume from 1 to 5000 µl, model dependent.

Simple and easy to use adjustable dispense volume.

Digital set volume selection.


Discard tip ejection for all models up to 1000 ul.

Comes complete with calibration tool.

MUDELS	Part No.	Capacity µl	Divisions µl	Precision %
Dispen-Var	5900010	1 - 10	0.1	< 1.5 - < 0.4
Dispen-Var	5900020	2 - 20	0.1	< 1.5 - < 0.8
Dispen-Var	5900100	10 - 100	1	< 1.0 - < 0.3
Dispen-Var	5900200	20 - 200	1	< 1.0 - < 0.5
Dispen-Var	5901000	100 - 1000	10	< 0.8 - < 0.28
Dispen-Var	5905006	1000 - 5000	10	< 0.4 - < 0.2

ACCESSORIES. PIPETTE TIPS AND RACKS. Accurate and precise, made of polypropylene.

Models Nº image	Part No.	Range µl	Colour	Qty. tips x part No.	Suitable for pipette part No.	Volumes pipette µl
	Tips 5490040	0.1- 10	Transparent	1000	5904009	5
1	Rack 5490035		Transparent	96	5904010	10
	11ack 0430000	0.1-10	папорагонс		5900010	1 - 10
2	Tips 5490041 Rack 5490036		Transparente Transparente	1000 96	5900020	2 - 20
					5904011	20
					5904012	25
	Tips 5490042 5 - 200	Transparent	1000	5904013	50	
3	Rack 5490037			5904014	100	
	11aux 3450037	J - 200	Transparent	ansparent 50	5904015	200
					5900100	10 - 100
					5900200	20 - 200
	5490044	101 - 1000	Blue	1000	5904017	500
4	5490034	101 - 1000	Blue	96	5904018	1000
	J45UUJ4	101 - 1000	Diue	50	5901000	100 - 1000
5	5490047	1000 - 5000	Blue	500	5905006	1000 - 5000
J	5490039	1000 - 5000	Blue	500	J 8UJUU0	1000 - 3000





Made of AISI 304 stainless steel. No. of places: 6 dosing pipettes up to 28 cm in length. Part No. 1001215

COMECTA

Original Comecta design





COMECTA: Multichannel variable volume dispensing pipettes "Multi-Var 8"

NUMBER OF CHANNELS: 8.

FEATURES

Due to its excellent tips alignment, they offer a great advantage in microtitrer techniques such as: agglutination, immunological radio assays, hemagglutination, antibiotic test, etc. (All models are calibrated according to DIN 12600 recommendations)

MODELS	Part No.	Tips Range μI	Volume A < :		Nomina CV	nt variation al volume <%	Increments µI	Tips types µl
			Min %	Max %	Min %	Max %		
Multi-Var 8	5905027	5-50	0.8	6	0.4	3	0.1	200
Multi-Var 8	5905028	10-100	0.8	4	0.3	2	0.2	200/300
Multi-Var 8	5905029	20-200	0.8	4	0.3	1.5	0.2	200/300



ACCESSORIES

Pipette tips rack. Part No. 5490037

Carousel support for a maximum 3 multichannel pipettes. Part No. 5490034



COMECTA Repeat multi dispenser "Repette-HW"

UP TO 49 SUCCESSIVE DISPENSES WITHOUT THE NEED TO REFILL. DISPENSING VOLUME FROM 2 µI UP TO 5000 µI.



FEATURES

Multi dispenser pipette.

Precise repeat dispenser.

Ideal for dispensing low and high viscosity liquids.

Designed for rapid dispensing.

Depending on the setting, the dispenser can discharge up to a maximum of 49 times successively, without the need to be refilled. Accepts 8 different syringes.

Precision table using COMECTA syringes. (These values are based on a maximum volume for each pipette using distilled water at ambient temperature):

Type syringe ml	Volume doses µl	Precision nominal volume A ≤ ±%		Variation coefficient nominal volume CV ≤ %	
		10%	2%	-10%	2%
0.5	10 - 50	0.8	4.0	0.6	1.4
1.25	25 - 125	0.8	4.0	0.3	0.8
2.5	50 - 250	0.7	3.5	0.2	0.8
5	100 - 500	0.5	2.5	0.2	0.6
12.5	250 - 1250	0.3	1.5	0.2	0.4
25	500 - 2500	0.3	1.5	0.2	0.4
50	1000 - 5000	0.3	1.5	0.15	0.4

ACCESSORIES

Steril syringes. Supplied in individual bag.

Part No.	Volume ml	Syringes x part No.
5490053	0,5	100
5490054	1,25	100
5490055	2,5	100
5490056	5	100
5490057	12,5	50
5490058	25	50
5490059	50	50

Dispense volumes from 2 µl up to 5000 µl, depending on syringe used.

Maximum viscosity per syringe is:

50 ml = 20 mPas, 5 ml = 260 mPas and 1.25 ml = 977 mPas.

Maximum pressure: 500 mbar.

Operating temperature: from 15 to 40 °C.


At the back of the dispenser shows a quick table reference of the volume and dose, that shows the repeatable dose from the syringe to be used and the correct slide position.

Includes an easy to locate wall or bench support.

Supplied complete with a 5 ml syringe and an autoclavable PP adapter for 25 and 50 ml syringes.

Repeat multi dispenser "Repette HW".

Part No. 5490050



SPAR	E	
8	5490060	Adapter made from PP for 25 and 50 ml syringes



Motorized pipette dispenser "Rota-Filler 3000®"

FOR ANY KIND OF PIPETTES FROM 1 TO 100 ML.

FEATURES

Three positions switch for fast or slow aspiration and dispensation or gravitatorial dispensation for free dispensing.

The pressure over the ergonomic bulbs allows modulation, speed control of aspiration and dispensation.

It is supplied with a couple of 0,45 µm filter by default for and effective protection against contaminated air either of the equipment during aspiration or of the sample during dispensation.

Two rechargeable and interchangeable batteries.

Low-battery light indicator signal.

Can be used while recharging.

TECHNICAL DATA

Charger of 65/12 mA. Charge time: 14 hours.

Weight: 215 g.

Dimensions: 170 height x 47 width x 143 length mm. Pipettes: from 1 to 100 ml (glass and plastic).

Motorized pipette dispenser "Rota-Filler 300®".

Part no: **5900700**

COMPLEMENT

5 hydrophobe filters pack of 0,45 µm.

Part no: 5900701



RF3000



Manual pipette dispenser "Pipet"

FOR ANY KIND OF PIPETTES FROM 1 TO 100 ML.

FEATURES

Ergonomic, comfortable, solid, lightweight and ease of use.

It has been designed with an ultracompressor silicone bulb that provides a soft manual control of the pipette, either in aspiration or dispensation.

Once the bulb is squeezed, control with the toggle both speeds, the aspiration and dispensation one. Push the middle button to expel blowing-out any liquid remaining in the pipette. $0.45 \mu m$.

Measure: 25 mm Ø. Weight: 130 g.

Manual pipette dispenser "Pipet". Part no 5900800

COMPLEMENT

5 hydrophobe filters pack of 0,45 µm.

Part no: 5900801





Measuring range "COM-PIP"

FEATURES

Made of polypropylene with silicon joints, for volumes up to 25 ml. Precise and safe, ergonomic design that allows filling and emptying without having to change the setting. Flexible thermoplastic seal against any type of burette. A button on the side empties the burette automatically. Easily dismantled and cleaned.

Part No. Capacity 1001131 0 to 2 ml 1001132

0 to 10 ml 0 to 25 ml

Identification colour

Blue

Green

Red

1001133



Made from AISI 304 stainless steel. Capacity 6 pipettes of differing sizes. Part No. 1001702







Digital cell counter for blood components "CC-1431"



FEATURES

Designed for any type of cell or type counting and blood analysis. Languages: Spanish, English and French.

Operating menu displayed on the screen, 4 LCD lines of 20 characters per line.

Visualization of data can be displayed as partial, total, and absolute or percentage of key strokes.

Possibility to correct the last count.

Memorizes the last pulsed count and identifies it on the screen.

Audible pulse control when the key is pulsed, advises when the count limit is reached. Up to 14 separate key inputs per program.

Power: 12 V. Includes mains adapter.

USB interface for download of results to a computer.

Includes software for PC with the following parameters: Patient Number, Laboratory name, Date, Hour, Comments and comparison of previous results after treatment.

FOUR PROGRAMS INCLUDED:

LEUC (leukocytes)

- Leukocyte count.
- 14 different cell types.
- 3 Keys for counts that are not used in the final count: («11. Normo» «12. Tox Gr» «14. Ery»).
- Total cell count (at the end of counting) selectable from 100 to 1000 in steps of 100. MYEL (myelocytes).
- Myelocytes count.
- 14 different cell types that accumulate to the total count.
- Total cell count selectable from 100 to 1000 in steps of 100.

RETI (reticuloecytes/erythrocytes)

- Reticuloecytes count.
- No. of cells by configurable view from 50 to 250.
- 1 Key for reticulocyte fields.
- 1 Key for erythrocytes that are not included in the final count.
- Total (end of count) fixed at 1000.

FREE (program, customer set)

- Configurable count.
- Each key can be activated or deactivated for the count.
- Total cell count (at the end of counting) selectable from 100 to 1000 in steps of 100.

KEYS

Two control keys:

- SELECT: choose one program or one option.
- ENTER: accept the option choose.
- 14 Keys for different cell types:


- Seg:	Segments or granulocytes.	- Myelo:	MieloBlasts.
- Eos:	Eosinocytes.	- At Ly:	Linfocytes untypical.
- Lym:	Linfocytes.	- My Bi:	Mielococytes Binucleads.

- Normo: NormoBlasts. - Mono: Monocytes. - Tox. Gr.: Granulocytes toxic. - Stab: Bands. - Baso: Basocytes. - Reti: Reticulocytes.


- Met My: Metamielicocytes. - Ery: Eritrocytes / EritoBlasts.



MODEL			
Part No.	Height / Width / Depth cm	Power W	Weight Kg
7121431	7 15 22	6	0.6



Report.



Software download and report.



APPLICATIONS

Molecular Biology, Genetics and microbiology.

FEATURES

The unit records each pulsed colony count on the digital screen to a total of 999 counts. Pressure to count can be made by any way (pencil, or pen tip etc.).

Epoxy coated casing. Reversible interior background part in white or black.

Steady, (glare free) fluorescent illumination, clear and wide view visibility.

Removable plates for housing varying sizes of Petri dishes.

Location for reversible contrast plates.

Magnifying glass X2 with adjustable height and rotation support

CONTROL PANEL

Mains switch.

Numerical digital display: 3 digits. Audible confirmation beep for each count. Zero reset button. RESET.





MODEL

Part No.	Ø of magnifying glass cm	Magnification	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4905000	12.5	2x	25 23 30	20	3.5

ACCESSORY 5 Wolffhuegel discs. Part No. 4905001

Supplied complete with Wolffhuegel disc.



Bench top ultraviolet radiation sterilizer

Fast destruction of germs by U.V. radiation based on the principle that albuminoid bodies absorb energy that which prevents germination and reproduction of the bacteria. Suitable for destroying: bacterias, virus, spores, yeast and mildew.

FEATURES

Main body made of expanded polystyrene, with a polished AISI 304 stainless steel tray that reflects the U.V radiation and a protective smoke coloured methacrylate top that allows the samples to be seen and at the same time, it is opaque to U.V. light. U.V. Wave length: U.V. 2537 Angstroms.

U.V. lamp life span: 5000 hours approx.


CONTROL PANEL

Illuminated Mains switch.

To prevent accidental sight or skin damage a safety micro-switch that cuts off the power if the U.V. protective cover is removed.



Part No.	Height/Width/Depth	Height/Width/Depth	Power	Weight
	(interior) cm	(exterior) cm	W	Kg
3000416	7 to 9 50 20	17 57 38	15	7.5



SPARE

Germicidal tube. Part No. 0022113



Westergren sedimentation pipette support stand

Sedimentation tube support stand.

Made of AISI 304 stainless steel.

Supports pipettes with a replaceable vial type rubber plug. Easy to clean.

Sectioned, tilted, 45 °C front panel scribble board for notes.

Designed so that the tubes can be mounted only with one hand.

Part No.

1000093 Three Place. 1000096 Six place.

1000912 Twelve place.





HEATING MANTLES

Two position heating Heating mantles with digital electronic regulation Heating mantels made with a double woven mineral wool Banks of heating mantles for Kjeldahl and soxhlet	page page page	227 228 229 230
	page	
Flexible heatings	page	231


HOT PLATES

Ceramic top hotplates	page	232
Round hotplates	pages	233 to 234
Rectangular hotplates	pages	234 to 235
Bank of circular hotplates and for kjeldahl and soxhlet	page	236

HEATING MANTLES FOR FLASKS

SAFFTV

FOLLOWING CE STANDARDS, IN THE CASE OF LIQUID SPILLS AN EARTH WIRE IS WOVEN INTO THE MANTLE TO PREVENT ELECTRIC SHOCKS.



COMMON FEATURES

- 1. Hand woven mineral wool with flexible braiding, which adapts snugly around the flask, with drain hole, replaceable .
- 2. The heating element is distributed homogeneously throughout the mantle. Connectors and cables are protected by mineral fibre tubes. Efficient low power consumption of 1 w/cm². Maximum mantle temperature of 450 °C.
- Double walled fibre wool thermally insulated casing.
 Pure Nickel connection terminals protected by mineral fibre tubing that prevent the deterioration of the contact due to high temperature.

Support for retort bar.

APPLICATIONS

For digestion distillation, extraction, evaporation and boiling processes.



MODELS SELECTION SUMMARY TABLE

Capacity ml	Ø flask mm	Fibroman C 2 positions	Fibroman D Woven	Fibroman O Funnel hole	Agiman Woven	Fibroman HT-W Electronic regulator	Micro for Kjeldahl banks banks	No. Places	For Kjeldahl/Soxhlet banks	No. Places
100	65	3003141	-	-	-	3031470	6014204	4	-	-
100	65	-	-	-	-	-	6014206	6	-	-
250	83	3003142	3004192	-	3000447	3031471	-	-	-	-
500	102	3003145	3004195	-	3000448	3031472	-	-	6003294	4
500	102	-	-	-	-	-	-	-	6003296	6
1000	132	3031410	3041910	3049110	3000449	3031473	-	-	-	-
2000	164	3031420	3041920	3049120	3000450	3031474	-	-	-	-
3000	190	3031430	-	-	-	-	-	-	-	-
5000	222	3031450	-	3049150	-	-	-	-	-	-
6000	235	3031460	-	-	-	-	-	-	-	-
10000	290	3314100	-	3491100	-	-	-	-	-	-
20000	365	3314200	-	3491200	-	-	-	-	-	-



With power heating switch selector

FEATURES

Two position isolated rocker switch heating control with power on indicator lamp.

TECHNICAL DESCRIPTION


- DUR ALLOY aluminium exterior case with epoxy coating.
- Suspended mineral wool woven heating mantle.
- Heating element homogeneously distributed throughout the woven mantle.
- Thermally insulated case, mineral wool fibre insulation.
- Nickel connectors.
- Retort stand clamp at the back of the unit.
- To prevent breakage or spills a drain hole is provided at the bottom of the mantle.
- Safety earth connection.

MODELS

Part No.	Flask capacity ml	Flask Ø approx. mm	Ø / H (exteri	•	Power W	Weight Kg
3003141	100	65	16	11	130	1
3003142	250	83	18	11	130	1,1
3003145	500	102	20	12	270	1,2
3031410	1000	132	22	13,5	410	1,4
3031420	2000	170	26	14	530	2
3031430	3000	190	29	18	620	2,2
3031450	5000	222	33	19	840	3,2
3314100	10000	290	38	22	1400	4,7
3314200*	20000	365	48	26	2300	11

^{*} The mantle part no. 3314200 supplied without the 2 position controller and switch. See Accessories on (page 230) and Regulators and controllers (pages 293 and 294).











Heating mantles for flasks "Fibroman HT-W" FOR ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C TO 400 °C.

STABILITY ± 0.5 to 1 °C according to temperature.

ELECTRONIC REGULATION OF TEMPERATURE WITH DOUBLE SCREEN TO SHOW THE PRESENT TEMPERATURE AND THE PROGRAMMED ONE.



FEATURES

External case made of AISI 304 stainless steel. Heating elements homogeneously distributed inside of the mineral wool woven and thermically insulated. Pure nickel connection terminals.

Temperature control system by inner thermocouple. TC probe connector included in the equipment for the liquid temperature control inside the flask.

Device to hold support bar.

CONTROL PANEL

On / Off switch.

Digital display indicating set and real temperature.

Heating operating led.

Push button A increase temperature.

Push button ▼ decrease temperature.

Side: TC probe connector for the liquid real tempera-



Device that allows to hold the rod has supported.

MODELS

Part No.	Flask Capacity ml	Ø Flask mm	Height / Width / Depth (exterior) cm	Ø Lenght Sensor TC mm	Stability °C	Power W	Weight Kg
3031470	100	65 to 67	13 22 26	4 170	±0.5 to 1	130	1,5
3031471	250	83 to 87	13 26 29	4 170	±0.5 to 1	216	3
3031472	500	102 to 105	13 26 29	4 170	±0.5 to 1	220	3,5
3031473	1000	132 to 139	14 30 33	4 170	±0.5 to 1	530	4
3031474	2000	164 to 170	16 31 34	4 200	±0.5 to 1	650	5

Note: the heating mantles are supplied with TC probe connected to the case.

ACCESSORIES



Support bar Dur-Aluminium, 12 mm Ø x 700 mm long. Part No. 6000270



Double bosshead. Part No. 1000069



Spring clamp to hold thermometers. Part No. 7000025



Clamp to support and hold condensers or other glassware. Part No. 70001161



Heating mantle for flasks "Fibroman-D"

FEATURES

Soft mantle made with a double woven mineral wool. A mains switch and safety earth are also fitted. The mantle can be supported by a retort ring or tripod stand.

MODELS

Part No.	Flask capacity ml		leight ior) cm	Power W	Weight Kg
3004192	250	15	6.5	130	0.8
3004195	500	17	7	270	0.9
3041910	1000	20	8.5	410	1
3041920	2000	24	11	530	1.5

See accessories (page 230) and controller (pages 293 and 294).

ACCESSORIES							
Retort Part No.	Tripod Ring Part No.	Mantle Part No.					
7041921	7041922	3004192					
7041951	7041952	3004195					
7419101	7419102	3041910					
7419201	7419202	3041920					





Heating mantle for flasks and funnels "Fibroman-0"

WITH TRIPOD SUPPORT STAND.

APPLICATIONS

Suitable for round bottom flasks, decanting funnels and reaction flasks with a bottom drain.

FEATURES

- Woven mineral fibre mantle.
- Epoxy coated exterior.
- On/Off switch.
- Safety earth.

MODELS

Part No.	Flask capacity ml	Ø / He (Exterior d	•	Ø cm Bottom hole	Tripod Height cm	Power W	Weight Kg
3049110	1000	24	10.5	8	34	410	1.5
3049120	2000	28.5	12	8	36	530	2.5
3049150	5000	37.5	17.5	9	38	840	3.5
3491100	10000	43.5	20	9	43	1400	5
3491200	20000	48	25	12	36	2300	13

Due to its capacity the part no. 3491200, 20000 ml capacity is manufactured in a cylindrical case. See accessories (page 230) and controller (pages 293 and 294).





Heating mantle for flasks

ADAPTABLE FOR THE "AGIMAN" STIRRER UNIT.

MODELS

Part No.	Flask capacity ml	Ø / H (Exterior d		Power W	Weight Kg
3000447	250	15	6.5	130	0.8
3000448	500	17	7	270	0.9
3000449	1000	20	8.5	410	1
3000450	2000	24	11	530	1.5

ACCESSORIES

Mantle magnetic stirrer "Agiman". With interchangeable mantle adapters. Electronic control of stir speed (up to 1600 r.p.m.).

Connection at the back for an electronic contact thermometer "Sensoterm" or digital electronic regulator "Electemp-TFT".

See controller (pages 293 and 294).

Part No. 7000446



Stirring mantle with adapter device.



Stirring mantle, complete.



Banks of heating mantles for digestion, Kjeldahl, Soxhlet extraction and round bottle flasks

CONTROLLABLE MANTLE TEMPERATURE UP TO 550 °C.

FEATURES

Independent heating control for each position.

2 height and angle adjustable bars to support condensers or other glassware for diquestion or extraction.

External case made of AISI 304 stainless steel.

Safety earth strip fitted.

CONTROL PANEL

Temperature control with a pulsed energy regulator. An indicator lamp shows when the mantle is on.

MODELS

Part No.	For flasks	Flask volume ml	No. H of places		Width terior)		Power W	Weight Kg
6003294	Kjeldahl / Soxhlet	300-500	4	18	76	25	1500	12
6003296	Kjeldahl / Soxhlet	300-500	6	18	105	25	2200	16
6014204	Micro-Kjeldahl	100	4	16	59	21	510	9
6014206	Micro-Kjeldahl	100	6	16	81	21	900	12

Supplied with:

2 bars and 4 clamps for Kjeldahl/Soxhlet models.

1 bar and 2 clamps for Micro-Kjeldahl models.

ACCESSORIES (DIGESTION)



Fume extractor tube clips made of AISI 304 stainless steel.

Part No.

6000008 For bank of Macro-Kjeldahl. 6000009 For bank of Micro- Kieldahl



Part No. **1000104**

For: Bank Macro-Kjeldahl,

4 positions.

1000106 For: Bank Macro-Kjeldahl,

6 positions.

1042044 For: Bank Micro-Kjeldahl,

4 positions.

1042066 For: Bank Micro-Kjeldahl,

6 positions.



Kjeldahl flask support made of AISI304 stainless steel. Flasks 300-500 ml.

Part No. 1000002


ACCESSORIES EXTRACTION



Double nut boss head. Part No. 1000069



Clamp for holding condensers. Part No. 7000403



EXTRACTION AND NEUTRALIZATION OF GASSES.

Reduces water wastage.



"Scrubber" Re-circulation pump

Designed to neutralize and absorb fumes generated in the Kjeldahl digestion process

Composed of a "Scrubber unit" that blocks the passage and neutralizes the condensed acids, and a re-circulating water pump that produces sufficient vacuum to aspirate the vapours.

No mains water connection required. Prevents waste contamination. Low noise level (<65dBA).

Re-circulating pump made from chemically resistant materials.

It is essential to put the "Scrubber" unit together with the neutralizer solution between the digestion block and the re-circulating pump.

"Scrubber" unit

Part No.	Height / (ex	Width terior)		Weight Kg
4001611	32	31	16	2

It is supplied with 3 Kg. acid vapours neutralizer solution.

Re-circulating water vacuum pump.

Part No.	Height / \((exte	Width . erior) o		Vacuum level bar	Pump rate litre/minute	Weight Kg
4001612	44	39	28	0.98	10	10

Spare part

3 Kg. acid vapours neutralizer solution. Part No. 4001610



Heating tapesHEATING TAPE SURFACE TEMPERATURE UP TO 210 °C. CONFORMS TO VDE EN 60335-A.



FEATURES

Heating elements distributed homogeneously between flexible moulded sheets of

APPLICATIONS

Used for heating small areas such as heating ducts for distillation, valves and pipes etc.

MODELS

Part No.	Length m	Width mm	Power W
3003178	0.80	25	180
3031710	1	25	210
3031715	1.5	25	265



Flexible heating cord



FEATURES

For tape surface temperatures up to 180 °C. Thermally insulated using flexible silicon tube.

APPLICATIONS

For heating small areas such as heating ducts, distillation, valves and pipes etc.

MODELS

Part No.	Length m	Ø Cord mm	Power W
1000781	1.5	3	45
1000782	3	3	90
1000783	6	3	180
1000784	12	3	360



Large surface area heating mats

MAT SURFACE TEMPERATURES UP TO 210 °C. CONFORMS TO VDE STANDARD EN 60335-A.



FEATURES

The heating element is insulated and distributed homogeneously between two sheets of flexible moulded silicon.

Heating Power: 0.39 W cm².

Mat surface temperature up to 210 °C. Held in place by hooked spring fittings.

APPLICATIONS

Suitable for heating drums up to 200 litres \emptyset 60 cm. Total: 3 bands of 171 cm in length x 15 cm wide.

Power consumption 3000 W.

Part No. 4000200

ACCESSORIES

Controller (pages 293 and 294).

HOTPLATES





Ceramic top hotplates "Ceramic-plac" ADJUSTABLE HOT PLATE TEMPERATURE UP TO 400 °C. MAXIMUM SUPPORTABLE WEIGHT: 40 KG.

SAFETY

CONFORMS TO THE EN 61010 STANDARD, OVER TEMPERATURE CUT-OUT FITTED.

Excellent resistance to chemicals. Unaffected by thermal shock



FEATURES AND ADVANTAGES

The "CERAMIC-PLAC" has the following features and advantages:

- Flat, smooth surface.
- Excellent resistance to chemical agents.
- High heating power within the heated area that permits extreme temperature within a short heat up time.
- Excellent thermal qualities and penetration of infrared radiation.
- Easy to clean.

FEATURES

- External case of AISI 304 stainless steel.
- Heating elements mounted on ceramic fibre thus making an insulated compact block.
- Temperature energy control from 6 % to 70 % of nominal power. (Model Part No. 3000919 and 3000920).
- Hydraulic thermostat temperature control up to 400 °C (Model Part No. 3000921).

MODELS

Part No.	Heating Zones	Heating area in cm	•	/ Width exterior)		Power W	Weight Kg	Control system
3000919	1	18 Ø	10	28	33	1800	5.5	Energy controller
3000920	2	14 Ø + 14 x 14 🗹	10	27	47	2200	7.5	With individual energy controllers for each zone
3000921	1	30 x 30 ⊄	10	49	49	3700	11	With hydraulic thermostat



Round Hotplates

HOTPLATE TEMPERATURE UP TO 400 °C. WITH 6 PLACE HEATING POWER CONTROLLER.

FEATURES

Cast steel plate with embedded heating elements throughout the plate surface. 🖎

Fast heating rate at maximum power setting.

Enamelled outer casing that is resistant to high temperatures and corrosive agents.

CONTROL PANEL

Controller 6 position power settings.

MODELS

Part No.	Plate Ø cm		/ Width exterior)	n / Depth cm	Power W	Weight Kg
1000718	12	8	20	20	1000	2.1
1000442	15	8	22	22	1600	3
1000443	19	8	26	26	2100	4.5



Hotplate - Circular "Selectemp" ADJUSTABLE HOT PLATE TEMPERATURE FROM 45 °C TO

ADJUSTABLE HOT PLATE TEMPERATURE FROM 45 °C TO 250 °C. Stability ± 2 °C.

FEATURES

Temperature control using a thermostat located in the center of the plate. Non warping steel plate with embedded heating element throughout the plate's surface.

Enamelled outer casing, resistant to high temperatures and corrosive agents.

CONTROL PANEL

Hydraulic thermostat temperature control synchronised with an indicator lamp. 12 levels regulable button.

MODEL

Part No.	Plate Ø	Height / Width / Depth	Power	Weight
	cm	(exterior) cm	W	Kg
1000444	15	8 22 22	1500	3



ACCESSORY

Circular sand, oil or water bath made of AISI 304 stainless steel, suitable for the "Selectemp" hotplate. Dimensions 200 mm dia x 115 mm high.

Comes complete with a lid with concentric reduction rings.

Part No. **7000008**





Hotplate double circular "Duplac"

FOR PLATE TEMPERATURES UP TO 400 °C.

FEATURES

Hot plate with non warping cast steel plate with embedded heating elements throughout the hot plate surface.

Independent plate temperature control is maintained by pulsed energy controllers. The exterior is epoxy covered steel.

CONTROL PANEL

Two temperature energy controllers.

Two indicator lamps function when the heaters are on.

MODEL

Part No.	Plate Ø	Height / Width / Depth	Power	Weight
	cm	(exterior) cm	W	Kg
4000634	19	15 47 24	1900	8.5





Hotplate-Rectangular "Combiplac"

PLATE TEMPERAURE ADJUSTABLE UP TO 400 °C.

FEATURES

Electronic energy temperature control.

The non warping steel hotplate has heating elements embedded throughout the hotplate surface maintaining the heat distribution. The main case is made of AISI 304 stainless steel.

CONTROL PANEL

Main power switch with "ON" indicator lamp. Electronic energy control from 0 to 100 %.

An indicator lamp shows when the unit is hot.

Connector at the back for the external digital thermometer, "Sensoterm" or "Electemp-TFT".

See Regulators and Controllers (pages 286 and 287).

MODELS

Part No.	Width / Depth (plate) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3000156	20 40	12 22 55	2300	9
3000718	30 30	13 33 46	2800	13



فنننننننننن

Sand bath adapters. Made entirely of AISI 310 stainless steel, removable for easy cleaning.

Part No.	Hotplate Part No.	Height/ (Sand 1		
6000719	3000156	5.5	20	40
6000720	3000718	5.5	31	31





Hotplate-Circular "Combiplac"

PLATE TEMPERATURE ADJUSTABLE UP TO 400 °C.

FEATURES

Electronic energy controller and embedded heating element throughout the non warping steel hotplate surface.

A retort bar can be located in a hole at the back of the unit.

The main case is made of AISI 304 stainless steel.

CONTROL PANEL

Main power switch with "ON" indicator lamp.

Temperature controlled by a 0 - 100 % energy controller.

A lamp indicates when the hot plate is "ON".

A Thermometer "Sensoterm" or "Electemp-TFT" that can be added and are listed as accessories.

See section Regulators and Controllers (pages 293 and 294).



MODELS

Part No.	Plate Ø cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3002301	12	10 16 27	550	2
3002302	19	13 21 30	950	3

ACCESSORIES

Suitable for the 12 cm \emptyset hotplates only Part No.3002301:

Bath of AISI 304 stainless steel of 140 mm Ø x 75mm high, with concentric reduction ring lid, Suitable for water, oil or sand. Part No. 7002510

Round bottom flask adapter made of aluminium and suitable for Kjeldahl flasks of 50 to 100 ml and 50 to 100 ml round bottom flasks. Part No. 1002290

Round bottom flask adapter made of aluminium and suitable for Kjeldahl flasks of 300 to 800 ml and round bottom flasks of 250 to 500 ml. Part No. 1002280

Support rod of 12mm Ø and 500 mm high, fits on to the main body of the hotplate. Suitable for the circular Combiplacs 3002301 and 3002302.

Part No.1000725



Electronic contact thermometers with digital readout "Sensoterm". With Pt100 sensor with stainless steel

AISI 304 sheath.

Part No. 1001009

Features. See Regulators and Controllers (pages 293 and 294).











Features. See Regulators and Controllers (pages 293 and 294).





Rectangular precision hotplate "Plactronic"ADJUSTABLE TEMPERATURE FROM AMBIENT +5 UP TO 200 °C STABILITY ±0.5 °C.

ADJUSTABLE TEMPERATURE FROM AMBIENT +5 UP TO 200 °C STABILITY ±0.5 °C WITH ELECTRONIC DIGITAL TEMPERATURE AND TIME CONTROL DISPLAY. SUITABLE FOR THERMAL TESTS (NOT SUITABLE FOR BOILING TESTS). SURFACE PLATE COATED WITH P.T.F.E. FOR PROTECTION AGAINST CHEMICALS.

FEATURES

Digitalelectronic temperature and timer control. and Digital display (3 digits).

Display resolution: 1°C. Temperature probe Pt 100.

DUR aluminium alloy hotplate's coated with P.T.F.E., with heating elements embedded in the plate surface. Main case, heater element and protective grill, made of AISI 304 stainless steel.

SAFETY:

CONFORMS TO EN61010. OVER TEMPERATURE SAFETY CUTOUT FITTED.

Maximum precision





CONTROL PANEL

- 1. Digital display of temperature and time.
- 2. Indicator lamp to show the parameters (temperature or time)
- 3. Push button to select visualised parameter.
- 4. Push button to increase parameter value.
- 5. Push button to decrease parameter value.
- 6. Push button start and stop.
- 7. Indicator alarm if the plate is over temperature.

MODELS

Part No.	Width / Depth (plate) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
6155100	24 24	13 25 42	700	7.5
6156100	20 40	13 22 57	800	8
6157100	25 60	13 27 77	1000	16.5



Rectangular hotplate "Placfin"

ADJUSTABLE HOT PLATE TEMPERATURE FROM AMBIENT +5 °C UP TO 130 °C. STABILITY ±3 °C. SUITABLE FOR THERMAL TESTS (BUT NOT BOILING).

FEATURES

Temperature control by a hydraulic thermostat.

Protected thermometer from 0 °C to 130 °C.

Hotplate surface protected by chemically treated aluminium, with heating elements distributed throughout the hotplate's surface.

Enamelled external case, resistant against corrosive chemicals.

CONTROL PANEL

Hydraulic thermostat control, synchronised with an indicator lamp when in operation.

MODEL

Part No.	Width / Depth	Height / Width / Depth	Power	Weight
	(plate) cm	(exterior) cm	W	Kg
7000497	23 22	11 22 30	400	3.7



ACCESSORY

Slide support adapter made of AISI 304 stainless steel. Capacity: up to 18 slide of 76 x 255 mm (in a horizontal position).
Part No. 7000498





Bank of circular hotplates "Uniplac" ADJUSTABLE TEMPERATURE HOTPLATE UP TO 400 °C.

SUITABLE FOR KJELDAHL DIGESTION ASSAYS AND SOXHLET EXTRACTION.

FEATURES

Independent hotplate control using pulsed energy regulators.

Cast steel hotplates with embedded heating elements distributed throughout the hot-

Adjustable height and angle stand bars for supporting condensers and other glassware etc.

External case made of AISI 304 stainless steel.

CONTROL PANEL

Independent temperature control for each position.

- Temperature control by pulsed energy regulators.
- Indicator lamps showing when the hotplate is on.

MODELS

Part No.	No. of places	Plate Ø cm	Height / Width / Depth (exterior) cm			Power W	Weight Kg
6002314	4	12	23	77	27	2200	15
6002316	6	12	23	105	27	3300	22

ACCESSORIES

Bath of AISI 304 stainless steel of 140 mm Ø x 75mm high, with concentric reduction ring lid, Suitable for water, oil or sand.

Part No. 7002510



Round bottom flask adapter made of aluminium and suitable for Kieldahl flasks of 50 to 100 ml and 50 to 100 ml round bottom flasks.

Part No. 1002290



Round bottom flask adapter made of aluminium and suitable for Kjeldahl flasks of 300 to 800 ml and round bottom flasks of 250 to 500 ml.

Part No. 1002280



ACCESSORIES FOR KJELDAHL DIGESTION

Kjeldahl support Bar for Kjeldahl flasks, made of AISI 304 stainless steel. Part No. 6023141 4 places Part No. 6023161 6 places

Clamp for Kjeldahl condenser tubes made of AISI 304 stainless steel.

Part No. 6000008



Borosilicate glass fume tube. Part No. 1000104 4 place

Part No. 1000106 6 place



Boss head for supporting clamps. Part No. 1000069



"Scrubber" Recirculating pump.

Features (see page 251).

Prevents water wastage.

Models

Unit "Scrubber"

Part No. Height / Width / Depth (exterior) cm Kg 32 31 16 2

It is supplied with 3 Kg. acid vapours neutralizer solution.

Water recirculating vacuum pump

Part No.	Height / (ex	Width terior)		Vacuum level bar	Vacuum pump rate litres/minute	Weight Kg
4001612	44	39	28	0.98	10	10

ACCESSORIES SOXHLET EXTRACTION



Clamp for flasks or condenser tubes. Part No. 7000403



Thermocycler for thermal cycles
Spectrophotometers
Colorimeters
Microtiter strip reader and water quality meter
Microplate reader and washer
pH-meters
Conductivity meter

page 238
pages 239 to 242
pages 243
page 244
page 245
pages 246 to 247
page 248



Thermocycler for thermal cycles "B960"

WITH FOUR ADAPTER MODULES FOR 0.2ML VIALS, 0.5ML VIALS AND WELLS. ADJUSTABLE TEMPERATURES FROM 0°C TO 99 °C.



Advanced technology. Peltier effect.

APPLICATIONS

It is used in molecular biology to amplify DNA via a method based in the polymerase chain reaction process.

BASIC PRINCIPLE

The equipment performs a certain quantity of thermal cycles depending on the method used and repeat them many times for a while, and at the end, the initial DNA fragment Fiber have been replicated thousands times.

For a better process output, changes between temperature levels must be made with minimum time. With thermocycler K96, the cycle temperature can be reached in seconds, even if beginning from remote positions of the last set point. These changes take place by keeping a perfect uniformity between different block points.

The system can also be programmed to produce a linear gradient of temperature widthways the block. This achieves that the highest level of productivity points of the process is optimized and located.

FEATURES


The thermocycler is composed of an inner lid system with heater and with adjustable height in order to perfectly adapt to the samples size. This prevents condensations in the upper side of the samples.

The equipment is based on a heat bomb controlled by continuous electric current and composed by several thermo electrical modules of Peltier effect, a low thermal resistance radiator and a fan convention system.


As it is integrated in the same block, this system allows increasing the process output and also fast transferring and extracting the block's temperature, going from the highest temperature level to the lowest in a minimum time.

The powerful control microprocessor allows monitoring the point where the process is at any moment, and showing it on the screen with real-time graphic images.

There is a useful elaborated software for processes organization, accessible by means of the keyboard and a high resolution LCD screen.



384 wells block, with red, yellow, green and black identification in diferent positions.



Graphic detail that shows temperature perfect uniformity in the different block wells, distinguished by different colours: red, yellow, green, black. It also shows how uniformity is maintained during all the thermal cycle periods, previously programmed in the display.



INTERCHANGEABLE BLOCKS

Each block incorporates a connector that identifies it and allows the thermocycler to recognize it. There's also an extractor handle which makes easy the block support use.



Thermal cycler door open with a heater system on the back side of the door to avoid condensations.

TECHNICAL SPECIFICATIONS

Temperature range: 0°C to 99°C.
Durability: 99min. 99sec.
Warming-up speed: 4°C/sec.
Cooling down speed: 4°C/sec.
Uniformity: At 95°C ±0.4°C.

From 20°C to 75°C ±0.2°C.

Precision: ±0.2°C.

Programmed gradient: From 2°C to 30°C according to the program.

Temperature range gradient: 30-99 °C.

Heater lid: From 70 to 115°C according

to the program.

Max. cycles number: 299.

Stored programs: Up to 1000.

Graphic display of 14.5cm, 320x240 pixels.

USB 2.0, LAN and RS232 output.

CONTROL PANEL

START switch.
Interactive graphic display.
Numeric and functions keyboard.
Height levelling control of the inner lid.

MODEL

Part No.	Height / Width / Depth (exterior) cm	Power supply	Power W	Weight Kg
5109000	25 27 38	220V/50-60Hz	780	7,8

ACCESSORIES

Interchangeable modules:

A. 96 vials of 0,2 ml. Code: **5109001**C. 96 vials of 0,2 ml + 77 vials of 0,5 ml. Code: **5109003**D. 384 wells Code: **5109004**



Visible range spectrophotometers "V-1100" and "VR-2000"

"V-1100" MODEL WITH MANUAL WAVELENGTH SETTINGS AND AUTOMATIC BLANK.

"VR-2000" MODEL WITH AUTOMATIC WAVELENGTH SETTINGS AND BLANK.





They are widely used in colleges and enterprises for general quantitatives analysis and experiments based in absorbance measurements.

COMMON FEATURES

High quality silicon photometric diode detector and 1200 lines/mm diffraction grating ensure the high quality accuracy and precision.

Digital display for an easy readout.

Automatic zero and blank (easy to use). Easy switching of transmittance, absorbance and concentration modes, just by pressing one key.

Large sample compartment, which can accommodate 5 to 100mm path length cuvettes with optional holders.

Its pre-aligned design makes it possible to change the halogen lamps by the user himself.

Optional:

Optional software based on Windows® which can expand the applications to standardcCurve and kinetics.

V-1100 MODEL

Easy of use.

Ergonomic and solid design for a continuous suitable use, for students and workers.

VR-2000 MODEL

Large LCD screen (128x64bits).

It can display a total of 50 groups of data (3 groups per screen).

It can display standard curve and kinetics curve graphic.

The system can also save the test results.

A total of 50 data groups and 10 standard curves can be saved in the RAM memory.

At most 9 standard samples can be used to establish a standard curve. The curve and the curve equation will be displayed simultaneously on the screen. The unknown concentration solutions, can be measured by the curve. If one knows the coefficient \mathbf{k} and \mathbf{b} of the formula: $\mathbf{C} = \mathbf{k} \mathbf{A} + \mathbf{b}$, one can input the value directly.


Data is stored in the memory in case of power cut.

ACCSESORIES

(see page 241).



MODEL	V-1100	VR-2000						
Part No.	4120025	4120026						
Wavelength range	325-1000 nm							
Spectral bandwidth	5 nm 4 nm							
Optical system	Single beam, diffraction grating 1200 lines/mm							
Wavelength adjustment	Manual	Automatic						
Wavelength accuracy	±2 nm	±1 nm						
Wavelength repeatability	1 nm	0,5 nm						
Photometric accuracy	±0,3% T	±0,5% T						
Photometric repeatability	±0,3% T							
Photometric range	-0,3-3 A, 0-200% T. 0-9999 Concentration							
Stray light	0,5% T	0,3% T						
Stability ±	0,004 A/h @ 500 nm	± 0,002 A/h @ 500 nm						
LCD Display	3,5 Digits	128x64 pixels						
Detector	Silico	on						
Sample compartment 10) mm path length standa	rd cuvette (100mm optional)						
Light source	Tungste	n lamp						
Output	USB, Parallel	port (printer)						
Mains supply	220 V / 50Hz AC or 110 V / 60 Hz AC							
External dimensions (HxWxD)	160 x 480 x 360 mm	180 x 470 x 370 mm						
Weight (Kg)	8	12						
SPARE Halogen lamp (visible)	V-1100 4312004	VR-2000 4312007						



Graphic display visualized detail of a regression line for concentration calculation. (VR-2000)



Ultraviolet and visible range spectrophotometers "UV-2005" and "UV-3100"

AUTOMATIC WALVELENGTH POSITIONING AND BLANK SETTING



APPLICATIONS

Suitable for pharmaceutical and biochemical laboratories, kinetics, quantitative analysis, wavelength scanning, multiple components and DNA/protein analysis.

COMMON FEATURES

Spectrophotometers UV-2005 and UV-3100 have been developed for accurate tests; Its stray light is only of 0.05% T. They are flexible, easy-to-use and maximize value. Value is provided from every day perfomance. The independent local software provides the following functions: Basic Mode, Quantitative Test, Kinetics and System Utilities.

Data is stored in the memory in case of power cut.

The Special application software provides a complete control of the spectrophotometer from a PC, through the built-in USB port.

Automatic wavelength settings.

Halogen and deuterium lamps can be switch on/off individually to extend lamp lifetime. Its pre-aligned design makes it possible to change lamps for the user himself.

Large sample compartment for 5-100mm path length cuvettes, with optional holders. A wide range of optional accessories can also be selected.

UV-2005 MODEL

Large LCD screen (128x64bits).

It can display a total of 200 groups of data (5 groups per screen).

It can display standard curve and kinetics curve graphic.

The system can also save the test results.

A total of 200 data groups and 200 standard curves can be saved in the RAM memory.

UV-3100 MODEL

Absorbance, transmittance and concentration measurements.

It establishes or uses stored calibration equations to measure solutions of unknown concentration.

Spectrum scan of sample at any selected wavelength range with choice of scanning speed and wavelength interval.

Measurement of absorbance changing versus time with reaction rate calculation function.

Measurement at multiple wavelengths to analyse and determine the composition of the mixtures.

Calculation of concentration and DNA purity.


Note: ratio at other wavelengths can be measured.

ACCSESORIES

(see page 241).



MODELS	UV-2005	UV-3100				
Part No.	4120020	4120021				
Wavelength range	190-	1100 nm				
Spectral bandwidth	2	2 nm				
Optical system	Single beam, diffraction	n grating 1200 lines/nm				
Wavelength accuracy	±0	,5 nm				
Wavelength repeatability	0,	3 nm				
Photometric accuracy	±0,	3% T				
Photometric repeatability	±0,	2% T				
Photometric range	-0,3-3 A, 0-200% T. 0-9999 Concentration					
Stray light	0,05% T @220 nm, 340 nm					
Stability	± 0,002 A	/h @ 500 nm				
LCD Display	Graphic(128x64)	Graphic(320x240)				
Scanning velocity	- High,Med,	Low.Max 3000 nm/min.				
Uniformity	-	± 0,002 A (200-1000nm)				
Sample compartment	10mm path length standard c	uvette (100mm optional)				
Light source	Halogen and deuteriu	m lamps (pre-aligned)				
Output	USB, parallel port (printer)					
Power requirements	220 V / 50Hz AC or 110 V / 60 Hz AC					
External dimensions(HxWxD	180 x 470 x 370 mm	160 x 480 x 360 mm				
Weight (Kg)	14	16				



Graphic display visualized detail of holmium crystal scanning, from a UV3100 spectrophotometer.

 SPARE PARTS
 UV-2005
 UV-3100

 Halogen lamp (visible)
 4312007

 Deuterium lamp (UV)
 4312006

ACCESSORIES

Thermostatted flow Peltier cuvette holder.

Temperature range: from 15 °C. to 40 °C. in steps of 0.1 °C (ambient temperature \pm 22°C.)

Precision: ± 0,2 °C.

Sample suction volume: 1, 2, 3, 4, 5, 7.5, 10, 12.5, y 15 ml/min.

LCD Display: 128x64 pixels.

Model

Part No.	Height / Width / Depth	Height / Width / Depth	Power	Weight
	(controller) cm	(cuvette holder) cm	W	Kg
4120018	11,5 27 17	10 12 7	30	3,5

Cuvette holder

For long cuvettes, pathlength between 1 and 50 nm. Part No. 4120030 For path leength100 mm. cuvettes. Part No. 4120031

Test tubes holder

For tubes from 10 to 20 mm. Ø Part No. 4120032

Ink printer (not thermal). Paper 2 1/4" (56 mm) wide roll. Includes power transformer and interface cables. (5V-3A)

Model

Part No.	Height A	Height / Width / Depth cm		
4120117	4	16	10	1













FEATURES

Range of polystyrene, glass and quartz cuvettes. Standard size pathlength 10mm x 45mm high. Special cuvette pathlength 40mm x 45 mm high (fig 6). Special cuvette pathlength 4mm x 45 mm high (fig 1).



MODELS

MODELO					
Standard cuvette	s				
Part No.	Figure	Material	Optical path length	Description	Presentation
5100020	9	Quartz	10 mm	Standard square, "macro"	Pack of 2 units
5100021	9	Glass	10 mm	Standard square, "macro"	Pack of 2 units
5100022	7	Polystyrene	10 mm	Standard square, "macro", disposable	Box of 100 units
5100025	8	Glass	10 mm	Flow through cell, square "micro"	Pack of 2 units

Special cuvettes					
Part No.	Figure	Material	Optical path length	Description	Presentation
5100014	6	Glass	40 mm	Rectangular	Pack of 1 unit
5100015	5	Glass	10 mm	Flow through cell. Window 5 x 10 mm	Pack of 1 unit
5100016	4	Glass	10 mm	Standard square, "micro"	Pack of 1 unit
5100017	3	Glass	10 mm	Standard square, "macro" with round lid	Pack of 1 unit
5100018	2	Quartz	10 mm	Semi-micro	Pack of 2 units
5100019	2	Glass	10 mm	Semi-micro	Pack of 2 units
5100023	7	Polystyrene	10 mm	Square Semi-Micro, disposable	Box of 100 units
5100024	1	Glass	5 mm	Square Semi-Micro	Pack of 1 unit







NECTA: Double Beam Spectrophotometer "UV-2300 II"

BAND WIDTH 1.5 nm. TOTAL STABILITY DOUBLE BEAM OPTICS. AUTOMATIC WAVELENGTH CORRECTION CALIBRATION. USB PORT FOR DATA STORAGE. WIDE RANGE OF ACCESSORIES.





APPLICATIONS

Research, chemistry, biotechnology, general spectroscopy analysis applications, environmental applications.

FEATURES

Monochromator high resolution optics that eliminate any optical aberrations, monochromator "Seya-Namioka", manufacturer of exclusive beam technologies and diffraction gratings in Japan.

Bandwidth 1.5 nm in accordance with European pharmacopoeia recommendations. (the relation between the maximum and minimum absorbence in Toluene and Hexane at 0.02% (V/V) should be more than 1.5T).

Several modes of operation including spectral scanning, time base scanning, multi wavelength determinations, peak and trough detection, etc.

Fast spectra scan displayed on a screen covering the whole spectral range: quick scan 3600 nm/minute, range 190 to 1100 nm

Validation function for GLP/GMP:

This function maintains and assures the optimum working parameters of the instrument. Parameter such as the wavelength precision and noise are monitored.

Memory facilities of analytical results:

The analytical parameters and results can be stored in the "flash" memory, connected by the USB. The stored information can be stored as text and can be transferred to a computer for reporting using MS WORD/EXCELL.

A DNA/RNA function is fitted to quantify the ratio at 260/280 nm.

Controllable from a computer with application specific optional software: "UV-Analyst Spectrum", (see accessories).



SPECIFICATIONS

Optical system: Double beam optics. Wavelength range 190 nm to 1100 nm.

Band pass: 1.5 nm

ABS range: -2.000 to 3.000 A or 0 to 300% T.

Stray light: less than 0.05% (220 nm Nal, 340 nm NaNO₂). Scan speed: 10, 100, 200, 400, 800, 1200, 2400 and 3600 nm/minute.

Wavelength accuracy: ±0.3 nm.

Photometric accuracy: ± 0.002 A from 0 to 0.5 A.

± 0.004 A from 0.5 to 1.0 A.

Baseline stability: 0.003A/hr (500 nm after 2 hours of use.).

Noise level: 0.003 A (at 500 nm).

Light source: Deuterium D2 and Halogen lamps.

Built in screen: LCD back light of 190 x 138 mm.

Connections: RS232 and parallel port.





Constant temperature 6 place cell changer, 10mm cuvettes, ideal for Kinetics.

Part No. 5110029



USB port. "Flash" memory not included.

SPARES

Tungsten halogen lamp. Part No. **5110021**

Deuterium lamp (UV).Part No. 5110022

AU	U	Eg	90	K	13
_					

Part No. **5110023** Flow cell 10 mm path length.

Part No. 5110024 Constant temperature single cell holder for 10mm flow cell.

Part No. 5110025 Micro cell holder for 10mm (50 µl) cell.
Part No. 5110026 Long path length cell holder for 100mm cells.

Part No. 5110027 5 Position cell changer for 10mm cuvettes.
Part No. 5110028 Constant temperature cell holder for 10mm

Part No. 5110028 Constant temperature cell holder for 10mm cuvettes
Part No. 5110029 6 place constant temperature cell changer for 10mm cuvettes.
Part No. 5110033 Software UV-Analyst Spectrum. Simple and logical, enhanced

software for multiple applications and results and manipulation of

data, calculation such as DNA/RNA ratios.



Semi Automatic clinical analyser "Photometer S-2000"

LOW PROFILE FUNCTIONAL CLINICAL ANALYSER.

APPLICATIONS

Clinical laboratories, Clinical Biochemistry, Haematology, E.I.A. Electrolytes and lones.

FEATURES

47 most commonly used Pre-programmed clinical methods: cholesterol, glucose, blood iron, LHD, AST, UREA etc.

Storage for up to 134 methods and up to 1000 stored results.

6 analytical fitting routines: final point, bio-chromatic, regression fit, multi-standard curve, best fit and kinetics.

Graphic display LCD (640 x 240 pixels, 256 colours).

Lamp life save function.

Graphic display.

Quality control of analysis function.

Built in printer.

Control by "roller ball mouse" or "external mouse."

Temperature conditions: 15 °C to 32 °C: Maximum humidity 85%

RS232 out put interface.

SPECIFICATIONS

Solid state photodiode detector.

Bichromatic optical system with 7 filters, 340, 405, 500, 546, 578, 620 and 670 nm.

Band pass: 10 nm. Flow cell: 30 ul.

Measuring range: from 0.000 to 2.500 ABS.

Resolution: display 0.001 ABS, calculations to 0.0001ABS.

Thermostatted Peltier cuvette holder: 25, 30 and 37 °C Precision ±0.1 °C.

Tungsten halogen lamp of 6V 10 W.

Comes complete with protective cover, spare lamp and two rolls of printer paper.



MODEL

Part No.	Height	Height / Width / Depth cm		Power rating	Power W	Weight Kg
4120010	14	45	33	110-220V/50-60Hz	150	10

ACCESSORY

Metallic thermostat dry block "BIO-BATH"

Small size with a fixed temperature of 37 $^{\rm o}{\rm C}$ with check thermometer, 8 tube positions for 75x13 VAC tubes.

Note: Accelerates the extraction of plasma for faster results and waiting time for the patient.

MODEL

Part No.	Height / Width / Depht			Stability	Power	Weight
	(exterior) cm			°C	W	Kg
7001561	8.8	9.6	7	0.5	8	0.5

SPARE. Halogen lamp 6 V 10 W. Part No. 4312015



Digital colorimeter "Clormic"

MICROPROCESSOR CONTROLLED.
AUTOMATIC ZERO ABSORBANCE AND 100% TRANSMITTANCE.
ALPHANUMERIC 20 CHARACTER 2 LINE L.C.D. DISPLAY.

APPLICATIONS

Chemistry Laboratory, quality control, environmental control.

FEATURES

Wavelength range: 400 to 800 nm, by using special filters with a 30 nm band pass.

12 position Filter disc wheel selection.

Standard filters: 420, 440, 490, 520, 550, 580, 620 and 680 nm.

Expanded Absorption range: -0.3 to 3.5 O.D. in real time.

Transmission: 0 to 100 T %.
Photometric accuracy: >1 %.
Photometric precision: ±1 %.
Photometric stability: 0.004 A/hr.
Light source: Long life tungsten lamp.

Detector: Solid state.

Sample chamber: 10 mm cuvettes or test tubes of 12 and 16 mm \emptyset .

Minimum volume: 1 ml.

Display: Alphanumeric LCD display of 2 lines of 20 characters.

Calculation functions: Transmission T %. Absorbence, Concentration by factor or standard concentrations.

Calibration: Self adjusting by software. **RS-232 interface.**

CONTROL PANEL

ON/OFF switch.

Interactive LCD display. Numeric and function keypad.



MODEL

Part No.	Built	Height / Width / Depth	Power	Weight
	in printer	cm	W	Kg
4120009	NO	11 18 28	10	4.5

SPARE

Lamp of 6 V / 6 mm. Part No. **4512009**



Microtiter Strip Reader "Reader M-2000"

STORAGE MEMORY FOR UP TO 59 METHODES INCLUDING: FACTORS, STANDARDS OR CALIBRATION CURVES WITH UP TO 10 STANDARDS. ADDITIONAL CALCULATION FEATURES.

APPLICATIONS

Clinical diagnostics, E.I.A., Veterinary, quality control.

FEATURES

Spectral range: 400 to 700 nm, using interference filters of 10 nm band pass.

Supplied complete with a filter of 450 nm (other wave-

lengths available consult for options.)
Measuring range: -0.2 to 3.5 0.D.
Photometric accuracy: >1%.
Photometric precision: ±1%.
Light source: Long life LED.

Detector: solid state. 59 method memories. Curve point editing feature.

Software: ABS. Concentration and polynomial curve fitting with up to 10 standards.

RS-232 output to connect a printer or to a PC.
Interface RS-232.



MODEL

Part No.	Height / Width / Depth		Power	Weight	
	cm		W	Kg	
4120013	11	18	28	7	2

ACCESSORY

Ink printer (not thermal). Paper 2 1/4" (56 mm) wide roll. Interface RS232. Includes power transformer and interface cables. (5V-3A)

Model

Part No.	Height / Width / Depth cm	Weight Kg
4120113	4 16 10	1



SPARE

Lamp of 4 mm Ø 5 V. Part No. **4313030**



Multi-parameter water quality meter "MP-2006"

LARGE TOUCH SCREEN. TO MEASURE PH, ORP, CONDUCTIVITY, TDS, SALINITY, DISSOLVED OXYGEN, IONS & TEMPERATURE

FEATURES

5.7 inches touch screen.

Capacity to store up to 2200 sets of data. GLP method. Possibility of storing, printing and deleting measuring data. USB and RS-232 ports for PC connection.

Measuring modes: continuous, timing and balance measuring modes.



	Range	Accuracy	Resolution
pH:	-2 to +19,999	±0,002 FS	0,001
ORP:	±1999,99 mV	±0,03 FS	0,01 mV
Conductivity:	0,001 µS to +199,9mS	±0,5 FS	0,001µS
TDS:	0,000 mg/L to 19.99g/L	±0,5 FS	0,001 mg/L
Salinity:	(0,0 to 8,00)%	±0,1%	0,1%
Dissolved oxygen:	0,00 to 19.99 mg/L	±0,5 mg/L	0,01 mg/L
lons:	0 to 19990(ug/L, mol/L)	±0,5%	±1
Ta:	-0,5 to +135	±0,5 FS	0,1



MODEL					
Code	Height	/ Width cm	/ Depth	Power	Weight Kg
4100700	20	20	7	OOON/EO COU-	4

ACCESSORIES

Standard electrodes of pH, conductivity and temperature. See page 247.

Special electrodes:

oheeiai eid	outouto.								
Code	Electrode type	Working temperature	Body material	Measure mm	Code	Electrode type	Working temperature	Body material	Measure mm
4120183	Platinum	0-50⁰	Glass	Ø 12 x 120	4120189	Calcium	5-60⁰	Glass	Ø 12 x 120
4120184	Platinum-2	9 0-50°	Glass	Ø 12 x 120	4120190	Nitrate	5-60⁰	Glass	Ø 12 x 120
4120185	Antimony	0-50⁰	ABS	Ø 12 x 120	4120191	Copper	5-60⁰	Glass	Ø 12 x 120
4120186	Tungsten	0-50⁰	Glass	Ø 12 x 120	4120192	Chloride	5-60⁰	Glass	Ø 12 x 120
4120187	Silver	0-50⁰	ABS	Ø 12 x 120	4120193	Oxygen	5-40⁰	ABS	Ø 12 x 120
					4120194	Reference	5-55⁰	Glass	Ø 12 x 120



Microplate reader 2100-C

TOTALLY AUTOMATIC TOUCH SCREEN CONTROL.



APPLICATIONS

Clinical diagnostics. Quality control in food analysis, Haematology.

FEATURES

Bi-chromatic optical system with 4 wavelength selection:

 $405,\,450,\,492,\,630$ nm, (other filters from 405 to 700 nm available). Memory for up to 500 programs, 10000 sample results and 1000 patients

Patient data software, laboratory personnel and analysis reports.

100 typical analysis methods pre programmed.

Capable of performing up to 12 different tests per plate.

 ${\it Calculations: ABS, Cut-Off, \ Curve \ fitting: regression, \ linear, \ logarithment \ logar$

mic, exponential and potential.

Measuring absorbance range: 0-2.5A

Displayed absorbance range: 0-3.5A

Accuracy:±1.0% or ±0.007 A.

Precision: ±0.5% or ±0.005 A.

Linearity: r > 0.995.

Resolution:0.001 ABS (display); 0.0001 ABS (calculated).

Measuring time: Continuous mode < 5 s, step by step mode < 15 s. Output RS-232 for download to a computer or a printer.

Parallel printer interface.

Programmable plate vibrating for mixing with control of time and speed.

CONTROL PANEL

LCD display (320 x 240 pixels) of 5.7 inches.

Operating system WINDOWS CE.

Operation by touch screen or optional "Mouse" control (mouse not in-

cluded).



MODEL	Part No.		Width / Depth erior) cm	Voltage	Weight Kg
2100-C	5109999	19	45 33	AC 110-250 V / 50-60 Hz	10

ACCESSORIES

Software for data transfer to computer. Part No. 5109997



COMECTA Microplate washer 2600-C

COMPLETELY AUTOMATIC. EASY TO USE. 50 WASH PROGRAMS.



APPLICATIONS

Clinical laboratories.

Blood banks.

Quality control food laboratories.

FEATURES

Wash head of 8 x 1. (12 x 1 optional).

Plate wash by 12 x 8 strips or 8 wells/strip.

Plate clean for flat, round and V bottom plates.

Automatic wash and rinse processes.

Reservoirs: 2 litre, wash solution.

2 litre, waste, with liquid level sensor.

Memory for up to 50 wash programs.

Software controll of head's depth adjust and horitzontal position.

Wash cycle times from 1 sec. to 2 hrs. Wash dosing range: from 50 to 2000ul.

Dosing resolution: 5% CV (based on 350ul of distilled water)

Residual volume (after aspiration):

< 1 µl (Cells, bottom U and V well plates.).

< 5 μ l (Cells, bottom of flat well plate.)

CONTROL PANEL

LCD graphic display of 90 x 50mm Membrane key pad of 9 keys. Emergency stop button.



MODEL	Part No.	Height / Width / Depth (exterior) cm	Voltage	Weight Kg
2600-C	5110000	14 45 33	AC 110-250 V / 50-60 Hz	12



Digital pH-meter "pH-2003" PORTABLE METER.

FEATURES

Manual temperature compensation from 0 to 60 $^{\circ}$ C with sample temperature key pad entry (does not read temperature).

Automatic buffer recognition of two of those three buffers: pH4, pH7 and pH9.

Electrode connector: BNC type. Includes: carry case and pH electrode.

SPECIFICATIONS

	Range	Precision	Resolution
рН:	0 to 14	±0.03	0.01
mV:	±1400	±2	1

Impedance: 5 x 10¹¹ Ohm.

MODEL

Part No.	Height / Width / Depth cm	Power source	Weight Kg
4120300	17 7.5 3	2 x 1.5V batteries	0.5



Digital pH meter "pH-2004"

PORTABLE METER WITH 1P65 PROTECTION AGAINST WATER AND DUST.

FEATURES

Automatic temperature compensation and read out from 0 to 60 °C with temperature sensor Rt-2252. Results memory for up to 250 readings.

Protection against water and dust IP65 (Impermeable to water splashes from any direction.) Automatic buffer recognition and calibration for 1 or 2 buffers, recognition of up to 5 selectable buffers.

RS-232 output to connect a printer or a computer. Electrode connector: BNC type. Includes: carry case, pH electrode and temperature sensor Rt-2252.

Direct connection to printer Part No. 4120113 (see page 247).

Universal arm support not included. Part No. 1001552 (see page 248).

SPECIFICATIONS

	Range	Precision	Resolution
pH:	0 a 14	±0.02	0.01
mV:	±1800	±2	1
Ta:	0 to 60	±0.5	0.1
Impedance:	1 x 10 ¹² Ohm.		

MODEL

Part No.	Height / Wid	•	Power source	Weight Kg
4120400	20 10	4.5	4 x 1.5V AA Batteries	0.5



Digital pH meter "pH-2005"

BENCH TOP MODEL WITH LARGE SCREEN.

FEATURES

Automatic or manual temperature compensation from 0 to 99.9 $^{\circ}$ C Automatic buffer recognition and calibration of 2 buffers, recognition of up to 3 buffers: pH 4, pH 7 and pH9.21.

RS-232 output to connect a printer or a computer.

Electrode connection BNC.

Temperature sensor connection Rt-2252.

Includes electrode's support stand.

SPECIFICATIONS

	Range	Precision	Resolution
pH:	0 to 14	±0.02	0.01
mV:	±1999	±1	1
Ta (°C):	0 to 99.9	±0.5	0.1
Impedance:	1 x 10 ¹² Ohm.		



See "Accessories for pH meters" for selection of electrodes (Page 242).



ACCESSORY

pH probe and additional temperature sensor, gel filled. Temperature sensor Rt 2252. With IP65 protective connectors.

Part No. 4120401



ACCESSORY

Portable magnetic stirrer, battery operated and submersible. Part No. 7001731 for "pH-2003" and "pH-2004" (see page 27).







MODEL

Part No.	Height / Width / Depth cm	Weight Kg
4120500	30 20 7.2	1.5

See accessories for electrode, temperature sensor and buffer solutions. Not compatible with printer Part No. 4120113 (see page 247).



Digital pH meter "pH-2006"

BENCH TOP, LARGE GRAPHIC TOUCH SCREEN CONTROL.

FEATURES

Ergonomic touch screen control functions. Includes support arm.

Automatic temperature compensation from -5 to 105 °C.

Up to 5 buffer calibration with 10 buffer recognition.

Stores up to 200 results.

RS-232 output to connect a printer or a computer.

Electrode connector: BNC type. Rt-2252 Temperature port. Selectable resolution range.

Direct interface for printer Part No. 4120113 (see accessories).

SPECIFICATIONS

	Range	Precision	Resolution
pH:	-2 to 18	±0.01	0.001
mV:	±1999.9	±0.1	0.01
T ^a (°C):	-5 to 105.0	±0.3	0.1

Impedance: 1 x 10¹² Ohm.

See accessories for electrodes, temperature probes, buffers and printer.



	_	_	÷	÷
IVI	Ш	Ш	H	ı
	_	_	-	•

Part No.	Height / Width / Depth cm	Weight Kg
4120600	29 20 7	1

ACCESSORIES for pH meters.





Universal arm support.

With holder for up to 4 sensors or electrodes. Very stable heavy base.

Ideal for working at various heights.

Part No. 1001552

General laboratory use.



μ	П	Electiones
F	ig.	Part No.

Fig.	Part No.	Range pH	Sensor T ^a	Temperature Range	Electrolyte	Dimensions mm	Material Body
1	4120102	0 - 14	No	0 - 80°	Liquid	Ø 12 x 120	Glass
2	4120125	0 - 14	Yes	0 - 60°	Gel	Ø 12 x 120	Ероху
3	4120101	0 - 14	No	0 - 60°	Gel	Ø 12 x 120	Ероху
4	4120114	0 - 10	No	0 - 60°	Gel	Ø 12 x 120	Ероху
5	4120104	0 - 14	No	0 - 60°	Gel	Ø 12 x 120	Glass
6	4120103	0 - 14	No	0 - 80°	Gel	Ø 8 x 60	Glass
7	4120126	0 - 14	Yes	0 - 80°	Liquid	Ø 12 x 120	Glass

у	Protective
у	Protective

Application

cover. General laboratory and external use. cover. General laboratory and outside use.

Surface analysis, Skin, meat, moist semisolids, culture media. Pointed sensor. Suitable for food, semi-solids, fruit meat etc.

For measuring in tubes. Liquid solutions in general.

Fast response, liquids in general.

Temperature sensors

Fig.	Part No.	Range	Туре	Dimensions mm	Material
8	4120121	0-100	Rt-2252	Ø 5 x 120	Stainless steel
9	4120120	0-100	Rt-2252	Ø 12 x 120	Glass

Ion selective electrodes

IUII UUIUUL	ivo cicotioaco				
Part No.	lon	Range	Limits	Range °C	Range pH
4120168	Chloride Cl ⁻	1-3 x 10 ⁻⁶	35.500 - 1	5-50	1-12
4120174	Nitrate NO⁻₃	1-7 x 10 ⁻⁶	62.000 - 0.4	0-50	2-11
4120177	Potassium K⁺	1-10-6	39.000 - 0.04	0-50	1-9

pH Buffer solutions.

4120107 Flask 250 ml Buffer ±0.02 pH 4.00 at 20 °C; 4.01 at 25 °C Flask 250 ml Buffer ±0.02 pH 7.02 at 20 °C; 7.001 at 25 °C 4120108 Flask 250 ml Buffer ±0.02 pH 9.23 at 20 °C; 9.18 at 25 °C 4120127 Frasco 250 ml KCI 3M solution (Electrolyte for pH electrodes) 4120109

MAGNETIC STIRRER "AGIMICRO" Features

Maximum stir volume: 2 litres. Adjustable stir speed from 40 to 1400 r.p.m.

Excellent chemical resistance.

Location hole at the back for a retort rod. (Accessory)

Ink printer (Not thermal).

interface cables. (5V-3A)

For use with pH and conductivity meters . Paper 2 1/4" (56 mm) wide roll. Interface RS232. Includes power transformer and



Control panel

Mains switch with luminous on indicator. Analogue control of speed in r.p.m.

Model

Part No.	Medidas cm	Power source	Weight Kg
7001638	Ø 12 x 5	230 V 1 W	0.5

Supplied complete with a 8 Ø x 32 mm P.T.F.E. coated stir bar.



ACCESSORY Support bar and electrode clamp. Part No. 7001639

Model			
Part No.	Height / Width / Depth cm	Power source	Weight Kg
4120113	4 16 10	DC 5 V, 3 A	1



Conductivity meter "CD-2004"

PORTABLE EQUIPMENT WITH 1P65 PROTECTION AGAINST WATER AND DUST.

FEATURES

Measures conductivity, TDS and salinity.

Automatic range switching.

Automatic temperature compensation.

Calibration using standards for conductivity, TDS and salinity.

RS-232 output to connect a printer or a computer.

Direct printer connection, see part no. 4120113.

Temperature coefficient programable.

Conductivity cell constant adjustment.

Memory for 250 results.

Ambient protection IP65 (protection against water splashes from any direction).

Includes Carry case, conductivity cell K= 1 and temperature sensor Rt-2252.

Conductivity cell with epoxy body. Temperature sensor Rt-2252.

Conductivity cell and temperature sensor.

IP65 protection connectors. Part No. 4120411

SPECIFICATIONS

Reading ranges: (automatic)
Conductivity TDS

0.00 to 19.99 µS/cm 0.00 to 10.00 mg/L

Salinity range: 0.00 to 8.00 Temperature range: 0.0 a 40.0 °C

Precision:

Conductivity: 1.5% base scale

Salinity: 0.20% Temperature: 0.3 °C



MODEL

HIODE	•		
Part No.	Width / Height / Depth	Power	Weight
	cm	source	Kg



ACCESSORY

Conductivity meter "CD-2005"

BENCH TOP METER WITH LARGE SCREEN.



Measures conductivity and TDS.

Automatic or manual setting temperature compensation with fixed coefficient of 2%.

Calibration with standards of conductivity and TDS. RS-232 output to connect a printer or a computer.

Cell constant adjustment.

Temperature sensor connection Rt-2252. Temperature measuring range: 0 to 60 °C.

SPECIFICATIONS

Conductivity ranges:

(Manual range change)

0 to 20.00 $\mu\text{S/cm}$

0 to 200.0 μ S/cm 0 to 2000 μ S/cm

0 to 10.00 mS/cm

Precision: 1.5% base scale + 1 digit. TDS ranges: 0 to 1000 mg/L.



Part No.	Width / Height / Depth cm	Weight Kg
4120510	30 20 7.2	0.5

Accessories to complement the meter, conductivity cells, temperature sensors and standards (see accessories). Not compatible with the printer 4120113.

Includes sensor stand.



Conductivity cell

Glass body. Working temperature range 0-100 °C

Fig.	Part No.	K	Sensor I" ("C)	Dimensions mm	Applications
1	4120220	1	Yes	Ø 12 x 120	General laboratory use up to 2000 µs
3	41202220.	1	Yes	Ø 12 x 120	Pure water up to 20 µS
4	4120223 1	0	Yes	Ø 12 x 120	High conductivity up to 200 ms

Temperature sensors

Fig.	Part No.	Range	Туре	Dimensions mm	Material
5	4120121	0-100	Rt-2252	Ø 5 x 120	Stainless steel
6	4120129	0-100	Rt-2252	Ø 12 x 120	Glass

Conductivity standards.

4120160 Flask 250 ml Standard ± 0.05 of 1278 mS/cm at 20 °C; 1413 mS/cm at 25 °C **4120161** Flask 250 ml Standard ± 0.025 of 4.915 mS/cm at 20 °C; 5.446 mS/cm at 25 °C **4120162** Flask 250 ml Standard ± 0.06 of 11.67 mS/cm at 20 °C; 12.88 mS/cm at 25 °C

MAGNETIC STIRRER "AGIMICRO" Features

Maximum stir volume: 2 litres.

Adjustable stir speed from 40 to 1400 r.p.m.

Excellent chemical resistance.

Location hole at the back for a retort rod. (Accessory)

Control panel

Mains switch with luminous on indicator. Analogue control of speed in r.p.m.

Model

Part No.	Medidas cm	Power source	Weight Kg
7001638	Ø 12 x 5	230 V 1 W	0.5

Supplied complete with a 8 \emptyset x 32 mm P.T.F.E. coated stir bar.



ACCESSORY Support bar and electrode clamp. Part No. 7001639



Universal arm support.

With holder for up to 4 sensors or electrodes. Very stable heavy base. Ideal for working at varying heights.

Part No. 1001552



-	251 to 252
paye	203
page	254
page	255 to 256
. 3	
nogo	257
paye	237
page	258
page	259 to 260
-	
paye	261 to 262
page	263 to 265
	page page page page page



NUTRITIONAL AND WATER ANALYSIS APPARATUS

EXTRACTION

Organic nitrogen determination by Kjeldahl "BLOC-DIGEST" Steam distillation of protein. PRO-NITRO "M"

Steam distillation of protein semiautomatic, PRO-NITRO "S"

Kjeldahl distillation automatic PRO-NITRO "A"

Extractor for the determination of cellulose and fibre, "DOSI-FIBER"

Cold extraction unit, "EF-1425"

Solvent extractor unit for the determination of residues, oils and fats

in nutritional and other materials DET-GRAS "N"

Sample Hydrolysis unit, "HI-1427"

Wine (Alcohol), distillation unit "DE 1626"

WATER ANALYSIS

Chemical Oxygen Demand in residual water "C.O.D."

Constant temperature refrigerated incubator cabinets, B.O.D. "MEDILOW S. M. I. 16"

Flocculator for water analysis, laboratory "FLOCUMATIC"

Flocculator for water analysis, portable "JARTEST"

De-mineraliser "LAB-ION"

Water Distillation unit, "AQUASEL", "L-3" and "AC-L8"

RECOMMENDED METHODS AND EQUIPMENT: PRO-NITRO M, S and A, DOSI-FIBER, EF-1425, HI-1427, DET-GRAS N, BLOC-DIGEST, DE-1626 and C.O.D.*

ANALYSIS OF CEREAL AND DERIVATIVES	Reference	Pro-Nitro M,S and A	Dosi-Fiber	EF-1425	Det-Gras N	Bloc-Digest
Determination index for cellulose	Method Wladesco		YES	YES		
Insoluble fibre in food	Method Van Soest		YES	YES		
Crude Fibre	Method Weende & Wijkströn		YES	YES		
Proteins	Method Kjeldahl	YES				YES
Soxhlet extraction for fat identification	Soxhlet Extraction				YES	
Crude Fat	Soxhlet Extraction				YES	
Arsenic	Determination A. A.					YES
Mercury	Determination A. A.					YES
ANALYSIS OF MILK AND DERIVATIVES	Reference	Pro-Nitro M,S and A	Dosi-Fiber	EF-1425	Det-Gras N	Bloc-Digest
Crude Fat	Soxhlet Extraction				YES	
Crude Protein	Method Kjeldahl	YES				YES
Casein N	Nethod Kjeldahl through precipitation of ca	asein YES				YES
ANALYSIS OF ALCOHOLIC BEVERAGES	Reference	DE-1626	Dosi-Fiber	EF-1425	Det-Gras N	Bloc-Digest
Degree of alcohol	Method Volumetric	YES				
Volatile Acid	Method Volumetric	YES				
Iron	Method Volumetric					YES
ANALYSIS OF FODDER AND RAW MATERIAS	Reference	Pro-Nitro M,S and A	Dosi-Fiber	EF-1425	Det-Gras N	Bloc-Digest
Raw Cellulose	Method Weende & Wijkströn		YES	YES		
Raw Protein	Method Kjeldahl	YES				YES
Raw Fat	Soxhlet estraction				YES	
Volatile base nitrates	Distillation Kjeldahl	YES				
Soluble Raw Fat en hydrochloric acid and Pipsin	Method Kjeldahl	YES				YES
Casein total	Method Kjeldahl in casein precipitate	YES				
Phosphorus	Method Photometric					YES
Hydrocyanic Acid	Distillation in Silver Nitrate	YES				
Pure zolidine	Soxhlet Extraction				YES	
ANALYSIS OF FRUIT BEVERAGES AND DERIV	VATIVES Reference	Pro-Nitro M,S and A	Dosi-Fil	ber [DE-1626	Bloc-Digest
Total Nitrogen	Method Kjeldahl	YES				YES
Volatile Acids	Distillation measurement				YES	
Arsenic	Dichromate Determination					YES
	Deference	Pro-Nitro M,S and A	Dosi-Fiber	D.Q.O./	ECO-8/ECO16	Bloc-Digest
ANALYSIS OF WATERS	Reference					_
	Method Kjeldahl	YES		-		YES
Total Nitrogen Total Phosphorus						YES YES

^{*} Reference: Methods are official publications issued by the MAPA (Ministry of Agriculture, Fisheries and Nutrition) 1993.



Dry block for Determination of Organic Nitrogen by the Kjeldahl method

Models Macro and Micro

THE EQUIPMENT FOR THE DETERMINATION OF ORGANIC NITROGEN IS MADE OF TWO BASIC ELEMENTS:

- BLOCK DIGESTER (MINERALISATION) WITH PROGRAMMABLE TEMPERATURE CONTROL AND GLASSWARE (MACRO 250 mi and Micro 100 mi).
- DISTILLATION UNIT "PRO-NITRO M", "PRO-NITRO S" (SEMIAUTOMATIC) AND "PRO-NITRO A" (AUTOMATIC).





Digestion Block "Bloc-digest"

FEATURES

Minimal sample manipulation

Uniform heating.

Capacity to store 20 programs of 4 steps for temperature and time.

RS-232 port for temperature register and digestion programation from a PC.

Gas collection system that does not require special water jet pumps.

Supplied complete with:

1 metal heater block.

1 programmer for time and temperature

1 tube support rack

1 gas collector

Digestion tubes.

Complete unit with: Dry-block connected to a programmable process unit, (time and temperature) rack support for tubes and fume extractor.



MODELS - COMPLETE UNIT MICRO

MICRO	Part No.	No. of positions
Bloc Digest m 12	4001047	12
Bloc Digest m 24	4001048	24
Bloc Digest m 40	4001049	40

MODELS - COMPLETE UNIT MACRO

MIODELO COMILECTE O	III IIIAUIIU	
MACRO	Part No.	No. of positions
Bloc Digest 6	4000629	6
Bloc Digest 12	4000630	12
Bloc Digest 20	4000631	20

ACCESSORIES: EXTRACTION SYSTEM AND NEUTRALIZATION OF ACID VAPOURS



Specially designed to absorb and neutralize vapours generated in the Kjedahl digestion processes.

Composed of a "Scrubber unit" that blocks the passage and neutralizes the condensed acids, and a re-circulating water pump that produces sufficient vacuum to aspirate the vapours.

It is essential to put the "Scrubber" unit together with the neutralizer solution between the digestion block and the re-circulating pump. Digester "Bloc Digest 20" requires 2 "Scrubber" units for a proper neutralization of fumes. Note: Pump with part number 4001612 has 2 suction inlets.

Reduces water consumption, no requirement for a constant connection for mains water.

Prevents emission of gas vapours and acidic water to waste.

Low noise level (<65dBA)

Re circulating pump made from chemically resistant materials

"Scrubber" unit

Part No.	Height / Width / Depth (exterior) cm	Weight Kg
4001611	32 31 16	2

It is supplied with 3 Kg. acid vapours neutralizer solution. **Spare part:**

3 Kg. acid vapours neutralizer solution.

Part No. 4001610

Water recirculation vacuum pump

Part No.	Height / Width / Depth	Vacuum level	Pump rate	Weight
	(exterior) cm	bar	Litres/minute	Kg
4001612	44 39 28	0.98	10	10



PARTS LIST AND ACCESSORIES

Heating blocks

To function correctly it is important to include in the purchase the controller unit for time and temperature RAT-2, The block should not be connected directly to the mains supply.



Models	Part No.	No. of positions	Ø tube mm	Height / (ex	Width / terior) o		Temperature °C	Power W	Weight Kg
MACRO	4000507	6	42	18	33	28	45 to 450	1500	18
MACRO	4000508	12	42	18	39	33	45 to 450	2100	25
MACRO	4000509	20	42	18	44	39	45 to 450	2500	31
MICRO	4001050	12	26	18	33	28	45 to 450	1500	16
MICRO	4001051	24	26	18	39	33	45 to 450	2100	22
MICRO	4001052	40	26	18	44	39	45 to 450	2500	27

Process programmer for time / temperature RAT-2. Part No. 4001538

Features

Temperature range from 45 °C to 450 °C. Memory for 20-4 steps programs. Maximum time per step: 600 minutes.

Acoustic indication for digestion end of program. Two selectable temperature gradients: kjeldahl/D.Q.O.

Temperature sense breakage alarm.

Independent control of safety maximum temperature. Bidirectional RS-232 serial connection. For temperature registration and digestion program edition with the RAT connected to a PC.


Software included.



Electronic controller RAT-2

Control Panel

The control panel and the RAT-2 display allow the creation and execution of the digestion program. During the digestion, it shows the block temperature, the elapsed time in the program step.




The software included makes easier the digestion program edition and the digester temperature register.

Tube support rack.

Made of chemically treated no corrosive Dural, with handles and side panel that prevents heat loss.



Models	Part No.	No. of places	•	/ Width .exterior) (
MACRO	4005071	6	15	17.5	12.5
MACRO	4005081	12	15	23	18
MACRO	4005091	20	15	28.5	23.5
MICRO	4001053	12	15	17.5	12.5
MICRO	4001054	24	15	23	18
MICRO	4001055	40	15	28.5	23.5





MACRO of 250 ml volume. Graduated to Series MICRO of 100 ml volume. 100 ml 42 mm Ø x 300 mm high. Part No. 4042300



26 mm Ø x 300 mm high. Part No. 4001045



Unit for collecting fumes supplied with a support rack. Made of borosilcate glass and stainless steel.



Models	Part No.	No. of places	Height / Wid (exterio	
MACRO	4005072	6	15 18	12
MACRO	4005082	12	15 23	18
MACRO	4005092	20	15 29	23
MICRO	4001056	12	15 18	12
MICRO	4001057	24	15 23	18
MICRO	4001058	40	15 29	23



Digestion and distillation tube Series Tube for digestion and distillation Water jet pump for vacuum extrac- Water Jet pump for vacuum extrac-Part No. 4000633



tion. Made of PVC designed for the fume tion. Metallic, suitable for the extraction extraction of the 20 and 40 place models. of fumes for the 6, 12 and 24 place models. Part No. 7000293



Steam distillation unit Kjeldahl "Pro-Nitro M"

DETERMINATION OF ORGANIC NITROGEN (KJELDAHL METHOD). Automatic Naoh dosage and temporize Stop.

Steam distillation Kieldahl unit.

Simple secure systematic analysis suitable for small to Medium throughput of samples.

FEATURES

Steam distillation system.

Compact steam generator with safety over temperature thermostat and over pressure device.

Safety door, the system will not operate if the door is open.

"Tube in place" sensor: if the tube is not located, the dosing process of NaOH will not take place.

Universal adapter for digestion/distillation tubes MACRO (Ø 42 mm) and MICRO (Ø 26mm)

Small footprint, saves bench top space: The H_2O and NaOH reservoirs are placed within the unit.

Stainless steel case with reinforced ABS plastic front. Automatic distillate titration kit. (See accessories).

SPECIFICATIONS

Measuring range: from 0.2 to 200 mg of Kjeldahl Nitrogen.

Programmable distillation time. Nitrogen recovery >99.5%

Distillation speed: from 35-40ml/minute Typical distillation time: from 7-10 minutes. Water consumption rate: from 80-100 litres/Hr. Steam generator water consumption: 2,5 Litres/Hr. Water reservoir for steam generator: 6 litres

NaOH reservoir: 2 Litres.

ALARMS

Low water level for the steam generator.

Safety door open or no distillation/digestion tube in place.

Steam generator over temperature.

AUTOMATED SEQUENCES

Open and closure of cooling water to the cooling coil. Automatic load of NaOH once the distillation has started. Select NaOH volume.

Stop at the end of the pre-set programmed time.

ADDITIONAL REQUIREMENTS

To complete Kjeldahl Nitrogen analysis a digestion block is also required. (See Bloc Digest pages 251 and 252).

MODEL

MODEL			
Part No.	Height / Width / Depth	Power	Weight
	cm	W	Kg
4002627	75 50 50	1800	30

Supplied with a MACRO Ø 42mm tube, set of reagent containers and tubing.





CONTROL PANEL

- 1. Low water indicator.
- 2. Door open or no tube presence indicator.
- 3. Over temperature indicator.
- 4. Mains on indicator.
- 5. Push button and indicator start/stop distillation.
- 6. NaOH volume selection.

ACCESSORIES

Tube for digestion and distillation Series MICRO of 100 ml volume.

Part No. 4001045

Digestion and distillation tube Series MACRO of 250 ml volume.

Part No. 4042300



Adapter kit for automatic determinations. Vessel with location positions for pH electrode, stirrer and reagents and distillate sample.

Part No. 4001724



Polycarbonate conical flasks durable. CAPACITY 250 ml

Part No. **5310100** with cap. Part No. **5310101** without cap.





Steam distillation unit Kjeldahl Semi-Automatic "Pro-Nitro S"

DETERMINATION OF ORGANIC NITROGEN (KJELDAHL METHOD). AUTOMATIC BORACIC AND NaOH DOSAGE, SAMPLE DRAINAGE AND TEMPORIZED STOP.

Semi-automatic steam distillation Kieldahl unit. Simple secure systematic analysis suitable for medium to large throughput of samples.



Steam distillation system. Compact steam generator with safety over temperature thermostat and over pressure device. Safety door, the system will not operate if the door is open.

"Tube in place" sensor: if the tube is not located, the dosing process of NaOH will not take place.

Universal adapter for digestion/distillation tubes MACRO (Ø 42mm) and MICRO (Ø 26 mm). Small footprint, saves bench top space: The H₂O, NaOH and H₃BO₃ reservoirs are placed within the unit.

Empty Digestion/Distillation tube system.

Stainless steel case with reinforced ABS plastic front.

Green LED 2 digit display.

Distillation program: (Add NaOH, Add Boric Acid, Distillation time, Empty tube.) Automatic distillate titration kit. (See accessories).

SPECIFICATIONS

Measuring range: from 0.2 to 200 mg Nitrogen.

Programmable distillation time.

Nitrogen recovery >99.5%

Distillation speed: from 35-40ml/minute Typical distillation time: from 7-10 minutes. Water consumption rate: from 80-100 litres/Hr. Steam generator water consumption: 2.5 Litres/ Hr. Water reservoir for steam generator: 6 litres

NaOH reservoir: 2 Litres. Boric Acid reservoir: 2 Litres

ALARMS

Low water level for the steam generator.

Safety door open or no distillation/digestion tube in place.

Steam generator over temperature.

AUTOMATIC

Single push button to start the distillation cycle:

- Boric acid dosing
- Start distillation.
- NaOH dosing
- Stop Distillation (Programmed time transpired.)
- Acoustic indicator at the end of the cycle.

ADDITIONAL REQUIREMENTS

To complete Kjeldahl Nitrogen analysis a digestion block is also required. (See Bloc Digest pages 251and 252).

MODEL

Part No.	Height	Height / Width / Depth cm		Power W	Weight Kg
4002851	75	50	50	1800	32

Supplied complete with a MACRO Ø 42 mm tube, set of connection tubes, set of reservoirs.

ACCESSORIES

Tube for digestion and distillation Series MICRO of 100 ml volume.

Part No. 4001045

Digestion and distillation tube Series MACRO of 250 ml volume.

Part No. 4042300



Adapter kit for automatic determinations. Vessel with location positions for pH electrode, stirrer and reagents and distillate sample.

Part No. 4001724



Polycarbonate conical flasks durable. CAPACITY 250 ml

6

09

1

000


Part No. 5310100 with cap. Part No. 5310101 without cap.






CONTROL PANEL

- 1. Illuminated indicator. Steam generator
- 2. Low water in the steam generator
- 3. Door open or no tube present indicator.
- 4. Over temperature indicator.
- 5. Push button and display to select parameters.
- Mode push button, Manual









10. Push Button, empty sample tube.



Automatic steam distillation unit Kjeldahl "Pro-Nitro A"



DETERMINATION OF ORGANIC NITROGEN (KJELDAHL METHOD)
FULLY AUTOMATIC OPERATION. FROM THE REAGENT DOSAGE TO THE TITRATION.

Steam distillation system Kjeldahl, complete with automated "ON-LINE" analysis (evaluation in real time). For systematic precise analysis, with minimum personnel intervention, simple and safe. Adequate for a laboratory with a medium to large throughput of samples.

The Kjeldahl steam distillation unit «PRO-NITRO A» evaluates the distillate at the same time as it is produced (evaluation «On-Line»), the evaluation and distillation are completed as one operation, reducing drastically the analysis time. This type of evaluation offers the following additional advantages: detects the point where the sample no longer produces Nitrogen, which means that, the distillation stops at the optimum maximum Nitrogen recovery and does not prolong the analysis longer than necessary.

The titration is a colorimetric method and is accepted by AOAC and does not require any periodic calibration.

FEATURES

Distillation by steam generation.

Automatic «On-line» colorimetric evaluation.

Steam generator with safety thermostat, over temperature and over pressure device. Safety, door closed, that prevents distillation if open.

Detects that a digestion/distillation tube is present. This prevents the dosing of NaOH if there is no tube located.

Universal adapter for MACRO (Ø 42 mm) and MICRO (Ø 26 mm) distillation tubes. **Space saving in the laboratory:** the reservoirs for the H2O, NaOH, Boric Acid and HCl are located inside the unit.

Empties the digestion/distillation tubes and the collector automatically.

Automatic stop when distillation is complete.

Large LCD display of 20 x 4 characters.

RS232 output to results printer.

Main system made from stainless steel with an ABS plastic front.

SPECIFICATIONS

Measuring range: 0.2 to 200 mg Nitrogen.

Nitrogen recovery:> 99.5%

Distillation speed: from 35 to 45 ml/minute

Coolant water consumption: 80 to 100 litres per hour. Steam generator water consumption: 2.5 Litre/Hr. Steam generator water reservoir capacity: 6 litres.

NaOH reservoir capacity: 2 Litres. Boric Acid reservoir capacity: 2 Litres. Titrant reagent reservoir capacity: 2 Litres.

Evaluation precision: 1.5% Minimum reagent dose 0.01ml.

ALARMS

Low water level for the steam generator.

Safety door open or no distillation/digestion tube in place.

Steam generator over temperature.

ADDITIONAL REQUIREMENTS

To complete Kjeldahl Nitrogen analysis a digestion block is also required. (See Bloc Digest pages 251 and 252).



MODEL

MODEL					
Part No.	Height /	Width	/ Depth	Power	Weight
		cm		W	Kg
4002430	75	50	50	1800	38

Supplied complete with a MACRO \emptyset 42 mm tube, set of connection tubes, set of reservoirs, 250 ml. of mixed indicator 4.8 and 100 gr. of sulphate ammonium.

AUTOMATION

Closing and opening of the condenser cooling water in line with the distillation process. Dosing of Boric Acid.

Dosing of NaOH once the distillation has started.

Select NaOH and Boric Acid volume.

«On-line» evaluation of distillate.

Auto detection of the end of the distillation process.

Special functions to maximise performance.

Special functions for maintenance.

REAGENTS

All the reagents used in the «PRO-NITRO A» are easily located:

- Solution of 30-40% NaOH.
- Solution of Boric Acid at 1% concentration (approx.) with mixed indicators (Bromo-cresol green and methyl red).
- Reagent for titration: HCl or H2SO4 from 0.05N or 0.25N adjusted to 0.001 Normal.



CONTROL PANEL

- 1. Menu to configure the date, time and selectable parameters.
- 2. Print the analysis information using the optional printer 4120113, purchased as an accessory.
- 3. <<ESC>> to cancel changes and escape from the menu.
- 4. Increase values and navigation through the menu.
- 5. Decrease values and navigate through the menu.
- **6.** <<ENTER>> to accept changes to parameters and navigation through the menu.
- 7. LCD display to visualise parameters and results

ADVANTAGES

Excellent precision on results.

Complete Nitrogen recovery from the sample.

Minimum operator intervention.

No calibration required.

RESULTS

The results can be downloaded to a printer (Optional), required for GLP, and includes the following data:

Minimum analysis time.

- Consecutive unrepeatable I.D. number of analysis.
- Date and time.
- Volume of NaOH.
- Volume of Boric acid.
- Reagent normality.
- Nitrogen detected.

15/10/05 12:16:08

Analysis Nr: 087598
NaOH: 75ml.
Boric: 25ml.

Normality: 0.1503

Results:

Reagent: 10.521ml

Nitrogen: 22.1382mg

ACCESSORIES



Ink printer (not thermal paper), size (4/6/10 cm) suitable for use with the PRONITRO A.

Paper 2 1/4" (56 mm) wide. Interface RS232.

Includes interface and mains cables.

Part No. 4120113

Digestion and distillation tube. Series MACRO of 250 ml volume. Graduated to 100 ml 42 mm \emptyset x 300 mm high. Part No. 4042300



Tube for digestion and distillation. Series MICRO of 100 ml volume. 26 mm Ø x 300 mm high.

Part No. 4001045





QUALITY CONTROL INFORMATION

ALL OF THE KJELDAHL DISTILLATION UNITS 4002430 REQUIRE A PROTOCOL ASSAY FOR THE RECOVERY OF NITROGEN WHEN MANUFACTURED.

THESE RESULTS COME WITH THE EQUIPMENT AND ARE VALID FOR IQ AND OQ CLARIFICATION.

COMPLEMENT



Digital colorimeter "Pro-A 520"

MICROPROCESSOR CONTROLLED.
AUTOMATIC ZERO ABSORBANCE AND 100% TRANSMITTANCE.
ALPHANUMERIC 20 CHARACTER 2 LINE L.C.D. DISPLAY.

APPLICATIONS

Reagent preparation for Pro-Nitro A.

FEATURES

Wavelength range: 400 to 800 nm, by using special fil-

ters.

Standard filter: 520 nm.

Expanded Absorption range: -0.3 to 3.5 O.D.

Transmission: 0 to 100 T %.

Photometric accuracy: >1 %. @ 1.000 A.
Photometric precision: ±1 %. @ 1.000 A.
Photometric stability: 0.004 A/hr. @ 0.000 A.
Light source: Long life tungsten lamp.

Detector: Solid state.

Sample chamber: 10 mm cuvettes.

Minimum volume: 1 ml.

Display: Alphanumeric LCD display of 2 lines of 20 cha- ON/OFF switch.

racters.

Calculation functions: Transmission T %.

Absorbence, Concentration by factor or standard concentrations.

Calibration: Self adjusting by software.

RS-232 interface.

CONTROL PANEL

Interactive LCD display.

Numeric and function keypad.

SPARE

Lamp of 6 V / 6 mm. Part No. **4512009**

MODEL

Part No.	Built in printer	Height / Width / Depth cm	, , , , , , , , , , , , , , , , , , , ,	
4120029	NO	11 18 28	10	4.5



Solvent recovery extractor for fats and oils "Det-gras N"

FOR THE DETERMINATION OF FATS BY SOXHLET METHOD. SAFETY PROTECTION 1P65.

Equipment for the extraction of fats and soluble substances that will dissolve in solvents. Based on the Soxhlet methods by Randall, the DET-GRAS-N extracts fats (and or other substances) from the sample by dissolving them in a solvent. If compared to the classic standard soxhlet method, a time saving of between 30-70% can be obtained due to a two stage extraction process and high speed in obtaining high temperatures.

APPLICATIONS

Practically all methods of classic soxhlet extraction can be substituted by the DET-GRAS N. Typical fat determinations are: in meat products, animal fodder, pre-prepared foods, fish etc. Also suitable for the extraction of soluble components in paper paste, textile fibres etc.

FEATURES

Heating by protected electric heater, distributed throughout the radiator surface that provides homogeneous heating throughout the samples.

Electrical safety according to IP65.

Robust temperature control by a durable thermocouple.

Safety over temperature device.

Outer case made of solvents resistant epoxy coated steel.

Extractions can be completed using glass or cellulose thimbles.

Compatible with most common solvents: Petroleum ether. Diethyl ether, hexane, acetone acetonitrile etc.

Supplied complete with 2 sets of seals, adapters for different solvent types.

Typical extraction time (food fats) 50 minutes.

SPECIFICATIONS

Cellulose extraction thimbles $\theta_{\rm ext}$ 26 x 60 mm. Reusable glass thimbles $\theta_{\rm ext}$ 34 x 80 mm. Working temperature of 90 to 240 °C.

Solvent recovery 60 to 80 %.

Solvent volume (by sample): up to 50 ml.

Program memories: 16.

Extraction time <
boiling>>: from 0 to 99 minutes.
Extraction time <<ri>from 0 to 99 minutes.
Time to recuperate solvents: from 0 to 99 minutes.




CONTROL PANEL

LCD display of 2 lines of 20 characters.

Keypad with 4 push buttons to select temperature, time and programs.

MODELS

MODELO							
Part No.	No. of places	Height A	Height / Width / Depth Cm		Water consumption litres/minute	Power W	Weight Kg
4002841	2	70	45	40	1	200	19
4002842	6	70	75	40	2	600	25





Accessories supplied with the 6 place model.

Extractor for the determination of cellulose and fibre "Dosi-fiber"

APPLICATIONS

Total Fibre (WEENDE, VAN SOEST or similar). Dietary fibre. Neutral detergent fibre. Acid detergent fibre. Other extraction processes which do not use acetic acid, acetic trichloric acid or nitric acid. Textile fibre. Wood and paper fibre.

EQUIPMENT DESCRIPTION

Integral extraction and filtration.

No sample transfer reduces the risk of sample loss, since the tubes, crucibles and filtre are transferred with the samples in place.

Excellent results due to reproducible operating conditions.

Easy sample handling due to the special crucible support stand.

Versatile and precise test procedure that allows the samples to be weighed at various stages in the extraction process.

FEATURES

Rugged external case with a "RILSAN" protective coating. All equipment components; condenser, valves, heater, compressed air valves and controls are all protected within the main case. Infrared heating elements.

CONTROL PANEL

Mains illuminated ON/OFF switch. Switch for compressed air pump. Heating element with electronic control.

TECHNICAL DATA

Sample size: 0.5 to 3 g (Normally 1 g). Reproducibility: approximately ± 1 % for fibre

level between 5-30 %. Measuring range: 0.1 to 100 %.

Cooling water consumption: 1 litre/minute.

EQUIPMENT

Comes complete with: crucible with a porosity P-2, Crucible support rack, heater lid and holder manipulator.

MODELS

Part No.	No. of places	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4000599	4	56 43 32	1000	19
4000623	6	56 57 32	1500	25



6 places extractor, part no. 4000623.

ACCESSORIES

Double hotplate for reagents. Power consumption 1750 W.

Part No. 4000634 Beaker for reagents. Part No. 1000635



SPARES

Support rackSupport rack4 crucibles.for 6 crucibles.Part No. 4000600Part No. 4000624

Crucibleswith a P - 2 porosity.
Part No. 4000601

EQUIPMENT FOR THE DETERMINATION OF TOTAL DIETARY FIBRE, ENZYMATIC METHOD

THE DETERMINATION OF TOTAL DIETARY FIBRE BY THE ENZYMATIC METHOD (AOAC, AACC) IS VERY DIFFERENT TO THE WEENDE AND THE DE VAN SOEST METHODS. THE EQUIPMENT FOR THIS ANALYSIS DEPENDS ON THE ENZYME KIT THAT YOU ARE GOING TO USE. ONCE THE TYPE OF ENZYME KIT IS SELECTED THE KIT INSTRUCTIONS WILL INDICATE THE NECESSARY EQUIPMENT TO COMPLETE THE ANALYSIS. THE COLD EXTRACTION UNIT EF-1425 PART NO 4001425 AND THE RECIPROCAL BATH PART NO 6032011 ARE FREQUENTLY USED (SEE PAGE 93).



Cold extraction unit "EF-1425"

DESCRIPTION

Specifically designed unit for cold extraction of up to 6 samples simultaneously, using acetone or other solvents. Complements the Dosi-Fiber for the analysis of crude fibre.

Components: 1 litre receptacle for solvent recovery, suction tube and 6 adapters to hold the filtration crucibles Part No. 4000601.

A pump or jet pump will be required to connect to the suction tube.

Comes complete with: blank stoppers in the case that not all of the 6 crucible places are used and a 1 litre solvent recovery vessel.

FEATURES

Stainless steel AISI 304, extractor, modular construction, easy to use and fast removal for efficient cleaning. Support and crucible rack, made from hardened glass.

MODEL

Part No.	•	/ Width erior) cı	Weight Kg	
4001425	45	20	34	3

ACCESSORY Jet pump. Part No. 7000293





Oenological Distillation unit "DE-1626"

APPARATUS TO OBTAIN THE DISTILLATION OF ALCOHOL FROM ALCOHOLIC BEVERAGES IN ACCORDANCE WITH CEE NO. 2676/90 (WINE ANALYSIS) AND CEE NO. 2870/2000 (ANALYSIS OF SPIRITS.)

APPLICATIONS

Degree of alcohol in wines, Degree of alcohol in spirits, Volatile acids and Sorbic acid.

FEATURES

Vapour distillation.

Designed especially to obtain the distillation of wines and spirits to determine the alcohol grade by volume, volatile acid, ascorbic acid and other components.

Distillation volume up to 400 ml.

Auto stop at the end of the distillation.

Simple operation with minimum service.

Vapour generator protection system.

- Tube Ø 42 mm for volatile acid samples.
- Tube Ø 52 mm for volatile acid samples and 100 ml proof alcohol.
- Tube Ø 80 mm for samples of 200 ml for proof alcohol.
- 200 ml Matrix flask to collect distillate.

TECHNICAL SPECIFICATIONS

Validated to the criteria of CEE No. 2676/90 and CEE No. 2870/2000 Distillation speed: 30-40 ml/min (Recover of 200 ml of distillate in 5/6 minutes).

Condenser coolant water rate: 80-100 l/hr.

Generator water consumption: 1 and 1.25 ml for each ml of distillate.

Vapour generator power: 2400 W.

CONTROL PANEL

Control Panel.

Electronic controller.

Mains switch.

START Push button.

STOP Push button.

ADD Push button.

By the use of the ADD button allows the distillate to set burette level, then, there is no need to fill up with water to a level of 200 ml.

VALIDATION

The oenological DE-1626 has been validated by the oenological laboratory INCAVI in Vilafranca del Penedés. This laboratory has completed analysis from the distillation unit DE-1626 for the determination of alcohol proof by volumetric analysis, following the criteria of CEE No 2676/90 and CEE 2870/2000. They have analysed samples of hydrated alcohol, different wines and spirits with a high alcohol content. The unit was successfully appraised by the laboratory.

MODEL

Part No.	Height .	/ Width cm	/ Depth	Power W	Weight Kg
4001626	90	30	30	2400	20

Supplied with tubes and flasks \emptyset 42, 52 and 80 mm and a 200 ml matrix flask.



Automatic termination of final distillation at 200 ml of distillate. By using the ADD function small quantities of distillate can be added to accurately fill the matrix flask to the mark without having to add water to the distillation.





ACCESSORIES



 Glass tube adapters.

 Part No.
 Dimensions

 4042300
 Ø 42 x 300 mm

 1001422
 Ø 52 x 300 mm

1000646 Ø 80 x 300 mm



Support rack for 6 tubes, made of stainless steel. Part No. Suitable tubes dimensions 4000648 Ø 42 mm

 4000648
 Ø 42 mm

 4001473
 Ø 52 mm

 4001613
 Ø 80 mm



Matrix flask 200 ml. Part No. 1001636

Alcoholmeters, graduations in 0.1 grade alcohol. Part No. Graduations

90 - 100

iii o.i yiau	c alcollol.
Part No.	Graduation
1001629	0 - 10
1001630	10 - 20
1001631	20 - 30
1001632	30 - 40
1001633	40 - 50
1001634	50 - 60

1001635



Alcoholmeters stand 6 places. Part No. 1000015

COMPLEMENTS FOR THE DENOLOGICAL DISTILLER "DE-1626"

ALL BASIC WINE ANALYSIS WITH A SIMPLE MANUAL EQUIPMENT, COMPOSED OF: OENOLOGICAL ANALYSER PHOTOMETER "M-3000", MICRO-CEN-TRIFUGE "CENCOM-1" AND DRY-BLOCK THERMOSTAT "ENOBLOC".



Oenological analyser photometer "M-3000"

SMALL BENCHTOP ANALYSER.

FOR COLOR MEASUREMENT, FOLIN-CIOCALTEAU INDEX AND FOR ENZYMATIC ANALYSIS.

APPLICATIONS

With the oenological application photometer, the following determinations can be performed:

- 1. Wine colouring intensity (absorbances addition at: 420, 520 and 620 nm.)
- 2. Folin index (polyphenols).
- 3. Grapes phenolic colour and maturity by means of Cromoenos® method.
- 4. Enzymatic analysis for grape juice, grape juice-wine and wine. Acids: acetic, malic, gluconic and lactic, glucose, fructose and glycerine.

FEATURES

Wavelength range: from 340 to 750 nm, with interference filters with a 10 nm band pass.

Manual filter selection (12 positions filter disc).

Supplied with filters: 340, 420, 520, 620 and 750nm. Other optional filters posible.

ABS, T% and concentration reading. ABS reading range: from -0.3 to 3.5 OD. T% reading range from 0 to 100%T.

Concentration range: from 0.001 to 9.999

Photometric accuracy: 1%.

Precision: ±1 %.

Photometric stability: better than 0.004 A/h. Light source: halogen lamp 6V/12W. Cell holder for 10 mm. path length cuvettes. Internal calibration by software.

CONTROL PANEL

Filter change knob.

Numeric keypad and for functions.

Alphanumeric display with interactive messages. RS-232 interface for connection to printer or computer.

MODEL

Part No.	Height .	/ Width mm	/ Depth	Power W	Weight Kg
4120019	11	18	28	35	3



SPARE

Halogen lamp 6 V 12 W. Part No. 4313040

Note: A visible ultraviolet range spectrophotometer UV-2005 is needed for a total polyphenols index determination (absorbance 280 nm) and the Cromoenos® method exact and accurate measurement. Part no. 4120020. Page: 235

ACCESSORIES

Glass cuvettes 10 mm. path length. 2 units pack. Part No. 5100021 Glass cuvettes 1 mm. path length, 2 units pack, Part No. 4120034

Polystyrene cuvette, semi-micro, 10 mm. path length. 100 units box. Part No. 5100023

Cuvettes adapter of 1 mm. Part No. 4120033



Microcentrifuge "Cencom I"

FOR MICROTUBES OF 2,2 / 1,5 - 0,5 AND 0,2 ML.

APPLICATIONS

To guarantee an excellent reproductibility in the enzymatic analysis, it is recommended to centrifuge at 14.500 rpm with 1,5 to 2 ml. tubes, to remove interferences (yeast, bacteria and colouring matter).

The centrifuge is essential for the grape extract cleaning by means of *Cromoenos*® method.

FEATURES

see page 114.

MODEL

Part No.	Tubes capacity	Height/ Width/ Depth (exterior) cm	Max. speed r.p.m.	R.c.f.(xg)	Power W	Time selected	Weight Kg
5022000	12	15 21 25	14500	14000	105	till 99' 99"	4

Short/Spin.

Comes complete with 12 adapters of 0.5 and 0.2ml.



Part No.	Tubes capacity	Height/ Width/ Depth (exterior) cm	Max. speed r.p.m.	R.c.f.(xg)	Power W	Time selected	Weight Kg
5022000	12	15 21 25	14500	14000	105	till 99' 99"	4

CONTROL PANEL

Push button Start/Stop.

Time selector push button.

Speed selector push button.

LCD display indicating time and speed.

Push button \(\rightarrow \) increase different parameters. Push button ▼ decrease different parameters.

Push button open lid when cycle has ended.

Push button motor drive rotation for quick acceleration

Start switch.

Metallic thermostat dry block "Enobloc"

APPLICATIONS

To use the enzymatic reagents performing at full capacity, it is recommended to incorporate the thermostat dry-block incubator for tubes/cuvettes at a fixed temperature of 35 °C. Part No. 7001570

Capacity: 15 cuvettes pathlength 10 mm, which can be macro or semi-micro of 1.5 ml..

Cromoenos® is a registered trademark by Bioenos S.L.



Part No.	. 3	/ Widt terior)	h / Depth cm	Stability °C	Power W	Weight Kg
7001570	11	18	28	±0,5	10	2,4



http://www.bioenos.com

(for analytical applications, please check out the website).

APPARATUS FOR WATER ANALYSIS



Thermo reactor for determining COD Chemical Oxygen Demand of waste water plants

STANDARD TO C.E.E. AND UNE 77-004. REFLUX METHOD.

FEATURES


Composed of a Hot block for 6,12 or 20 samples, process programmer and lifting racks. Uniform heating throughout the block to all samples. Automatic temperature and time control during digestion.

C.O.D. unit comprises of:

- 1 metallic hot block.
- 1 Processor for temperature and time.
- 1 Support rack for tubes.
- 1 Support rack for condensers.
- C.O.D. tubes with an aperture of 29/32.
- Condensers for C.O.D.

MODELS - EQUIPMENT DESCRIPTION

	Part No.	No. of places
DQO. - 6	4000638	6
DQO. - 12	4000639	12
DQO. - 20	4000640	20



INDIVIDUAL PARTS LIST AND SPARES

Hot block. This block cannot be used without the RAT controller of temperature and time.



Part No.	No. of places	Ø tube mm	Height / Width / Depth (exterior) cm	Temperature °C	Power W	Weight Kg
4000507	6	42	18 33 28	45 - 450	1500	18
4000508	12	42	18 39 33	45 - 450	2100	25
4000509	20	42	18 44 39	45 - 450	2500	31



System programmer for temperature/time RAT-2. Part No. 4001538

Display of time, temperature and program.

Selectable temperature from 45 to 450 °C. Memory

capacity: 10 programs of 5 stages of each program.

Maximum time per stage: 600 minutes. Visual and

acoustic alarm at the end of the cycle. Alarm if the

Condenser tube support. Part No.

4000643 6 places. **4000644** 12 places. **4000645** 20 places.



Tube support rack.

Made of chemically treated dural sheet metal, with handles and removable panels that maintain the heat around the tubes and allow the user to examine the samples during the process.

Part No.

4005071 for C.O.D. 6 (block 4000507). **4005081** for C.O.D. 12 (block 4000508). **4005091** for C.O.D. 20 (block 4000509).



Mains switch. Temperature display.

temperature probe fails.

Push button to select temperature.

Push button to increase value.

Push button to start.

Push button to stop.

Push button to decrease value.

Push button to select time.

Display time and program number Push button to select program.

Tubes for C.O.D. Part No. **1000641**









Laboratory flocculator"Flocumatic"

FOR THE DETERMINATION OF THE NECESSARY AGENTS REQUIRED FOR SEDIMENTATION. DIGITAL ELECTRONIC CONTROL OF SPEED AND TIME FUNCTIONS.

Maximum torque per stirring unit is 40 Ncm. Maximum viscosity 30.000 mPas.

APPLICATIONS

Optimization of coagulant additives and poly-electrolytes for surface and residual water treatment.

Evaluation of the efficiency of an absorbent in toxic agents.

Polyvalent. Without transmissions by belt. Due to its powerful torque it allows stirring and mixing high viscosity substances.

FEATURES

Stirring equipment for 4 or 6 places that accommodates beakers up to 1000 ml tall shape or 2000 ml short shape.

Stirring speed from 15 to 200 r.p.m.

Silent running.

Adjustable timer from 1 to 999 minutes or continuous operation.

Easily adjusted height over head stir, bars made of AISI 304 stainless steel are easily fixed in position.

Two vertical and horizontal adjustable illumination units can be placed in the unit. This adjustable accessory allows the light angle to be varied according to the samples being run. An original Selecta design. See accessories part numbers 3000834 and 3000915, 4 and 6 place models that obtain simultaneous front and back lighting. The main case of the unit is made of epoxy coated steel, with AISI 304 stainless steel supports with anti slip rubber feet.



- 1. Alphanumeric display. Indicates rpm and time.
- 2. Parameter invalidate button.
- 3. Parameter value decrease button.
- 4. Parameter value increase button.
- 5. Parameter validation button.





Flocculator "Flocumatic" with horizontal and vertical illumination (6 places).

Flocculator "Flocumatic" with horizontal illumination (4 places).



MODELS

Part No.	No. of places	Speed r.p.m.	•	3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3		Illumination power W	Power W	Weight Kg
3000833	4	15 - 200	52	70	23	40	160	24
3000914	6	15 - 200	52	99	23	50	170	36

ACCESSORIES

Illumination units, simultaneous lighting horizontal and vertical. Part No. Model part No.

3000834 3000833 4 places. **3000915** 3000914 6 places.





Portable flocculator "4P-M"

FOUR POSITIONS. WITH ILLUMINATION.

FEATURES

Mains voltage or portable connection to a cigarette lighter within a car (12 V DC), for achieving analysis for a site chemical dose measurements of waste water treatment plants.

Head stirrers made of AISI 304 stainless steel with variable height.

Suitable for vessels up to 1 litre capacity.

CONTROL PANEL

- Alphanumeric display. Indicates rpm and time.
- 2. Parameter invalidate button.
- 3. Parameter value decrease button.
- 4. Parameter value increase button.
- 5. Parameter validation button.



MODEL

Part No.	·	Width terior)	/ Depth cm	Power W	Power supply	Weight Kg
3000835	37	25	27	20	220-240V/50-60Hz	6



Made from stainless steel and glass, easy disassembly for cleaning and storage.

COMMON FEATURES

Easy to use, suitable for continuous distilled water requirements.

Boiler unit and heater element made of AISI 304 stainless steel, with lid and condenser made of orosilicate glass. The glassware pressure seal connection is an original Selecta design.

Safety system cuts off power to the heater element if there is not enough water. Distilled water produced is pyrogen free.

MODELS	Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20 °C	Ø / Height (ext.) cm	Power W	Weight Kg
Aquasel 4	4903004	4	1	3 µs/cm	18 43	2400	5
Aquasel 6	4903006	6	2	3 µs/cm	23 68	4800	8

Supplied complete with silicon seals.





Distillation still "L-3"

DISTILLATION CAPACITY: 3 LITRE PER HOUR.

FEATURES

Made entirely of borosilicate glass, with protective case support, designed for functionality, that permits easy access to the quartz heater elements.


Functions: automatic with continuous water production. Produces distilled water of a high purity, free of pyrogens and metal ions.

Equipped with quartz heater elements and protection device using a safety thermostat that is activated in the case of low water pressure to the condenser or boiler unit. The unit automatically reactivates when the water pressure is restored.

Suitable for wall mounting.

MODEL

Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20 °C		Height / Width / Depth (exterior) cm		Power W	Weight Kg
4903000	3	1	3 µs/cm	34.5	56	19	2200	4.5





Water distiller "D-4 Large"

DISTILLATION CAPACITY: 4 LITRES/HOUR.

FEATURES

Equipment entirely made of borosilicate glass, with protective case designed for functionality that allows an easy access to the glass and heating elements.

Automatic operation.

Continuous production of high purity distilled water, free of dissolved gases, metals, salts and pyrogens. Spiral heating system.

Double safety system:

By a presostat that turns the heater off in case of lack of water supply.
 Water level device that avoids the computer to function until the heating element is covered with water.

CONTROL PANEL

Power switch.

Lack of water flow indicator lamp.

MODEL

Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20° C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4903001	4	1	2 μs/cm	69 28 28	3000	6,3





Water distiller, "AC-L4" and "AC-L8"

DISTILLATION CAPACITY: 3.5 AND 8 LITRES/HOUR.

FEATURES

Tank and heating element made from stainless steel.

External case made from epoxy covered steel.

Easily dismantled for cleaning and service. All seals made from silicon.

In the case of over temperature or low water supply, a safety thermostat cuts power to the heater.

Water supply tube Ø 10 to 11 mm.

Distilled water tube connection Ø 12 mm.

Voltage: 230 V / 50/60 Hz

CONTROL PANEL

Illuminated power switch.

Safety switch thermostat with rearm system.

MODELS	Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20 °C	Height	/ Width cm	/ Depth	Power W	Weight Kg
AC-L4	4903007	4	1	2,5 µs/cm	44	25	23	3000	12
AC-L8	4903008	8	1,4	2,5 µs/cm	61	26	26	6000	14



Water distiller "L-4 Cabinet"

DISTILLATION CAPACITY: 4 LITRES/HOUR.

FEATURES

Equipment entirely made of borosilicate glass, with protective case with transparent plastic door that allows an easy access to the glass and heating elements.

The water distiller incorporates a boiler, a large surface cooling and two heating elements with protective quartz. Desktop model with device to hang it on the wall.

Continuous production of high purity distilled water free of dissolved gases, metals, salts and pyrogens. Easy to dismantle for a complete cleaning.

Automatic operation.

Totally automatic operation:

- · The equipment disconnects when the distilled water tank is full.
- The equipment reconnects when emptied.
- Automatic operation of the distiller.

Double safety system: 1. By a thermostat that turns the heater off in case of lack of water supply. 2. Water level device that avoids the computer to function until the heating element is covered with water.

CONTROL PANEL

Start up main switch.

Light indicators of heating, cleaning, of lack of water and tank full.

MODEL

Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20° C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4903002	4	1	2 µs/cm	40 55 22	1x3000	12



Water distiller "R-4 Reser"

DISTILLATION CAPACITY: 4 LITRES/HOUR. STORAGE TANK: 8 LITRES.

FEATURES

External case epoxy covered.

Cooling and heating elements made of stainless steel.

Easy to dismantle for cleaning.

Thermostat that disconnects the heater when there's lack of water supply.

Monitoring function that disconnects the distiller when the storage tank is complete

MODEL

full, and resumes as it empties.

Distilled water dispenser located in the front part.

Distilled water and water supply racords of 3/4 inch.

Desktop model with device to hang it on the wall.

CONTROL PANEL

Illuminated power switch.

Lack of water flow indicator lamp.





Note on the outlet water: Depending on the inlet water quality, it can be of type II or III (according to ASTM). All parameters have been tested under standard laboratory conditions.





Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20° C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4903005	4	1	2.5 us/cm	47 54 31	3000	24



COMECTA Water double distillation "L-4B"

DISTILLATION CAPACITY: 4 LITRES/HOUR.



FEATURES

Made completely from borosilicate glass with protective case that has easy access to the heater elements and glass parts.

Automatic function and continuous production.

Produces high quality double distilled water free from pyrogen and metallic ions. It functions through a dual distillation process:

The distilled water from the first stage is collected by the heater reservoir of the second stage, that starts the double distiller function when when the water's level is optimum, thus providing double distilled process.

Protection with a safety thermostat is activated if there is insufficient cooling water. Reactivated automatically when the water supply is re-established.

Water flow stop if there is an electrical fault.

If it is required to store the distilled water and to prevent over filling, the system stops, switches off power and disconnects the water flow when is full.

Voltage: 400 V / III+N, 50/60 Hz = 8.67 A



Part No.	Distillation capacity I/h	Water consumption I/min	Conductivity at 20 °C	Height/	Width cm	/Depth	Power W	Weight Kg
4903010	4	1.3	0.5 us/cm	46	66	50	2x3000	25



SPARES

Heater element made of quartz. Part No 4903011 Glass still, borosilicate glass (front -1). Part No. 4903012 Glass still, borosilicate glass (interior-2). Part No 4903013 Condenser, borosilicate glass (the equipment comes with 2 condensers. Price provided is for 1 spare condenser). Part No. 490301



Water distiller, specially for autoclaves "Dest-4"

DISTILLATION CAPACITY: 1.5 LITRES/HOUR. INTERNAL CAPACITY 4 LITRES.

Case made from thermo resistant plastic. Internal stainless steel. Cooling by forced air though a condenser. Without glass heater elements, Easy to use. Safety cut out if over temperature due to insufficient water.

Supplied complete with a 4 litres plastic bottle and three dispenser tap with filter. Note: for daily use it is recommended to change the filter every 3 months.

No external additional water required.



SPARE. Dispensing filter tip. Part No. 4001730

MODEL

Part No.	Distillation capacity I/h	Conductivity µs/cm	Ø Height cm	Power W	Weight Kg
4001729	1.5	5	29 39	750	3.5

COMECTA Descaler "C-3"

APPLICATIONS

Pre-treated water that contains large quantities of CaCO3 (lime). For general laboratory use.

Recommended for use with distillation units where water hardness is above 25° french.

TECHNICAL SPECIFICATIONS

Metal container of AISI 304 stainless steel with regeneration of salt by using a two way reverse valve system. Resin capacity: 12 litres. Regeneration salt: 2 Kg. Regeneration cycle: 1200 litres to 60° french / 4800 litres to 35° french.

Regeneration cycle for distillation of: 300 to 800 litres. Maximum mains water pressure: 4.5 Kg/cm². Maximum supplied admissible hardness: 60° french.

Out put hardness: 1° french.

External measurement: 62 cm high x 19 cm Ø.

Weight: 20 Kg. Part No. 0703052



Comes complete with in and out 3/4" tubing.

COMECTA 50 Litre distilled water reservoir

FEATURES

Recommended for the storage of distilled water. Made of high density polyethylene. 50 litre capacity, with handles and dispenser tap. 15cm Ø top lid with stopper. Dimensions: 62 cm high x 38 cm Ø.

Part No. 0106006



Water reservoir trolley.

Robust steel trolley with epoxy coating. Lockable wheels, Inner section for holding jars and carboys etc.

Fits containers 38 cm Ø.

Height 48 cm.

Part No 5903049






ANTIVIBRATION BALANCE TABLE



MELTING POINT APPARATUS



PRECISE BALANCES



LIGHT SOURCE



MICROSCOPES



REFRACTOMETERS

Precise balances
Balance accessories and precision weights
Stereomicroscopes
Microscopes
Microcameras
Refractometers
Polarimeter
Melting point apparatus

page267pages268 to 269pages270 to 271pages272 to 278page278pages279 to 281page281

282

page

Precise analytical balance "FA-2204B"

FEATURES

Digital 7 segment screen. Time to stability: 8 seconds.

Adjustable time and sensitivity control.

ON/OFF automatic zero adjustment.

Automatic internal 200g calibration.

Tare range by subtraction up to the maximum capacity.

Piece counting function.

Density measurement in the base of the unit.

Units of measurement in g, ct, and oz.

Level adjustment (bubble).

Data interface: RS-232 for connection to a printer.

Operable ambient temperature 5 to 40 °C Power supply: 110/220 V-50/60Hz. Supplied complete with cover.



MODELS	Part No.	Maximum weight capacity g	Precision ±g	Linearity ±g	Height / Width / Depth (exterior) cm	Height / Width / Depth (Draft shield) cm	Ø plate cm	Weight Kg
FA-2204B	5830039	200	0.0001	0.0002	35 21 34	24 19 17.5	8	8.5

COMECTA Electronic Balance "Laborcom"



COMMON FEATURES

Easy to read seven segments display. Stabilization time: 3 seconds. ON/OFF. Automatic Zero adjustment External automatic calibration.

Tare range by subtraction up to maximum capacity. Piece couting function. Overload indicator. Operating temperature: from 5 to 35 °C. Voltage: 220 V ±10 V / 50 Hz.

MODELS	Part No.	Maximum weight capacity g	Precision ±g		/ Widtl exterior)		n Ø plate cm	Weight Kg
LABORCOM 110-L	5830025	110	0.01	7	17.5	24	12.5	1.6
LABORCOM 210-L	5830026	210	0.02	7	17.5	24	12.5	1.6
LABORCOM 600-L	5830031	600	0.1	7	17.5	24	12.5	1.6
LABORCOM 1100-L	5830032	1100	0.1	7	17.5	24	12.5	1.6
LABORCOM 2100-L	5830033	2100	0.1	7	17.5	24	12.5	1.6
LABORCOM 11000-L	5830027	11000	1	7	17.5	24	16.5 x 16.5	1.8



Shaking monitor for blood collection "HEDO2"

FOR BLOOD COLLECTION PROCESSES IN BAGS

FEATURES

HEDO2 is a digital blood collection monitor, created to control and monitoring blood collection processes. It has been designed to perform simple, precise and accurate operations.

It is composed of:

Easily removable tray.

Wired or wireless interface networking.

Automatic clamp closing at the end of the donation.

Automatic tare calculation.

Real time monitoring of the donation showed on the digital display.

Stand alone mode (operating with rechargeable battery). 8 hours life.

Sensor to detect the tube presence.

Power supply: 110/220 V-50/60Hz.

MODEL

Part No.	Max weight	Precision	Resolution	Height / Width / Depth	Power	Weight
	Capacity g	g	g	(external) cm	W	Kg
5010000	999	± 3	1	19 32 32	60	5

It is supplied with a transport suitcase.



ACCESSORIES

Wireless data transfer Hedo-2. Part No. 5010001 Wireless connection module Hedo-2, Part No. 5010002 Barcode reader Hedo-2. Part No. 5010003 Alarm lamp pool Hedo-2. Part No. 5010004



WEIGHING DISHES "ACCU-WEIGH"

Anti static. Reproducible tare weighing. Contamination free. Large selection of sizes.

Common features:

Made from high density polystyrene with flat base, round edge, contamination free surface and hydrophilic. Flexible in form for easy manipulation of the sample either liquid or solid.

Manufactured in two forms:

- Diamond form, colours white and black,
- Square form, colour white.

Part No.	Form	Capacity	Height /	Width	/ Depth	Weight	Colour
5726004	Diamond	5	3.8	31	54	0.25	White
5726005	Diamond	5	3.8	31	54	0.25	Black
5726006	Diamond	30	13.5	55	85	1.00	White
5726007	Diamond	30	13.5	55	85	1.00	Black
5726008	Diamond	100	20	92	132	2.6	White
5726009	Diamond	100	20	92	132	2.6	Black
5726010	Square	7	8.5	44	44	0.65	White
5726011	Square	100	24	80	80	2.2	White
5726012	Square	250	24	134	134	6.1	White

Comes in packets of 100 units.



ALUMINIUM WEIGHING BOATS

Static free, pressed dishes with rolled edge.

Withstands temperatures up to 450 °C.

Due to their form the dishes can be used one on top of another to protect samples. Suitable for testing, solid samples and mixtures of emulsions for moisture content, by differential weight loss method.

Useful for weighing, comparing, freezing or evaporation.

Consumable, single use.

Part No.	Thickness mm	Height / Ø (usable) mm	Comes in
5726000	0.2	5 42	packets of 112 units
5726001	0.2	10 57	packets of 110 units



ANTIVIBRATION BALANCE TABLE

Formica covered metal table with central polished granite plate to support the balance, mounted on silent blocks.

Adjustable shock device. Maximum weight 40 Kg.

Part No.	Height / Width / Depth	Height / Depth	Weight
	(exterior) cm	(balance support) cm	Kg
5838101	80 90 70	45 39	70





BALANCE DENSITY MEASUREMENT "PR-SE" (PATENTED)

Applications: System for reading the density of liquids using the immersion of a glass sinker of known volume. Can be used for opaque and viscous liquids for all classes of liquid.

Features: Made of borosilicate glass of low dilatation coefficient and a constant volume. Mode of operation:

- Locate the tube with liquid on an electronic balance.
- Push the TARE to achieve a zero stable reading
- Introduce and submerse the measuring suspended sinker centrally without tilting the tube, keep suspended.

At this moment the screen will indicate the DIFFERENTIAL WEIGHT, that is exactly the same as the density of the liquids multiplied by 10.

To obtain a density value in mS, it is necessary to have a balance with a minimum resolution of $0.01\ \mathrm{g}$.

Part No.	Density interval	Maximum tolerance	Temperature
	g/ml	error g/ml	assay/reference °C
5600000	0.600 - 2.000	±0.0005	20 / 20



Calibrated precision weights class E2 and F1. Certification W.E.C.C. (Western European Calibration Co-operation)

COMMON FEATURES

In the manufacturing of these weights the meteorological requirements are applied to conform with the OIML verification for accuracy, raw material, geometric form. identification and conservation. To quarantee their accuracy, the individual weights are supplied in a partitioned wooden case with a clasp lid and a small tray housing for locating the mg weights. The wooden boxes are lined with velvet for the g and Kg weights.

Weight	Part No. Weight Class E2	Part No. Certificate E2	Tolerance Max. error ± mg	Part No. Weight Class F1	Part No. Certificate F1	Tolerance Max. error ± mg
1 mg	5838102	5838124	0.006	5838146	5838169	0.020
2 mg	5838103	5838125	0.006	5838147	5838170	0.020
5 mg	5838104	5838126	0.006	5838148	5838171	0.020
10 mg	5838105	5838127	0.008	5838150	5838172	0.025
20 mg	5838106	5838128	0.010	5838151	5838173	0.03
50 mg	5838107	5838129	0.012	5838152	5838174	0.04
100 mg	5838108	5838130	0.015	5838153	5838175	0.05
200 mg	5838109	5838131	0.020	5838154	5838176	0.06
500 mg	5838110	5838132	0.025	5838155	5838177	0.08
1 g	5838111	5838133	0.030	5838156	5838178	0.10
2 g	5838112	5838134	0.040	5838157	5838179	0.12
5 g	5838113	5838135	0.050	5838158	5838180	0.15
10 g	5838114	5838136	0.060	5838159	5838181	0.20
20 g	5838115	5838137	0.080	5838160	5838182	0.25
50 g	5838116	5838138	0.10	5838161	5838183	0.30
100 g	5838117	5838139	0.15	5838162	5838184	0.5
200 g	5838118	5838140	0.30	5838163	5838185	1
500 g	5838119	5838141	0.75	5838164	5838186	2.5
1 Kg	5838120	5838142	1.5	5838165	5838187	5
2 Kg	5838121	5838143	3.0	5838166	5838188	10
5 Kg	5838122	5838144	7.5	5838167	5838189	25
10 Kg	5838123	5838145	15	5838168	5838190	50

Weight material

From 1 to 5 mg = aluminium, density approx. 2.7 g/cm^3 .

From 10 to 500 mg = alpaca, density 8.6 g/cm^3 .

From a and Ka = antimagnetic polished stainless steel, density approx. 8.0 a/cm³. Without adjustable chamber E2 and with adjustable chamber F1. The calibration certificate is expedited by DKD (member of the W.E.C.C.), which contains the determination and deviation between the nominal mass values, the actual measured value and the tolerance measurement.



SET OF STANDARD WEIGHTS WITH BOX AND MATCHING CERTIFICATE

Set of weights E2 of 1 mg to 500 mg Part No. 5831822 Set of weights F1 of 1 mg to 500 mg Part No. 5832822

Material: Wooden box (E2), or resin box (F1), with clasp, tweezers and a small brush.

Weight details: 1 mg, 2 mg, 2 mg, 5 mg, 10 mg, 20 mg, 20 mg, 50 mg, 100 mg, 200 mg, 200 mg and 500 mg. Total 12 weights = 1.11 g

Set of weights E2 of 1 mg to 50 g Part No. 5831302 Set of weights F1 of 1 mg to 50 g Part No. **5832302** Certificate for set of weights E2 of 1 mg to 50 g Certificate for set of weights F1 of 1 mg to 50 g

Part No. 5862301 Part No. 5862401

Material: Wooden box with clasp, lined with velvet, a special glove, housing for extractable fractional weights, tweezers and a small brush.

Weight details: 1 mg, 2 mg, 2 mg, 5 mg, 10 mg, 20 mg, 20 mg, 50 mg, 100 mg, 200 mg, 200 mg, 500 mg, 1 g, 2 g, 2 g, 5 g, 10 g, 20 g, 20 g y 50 g Total 20 weights = 111.11 g

Set of weights E2 of 1 mg to 100 g Part No. **5831303** Set of weights F1 of 1 mg to 100 g Part No. 5832303 Certificate for set of weights E2 of 1 mg to 100 g Part No. 5862302 Certificate for set of weights F1 of 1 mg to 100 g Part No. 5862402

Certificate for set of weights E2 of 1 mg to 500 mg Part No. 5862350

Certificate for set of weights F1 of 1 mg to 500 mg Part No. 5862450

Material: Wooden box with clasp, lined with velvet, a special glove, housing for extractable fractional weights, tweezers and a small brush.

Weight details: 1 mg, 2 mg, 2 mg, 5 mg, 10 mg, 20 mg, 20 mg, 50 mg, 100 mg, 200 mg, 200 mg, 500 mg, 1 g, 2 g, 2 g, 5 g, 10 g, 20 g, 20 g, 50 g and 100 g Total 21 weights = 211.11 g

Set of weights E2 of 1 mg to 200 g Part No. 5831304 Set of weights F1 of 1 mg to 200 g Part No. 5832304 Certificate for set of weights E2 of 1 mg to 200 g Part No. 5862303 Certificate for set of weights F1 of 1 mg to 200 g Part No. 5862403

Material: Wooden box with clasp, lined with velvet, a special glove, housing for extractable fractional weights, tweezers and a small brush.

Weight details: 1 mg, 2 mg, 2 mg, 5 mg, 10 mg, 20 mg, 20 mg, 50 mg, 100 mg, 200 mg, 200 mg, 500 mg, 1 g, 2 g, 2 g, 5 g, 10 g, 20 g, 20 g, 50 g, 100 g, 200 g y 200 g Total 23 weights = 611.11 g

Stereomicroscope "XTX-3C LED and ZTX-20 LED"







APPLICATIONS

By design and quality optics, these microscopes are recommended especially for long duration work. Ideally suited for biology, botany, archaeology, geology, teaching and electronics.

COMMON FEATURES

Three dimensional optical resolution, wide field of view without side aberrations, excellent crisp contrast and true colour.

Rechargeable battery 3x1800Ah and external charger.

Binocular head inclined at 45° with adjustable interpupillary distance from 55 to 75 mm.

Rotating head of 360° with easy access for large items.

Dioptric correction in the left eyepiece of ±5dp.

Ocular head held with secure fixing screws.

Rubber eyecups.

Pair of eyepieces, wide field of view WF 10x (23 mm \emptyset). In the model ZTX-3C. Pair of eyepieces, wide field of view WF 20x (20 mm \emptyset). In the model ZTX-20. Objectives of 2 x and 4 x, located in a double objectives support turret with rotation in both directions.

Double focus: by means of a vertical slip clutch of the head up and down a column, with fine adjustment using a bilateral macrometer with rack and pinion tension and control.

Solid base consisting of a platform with double contrast (black and white) circular plate and translucent etched light diffuser, with fixing screw and holding clips.

Model XTX-3C LED: System with double illumination of low voltage 12 V/10 W, with incident light (episcopic) in the head and transmitted (diascopic) in the circular platform.

Model ZTX-20 LED: system with triple illumination of low voltage, with incident light (episcopic) in the head and transmitted (diascopic) in the circular platform, Activated individually or simultaneously, with controllable intensity. Also the incident light incorporates a condenser lens with an extractable blue filter.

Dimensions: XTX-3C LED: 270 high x 117 wide x 190 mm deep. Dimensions: ZTX-20 LED: 350 high x 320 wide x 190 mm deep.

Stereomicroscope model XTX-3C LED. Part No. 5313113
Stereomicroscope model ZTX-20 LED. Part No. 5313090

ACCESSORIES

Eyepiece with wide field of view WF 5x (single).

Eyepiece with wide field of view WF 10x (single).

Eyepiece with wide field of view WF 10x (single).

Eyepiece with wide field of view WF 15x (single).

Eyepiece with wide field of view WF 20x (single).

Eyepiece with wide field of view WF 20x (single).

Eyepiece with wide field of view WF 20x (single).

Eyepiece with wide field of view WF 20x (single).

Eyepiece with wide field of view WF 15x (single).

Eyepiece with wide field of view WF 20x (single).

Part No. 5313091

Part No. 5313092

Part No. 5313093

Part No. 5313093

Part No. 5313093

Part No. 5313093

Part No. 5313095

SPARES

Set of LED lamps (for XTX-3C LED).Part No. 5313013Set of LED lamps (for XTX-20 LED).Part No. 5313015

MODELS	Part No.	Objectives	Eyepieces	Magnification	Ø Field of view mm	Working Distance mm	Illumination	Weight kg	
VTV acien	5313113	2x	WE 10v	20x	10.0	E 7	Transmitted or incident	0.0	
XTX-3C LED	3313113	4x	W F-IUX	WF-10x 57 Transmitted or incide 40x 5.0	ITAIISIIIILLEU OF IIICIUEIIL	2.8			
ZTX-20 LED	5313090	2x	WF-20x	40x	5	40	Transmitted or incident	3.6	
TIV-50 FED	4x			80x	2.5	40	independent or simultaneous	J.U	

Comes complete with key for tightening the axis tension, spare lamp and plastic dust cover.

COMECTA Stereomicroscopes with Zoom "SQF-L LED, SQF-D LED and SQF-E LED" IN THE SQF-





APPLICATIONS

By design and quality optics, these microscopes are recommended especially for long duration work. Ideally suited for biology, botany, archaeology, geology, teaching and electronics.

Model SQF-L LED.

COMMON FEATURES

Variable magnification using a continuous zoom control, with control knobs on either side to enable parfocal centralised focus.

Three dimensional optical resolution, wide field of view without side aberrations, excellent crisp contrast and true colour.

Optics are protected against knocks.

Binocular or Trinocular head (Model dependent), inclined at 45°.

Adjustable interpupillary distance from 54 to 75 mm.

Rotating head of 360° with easy access for large items.

Dioptric correction of ±5 dp in both evenieces.

Ocular head fixed with secure screws.

Rubber eyecups.

Pair of eyepieces wide field of view WF20x (12 \emptyset) or WF10x (22 \emptyset) model dependent.

Macrometric focus: by means of a vertical slip clutch with a fine adjustment using rack and pinion tension control.

Solid base consisting of a platform with double contrast (black and white) circular plates and translucent etched light diffuser, with fixing screw and holding clips.

Triple illumination LED system, low voltage with incident light (episcopic) in the head and transmitted (diascopic) in the circular platform: activated individually or simultaneously, with controllable intensity.

Model SOF-E LED trinocular

Same characteristics to the SQF-D LED model, except that the head is trinocular and it is supplied with a photography adapter and eyepiece WF 10x with an adjustable focus visor. The photographic eyepiece and camera with support T2 PK and trigger cable have to be ordered as accessories.

ACCESSORIES

1100200011120	
Supplementary objective 0.5x. SQF-L	Part No. 5313096
Supplementary objective 0.75x. SQF-L	Part No. 5313097
Supplementary objective 1.5x. SQF-L	Part No. 5313098
Supplementary objective 2x. SQF-L	Part No. 5313099
Eyepiece with wide field of view HWF 10x (single).SQF-D & E	Part No. 5313081
Eyepiece with wide field of view HWF 15x (single).SQF-D & E	Part No. 5313082
Eyepiece with wide field of view HWF 20x (single).SQF-D & E	Part No. 5313083
Photographic adapter 4x (single). SQF-E	Part No. 5313078
Video camera for eyepiece (See page 278)	Part No. 5435030
Adapter CCD of 0.5x suitable for a video camera.	Part No. 5313089
Eyepiece WF10x with grid (single). (Model SQF-D and E)	Part No. 5313084
Wooden case with handle and lock.	Part No. 5313016

SPARES

Set of LED lamps SQF-L Part No. 5313015
LED lamp 3 W SQF-D LED and SQF-E LED (transmitted). Part No. 5313044
LED lamp 3 W SQF-D LED and SQF-E LED (incident). Part No. 5313045

Technical table with the different magnifications, working distance and Ø of the field of view resulting in the configuration of the diverse eyepieces and supplementary objectives available.

MODELS	Part No.	Head	Zoom Range	Eyepieces	Standard magnification	Ø Field of view mm	Working distance mm	Illumination system	Weight kg
SQF-L LED	5313311	Binocular	1x - 4x	10x	1x - 4x	26 - 5.8	90	Transmitted	5.5
SQF-D LED	5313309	Binocular	0.67x - 4.5x	20x	13.4x - 90x	17.9 - 2.7	105	or incident or both	6.1
SQF-E LED	5313308	Triocular	0.67x - 4.5x	20x	13.4x - 90x	17.9 - 2.7	105	simultaneously	6.4

Comes with a dust cover.

COMECTA Microscope Monocular "M-100 FL LED"

APPLICATIONS

Directed at the requirements of higher education, universities and suitable for: veterinary, dentistry, biology, industry, etc.

FEATURES

Monocular head inclined angle 45° and rotatable through 360°.

Eyepiece, with wide field of view WF 10x (18 Ø) fixed in position using secure screw.

Objectives, parfocal DIN 45 mm achromatic brightfield of 4x (A.N. 0.10), 10x (A.N. 0.25) and 40x retractable (A.N. 0.65).

Total magnification: 40x, 100x and 400x,

Revolving triple objectives turret, rotation in either direction with a fixed axis stop locator.

Focus, vertical slip clutch, bilateral micrometer adiustment of 0.002 mm, separate macrometer with adjustable tension. Safety fixing screw for tension adjustment and end of focus limit stop.

Stage with spring clip slide holder of 110 x 126 mm. Condenser, fixed lens and aperture (A.N. 0.65), with diaphragm disk of 5 apertures.

Illumination, transmitted light (diascopic), LED 0.5 W, with collimator lens and blue filter for chromatic corrections.

Robust construction, with power switch and illuminator built into the base.

Dimensions: 360 height x 130 width x 175 mm depth. Weight 2.8 kg.

Supplied complete with key to adjust the tension of the slip clutch and cover.

Microscope monocular M-100 FL LED (400x).

Part No. 5313119

ACCESSORIES

Eyepiece Widefield WF-15x with cross hairs (Pointer). Part No. 5313076

Eyepiece Widefield WF-20x with reticule.

Part No. 5313077 Condenser ABBE 1.25 A.N. with filter and iris diaphragm. Part No. 5313003

Graduated mechanical stage.

Part No. 5313009

Video Micro-camera, 3.1 MP-:

- with connection USB for direct connection to a PC, including software.

Part No. 5435034

Wooden case with carry handle and lock.

Part No. 5313022

SPARE

LED lamp 0.5 W. Part No. 5313028



ACCESSORY

Set of polarisers with circular graduated base and filter. Part No. 5313075



COMECTA Microscope Binocular "701 LED"

FEATURES

Binocular head with 45° tilt and full 360° rotation. **Interpupillary** distance: 52 to 72 mm.

Dioptric adjustment of ±5 dp in both eyepieces, sharpness compensation for visual differences.

Pair of eyepieces DIN standard WF10x (18 mm Ø) and WF16x(11mm \emptyset).

Objectives DIN standard, achromatic, brightfield of 4x (A.N. 0.10), 10x (A.N. 0.25), 40x (A.N. 0.65) retractable and 100x (A.N. 1.25) immersion oil and retractable. Colour coded objectives for easy identification.

Total magnification: 1600x.

Four position objectives turret rotating in both directions with stop locator when it is in the correct position. Large sub stage 135x145 mm, variable height with a rack and pinion mechanism with bilateral coaxial control. Micrometer adjustable stage graduated to 0.002 mm and macrometer adjustment tension, Movement up to 30 mm.

Safety stop with pre-set adjustable safety level.

Adjustable slide holder, with 0.1 mm scale with a smooth X-Y 40 mm longitudinal and 80 mm traversal movement.

Condenser Abbe of 1.25 A.N. with diaphragm iris and exchangeable filters. Rack and pinion height, adjustable.

Illumination transmitted (diascopic), with removable condenser lens and pre-centred lamp LED 3W, with variable intensity control. Suitable for use with an external light source using a plan or concave reflector mi-

Robust construction that incorporates the illumination source and mains switch in the base.

Size: height 370 x width 200 x depth 225 mm.

Weight 6 Kg.

Supplied complete with blue set filters, matt glass filter, fuse, oil, dust cover.

Power supply: 110/220 V-50/60Hz.

Microscope Binocular 701 LEd. Part No. 5313111



ACCESSORIES

Eyepiece with wide angle lens WF-5x (single). Part No. 5313034

Eyepiece with wide angle lens WF-15x (single). Part No. 5313004

Eyepiece with wide angle WF-10x with reticule (single).

Part No. **5313035**

Eyepiece with wide angle WF-10x with pointer (single).

Part No. 5313007

Video camera for eyepiece (see page. 278) Part No. 5435030

SPARE

LED lamp 3 W. Part No. 5313027

COMECTA Polarising Binocular Microscope 206 LED



APPLICATIONS

Wide range of applications, includes bright-field, dark-field, polarization and phase contrast. Ideal for universities, clinics, biology, bacteriology, agriculture and industry.

FEATURES

Binocular head inclined 30° rotatable through 360°.

Adjustable interpupillary distance of 55 to 75 mm.

Diopter adjustment of ± 5 dp. in the right eyepiece, compensation for any differen-

Paired eyepieces DIN, wide field of view WF10x (180) and a pair of WF16x (0). Objectives DIN, achromatic brightfield of 4x (A.N. 0.10), 10x (A.N. 0.25), 40x (A.N. 0.65) retractable and 100x (A.N.1.25) oil immersion retractable; colour coded for rapid identification.

Total magnification: 1600x

Quadruple revolving objectives turret, rotates in both directions with notched location positions.

Large stage of 160 x 140 mm, height adjustable by mechanical rack and pinion bilateral coaxial focus; 0.002 mm graduated micrometer and macrometer with adjustable tension (distance 30 mm).

Pre-adjusted, Safety stop limit.

Mechanical stage. (slide holder), graduated scale of 0.1 mm, adjustable in both (x-y) using a smooth rack and pinion movement of 50 mm longitudinal and 70 mm


Condenser Abbe of 1.25 A.N. with iris diaphragm and replaceable filters.

Adjustable height by rack and pinion.

Illumination: transmitted (diascopic), with removable condenser lens and LED 3 W lighting with regulable intensity.

Strong and robust, illuminator incorporated into the base.

Set of simple polarisers: controllable analysis polariser.





Set of simple polarisers.

Comes complete with a set of blue, green and yellow filters, set of simple polarisers, lamp, 2 spare fuses, bottle of immersion oil and plastic cover

MODEL	Part No.	Height / Width / Depth mm	Power supply	Weight Kg
206	5313115	425 245 335	220 V / 50 Hz	7.5

SPARE

3W LED lamp. Part No. 5313046

ACCESSORIES

Eyepiece with wide angle lens WF-5x Part No. 5313034 (single) Eyepiece with wide angle lens WF-15x (single) Part No. 5313004 **Eyepiece with wide angle lens WF-10x** with reticule (single) Part No. 5313035 **Eveniece with wide angle lens WF-10x** with pointer (single) Part No. 5313007

Complete kit for phase contrast, made of: Objectives DIN plan negative phase contrast: 10x (A.N. 0.25), 40x (A.N. 0.65) and 100x (A.N. 1.25) retractable oil immersion, condenser turret phase A.N. 1.25, central telescope and blue filter.

Part No. 5313017

Darkfield condenser: Permits the observation of low contrast and small size, dispersed in a homogeneous substrate samples which they appear luminescent under a uniform darkfield.

Darkfield condenser A.N. 0.83 - 0.91 dry (operable under low magnification).

Part No. 5313018

Darkfield condenser A.N. 1.25 - 1.36 oil immersion (operable under high Part No. 5313019 magnification).

Video camera for eyepiece (see page 278)

Part No. 5435030



5313017



5313018





5435030



Microscopes "3000-A LED, 3000-B LED and 3000-C LED"



EXCELLENT RESOLUTION PLANACHROMATIC OPTICS. HIGH SPECIFICATION OPTICAL MICROSCOPES. MAXIMUM VERSATILITY WITH A WIDE RANGE OF CONFIGURATIONS, MODELS SHOWN ARE BINOCULAR AND TRINOCULAR WITH PHASE CONTRAST, BRIGHT AND DARK FIELD OPTICS, THIS INSTRUMENT IS IDEAL FOR A WIDE RANGE OF APPLICATIONS.

Optics Anti mould / Anti fungal



APPLICATIONS

Biology, bacteriology, cytology, immunology, diagnostics, pathology, agriculture,industry and university laboratories.

3000-A LED (BINOCULAR) AND 3000-B LED (TRIOCULAR).

FEATURES

 $\mbox{\bf Head, binocular or trinocular}$ (model dependant), inclined $30^{\rm o}$ with $360^{\rm o}$ rotation.

Interpupillary distance adjustable from 53 to 75 mm.

Dioptric adjustment of ± 5 dp in the left eyepiece, to compensate any visual differences.

Pair of DIN eyepieces wide field of view WF 10x (20 mm \emptyset).

DIN objectives, planachromatic bright field of **4x** (A.N. 0.10), **10x** (A.N. 0.25), 40x (A.N. 0.65) retractable and 100x (A.N. 1.25) oil immersion retractable, colour coded for rapid identification.

The planachromatic objectives provide an excellent image to the edge of the optical field of view.

Quadruple revolving objectives turret, rotates in both directions with spring ball stop locator when it is in the correct position.

Large substage of **180x140 mm**, rack and pinion adjustable height. Coaxial focus using controls located on both sides, graduated **micrometer** in steps of 0.002 mm and **macrometer** (with adjustable tension of up to 30 mm).

Safety stop with present adjustment.

Adjustable slide holder with a graduated scale of 0.1 mm, made up of a smooth (X-Y) rack and pinion mechanism, longitudinal and transversal movement of 50 mm and 75 mm respectively.

Abbe condenser of 1.25 A.N. precentred, with double lens, iris diaphragm and retractable from the front. Adjustable height using a rack and pinion mechanism.

The Abbe condenser can be easily replaced for special applications with a phase contrast condenser or a darkfield condenser (see accessories).

Transmitted illumination with LED that is formed with a condenser lens, viewing diaphragm removable filters and a precentred, low voltage, adjustable intensity.

Robust stable structure that includes the illumination system in the base, adjustable intensity potentiometer and On/Off switch.

Dimensions: 410 Height x 200 Width x 310 mm Depth.

Microscope Binocular Part No. 3000-A LED. Part No. 5901981 Microscope Triocular Part No. 3000-B LED. Part No. 5901982



Model "3000-B". Part No. 5901982.

MODEL 3000-C LED PLAN PHASE CONTRAST - BINOCULAR.

FEATURES

Similar characteristics to the model 3000-A LED except for a special configuration of Plan Phase Contrast objectives:

The objectives, **DIN quality, are planachromatic with positive phase contrast** of 10x (A.N. 0.25), 25x (A.N. 0.40) 40x (A.N. 0.65) retractable and 100x (A.N. 1.25) retractable and oil immersion.

The condenser is a phase turret A.N. 1.25, with central telescope and a set of white, blue and green filters. With sharp detail for examination of transparent thin sections that can differentiate the refraction index that the human eye cannot differentiate. Excellent for applications in biology, medicine, crystallography, agriculture, industrial plastics and rubber etc.

Microscope Binocular model 3000-C LED. Part No. 5901983

MODELS	Part No.	Head	Eyepieces	Objectives	Technique	Illumination	Weight kg
3000-A LED	5901981	Binocular		Planachromatic (PL) 4x - 10x - 40x - 100x	Brightfield		8.4
3000-B LED	5901982	Triocular	WF 10x	r ialiaciliulilatic (i L) 4x - iux - 4ux - iuux	brightheiu	3 W	8.7
3000-C LED	5901983	Binocular	(20 Ø)	Planachromatic (PL) phase contrast 4x - 25x - 40x - 100x	Phase contrast	LED	8.4

Supplied complete with matt glass filter, fuse and a dust cover.

SPARE	ACCESSORIES
LED lamp 3. Part No. 5313026	See page: 278



Trinocular fluorescence microscope "3002-F"



EXCELLENT RESOLUTION. POLYVALENT. CAN WORK IN EPIFLUORESCENT OR TRANSMITTED LIGHT MODE. WITH TRINOCULAR ADAPTER FOR CAMERA AND VIDEO CAMERA. PHASE CONTRAST AND DARK FIELD MODES AS ACCESSORIES.

OPTICS ANTI MOULD / ANTI FUNGAL

APPLICATIONS

Biology, cytology, oncology, genetics, scientific investigation, clinical diagnostics, epidemic prevention and embryology.

FEATURES

Trinocular head. inclined 30° and rotatable through 360°.

Interpupillary distance adjustable from 55 to 75 mm.

Dioptric adjustment of ±5 dp in the left eyepiece, to compensate for any visual differences.

Pair of wide field DIN eyepieces WF 10x (22 Ø).

Bright field planachromatic DIN objectives of 4x (A.N. 0.10), 10x (A.N. 0.25), 40x (A.N. 0.65) retactable, and an oil immersion retractable 100x (A.N. 1.25) they are colour codified for a quick identification.

Planachromatic objectives provide an excellent image to the edge of the optical field of view.

Quadruple revolving objective turret, with both sides rotation and with a roller baring and a notch stop.

Large sub stage of 180x150 mm, rack and pinion height adjustable, by coaxial focus bilateral controls, graduated micrometer in steps of 0,002 mm and macrometer with adjustable tension (30 mm path). With transmitted light lock plate and **UV** radiation protection screen.

Safety stop with preset adjustable safety level.

Mechanical slide holder (moveable), with graduated scale of 0,1 mm, made up of a smooth movement mechanism (x-y), by a 50mm longitudinal and a 75mm transversal rack and pinion movement.

Abbe condenser of 1.25 A.N., precentred, with double lens, iris diaphragm and retractable frontal lens. Rack and pinion height, adjustable.

Abbe condenser can be easily replaced by a phase turret condenser (for contrast) or dark field condenser, for special observation conditions (see accessories).

Transmitted illumination (diascopic), according to Kölher, that consists of a condenser lens, a field diaphragm, removable filters and low voltage halogen precentred lamp of 12V/30W with adjustable intensity.

Robust stable structure that includes the illumination system, the intensity potentiometer and the On/Off switch on its base.



Lamp power supply connected to 110 V/220 V with intensity, voltage and power adjustable by a knob.

Illumination module for reflected fluorescence (Waveband B and G, N free).

The unit includes diaphragm, selectable attenuation filter and 5 positions filter disk: ultraviolet, violet, blue and green filters, and a empty position for transmitted light.

ACCESSORIES

Waveband V and UV. Part No. 5901986

MODEL	Part no.	Eyepieces	Objectives	Technique	Illumination	Weight kg
3002-F	5901984	WF 10x (22 Ø)			Transmitted 12v / 30w & Mercury fluorescence 100w / DC	10,5

Supplied complete with a set of blue, green and white filters, bottle of immersion oil ηD 1.515, mercury lamp 100 w, spare fuse and dust cover.

SPARE PARTS Halogen lamp 12 V / 30 W. Part No. 5313021. Mercury lamp 100 W DC. Part No. 5313055

Fluorescense omega filters: U V. Part No. 5313060 Violet. Part No. 5313061 Blue, Part No. 5313062 Green. Part No. 5313063

ACCESSORIES FOR MICROSCOPES MODELS "3000-A LED", "3000-B LED", "3000-C LED" AND "3002-F"

Eyepiece with wide angle lens WF-15x (single). Part No. 5313073. Eyepiece with wide angle lens WF-20x (single). Part No. 5902074

Contrast phase kit: Objectives DIN, planachromatic positive phase contrast of 10x (A.N. 0.25), 25x (A.N. 0.4), 40x (A.N.0.65) retractable and 100x (A.N.1.25) retractable oil immersion and phase condenser turret, A.N. 1.25, with central telescope and green filter. Part No. 5901985

Dark Field condenser: Permits the observation of low contrast particles dispersed in a homogenous substrate, that appear luminous in a dark field.

- Dark field condenser: A.N. 0.77 ~ 0.91 dry (low magnification). Part No. 5313071
- Dark field condenser: A.N. 1.25 ~ 1.40 immersion (high magnification). Part No. 5313072

Micrograph Kit, camera mechanism WF 10x with a visor for focus correction and a photographic objective of 2.5 -4x. (Adaptable for Digital reflex photographic kit). Part No. 5902008

Digital reflex photographic kit Composed of a mount, a shutter and a SD memory card of 2Gb. Part No. 55313222

Adapter for a CCD video camera: attachment. Part No. 5313070

Wooden case. Part No. 5313016







5901985

5313222



COMECTA Microscope biological inverted, trinocular "XDS-1B LED"







APPLICATIONS

Optical microscope with many benefits, specially designed for biomedicine, biology and for culture analysis in matrix flasks and in microtitre plates, flasks, culture vessels. tubes etc. through the use of its objectives and long working distance condenser, technique uses brightfield or phase contrast.

FEATURES

Robust stable design with lateral arm wrist supports for maximum comfort and easy manipulation.

Trinocular head, inclined 45°.

Interpupillary distance adjustable from 53 to 77 mm.

Diopter correction of ±5dp in the left eyepiece, used to correct differences in sharpness.

Pair of DIN wide field WF 10x (22 Ø), telescopic central evepiece phase contrast DZ11x DZ11x and trinocular photography eyepiece \mathbf{S} 5x (9.5 \emptyset).

Objectives DIN planachromatic brightfield,long working distance LWD PL of 10x (A.N. 0.25/working distance W.D.=4.3 mm), 20x (A.N. 0.40/W.D.=8 mm), 40x (A.N. 0.40/W.D.=3 mm) and on negative phase contrast LWD PL ph of 10x (A.N. 0.25/W.D.=5 mm), colour coded for rapid identification. Maximum thickness acceptable for receptacle bases 1.2 mm.

Quadruple objectives turret with rotation in both directions with a roller baring and notch stop. Rack and pinion height adjustment mechanism with bilateral coaxial focus controls: micrometer of 0.2 mm rotation graduated at 0.002 mm increments and macrometer with adjustable tension.

Total working distance 12 mm.

Condenser with long working distance (55 mm), brightfield and phase contrast. Precentred, with iris diaphragm and flip out phase contrast. **Including a phase ring** of 25x. Rack and pinions mechanism for height regulation or lateral removal for more height receptacles of up to 145 mm.

Large stage: dimensions 227x208 mm. Includes 3 stages for optimum accommodation of differing sample vessels: plates and or culture flasks, microtitre plates, etc.

Mechanical sub-stage with graduated scale of 1 mm and vernier of 0.1 mm, made up of a smooth mechanism (x-y) rack and pinion drive with 77 mm longitudinal and 114 mm transversal.

Illumination system, diascopic transmitted light by a halogen lamp of 6 V / 20W, precentred by external knobs and controllable intensity.

The unit includes a field diaphragm with filter carousel of green and grey.

MODEL	Part No.	Height / Width / Depth exterior cm	Weight Kg
XDS-1B LED	5901988	58 30 53	17



Comes complete with 2 flasks, 2 culture plates, wooden holder for eyepieces, objectives and accessories, 2 lamps, 2 spare fuses and a protective cover.

ACCESSORIES

Objectives planachromatic long working distance

Negative phase contrast:

Objective PLL ph 20x

(A.N. 0.40 W.D.=8 mm)

Part No. 5313065

Objective PLL ph 40x

(A.N. 0.60 W.D.=3,5 mm)

Part No. 5313066

Adapter thread CCD 0.5 x

Part No. 5313067

Adapter thread 1 x

Part No. 5313068

Eyepiece WF10x (22 mm) With reticule

Part No. 5313064



5313065

5313066

5313064

Digital reflex photographic kit, composed of a mount, a shutter and a SD memory card of 2Gb. Part No. 5313222



5313222

3W LED lamp. Part No. 5313029

SPARES





COMECTA Microscope metallographical inverted "18AT"

EXCELLENT RESOLUTION. PLANACHROMATIC OPTIC.



APPLICATIONS

Microscope with high capabilities.

Made for large sample observations. Samples are placed on top of an open obstruction free stage.

Binocular or Trinocular configuration with a range of accessories, that convert the instrument into an versatile microscope ideal for the identification and analysis of different materials and alloys. Suitable for a wide range of scientific research applications in metallurgy, physical properties of metals, heat treatments and quality control.

FEATURES

Head, binocular or trinocular (model dependent), 45° inclined angle.

Adjustable interpupillary distance from 53 to 77 mm.

Dioptric adjustment of ±5 dp in the left eyepiece, compensating for any visual differences.

Paired eyepieces, DIN quality. Wide field of view WF 10x (22 mm Ø).

Objectives DIN grade, planachromatic brightfield of **10x** (A.N. 0.25), **20x** (A.N. 0.50), **50x** (A.N. 0.70), retractable, and **100x** (A.N. 0.85), retractable, using immersion oil. Colour coded for easy and rapid identification.

Magnification range: 100x, 200x, 500x and 1000x.

The planachromatic objectives provide an excellent image, free from aberrations at the edge of the field of view.

Quadruple rotating objectives turret, rotates in both directions with stop locator when in the correct position.

Large sub stage of **242x200 mm**, rack and pinion adjustable height, with bilateral coaxial focus. Micrometer graduated in 0.002 mm steps and **macrometer**, with adjustable tension up to 19 mm.

Safety stop with a pre-set safety level.

Slide holder with scale of 0.1 mm, with a smooth X and Y rack and pinion movement of 15 mm longitudinal and 15 mm traversal movement.

Illumination: incidental (episcopic), Kölher, that consists of a condenser lens, condenser Abbe of 1.25 A.N., pre-centred with an iris diaphragm,



filters and a pre-centred halogen 6 V / 30 W low voltage controllable light source. **Structure, robust** with an illumination system in the base, intensity potentiometer and illumination switch.

Dimensions: 40 height x 25 Width x 61 cm depth.

Weight: 15 Kg.

Comes complete with a set of filters: blue, yellow, green and matt glass filter.

Microscope Triocular Part No. 18AT. Part No. 5901990

SPARE

Halogen lamp 6 V / 30 W. Part No. 5313069

ACCESSORIES



5313222

Digital reflex photographic kit composed of a mount, a shutter and a SD memory card of 2Gb. CMOS sensor higher than 10 megapixels. Speed up to 3 pictures per second. LCD display of 2,5 inches. Processed with total image control. Large area viewfinder and high luminosity.

Part No. 5313222

Eyepiece with wide field, of view WF-10x (22 mm) with reticule. Part No. 5313064

Adapter thread CCD 0.5x.

Part No. 5313067

Adapter thread CCD 1x.

Part No. 5313068

Micrograph Kit comprising of:

Photographic adapter mechanism with eyepiece WF-10x with focus adjustment visor and objective of 4x (adaptable to **Digital reflex photographic kit**).

Part No. **5902014**



COMECTA: Microcameras "R-350"





FEATURES

Revolutionary for its small size. Quality image, digital technology. Can be directly used with a monitor or TV via an a euro connector. Can be input to a video for recording images. (Suitable for a PC with a PCI card for image processing and storage on diskettes or a hard drive). Designed for different applications, especially for diagnostics of dermatological illnesses of the skin scalp or capillary fibre. Other magnifications available: 10x, 100x and 200x.

Sensor	2
Resolution lines	U
Pixels (H x V)	0
SystemPAL/NTS0	C
Minimum illumination	Χ
S/R better than	b
Horizontal line frequency15625 H	Z
Vertical line frequency	
Light source for model	V
Power consumption2V	V
Voltage6V (220-50Hz)
Dimensions: Model R-350-350x	n
Dimensions: Model R-350-2x	n
Weight 0.3 Kg	n

MODELS

R-350-2x. 2 magnification, with euro connector.

Part No. 5435038 System PAL Part No. 5435032 System NTSC

Suitable for microscopes:

3000B Part No. *5901982* **3001F** Part No. *5901984* **SQF-E** Part No. 5313308 **XDS-1B** Part No. 5901988 18AT Part No. 5901990

R-350-350x. 350 magnification, with euro connector, suitable for monitor or TV monitor.

Part No. 5435035 System PAL Part No. 5435036 System NTSC



Microcamera model R-350 2x, suitable for microscopes.



COMECTA: Evepiece video camera

ADAPTS DIRECTLY ONTO THE EYEPIECE OF THE MICROSCOPE OR TRINOCULAR HEAD.

Micro video camera 3.1 MP with a USB connection, direct to a PC, includes software, Part No. 5435030





COMECTA: Light source, cold source, fibre optic "Z-150" **CONSTANT COLOR**

APPLICATIONS

Microscopy illumination for small samples free from colour alteration.

Especially recommended for microphotography of live biological samples, cultures, micro-organisms, process control, electronic components, dental studies, etc.

Halogen lamp 150 W. Adjustable illumination intensity control. Forced air circulation. Heat absorption filter. Safety transformer fitted. Primary fuse protection.

CONTROL PANEL

Mains switch. Potentiometer controlled light intensity.

MODEL

Part No.		/ Width terior c	n / Depth cm	Power W	Weight Kg
5901016	15	19	25	150	4

SPARE Lamp 150 W. Part No. 5901017

Supplied with double flexible arm 50 cm long.





COMECTA: Microscope accessories and consumables



GLASS SLIDES

Glass slides (smooth): 26 x 76 x 1.5 mm.

Part No. 5600400 Box of 72.

Glass slides (matt end for notes): 26 x 76 x 1.5 mm.

Part No. 5600399 Box of 72. Cover slip: 18 x 18 mm. Part No. 5600396 Box of 100.

Cover slip: 22 x 22 mm. Part No. 5600397 Box of 100.

Cover slip: 24 x 24 mm. Part No. 5600398 Box of 100.



SLIDE DISPENSER

Base made of ABS with a transparent Polycarbonate lid, with rotating manual dispenser of individual slides. Capacity 50 slides of 76 x 26 mm.

Part No. 1001366



CLEANING KIT

Suitable for optical microscopes.

Part No. 5313000



WOODEN MICROSCOPE BOX

Solid wood, microscope box, of 3 mm thickness, treated with varnish, fitted with a metallic door, clasp lock and a carry handle.

Microscope Box suitable for Model M-100FL, XTX-3C and ZTX-20.

420 high x 255 wide x 285 mm deep. Part No. 5313022 Microscope Box suitable for Model SQF-L, SOF-E. 2000-A. 2000-B and 2000-C.

Internal dimensions: 420 high x 265 wide x 317 mm deep. Part No. 5313016

Slide support, Slide box for storing slides, slide staining stands, etc. (See pages 215 to 219).



COMECTA Refractometer ABBE "2WAJ"

LCD DIGITAL THERMOMETER BUILT-IN.

APPLICATIONS

Ideal for measuring the refractive index ηD of liquids, solids and sugar concentrations (0 - 95% = η D 1.333-

Indispensable instrument for industrial chemistry, petroleum, pharmaceutical, food oils and fats, sugars, universities etc.

FEATURES

Telescopic system X2 magnification.

Eyepiece with dioptric correction.

Connection for external water bath.

Thermometer 0-70 °C. To see the sample temperature. Robust metallic construction.

Weight: 2.6 Kg.

Measurements: 240 mm high x 100 mm wide x 200 mm depth.

Measuring range:

Refraction index η D: 1300 - 1700. Brix: 0 - 95%.


Measuring definition:

Refraction index η D: 0.0005. Precision ηD : ± 0.0002 . Brix: 0.25%. Precision: ±0.1%.

Part No. 5901008

Supplied complete with

- Aluminium carry case,
- 1 bottle of Naphthalene Bromide,
- 1 calibration standard,
- 1 screwdriver and dust cover.





Digital ABBE refractometer "WYA-2S"

APPLICATIONS

Ideal for measuring the refractive index of liquids, solids and sugar concentrations (° Brix), industry: chemical, petroleum, food, pharmaceutical, universities and research etc.

FEATURES

Apparatus equipped with an optical system made of a telescope with dispersion correction, used to define the division line between dark and bright areas.

Precise optics with automatic temperature compensation displayed as Brix and sample temperature.

Stabilized temperature connection for a water bath. RS-232 output with software included.

Measuring range:

Refraction index η D: 1.3000 - 1.7000

Height / Width / Depth

(exterior) cm

Brix Bx - TC: 0 - 95% Brix Bx: 0 - 95%

MODEL Part No.

CONTR	OL PANEL				
1		P		12 _/ 11	
2				Ϋ́	10
3					_
4		ILE	7///	// /c	7
5		A 500	7//	116	///
6				<u>3</u> (4
7				る \	\mathcal{N}
8		1	1 813		//
	- 1	1	a		ال

Measurement scale definition:

Refraction index ηD : ± 0.0002 Brix Bx - TC: ± 0.1%

Brix Bx: ± 0.1%

Temperature:

Working range: 0 - 50 °C.

Correction range of Bx, in respect to tempera-

ture: 15 - 45 °C.



Supplied complete with dust cover and maintenance box that includes:

- 1 bottle of Naphthalene Bromide,
- 1 calibration standard,
- 1 Hex key,

Weiaht

Kg

supply

38 18 33 220V ±20V 50Hz 10

- 3 lamps 6.3V, 0.25 A.
- 3 fuses 1 A.

Optimum precision



Technical description

- 1. Monocular system.
- 2. Telescope adjustment.
- 3. Digital display.
- 4. Mains switch "POWER".
- 5. Push button "READ".
- 6. Push button for Brix with temperature correction "Bx TC".
- Push button for the refractive index "ηD".
- Push button for Brix without temperature correction "Bx".
- Optical control.
- 10. Light intensity unit.
- Refractometer prism.
- Push button for temperature indicating. "TEMP".

COMECTA: Refractometers



COMMON FEATURES

Robust metal construction.

Excellent optical quality.

Eyepiece with adjustable dioptric ring to compensate any visual sharpness.

Easy to use, with direct easy to read scales.

Adjustable calibration screw.

Supplied complete with pouch, Pasteur pipette and screwdriver.



REFRACTOMETER FOR SALINITY \(\infty \) AND DENSITY "C-1" **APPLICATIONS**


For the investigation in clinical laboratories (solutions and dissolutions), Industry. food provisions (vegetables, fruit, fish and shell fish).

TECHNICAL SPECIFICATION

Automatic compensation

Temperature:10 - 30 °C. Salinity range: 0 to 100 ‰. Division Scale: 1 ‰.

Part No. **5901001**






REFRACTOMETER BRIX% INDEX "C-2", "C-3", "C-4" AND "C-5" **APPLICATIONS**

Concentrates of sugar in fruit, juices, gaseous drinks, milk, wine, industrial oils etc.

TECHNICAL SPECIFICATION

"C-2" Range Brix: 0 - 32 %. Precision±0.2 %.

"C-4" Rango Brix - 10 %. "C-5" Rango Brix 0 - 18 %. Division scale:0.1 %.





Part No.

5901002 Model C-2

5901020 Model C-2 with temperature compensation

5901003 Model C-3


5901021 Model C-3 with temperature compensation


5901022 Model C-4


5901023 Model C-4 with temperature compensation

5901024 Model C-5

5901025 Model C-5 with temperature compensation








CLINICAL REFRACTOMETER 2 SCALES "C-6" APPLICATIONS

Suitable for the determination of proteins, serum and specific gravity (density) of urine.

TECHNICAL SPECIFICATION

Proteins in serum: 1 - 12 g/100 ml.

Part No. **5901006**





COMECTA

COMECTA Digital portable refractometers "NR-101" and "NR-151"



APPLICATIONS

Model NR-151. Measurement of % Brix and refractive index η D.

Ideal for concentrations of sugars in % Brix and the measurement of refractive index ηD . For fruit, juices, drinks, milk, wine, edible oils, food industry, chemicals, pharmaceutical, science, petrochemical, oils, lubricants etc.

Model NR-101. % Brix, refractive index η D and % salinity.

Ideal for concentrations of sugar % Brix, measuring refractive index ηD and NaCl in % of salinity.

Ideal for fruit, juices, drinks, milk, wine, edible oils, food industry, chemicals, pharmaceutical, science, petrochemical, oils, lubricants etc.



Portable and easy to use.

Alphanumeric display.

Determinations using a high precision photo cell integrated to a CCD camera and digital image processor.

Chamber prism of anti corrosive stainless steel and lid.

Automatic temperature compensation (ATC) to 20 °C, with a sample temperature range from 5 to 40 °C.

Temperature selection °C or °F.

Calibration using distilled water.

Function selection limits and tolerable alarms.

Microprocessor with instant response.

Disconnects automatically if inoperable for 3 minutes.

Low voltage alarm.

Power supply 9V battery.

Operating temperature from 0 to 40 $^{\circ}$ C at up to 90% relative humidity and a maximum height of 2000 m above sea level.





MODELS	Part No.	Reading and Measurement	Measurement range	Measurement Resolution	Measurement Precision	Temperature Range	Temperature Resolution	Temperature Precision	Height/Width/ cm	Depth	Weight Kg
ND 454	5004040	Brix %	0.0 - 95.0%	0.1%	±0.1%						
NR-151	5901012	Refraction index η D	1.3330 - 1.5318 17 D	0.0001 Ŋ D	±0.0002 17 D	0 +40 °C	0.1 °C	±1 ℃ ±2 ℉			
		Brix %	0.0 - 60.0%	0.1%	±0.1%		32 +140 °F 0.2 °F		6 18 10	10	0.9
NR-101	5901013	Refraction index η D	1.3330 - 1.4419 17 D	0.0001 Ŋ D	±0.0002 17 D	02 · 140 1					
		Salinity %	0.0 - 28.0 %	0.1 %	±0.1 %						

Comes complete with protective carry case, dosing burette, bottle of distilled water and a 9V battery.



COMECTA Disk polarimeter "polar-2010"

LAURENT TYPE MODEL
SODIUM LAMP ILLUMINATION



APPLICATIONS

Designed for colleges and universities, chemical laboratories, pharmaceutical, clinics, cosmetics and food industry.

FEATURES

Vernier reading definition. Angle \pm 0,05 °.

Measuring range of optical rotation: \pm 180 °.

Rotation scale division: 1º.

Measuring lens 4x.

Ocular with dioptres adjustment.

Slide rack compartment length: for tubes up to 220 mm.

Monochromatic lamp wave length: 589,3 nm.

MODEL

Part No	External measurements cm.	Consumption	Weight
	Height / Width / Depth	W	Kg
5120210	32 43 14	20	5



SPARE	PART	S
Cadium	lamn	Day

Sodium lamp. Part No. 5120201



COMECTA Automatic Melting Point Apparatus "WRS-1B" "WRS-2A"



AUTOMATIC DETERMINATION NO VISUAL CONTROL REQUIRED.

A different conception.

Totally automatic, without the need to visually determine the melting point process

APPLICATIONS

Research laboratories, quality control. Test and control of purity after synthesis in pharmaceutical laboratories.

AUTOMATIC MELTING POINT WRS-1B

Single sample analysis

LCD display, simple alphanumeric keypad.

Melting point determination by photoelectric detection.

Automatic functions. Reads and registers the starting point, fusion point and end point.

Part No. 5609401

AUTOMATIC MELTING POINT WRS-2A

Simultaneous analysis, three samples.


LCD graphic display and 16 keys numeric keypad.

Melting point determination detector photoelectric.

Automatic functions. Reads and registers the starting point fusion point, end point,

melting point curve and average value determination.

Part No. 5609402



Details of graphic display, calculated average of melting point curves.



MODELS		WRS-1B	WRS-2A	
Part No.		5609401	5609402	
Digital LCD screen		Numeric	Numeric and graphical	
Maximum capillary tube capacity, 1.4 mm Ø ext.		1	3	
Temperature range		from 50 °C	to 300 °C	
Readable temperature resolution		0.1	°C	
Melting point measurement precision		Upto 200 °C ±0.5 °C a	nd from 200 °C to 300 °C ±0.8 °C	
Selectable ramp rate in °C/minute		0.2 - 0.5 - 1 - 1.5 - 2	- 3 - 4 and 5 °C/min.	
		0.4 °C	0.3 °C ramp rate of 1 °C/min.	
Repeatability		ramp rate of 1 °C/min	0.2 °C ramp rate of 0.2 °C/min.	
	increasing	from 50 °C of 300 °C 3 minutes		
Decending time to reach the start temperature	decreasing	from 300 °C of	50 °C 5 minutes	
Sample holder		for capillary tubes of 1.4 mm Ø ext. / 1 mm Ø int.		
Sample fill height		3 mm		
Output interface RS-232		9600 baudrate, 1 bit stop and 8 data bits		
Power		110 W		
Power supply		220V-50/60Hz		
Weight		9.5 Kg	12.5 Kg	
Weight / Dimensions, height x width x depth in mm		380 x 315 x 200	398 x 278 x 210	

Supplied complete with 1000 capillary tubes of 1.4 mm Ø ext set of maintenance tools comprising of star and flat screwdrivers, drill bit of 1.5 mm Ø, allen key 3 mm, lamp and fuse, software and protective cover.

ACCESSORY Capillary tubes of 1.4 mm Ø ext. 1.0 mm Ø int. Box of 5000 tubes Part No. 5609404


SPARE Lamp 6.5 V 0.3 A. Part No. **5609405**



ROTARY VISCOMETERS



FLOW CUPS



GLASS VISCOMETERS



MONITOR AND CONTROL



GLASS THERMOMETERS



REACTORS



LABORATORY FUME CUPBOARDS

Laboratory fume cupboards Vertical laminar flow cabinets Rotary viscometers Flow cups for measuring liquid viscosity **Viscometer precision bath Glass viscometers Monitor and control apparatus Glass thermometers Laboratory reactors**

Filters

Analytical consumables

page 284 page 285

pages 286 to 288 pages 289 to 290

page 290

page 291 to 292

pages 293 to 297

pages 296 to 297

pages 298 to 299

pages 300 to 301

pages 302 to 311



Laboratory fume cupboards "VG-BS"

FOR STEAM, ACIDS, BASES AND SOLVENTS EXTRACTING.



SAFETY: DESIGNED AND BUILT MEETING EUROPEAN STANDARD UNE-EN 14175-2

COMMON FEATURES

Antiacid continuous ceramic worktop with anti-dumping peripheral rim and evacuation

Inner working chamber built in thermo stable compact resins, antiacids and anti-moisture. M1 fire resistance.

Extraction system with perfectly balanced injected polypropylene turbine. Coupled and protected engine outside the fumes way.

Vertical guillotine door with adjustable height in any position and three sections of 3 + 3mm laminated glass slide which creates a high resistance protection shield.

Self-supporting structure, made of structural elements with 2 mm steel tube, degreased, phosphated and coated with a 70 μ layer of epoxy paint vitrified in the oven. This structure is completely isolated from the working cabin and from the gases contact.

External cover in white melamine board or thermo stable resins in high humidity environments.

ASSEMBLY SYSTEM INTO TWO BODIES

Upper body with cabinet and embedded extractor.

Lower body with a structure, ceramic worktop and remote control panel for fluids and gases in an horizontal housing.

Optional: working system (Walking) at ground level for the input and working with machinery, reactors and organic synthesis appliances. (Request information).

TECHNICAL FEATURES

Lighting level: More than 500 Lux in working area IP 55 or ATEX*.

Extraction volume: from 420 to 4.500 m3/h. Engines: from 0,75 to 2 CV - IP55 or ATEX*. Voltage: 220V/50Hz. Single phase / three phase.

*ATEX: if we work with flammable, solvents or hydrocarbons.

On request: Automatic control unit for variable flow (motor) depending on door's opening.

Failure alarm.

On request: Perimeter and ceiling in continuous tiles, for working with sulfuric acid. On request: Perimeter and ceiling in continuous tiles, for working with perchloric acid. It comes with gas washing.

CONTROL PANEL

ON / OFF push button with green LED and electrical sockets 230V/16A. Fuel regulator knob (Yellow).

Nitrogen regulator knob (Blue/ green button).

Vacuum regulator knob (Grey/ black button).

Air regulator knob (Blue/ yellow button).

Water regulator knob (Green/ blue button)

Safety locking pin on door for recommended working position.



Can be supplied with custom measurements on demand.

MODELS	Part No.	Height / Weight / Depth Upper body (exterior) cm	Height / Weight / Depth Lower body (exterior) cm	Height / Weight / Depth Upper body (inner) cm	Height / Weight / Depth Lower body (inner) cm	Upper Weight Kg	Lower Weight Kg
VG-BS-1200	5001613	170 122 90	90 122 90	120 96 62	72 96 55	124	84
VG-BS-1420	5001614	170 142 90	90 142 90	120 116 62	72 116 55	162	116
VG-BS-1720	5001615	170 172 90	90 172 90	120 146 62	72 146 55	200	148
VG-BS-2002	5001616	170 202 90	90 202 90	120 176 62	72 176 55	238	180

ACCESSORIES PVC Tubes			PVC Elbow joint		Safety cabir	iets		
Part No.	Туре	Ø / Length cm	Part No.	Ø / Degrees	Part No.	Acids & base	Flammables	Doors
5001617	Rigide	20 200	5001621	20 45°	5001625	Yes	No	1
5001618	Rigide	25 200	5001622	25 45⁰	5001626	Yes	No	2
5001619	Flexible	20 100	5001623	20 90°	5001627	No	Yes	1
5001620	Flexible	25 100	5001624	25 90⁰	5001628	No	Yes	2



Vertical laminar flow cabinets "H-1" and "H-2"

TFT TOUCH SCREEN CLASS 1



SAFETY:

US STANDARD ST 209 (CLASS 100)

EUROPEAN STANDARD ISO 14644 (CLASS 5) HEPA FILTER (HIGH EFFICIENCY PARTICULATE AIR)

APPLICATIONS

Water quality control, laboratories and food industries, pharmacies services, haematology and clinical analysis, drug-injecting antibiotics, microscopic analysis, assisted fertilization, cell cultures, etc...

FEATURES

Lower setting made of PMMA transparent methacrylate.

Working surface made of stainless steel AISI 304.

Ultraviolet lamp.

Metallic parts coated with protective polymerization.

Extraction system with silent turbine.

Operating hours counter for HEPA filter replacement.

TECHNICAL SPECIFICATIONS

Illumination level (W-Lux): 500. Airflow speed (M/S): 0.7 - 1.0.

Noise level (dB): < 67.

Power supply (V-Hz): 220/50.

HEPA filter class H14: 99,999% efficiency for particles of 0.5 microns. Working

temperature in °C: between 15 and 35.



TFT touch screen

Stop/start push button.

Light push button.

Sterilization push button.

Ventilation push button.

Operating hours counter push button.



MODELS	Part No.	Height / Width / Depth (exterior) cm	Height / Width / Depth (interior) cm	Fluorescent Lamp (W)	UV lamp W	Power W	Weight Kg
H1	5609511	55 46 78	54 38 40	15x1	15x1	150	40
H2	5609510	70 55 90	68 45 50	18x1	18x1	160	60

ACCESSORIES	Supporting		
FOR MODEL	Part No.	Height / Width / Depth (exterior) cm	Weight Kg
H1	5609512	78 55 47	5
H2	5609513	78 70 55	8





Shoe soles coater "ZA-1"

AUTOMATIC.

DIGITAL TEMPERATURE CONTROL OF AIR.

FEATURES

Shoe soles coated with a PVC film.

It cuts the film by thermal contraction providing hot air with temperature control. Different sizes shoes coated in seconds.

It saves on consumables. A roll of film can be coated up to 1000 shoes.

According with CE and ROHS regulations. Consumables 100% recyclable.

MODEL

Part No.		mension : / Width	s cm n / Depth	Voltage	Power W	Weight Kg
5000051	33	45	80	220V-50Hz	1300	27





SPARE PARTS

PVC coil for 1000 coatings. Part No. **5000052**





Rotary Viscometers

Rotary viscosimeters "ST-2020"

INTRODUCTION

Rheology is the study of the effects experimented in a substance when a mechanical force is applied on a it (flow and deformation) under different external conditions. It is used to describe the consistency of different products and is normally defined by the components: viscosity and elasticity.

Measuring viscosity is determined by the tangible force required to displace the materials particles with a specific deformation-flow i.e. velocity. The relationship between the tangible force and the deformation flow obtains the viscosity result. Ambient conditions such as temperature and pressure also have an effect on viscosity. The measurement of viscosity is not just limited to the research laboratory, it has progressively entered the field of industrial quality control.

PRINCIPLES OF VISCOSITY

These instruments operate by means of a cylinder or disk (spindle) that is submerged into the material to be analysed and by measuring the resistance of the substance at a selected known speed. This resistance results is the measurement of the viscosity according to the flow characteristics of the reference spindle; the instrument calculates the result and directly displays the viscosity that is reported in **cP** (**CGS**) or **mPa-s** (**SI**).

A wide range of viscosity can be measured using viscometers that are equipped with different types of spindles and speed ranges. The design of the spindles and the principals of measurement principles are regulated by **ISO 2555** and **ISO 1652** standards. All spindles are made of AISI 316 stainless steel. Each spindle can be identified by a letter and a number.

SELECTION TABLE

Standard measuring range of the viscometers, without additional accessories						
Part no	1001616	1001617				
Model	ST-2020 L	ST-2020 R				
Units	centiPoise (cP)	centiPoise (cP)				
Standard spindle	L1 to L4	R2 to R7				
Speed range r.p.m.	1 to 60	0,1 to 100				
Measuring range	20 to 600.000 c P	20 to 40.000.000 cP				
Temperature range °C	0,0 a 100,0	0,0 a 100,0				
Power requirement	115/230V to 12VDC 1.2A	115/230V to 12VDC 1.2A				
Power	15 W	15 W				
Weight	5 K g	5 Kg				

FEATURES

L.C.D. display of parameters and results:

- Selected speedr.p.m.
- Selected spindle
- Viscosity result
- Base scale percentage%.
- Sample temperature: °C or ° F.
Auto alarm in the case of any fault being detected.
Off scale detection and indication by an audible
and visual signal.

Step controlled speed to prevent spindle vibrations.

Velocity from 0,1 to 100 r.p.m.

Mains power surge protection.

RS 232 unidirectional interface, download to a computer.

TECHNICAL DATA

Precision: $\pm 1\%$ base scale.

Repeatability: 0.2%.

Supplied complete with:

- Anti shock carry case.
- Main unit.
- Support base.
- Spindle protector.
- Spindle support.
- Set of spindles (model dependant)
- Temperature prove

DIGITAL THERMOMETER

Temp range:- from 0.0 °C to + 100.0 °C (+32.0 °F to + 212.0 °F).

- Resolution: 0.1 °C (0.1722 °F).
- Precision: ± 0.1 °C.

Model ST2020L

ACCESSORIES

Standard spindles for L model.

L1 Part No. 1000998

L2 Part No. 1000999

L3 Part No. 1001000

L4 Part No. 1001001

Standard **spindles** for R model.

R1 Part No. 1000995 (for low viscosity samples)

R2 Part No. 1001030

R3 Part No. **100 1031**

R4 Part No. 1001032

R5 Part No. 1001033 R6 Part No. 1001034

R7 Part No. 1001035



L3,L4 : Suitable for models L.





Rack of standard spindles R2, R3, R4, R5, R6 and R7: Suitable for model R.



Wide range rotary viscosimeters "STS-2011"

MODEL WITH TEMPERATURE READOUT AND SHEAR RATE MEASUREMENT.

FEATURES

10 different options of language.

L.C.D. display of parameters and results.

Parameters display:

- Selected speedr.p.m.

- Selected spindle

- Viscosity result cP (mPa-s) or cSt.

- Base scale percentage%.

- Sample temperature: °C or ° F.

- Deformation ratio (with special spindle)

- Torsion force (with special spindles)

This instrument determines both relative and absolute viscosity.

Data can be changed between S.I. and C.G.S.

Automatically checks for correct operation point by scanning at different speeds.

Auto alarm in the case of any fault being detected.

Off scale detection and indication by an audible and visual signal.

Step controlled speed to prevent spindle vibrations.

Calibration by the user himself.

18 preselected speeds from 0.3 a 100 r.p.m.

The operator can select any speed within this range. USB port.

Mains power surge protection. 10 memoirs of working programs.

TECHNICAL DATA

Temperature range:- from 0.0 °C to + 100.0 °C

(+ 32.0 °F to + 212.0 °F).

- Resolution °C: 0.1 °C (0.1722 °F).

- Precision °C: ± 0.1 °C.

- Probe type: Pt 100.

Direct results in cP(mPa-s) or cSt.: models STS-2011 L & R. Direct results in Poise (Pa-s) or St.: model STS-2011 H.

Precision: ±1% base scale. Supplied complete with:

- Anti shock carry case.

- Main unit.

- Support base.

- Spindle protector.

- Spindle support.

- Set of spindles (model dependant) (see table).



I	MODELS	Part No.	Measuring	Standard	Power	Power	Weight
			range	spindles	requirement	W	Kg
9	STS-2011 L 1	001611	20 to 2.000.000 cP	L1, L2, L3, L4	100-240 V 50/60 Hz	15	5
9	STS-2011 R 1	001612	100 to 13.000.000 cP	R2, R3, R4, R5, R6, R7	100-240 V 50/60 Hz	15	5
	STS-2011 H 1	001613	200 to 106.000.000 cP	R2. R3. R4. R5. R6. R7	100-240 V 50/60 Hz	15	5

ACCESSORIES

Standard **spindles** for L model.

L1 Part No. 1001036

L2 Part No. 1001037

L3 Part No. 1001038

L4 Part No. 1001039

Standard spindles for R and H model.

R1 Part No. 1000990 (for low viscosity samples)

R2 Part No. 1001070

R3 Part No. 1001071

R4 Part No. 1001072

R5 Part No. 1001073

R6 Part No. 1001074

R7 Part No. 1001075

Small sample volume adapters (APM).

Suitable for sample volumes from 8 to 13 ml.

Requires the "TL" or "TR" set of spindles.

Unit suitable for temperatures from 0 °C to +100 °C.

Part No. 1000986 Small sample adapter APM with water jacket.

Part No. 1001623 Temperature sensor suitable for APM.

Special Spindles

Part No. 1001224 set of spindles TL5 - TL6- TL7 for models L.

Part No. 1001225 set of spindles TR8 - TR9 - TR10 - TR11 for models R and H.

Low viscosity adapters. (LCP).

Required if low viscosity measurements are necessary.

Reproducible results and measurements of viscosity from 1 cP.

Suitable for models L and R. Sample volume: 16 to 18 ml. Includes special spindle LCP.

Part No. 1000985. Adapter LCP with water jacket. Part No. 1001624 Temperature sensor suitable for LCP.

Displacement spiral helix adapter.

Required for low fluidity samples.

Part No. 1000988



Part No. 1000985

Spindle R1

Part No. 1000986

Part No. 1000988





Temperature control equipment for viscosity measurements for viscome-

ters "ST-2020" and "STS-2011"

The influence of temperature while measuring viscosity is considerable, lower temperatures increase viscosity. Therefore it is indispensable to control temperature when precise viscosity measurements are required. The control of temperature by using a thermostatically controlled bath is the most efficient, because of the recirculation of liquid produces a rapid and stable temperature.

We recommend our range of immersion thermostats for this application.

TEMPERATURE CONTROL AND CONFIGURATION FOR VISCOSITY MEASUREMENTS:

FRIGITERM-TFT-10 Part No. 3000546, FRIGTERM-TFT-30 Part No. 3000547 Suitable for below ambient working environments (see page 103 for specifications).

For temperatures above ambient see the **DIGITERM-S-150** Part No. **3000543**, or **DIGITEM-TFT-200** Part No. **3000544** complete with 12 litre tank Part No **6000391**, (see pages 101 and 102 for specifications).

An adapter kit for the thermostat bath is required:

Part No. 1001625 for STS-2011 Part No. 1001626 for ST-2020

Adapter for thermostat bath comprising of an extension spindle and 4 leg adjusters for the bath. (Kit can be fitted by the user.)

This kit adjusts the bath measuring height to enable samples to be placed inside.

Part No. 1001627. Support for 2 500 ml beakers, to be placed inside the bath.

Part No. 1001628 Support base for the **DIGITERM** immersion thermostats. This accessory is recommended to close the bath and maintain a constant temperature.

CONSTANT TEMPERATURE CONTROL ACCESSORY FOR SMALL SAMPLE VOLUME ADAPTERS. USE WITH THE RE-CIRCULATING JACKET ACCESSORIES. PART NO. 1000996

To work below ambient temperatures, we recommend the use of the **FRIGITERM-TFT-10** Part No. **3000546**, or **FRIGITERM-TFT-30** Part No. **3000547** (see page 103 for specifications.)

For temperatures above ambient see the **DIGITERM-S-150** Part No. **3000543**, or **DIGITEM-TFT-200** Part No. **3000544** complete with 12 litre tank Part No. **6000391**, (see pages 101 and 102 for specifications).

It is necessary to configure the immersion thermostats for "external re-circulation." The Bath adapter kit is not required.





Digiterm thermostat bath with support base Part No. 1001628. 12 litre bath Part No. 6000391 and adapter kit Part No. 1001625 or 1001626 To be used with a rotary viscometer and beaker support Part No.



Immersion thermostat Digiterm controlling temperature of viscometer.



COMECTA Rotary Viscometer "NDJ-1"

APPLICATIONS

The instrument operates by rotating a disk or cylinder, (spindle), that is submerged in the liquid or fluid to be analysed. A pre-selected speed is set, the unit measures the absolute resistance from viscosity of the fluid being analysed. Suitable for samples such as: foods, cosmetics, fats and oils, pharmaceutical products, paints and plastics, etc.

FEATURES

The asynchronous motor is connected to a graduated disk with 4 different speeds that propel the spindle via a spiral and die.

Supplied complete with a set of 4 spindles in a box, numbered 1 to 4 with spindle stand. The viscometer includes a level and adjustable screw feet support base and protective case.

NDJ-1	5120230	10-100.000 mPaS	±5% Liquid Newtons	6 - 12 - 30 - 60	15	6
MODEL	Part No.	Measuring range	Tolerance	Spindle Speed r.p.m. for 1 to 4	Power W	Weight Kg





Flow cups for measuring liquid viscosity

CUP N° 4 STANDARD DIN 53211. CUPS STANDARD ISO 2431. CUP FORD STANDARD ASTM D-1200. CAPACITY: 100ML.

FEATURES

Suitable for measuring kinematic viscosities from 5 to 700 cSt, model dependent. Metallic cups, chrome finished and calibrated.



Flow cups with handle. Models DIN 53211 N° 4 and Ford ASTM D-1200.



Flow cups. Standard models.

ACCESSORY

Support stand with adjustable level. Part No. 7001021



Part No.	Standard	Bore Ø mm	Format	Range cSt	Admisable fall times
1000123	DIN 53211	4	Standard	90 to 700	25" to 100"
7001239	DIN 53211	4	Heated	90 to 700	25" to 100"
1000347	DIN 53211	4	With handle	90 to 700	25" to 100"
1001013	ISO 2431	3	Standard	5 to 42	30" to 100"
7001017	ISO 2431	3	Heated	5 to 42	30" to 100"
1001014	ISO 2431	4	Standard	35 to 135	30" to 100"
7001018	ISO 2431	4	Heated	35 to 135	30" to 100"
1001015	ISO 2431	5	Standard	100 to 350	30" to 100"
7001019	ISO 2431	5	Heated	100 to 350	30" to 100"
1001016	ISO 2431	6	Standard	190 to 680	30" to 100"
7001020	ISO 2431	6	Heated	190 to 680	30" to 100"
		2.53		25 to 120	20" to 100"
1000705	ASTM D-1200	3.40	Standard	40 to 220	20" to 100"
		4.12		70 to 370	20" to 100"
		2.53		25 to 120	20" to 100"
7000706	ASTM D-1200	3.40	Heated	40 to 220	20" to 100"
		4.12		70 to 370	20" to 100"
		2.53		25 to 120	20" to 100"
1000707	ASTM D-1200	3.40	With handle	40 to 220	20" to 100"
		4.12		70 to 370	20" to 100"




Heated cups with threaded base, can be connected to a water bath or to a temperature regulator Electemp.

ACCESSORIES

Water bath with heater and screw on disk, with level. Part No. 7001022 Temperature controller Electemp-TFT. Part No. 3000897 Pt 100 sensor probe for the Electemp-TFT and water bath. Part No. 7001496 (See page 293).

Comes with a calibration certificate with results traceable to know standards and equipment





Ford cup thermostat for viscosity measurement "TV-1452"

CUP MODELS THAT CONFORM TO THE FOLLOWING STANDARDS: DIN 53211, ISO 2431 AND ASTM 1200. ELECTRONIC DIGITAL CONTROL FROM 10°C TO 60°C.

Heating and cooling by Peltier effect

APPLICATIONS

Thermostat Ford cups for viscosity measurement need to be maintained at a precise temperature of 10 and 60 $^{\circ}\text{C}.$

FEATURES

Made up of an independent Peltier thermostat control system that heats and cools. Made from AISI 304 stainless steel, with adjustable levelling feet supports, and central orifice for locating the cup.

CONTROL PANEL

Illuminate ON/OFF switch
Digital temperature display
Key pad to select readout and temperature.



Part No.	Temperature °C	Stability °C	Homogeneity °C	Read error	Resolution	Height	/ Width	/ Depth	Power W	Weight Kg
3001452	10 to 60 °C	±0.1	±0.2	±0.5	0.1	23	34	30	130	5

See the different standards, models and cups (see page 289.)



Viscometer Precision Bath "VB-1423"

FOR CONTROLLABLE TEMPERATURES FROM AMB.+5 °C TO 100 °C.
STABILITY ±0.05 °C. HOMOGENEITY ±0.05 °C. RESOLUTION 0.1 °C.

SAFETY:

STANDARD DIN 12879.2 CONTROLLABLE SAFETY THERMOSTAT WITH MANUAL RESET.

FEATURES

Temperature sensor Pt100 thermo-resistor, stainless steel AISI 304 lid with three viscometer locations ports, three independent lids and an additional location port for the control thermometer. The main body of the bath is made of a 20 litre borosilicate glass tank. A white plate is located at the back to help optimise and read the viscometers.

CONTROL PANEL

- 1. Mains power iluminated switch.
- 2. Temperature regulator:
 - 3. Real time temperature display.
 - 4. Push button to increase value.
 - **5**. Push button to decrease value.
 - **6**. Push button to configure operation.
- 7. Safety thermostat safety lamp.

Used for measurements with glass viscometers

Made for the calibration of viscometers according to the following standards UNE 400313, ISO 3105, ASTMD 445 and 2515





MODEL

MODEL						
Part No.	Control range °C	Capacity litres	Height / Ø (tank) cm	Height / Ø (total) cm	Power W	Weight Kg
3001423	amb.+5 up to 100	20	32 30	47 30	1000	8

ACCESSORIES

Universal viscometer support made from PTFE with stainless steel AISI 304 support. Suitable for the following viscometers:



- Cannon-Fenske for transparent liquids.
- Cannon-Fenske for opaque liquids.
- Ubbelohde.
- Ostwald.
- BS U Tube.

-Cannon-Manning semi-micro.

- Ubbelohde type BS/IP/SL, BS/IP/SL(S) & type BS/IP/MSL.
- DIN Ubbelohde.

Part No. 1001453

Calibration Chronometers (see page 295).

Thermometers for viscometer baths.

Part No.

 1001454
 Thermometer ASTM 120C at 38.6 to
 41.4°C divisions 0.05 °C.

 1001455
 Thermometer ASTM 121C at 98.6 to 101.4 °C divisions 0.05 °C.

 1001456
 Thermometer ASTM 91C at 20.0 to 50.0 °C divisions 0.1 °C.

 1001457
 Thermometer ASTM 92C at 40.0 to 70.0 °C divisions 0.1 °C.

 1001458
 Thermometer ASTM 93C at 60.0 to 90.0 °C divisions 0.1 °C.

1001459 Thermometer ASTM 94C at 80.0 to 110.0 °C divisions 0.1 °C.





VISCOMETER "UBBELOHDE"

ASTM D445 - ASTM D446 - ISO 3104 - ISO 3105. Suitable for transparent liquids. Complete with calibration certificate. Total length 283 mm. Permanent amber markings.

ACCESSORY

Viscometers support rack. Part No. 1025812 Capacity: 6 viscometers.





Ubbelohde

Part No.	Size	Nominal constant	Viscosity	y ran	ge in cSt
5600001	0	0.001	0.3	to	1
5600002	OC	0.003	0.6	to	3
5600003	OB	0.005	1	to	5
5600004	1	0.01	2	to	10
5600005	1 C	0.03	6	to	30
5600006	1B	0.05	10	to	50
5600007	2	0.1	20	to	100
5600008	2 C	0.3	60	to	300
5600009	2 B	0.5	100	to	500
5600010	3	1.0	200	to	1000
5600011	3 C	3.0	600	to	3000
5600012	3 B	5.0	1000	to	5000
5600013	4	10.0	2000	to	10000
5600014	4C	30.0	6000	to	30000

VISCOMETERS "CANNON-FENSKE" TRANSPARENT

ASTM D445 - ASTM D446 - ISO 3104 - ISO 3105 - IP 71. Suitable for transparent liquids. Complete with calibration certificate. Constant at 40 °C to 100 °C. Total length 250 mm. Permanent amber markings.



Transparent

Part No.	Size	Nominal constant	Viscosi	ty ra	nge in cSt
5600050	25	0,002	0,4	a	1,6
5600051	50	0,004	0,8	a	3,2
5600052	75	0,008	1,6	a	6,4
5600053	100	0,015	3	a	15
5600054	150	0,035	7	a	35
5600055	200	0,1	20	a	100
5600056	300	0,25	50	a	200
5600057	350	0,5	100	a	500
5600058	400	1,2	240	a	1200
5600059	450	2,5	500	a	2500
5600060	500	8	1600	a	8000
5600061	600	20	4000	a	20000

VISCOMETERS "CANNON-FENSKE" OPAQUE

ASTM D445 - ASTM D446 - ISO 3104 - ISO 3105. Suitable for opaque and transparent liquids. Complete with calibration certificate. Constant at 40 °C and 100 °C. Total length 295 mm. Permanent amber markings.



Opaque

Part No.	Size	Nominal constant	Viscosi	ty ra	nge in cSt
5600065	25	0,002	0,4	a	1,6
5600066	50	0,004	0,8	a	3,2
5600067	75	0,008	1,6	a	6,4
5600068	100	0,015	3	a	15
5600069	150	0,035	7	a	35
5600070	200	0,1	20	a	100
5600071	300	0,25	50	a	200
5600072	350	0,5	100	a	500
5600073	400	1,2	240	a	1200
5600074	450	2,5	500	a	2500
5600075	500	8	1600	a	8000
5600076	600	20	4000	а	20000

Viscosity range in cSt

to

to

to

to

3

10

30

100

300

1000

3000

to 10000

0,9 to

2

6 to

20 to

60

200

600

2000



COMECTA Glass Viscometers

MICRO-UBBELOHDE VISCOMETER 3 ml

Suitable for transparent liquids. Calibration certificate. Nominal overall length 290 mm. Permanent ring marks in amber. Sample filling volume: 3 ml.



Micro-U	bbei	ohde
---------	------	------

Code	Size	Nominal constant	Viscosit	y rang	e in cSt
5600020	1	0,01	0,4	to	6
5600021	lc	0,03	1,2	to	18
5600022	II	0,1	4	to	60
5600023	llc	0,3	12	to	180
5600024	III	1,0	40	to	800

"U-TUBE REVERSE FLOW" VISCOMETER BS/IP/RF

Suitable for opaque liquids. Calibration certificate. Nominal overall length 275 mm. Permanent ring marks in amber.



U-TUbe reverse flow

Code	Size	Nominal constant	Viscosit	y rang	e in cSt
5600025	1	0,003	0,6	to	3
5600026	2	0,01	2	to	10
5600027	3	0,03	6	to	30
5600028	4	0,1	20	to	100
5600029	5	0,3	60	to	300
5600030	6	1,0	200	to	1000
5600031	7	3,0	600	to	3000
5600032	8	10	2000	to	10000
5600033	9	30	6000	to	30000
5600034	10	100	20000	to	100000
5600035	11	300	60000	to:	300000

Nominal constant

0,003

0,01

0,03

0,1

0,3

1,0

3,0

10,0



Code

5600036

5600037

5600038

5600039

5600040

5600041

5600042

5600043

Size

A

В

C

D

Ε

G

Н



Ostwald

BS/U TUBE VISCOMETER Suitable for transparent liquids. Nominal overall length 300 mm. Permanent ring marks in amber.

Code	Size	Nominal constant	Efflux time between ring marks Water at 20 °C. (seconds)
5600044	45	0,022	45
5600045	85	0,011	85
5600046	100	0,01	100
5600047	125	0,008	125

OSTWALD PATTERN VISCOMETER

Suitable for transparent liquids. Calibration certificate at 20°C. Nominal overall length 250 mm. Permanent ring marks in amber. Sample filling volume 10 ml.

Calibration certificate.

Monitor and Control Apparatus



Electronic controller "Electemp-TFT"



FEATURES

Temperature range: -150 °C to 1300 °C

Stability: ±0.1 °C.

Resolution: ±0.1 °C Pt 100 /±1 °C K Type probe.

Touch screen

- Main On/Off switch.
- Acoustic and visual alarm.
- Clock calendar.
- Up to 10 work programs.
- Up to 6 segments per program.
- Stability time in each segment (from 1 min to 99h).
- Ramps between segments.
- Alarms and events storage.
- Probe error detection.
- Self Diagnostics .
- Network failure detection and saving.
- Over temperature and low temperature alarms and memorization (date, start time, end time and temperature).
- Safety thermostat (TS) by software .
- RS-232 and USB output.
- Software for PC.
- Temperature control auto-tuning.
- Manual user on screen.
- Configurable parameters: Date / time, temperature correction , data collection interval, language (English, Spanish and French) , °C / °F selection , over temperature and low temperature limit.

CONTROL PANEL

Setpoint temperature and time indicator.

Temperature limit indicator.

Timer virtual pushbutton.

Reverse virtual pushbutton.

Start virtual pushbutton.

Stop virtual pushbutton.

Programming virtual pushbutton.

BACK PANEL

Power output (SCHUKO plug) 3100 W at 230 V.

Protection fuse.

Pt100 probe connection.

K type probe connection.

Data output.

APPLICATIONS

Digital thermometer.

To control and display temperature in any equipment with temperature control system, that uses a contact thermometer, such as: Baths, Immersion thermostats, Hotplates and Hotplate stirrers, Mantles, etc.

Can also be used by a immersion thermostat controlling a heater element. (maximum power 3100 W).

As a regulator or temperature reader in equipment that are labelled with .



MODEL

Part No.		Width / (exterio		Max. power W (at 230 V)	Weight Kg
3000897	10.5	13.5	21	3100	1.250

ACCESSORIES

Pt 100 probe. Immersion dimensions: 4 mm \emptyset , 130 mm length. Cable: 150 cm length.

Part No. 1000893

Type K probe. Immersion dimensions: $2 \text{ mm } \emptyset$, 150 mm length. Cable: 150 cm length.

Part No. 1000895



Electronic contact thermometer "Sensoterm II"

DIGITAL DISPLAY, PT 1000 PROBE, AISI 304 STAINLESS STEEL. STABILITY ± 1 °C, ± 1 Digit. Resolution ± 1 °C.

ONLY SUITABLE FOR INSTRUMENTS THAT CAN USE AN ELECTRICAL CONTACT THERMOMETER.

FEATURES

Advantages include a separate temperature probe that maintain the main body away from steam and other vapours. Safer than a glass contact thermometer, eliminating the risk of broken it and a mercury spillage. Main body made of ABS plastic.

TECHNICAL DATA

Temperature range: from -50 °C to 300 °C.
Temperature probe: 3.5 mm Ø x 250 mm immersion.
Temperature Display: 8 mm high LCD display.
Ambient working range: 0 to 60 °C and a relative humidity of 80 %.



Part No.	•	eight / Width / Depth cm (exterior)		Weight Kg
1001009	13	7	2.4	0.250





Instruments labelled with labelled with that can use the electronic controller "Electemp-TFT" and the electronic contact thermometer "Sensoterm".



Sand Bath "Combiplac-Sand" Parts No. 6000708 and 6000709.



Hotplate - Rectangular "Combiplac" Parts No. 3000156 and 3000718.



Hotplate - Circular "Combiplac" Parts No. 3002301 and 3002302.



Magnetic Stirrers "Agimatic" Parts No. 7000243, 7002431, 7000446, 7001511, 7001609, 7000837 and 7000838.

NOTE: Any equipment that can be connected to an electrical contact thermometer can also be used with an "Electemp-TFT" or "Sensoterm" unit.



Electronic energy controller "Enertron"

FEATURES

Regulates the power to the heating element by pulsing electronic circuit, free of noise and interferences

CONTROL PANEL

Main ON/OFF switch. Mains "ON" indicator lamp. "ON" indicator power out. Controllable from 4% to 100 %.

BACK OF CONTROLLER

Output of power (Schuko plug system).

Protection fuse.

Connection for a contact thermometer.

APPLICATIONS

Regulates the power of a heating element in ovens baths, radiators, hotplates, etc.

MODEL

Part No.		Width / (exterior		Max. power W (at 230 V)	Weight Kg
3000886	10.5	13.5	21	2500	1.2



"Enertron" application to a hotplate, showing the unit connected with an electronic contact thermometer.

COMECTA Digital thermometer "TC-9226-A"

FEATURES

Main body made of PVC with anti-slide for better stability.

Temperature range from -200 °C to 1370 °C (-328 to 2498 °F).

Reads °C or °F, reading within 1 second.

Memorizes the last 10 temperature readings. Displays the reading for 3 seconds. It is posibles to hold the reading on the screen.

Switches off automatically after 1 hour (if it is in a continuous mode).

LCD with back light illumination.

Low battery power indicator.

Powered by 3 x 1.5V "AAA" batteries".

Compatible with most Type K sensors.

Operating conditions: from -10 to 50 °C.

Dimensions: (height x width x depth): 17 x 7 x 3 cm. Weight 150g.

MODEL	Part No.	Temp. range °C	Resolution °C	Precision °C
CT-9226-A	5903030	-200 to 1370	0.1	from -200 to -50: ±2 from -50 to 740: ±1 from 740 to 1370: ±2

The precision figures supplied are based on a measurement range from 18 to 28 °C (does not include the sensor error).



ACCESSORY

Sensor type K of 8 mm Ø x 30 cm length. Part No. 5903025

COMECTA Digital Humidity meter with Thermometer "HT-305E"

protective cover and carry case.

Supplied complete with humidity / temperature sensor,

Part No. 5903023

FEATURES

Temperature range from 0 °C to 60 °C. Humidity range from 10% to 95% RH. Digital Display °C or in °F Fast response: 0.4".

Power: battery 9 V.

Dimensions: (height x width x depth): 3.2 x 7 x 15.5 cm.

Weight: 370 g.

MODEL	Part No.	Temp. range °C	Resolution °C	Precision °C	Humidity range %	Humidity resolution %	Humidity precision %
HT-305E	5903023	from 0 to 60	0.1	±0.8	from 10 to 95	0.1	> 70%: ±3% +1% < 70%: +3%







Ink Printer "Microprinter"

APPLICATIONS

Small ribbon printer (not thermal) with paper roll 2 1/4" (56 mm) wide.

Interface connection: RS232.

Applications: Measuring equipment, pH meters, conductivity meters, spectrophotometers, etc.

Includes interface cable and power transformer 5V - 3A and cable.



Part No.	Height	/ Width / Depth cm		Power supply	Weight Kg
4120113	4	16	10	DC 5V 3A	1



Automatic programmer "Promat" 24 H

Connects and disconnects electrical power by a pre-programmed timer switch.

Adjustable timer: Cycles up to 24 hours, with controllable connection disconnection up to every 15 minutes. Manual interrupt setting to disconnect the program, luminous indicator shows when power is being applied.

Maximum current: 16 Amp.

Part No. 1000500







COMECTA Timer alarm "Relavis"

Audible sound when the program time has lapsed. Time setting from 0 to 60 minutes in 1 minute intervals.

Part No. 5900499





COMECTA Digital alarm timer "Count Down-Up"

Maximum count time 99 min. 59 sec. Timer: count-down and count-up. Portable: table top with magnetic holder or clip if used on a pocket or overall.

Audible alarm fitted. Part No. 5900500





COMECTA Digital Chronometer "Diamond"

Maximum timer: 23 hr 59 min 59 sec. Division 1/100 seconds. Control: LAP/SPLIT (timer: split time and accumulative time). Date incorporated.

Audible Alarm (duration 60 seconds).

Audible Alarm "snooze" function (once pressed the stop alarm goes off every 5 minutes until it is completely shut off).

Part No. 5900501





COMECTA Thermometer reader "Micro-Temp"

SUITABLE FOR MEASURING TEMPERATURE FROM -50 °C TO +300 °C OR FROM -58 °F TO +527 °F.

Applications: process control, laboratories, pharmacies, food processes and agriculture in general. Digital screen readout in °C or °F.

Stainless steel penetration probe of 3.5 mm Ø x 120 mm long. Retention memory of minimum and maximum data readout.

Disconnects automatically after 1 hour.

Powered by a 5V - 3A battery.



MUDEL		
Part No.	Precision	Resolution
5903029	±1 °C / 1 °F	0.1 °C / 0.1 °F

Comes complete with protective cover and support clip.





COMMON FEATURES

All thermometers are made of glass rod or with an opal scale, with a reinforced bulb. A wide selection of scales and temperature ranges is available.

Scales etched in acid and graduated with indelible ink.

Resistant to most chemical agents.

All thermometers are supplied with a protective plastic cover.






Solid stem Themometers

Yellow backed, mercury filled.

Part No.	Range °C	Division scale °C	Total length cm
1090100	-10 +100	1/1	30
1090150	-10 +150	1/1	30
1090250	-10 +250	1/1	30
1090360	-10 +360	1/1	30
1093100	-10 +100	0.5	30
1092050*	0 +50	0.1	45

* White back.



-10 +100

1/1

cm

30

Yellow backed with suspension ring, mercury filled.

1001257

Part No.	Range °C	Division scale °C	Immersion length cm	Total length cm
1001299	0 +200	1/1	5	25
3009100	-20 +100	0.5	15	34
3009200	0 +200	1/1	15	34
1009100	0 +60	0.5	15	30



Enclosed scale thermometers

With prismatic capillary, with easy to read scale. Contains mercury.

	ζ	7		
	10 The state of th	2 10 18		
		2 4		
,			,	

·	,	•		
Part No.	Range °C	Division scale °C	Immersion length cm	Total length cm
1095060	-10 +50	1/1	3	30
1095100	-10 +100	1/1	3	30
1095150	-10 +150	1/1	3	30
1095250	-10 +250	1/1	3	30
1095360	-10 +360	1/1	3	30
1096050	-10 +50	0.5	3	30
1096100	-10 +100	0.5	3	30
1097050	-10 +50	0.1	3	45
1000120	-10 +120	1/1	10	25
1000100	-10 +100	1/1	13	25
1000872	-10 +120	1/1	15	33
1000873	0 +220	1/1	15	33



Precision thermometers

Mercury thermometers.

Complete with **calibration certificate** PROTON registration control no. 02-C.37 and 02-V.06 (Meteorological control).

Traceable certificate indicating the proof of continued quality of work-manship against different temperatures, with method and conditions of measurement against known traceable standards.



Part No.	Range °C	Division scale °C	Туре	Total length cm
1001258	-10 +100	1/1	solid stem	30
1001259	-10 +250	1/1	solid stem	30
1001260	-10 +50	0.5	solid stem	30
1001261	-10 +100	0.5	solid stem	30
1001262	0 +50	0.1	solid stem	45
1001263	-10 +100	1/1	enclosed scale	30
1001264	-10 +250	1/1	enclosed scale	30
1001265	-10 +50	0.5	enclosed scale	30
1001266	-10 +100	0.5	enclosed scale	30
1001267	-10 +50	0.1	enclosed scale	45

Case with 7 glass rod mercury thermometers, white back, scale in 0.1 °C.

Range °C*	Total length
-30 +20	45 cm
0 +50	45 cm
+50 +100	45 cm
+100 +150	45 cm
+150 +200	45 cm
+200 +250	45 cm
+250 +300	45 cm
*One thermom	

*One thermometer per range.

Part No.

1001268 with calibration certificate PROTON.

1001298 without calibration certificate.





Electric contact thermometers

Adjustable temperature setting with the use of a dial on the top of the thermometer. To provide better mercury flow, the capillary interior is filled with an inert gas. Supplied complete with a 100 cm cable with the appropriate connection on the end.

Part No.	Range °C	Division scale °C	Immersion length cm	Total length cm
3521015	0 +100	1/1	15	39
3522015	0 +200	2/1	15	39
3523015	0 +300	5/1	15	39
3352604	0 +60	1/1	4	28
3352104	0 +100	1/1	4	28
3352204	0 +200	2/1	4	28





Universal reactor for laboratory "RV-12"

ADAPTABLE TO ALL THE STANDARDIZED GLASS ELEMENTS WITH A MOUTH OF 100 MM Ø. REACTOR CAPACITY: 12 LITRES.

APPLICATIONS

Universal reactor suitable for all laboratory and pilot plant's operations.

FEATURES

Standard mouth for connection to DN100 glass reactor's covers.

Reactor in stainless steel AISI 316, with a capacity of 12 litres and with a drain tap.

STANDARD EQUIPMENT

- 1. Standard support. Part No. 1001748
- 2. Clamp with support ring for the reactor. Part No. 1001747
- 3. **Reactor** in stainless steel AISI 316, capacity 12 litres with drain tap. External sizes: Height 38 cm. x 25 cm. Ø. Weight: 4,5 Kg Part No. 1001746

COMPLEMENTS

- 4. Clasp for fastening normalized flat gaskets for reactor "RV-12" or any glass reactor with normalized mouth of 100 mm. Ø. Part No. 1001758
- 5. Stirrer RZR-1 with adjustable speed. Part No. 5050111
- 6. Double bosshead MAXIM. Part No. 7000408 or Clamp Part No. 1001765
- 7. Flexible shaft holder for 8 mm. Ø stirring shaft. Part No. 5050903
- **8. Half-moon shape stirring shaft**, 8 mm. Ø stainless steel AISI 304 (length 90 x shovel width 12 cm). Part No. 1001754
- **9. Helix shape stirring shaft,** 8 mm. Ø stainless steel AISI 304 and PTFE shovel (length 90 x shovel width 4 cm). Part No. 1001757
- 10. Electronic controller ELECTEMP-TFT. Part No. 3000897 with connector for inner temperature control, if the heating mantle is used.
- 11. Support plate with bosshead for Electemp-TFT regulator. Part No. 1001750
- **12. Pt100 probe** Part No. **1001749** 200 mm length x 4 mm. \emptyset , with connection cable to Electemp-TFT regulator.
- 13. Ring holder with double bosshead for heating mantle. Part No. 1001751
- 14. Heating mantle with doble heating elements for temperatures up to 400 °C. Power consumption 1400W. Part No. 3001752

OPTIONAL


Aisi 316 stainless steel device, with standard double mouth and coll for liquid heating or refrigeration inside the reactor.

This device can be connected to the universal reactor part no. 1001746 or to glass reactors with a 10 l. capacity and with DN-100 standard mouth.

Coil measurement: A: 21 cm. B: 15 cm. C: 10 cm. Part No. 1001755

For reactors higher than 40 cm., it must be manufactured according to drawing.

Note: Several tests have been made in the laboratory with the silicone oil "47V50" and with the **Ultraterm-TFT-300 thermostat**. Part No. 3000545 (See page 103), reaching temperatures up to 170 °C in 35' inside the reactor properly isolated with "durablanket120". With methanol fluid transmitter liquid and with **Frigiterm-TFT-30 thermostat**. Part No. 3000547 (see page 103), reaching temperatures up to -21 °C in 57' inside the reactor properly isolated with "durablanket120".



Tapa de vidrio



Customizable reactors on demand

FOR LABORATORY AND PILOT PLANT WITH HEATING ANS COOLING SYSTEMS BOROSILICATE GLASS 3.3. CAPACITY: FROM 0,5 UP TO 20 LITRES.
AISI 304 - 316 - 316L STAINLESS STEEL LITRES CAPACITY: FROM 20 UP TO 100 LITRES.


Borosilicate glass 3.3: Spherical, cylindrical and thermostatic with or without PTFE release valve with temperature control accessories, PH, oxygen, turbidity, etc. Heidolph stirring systems and variable PTFE blades.

Aisi 304 - 306 - 306 stainless steel: Cylindrical and thermostatic with or without release valve. With temperature control accessories, PH, oxygen, turbidity, etc.

Fixed and variable stirring systems of different powers according to viscosity and density, combined with stirring elements, dispersers and emulsifiers micronizers and helical.

REQUEST INFORMATION.







Reactor for pilot plant "R-605" and "R-620"

CAPACITY UP TO 20 LITRES

APPLICATIONS

This equipment is designed for reflux and distillation operations with a stirring element, for all laboratory's pilot plant operations and small productions.

FEATURES

Reactor module built-in, made of glazed steel in epoxy, composed of an engine group stirring and a reaction borosilicate glass 3.3 with graduation and dual Camera to temper, for a volume of 4-20 liters in "R-620" model and 1-5 liters in "R-605" model. Stirring anchor blade and elements in contact with the liquid or reaction steam, composed of borosilicate glass and PTFE.

R-605 model addition and reflux. Distillation system on demand.

Reactor for multiple reactions: add PTFE tap Part No. 1001752, support Part No. 1001744 and digital thermometer Part No. 5903030 or pH analyser/Redox and temperature Part No. 4120800.

CONTROL PANEL

Main switch ON / OFF.

Regulator knob for stirring speed from 0 to 800 rpm, and with digital reading.

MODELS

Part No.			nsions h / Depth	Reactor capacity	Operating temperature	Power W	Weight Kg
1001705	130	70	35	5	-30 °C to 200 °C	90	10
1001745	180	62	50	20	-30 °C to 200 °C	250	68



ACCESORIES

pH analyser/Redox and temperature from -10°C up to 130°C. Backlit graphic display. For pH/ ORP measurement, temperature by probe measurement, temperature automatic compensation. 1 relay output for alarm or temperature, 1 relay output for electrode wash or temperature. pH measurement: 00.00÷14.00Ph.

Redox ±1500Mv± 1Mv. Part number **4120800**

HT glass sensor pH 0-14 +Temp -10 to 100°C. Part number 4120801

Glass sensor adapter for R-605. Part number 4120802 Glass sensor adapter for R-620. Part number 4120803

Digital thermometer "TC-9226-A"

Temperature range: -200 °C. to 1370 °C.

Part No. 5903030


Universal drain tap in PTFE and bellow for disposal of different density and viscosity liquid / liquid and liquid / solid suspensions.

Part No. 1001753 for model RV-605. Part No. 1001752 for model RV-620.

Over elevation support base.

Part No. 1001773 for model RV-605. Height 38 cm. Part No. 1001744 for model RV-620. Height 65 cm.





Vacuum inlet

Drip bottle feeding inlet

Drip inlet valve

inlet

Cooling cycle water

Peristaltic dosing pump "Peristaltic PR-2003."

For PH stabilisation (acid-base).

Also applicable to dosages in process.

Features: See page 49. Part No. 4129050

FOR REACTOR R-605:

Recirculation thermostat "Ultraterm-TFT-200".

For adjustable temperatures from ambient +5 °C. to 200 °C.

Part No. 3000545 (See page: 103).

Recirculation cryothermostat "Frigiterm-TFT-30"

For adjustable temperatures from -30 °C. to 100 °C.

Part No. 3000547 (See page: 103).

FOR REACTOR R-620:

High power recirculation thermostat "Interheater

For adjustable temperatures from ambient +5 °C. up to 90 °C.

Part No. **6001491** (See page: 107).

Water refrigerador and recirculator

"Unichiller 015"

For temperatures from -20°C to +40 °C.

Part No. **5910027**





FILTRATION. FILTER PAPER FOR LABORATORY





Quantitative filter papers

ASH CONTENT LESS THAN 0.01%. CRITICAL GRAVIMETRÍES AND ANALYSIS



Filtration	Basis weight g/m ²	Thickness mm	Ashes %	Туре	Units box
Fast	97	0,200	< 0,01	Flat	100
Medium	97	0,200	< 0,01	Flat	100
Slow	97	0,170	< 0,01	Flat	100

Fast: Ø 70: 5301070. Ø 90: 5301090. Ø 110: 5301110. Ø 125: 5301125. Ø 150: 5301150. Ø 180: 5301180.

Medium: Ø 70: 5302070. Ø 90: 5302090. Ø 110: 5302110. Ø 125: 5302125. Ø 150: 5302150. Ø 180: 5302180.

Slow: Ø 70: 5303070. Ø 90: 5303090. Ø 110: 5303110. Ø 125: 5303125. Ø 150: 5303150. Ø 180: 5303180.






General use qualitative filter papers

FILTRATIONS FOR ROUTINE ANALYSIS IN ALL KINDS OF LABORATORIES





Fast: Ø 70: 5304070. Ø 90: 5304090. Ø 110: 5304110. Ø 125: 5304125. Ø 150: 5304150. Ø 180: 5304180.

Medium: Ø 70: 5305070. Ø 90: 5305090. Ø 110: 5305110. Ø 125: 5305125. Ø 150: 5305150. Ø 180: 5305180.

Slow: Ø 70; 5306070, Ø 90; 5306090, Ø 110; 5306110, Ø 125; 5306125, Ø 150; 5306150, Ø 180; 5306180,





Filter paper reams

PROTECTION AND ABSORPTION IN WORKBENCHES IN ALL KIND OF LABORATORIES

Grade	Basis weight g/m2	Thickness mm	Absorption speed* mm/10 min	Moisture resistance kN/m	Units box
1305	73	0.170	75/70	0.290/0.260	500
1510	60	0.130	60/55	0.280/0.230	500

1305: 42 x 52 (cm): 5130502. 32 x 42: 5130501. 1510: 42 x 52(cm): 5151002. 32 x 42: 5151001.





Absorbent paper with polyethylene

PROTECTION AND ABSORPTION THANKS TO ITS DOUBLE LAYER: ONE OF ABSORBENT PAPER AND THE OTHER OF WATERPROOF **POLYETHYLENE**

Grade Basis weight Thickness Area absorption g water/m² 0.250 1505 125

1505: 42 x 52 (cm): 5150502. 32 x 42: 5150501. 50 (cm) Width x 50 (meters): 5150503



Glass microfiber filters without binders

RECOMMENDED IN AIR AND WATER ANALYSIS ACCORDING TO INTERNATIONAL STANDARDS

Grade	Retention µm	Basis weight g/m ²	Thickness mm	Retention DOP %	Binders	Units Box
MFV1	1.6	52	0,26	99.998	No	100
MFV3	1.2	52	0,26	99.998	No	100

 $MFV1: \hspace{0.2cm} mm \hspace{0.1cm} \emptyset \hspace{0.1cm} 47; \hspace{0.1cm} \textcolor{red}{5291047}. \hspace{0.1cm} \emptyset \hspace{0.1cm} 70; \hspace{0.1cm} \textcolor{red}{5291070}. \hspace{0.1cm} \emptyset 90; \hspace{0.1cm} \textcolor{red}{5291090}. \hspace{0.1cm} \emptyset \hspace{0.1cm} 110; \hspace{0.1cm} \textcolor{red}{5291110}. \hspace{0.1cm} \emptyset \hspace{0.1cm} 125; \hspace{0.1cm} \textcolor{red}{5291125}. \hspace{0.1cm} \emptyset \hspace{0.1cm} 150; \hspace{0.1cm} \textcolor{red}{5291150}.$

MFV3: mm Ø 47: 5293047. Ø 70: 5293070. Ø90: 5293090. Ø 110: 5293110. Ø 125: 5293125. Ø 150: 5293150.





Cellulose extraction cartridges

NOBLE CELLULOSES AND COTTON LINERS TOTALLY FREE FROM IMPURITIES.

FOR REMOVING PAINTINGS, ANIMAL FAT, COSMETICS, VARNISHES, LACQUERS, OIL SEEDS, RESINS AND WOOD FATS, ETC... WALL THICKNESS: BETWEEN 1 AND 2 MM. UNITS PER BOX: 25

Part No.	Interior Dimensions mm Ø/ L	Part No.	Interior Dimensions mm Ø/ L	Part No.	Interior Dimensions mm Ø/ L
5297001	19 x 90	5297005	28 x 100	5297009	33 x 94
5297002	22 x 80	5297006	30 x 80	5297010	33 x 100
5297003	22 x 100	5297007	30 x 100	5297011	35 x 100
5297004	25 x 100	5297008	33 x 80	5297012	43 x 123





Cellulose Acetate Membrane Filters

LOW PROTEIN ABSORPTION FOR FILTERING DIFFERENT SAMPLES: BIOLOGICAL, AQUEOUS, ALCOHOLS, OILS AND HIDROCARBONS. AVAILABLE IN PORES FROM 0.22 UP TO 5 µM AND IN DIAMETERS FROM 13 UP TO 293 MM.

AVAILABLE ON DEMAND IN Ø 13 - 25 - 142 Y 293 MM.

PH COMPATIBILITY: 4-8.

Part No.	Membrane	Diameter	Pore	Colour	Surface	Format	Units x
	Material	mm	μm				box
5297013	Cellulose acetate	47	0.22	White	Flat	Non-Sterile	100
5297014	Cellulose acetate	47	0.45	White	Flat	Non-Sterile	100
5297015	Cellulose acetate	47	0.8	White	Flat	Non-Sterile	100
5297016	Cellulose acetate	47	3.0	White	Flat	Non-Sterile	100
5297017	Cellulose acetate	47	5.0	White	Flat	Non-Sterile	100
5297018	Cellulose acetate	90	0.22	White	Flat	Non-Sterile	100
5297019	Cellulose acetate	90	0.45	White	Flat	Non-Sterile	100
5297020	Cellulose acetate	90	0.8	White	Flat	Non-Sterile	100
5297021	Cellulose acetate	90	3.0	White	Flat	Non-Sterile	100
5297022	Cellulose acetate	90	5.0	White	Flat	Non-Sterile	100





TECHNICAL SPECIFICATIONS

Materials:

100% borosilicate glass 3.3 except aluminium clamps.

Capacities:

Funnel 300 ml. Filtration flasks 1000 ml.

Plate porosity: Grade 2/3

Sterilization: Autoclave (max. 134°C), dry heat (max. 180°C).

Diameter filters: 47 mm

Output connections to vacuum: Tube 8 - 10 mm. diameter.

Complete equipment. Part No. 5297023



Vacuum pumps. See page 50 and 51.



AQUISEL

MEDICAL DEVICES PRODUCTOS SANITARIOS

MEDICAL DEVICES MANUFACTURER FOR DISPOSABLE IN VITRO DIAGNOSTIC.

FABRICANTE DE PRODUCTOS SANITARIOS PARA DIAGNÓSTICO IN VITRO. UN SOLO USO.





Eliminate unnecessary risks with the AQUISEL Valve-Stopper system.

Evite riesgos innecesarios con el sistema Tapón-Válvula AQUISEL.



Container with AQUISEL Valve-Stopper.

Contenedor con Tapón-Válvula AQUISEL.



Fixation and concentration of faecal parasites "Aquisfek" Fijación y concentración de parásitos fecales "Aquisfek"





Pipettes for E.S.R. with zero levelling and automatic closure system for conventional tubes or AOUISEL Valve-Stopper.

Pipetas VSG de enrase a cero. Con émbolo de presión para tubos convencionales o con Tapón-Válvula AQUISEL.











B AQUISEL VALVE-STOPPER

System available for sample manipulation without the need of uncorking the tube, maintaining it hermetically sealed and avoiding the risk of contamination.

TAPÓN- VÁLVULA AQUISEL

Sistema que permite su manipulación sin necesidad de destaponar, permaneciendo hermético y evitando el riesgo de contaminación.

The AQUISEL system consist of a tube, valve-stopper and an appropriate reagent for specific use.

The AQUISEL tube is made of a high resistance transparent polypropylene. The valvestopper is a very elastic material and has a special cross-cut design which allows an automatic recover to its initial closed position, after insertion or extraction of the sample.

The tube's reagent contents is offered through different possibilities identified by colours in accordance with International Standards.

It can be used by most of the current analytical procedures, even if they are manual, semi-automatic or automatic procedures.

It allows inter-hospital road and postal transport thanks to its hermetic seal and high resistance to uncorking due to vibrations, hits, cold, heat, etc.

USE AND APPLICATIONS

- 1. Inserting (blood) samples with syringe (without needle).
- 2. Extraction and dosage with semi-automatic pipette.
- 3. Capillary blood extraction.
- 4. Capillary filling for microhematocrite.
- 5. Automatic sampling for haematology in primary tubes.
- 6. Sample road transport.
- 7. Centrifuge heads adaptable.

El sistema AQUISEL se compone de tubo, tapón-válvula y del reactivo apropiado para su uso específico.

El tubo del sistema AQUISEL está fabricado en polipropileno de alta resistencia y de gran transparencia. El tapón-válvula es de material muy elástico y **posee unos cortes que por su diseño especial recuperan automáticamente su hermeticidad, después de la introducción o extracción de la muestra.**

Su contenido en reactivo ofrece diferentes posibilidades que se distinguen por su color según Normas Internacionales.

Puede ser utilizado con la mayoría de los procedimientos analíticos actuales, ya sean manuales, semiautomáticos o automáticos.

Permite el transporte neumático intrahospitalario e incluso por correo por su hermeticidad y gran resistencia a destaponarse con movimientos, golpes, frío, calor, etc

USO Y APLICACIONES

- 1. Introducción de muestra con jeringuilla (sin aguja).
- 2. Extracción y dosificación con pipeta semiautomática.
- 3. Recogida de sangre capilar.
- 4. Llenado de capilar para microhematocrito.
- 5. Toma de muestras automática para hematología en tubos primarios.
- 6. Transporte neumático de muestras.
- Adaptable a cabezales de centrífuga.











Tubes with AQUISEL Valve-Stopper Tubos con Tapón-Válvula AQUISEL

There are different size and volume tubes with additives: anticoagulant, antiglicolytics, clotting accelerators or empty tubes.


Diferentes formatos y volúmenes de tubos con aditivos: anticoagulantes, antiglicólicos, aceleradores de la coagulación o vacíos.



Tubes for small blood volume, size 55x12Tubos para pequeños volúmenes de sangre de 55x12.

Suitable for paediatrics or small sample volumes.

Indicados para pediatría, o para pequeños volúmenes de muestra.



Product presentationPresentación del producto

The tubes are presented in polystyrene racks of 100 units. Every 50 units are covered with a transparent lid for protection.

Los tubos van colocados en gradillas de poliestireno de 100 unidades. Cada 50 unidades se cubren con una tapa transparente de protección.







Accesory: CAP-GALET Accesorio: CAP-GALET

Funnel adapter which can be used together with the primary AQUISEL tubes and allows sample blood transfer to secondary tubes, useful for capillary blood extractions in paediatrics, for low blood volumes, can be used in all auto-analysers for primary tube sampling, etc. All these procedures could be completed without uncorking the primary AQUISEL tubes. Once the CAP-GALET is removed, the tube remains sealed again for conservation or disposal.

Adaptador en forma de embudo que juntamente con los tubos primarios AQUISEL permite trasvasar las muestras a tubos secundarios, realizar extracciones de sangre capilar en pediatría, útil para bajos volúmenes de sangre, ser utilizado en todos los autoanalizadores con toma-muestra de tubo primario, etc. Todos los procesos se realizan sin necesidad de destaponar los tubos primarios AQUISEL. Una vez extraído el CAP-GALET, el tubo permanece de nuevo cerrado para su conservación o desecho.

Product table - Cuadro de productos



AQUISEL TUBES WITH VALVE-STOPPERS - TUBOS AQUISEL CON TAPÓN-VÁLVULA

Nº D	TUBES QUANTITY IN THE SHIPPII E TUBOS EN LAS CAJAS DE EXPEDICI	•		,	600)	800)	1.00	00	800	0	800)	1.00	00
	Nominal size of the tube: le Tamaño nominal del tubo: lo				100 x	:16	75 x16		55 x16 FLAT BOTTOM CON FALDON		80 x13		70 x12		55 x	
	CONTENT / REAGENTS CONTENIDO / REACTIVO		DOSAGE TYPE TIPO DE DOSIFICACION	STOPPER COLOUR COLOR DEL TAPON	Code Código	ml blood sangre	Code Código	ml blood sangre	Code Código	ml blood sangre	Code Código	ml blood sangre	Code Código	ml blood sangre	Code Código	ml blood sangr
K3E	EDTA 3K		LIQUID LIQUIDO	MAUVE Malva					1501103 1501113	2	1581104 1581114	2,5	1581105	3	1501106 1501116	2
K3E	EDTA 3K		LIQUID LIQUIDO	MAUVE MALVA					1501123 1501133	4* 2*	1581124 1581134	4* 2,5*			1501126 1501136	0,5
4 NC	CITRATE 3 Na - CITRATO 3 Na 0,106 M ~ 3,8% 1:4 VSG / ESR		LIQUID LIQUIDO	BLACK NEGRO					1501223	2	1581204	1,6 1	1581205	1,2	1501206 1501226	1 0,32
9 NC	CITRATE 3 Na - CITRATO 3 Na 0,106 M∼3,8% 1:9 CLOTTING FACTORS FACTORES DE COAGULACIO		LIQUID LIQUIDO	LIGHT BLUE AZUL CLARO							1581304	4	1581305	2	1501306 1501326	1 0,5
LH	Li HEPARIN - HEPARINA Li	O.N.	LIQUID LIQUIDO	GREEN VERDE	1561401	10					1581404	4	1581405	2	1501406	1
LH	LI HEPARIN WITH SEPARATOR GE Heparina Li con gel separadi		LIQUID AND GEL Liquido y gel	GREEN VERDE											1501416	1
LH	Li HEPARIN + MONOIODO ACETAT HEPARINA Li + MONOIODO ACETAT		LIQUID LIQUIDO	GREEN VERDE							1581604	4			1501606	1
	MONOIODO ACETATE LI + CLOTTING GRA Monoiodo acetato Li + granulos act		GRANULES GRANULOS	LIGHT GREY GRIS CLARO	1561521	8					1581524	4	1581525	2		
FX	FLUORIDE Na + OXALATE K FLUORURO Na + OXALATO K		LIQUID LIQUIDO	DARK GREY GRIS OSCURO							1581704	4	1581705	2	1501706	1
Z	EMPTY AND UNLABELLED VACIO SIN ETIQUETA			WHITE BLANCO	1561961	10	1581962	7	1501963	5	1581964	5	1581965	3	1501966	2
Z	EMPTY WITH BLANKLABEL FOR SERU VACIOS CON ETIQUETA EN BLANCO PARA SEPARATOR GEL FOR SERUM		POWDER PULVERIZADO	WHITE Blanco White	1561971	10					1581974	5	1581975	3	1501976	
Z	GEL SEPARADOR PARA SUERO PS GRANIII ES FOR SERIIM		GEL GRANULES	BLANCO WHITE	1561981	9					1581984	4			1501986	<u> </u>
Z	GRANULOS PS PARA SUERO CLOTTING ACCELERATOR ACTIVADOR DE LA COAGULACIO		GRANULOS LIQUID LIQUIDO	BLANCO RED ROJO	1561951	9	1581952	6			1581954	4			1501956 1501806	1,5
Z	ACTIVADOR DE LA COAGULACIO CLOTTING ACCELERATOR SEPARATOI ACTIVADOR DE LA COAGULACION Y GEL SEPA	R GEL	LIQUID AND GEL Liquido y Gel	RED ROJO	1561811	9					1581814	4			1501816	
Z	CLOTTING ACCELERATOR PS GRANULES FOR ACTIVADOR DE LA COAGULACION Y GRANULOS PS PA	R SERUM	GRANULES GRANULOS	RED ROJO	1561941	9	1581942	6			1581944	4	1581945	3	1501946	1,5
*;	Special for Haematology	Shippi	ng Lo	ng Largo (cm)	47		47		47		39		36		36	
au	rtomatic machines pecial cut 5 mm)	boxe s		Ith Ancho (cm)	23		23		23		22		21		21	
*E	special para máquinas	Medid Caja	1	igh Álto (cm)	34		34		34		36		34		34	
	ntomáticas de hematología orte especial del tapón 5 mm) 75v13 eize tube for ADUISEL	Expedia	ción vve	ight Peso (Kg)	4,3		4,3		4,0		3,5		3		3,5	

75x13 size tube for AQUISEL P-2 pipette. Tubo 75x13 para pipeta AQUISEL P-2.

Special tube for AQUISEL P-3 pipette.
Tubo especial para pipeta AQUISEL P-3.

Special tube for AQUISEL P-4 MICRO pipette.
Tubo especial para pipeta AQUISEL P-4 MICRO.



ACCESORY: CAP-GALET - ACCESORIO: CAP-GALET

Code 1650000. Supplied in bags of 500 units.

Código **1650000**. Se suministra en bolsas de 500 unidades.

AQUISEL TUBES WITH NON-CUT STOPPER - TUBOS AQUISEL CON TAPÓN SIN CORTES

CONTENT/REAGENTS CONTENIDO/REACTIVO	DOSAGE TYPE TIPO DOSIFICACION	STOPPER COLOUR COLOR TAPON	Nº TUBES PER BOX Nº TUBOS POR CAJA	Code Código	ml blood ml sangre	LENGTH x Ø Tube mm LONG. x Ø Tube mm.	Width / Depth / Height Box cm Ancho / Fondo / Alto Caja cm	Weight Peso Kg
K3E EDTA 3K	LIQUID Liquido	MAUVE MALVA	800	1681154	3	80 x 13	39 22 36	3

AQUISEL PIPETTE

PIPETTES FOR ERYTHROCYTE SEDIMENTATION READ E.S.R. WESTERGREN METHOD.

Pipettes made in polystyrene with similar properties as glass ones.

The ZERO-CLOSURE system allows the pipette filling, preventing that the blood exceeds the zero point and also avoiding total or partial loss. Automatic sealing makes it a disposable product.

PIPETA AQUISEL

PIPETAS PARA LA DETERMINACIÓN DE LA VELOCIDAD DE SEDIMENTACIÓN GLOBULAR V.S.G. MÉTODO WESTERGREN.

Pipetas fabricadas en poliestireno con un comportamiento equivalente al vidrio.

El sistema de CIERRE-CERO permite el llenado de la pipeta impidiendo que la sangre supere el punto cero y evitando el vaciado total o parcial. Al quedar automáticamente sellado se convierte en un artículo de un solo uso.




For any kind of tube use, even with AQUISEL or with a conventional stopper.

Suitable for mechanical or manual aspiration.

Pipeta AQUISEL A-2 sin émbolo.

Para cualquier tipo de tubo, con tapon AQUISEL o convencional.

Llenado por aspiración mecánica o manual.



AQUISEL Pipette P-2 with rubber.

For 75x13 tubes with AQUISEL or with a conventional stopper.

Pipette filling by pressure till the bottom of the tube.

Pipeta AQUISEL P-2 con émbolo.

Para tubos de 75x13 con tapon AQUISEL o convencional.

Llenado de la pipeta por presión hasta el fondo del tubo.








AQUISEL Pipette P-3 and P-4 MICRO systems

Set of AQUISEL pipettes and tubes. Tubes come with the valve-stopper which allows the pipette insertion without the need of uncorking. Pipette filling by pressure.

Sistemas Pipeta AQUISEL P-3 y P-4 MICRO

Conjuntos formados por pipetas y tubos AQUISEL. Los tubos llevan incorporado el TAPÓN-VALVULA lo que permite la introducción de la pipeta sin necesidad de destaponar. Llenado por presión.














PRODUCT PRESENTATION - PRESENTACIÓN DEL PRODUCTO





AQUISEL Pipettes A-2

Expanded polystyrene packaging of 500 units.

Pipetas AQUISEL A-2

Envase de P.S. expandido. De 500 pipetas de capacidad.



AQUISEL Pipettes P-2

200 units cardboard package. Supplied in 800 units boxes.

Pipetas AQUISEL P-2

Envase de cartón de 200 unidades. Se suministran en cajas de 800 unidades.



Set P-3. AQUISEL Pipette and tube

200 units cardboard package and 100 tubes rack made of expanded polystyrene, supplied in 600 sets boxes.

Conjunto P-3. Pipeta con tubo AQUISEL

Envase de cartón de 200 pipetas y de P.S. expandido de 100 tubos. Se suministra en cajas de 600 conjuntos.



Set P-4 Micro. AQUISEL Pipette and tube

100 pipettes cardboard package and 100 tubes rack made of expanded polystyrene, supplied in 400 sets box.

Conjunto P-4 Micro. Pipeta con tubo AQUISEL

Envase de cartón de 100 pipetas y de P.S. expandido de 100 tubos. Se suministra en cajas de 400 conjuntos.

MODEL Modelo	CODE Código	UNITS / BOX N° UNIDADES CAJA	BOX WEIGHT Kg Peso Caja Kg	BOX DIMENSIO Width Ancho	NS cm - MEDI Depth Fondo	IDAS CAJA cm Height Alto
A-2	1750224	500	1,5	27	14	11
P-2	1780225	800	4	55	27	26
Set P-3 Conjunto P-3	1760200	600 pipettes + 600 tubes 600 pipetas + 600 tubos	6	47	43	26
Set P-4 Micro Conjunto P-4 Micro	1740100	400 pipettes + 400 tubes 400 pipetas + 400 tubos	3,6	20	40	37

COMPLEMENTS - COMPLEMENTOS



Methacrylate rack

10 place rack. For P-3 / P-4 set. Code

1770003 For tube Ø 13 mm - set P-3 1770004 For tube Ø 12 mm - set P-4 Micro

Gradilla en metacrilato

10 plazas numeradas. Para conjunto P-3 / P-4. Código

1770003 Para tubo Ø 13 mm - conjunto P-3 1770004 Para tubo Ø 12 mm - conjunto P-4 Micro



Metallic Westergren stand

10 place autoclavable stand.

For any kind of AQUISEL pipettes with Ø 13 mm tubes and A-2 pipettes with Ø 16 mm tubes. Easy and functional system that allows the total set disposal.

Code

1770001 For Ø 13 mm tubes 1770002 For Ø 16 mm tubes Loose tubes for pipettes:

1581209 for AQUISEL P-3 (800 units/box) 1501226 for AQUISEL P-4 Micro (1000 units/box)



Soporte Westergren metálico

10 plazas numeradas. Autoclavable.

Para todo tipo de pipetas AQUISEL con tubos de Ø 13 mm y para pipetas A-2 con tubos de Ø 16 mm. Sistema simple y funcional que permite el vertido total del conjunto.

Código

1770001 Para tubo Ø 13 mm 1770002 Para tubo Ø 16 mm

Tubos sueltos para pipetas:

1581209 para AQUISEL P-3 (800 uds/caja)

1501226 para AQUISEL P-4 Micro (1000 uds/caja)

AQUISEL URINE CONTAINER

CONTAINER FOR URINE SAMPLES COLLECTION WITH AQUISEL VALVESTOPPER.


The AQUISEL VALVE-STOPPER allows the urine manipulation without uncorking, eliminating the contamination risks, samples loss, spillage, odours, etc.

All the AQUISEL urine container presentations will be supplied with the GALLET accessory, which makes easy the sample dosage to a tube for sedimentation analysis, or the urine test strip insertion.



CONTENEDOR PARA LA RECOGIDA DE MUESTRAS DE ORINA CON TAPÓN-VÁLVULA AQUISEL.

El TAPÓN-VÁLVULA AQUISEL permite la manipulación de la orina sin destaponar. Evita los riesgos de contaminación, pérdidas de muestra, derrames, malos olores, etc. Todas las presentaciones del contenedor de orina AQUISEL incorporan el accesorio GALET que facilita la dosificación de la muestra a un tubo para el análisis de sedimentación, o la introducción de la tira reactiva.





Insertion of the GALET in the AQUISEL VALVE-STOPPER.

Introducción del GALET en el TAPÓN-VÁLVULA AQUISEL.



Urine sample extraction through the GALET.

Extracción de la muestra de orina a través del GALET.



Urine test strip insertion.

Introducción de la tira reactiva.



Dosage through automatic or conventional pipette.

Dosificación por medio de pipeta automática o convencional.

CODE	BAG	ASSAMBLED	UNITS BOX	WEIGHT Kg	DIMENSI	ONS cm - MED	IINAS em
CÓDIGO	BOLSA	MONTADO	UDS. CAJA	PESO Kg	Width Ancho	Depth Fondo	Height Alto
1660223	NO NO	NO NO	600	12	60	40	43
1640234	NO NO	YES SI	400	8	60	45	43
1630233	YES SI	YES SI	300	7	60	45	43



Code 1600230 GALET CONTAINER (1000 units package).

Código 1600230 GALET CONTENEDOR (paquetes de 1000 unidades).








FIXATION AND CONCENTRATION OF FAECAL PARASITES "AQUISFEK"


NEW

New compact design of Aquisfek device, which allows an easy and convenient sample collection, a safe transportation with an easier handling system and a more effective filtration.

FIJACIÓN Y CONCENTRACIÓN DE PARÁSITOS FECALES "AQUISFEK"

Nuevo diseño compacto del dispositivo Aquisfek, el cual permite una recogida de la muestra fácil y cómoda, un transporte seguro con un sistema de manipulación más sencillo y una filtración más eficaz.







Procedure

SAMPLE COLLECTION.

EASY TO HANDLE BY THE PATIENT.

Separate the two parts (A and B) of the device and with the aid of the spoon built-in at the end of the filter, collect the sample.

Introduce the spoon with the sample back into the liquid and screw to close.

PROCESSING THE SAMPLES IN THE LABORATORY.

FAST AND EASY.

Once the sample has been homogenized, the whole device is centrifuged. After centrifugation, discard the part A (spoon) and with the aid of a pipette, remove the supernatant in the tube (Part B).

The concentrated sample residue is transferred to a slide for a further microscopic observation.

TWO DIFFERENT FORMATS.

KIT 1. WITH FORMALIN / ACETATE. PARASITE COMMONL FIXER.
KIT 2. WITHOUT FORMALIN TO MINIMIZE THE RISK OF ENVIRONMENTAL
CONTAMINATION.

KIT 1.

Part number

1800402 72 units supplied (24 units with 3 units in each bag).

1800404 100 individual units supplied.

KIT 2. Without formalin.

Part number

1800400 72 units supplied (24 units with 3 units in each bag).

1800401 100 individual units supplied.

KIT 1.

Código

1800402 Se suministran 72 unidades (24 bolsas con 3 unidades en cada bolsa)

KIT 2. Sin formol para minimizar el riesgo de contaminación ambiental.

1800404 Se suministran 100 unidades individuales.

KIT 1. Con Formol / Acetato. Fijador habitual de parásitos.

KIT 2. Sin formol.

RECOGIDA DE LA MUESTRA.

da al filtro coger la muestra.

Introducir y cerrar el sistema.

sobrenadante del tubo (Parte B).

DOS FORMATOS DIFERENTES.

FÁCIL Y RÁPIDO.

ga el dispositivo.

al microscopio.

GRAN SENCILLEZ DE MANIPULACIÓN POR EL PACIENTE.

PROCESO DE LAS MUESTRAS EN EL LABORATORIO.

Separar las dos partes (A y B) del dispositivo y con ayuda de la cucharilla incorpora-

Después de homogenizar la muestra, sin extraerla del tubo de recogida, se centrifu-

Desechar la parte de la cuchara del dispositivo (Parte A) y separar, con una pipeta, el

El residuo centrifugado se transifere a un porta para la observación de los parásitos

Código

1800400 Se suministran 72 unidades (24 bolsas con 3 unidades en cada bolsa)

1800401 Se suministran 100 unidades individuales.


Note: For centrifugation, Cencom 2 is recommended. Part number 7002240. See page 114

Nota: Para su centrifugación se recomienda la centrifuga Cencom 2. Código 7002240. Ver pág 142



Tubes with screw cap





1. 10 ml round bottom tubes Ø16 x 100mm.

Polystyrene, screenprinted and with screw cap. Packed 50 units.

Part No. 5004000 Non-sterile. Part No. 5004001 Sterile.

2. 10 ml round bottom tubes Ø16 x 100mm.

Polystyrene and with screw cap.

Part No. 5004002 Non-sterile. Packed 1000 units. Part No. 5004003 Sterile. Packed 200 units.

3. 13.5 ml conical bottom tubes. Ø16 x 110mm.

Polystyrene and with screw cap.

Rcf: 8,000 g.

Part No. 5004004 Non-sterile. Packed 1000 units. Part No. 5004005 Sterile. Packed 200 units.

4. 15 ml round bottom tubes Ø16 x 120mm.

Polystyrene and with screw cap.

Part No. 5004006 Non-sterile. Packed 750 units. Part No. 5004007 Sterile. Packed 150 units.

5. 15 ml conical bottom tubes. Ø17 x 120mm.

Polystyrene, screenprinted and with screw cap. Rcf: 3,500 $\ensuremath{\text{q}}.$

Packed 150 units.

Part No. 5004008 Non-sterile. Part No. 5004009 Sterile.

6. 15 ml conical bottom tubes. Ø17 x 120mm.

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables.

Packed 150 units.

Part No. 5004010 Non-sterile.. Part No. 5004011 Sterile.

7.50 ml flat bottom conical tubes. Ø30 x 115mm.

Polystyrene, screenprinted and with screw cap. Autoclavables.

Packed 100 units.

Part no. 5004012 Non-sterile. Part no. 5004013 Sterile.

8. 50 ml conical bottom tubes. Ø30 x 115mm.

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables.

Packed 100 units.

Part no. **5004014** Non-sterile. Part no. **5004015** Sterile.

9. 50 ml conical bottom tubes. Ø30 x 115mm.

For photosensitive materials.

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables at 121 °c. Resistant to temperatures -80°C.

Packed 25 units.

Part no. 5004016 Sterile..

10. Racks for cylindrical and conical tubes

Part no.	Made of	Tubes	For tubes
		number	up to Ø
1001715	inox. AISI 304	15 (5 x 3)	17 mm
1025524	inox. AISI 304	24 (8 x 3)	17 mm
1001716	inox. AISI 304	8 (4 x 2)	30 mm





Petri dishes In polystyrene.

1. Plates without compartments. \emptyset 90 mm. 3 winds

Part No. 5004017 Aseptic. Packed 20 units. Part No. 5004018 Sterile. Packed 20 units.

2. Two compartment plates. Ø 90 mm. 3 winds.

Part No. 5004019 Aseptic. Packed 20 units. Part No. 5004020 Sterile. Packed 20 units.

3. Three compartment plates. Ø 90 mm. 3 winds.

Part No. 5004021 Aseptic. Packed 20 units. Part No. 5004022 Sterile. Packed 20 units.

4. Four compartments plates. Ø 90 mm. 3 winds.

Part No. 5004023 Aseptic. Packed 20 units. Part No. 5004024 Sterile. Packed 20 units.

5. Contact Petri dishes crosslinked for counting.


Part No. 5004025 Sterile. Ø 55 mm. Packed 10 units.

Part No. 5004026 Sterile. Ø 90 mm. Packed 20 units. 6. Square Petri dishes. 120x120 mm. 3 winds.

Part No. 5004027 Sterile. Packed 10 units.

7. Transport and storage support for Petri dishes of Ø 90 mm.

Capacity: 48 plates. Part No. 5004028









Swab in tube, sterile

CE CLASS LLA (MDD)

Sterile swab in polypropylene round bottom tube with sticker-sealed, sterility assurance

Dimensions: Ø13 X 165MM. Sterilized by ethylene oxide.

Part No.	Features	Quantity box units	Weight Box kg
1001660	Wood + cotton	4x500	14
1001661	Snappable PS +cotton	4x500	14,20
1001662	Snappable PS +viscose	4x500	14,20



SWAB: CE CLASS LLA (MDD). TUBE WITH MEDIA: CE (IVD)

AMIES media is a modification of Cary Blair media, which is at the same time a modification of Stuart media. Basically the glycerophosphate is replaced by an inorganic phosphate, and the methylene blue by pharmacological charcoal. Calcium and magnesium are also added, thus maintaining the permeability of the bacterial cell. This media assures the viability of organisms such as: Neisseria sp., Haemophilus sp., Corynebacterium sp., Trichomonas vaginalis Streptococcus pyogenes, Streptococcus pneumoniae, Shigella flexneri, Salmonella typhi, Brucella abortus, Enterobacterias, etc. Some microorganisms can resist on the media for up to 3 days, although the recovery of microorganisms is better if cultured in the first 24 hours. Amies transport swabs are available with or without charcoal.

Swabs are sterilised by radiation.

Expiry date: 30 months after sterilisation date

Peel-pack dimensions: 38 x 210 mm.

Part No.	Features	Quantity box units	Weight box kg
1001663	Snappable polystyrene + viscose	6x100	8,5
1001664	Snappable polystyrene + viscose (with charcoal)	6x100	8,5





SWAB: CE CLASS LLA (MDD). TUBE WITH MEDIA: CE (IVD)

The modified Stuart media allows the conservation and transportation of a large number of pathological microorganisms, such as: Neisseria gonorrhoeae, Heamophilus influenzae, Neisseria meningiditis, Bordetella pertusis, Corynebacterium diphteriae, Trichomonas vaginalis, Staphylococcus aeureus, Streptococcus sp., Salmonella sp., Shigellas sp., etc. Once used, it can be for up to 24 hours. It is very reduced due to the presence of thioglycolate, which difficults the enzymatic reactions of the bacteria. The multiplication of the bacteria is prevented due to the lack of nitrogen in the media. Swabs are sterilised by radiation.

Expiry date: 30 months from sterilisation date.

Dimensions of the peel-pack: 38 x 210 mm.

1001665	Wood + cotton	6x100	8,5
Part No.	Features	Quantity box units	Wight box kg



SWAB: CE CLASS LLA (MDD). TUBE WITH MEDIA: CE (IVD)

Specially formulated for the transport of virus, very useful for the preservation of the sample during transportation to the laboratory. The media maintains the sample viable for up to 72 hours and in some cases for longer time. The addition of antimicrobial substances inhibit the growth of bacteria and fungi. Ideal for nasal, pharyngeal samples and also ocular, skin and mucus samples.

Swabs are sterilized by radiation.

Expiry date: 18 months from sterilisation date.

Part No.	Features	Quantity box units	Wight box kg
1001666	Snappable polystyrene + dacon	6x100	7,8





