


- ① **RV 10 digital V**
EUR 4.517,00
Ident. No. 0010004799
- ② **RV 10 digital V-C**
EUR 4.896,00
Ident. No. 0010004802
- ③ **RV 10 digital FLEX**
EUR 4.336,00
Ident. No. 0010004805
- ④ **RV 10 digital dry ice condensor**
EUR 5.012,00
Ident. No. 0010006365
- ⑤ **RV 10 digital dry ice condensor C**
EUR 5.390,00
Ident. No. 0010006364




TECHNICAL DATA	RV 10 DIGITAL
Motor type	Brushless DC drive motor
Speed range	5 – 280 rpm
Speed display	Digital
Clockwise and counter-clockwise	Yes
Heating temperature range	RT – 180 °C
Type of cooling	Vertical
Cooling surface	1500 cm ²
Filling volume max.	4 l
Soft start	Yes
Head angle adjustable	0 to 45°
Stroke displacement	140 mm, motorized
Setting of lower end stop	60 mm, contact free
Vacuum controller	Optional
Dimensions (W × H × D) without glass parts	500 × 430 × 440 mm
Protection class acc. to DIN EN 60529	IP 20



- ① **RV 10 auto V**
EUR 7.693,00
Ident. No. 0010004835
- ② **RV 10 auto V-C**
EUR 8.150,00
Ident. No. 0010004838
- ③ **RV 10 auto FLEX**
EUR 7.516,00
Ident. No. 0010004841
- ④ **RV 10 auto pro V**
EUR 9.476,00
Ident. No. 0010004815*
- ⑤ **RV 10 auto pro V-C**
EUR 9.834,00
Ident. No. 0010004812*
- ⑥ **RV 10 auto pro FLEX**
EUR 9.536,00
Ident. No. 0010004818*
- ⑦ **RV 10 auto pro V-C Complete**
EUR 10.889,00
Ident. No. 0010007353



TECHNICAL DATA	RV 10 AUTO
Motor type	Brushless DC drive motor
Speed range	5 – 300 rpm
Speed display	Digital
Clockwise and counter-clockwise	Yes
Heating temperature range	RT – 180 °C
Type of cooling	Vertical condenser
Cooling surface	1500 cm ²
Filling volume max.	4 l
Soft start	Yes
Head angle adjustable	0 to 45°
Stroke displacement	140 mm, motorized
Setting of lower end stop	60 mm, contact free
Vacuum controller	Integrated
Dimensions (W × H × D) without glass parts	500 × 430 × 440 mm
Protection class acc. to DIN EN 60529	IP 20



*Pro packages include VACSTAR digital