



DB 1.1



DB 3.2



DB 4.6



DB 4.6



DB 5.2



DB 5.8



DB 8.1

Product	Description	Ident. No.	Price
MICRO CENTRIFUGE TUBES			
DB 1.1	Single block, PCR tubes, Ø 7,9 mm	0004467600	EUR 109,00
DB 1.2	Single block, Eppendorf tubes, Ø 11,5 mm	0004467700	EUR 109,00
DB 1.3	Single block, microtubes, Ø 11,1 mm	0004467800	EUR 109,00
DB 1.4	Single block, VWR Eppendorf tubes, Ø 10,9 mm	0004467900	EUR 109,00
DB 1.5	Single block, corning tubes, Ø 10,9 mm	0004468000	EUR 109,00
CONICAL BOTTOM CENTRIFUGE TUBES			
DB 2.1	Single block, conical tubes, Ø 17,1 mm	0004468100	EUR 109,00
DB 2.2	Single block, conical tubes, Ø 29 mm	0004468200	EUR 109,00
COMBINATION BLOCKS			
DB 3.1	Single block, test tube comb., 6 12 13 25 mm	0004468300	EUR 109,00
DB 3.2	Single block, centrifuge comb., 1,5 15 50 ml	0004468400	EUR 109,00
DB 3.3	Single block, microtube comb., 0,5 1,5 2 ml	0004468500	EUR 109,00
STANDARD TEST TUBES			
DB 4.1	Single block, 30 round tubes, Ø 8,3 mm	0004468600	EUR 116,00
DB 4.2	Single block, 24 round tubes, Ø 10,7 mm	0004468700	EUR 116,00
DB 4.3	Single block, 16 round tubes, Ø 13,9 mm	0004468800	EUR 116,00
DB 4.4	Single block, 20 round tubes, Ø 13,9 mm	0004468900	EUR 116,00
DB 4.5	Single block, 12 round tubes, Ø 17,5 mm	0004469000	EUR 116,00
DB 4.6	Single block, 8 round tubes, Ø 21 mm	0004469100	EUR 116,00
DB 4.7	Single block, 6 round tubes, Ø 26,2 mm	0004469200	EUR 116,00
DB 4.8	Single block, 4 round tubes, Ø 35 mm	0004469300	EUR 116,00
DB 4.9	Single block, 12 round tubes, Ø 19,1 mm	0004469400	EUR 116,00
VIALS			
DB 5.1	Single block, 12 mm vials, Ø 12,7 mm	0004469500	EUR 116,00
DB 5.2	Single block, 15 mm vials, Ø 15,8 mm	0004469600	EUR 116,00
DB 5.3	Single block, 17 mm vials, Ø 17,8 mm	0004469700	EUR 116,00
DB 5.4	Single block, 19 mm vials, Ø 19,7 mm	0004469800	EUR 116,00
DB 5.5	Single block, 21 mm vials, Ø 21,7 mm	0004469900	EUR 116,00
DB 5.6	Single block, 23 mm vials, Ø 23,8 mm	0004470000	EUR 116,00
DB 5.7	Single block, 25 mm vials, Ø 25,8 mm	0004470100	EUR 116,00
DB 5.8	Single block, 28 mm vials, Ø 28,8 mm	0004470200	EUR 116,00
DB 5.9	Single block, 16 mm vials, Ø 16,4 mm	0004470300	EUR 116,00
STANDARD TEST TUBES			
DB 6.1	Single block, 10 × 8 PCR tube strips, Ø 6,4 mm	0004470400	EUR 109,00
DB 6.2	Single block, 64 PCR tubes, Ø 6,4 mm	0004470500	EUR 109,00
DB 6.3	Double block, 1 × 96-Well PCR plate, Ø 6,4 mm	0004470600	EUR 267,00
DB 7.1	Double block, 96- or 384-Well-plate	0004470800	EUR 160,00
DB 8.1	Single block, 2 × 6 cuvettes, Ø 12,5 mm	0004470700	EUR 109,00
TEMPERATURE SENSORS			
PT 1000.60	Temp. sensor, stn. steel, Ø 3 mm, depth: 230 mm	0003516800	EUR 133,00
PT 1000.90	Temp. sensor, coated, Ø 3 mm, depth: 230 mm	0004480600	EUR 188,00

C-MAG HP 4

EUR 529,00

Ident. No. 0003581600

**C-MAG HP 7**

EUR 582,00

Ident. No. 0003581800

**C-MAG HP 10**

EUR 795,00

Ident. No. 0003582000



TECHNICAL DATA	C-MAG HP 4	C-MAG HP 7	C-MAG HP 10
Heat output	250 W	1000 W	1500 W
Heating temperature range	50 °C – 500 °C	50 °C – 500 °C	50 °C – 500 °C
Connection for ext. temperature sensor	–	ETS-D5	ETS-D5
Heating rate heating plate	2,5 K / min	5 K / min	5 K / min
Set-up plate material	Ceramic	Ceramic	Ceramic
Set-up plate dimensions	100 × 100 mm	180 × 180 mm	260 × 260 mm
Dimensions (W × H × D)	150 × 105 × 260 mm	220 × 105 × 330 mm	300 × 105 × 415 mm
Protection class according to DIN EN 60529	IP 21	IP 21	IP 21