

The right pipette tip for your Transferpette®

Optimal results are achieved in combination with the BRAND micropipettes. Quality tips from BRAND provide a perfect fit for precise analysis.

Pipette and filter tips from BRAND are tested for BRAND pipettes and for most pipette models from GILSON®, Thermo Fisher Scientific FINNPIPETTE®, Eppendorf®, and Sartorius®.

The 5 ml tip is tested exclusively for BRAND and Thermo Fisher Scientific FINNPIPETTE®. The 10 ml tip is tested for BRAND and Eppendorf®.

Note: Pipette shafts may be subject to modifications. Compatibility depends, among other things, on the manufacturer, pipette type, serial number, and manufacturing date.

We recommend using the free sample bag to check tip compatibility.




Transferpette® S



Transferpette® electronic

| Pipette tips Filter tips Volume range | Transferpette® single-channel Nominal volume | | | | | | | | | | | | | Transferpette® multi-channel Nominal volume | | | | |
|---|--|-----------------|-----------------|-------|-------|-------|--------|--------|--------|---------|------|------|-------|---|-------|--------|--------|--------|
| | 1 µl | 2,5 µl | 10 µl | 20 µl | 25 µl | 50 µl | 100 µl | 200 µl | 500 µl | 1000 µl | 2 ml | 5 ml | 10 ml | 10 µl | 50 µl | 100 µl | 200 µl | 300 µl |
| 0,1 - 20 µl | ✓ | ✓ ₂₎ | ✓ ₂₎ | | | | | | | | | | | ✓ | | | | |
| 0,5 - 20 µl | | ✓ | ✓ | | | | | | | | | | | ✓ | | | | |
| 1 - 50 µl | | ✓ ₂₎ | ✓ | | | | | | | | | | | ✓ | | | | |
| 2 - 200 µl ¹⁾ | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | | ✓ | ✓ | ✓ | ✓ |
| 5 - 300 µl | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | | ✓ | ✓ | ✓ | ✓ |
| 50 - 1000 µl ¹⁾ | | | | | | | | | ✓ | ✓ | | | | | | | | |
| 0,5 - 5 ml | | | | | | | | | | | ✓ | ✓ | | | | | | |
| 1 - 10 ml | | | | | | | | | | | | | ✓ | | | | | |
| 0,1 - 1 µl | ✓ | ✓ ₂₎ | ✓ ₂₎ | | | | | | | | | | | ✓ | | | | |
| 0,5 - 10 µl | | ✓ | ✓ | | | | | | | | | | | ✓ | | | | |
| 1 - 20 µl | | | ✓ | | | | | | | | | | | ✓ | | | | |
| 2 - 20 µl | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | | ✓ | ✓ | ✓ | |
| 5 - 100 µl | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | | ✓ | ✓ | ✓ | ✓ |
| 5 - 200 µl | | | | | | | | ✓ | | | | | | | ✓ | ✓ | ✓ | ✓ |
| 50 - 1000 µl | | | | | | | | | ✓ | ✓ | | | | | | | | |



✓ = Tip volume less than pipette's nominal volume

¹⁾ Bulk tips are yellow or blue colored; racked tips are colorless in a yellow-colored or blue-colored tip tray

²⁾ Calibration and possible adjustment needed