

Accuracy table



Offers 59 different partial volumes with different numbers of dispensing steps, depending on the PD-Tip size and stroke setting number used.

| Setting | Tip size (ml) | | | | | | | | | | Steps |
|-------------|---------------|-----|-----|-------|-----|-----|------|------|------|------|-------|
| | 0.1 | 0.5 | 1 | 1.25 | 2.5 | 5 | 10 | 12.5 | 25 | 50 | |
| 1 | 2 | 10 | 20 | 25 | 50 | 100 | 200 | 250 | 500 | 1000 | 49 |
| 1.5 | 3 | 15 | 30 | 37.5 | 75 | 150 | 300 | 375 | 750 | 1500 | 32 |
| 2 | 4 | 20 | 40 | 50 | 100 | 200 | 400 | 500 | 1000 | 2000 | 24 |
| 2.5 | 5 | 25 | 50 | 62.5 | 125 | 250 | 500 | 625 | 1250 | 2500 | 19 |
| 3 | 6 | 30 | 60 | 75 | 150 | 300 | 600 | 750 | 1500 | 3000 | 15 |
| 3.5 | 7 | 35 | 70 | 87.5 | 175 | 350 | 700 | 875 | 1750 | 3500 | 13 |
| 4 | 8 | 40 | 80 | 100 | 200 | 400 | 800 | 1000 | 2000 | 4000 | 11 |
| 4.5 | 9 | 45 | 90 | 112.5 | 225 | 450 | 900 | 1125 | 2250 | 4500 | 10 |
| 5 | 10 | 50 | 100 | 125 | 250 | 500 | 1000 | 1250 | 2500 | 5000 | 9 |
| Volume (µl) | | | | | | | | | | | |

