

**New Edition**  
**Improved hardware**  
**Compact design**

- Gold Plated BNC
- Multichannel
- Pre-programmed standards
  - pH: 11 preprogrammed + 5 user editable
  - Conductivity: 3 preprogrammed + 3 user editable
- No interference between electrodes
- Stability algorithm with intuitive indicator
- Hold function
- Selectable resolution
- Range lock
- Capacitive compensation
- Galvanic isolated USB interface
- High Accuracy
- Free software and firmware updates



The C3210 and C3230 are all full-parameter multi-channel instruments. 2 channels with each channel having it's own dedicated measurement hardware without interference between the channels. All values can be displayed simultaneously on the screen.

New SZ20T: maintenance free DO probe



SH300

### Specifications depending on model

Measurement Channels	2
Temperature Channels	2
pH	-2.000...+16.000 pH
mV	±2000.0 mV
Ion	0.01 ng/l...100 g/l (C3230 only)
Conductivity	0...2000 mS/cm
Resistivity	0...200 MΩ.cm
Salinity	0.0...70.0
TDS	0...100.0 g/l
Dissolved oxygen	0...60.00 mg/l 0...600%
Air pressure	600...1300 hPa
Temperature	-5.0...+105.0°C
Warranty	36 months
Made in Belgium	

Code	Description
<b>C3210</b>	2 channel pH/conductivity/DO meter
<b>C3230</b>	2 channel pH/Ion/conductivity/DO meter

Code	Description	Content
<b>C32xxP</b>	pH meter kit, glass electrode	meter + pH electrode SP20T + 2x50ml pH buffers + 50ml electrolyte
<b>C32xxPE</b>	pH meter kit, epoxy electrode	meter + pH electrode SP10T + 2x50ml pH buffers + 50ml electrolyte
<b>C32xxK</b>	EC meter kit, glass electrode	meter + EC electrode SK20T + 50ml EC standard
<b>C32xxKE</b>	EC meter kit, epoxy electrode	meter + EC electrode SK10T + 50ml EC standard
<b>C32xxPK</b>	pH/EC meter kit, glass electrodes	meter + pH electrode SP20T + EC electrode SK20T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard
<b>C32xxPKE</b>	pH/EC meter kit, epoxy electrodes	meter + pH electrode SP10T + EC electrode SK10T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard
<b>C32xxZ</b>	Oxygen meter kit	meter + DO electrode SZ20T
<b>C32xxT</b>	Complete meter kit, glass electrodes (DO epoxy)	meter + pH electrode SP20T + EC electrode SK20T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard + DO electrode SZ20T
<b>C32xxTE</b>	Complete meter kit, epoxy electrode	meter + pH electrode SP10T + EC electrode SK10T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard + DO electrode SZ20T
<b>C32xxX</b>	Meter kit without electrodes	meter + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard
<b>SH300</b>	Flexible electrode holder (optional)	

Kits are available for each meter. Replace **xx** with correct meter number.



Artistic in science

Consort bvba  
Parklaan 36  
B-2300 Turnhout  
Belgium

Telephone [\[+32\] \[0\]14 41 12 79](tel:+322014411279)  
Email [info@consort.be](mailto:info@consort.be)  
Sales [sales@consort.be](mailto:sales@consort.be)  
Support [support@consort.be](mailto:support@consort.be)  
Website [www.consort.be](http://www.consort.be)