

Product name:

Methanol**Synonyms:** MeOH, Hydroxymethane, Methyl alcohol, Carbinol**CAS number:** 67-56-1

<i>General information:</i>		
	Chemical formula: CH ₃ OH	Molar Mass: 32.04 g/mol
	Appearance: Colorless liquid	Density: 0.792 g/cm ³ (20 °C)
	Solub. in water: Soluble (20 °C)	Melting point: -98 °C
	Boiling point: 64.5 °C	Flash point: 11 °C
	Storage temp.: below +30 °C	Cautions: Highly flammable, Toxic
	R: 11-23/24/25-39/23/24/25	S: 7-16-36/37-45
	HS-No.: 2905 11 00	EC-Index-No.: 603-001-00-X

Hazard Symbol



Safety Information

Grade

Extra pure

Specification:				
	Assay	≥	99	%
	Identification A (IR)		Conforms	
	Identification B (GC)		Conforms	
	Acidity	≤	0.45	mL
	Alkalinity	≤	0.2	mL
	Readily carbonizable substances		Conforms	
	Readily oxidizable substances		Conforms	
	Nonvolatiale residue	≤	0.01	%
	Acetone and aldehydes	≤	0.03	%
	Water	≤	0.5	%

Ordering code	1 L	2.5 L	5 L	10 L	20 L
P.E.	211.06.01000.02	211.06.02500.02	211.06.05000.02	211.06.10000.02	211.06.20000.02
Glass		211.06.02500.01			