evosafe ecologic



SNIJDERS ULTRA-LOW TEMPERATURE -86°C FREEZERS VF475-86G



Туре	Upright		
Volume total	480 litre		
Usable volume 2/4 compartments	447/440 litre		
Storage capacity 2" boxes	· ·		
External dimensions	352		
	905/850 x 888 x 1992 mm		
Internal dimensions	600 x 626,5/592 x 1276 mm		
Compartments/inner doors/platforms	4/2/3 in height adjustable		
Main door	With lock and keys		
Inner doors	Stainless steel, insulated		
Material internal/external	Stainless Steel/Epoxy Coated steel		
Placed on	Heavy duty castors, in front with brakes		
Insulation	Polyurethane		
Cooling	Two (2) cascade connected compressors		
Refrigerants	First stage: R455A/Second stage: R170		
Condenser	Air-cooled		
Temperature adjustments	-50°C till -86°C (other range; on request)		
Controller	Web-based		
Display	7" Touchscreen in main door at eye level		
Temperature history	On display with scroll function		
Data storage	At least 12 months		
Temperature sensors	PT1000		
Backup battery	Alarming and monitoring during power failure		
Alarming	Optical and acoustic, but also possible via E-Mail		
Communication	Telephone dialling activation system		
	Potential (Voltage) free contact		
	RJ45 (Ethernet) connection		
	USB connection		
Software	I-freeze		
CO ₂ -backup system	Partially prepared		
Transit/vacuum release	Ø 26 mm/vacuum release situated in transit		
Cable entry ports	Two (2)		
Compressor room	Well-arranged layout and very easy filter exchange		
Requirements	220/240 VAC, 50 Hz, fuse 16 Amp. slow		
Power and current rating	1015 VA/950 W/4,7 Amp.		
Energy consumption	11 kWh/24h		
Heat exchange	0,43 kW/h		
Noise level	53 dBA		
Net weight	315 kg		
CE certified	EN61010-1:2010, EN60601-1-2:2007		



evosafe ecologic



SNIJDERS ULTRA-LOW TEMPERATURE -86°C FREEZERS VF475-86G

OPTIONS

Hybride Cooling System

Water cooling

CO₂-backup system

LN₂-backup system

Extended backup system

Temperature chart recorder

Voltage stabilizer

Transformer

IQ/OQ/PQ

Electronic lock main door

Padlock on main door

Left hinged door(s)

3, 4, 5 or 6 compartments/inner doors

Inner doors lockable by key

Levelling feet

Alarming by texting and/or phone

Sliding shelves with hinged front

Drawers in Aluminium or Stainless Steel

Racking systems in Aluminum or Stainless Steel

Cryo Boxes

Cryo Gloves




Laurent Janssensstraat 105 5048 AR Tilburg The Netherlands T +31 13 750 15 55 F +31 13 463 86 35 www.snijderslabs.com









DRAWER RACKING SYSTEM*1

Model #	Comp.	# Racks	# Boxes per rack	Total boxes 133x133	Total plates 128x85	Max. height box
DSA-7336-136G	4	16	42		672/784*2	33 mm
DSA-5504-136G	4	8	20	352*²	528/616*2	50 mm
DSA-6504-136G	4	8	24	35212		
DSA-11524-136G	2	8	44	352	528/616*2	52 mm
DSA-D211524-136G	2	4	88	352	528/616*2	52 mm
DSA-5534-136G	4	16	20	320	480/560*2	53 mm
DSA-7844-136G	2	8	28	224		84 mm
DSA-D27844-136G	2	4	56	224		84 mm
DSA-3914-136G	4	16	12	192		91 mm
DSA-21284-136G	4	16	8	128		128 mm

*1 other dimensions on request.

*2 only on request.

SIDE ACCESS RACKS

Model #	Comp.	# Racks	# Boxes per rack	Total boxes 133x133	Maximum height box
HR-4-6-45	4	16	24	384	45 mm
HR-4-5-53	4	16	20	320	53 mm
HR-4-3-76	4	16	12	192	76 mm
HR-5-FLEX	4	16	20/12	320/192	55/95 mm















