

## اتوکلاو ۲۵ لیتری Autoclave 25Lit



### مشخصات فنی :

- تمام استیل ( درب ، طوقه ، مخزن ، کاور دستگاه و لوله کشی )
- دارای برد دیجیتال با تکنولوژی PID جهت همدمایی بالاتر
- قابلیت تنظیم دما با گستره دمایی ۱۰۵ تا ۱۲۵ درجه سانتیگراد و تنظیم زمان ۱ تا ۹۹ دقیقه
- سیستم مصرف بهینه انرژی با راندمان بالا
- سیستم کنترل پسماند دمای المنت
- دارای کنکتور الکترونیکی جهت قطع کامل برق المنت ها در زمان بروز خطاها
- کنترل سطح آب الکترونیکی ( با سیستم قطع برق المنت در زمان کمبود آب )
- دارای سطل آلومینیومی برای بارگذاری وسایل
- دارای شیر اطمینان طرح ایتالیایی مخصوص اتوکلاو ( با بدنه برنجی و فنر استیل )
- گیره های نگهدارنده درب نشکن و آنتی باکتریال ( ABS تزریقی )
- دارای یک سال گارانتی و ده سال خدمات پس از فروش
- دارای مهر استاندارد، ایزو ۱۳۴۸۵ (استاندارد مدیریت کیفیت تجهیزات پزشکی) و ایزو ۹۰۰۱

### Specifications

- All steel (door, crown, tank and cover)
- Has a digital range with P.W.M technology
- High efficiency energy consumption system
- Element temperature residual control system
- Electronic water level control (Elemental power cut during water scarcity)
- Anti-bacterial and antiviral retaining clips
- Standard Seal, ISO 13485, ISO 9001 and Approvals of the General Office of Medical Equipment